

General Disclaimer

One or more of the Following Statements may affect this Document

- This document has been reproduced from the best copy furnished by the organizational source. It is being released in the interest of making available as much information as possible.
- This document may contain data, which exceeds the sheet parameters. It was furnished in this condition by the organizational source and is the best copy available.
- This document may contain tone-on-tone or color graphs, charts and/or pictures, which have been reproduced in black and white.
- This document is paginated as submitted by the original source.
- Portions of this document are not fully legible due to the historical nature of some of the material. However, it is the best reproduction available from the original submission.

DR NO. MA-04

SD 74-SA-0047-2

ANALYSIS OF COMMERCIAL EQUIPMENT AND
INSTRUMENTATION FOR SPACELAB PAYLOADS
(CONTRACT NAS8-30541)

VOLUME II. TECHNICAL REPORT

SEPTEMBER 16, 1974

Approved:

R.P. Arras

R. P. Arras,
Study Manager

(NASA-CR-120779) ANALYSIS OF COMMERCIAL
EQUIPMENT AND INSTRUMENTATION FOR SPACELAB
PAYLOADS, VOLUME 2 (Rockwell International
Corp., Downey, Calif.) 383 p HC \$10.25

N75-24808

Unclas

CSCL 22B G3/18 24200



Space Division
Rockwell International



Space Division
Rockwell International

FOREWORD

The Analysis of Commercial Equipment and Instrumentation for Spacelab Payloads was performed by the Space Division of Rockwell International Corporation under Contract NAS8-30541 for the George C. Marshall Space Flight Center of the National Aeronautics and Space Administration. The study explores the feasibility of using commercially available laboratory equipment and instrumentation in the Spacelab, a Shuttle payload, in support of various sortie-mode experiments. The work was managed by Richard P. Arras (Telephone (213) 594-3807) of the Applications Programs area of the Space Division of Rockwell International. The study was administered under the technical direction of Mr. Charles W. Quantock (Telephone (205) 453-3425) of the Payload Studies group of Program Development of the Marshall Space Flight Center.



ACKNOWLEDGEMENTS

The continuing guidance and support given by the Contracting Officer's representative, Mr. C. W. Quantock of MSFC, is gratefully acknowledged. In addition, acknowledgement is made to Messrs. J. Hendrick, P. Fagan, P. Hinderer, E. Kraly, T. Meinhardt, M. Pope and E. Tivenan of Rockwell International, and G. Bordeaux, J. Hammond, and C. Stata of Beckman Instruments for their technical assistance. Without their technical support and the cooperation of many commercial equipment manufacturers this study could not have been accomplished. Finally, the study manager wishes to acknowledge Ms. L. Aldrich and Ms. C. Perry who typed all the material in these volumes.

PRECEDING PAGE BLANK NOT FILMED

CONTENTS

Section		Page
1.0	INTRODUCTION	1-1
1.1	STUDY OBJECTIVES	1-1
1.2	OVERALL APPROACH	1-2
1.3	SUMMARY	1-5
	1.3.1 Selection of Equipment for Analysis	1-5
	1.3.2 Analysis of Selected Equipment	1-5
	1.3.3 Analysis of Specifications	1-6
	1.3.4 Modification Source Tradeoff	1-7
	1.3.5 Recommended Equipment Procurements	1-7
	1.3.6 Conclusions	1-7
2.0	SELECTION OF EQUIPMENT FOR ANALYSIS	2-1
2.1	INTRODUCTION	2-1
2.2	FUNCTIONAL REQUIREMENTS DATA BANK	2-2
	2.2.1 Explanation of Tabulation	2-3
	2.2.2 Functional Modifications	2-5
2.3	SHOPPING LIST	2-7
	2.3.1 Explanation of Tabulation	2-7
2.4	SELECTION OF CANDIDATE EQUIPMENT FOR FURTHER ANALYSIS	2-7
	2.4.1 Equipment Selection Approach	2-9
	2.4.1.1 Equipment-Type Comparison	2-11
	2.4.1.2 Specific Package Comparison	2-16
2.5	CONCLUSIONS FROM EQUIPMENT SELECTION ACTIVITY	2-19
3.0	ANALYSIS OF SELECTED EQUIPMENT	3-1
3.1	INTRODUCTION	3-1
3.2	SUITABILITY ANALYSIS	3-2
	3.2.1 Approach	3-2
	3.2.2 Spacelab/Experiment Equipment Interface Requirements (SEEIR Document)	3-6
	3.2.2.1 Document Use	3-6
	3.2.2.2 SEEIR Document Basic Philosophy	3-6
	3.2.2.3 General SEEIR Contents and Organization	3-7
	3.2.2.4 Explanation of Specified Requirements	3-8
	3.2.3 Results of Analysis	3-17
	3.2.3.1 Examples of Suitable and Unsuitable Equipment Features	3-17
	3.2.3.2 Modification Action Summary	3-37
	3.2.4 Suitability Analysis Conclusions	3-38
3.3	EQUIPMENT DESIGN MODIFICATIONS	3-38
	3.3.1 Standard Design Guidelines	3-39
	3.3.1.1 ARINC To-NIM Conversion	3-39
	3.3.1.2 Cable Clamps and Tie-Wraps	3-40
	3.3.1.3 Card Guides - Printed Circuit Boards	3-40

3.3.1.4	Component Clips	3-40
3.3.1.5	Fasteners (Vibration-Proof)	3-41
3.3.1.6	Glass Viewing Covers	3-41
3.3.1.7	Material/Hardware Replacement	3-41
3.3.1.8	Non-Rack-Mounted Instruments	3-41
3.3.1.9	Mechanical Securing of Rack-Mounted Instruments	3-41
3.3.1.10	Printed Circuit Board Securing	3-41
3.3.1.11	Protrusion Protection	3-44
3.3.1.12	Venting of Chassis, Covers, or Bulkheads	3-45
3.3.1.13	Modification Certification	3-45
3.3.2	Summary of Recommended Modifications	3-45
3.3.2.1	Staking	3-45
3.3.2.2	Circuit Board and Parts Stiffening	3-47
3.3.2.3	Other Supports	3-47
3.3.2.4	Electromagnetic Interference (EMI)	3-47
3.3.2.5	Material Replacement	3-47
3.3.2.6	Zero-G Compatibility	3-48
3.3.2.7	Item Mounting and Integrity	3-49
3.3.2.8	Audible Noise	3-49
3.3.2.9	Particle Containment	3-49
3.3.3	Summary of Recommended Processes	3-50
3.3.3.1	Flame Retardent	3-50
3.3.3.2	Vibration Encapsulation	3-51
3.3.3.3	Overspray (Painted Surfaces)	3-51
3.3.3.4	Expelling Trace Contaminants	3-51
3.3.3.5	Testing of Modified Hardware	3-52
3.4	COST ANALYSIS	3-52
3.4.1	Cost Analysis Approach	3-52
3.4.1.1	Modified Equipment Costs	3-53
3.4.1.2	Custom-Built Equipment	3-58
3.4.2	Modification Cost Summary	3-60
3.4.3	Cost of Custom Equipment Summary	3-65
3.4.4	Comparison of Modified Hardware Cost to Custom-Built Hardware Cost	3-67
3.5	SELECTED EQUIPMENT ANALYSIS CONCLUSIONS	3-70
4.0	SPECIFICATION ANALYSIS	
4.1	INTRODUCTION	4-1
4.2	NASA SPECIFICATION VERSUS INTERFACE SPECIFICATION	4-1
4.2.1	Areas of Disagreement	4-3
4.2.1.1	Requirements Increasing Engineering Costs	4-3
4.2.1.2	Requirements Causing Additional Equipment Modification	4-5
4.2.2	Impact of EC006M00000A on Selected Equipment	4-7
4.2.2.1	Modifications	4-7
4.2.2.2	Delta Modification Cost Analysis	4-10
4.3	ALTERNATIVE SPECIFICATION REQUIREMENTS	4-14
4.3.1	CAM Equipment Rack Vibration Isolation	4-14
4.3.1.1	Objective	4-14



	4.3.1.2	Comparison of Alternatives	4-15
	4.3.1.3	Conclusions	4-16
4.3.2		Material Control Tradeoff	4-16
	4.3.2.1	Objective	4-16
	4.3.2.2	Discussion of Results	4-16
	4.3.2.3	Conclusion	4-18
4.3.3		Moisture and Fungus and Corrosion Requirements Tradeoff	4-19
	4.3.3.1	Objective	4-19
	4.3.3.2	Tradeoff Conclusions	4-19
4.3.4		Power Distribution for Experiment Equipment Study	4-19
	4.3.4.1	Objective	4-19
	4.3.4.2	Discussion of Results	4-20
	4.3.4.3	Conclusions	4-23
5.0		MODIFICATION SOURCE TRADEOFF	5-1
5.1		GENERAL APPROACH TO THE TRADEOFF	5-1
5.2		EQUIPMENT MODIFICATION PROCESS FLOWS	5-2
5.3		FACILITIES AND OTHER RESOURCE REQUIREMENTS	5-7
5.4		EVALUATION CRITERIA ANALYSIS	5-11
	5.4.1	Modification Costs	5-11
	5.4.2	Modification Time	5-13
	5.4.3	Modification Ability	5-13
	5.4.4	Modified Equipment Performance Risk	5-15
5.5		RECOMMENDED APPROACH	5-15
6.0		EQUIPMENT RECOMMENDED FOR PROCUREMENT	6-1
6.1		CHARACTERISTICS OF EQUIPMENT SELECTED FOR PROCUREMENT	6-1
6.2		TECHNICAL ASSESSMENT OF EQUIPMENT TYPES	6-2
6.3		PRELIMINARY EQUIPMENT-TYPE ASSESSMENT	6-2
6.4		PROCUREMENT RECOMMENDATIONS	6-6
7.0		PROGRAMMATIC SAVINGS	7-1
7.1		PROGRAMMATIC EXTENSION SCENARIOS	7-1
7.2		SAVINGS PROJECTIONS	7-2
8.0		CONCLUSIONS	8-1

ILLUSTRATIONS

Figure		Page
1-1	Overall Study Flow	1-3
2-1	Commercial Equipment Identification and Selection Logic .	2-1
2-2	Functional Requirements Sorted According to Science and Application Disciplines	2-4
2-3	Function Modification Summary	2-6
2-4	Sample from Commercial Equipment Shopping List	2-8
2-5	Equipment Selection Approach	2-10
2-6	Estimation of Relative Equipment Procurement Costs	2-12
2-7	Equipment-Type Primary Operating Principle	2-15
2-8	Discipline Representation	2-15
2-9	Design Specification Representation	2-18
2-10	Package Installation Configuration Representation	2-18
2-11	State-Of-The-Art Representation	2-18
3-1	Logic Flow Diagram for the Analysis of Selected Equipment .	3-1
3-2	Typical Evaluation Summary Sheet	3-4
3-3	Example Copy Of A Typical SEEIR Modification Summary Sheet	3-5
3-4	Interface Specification Description Key Contents and Organization	3-7
3-5	Shield Effectiveness Versus Frequency	3-12
3-6	Collins Transceiver	3-18
3-7	Collins Transceiver - Chassis View	3-20
3-8	Hewlett-Packard Spectrum Analyzer	3-21
3-9	Interior View of H/P Spectrum Analyzer Display Section .	3-22
3-10	H/P Spectrum Analyzer - RF Section	3-23
3-11	Interior View of H/P Spectrum Analyzer - IF Section	3-25
3-12	Interior View of H/P Spectrum Analyzer Automatic Preselctor	3-26
3-13	NIM Bin With Modules - Frontal View	3-27
3-14	NIM Bin - Rear View Showing Power Supply	3-28
3-15	Tennelec NIMs - EMI Covers Removed	3-29
3-16	Beckman Gas Chromatograph	3-31
3-17	Beckman pH Meter	3-31
3-18	Sylvania CO ₂ Laser Interior	3-32
3-19	Research, Inc. Display Terminal	3-33
3-20	Research, Inc. Display Terminal - Interior View	3-35
3-21	Coulter Blood Cell Counter	3-36
3-22	ARINC-To-NIM Conversion	3-39
3-23	Card Guide, Printed Circuit - Positive Securing	3-40
3-24	Component Clips	3-40
3-25	Mechanical Securing of Heavy Rack-Mounted Instruments	3-42
3-26	Mechanical Securing - Printed Circuit Boards	3-43
3-27	Protrusion Protection	3-44
3-28	Cost Analysis Approach	3-52
3-29	Work Breakdown Structure	3-53
3-30	Engineering Cost Estimate (Modified Equipment)	3-55
3-31	Sample Custom Cost Estimate (Tennelec Counter/Timer)	3-59
3-32	Modification Versus Retail Cost	3-62



3-33	Design Specification Impact on Modification Costs . . .	3-63
3-34	Gravity Dependence Impact on Modification Costs . . .	3-64
3-35	Operating Principle Impact on Modification Costs . . .	3-65
3-36	Modified Versus Custom Cost	3-69
3-37	Electronic Instrumentation Most Cost Effective to Modify .	3-70
4-1	Specification Comparison - Approach	4-2
4-2	EC006M00000A Versus SEEIR Modification Costs	4-13
4-3	Modified Versus Custom Costs	4-13
4-4	Risk Analysis for Relaxing Material Control Requirements from NASA Spec to Study Spec	4-17
4-5	Equivalent Circuits for Power Distribution Alternatives .	4-21
5-1	Original Equipment Manufacturers' Modification Process Flow	5-3
5-2	Process Flows for Centralized Contractor and NASA . . .	5-5
5-3	Sample of Commercial Equipment Questionnaire	5-8
5-4	Rate of Modified Equipment Acquisition	5-14
7-1	Programmatic Extension Scenarios	7-2
7-2	Programmatic Cost Comparison	7-4

TABLES

Table		Page
2-1	Functional Requirements Data Bank Source Documents	2-3
2-2	Functional Modification Category Definitions	2-6
2-3	Equipment-Type Selection	2-14
2-4	Selected Equipment Representation Summary	2-17
3-1	Summary of Activities For Each Selected Item	3-3
3-2	Calculated Shield Effectivity Of 0.040-Inch (0.1 cm) Sheet Iron (Copper)	3-13
3-3	Modification Actions Summary	3-37
3-4	Detailed Summary of Item Modifications	3-46
3-5	Manufacturing Cost Estimate (Modified Equipment)	3-56
3-6	Modified Tennelec Timer Coscs	3-58
3-7	State-Of-The-Art Rationale	3-59
3-8	Modified Equipment Cost	3-61
3-9	New Hardware Development Cost	3-66
3-10	Comparison of Modification Cost to Custom-Built Cost	3-68
4-1	NASA and Study Specification Origins	4-2
4-2	Delta EC006M Requirements Over SEEIR	4-4
4-3	EC006M Delta Modification Summary	4-8
4-4	Delta Modification Cost Summary	4-12
4-5	Delta Weight Summary	4-22
4-6	Delta Cost Comparison of AC and DC Power Distribution	4-24
5-1	Results of Equipment Manufacturer Survey	5-9
5-2	Equipment Modification Trade	5-16
6-1	Equipment Listed by Technical Importance	6-3
6-2	Equipment Types Ranked According to Spacelab Experiment Program Savings	6-4
6-3	Equipment-Type Rejection Rationale	6-5
6-4	Procurement Recommendations	6-6

1.0 INTRODUCTION

To date, the major portion of the cost of space hardware has been incurred by the design, development, and qualification of custom-built equipment. This equipment is typically produced to meet stringent specifications associated with the high reliability required and the environment encountered on infrequent space missions using expendable launch vehicles. The Shuttle experiment program will be different in the following ways: (1) flights will be more frequent, (2) the Spacelab environment will be earth-like, and (3) on-board equipment will be returned to earth after each flight. These differences, together with the presence of man, increase the potential for use of commercial equipment in the Spacelab.

The use of available commercial, airborne or military (CAM) equipment could result in significant reduction in sortie payload costs, thus allowing many more payloads to be flown. This will be a particular benefit early in the Shuttle program. In addition to the cost savings provided by using CAM equipment, it is anticipated other benefits will result such as shortening equipment procurement lead time and providing flexibility for the principal investigator to employ some of the same equipment in space as used in his ground-based laboratory. These benefits have the potential of significantly improving the value of the space program.

The importance of these benefits has been recognized in the planning for Spacelab which includes the design objective that standard laboratory equipment be accommodated without extensive modification. This study is particularly timely since a better understanding of the role of CAM equipment is needed soon to effectively implement the above objectives and to plan long lead-time payload developments.

This study examines a representative group of experiment equipment and develops quantitative cost tradeoff data to extend earlier feasibility assessments. Data and recommendations have been generated to assist in developing NASA specifications that do not unnecessarily restrict the use of CAM equipment in the Spacelab. The study further recommends equipment for NASA procurement for testing and identifies the preferred source of equipment modifications.

1.1 STUDY OBJECTIVES

The objective of this study is to investigate analytically the feasibility of using commercially available laboratory equipment and instrumentation in the Spacelab in support of various experiments. "Commercially available" has been interpreted to mean any equipment identified in a catalog including commercial, airborne and military (CAM) equipments. Feasibility is to be demonstrated by the breadth of application of CAM equipment to experiment equipment requirements in the Spacelab, and the cost effectiveness of utilizing this class of equipment instead of custom-built aerospace equipment typical of past designs.

Six specific objectives were satisfied in the analysis of the feasibility of CAM equipment:

1. Compilation of a list of representative equipment and instrumentation meeting the requirements of Spacelab experiments. This list is indicative of the available hardware that could be flown aboard the Spacelab.
2. Selection of a limited number of these items for analysis. These items were representative of the full equipment list to assure coverage of as many design variations as possible.
3. Collection of engineering specifications and data relative to these items. A library of manuals and specifications has been accumulated for the selected equipment.
4. Analysis of these data with respect to applicability for use in the Spacelab. The suitability of the selected equipment to space operations was assessed from specifications and visual examination. Modifications were identified where necessary and the cost of the modified hardware was compared to similar custom-built hardware.
5. Tradeoff analyses of specification requirements with respect to cost. Alternative requirements, compatible with the use of commercial hardware, were evaluated and specification changes were recommended which will reduce Spacelab equipment costs.
6. Recommendations for procurement of specific items for more extensive future test and analyses.

1.2 OVERALL APPROACH

This section presents a task-by-task discussion of the technical approach to the study. The study approach is discussed with reference to the study flow diagram of Figure 1-1.

Task 1. Identification and Selection of Potential Commercial Equipment and Instrumentation for Analysis

The Shuttle System Payload Data Activity (SSPDA) and other studies have identified types of instruments and equipment required to accomplish Spacelab experiments. Candidate commercially available equipment and instrumentation were identified, the degree of modification required to meet experiment functional requirements was assessed for each of them, and items were selected for detailed analyses throughout the study.

A listing of Spacelab instruments and equipment items and their functional requirements was developed, based on SSPDA data and other studies. Commonality of instrument usage across various disciplines was assessed. The instrument and equipment items identified were then categorized by their commercial

availability and extent of modification required to meet experiment functional requirements. Vendor catalogs were searched to identify available equipment meeting these requirements. In parallel with these activities a selection logic was developed to identify items representing a cross-section of Spacelab experiment equipment for further analysis in the study. A limited number of equipment items were selected, based on this selection logic, for analysis in Task 2.




Figure 1-1. Overall Study Flow

Task 2. Analysis of Selected Items of Equipment and Instrumentation

The most cost-effective approach for Spacelab equipment and instrumentation development was determined from a cost analysis which compared commercially available equipment modified to be suitable for space flight with equipment custom built for the Spacelab.

The suitability of commercially available equipment for installation and operation in the Spacelab was determined first. Physical and performance data were collected for each of the pieces of selected equipment. This activity spanned a number of weeks and occurred in parallel with all other Task 2 sub-tasks because of the time required to obtain comprehensive information. Available information from the Spacelab and Shuttle programs provided the basis for the development of a Spacelab/Equipment Interface Requirements document which served as the specification for all equipment to be placed in the Spacelab. Equipment characteristics were compared with this interface requirements document to identify discrepancies which make the existing equipment unsuitable for space flight.

The cost of equipment modifications to eliminate the previously identified discrepancies were compared with custom-built equipment which meet space flight requirements. Modifications required by the suitability analysis were developed



conceptually. Descriptions of these modifications provided information to manufacturing planning and estimating personnel for determination of the shop cost for the modification. Supplemental costs to be added to the shop costs were identified for a complete determination of the cost of modification. In a parallel activity, the cost of custom building the equipment was determined. The two costs were compared and the cost-effective approach identified.

Task 3. Compare the Specifications of the (Modified) Experiment Equipment and Instrumentation with NASA Specification

The objectives of this task were to (1) evaluate the impact of imposing NASA Specification EC006M00000A requirements on the selected commercial equipment requirements, and (2) define and evaluate feasible alternative approaches to the accommodation of commercial equipment in the Spacelab. Task 2 established a set of Spacelab requirements and the needed item-by-item modifications and costs for compliance. Task 3 established the added cost for further modification of design, test, and processing for EC006M00000A compliance and then determined how the feasibility of using commercial hardware in space could be improved by changing this specification.

The initial activity consisted of a comparison of the EC006 and the study specifications. Delta costs that would accrue if the NASA specification were enforced were derived. Modified user equipment capabilities/requirements were compared to EC006M00000A for the purpose of identifying additional discrepancies to those identified in Task 2. Having identified the differences between modified equipment characteristics and EC006 specification requirements, the hardware and process changes to the modified items necessary to meet EC006 were identified, costed, and evaluated using the same procedures as in Task 2.

The characteristics of unmodified equipment were assembled and compared with corresponding requirements in the EC specification to identify requirements that must be cost effectively relaxed. Alternative Spacelab equipment specification requirements and the resultant system impacts were established. The advantages and disadvantages (from cost, performance, reliability, and safety points of view) of adopting each alternative requirement were assessed.

Task 4. Develop Recommendations

Task 4 integrated the outputs from all the previous study tasks to: (1) assess the overall feasibility of using commercial equipment in the Spacelab; (2) make specification change recommendations; and (3) develop a list of equipment for NASA procurement. In addition, a major trade was performed to determine the preferred source of the necessary equipment modifications.

Candidate specification changes were evaluated with respect to the entire Spacelab experiment program and change recommendations were made which would make the use of commercial equipment more cost effective. The modification source trade was performed defining the most cost-effective agency to perform the prescribed modifications. These outputs, as well as Tasks 2 and 3 results, were integrated to produce an overall feasibility assessment of the applicability of commercial hardware to the Spacelab program. The principal study findings were then used to select a list of hardware items for NASA procurement.

1.3 SUMMARY

The final report of this study consists of three volumes:

Volume I Executive Summary
Volume II Technical Report
Volume III Design Analyses and Trade Studies

Volume II presents the discussion of technical results of the study. It consists of two parts. The first part contains the technical discussion, and the second contains five appendixes to the discussion. This discussion is organized similar to the study task breakdown. A brief summary of the results discussed in each section follows.

1.3.1 Selection of Equipment for Analysis

Section 2 of this volume describes the approach employed to select specific available hardware items for analysis. The major outputs of this effort were a Spacelab experiment functional requirements data bank, a shopping list of available hardware that meets these functional requirements, and a list of equipment items selected for further analysis. The functional requirements data bank consists of a tabulation of equipment requirements for over 600 equipment items. These items were identified in SSPDA and other specific sortie payload studies. Functional modification classifications are also documented on this tabulation. The shopping list consists of a separate tabulation of approximately 500 equipment items which satisfy the requirements established in the functional requirements data bank. This tabulation and the data bank tabulation have been included in the appendixes. Thirty-four packages, consisting of 43 separate items, were selected from the shopping list for further analysis. The selection criteria emphasized higher cost equipment items, assuming that they offered the greatest potential for savings through the use of available hardware. The resultant group was also representative of the various types of equipment expected for Spacelab experiment operations.

1.3.2 Analysis of Selected Equipment

The suitability of the selected equipment items for installation in the Spacelab and the cost effectiveness of modifying available CAM equipment are presented in Section 3.0 of this volume. The major outputs of this effort were the Spacelab/Experiment Equipment Interface Requirements (SEEIR) specification, definition of the modifications required to make the selected equipment suitable for installation and operation in the Spacelab, and a comparison of the cost of this modified equipment to the cost of the same hardware custom built.

The SEEIR represents the first equipment specification derived with the use of commercial equipment in space in mind. It does not contain many of the typical requirements included in previous manned spacecraft specifications. The specification appears in Appendix D of this volume.



Unsuitable characteristics of the 34 selected equipment items were defined from visual inspection and a comparison of documented specifications for each item. The internal construction of 29 of 34 of the selected items was visually examined by the study team to assess structural integrity for acceleration and vibration loading, adequacy of cooling provisions, and the acceptability of materials used. Modifications ranged from stiffening printed circuit boards and replacement of unsecured screws to survive vibration, to a total redesign of gravity-dependent functions. Three units--a Coulter blood cell counter, a refrigerator/freezer and a dewar--required such extensive modification that they could only be custom built.

The comparison of modified and custom-built hardware costs indicates significant savings can be achieved by modifying available hardware. The average cost to custom-build the 34 selected items was estimated to be \$225,000. The average cost of similar available hardware including all modifications is \$31,000, an average savings of approximately 85 percent of the custom-built hardware cost. Detailed descriptions of the suitability analysis, modification details, and cost estimates are presented for each equipment item in Section 2.0 of Volume III.

1.3.3 Analysis of Specifications

Section 4.0 presents the results of the comparison of the study-developed interface specification (SEEIR) and a typical NASA equipment specification consistent with the previous manned spacecraft practice. Major output of the effort is the determination of the cost impact of imposing a rigorous manned spacecraft-type specification on the selected equipment items and identification of approaches to relaxing specifications to make them more suitable to the use of available CAM equipment.

The impact of the NASA specification (EC006M00000A) upon the selected hardware is two-fold. First, additional modifications to the selected hardware are required to meet more rigorous criteria such as worst case temperature limits. Secondly, additional systems engineering time is required to document additional analytical and test results required by the EC specification. Delta costs incurred by these additional activities averaged \$30,000 per unit. Descriptions of the delta modifications and costs are presented for each equipment item in Section 2.0 of Volume III.

Four tradeoffs are also presented in Section 4.0. These tradeoffs examine relaxation of the equipment specification. Approaches to relaxing the equipment vibration environment, reducing material control restraints, and eliminating fungus and corrosion requirements are examined. A fourth tradeoff was performed to determine whether elimination of the power supply sections of most equipment is feasible and cost effective. The results of the tradeoffs indicate possible cost savings with little increase in risk if specifications are relaxed for corrosion and fungus. Material control can be relaxed to be consistent with that defined in the Spacelab/Experiment Equipment Interface Requirements document without reducing the safety integrity of the Spacelab. Reduction of the vibration environment to a level eliminating all equipment modifications defined for vibration integrity could not be cost effectively accomplished. Finally,

elimination of the power supply sections to be compatible with dc power resulted in the addition of dc-to-dc converter sections to each unit which show no performance improvement over the original configuration but do cause increased modification costs.

1.3.4 Modification Source Tradeoff

The cost-effective modification agency was determined by a tradeoff analysis of three alternative modification sources. Section 5.0 compares the advantages and disadvantages of the original equipment manufacturer, a centralized industrial contractor, or NASA performing equipment modifications. The analysis of the three alternatives consisted of an assessment of the differences between hardware process flows through each facility, polling of a number of manufacturers to determine their capability and willingness to perform the types of modifications identified by this study, and a review of the production facilities of two instrumentation manufacturers--Beckman Instruments and Hewlett-Packard.

The results of the assessment show that the original equipment manufacturer is the most cost-effective modifier of its own equipment. This conclusion was based on the immediate availability of equipment for modification, the potential for modifying equipment during the existing assembly process, and the familiarity of the manufacturer with his own equipment.

1.3.5 Recommended Equipment Procurements

Section 6.0 presents a list of eight pieces of equipment recommended for procurement by NASA-MSFC for test evaluation. The CVT program will be evaluating CAM hardware characteristics with an extensive testing program. This study has recommended the purchase of eight units in addition to those currently in the CVT inventory. The objective of the selection process was to satisfy two goals. First, equipment was identified that could provide the largest dollar savings to the Spacelab experiment program if available hardware were used. The second goal was to recommend an equipment group that would maximize the data return for the minimum expenditure of funds.

1.3.6 Conclusions

Use of CAM equipment in the Spacelab will be cost effective. Extension of the savings estimated for the 34 selected equipment items to the entire Spacelab experiment program indicates potential savings ranging from \$230 million to \$1230 million, depending upon the equipment development and replacement scenario adopted. The most probable savings amount is believed to be approximately \$400 million.

2.0 SELECTION OF EQUIPMENT FOR ANALYSIS

2.1 INTRODUCTION

The objective of the first study task was to select a small group of available hardware representative of the types of hardware that could be used for Spacelab experimentation. This section describes the activities performed in the selection of this hardware.

Figure 2-1 illustrates the logic leading up to hardware selection. Results of the Shuttle System Payload Data Activity (SSPDA) and specific detailed payload studies provided package functional performance requirements data. These functional requirements have been assembled on a single computer tabulation. Common equipment requirements were identified across disciplines by collating the inputs by equipment type. Instruments and equipment types were categorized according to their availability and extent of modification required to meet experiment functional requirements. Catalogs from various suppliers were consulted to identify available equipment that meet the defined functional requirements. These equipments are tabulated on IBM cards forming a shopping list from which equipment items for further analysis were selected. A selection criterion stressing total procurement cost and representation of the design characteristics of available equipment was employed to select items from the shopping list.




Figure 2-1. Commercial Equipment Identification and Selection Logic

Thirty-four available packages were selected. In a few cases, the selected packages consisted of more than one hardware element, leading to a total of 43 specific pieces of hardware for analysis. The suitability to space operations and the modifications required to assure operation on board the Spacelab are discussed in Section 3.0.

2.2 FUNCTIONAL REQUIREMENTS DATA BANK

Functional requirements have been defined for 615 equipment items, located in the Spacelab, identified in the following science and application disciplines: Astronomy, Solar Physics, Atmospheric and Space Physics, Earth Observations, Earth and Ocean Physics, Space Processing Applications, Life Sciences, Space Technology, and Communication and Navigation. Four hundred forty-four different types of equipment were identified from the 615 data bank entries. These data have been recorded on IBM cards to allow sorting by equipment type and to facilitate any updating or revision of this information. The resultant tabulation is presented in Appendix A of this volume.

Many Shuttle sortie mode payload studies have documented the experiment equipment performance requirements necessary to achieve experiment objectives. The Shuttle System Payload Data Activity (SSPDA) has identified many of these equipment items. Information from SSPDA, supplemented by the data from individual payload studies, was gathered, forming a functional requirements data bank for use in this study.

The source information used to generate the data bank is shown in Table 2-1. Specific payload studies generally provided additional data to the SSPDA tabulation. In many cases, the SSPDA data did not report information at the package level, resulting in the use of information from these payload studies as the primary source of data. For example, in the Space Processing and Life Science disciplines, experiment modules containing families of equipment are identified on the SSPDA tabulation. In these cases, equipment items and their functional requirements were obtained entirely from the payload study. The Manned Earth Observatory (MEO) study posed an additional problem in that this sortie payload is no longer identified on the Shuttle mission model. SSPDA tabulations reference this study in both the Earth Observations and the Earth and Ocean Physics disciplines; however, the experiment objectives are not the same. Since the intent of this study is to examine hardware representative of all types of sortie payloads it was decided to include information from this study on the tabulation. The MEO information is included under the Earth Observations discipline.

Hardware identified in the functional requirements data bank are the only kinds of interest to the study. Spacelab subsystem hardware has been excluded. All hardware is located in the Spacelab. Applicability of commercial hardware to pallet installation exposed directly to the space environment appears limited. Therefore, pallet-mounted hardware was excluded. Consumables such as film, tape reels, etc., and low-cost components identified in some studies are also not included to keep the tabulation manageable.

Table 2-1. Functional Requirements Data Bank Source Documents

1. NAS8-29462, Shuttle System Payload Data Activity (GDCA/Rockwell/GE)
2. Payload Element Data Sheets and Final Report of the Space Shuttle Payload Planning Working Groups (NASA-GSFC)
3. NAS8-28938, Requirements and Concepts for Materials Science and Manufacturing in Space Payload Equipment Study (TRW)
4. NAS8-29150, Life Science Payload Definition and Integration Study
5. Advanced Technology Laboratory Phase A Study Final Report (NASA-Langley)
6. NAS1-11674, Manned Activity Scheduling Study Final Report (GDCA)
7. NAS8-27540, Definition of Experiments and Instruments for a Communication/Navigation Research Laboratory (TRW)
8. NAS8-28144, Astronomy Sortie Mission Study (MMC)
9. NAS8-28047, Plasma Physics and Environmental Perturbation Laboratory (TRW)
10. NAS9-12255, Preliminary Design Study for an Atmospheric Science Facility (MMC)
11. NAS9-12649, Scientific Design of a Shuttle Auroral Observatory System (U. of Alaska)
12. NAS8-28013, Mission Requirements for a Manned Earth Observatory (TRW)
13. NAS8-27861, Zero-Gravity Cloud Physics Laboratory (MDAC)

Functional requirements for each item were obtained from the referenced literature when available. If the data were not available, experts were consulted by the study team and their best judgment was used to define appropriate functional requirements. Commercial hardware could not have been identified for many pieces of experiment equipment if this had not been done. As experiments become better defined, functional requirements can be updated and applicable equipment can be identified.

2.2.1 Explanation of Tabulation

A sample from the functional requirements data bank tabulation is shown in Figure 2-2. The listing is organized as follows. Each piece of equipment is described by a title card, a narrative description of the function of that piece of equipment, and a listing of its key functional requirements. The title card contains additional information as shown by the headings on the figure. Numbers in the left-hand column identify the equipment-type category. This categorization is consistent with that defined by Science Magazine in its "Guide to Scientific Instrumentation, 1972-73." The title of the piece of equipment is followed by its SSPDA equipment identification number (i.e., SP001). The next group of numbers identifies the quantity of that item on the payload, its location (I, inside Spacelab; E, external to Spacelab), the state of the art of this hardware (1, current; 2, advanced), the source of the hardware (1, commercial; 2, aircraft; 3, manned space; and 4, unmanned space), and finally the functional modification required by the available hardware (1, none; 2, minor; 3, major; 4, assembled from piece parts; and 5, new development). The final columns sequentially order the data by scientific discipline, SSPDA payload, and equipment item.

		FUNCTIONAL MODIFICATIONS				DISPLAY SEQ.	
		SOURCE					
		STATE OF ART					
		LOCATION					
		QUANTITY					
EQUIP TYPE	EQUIPMENT TITLE	SSPDA EQUIP IDENTIFICATION					
186000	HOT WALL TURE FURNACE	SP002	31111	60230101			TITLE CARD
186000				60230102			
186000	THIS UNIT IS A GENERAL PURPOSE, HOT WALL HEATING DEVICE			60230103			NARRATIVE
186000	PROVIDING ACCURATE CONTROL OVER THE HOT ZONE TEMPERATURE WITH			60230104			
186000	RESPECT TO A FLAT PROFILE OR TO A SPECIFIED GRADIENT PROFILE.			60230105			
186000				60230106			
186000	VACUUM LEVEL.....2E-8 PSI (1.95E-4N/M2)			60230107			FUNCTIONAL REQUIREMENTS
186000	OPERATING TEMPERATURE.....2200F (1200C)			60230108			
186000	MAX EXTERNAL SURFACE TEMP.....110F (45C)			60230109			
186000	HOT ZONE DIMENSIONS.....2.5 IN (6CM) DIA. X 6IN(15CM) LONG			60230110			
186000	MAXIMUM DEVIATION FROM FLAT			60230111			
186000	PROFILE.....+/-4F (+/-2C)			60230112			
186000				60230113			
186000	HOT WALL FURNACE	SP002	21111	60230201			TITLE CARD
186000				60230202			
186000	THIS UNIT IS A GENERAL PURPOSE, HOT WALL HEATING DEVICE PROVIDING			60230203			NARRATIVE
186000	ACCURATE CONTROL OVER THE HOT ZONE TEMPERATURE WITH RESPECT TO			60230204			
186000	A FLAT PROFILE OR TO A SPECIFIED GRADIENT PROFILE.			60230205			
186000				60230206			
186000	VACUUM LEVEL.....2E-8 PSI (1.5E-4N/M2)			60230207			FUNCTIONAL REQUIREMENTS
186000	OPERATING TEMPERATURE.....3300F (1800C)			60230208			
186000	EXTERNAL SURFACE TEMPERATURE..110F (45C)			60230209			
186000	HOT ZONE DIMENSIONS(MIN).....1IN. (2.5CM) DIA. X 5IN. (13CM) LONG			60230210			
186000	MAXIMUM DEVIATION FROM FLAT			60230211			
186000	PROFILE.....+/-7F (+/-4C)			60230212			
186000	VIEWSPTS.....2			60230213			
186000				60230214			

Figure 2-2. Functional Requirements Sorted According to Science and Application Disciplines

Identification of each equipment item included its respective SSPDA identification number where possible. Exceptions to this practice were required whenever additional data were found in the supporting studies. Only external equipment is identified in the SSPDA tabulation in the Astronomy, High-Energy Astrophysics and Solar Physics disciplines. The Astronomy Module study by Martin was used as a model for the internal equipment of these three disciplines. This study indicated the basic support hardware required by the display and control console in the Spacelab. Study of these disciplines indicates that most of the experiment hardware is pallet mounted, and very little additional hardware except a control console would be located in the Spacelab. Listing of the Zero-G Cloud Physics Laboratory hardware is consistent with SSPDA, but the MEO is not (for reasons explained earlier). This hardware is clearly indicated by an MEO rather than an SSPDA part number in the listing. Space processing numbers relate to the SSPDA modular assembly that includes the subject package. In cases where a package is used in more than one module, only the first module was indicated. However, the total quantity was adjusted to include the items identified in the other modules. Changeout of hardware from one module to another was assumed not to be cost effective. Finally, Life Science and Communication/Navigation part numbers that were taken from the sortie payload reports were given 900 (or higher) series numbers.

2.2.2 Functional Modifications

Commercially available hardware may require modification to achieve an experiment performance objective, or it may require modification to be capable of operation on board the Spacelab. The study primarily addresses this latter type of incompatibility (Section 3.0). However, a gross assessment of the capability of available hardware to meet the functional requirements of Shuttle-sortie mode experiments was accomplished during this study task. This assessment was based on a similar assessment performed during the SSPDA study, available data from commercial equipment catalogs, sortie payload studies, and the study team's best judgment. The results of the functional modification assessment appear in the functional requirements data bank.

The functional modification categories are explained in Table 2-2. Figure 2-3 summarizes the results of the functional modification assessment. It is significant that the functional requirements of 282 out of the 444 different items identified could be satisfied by available hardware. Of the 282 items, 247 (55 percent of the total) did not require any modification to meet these requirements. These results show good correlation with experience in the ASSESS program.¹ The applicability of available equipment results from two factors. First, much of the equipment inside the Spacelab is experiment support hardware. Hardware such as oscilloscopes, frequency generators and gas chromatographs exist with diverse capabilities. Specialized hardware unique to space exploration is more frequently located on the Spacelab pallet than in the pressurized module. Secondly, many of the experiments are similar to those performed on the ground, differing only in that they are performed in a space environment. Since commercial laboratory instrumentation is typically used in the ground experiments, it is not surprising that the equipment would meet defined functional requirements.

¹NASA TMX 62,287. Study of Airborne Science Experiment Management Concepts for Application to Shuttle, Volume 2. Donald R. Mulholland, et al. Airborne Science Office, NASA-Ames Research Center, Moffett Field, California; July 1973.

Table 2-2. Functional Modification Category Definitions

Category	Definition
None	Package can accomplish all functional requirements in the laboratory.
Minor	All functional requirements can be met except one or two. Only minor modification need be done to meet these requirements. No advancement in the state of the art is required.
Major	Primary functional requirements can be met but many secondary requirements cannot, or more than two functional requirements require modification but basic unit is still applicable. No advance in state of the art is required.
Assembled from available piece parts	Required package does not exist in a catalog. No advancement in the state of the art is required. Package can be made from existing piece parts.
New development	Commercial package not available; requires advancement in state of the art.




Figure 2-3. Function Modification Summary

The remaining segments on the figure show the number of items that must be custom built. Units which must be custom built, but are not new developments, are included in the piece-part segment. Equipment items that fell into this category were assemblies and control electronics. This segment may be distorted by 63 control packages identified for the Ames payload. This quantity might be reduced if a central processor and data bus concept were employed. The percentage of available equipment meeting functional requirements without modification would then increase to 65 percent. Items categorized as new developments include items which have been identified as SRT items in the sortie payload studies such as the particle counters identified in the Zero-G Cloud Physics Laboratory or are typically custom-built units such as cloud chambers in the same payload, or combustion chambers employed in the Advanced Technology Laboratory.

2.3 SHOPPING LIST

Available CAM equipment which satisfies the functional requirements tabulated in the data bank has been identified by reviewing manufacturer catalogs. Over 425 different items have been identified. A tabulation of these items has been included in Appendix B. Manufacturer, part number, current cost, and key performance capabilities are listed for each item. The listing is representative of the equipment available today and is by no means exhaustive. It includes a large enough quantity of equipment to assure selection of representative hardware for further analysis in the study.

2.3.1 Explanation of Tabulation

A sample of the shopping list tabulation is shown in Figure 2-4. The equipment type is first identified. Then, the requirements in the functional data bank are tabulated. (These requirements are not shown on the figure.) Next, the applicable commercial items are identified by part number and cost for each equipment type. For example, the Astro Model A222 (costing \$1000), the Lindberg Model F8417 (costing \$355) and the Leco Model 521-300 (costing \$1370) hot wall tube furnaces satisfy the requirements for this type of furnace defined in the Space Processing sortie experiment payloads. The Astro Model 1100V (costing \$3100) satisfies the requirements for the hot wall furnace in this discipline. The three numbers following the price indicate the specification to which the hardware is built (1, commercial; 2, aircraft, 3, military; 4, unmanned space; and 5, manned space), the technology status (1, current; 2, advanced), and the modification category (same as explained in the functional requirements data bank). For example, the Lindberg and Leco furnaces require major modification (Category 3) to achieve a vacuum in the oven. Significant performance capabilities for each identified piece of hardware are listed below the identifier cards. Numbers in the right-hand columns are for sequencing.

2.4 SELECTION OF CANDIDATE EQUIPMENT FOR FURTHER ANALYSIS

A criterion for selection of hardware for further analysis, from the hundred of items in the shopping list, was developed. This criterion was used to select the items for in-depth analysis to determine their suitability for space operations and to identify the modifications necessary to make them suitable.

ORIGINAL PAGE IS
OF POOR QUALITY

2-8

HOT WALL TUBE FURNACE	PART NUMBER	FUNCTIONAL MODIFICATION		
		STATE OF THE ART		
		SPECIFICATION		
		COST		
ASTRO INDUSTRIES INC	#A222	\$1000		111
HOT ZONE DIMENSIONS		5 CM DIA. BY 28 CM LONG		
MAX SUSTAINED TEMP		1200 DEG C		
LINDBERG HEVI-DUTY	#F8417	\$355		113
MAX OPERATING TEMPERATURE		2200F (1200C)		
CHAMBER SIZE		1IN I.D.X 12IN LONG		
POWER		970W		
LECO	#521-300	\$1370		113
INDUCTION FURNACE				
TEMPERATURE (MAX)		>3000F (>1700C)		
POWER		1.5KW		

HOT WALL FURNACE

SP-02 302 2

ASTRO INDUSTRIES INC	#1100V - 1060-M1	\$3100	111
HOT ZONE DIMENSIONS	2.5 CM DIA. BY 15 CM LONG		
MAX SUSTAINED TEMP	1800 DEG C		
TIME TO TEMPERATURE	10 MINUTES MAX		
HEATING ELEMENT	MOLYBDENUM MESH WITH MOLYBDENUM SHEET RADIATION SHIELDS		
OPERATION ENVIRONMENTS	VACUUM, INERT GAS, AND REDUCING ATMOSPHERES		
OBSERVATION	SHUTTERED RADIAL VIEWPORT		
ACCESS	TOP AND/OR BOTTOM		
COOLING	WATER COOLED		

Figure 2-4. Sample from Commercial Equipment Shopping List

Four characteristics were considered to be significant for the equipment selected for analysis: (1) high procurement cost, (2) differing physical characteristics, (3) representation of many types of equipment, and (4) manufacture by a company willing to provide detailed data necessary for analysis.

The criterion stresses relatively more expensive experiment equipment items because they offer the greatest potential for savings to the Spacelab experiment program through the use of commercial equipment. Frequency of use and unit cost are combined in a single cost parameter to estimate a relative total program procurement cost for each item. For example, an item which costs five times more than another, but only used once, will have a smaller total procurement cost than a second item which is identified in ten different sortie payloads.

Representation of equipment built to different requirements is also desirable to assure identification of as many design issues as possible. At least one piece of equipment had to be selected from each SSPDA science and application discipline. In this way, hardware designs resulting from special functional requirements unique to a specific type of experimentation were examined. In addition to hardware built to commercial specifications, similar hardware built to aircraft and military specifications was selected to gain understanding of the design differences resulting from different specification requirements.

Equipment with differing physical characteristics was selected to assure examination of many design issues. Packages which are rack-mounted and others that are bench-mounted were selected. Integral and modular packaging techniques, such as those used by the nuclear instrumentation modules (NIM) and CAMAC, are represented. Equipment employing integrated circuits (IC's) and metal oxide on silicon-field effect transistors (MOS-FET) was identified to gain understanding of the suitability of hardware with technology more closely related to that of the Shuttle era.

Final selection of specific equipment items for analysis was based on the availability of manufacturer data. Similar studies have found that many manufacturers have not wished to provide data on their commercial equipment for either economic or competitive reasons. This selection criterion is a practical condition which assures the most complete information for the subsequent detailed analysis.

The approach employed to select hardware items for specific analysis is shown on Figure 2-5. The selection process had two phases. In the first phase equipment types were compared to determine which types would be best to analyze. During the second phase, specific available packages were identified for each selected equipment type which provided the broadest cross-section of design issues.

2.4.1 Equipment Selection Approach

Selection of the equipment for analysis was made from the 425 items in the shopping list by successively grouping and eliminating items with common characteristics. First, the list was condensed into 90 different types of hardware by grouping common equipment such as all receivers regardless of frequency band,

centrifuges regardless of speed, spectrometers and spectrophotometers regardless of spectrum band, etc. Next, these types were ranked according to relative procurement expense. Then, starting with the most expensive equipment type, the primary operating principle of successively less expensive types was compared. Less-expensive types with similar operating principles, such as electronic packages, were eliminated in favor of equipment with different operating principles such as optical or chemical analysis hardware. The representation of the resulting 25 types of equipment in each of the science and application disciplines was then validated to assure that at least one item from a payload in each discipline had been selected.




Figure 2-5. Equipment Selection Approach

The selected hardware group was increased to 34 to assure representation of different design characteristics. Additional items which were packaged as NIM (nuclear instrumentation modules) and ATR (air transport requirement or Austin Trumbull Radio) modules were added to commercial counterparts selected for a single equipment type to examine different installation approaches. Items which were representative of Mil-Spec and airborne specifications were also added to the commercially available units. Finally, selected advanced state-of-the-art hardware was added by the group.

The final group of selected hardware was changed slightly after soliciting information from the various manufacturers and determining which units were available for inspection. Twenty-nine out of 34 units were visually inspected, improving the depth of the suitability analysis.

2.4.1.1 Equipment-Type Comparison

Relative Procurement Cost. The quantity of each type of equipment estimated for the sortie experiment program and its retail cost was used to estimate the relative procurement cost for that item in the program. Equipment quantities identified in the functional requirements data bank were extended to the entire experiment program by applying the following ground rules and employing the payload flight schedule defined in the October 1973 mission model.

1. An equipment life of five years.
2. A single experiment module can be launched no more frequently than once every six months.
3. Hardware identified for one payload cannot be flown with a other payload unless it was identified as core equipment.
4. A single piece of equipment cannot be used in more than one discipline. A duplicate must be purchased if requirements are identical.

The five-year life is based on technological change, cyclic life considerations, and the duration of the sortie experiment program. Technology changes in any five-year period are expected to make replacement of experiment equipment with improved hardware desirable. For a few items that historically have experienced rates of slow technology evolution, a lifetime greater than five years was assumed. Greater scientific benefit will be achieved per flight with improvements in equipment capability justifying such changes. Since the equipment is not initially built for space, a maximum of 10 cycles (one launch every 6 months for 5 years) through the launch-space-reentry environment was judged to be a reasonable goal. A lower number of cycles would emphasize equipment of payloads which have repeat flights. Since equipment components of repeat flights are not totally defined, over-emphasis might lead to improper equipment selection. Finally, a few sortie payloads fly for more than 10 years. Therefore, this ground rule generally requires at lease one replacement during the sortie experiment program. Such a replacement schedule would be expected in most cases under any circumstances. The estimation of equipment quantities appears in Volume III, Section 5.0.

An average six-month turnaround duration for refurbishment of Spacelab experiment modules is based on in-house studies of experiment integration activities¹. This duration may be shorter for payloads which require little or no refurbishment, and longer for other payloads which have significant change-out or recalibration requirements (such as astronomy payloads). Six months is considered a fair average. A much shorter turnaround is envisioned for the Spacelab support module; however, it is assumed not to house experiment hardware.

¹NAS1-12933, Spacelab User Implementation Assessment Study

Equipment identified for a given experiment module was assumed to remain with that experiment module for many reasons: (1) refurbishment of experiment modules will probably occur at many locations--making the shuttling of equipment packages a complex, costly logistics problem; (2) the cost of additional checkout and calibration activities caused by removal and reinstallation of hardware will be greater than the purchase of an additional piece of hardware; and (3) the time spent for these additional integration activities would lengthen the turnaround time of experiment modules, possibly impacting the sortie experiment flight schedule. Similar logic applies to the use of the same hardware in different disciplines. Changeout complexities could be increased by minor modifications or accessories that may be particular to a discipline which are incompatible with another discipline.

Figure 2-6 shows an example of the determination of the relative procurement cost for a unit on board the Advanced Technology Laboratory (ATL). The October 1973 mission model shows that ATL flights occur throughout the 12 years of the mission model. A five-year equipment life study guideline requires that three sets of equipment be installed in the lab. Since a six-month turnaround time is required to refurbish each Spacelab experiment module, the ATL would require two experiment modules, doubling the equipment complement. A resultant total of six equipment sets are, therefore, needed in the ATL program.

OCTOBER 73 NASA MISSION MODEL

	CY	80	81	82	83	84	85	86	87	88	89	90	91	TOTAL
ADVANCED TECHNOLOGY LAB	2	4	4	4	4	4	4	4	4	4	4	4	4	46

5-YEAR LIFE				1					1				1	
6-MONTH TURNAROUND ADJUSTMENT				+1					+1				+1	
EQUIPMENT SETS							6							
OSCILLOSCOPE							x 1							
TOTAL EQUIPMENT QUANTITY							6							
RETAIL COST RANGE							x 1K TO 3K							
RELATIVE PROCUREMENT COST FOR ATL PROGRAM							6K TO 18K							

Figure 2-6. Estimation of Relative Equipment Procurement Costs

The cost range for oscilloscopes needed on the ATL is then determined by multiplying the equipment sets (6) by the number of units on the ATL payloads (1), and by the retail cost range of that equipment (\$1000 to \$3000). Similar calculations would be made in other disciplines also using oscilloscopes. The sum total of these amounts would define the relative programmatic procurement cost for that item.

These total costs should not be interpreted as actual programmatic costs. Procurement costs will probably be much higher because of modifications or, in



some cases, custom-building the hardware. However, this approach employs costs from a common reference, retail costs, and therefore the relation of the equipment in the ranking is valid. The ranking is only intended as a selection tool and not as an indicator of the absolute magnitude of Spacelab equipment costs.

Since the selection criterion uses relative cost as a discriminator, the ranking of hardware by cost to the experiment program would not change significantly if aerospace costs, rather than retail costs, were used. Retail costs are an indication of complexity, state of the art, and production rate. Aerospace costs generally are estimated relative to state of the art and complexity. Manufacturing lot sizes are always small. Savings through mass production drives the retail cost of some items down lower than they would be on a relative scale based on low-production aerospace development costs. However, not many scientific instruments are mass-produced. Those that are, are used in sufficient quantities in the sortie experiment program to raise their total procurement cost and--as a result--their relative ranking. Also, since a significant number of equipment items will be selected, location in the ranking relative to individual items is not important--only an item's location within a segment of the ranking.

Table 2-3 lists the common equipment types ranked according to relative procurement cost. Highest rankings went to the equipment types which had the highest average relative procurement cost.

Primary Operating Principle. The number of equipment types of interest was reduced from 90 to 25 by comparing operating principles of the equipment and eliminating those that were similar. Table 2-3 illustrates this process. Starting with the most expensive units, units are selected until lower cost units which have similar design characteristics are identified. These similar units are passed over in favor of other units that cost less, but are representative of different operating principles. (Cameras are exceptions to this logic and are discussed later.) For example, a TV monitor was not selected because it has a CRT tube and electronic components similar to the previously selected keyboard/display terminal or the oscilloscope. The holographic imager is a device using optics and a laser to generate three-dimensional images. It was rejected because a laser had previously been selected. Similar logic was used on other items.

Still, cine and TV cameras were excluded from selection because commercial units have already flown in space. Hasselblad and Maurer cameras were used on Apollo and Skylab. A vidicon developed by RCA for space operations was also used and is available. Selecting alternative sources to these units and performing a suitability analysis on them would not yield any new information. Also, cameras are small enough to store in a protective container when not in use, reducing the need for significant modifications. Although they were not analyzed, their successful use in space indicates that available cameras can cost-effectively be employed on the Spacelab.

The resultant breakdown of basic design approaches used for the selected types is shown on Figure 2-7. Total representation exceeds the number of equipment types (25) because of dual representation by some types such as an

Table 2-3. Equipment-Type Selection

Item	Average Relative Procurement		Rejection Rationale
	Cost (\$K)	Selection Action	
Tape Recorders	1560	S	Space-rated hardware available
Spectrometer/Spectrophotometer	840	S	
Lasers	720	S	
Computer	610	S	
TV Camera	570	R	
Receivers	530	S	
Blood Cell Counters	380	S	
Spectrum/Wave Analyzers	350	S	Space-rated hardware available
Frequency Generators			
Strip Chart Recorders/X-Y Plotters	290	S	
Cine Cameras	280	R	
Keyboard/Display Terminal	242	S	
Gas Chromatograph/Gas Analyzers	220	S	
Transmitter	156	S	Similar to gas analyzers
Metabolic Analyzer	122	R	
Microscope	114	S	
Computer Input/Output Devices (A/D converters, amplifier, etc.)	106	S	
Oscilloscope	102	S	
Timer	69	S	
Volt Ohm Meter	64	S	Similar to spectrum analyzer
Signal Generator	63	R	
Electrophoresis	54.4	S	
Symbol Generator	47	R	
Pyrometer	38	R	
Accelerometer	35	R	Similar to previously selected electronics
Scanner Programmer	35	R	
Laser Optics	32.8	R	
Furnace/Oven	30	S	
Radiation Sources	30	R	
Dewar	27.6	S	
Photo Processor	25	R	Similar to laser
Incubator	24	R	
TV Display	24	R	
Power Supply/Conditioner	22	S	
Refrigerator/Freezer	20	S	
Lyophilization Unit	19	R	
Acoustic Transducer	16	R	Similar to previously selected electronics
Optical Particle Detector	14.5	R	
Still Camera	14	R	
Bandpass Filter	13.7	R	
Centrifuge	13.5	S	
Blood Circulation Instrumentation	13.2	R	
Radiation Counter	13	S	Similar to previously selected electronics
Physiological Monitor	12.4	R	
Microfilm Viewer	12	R	
Power Meter	11	R	
Microtome	7.4	S	
Acoustic Generator	7.2	R	
pH Meter	5.9	S	Similar to power supply
Compressor/Pumps	5	R	

ORIGINAL PAGE IS
OF POOR QUALITY

electrophoresis device which is electro-chemical. The preponderance of electrical equipment is representative of the distribution of equipment types in the various sortie laboratories. Electronic instrumentation is prevalent in all the laboratories. Representation of the other technology areas is adequate to understand their associated design issues.




Figure 2-7. Equipment-Type Primary Operating Principle

Science and Application Discipline Representation. The selection criteria required selection of at least one unit from each discipline. The selection approach tends to emphasize those items which are needed in significant quantities and, as a result, the selected items are identified in more than one discipline. Figure 2-8 shows that each discipline is represented by at least five selected units reflecting this cross-discipline representation.




Figure 2-8. Discipline Representation

The quantities identified for each discipline reflect the general location of experiment equipment in that discipline. The experiment equipment of the Astronomy, High-Energy Astrophysics and Solar Physics disciplines are located on the pallet. The pressurized module acts only as a control center and as a result only contains core-type support hardware. Therefore, the representation of equipment from these disciplines is understandably low. Many experiments are performed in the Spacelab in the Space Processing, Space Technology and Life Sciences disciplines. In the case of Life Sciences, all are internal.

This discipline also has the greatest potential for using commercial equipment because its instrumentation requirements are frequently the same as those on the ground. Space Technology has high representation because its experiments are similar to those in Earth Observations, Life Sciences, and Communication/Navigation disciplines which are well represented.

2.4.1.2 Specific Package Comparison

Table 2-4 lists the selected equipment types and the specific available units that were analyzed for their suitability to space operations. Also shown is the representation of the selected units with respect to the selection criteria.

Design Characteristics Representation. Hardware was selected that provided broad representation of available hardware. Units were selected from firms that were leaders in their industry such as Beckman for biomedical instrumentation, Tektronix and Hewlett-Packard for electronic instrumentation, Clay-Adams for centrifuges, and American Optical for microscopes. Other units were selected because other sortie payload studies identified them as suitable for Spacelab operations, such as the Astro furnace. Still others were selected because they displayed design issues of interest. To achieve such representation, more than one available unit was selected for certain equipment types. For example, a NIM-type timer by Tennelec was selected along with the 19-inch rack-mounted Datatron timer. In another case, both commercial and mil-spec volt/ohm meters were selected. Additions such as these caused the final number of items for analysis to increase from 25 to 34.

Representation of different design specifications by the selected hardware is shown graphically in Figure 2-9. Available airborne equipment applicable to this study is limited to recorders, computers, signal conditioners, receivers and transmitters, leading to the low representation of this type of hardware. Three airborne units were selected. Selection of the Ampex AR700 occurred because it was not only an airborne unit, but it was the specific model that was modified by Martin Company for use on the Skylab. Selection of an airborne receiver and a transmitter was natural because of the availability of this type of hardware. The Mil-Spec units represent varying degrees of Mil-Spec design. The Fluke DVM and the Hewlett-Packard oscilloscope are more rugged than the Singer unit and the Fluke spectrum analyzer. The former units are designed for field use which accounts for the ruggedness of their designs. If the Army Field Medical Laboratory¹ effort ever moves into the hardware phase, additional Mil-Spec units in the biomedical field will be available to Spacelab.

Packaging approaches of the selected hardware include ATR, NIM, and conventional bench- and rack-mounted configurations. The distribution of the different configurations is shown in Figure 2-10. Many airborne units can be obtained as ATR modules. Of the three airborne units selected, the Collins receiver and the RHG transmitter were of this configuration. NIM configuration equipment appears limited to nuclear-related and computer-related hardware. In addition to the particle counter, the timer, and the analog-to-digital converter, a NIM bin power

¹SID 67-269, Phase II, Army Field Medical Laboratory Program, Phase II, Laboratory Feasibility Equipment Evaluation and Development Data Systems Study Report, Rockwell, 2/28/67



Table 2-4. Selected Equipment Representation Summary

EQUIPMENT-TYPE COMPARISON											SPECIFIC PACKAGE COMPARISON										
Equipment Type	Discipline Representation										Operating Principle Representation	Design Characteristics									
	Astronomy	High-Energy Physics	Solar Physics	Atmosphere and Space Physics	Earth Observations	Earth and Ocean Physics	Space Processing	Life Sciences	Space Technology	Communication/Navigation		Design Type	Package Config.	State of The Art	Commercial	Airborne	Nav-Spec	ATR	Special	Advanced	Current
Tape Recorder	X	X	X	X	X	X	X	X	X	X	X	Ampex	AR700								
Computer	X	X	X	X	X	X	X	X	X	X	X	Moneywell	5600								
Keyboard, Display Terminal	X	X	X	X	X	X	X	X	X	X	X	Digital Equipment	PDP-8a								
Blood Cell Counter												Research, Inc.	3311RN								
Timer	X	X	X	X	X	X	X	X	X	X	X	Coulter Electronics	En								
Spectrometer, Spectrophotometer												Datatron	3150								
Receiver-Transceiver												Tennelec	TC540								
Transmitter												Beckman	24								
Signal Conditioners (Amplifiers, A/D Converters, etc.)	X	X	X	X	X	X	X	X	X	X	X	Singer	NM37-57								
Optimum Wave Analyzers												Collins Radio	618M-2								
Frequency Generators												RHG	W1900								
Laser Assembly												Nell	126								
Gas Chromatograph, Gas Analyzer												Nuclear Data	ND100								
Microscope												Hewlett-Packard	1411 8552/8555 8445								
Strip Chart Recorder, XY Plotters												Singer	558 50-1L								
Volt/Ohm Meter												Fluke	645M								
Oscilloscope												Sylvania	948								
Electrophoresis Apparatus												Beckman	6700								
Furnace, Oven												American Optical	XH20TGQW								
De-war												Moneywell	185B								
Power Supply												Fluke	8700A								
Centrifuge												Fluke	B125A								
Refrigerator Freezer												Tektronix	4B5								
Particle Counter Assembly												Hewlett-Packard	AN/USM-281A								
Microtome												Beckman	R-100								
pH Meter												Astro	1000A								

ORIGINAL PAGE IS
OF POOR QUALITY

supply was selected so that the major elements of the NIM system could be better understood. CAMAC hardware was not selected, but its packaging approach is similar to the NIM approach. The remaining units were either bench- or rack-mounted. Most of the electronic packages can be obtained in a 19-inch EIA rack-compatible configuration. NIM hardware is also compatible with 19-inch racks. The ruggedized military units made for field use are exceptions to this generalization and are not available in a rack-mounted configuration.




Figure 2-9.
Design Specification
Representation




Figure 2-10.
Package Installation
Configuration Representation

The selected hardware displayed a good balance of current and advanced technology as shown in Figure 2-11. The term "advanced technology" requires interpretation since any available hardware could be classed as current by the mere fact that it is available. Components that are representative of the hardware expected to exist during the Shuttle era are advanced state-of-the-art units which utilized integrated circuits, MOS-FET devices, laser and fiberoptics for light transmission. Units selected for these advanced characteristics included the Research Incorporated display terminal which makes extensive use of integrated circuit chips, the Neff amplifier which employs MOS-FET devices, and the Honeywell strip chart recorder which replaces the typical pen/galvanometer marker with transmission of a CRT image by fiberoptics into photo-sensitive paper.




Figure 2-11. State-Of-The-Art Representation

Availability for Inspection. Final selection of model number also hinged upon the availability of hardware for inspection. Contact with suppliers was, in general, good. Many brought demonstrator units into Space Division and allowed the team to dismantle and photograph them. In other cases, the hardware was not available for inspection on the west coast. The model number was either changed or the suitability analysis had to rely only upon available specification data.

Hardware inspection was achieved on 29 of 34 items. Two activities contributed significantly to this high success ratio. First, the capital equipment tabulations of all Rockwell divisions in the Los Angeles area were consulted to identify selected hardware owned by this company which the team was assured of examining. Seven units were identified from this search and later examined. Second, a tour of the Beckman Instruments facilities yielded even better results. Nine more units were inspected at the Beckman plant. One unit in particular, the Coulter blood cell counter, could not have been inspected if Beckman had not had one in their laboratories. Contacts with Coulter Electronics had not been successful. When talking to the president of the company he recommended that we buy a unit if we wanted to inspect it. This type of response was not typical of the manufacturers and sales representatives contacted by the study team.

2.5 CONCLUSIONS FROM EQUIPMENT SELECTION ACTIVITY

The three major achievements of the equipment identification and selection activity were: (1) the generation of a functional requirements data bank documenting the equipment requirements for all internally located equipment identified in the sortie science and application disciplines; (2) development of a shopping list identifying available hardware that meets these functional requirements; and (3) selection of 34 items from the shopping list for suitability analysis.

During this effort, one conclusion was reached which is significant to the use of available hardware in space. Approximately 60 percent of the experiment equipment requirements of Spacelab could be satisfied by available hardware. Further study of available hardware catalogs might increase this percentage even higher. There is such a wide variety of equipment available that the requirements of the Spacelab hardware usually can be achieved by some supplier. This fact should encourage NASA and principal investigators to examine catalogs of available hardware to satisfy experiment requirements before initiating a program to custom build specialized hardware for Spacelab experiments.

**ORIGINAL PAGE IS
OF POOR QUALITY**

3.0 ANALYSIS OF SELECTED EQUIPMENT

3.1 INTRODUCTION

This section presents the suitability and cost analysis for the 34 selected equipment items identified in the previous section. Equipment characteristics were compared to Spacelab interface requirements to identify incompatibilities. Modifications were defined to correct these incompatibilities and the cost of these modifications were estimated. The costing of custom procurement of the same equipment is also described. The costs of the two alternatives were compared in order to determine if off-the-shelf experiment equipment is cost effective for Spacelab applications.

Figure 3-1 illustrates the analytical approach employed. The baseline criteria for the modification and costing are provided by a Spacelab/Experiment Equipment Interface Requirements (SEEIR) document specially generated for this study. This document reflects a policy of minimum NASA/Spacelab involvement with the experimenter prior to placement of the equipment on board Spacelab. The SEEIR document is included as Appendix D to this volume. Further discussion of rationale behind the document is presented in Section 3.2.2.




Figure 3-1. Logic Flow Diagram for the Analysis of Selected Equipment

The SEEIR criteria plus a set of standard design guidelines were established to a level of detail adequate to discern significant unsuitable design features and define acceptable solutions for the selected units. Suitability was determined from a combination of visual inspection of the equipment items and available data from catalogs, maintenance manuals, and responses by suppliers to questions. Basic suitability and specific unsuitabilities were uniformly determined by using the same engineering team for all items. Where logical, similar design modifications were defined for common unsuitable characteristics in order to arrive at consistent cost impacts. The suitability analysis and design guidelines are described in more detail in Section 3.2, and the modifications and results are covered in Section 3.3.

A "grass roots" cost estimating procedure was employed to establish modification costs. Experienced instrument manufacturing personnel from Beckman Instruments provided fabrication costs based upon current manufacturing practices. Rockwell and Beckman collaborated in developing standards for engineering and test costs. These costs, along with documentation and project management cost, defined the total modification cost for each item. The costing of custom-procured hardware employs traditional techniques of cost estimating. Cost estimating relationships (CER's), applicable to particular types of equipment, were utilized to estimate procurement costs. The cost analysis effort is described in detail in Section 3.4.

3.2 SUITABILITY ANALYSIS

3.2.1 Approach

The basic approach to evaluating the 34 selected items of equipment relied heavily upon visual inspection of the equipment items. Twenty-nine of the 34 items were visually inspected. All covers, subassemblies, etc., were removed to observe pertinent design features. Criteria checklists and reference to standard design guidelines that reflect acceptable solutions to various SEEIR criteria were used to evaluate each item for acceptance and shortcomings.

Table 3-1 summarizes the suitability and cost analysis activities for each equipment item. All items inspected, except the ND-100 multichannel analyzer, were photographed to show overall construction as well as to highlight features judged inadequate to meet SEEIR criteria. The photographs were referred to throughout the study for modification definition and costing. The RHG transmitter was not available for physical inspection, but the manufacturer provided internal pictures of it and similar units to illustrate its construction. Published data such as specifications, maintenance manuals and operations manuals were acquired to aid judgment on adequacy of design. Also, the counsel of other technical personnel engaged in testing or design on other hardware programs was sought to help verify judgmental standards. This was particularly true in the areas of vibration and material usage.

Information was collected on checklist work sheets similar to that shown by Figure 3-2. The acceptability, uncertain acceptability, or unacceptability of each equipment characteristic was checked on the work sheet and amplifying

ORIGINAL PAGE IS
OF POOR QUALITY

Table 3-1. Summary of Activities For Each Selected Item

Item		Pictures	Data	Inspection	Modification Definition	Modification Design	Fabrication Cost	Total Mod Cost	Custom Cost	Cost Comparison
AMPEX	AR700	X	X	X	X	X	X	X	X	X
HONEYWELL	5600	X	X	X	X	X	X	X	X	X
DIGITAL EQUIPMENT	PDP8	X	X	X	X	X	X	X	X	X
RESEARCH	3311RN	X	X	X	X	X	X	X	X	X
COULTER	Fn	X	X	X	X	O	Q	O	X	O
DATATRON	3150	X	X	X	X	X	X	X	X	X
TENNELEC	TC545	X	X	X	X	X	X	X	X	X
BECKMAN	24	X	X	X	X	X	X	X	X	X
SINGER	NM37/57	X	X	X	X	X	X	X	X	X
RHG	W1900	X	X	O	X	X	X	X	X	X
COLLINS	618M-2	X	X	X	X	X	X	X	X	X
NEFF	126	X	X	X	X	X	X	X	X	X
NUCLEAR DATA	100	O	X	X	X	X	X	X	X	X
HEWLETT-PACKARD	141T	X	X	X	X	X	X	X	X	X
SINGER	SSB 50-IL	(1)	X	O	X	X	X	X	X	X
FLUKE	645M	(1)	X	O	X	X	X	X	X	X
SYLVANIA	948	X	X	X	X	X	X	X	X	X
BECKMAN	6700	X	X	X	X	X	X	X	X	X
AMERICAN OPTICAL	H20	X	X	X	X	X	X	X	X	X
HONEYWELL	1858	X	X	X	X	X	X	X	X	X
FLUKE	8200	X	X	X	X	X	X	X	X	X
FLUKE	8125	(1)	X	O	X	X	X	X	X	X
TEKTRONIX	485	X	X	X	X	X	X	X	X	X
H/P	AN/USM	O	X	O	X	X	X	X	X	X
BECKMAN	R	X	X	X	X	X	X	X	X	X
ASTRO	1000	X	X	X	X	X	X	X	X	X
CRYO ASSOC.	LR90	X	X	X	O	O	O	O	X	O
SORENSEN	SLR40	X	X	X	X	X	X	X	X	X
POWER DESIGN	AEC320-3	X	X	X	X	X	X	X	X	X
CLAY ADAMS	CT1004	X	X	X	X	X	X	X	X	X
REVC0	ULT185	X	X	X	O	O	O	O	X	O
ORTEC	401	X	X	X	X	X	X	X	X	X
AMERICAN OPT	820	X	X	X	X	X	X	X	X	X
BECKMAN	Phasar I	X	X	X	X	X	X	X	X	X

(1) Pictures available in maintenance manuals or other data.
O No activity.

COMMERCIAL EQUIPMENT EVALUATION SUMMARY SHEET FOR INDIVIDUAL UTILIZATIONS IN SPACELAB						
Name: _____		Model: _____		Company: _____		
Feature	No Information	Source		Acceptable	Acceptable-- Verification Required	Unacceptable
		Literature	Inspection			
Comment						
I CONSTRUCTION A. Shatterables (with release restraints) B. 9-G mounting/integrity C. Protrusions and edges safety D. 19-inch rack mount capability E. Shock-vib-accel-acoustics resistance F. Depressur. hazard suppression G. EMI generation suppression H. EMI susceptibility protection II MATERIALS USAGE A. Flaking and peeling resistance B. Concentrations of flammable/ unidentified materials C. Non-prevalant commercial materials (or warnings on handling/usage of item) D. Resistance to combustion ignition III INPUT POWER COMPATIBILITY IV THERMAL COMPATIBILITY V DATA SYSTEM COMPATIBILITY VI ATMOS. CONTAM. RESISTANCE VII ZERO-G OPERABILITY A. Human factors B. Functional operation C. Loose parts/constraints VIII AUDIBLE NOISE						

Figure 3-2. Typical Evaluation Summary Sheet



remarks added as inspections and data evaluations were made. General modification requirement summary sheets were then prepared from the checklist data and from other data relating to peculiar requirements not covered by the general checklist. A sample modification requirements summary sheet is shown in Figure 3-3. In general, a conservative approach for costing purposes was taken by assuming the need for modifications where a reasonable doubt of suitability existed, as well as for design characteristics that were clearly unacceptable. Testing of these items would resolve the uncertainties. These modification summary sheets were used by design engineers to select standard modifications or to define any necessary unique designs as needed for realistic costing.

This same process was essentially repeated, using the same photographs and data accumulated in this task, to determine the delta modification requirements to implement the current version of NASA Equipment Specification EC006M00000A. The EC006M evaluation is described in Section 4.0.

MODIFICATION REQUIREMENTS SUMMARY

FOR

TENNELEC COUNTER (NIM ASSY) - TC545A Counter/Timer
(Excludes Nuclear Detector) - TC216 Linear Amp/Set
- TC909 Power Supply

- A. Provide 0.5-in. min. radii guard on NIM bin to prevent bodily contact with panel controls.*
- B. Vibration Protection (shockmount, worst axes random plus sine specs):
 - 1. Stiffen/support TC545 and TC216 circuit boards from excessive deflection (oil canning).
 - 2. Support/stiffen lead-mounted parts, such as ceramic capacitors, on circuit boards to prevent fatigue bending.
 - 3. Support end-mounted trimpots in TC545 and TC216.
 - 4. Support large lead-mounted capacitors in the TC909.
 - 5. Tie wires/bundles to prevent fatigue bending at connections (3 modules); see photos.
 - 6. Push-on transistor heat dissipators (2 in TC909) from loosening.
 - 7. Screws, knobs, thumbscrews from backing out (3 modules).
- C. Bake out for equivalent 150 hours operation to expel surface gases (3 modules).
- D. Provide flame retardant for unidentified organic electronic part materials on circuit board and circuit board itself (3 modules).**

*Assume that a 19-in. rack mount bin will be used in place of the portable unit inspected.

**Wiring is TFE per Tennelec rep. (Shuttle rated A for flammability and offgassing)

Figure 3-3. Example Copy Of A Typical SEEIR Modification Summary Sheet

3.2.2 Spacelab/Experiment Equipment Interface Requirements (SEEIR) Document

3.2.2.1 Document Use

In order to evaluate the selected equipment for suitability, it was necessary to have evaluation criteria available. Such criteria had to be at a level of detail to permit the analysis of specific designs of a wide range of functional equipment items. The Spacelab/Experiment Equipment Interface Requirements Specification was developed to fill this need. As the study progressed, the need for still more definition of acceptable design practice became necessary in order to assure costing consistency. The SEEIR design evaluation/acceptance interpretations were used to fill the need. Areas for further study have been identified in many sections. The results of these studies could be incorporated into the SEEIR at some future date.

The use of the current version of the SEEIR, for purposes other than this study, should be tempered by the fact that many of the cited values will probably change as a result of continuing Spacelab design definition. Also, a number of engineering judgments were needed to resolve source data that were incomplete or in conflict. More justification and analysis may be warranted to support these derived requirements.

3.2.2.2 SEEIR Document Basic Philosophy

The basic philosophy for the SEEIR document was to minimize the cost of experiment hardware by transferring the responsibility of successful equipment operation to the principal investigator (PI), recognizing however that some aspects of safety and experiment environment could not be readily compromised. The intent is to take advantage of the benign environment existent in the Spacelab and acknowledge exposure to typical, rather than worst case, ground environments in order to maximize the potential use of commercial equipment designs.

The PI would be responsible for assuring NASA that his equipment met safety requirements and his operational objectives. NASA would supervise nominal integrated tests and inspections, prior to taking the equipment on board Spacelab to verify that hazards or interferences to other Spacelab equipment do not exist. Savings could accrue to the PI and NASA by minimizing NASA involvement in such areas as program plans, design reviews, coordination/approvals, data generation and review and demonstration tests. Reliability, maintainability and other assurances would be consistent with the PI's success objective, the scheduled reuse of equipment and between-flight refurbishment capabilities.

The SEEIR is a first step toward a Spacelab experiment equipment specification. Its intended use was only as an aid to the suitability analysis performed in this study. If the SEEIR is used as a basis for the final NASA experiment equipment specification, further investigation is needed to define the mandatory assurance testing and inspections for various categories of equipment in order to achieve adequate safety confidence at low cost. Additional investigations in the following areas are also recommended: generation of simplified vibration requirements, development of a design evaluation "cook book" for experimenters/designers, and definition of precise criteria for acceptance of unidentified materials or conditional flammables/toxins.

ORIGINAL PAGE IS
OF POOR QUALITY

Data sources are listed in the Applicable Documents section (paragraph 2.0) of the SEEIR. The original SEEIR document was updated on the basis of the Spacelab Systems Requirements document of the ESRO RFP when it became available. A primary source of detailed environmental criteria was the ERNO Preliminary Spacelab Design Specifications (December 1973). The vibration specifications, in particular, are relatively severe for typical commercial equipment and should be reviewed and updated to the most recent design predictions in any subsequent uses of the SEEIR document. Shuttle program documents were used to define acceptable materials. A SEEIR update should include more direct reference to general application NASA documents. The role of source documents used for SEEIR is covered further in the following discussions of the detailed rationale behind some of the more critical criteria.

3.2.2.3 General SEEIR Contents and Organization

Figure 3-4 depicts the key contents and organization of the SEEIR. The specification mainly addresses the topics listed in the central triangular region of the diagram. These topics are covered in two sections. The first section defines the Spacelab capabilities and requirements, which can be categorized into three main groups: safety, environments and equipment support capability. The second section comprises the actual experiment equipment requirements. The latter can be categorized as those requirements that impact Spacelab safety and environmental preservation requirements, and those that are necessary only to contend with the Spacelab-provided environment and support.




Figure 3-4. Interface Specification Description
Key Contents and Organization

Equipment requirements that impact Spacelab requirements are considered mandatory. Other requirements are necessary for experiment success, but do not affect the integrity of the Spacelab.

The first section of SEEIR documents Spacelab characteristics. Experiment equipment requirements derived from these characteristics are placed in the latter part of the document. In many cases the equipment requirements, such as environments, are identical to Spacelab design capability, and hence, merely refer to applicable Spacelab requirements. In other cases, the equipment requirement may need more detail, or represent an allocation of the Spacelab total support capability.

3.2.2.4 Explanation of Specified Requirements

This section discusses the sources of data, problems, assumptions, rationale, use and selection of specific experiment equipment requirements that make up the SEEIR. Requirements can be conveniently categorized as safety provisions, equipment environments, and support capability. The discussion has been organized by key topics rather than by specification paragraphs. Pertinent subparagraphs under 3.1.1, Spacelab Design Characteristics and 3.2, Experimenter Equipment/Materials Requirements, of the SEEIR document are covered where significant issues may be involved.

3.2.2.4.1 Safety Provisions

Caution and Warning (C&W). The ESRO Spacelab RFP requires C&W for all potential hazards that may develop during the mission. The ERNO Spacelab studies and design data provide audio and visual display of caution and warning conditions. Independent emergency power for C&W is indicated. Each item will be evaluated for hazards and unique C&W needs. Functional failure of most experiment equipment items will not generally affect Spacelab operation. Hazards such as fire or leakage can be handled by "tack-on" sensors and integrated logic converters. In most cases sensors could be attached to the mounting rack or equipment exterior when prudent to augment overall Spacelab sensor systems. More definitive study is suggested to develop ground rules for providing necessary Spacelab C&W with a minimum of impact on experiment equipment and also to define those conditions where experiment equipment must provide specific design for C&W.

Flammability-Toxicity-Outgassing. The allowable materials listing from the Shuttle program was used to judge whether known materials were acceptable with respect to flammability, toxicity and outgassing for an air atmosphere. Mercury, cadmium and methyl chloroform were prohibited for any use based upon consultation with in-house materials experts.



All constituent materials could not be precisely identified for the equipment analyzed. However, it was concluded that the risk of dangerous quantities of toxics is minimal for off-the-shelf items in extensive use. Occupational Safety and Health Administration Act (OSHA) requirements, military specifications and accumulated experience make similar components acceptable for home, factory or laboratory. These same specifications reduce the risk to Spacelab inhabitants from materials in available equipment. However, since outgassing and flammability are rigid safety issues, all unidentifiable materials are assumed to be flammable. Threshold limits for chemical contaminants quoted in the specification were originally developed for the Space Station Phase B study and repeated in ERNO documentation. All wiring and PVC-appearing plastics are assumed to be PVC and flammable if not identified as a nonflammable such as TFE, phenolic epoxy, etc., which are approved by Shuttle. Circuit boards are coated with an impermeable film to prevent outgassing and fire. Typical electronic components are assumed acceptable without specific materials identification when in totally enclosed containers, preventing migration of air to support combustion, excepting certain parts/materials which are specifically prohibited such as mercury, cadmium, etc. The prohibited parts would be easily identified with the aid of cook-book rules to be provided (i.e., lamp bulbs and nixie tubes may contain mercury and must be replaced, etc.). Where deemed necessary, it is assumed that fire-retardant encapsulants could be applied. Cadmium-coated chassis are prohibited. It is assumed that unidentified materials could sputter flaming particles, therefore, screened or solid enclosures or fire-retardant encapsulation is required where ignition is possible. Combustible paths to points external to the case are to be eliminated.

One-hundred percent material identification is not required by the specification. Justification for not requiring 100-percent materials identification is strengthened where fused circuits and low-circuit voltages essentially eliminate short-circuit fire ignition. The precautions of conformal coating and 150-hour bakeout reduce hazards to an acceptable level. Metal enclosures, further enhanced by rack metal barriers, minimizing the potential for external ignition, provide added justification.

To reduce the hazard from trace gas contamination, it was decided, accordingly, to require a 150-hour equivalent operational bakeout to drive off surface gases and volatile heavy organics. The effectiveness of this process would be verified during Spacelab integration tests. Periodic gas samples from the bakeout chamber are recommended to assure that the volatiles driven off have dropped to an acceptable level after 150 hours. Longer durations may be necessary in a few cases. Elevated temperature and/or reduced atmospheric pressures could be utilized to expedite the process.

While the rationale for the SEEIR appears reasonable, a study effort for a more detailed analysis and documentation of rules and requirements is suggested. Cadmium-plated chassis (coated or uncoated) may, in fact, be acceptable--at least within enclosed/protected containers--and/or where credible ignition sources do not exist.

Electrical Grounds. Prevention of electrical shock from normally accessible surfaces, whether from static electricity or faulty circuits, is mandatory. Such control is consistent with 1970 OSHA standards. Compliance must also satisfy EMI control requirements for signal and case grounds.

Edges and Protrusions. The ERNO preliminary design data require 0.5-inch radii on edges and corners and no dangerous protrusions (presumably in travel lanes). This criterion for equivalent protective guards was adopted for SSEIR. Other source data suggested that a lesser radii may be tolerated. However, for this study the lesser radii would have little effect on cost estimates.

Mounting Integrity-Loose Items. Shuttle documentation indicates a maximum 9-g crash load on item mounting and integrity (i.e., no loose parts may become projectiles) in the +X direction. ERNO data also indicate up to 14 g rms (20 to 2000 Hz) random plus up to 3-g peak sinusoidal vibration for launch/reentry conditions for hard-mounted (to Spacelab) items. Up to 4.4 g rms (20 to 2000 Hz) random plus 3-g peak sinusoidal is stated for shock-mounted items. Since this was the only definitive data available it was accepted as baseline.

Loose items in Spacelab that might cause short circuits, injury, discomfort, inconvenience or damage were precluded. These include items that could come loose in a 9-g crash or "float" about in zero-g as a result of prior exposure to vibration. Chassis and part-mounting screws, thumbscrews, control panel knobs, door latches, and connector matings were reviewed to assure that they withstand both vibration and acceleration loading.

Touch-Temperature. The ESRO RFP and other data indicate that touch-surfaces should not exceed 115 F, which was accepted as baseline. Some MIL specifications, such as MIL-STD-1472A, allow 120 F (for handling) versus 140 F (inadvertent). The selected value was the most conservative and is consistent with other manned space programs.

3.2.2.4.2 Equipment Environments

Flaking-Peeling-Shattering. The ESRO RFP and ERNO documents site a Class 100,000 cleanliness environment for Spacelab on the ground, with 5-micron filtering subsequent to launch. Consistent with this requirement, internal and external surface coatings must be free of visual evidence of flaking or peeling. A more stringent requirement was unnecessary due to the relatively short mission duration in which degradation may take place.

A more critical requirement is to prevent shatterables (such as glass) from getting into the Spacelab atmosphere. This would include glass nixie and CRT tubes and fuses inside an item that could break for any reason and escape through holes in the case. An unresolved issue is the micron size and kinds of shattered particles and/or quantities that can be tolerated. Sealing or screening of case openings involves secondary impacts on thermal air-flow cooling. Within the range of uncertainty on allowable escape sizes, the modification costs would not vary significantly, precluding a need for detailed analysis of optimum mesh sizes.

EMI Generation. Extensive criteria in ERNO, ESRO RFP, MIL-I-6181D, MIL-STD-461 and -462, MSFC SL-E-001 and -002, and various spacecraft programs documents were reviewed. Satisfactory EMI control is historically difficult to achieve without fixes during verification tests, often because various shielding and circuit isolation rules are not followed.



Use of wave propagation theory resulted in the shielding attenuation characteristics shown in Figure 3-5 and Table 3-2. Attenuation as a function of shield material and thickness is given for E (electric), H (magnetic) and plane waves. Equipment items that are totally EMI enclosed (not necessarily physically totally enclosed) with 0.040-inch (0.1 cm) thickness ferromagnetic material attenuate EMI by about 200 dB or greater (factor of 10^{20}) for all plane waves and electric fields, and for magnetic fields over 10 kHz in frequency. Significant attenuation (factor of almost 10^6) occurs for H-fields down to 1 kHz. The 400 cps and 50 to 60 cps power frequencies are relatively less attenuated. Hence, high currents of these frequencies may need special or extra shielding such as multilayers including copper. Since magnetic fields attenuate rapidly as a function of distance, only circuits close to relatively large currents are usually affected (i.e., 12 inches/30 cm or less). Therefore, this intercoupling is more of a problem with external cabling, which is not considered directly in the item evaluations. However, the item output interface should allow good cabling practices. That is, interface connectors should facilitate routing of high-current, high-frequency signals or power supply wiring should be routed to special wire bundles and effectively shielded. Also separate signal returns, twisted and/or shielded, should be allowed by input/output circuits and connector interfaces when needed.


Finally, to prevent ground loops (that may induce input signal components caused by circulating currents due to differences in ground potentials) the dc signal returns and chassis/case grounds should be isolated so that they can be separately (tree-branch fashion) returned to single connecting points in Spacelab. Small filter capacitors connecting chassis to signal returns were deemed acceptable.

These criteria appear consistent with the specifications reviewed, but are more design definitive so as to provide a specific means with which to evaluate each item for EMI design compatibility. Again, the complex nature of EMI suppression warrants more detailed study in order to generate more complete cook-book design criteria, and exceptions, for actual use in selecting Spacelab equipment and needed modifications.

Considerations for item sensitivity to externally generated EMI is a non-mandatory Spacelab requirement that is treated in the next section. Also, note that internal-to-an-item EMI control is considered to be the experimenter's responsibility and therefore not addressed within SEEIR.

Audible Noise. The ESRO RFP referred to the NC50 standard noise criteria (NC) curves for general Spacelab environment. ERNO stated a 55-dB requirement, without reference to frequency spectrum, for Spacelab and also that 40 dB (no reference to distance or frequency was given) is the budget for each item in order to assure meeting the overall Spacelab requirement. The NC50 curve (MIL-STD-1472A, dated 15 May 1970) varies from about 73 dB at approximately 20 cps (center frequency of octave bands) to 47 dB at 4800 cps, with 55 dB at 300 cps. It was decided to retain the 40-dB item requirement for the 20 to 4800 cps audio range at two feet in front of the installed item.

ORIGINAL PAGE 11
OF POOR QUALITY



REFERENCE: ELECTRO TECHNOLOGY, NOVEMBER 1966, P. 85 - 94

Figure 3-5. Shield Effectiveness versus Frequency

**Table 3-2. Calculated Shield Effectiveness Of
0.040-Inch (0.1 cm) Sheet Iron (Copper)**

FREQUENCY	FIELD INTENSITY REDUCTION (DB)					
	H-FIELD		PLANE WAVE		E-FIELD	
60 Hz	12	(25)				
400 Hz	40	(30)	100	(150)	240	(300)
1 KHz	55	(34)	140	(150)	200	(240)
10 KHz	185	(45)	255	(160)	310	(200)
100 KHz	615	(90)	680	(200)	750	(240)
1 MHz	1,350	(150)	1,300	(220)	1,350	(250)
10 MHz	3,200	(470)	3,270	(480)	3,300	(550)
100 MHz	5,500	(1,250)	5,500	(1,250)	5,500	(1,300)
1 GHz	6,500	(4,000)	6,500	(4,000)	6,500	(4,000)
10 GHz	5,500	(12,000)	5,500	(12,000)	5,500	(12,000)

Some items have intermittent high-noise levels. A tape recorder rewinding a tape is a typical example. This type noise can be a beneficial warning signal and is generally of short duration. Therefore, to quiet such designs was considered unnecessary and probably not in the best interests of safe operation aboard the Spacelab. Hence, exception was made in the SEEIR for such short-term loud noises.

Again, further analysis--possibly, working models--is warranted to set more substantiated criteria. It was found that cooling fans in some equipment are relatively noisy. Cost of their modification was included in the study, but exceptions could possibly be taken in the final analysis.

Vibration-Shock-Acoustics. Vibration, shock and acoustics can also be considered a safety requirement to the extent that contribution to causing loose or shattered items in Spacelab could occur. However, this is only one possible cause of shattering. Of prime concern is the disturbance of electrical and mechanical adjustments, loosening of structure, and component failure due to fatigue or yield-limit stresses. Performance during the high vibration periods is assumed to be unnecessary since equipment will be "OFF" in general.

ERNO documentation states that shock will not be a driver. Since sources of shock are isolated through various attachments and structure, the energy will probably be converted to a vibration input at most item mounts, supporting this conclusion.



The ERNO acoustics-level information was included in SEEIR for completeness. The acoustics input from Shuttle is presumed to be transformed into the ERNO-specified acoustics and vibration as a result of structure and attachment characteristics. However, the acoustics impact on experiment equipment in general is considered insignificant compared to direct vibration inputs at mounts. Each item is, however, evaluated for any unusual susceptibility to acoustic inputs.

Of particular concern for item vibration acceptability are loosening of screws (chassis, parts, knobs, thumbscrews, etc.), circuit card deflection (oil-canning), parts and wiring fatigue breakage, and yield of solder joints or leads for "heavy" lead-mounted parts, secondary resonance drivers (transformers on "flat plates"), cathode ray tube (CRT) neck supports, cantilevered mounted parts, electronic and mechanical adjustments (decalibration), connector demating (harness, circuit cards, input/output) and other parts subject to resonance, yield or fatigue stress breakage. Each item was evaluated for such factors. However, there is great uncertainty as to the specific ability of an item to withstand vibration. The evaluation judgments made are felt to be representative for group costing purposes, but evaluation tests and/or refined evaluation rules are needed to better identify all problem modes on specific equipment.

Atmospherics. Specifications for relative humidity, particulate contamination, and the nominal equivalent air atmosphere were defined in the specification. These criteria were derived from the ESRO RFP and/or an interim NASA Sortie Lab reference document for experimenters. It was assumed that the experiment equipment would not be required to be functional after a Spacelab depressurization since depressurization is not a Spacelab design requirement and would undoubtedly result in an aborted mission in any event. The safety aspects of depressurization, however, must be considered regarding the release of fluids/gases or shatterables.

3.2.2.4.3 Item Performance

Mounting Locations. The ESRO RFP, ERNO and MBB designs indicate that 19-inch standard rack-mounted equipment items are desirable. Cooling, volumetric efficiency, integration, human factors, safety and other factors can presumably be more efficiently handled in this way. Racks can be made available for ceiling and under-floor as well as wall mounting. Therefore, the SEEIR provides for evaluating each item for suitability for rack-mounting. The requirement implies evaluating first for volumetric suitability. Where existing item design is for 19-inch racks (with or without catalog accessories), the capability for Spacelab installation (human factors, vibration and crash loads with respect to mass) is still necessary. NIM, CAMAC and other smaller size modules can fit into 19-inch rack-sized bins, crates or other substructures. Items not suitable for rack-mounting are assumed to be bench/shelf-mounted. The rack is assumed to have vibration/shock isolators with the resulting vibration characteristics cited in SEEIR. Benches were assumed hard-mounted so that any shock-mounting must be incorporated into the item or its mounting. In the absence of rack design criteria, it was assumed that cooling air would be made available in parallel to each item with needed ducts, deflectors or directors. Also, a rack-mount guide rail and support pin configuration was adopted for relatively heavy items.

EMI Susceptibility. It was necessary to convert EMI system requirements, such as a 6-dB margin between susceptibility and background EMI levels, or test levels/condition in MIL-STD-461 and -462, to design criteria that could be evaluated on each item. Internal generation-interference modes are assumed to be the responsibility of the experimenter. The requirement of the 0.040-inch (0.1 cm) metal case for generated EMI attenuation also provides incoming attenuation (shown by Figure 3-5 and Table 3-2), the same as for the EMI generation case discussed previously (3.2.2.5.2).

Isolation of input/output (I/O) circuit signal returns and common and power grounds avoids "reception" of common mode noise as well as the prevention of "transmission" discussed above. In addition, voltage levels and I/O circuit impedance requirements are specified in order to reduce the sensitivity to noise signals from high impedance coupled sources. Five to 12 volts output signal voltages were chosen since they are popular levels successfully used in industrial and avionics noise environments when proper design and installation rules are observed. The 10-K ohm minimum output impedance is easily obtained with transistorized circuits. A lower value could prove more desirable and would not appreciably affect the cost analysis for the study. The same voltages and impedance were assumed for full-scale analog signals. Theoretically, analog signals could run near zero volts although generally the signals of interest are made to run nominally at 2/3 to 3/4 full scale.

Input commands and other signals of an item must be compatible with the output voltage and impedance criteria of interface units. It was assumed that command logic would be in the 5- to 12-volt range with a 1-K ohm impedance. This would be a Spacelab design requirement.

The Spacelab remote acquisition units (RAU's) should provide for adjustable scaling within the noted voltage and impedance ranges, as a minimum, in order to avoid the general need to add circuitry to off-the-shelf items.

The voltage-impedance-connector requirements were not applied to front panel input/outputs intended for test probes or reconfigurable patch wiring to adjacent modules when evaluating items for suitability.

Input Power Form. The power forms established are consistent with the ESRO RFP. Regulation tolerances associated with each voltage level were selected from other documents and typical practice. It is assumed that the quantity of each Spacelab power type can be adjusted with parallel building-block converter units as long as the total available power, including conversion inefficiency, does not exceed the total allotted to experiment equipment.

A key assumption is that all ac power forms can be held to ± 1.0 cps in frequency. This appears reasonable, assuming electronic converters, and desirable since commercial equipment sometimes uses line frequency for clock/gate/trigger timing. Finally, the 50-Hz source was added since the ERNO RFP indicated the requirement. However, simplification in power converters appears possible by using only a 60-cycle source since 50-cycle equipment generally operates satisfactorily on 60 cycles. More effort is also justified to further refine Spacelab power supply requirements such as transient regulation, short-term frequency stability, etc.



The SEEIR document states that equipment items using 1.0 kw or more shall consider more efficient use of power. In retrospect, since this represents one fourth to one third of the available power, and since most experiments appear well below this need, 250 watts appear to be a better upper limit to require special efforts and integration coordination. Each mission will require a Spacelab management effort to coordinate total power usage timelines in order to assure adequacy of support without overloads.

More efficient means of power conditioning and distribution might be effectively implemented in the Spacelab; i.e., each rack could contain all power converters needed for that rack. The possibility of requiring the experiment packages to use +28 vdc only was considered. However, a trade study (see 4.3.4) showed such an approach to be undesirable. Most packages convert ac to many secondary voltages, few of which are near +28 vdc.

Data Handling. Data handling is assumed to be a Spacelab core function and, similar to power, not budgetable to individual items in a general SEEIR document. The capability stated in SEEIR generally reflects interpretations from ERNO documents as well as ESRO RFP requirements. This has little impact upon equipment items except to identify the need for interface capability. RAU capability to scale 5- to 12-volt logic or full-scale analog was assumed, consistent with the EMI requirements as discussed above. Where outputs to cabling do not meet requirements, buffer amplifiers were added.

Thermal. Values for ambient, wall, touch and ground storage temperatures were basically taken from ERNO documentation and/or the ESRO RFP. These temperatures, along with relative humidity, are key values in evaluating requirements. However, since they are relatively benign, small changes would have little impact on equipment modifications.

It was assumed in SEEIR that coldplates and rack air cooling would be provided by Spacelab. Although two units were water-cooled, the furnace and the laser, none of the units analyzed required coldplate cooling. Instead, the designs provided for internal water cooling. It is not clear at this time whether a liquid coolant interface to experiment equipment will be available in the Spacelab. Further studies of equipment planned for the Spacelab should be made to resolve this issue.

It was assumed that units with fans did not need additional cooling provisions unless obvious dead-air spaces containing significant power dissipation were present.

Zero-G Compatibility. The main areas of concern, as a consequence of placing 1-g designs into a zero-g environment, are functional capability and power dissipation. Human factors is a potential area for concern, however, it appears that usual operations on equipment such as adjustments, controls and switches can be operated effectively in space as long as shirtsleeve conditions and positioning restraints (foot/hand grasps, etc.) are adequate. Loose accessories, however, must be constrained from floating away in zero g.

The functional operation of each item requires careful scrutiny for subtle as well as obvious effects of gravity operations. A general test is to ascertain that operation is independent of orientation in a 1-g environment, which is not currently stated in SEEIR, but was used as a tool in evaluating each item.

The absence of natural convection in a zero-g environment requires evaluation of each unit for adequate thermal dissipation. For instance, the need for convective air currents in enclosed boxes is not obvious without access to design data. Therefore, it was assumed that moving air would be required over any power dissipation point unless it could be shown that direct conduction paths were adequate to dissipate the thermal energy to the external surface of the unit. External air circulation was assumed to be adequate in all cases.

3.2.3 Results of Analysis

A summary of results of the suitability analysis is presented in this section. The first portion includes illustrative examples of suitable and unsuitable conditions observed on the 34 items evaluated. A summary of the unsuitable conditions for all 34 items is presented in the last part of this section.

The type of modifications recommended for each item is presented in Section 3.3. Volume III provides the detailed evaluation of each item.

3.2.3.1 Examples of Suitable and Unsuitable Equipment Features

The following sequence of photographs were taken in conjunction with the suitability evaluation. Several dozen photos were taken of all but 6 of the 34 items. Those shown illustrate the range of suitability problems encountered.

Collins Transceiver. The Collins VHF transceiver, illustrated in Figure 3-6, is used in two-way communications and is potentially usable for the COM-NAV missions. It is illustrative of an ATR (Austin Trumbull Radio--commonly known as Air Transport [Association] Requirement) modular design with extensive, reliable commercial and private aircraft service. The design criteria are defined by ARINC (Air Radio, Inc.) for these units. It is an example of hardware that requires little if any design modifications for Spacelab experiment use. The close-to-the-board IC mountings and circuit board fastening indicate good resistance to vibration effects. Its totally enclosed case design eliminates any flammability problem since fire cannot propagate from this box, and the stagnant air inside the box would cause any flareup to be self-extinguished. The unit is qualified to operate up to 55,000 feet, which is an equivalent cabin pressure of approximately 1.5 psia. As a result, thermal dissipation by conduction to the outer case and mountings is designed into the unit. Even though the unit is enclosed, a 150-hour bakeout to drive off residual surface contaminants was recommended for costing purposes.




Figure 3-6. Collins Transceiver

ORIGINAL PAGE IS
OF POOR QUALITY



Views of the transceiver chassis, shown in Figure 3-7, illustrate mounting of heavier electronic components and internal EMI shielding. The heavy transformer is center-mounted rather than end-mounted. Electronic components are lead-mounted in such a way as to minimize the possibility of fatigue bending due to vibrations, and have relatively large lead cross-sections to support the component masses. This item could be a candidate for procurement later to confirm its apparent vibration integrity.

Hewlett-Packard Spectrum Analyzer. The H/P spectrum analysis equipment, pictured in Figures 3-8 through 3-12, illustrates typical high-frequency analyzer equipment studied. This system analyzes radio frequency inputs for various bandwidth energies and waveform characteristics, and is potentially usable in atmospheric physics, earth observations, life sciences, and space technology payloads.

The system illustrated consists of a mainframe display unit. Plug-in modules are RF and IF spectrum analyzer modules. These units are combined in a single package with built-in forced-air cooling. The automatic pre-selector on top of the mainframe unit is structurally independent with its own forced-air cooling fan.

Modifications recommended for this unit are not extensive. Control panel knobs were replaced because they are plastic. A different non-shatterable faceplate over the CRT display face is necessary to protect the astronaut from fragments resulting from implosion. The panel handles require modification or guards since the edge radii are less than 0.5 inch (1.25 cm). The panel knobs require shielding from bodily contact. The design guidelines discussed in Section 3.3.1 were developed to solve this typical guardrail requirement. The front panel signals going through the BNC-type coaxial cable input/output connectors are not required to meet the impedance and voltage level requirements of the EMI criteria of the SEEIR because these inputs cannot be modified without affecting performance.

Figure 3-9 displays CRT installation typical of electronic instrumentation. Modifications identified include closing all openings in the CRT shield with suitable mesh screening in order to trap glass particles should CRT breakage occur. The filter on the fan inlet would probably be adequate to prevent flaming particles from passing through. The outlets are in the form of open holes approximately 3/16-inch (0.5 cm) in size, which are typically assumed to need screens to contain shatterables and/or potential flaming particles. High heat dissipation power transistors are mounted close to the fan outlet, eliminating any cooling shortcoming in zero g. The large electrolytic capacitor bank needs support due to cantilevered mounting.

The module shown in Figure 3-10 illustrates ruggedness and EMI shielding typical of many well-designed commercial high-frequency instruments. Well-anchored cabling indicates good vibration integrity although the wiring bundles must be replaced because the wires are insulated with PVC. Polyvinyl chloride is not an acceptable material from flammability and outgassing standpoints. The large number of electrical calibration adjustments need staking for protection against vibration.




Figure 3-7. Collins Transceiver - Chassis View




Figure 3-8. Hewlett-Packard Spectrum Analyzer

ORIGINAL PAGE IS
OF POOR QUALITY




Figure 3-9. Interior View of H/P Spectrum Analyzer Display Section




Figure 3-10. H/P Spectrum Analyzer - RF Section

ORIGINAL PAGE IS
OF POOR QUALITY



Figure 3-11, the photograph of the IF unit, illustrates typical circuit board construction. The boards with relatively heavy parts in particular may need added support/dampening to prevent "buzz" and "oil-canning." Conformal coating to better anchor the relatively high-mounted, thin-leaded ceramic capacitors and reduce fire and trace contaminant hazards is recommended.

Figure 3-12, showing the interior of the preselector unit, indicates relative ruggedness. Modifications include support for the tubular waveguides and the large end-mounted component on the circuit board. The large transformer on the large flat plate bottom mounting surface is questionable. However, its location near the sides of the unit was assumed sufficiently stiff to preclude modification.

NIM Bin. Figure 3-13 illustrates typical NIM modules and their 19-inch rack-compatible bin which contains these modules. The modules shown include signal conditioners set up as a single-channel wave analyzer for nuclear instrumentation inputs.

NIM components have potential use on all Spacelab missions. They offer the advantage of interchangeable modules for easy maintenance, relatively rugged design for vibration integrity (although circuit boards do require additional stiffening), and venting provisions compatible with forced-air cooling.

A major disadvantage of this type of design is the poor volumetric and weight efficiency, which is a penalty for the convenience of a modular design approach that allows a wide range of custom functions to be mechanized.

A typical NIM power supply attached to the rear of the NIM bin is shown in Figure 3-14. The resultant cantilevered weight will require either special end support for the bin, or moving the power supply to an independent rack mount location. Depending on the voltage levels of available power in the Spacelab, the power supply could be eliminated. NIM modules are compatible with 12-volt and 24-volt dc power. Also, off-the-shelf power supplies of various types that mount in the bin are available.

Tennelec NIM. Figure 3-15 shows typical NIM design practices. TFE wiring is often used. Circuit boards tend to be rather large and may need stiffeners to prevent excessive deflections (oil-canning). Fiberglass circuit boards are also often used for excellent fire resistance. Conformal coating is recommended to support lead-mounted ceramic capacitors.




Figure 3-11. Interior View of H/P Spectrum Analyzer IF Section

ORIGINAL PAGE IS
OF POOR QUALITY




Figure 3-12. Interior View of H/P Spectrum Analyzer Automatic Preselector




Figure 3-13. NIM Bin With Modules - Frontal View

ORIGINAL PAGE IS
OF POOR QUALITY

3-27

SD 74-SA-0047-2




Figure 3-14. NIM Bin - Rear View Showing Power Supply

ORIGINAL PAGE IS
OF POOR QUALITY





Figure 3-15. Tennelec NIMS - EMI Covers Removed

ORIGINAL PAGE IS
OF POOR QUALITY

Beckman Gas Chromatograph. The gas chromatograph, shown in Figure 3-16, identifies gaseous constituents for petroleum industry process control. It is capable of the same scientific accuracies as laboratory instruments, and can be used on the space technology missions as well as other missions.

Of particular interest is the explosion-proof construction of this unit. Integral electronics are contained in the three heavy (white) cans, which are mounted to thick steel plates. TFE-insulated wiring is used for low flammability. The gas routing switch and plumbing, heating unit and the sensor (right side) are in a ruggedly constructed oven.

The insert picture shows a unique built-in troubleshooting panel located in the process controller electronics. This controller is separately located from the sensor unit and can be rack-mounted in the Spacelab.

The analyzer plumbing needs support and the door (not shown) needs positive latching for vibration and 9-g loads. Some exposed asbestos fiber insulation needs replacement or coating. The unit currently consumes 1500 watts for the test chamber heater. This heater can be replaced by a component dissipating 250 watts for Spacelab application because the volume of air currently flowing through the oven can be significantly reduced.

Beckman pH Meter. The Beckman pH meter, shown in Figure 3-17, measures the hydrogen ion activity in liquids and has potential application to the space processing and life science payloads. The unit is basically a high-gain amplifier that measures impedance, with a liquid crystal display unit.

The following changes are proposed. The plastic case, typical of some commercial units, is not an acceptable material. The case should be replaced with metal. No cooling provisions other than natural convection are provided by this design. However, the total power requirement for the unit is only 3 watts. Since there are no single high-dissipation components and the total dissipation is low, modifications for thermal control were considered unnecessary. The glass liquid crystal display unit requires a lexan plastic shield to protect the astronauts from glass particles in case of breakage. A rack-mounted installation is recommended. The pH sensors (not shown) currently consist of open-probed containers and must be redesigned for zero-g operation.

Sylvania CO₂ Laser Interior. The laser unit, shown in Figure 3-18 with the cover removed, consists of glass tubes up to about one yard long which contain CO₂ with integral water-cooling jackets and adjustable optics to obtain lasing action and to direct the light energy out. The unit has potential application to the space processing missions where a concentrated energy source is needed for heating material samples and other experiments needing a monochromatic light source.

The unit will require special mounting with integral shock mount provisions. Alignment in orbit may be required after exposure to boost vibration. Of interest is a completely silicone-encapsulated high-voltage power module, which is acceptable from a materials usage viewpoint. The use of external water cooling requires adequate plumbing support, particular at laser connections. Also, the availability of a liquid-coolant interface with the Spacelab thermal control subsystem must be verified.





Figure 3-16. Beckman Gas Chromatograph



Figure 3-17. Beckman pH Meter




Figure 3-18. Sylvania CO2 Laser Interior

ORIGINAL PAGE IS
OF POOR QUALITY



Research, Inc. Display Terminal. The display unit, shown in Figure 3-19 without enclosure cover and plastic CRT shield, provides entry and display of data associated with computer storage and has potential use on all Spacelab missions. Although an advanced functional design, it is judged to represent an extreme of the modification spectrum for the selected electronic equipment. In this view, it can be seen that the large vertical circuit boards and the long keyboard circuit board require support to survive the vibration environment. Socket-mounted IC's and some other components may also require cementing. Ribbon cabling must be replaced because of flammable insulation. Wire runs need tie-downs. The large transformer is end-mounted on the chassis and likely to generate secondary vibration inputs. The CRT plastic face shield (not shown) needs to be better mounted and secured.




Figure 3-19. Research, Inc. Display Terminal

Additional nonsuitable design features are shown in Figure 3-20, and require the following additional modifications:

- . Support large electrolytic capacitors on right circuit board. One capacitor is attached to two surfaces not rigidly attached to each other.
- . Support bottom sweep generator circuit board and its components, particularly the horizontal oscillator coil.
- . Isolate the powdered iron core flyback transformer from vibration to prevent fracturing.
- . Provide EMI shielding for sweep and high-voltage generator circuits since these type circuits are especially good EMI generators.
- . Provide support for CRT neck deflection coils and focus components; cement adjusted positions.
- . Provide CRT enclosure that would contain broken glass.
- . Provide porting for forced-air cooling to heat-sinked transistors, transformers, and circuit boards.
- . Provide connectors compatible with EMI design control practices.
- . Remove flaking high-voltage corona spray material at the high-voltage connection to tube and replace.

As with most other equipment items, all screws and adjustments not having a vibration-proof feature should be cemented against vibration.

Coulter Blood Cell Counter. The unit pictured in Figure 3-21 counts small particles, typically blood cell or other biological components, in a liquid suspension solution. The function of this unit is applicable to the life science and space technology payloads. A precise volume of liquid which is sucked through a micron-sized orifice is measured. A falling mercury manometer column, prepositioned with an integral vacuum pump, provides the calibrated volume and secondary vacuum to intake the liquid sample. Electronic impedance-type sensors and associated circuits count and display the count of particles passing through the orifice for the measured volume. Optics and CRT displays allow direct visual monitoring of particle size and orifice clogging.

This unit requires a complete redesign. Its function is gravity-dependent. It contains mercury, has little structural integrity, and is full of glass. It is the least suitable piece of equipment for spaceflight that was examined.




Figure 3-20. Research Inc. Display Terminal - Interior View





Figure 3-21. Coulter Blood Cell Counter





3.2.3.2 Modification Action Summary

The summarized results of the evaluation of 34 units are presented in Table 3-3. The table shows the number of items found to have unsuitable characteristics requiring modification. It should be emphasized that a conservative approach was taken in defining modifications in order to not undervalue modification costs nor over-simplify the degree of modification. Most modifications appear straightforward and have minor impact on a given package. Section 3.3 provides more detail on specific modifications for each item.

Table 3-3. Modification Actions Summary

Action	Number Required (on 34 Units)
Vibration protection	30
Improved cooling	9
Substitute materials	27
Toxics bakeout	30
9-g integrity measures	10
EMI attenuation	10
Shatterables containment	18
Edges-protrusions guards/rounding	32
Functional redesign (zero-g)	9*
Audible noise reduction	4
*3 units require complete new design	

Some vibration protection was needed for almost all units. Testing could show a greater or lesser need for redesigns. Final Spacelab definition and shockmount studies may also change the severity of the problem. Without providing much greater depth of analysis on each unit, the modifications defined can only be categorized as "representative" for cost study purposes. Resonance and cross-axis effects, etc., could result in some rather sturdy-appearing mounts to be insufficient.

Surprisingly few cooling modifications were identified. Many units had self-contained forced-air fans and good air flow paths. These same cooling fans generated noise in excess of the allowable 40 dB noise criteria in four cases. Airborne units are designed for conductance of heat to the case so that they can operate at altitudes up to 55,000 feet without internal air circulation. Most rack-mounted equipment items are vented to allow forced-air cooling in the rack installation. It was assumed that rack air could be supplied in parallel to each unit where needed.

Material substitutions were high due to conservative replacement of wiring, plastic panel knobs and fixtures if not positively identified as acceptable. By excluding knob replacement, only 20 items need material replacement. Only 13 units need other than knob/wiring-related replacements. These include 5 nixie tubes which contain trace amounts of mercury. Likewise, toxic bakeout

was assumed to be necessary for all units containing organics in order to preclude heavy initial outgassing from components, solvent cleaners, process catalyzers, etc.

EMI provisions were generally satisfactory based on a ferromagnetic case and input/output circuit compatibility with assumed Spacelab EMI grounding and cabling practices. Modifications were required for 10 units.

Slightly more than one-half the units contained glass shatterables requiring added containment or replacement. Most units had exposed sharp edges and protrusions. Control panel knobs were considered to be protrusions. Handles and guards were generally judged to be inadequate. Redesigned guard rails were assigned to the units, but could be part of the rack design.

Finally, several biomedical systems use open vials in sample handling (pH meter, electrophoresis system, blood cell counter, and spectrophotometer) or mechanical gravity seating (microscope) and require functional redesign. The freezer and dewar depend upon gravity for vapor phase separation processes. These units and the Coulter counter were not amenable to modification; i.e., a different operating principle must be used, requiring total new developments.

3.2.4 Suitability Analysis Conclusions

As might be expected, the most significant modifications appear to be for vibration integrity and replacement of unsuitable materials. However, both areas require test verification and further in-depth study because:

1. The final criteria to be met have not been established. Trade studies discussed in Section 4.3 indicate a benefit for Spacelab to provide a better environment than given by SEEIR.
2. Both characteristics are difficult to evaluate by inspection. Vibration effects are complex; resonances could exist at locations which are difficult to predict. Materials could not frequently be identified because of inadequate supplier information, which will be the case for the actual selections of equipment items. A better definition of allowable risks due to Spacelab safety features should be pursued.

Rockwell experience indicates that ATR, MIL or NIM-type equipment items have inherent design characteristics needing few modifications. Chemical and biomedical equipments tend to be less suitable because most have been designed for sensitive laboratory work in a benign environment. Also, test samples are commonly liquid, utilizing gravity during processing.

3.3 EQUIPMENT DESIGN MODIFICATIONS

This section summarizes the efforts and results related to specific modification designs and processes required to make suitable those equipment items found to be unsuitable. As the study progressed, certain unsuitabilities and redesign situations became repetitious. A set of standard design guidelines evolved and were used to help assure consistent design changes which, in turn,



led to a more consistent cost analysis. The design guidelines are presented in Section 3.3.1. A summary of recommended processes is presented in Section 3.3.2.

3.3.1 Standard Design Guidelines

The following guidelines were established to enable consistent costing for similar equipment unsuitabilities. The application of specific guidelines may not necessarily be the best redesign approach in some cases. However, they are adequate to assure proper cost estimates for modified hardware. Equipment also required non-typical modifications not covered by the guidelines. Such unique designs are found in Volume III under the modification description section for each item.

3.3.1.1 ARINC-To-NIM Conversion

Since all ATR (ARINC) cases have a standard 7.625 in. (19.4 cm) maximum height, they will fit inside a NIM (nuclear instrumentation module) bin of the 8-3/4 in. (22.2 cm) type whose front opening height is 7.785 in. (20 cm). The ATR case is modified by the addition of a front panel and side panels which provide top and bottom runners which mate with the guides in the NIM bin as shown in Figure 3-22. The connector mounting plate in the rear of the NIM bin is removed in the area where the over-long ATR case protrudes.




Figure 3-22. ARINC-To-NIM Conversion

3.3.1.2 Cable Clamps and Tie-Wraps

Cable clamps and tie-wraps must be nylon. Unsupported lengths of cable bundles and single wires must be secured. Wires going through holes in metal bulkheads must be protected by nylon grommets.

3.3.1.3 Card Guides - Printed Circuit Boards

Card guides which do not provide vibration protection must be replaced with spring clip guides which provide support as well as high-frequency isolation and good thermal conductivity. See Figure 3-23.




Figure 3-23. Card Guide, Printed Circuit - Positive Securing

3.3.1.4 Component Clips

Lead-mounted components over 4.0 grams shall be supported by rivet-mounted metallic clips as shown in Figure 3-24. Where clips may short out conductor lines of the printed circuit, a thin glass epoxy insulator must be placed between the clip and the board surface.




Figure 3-24. Component Clips

3.3.1.5 Fasteners (Vibration-Proof)

All structural fasteners must be those with thread-locking features. An epoxy formulation may be used, such as "Glyptol," in nonstructural situations such as potentiometer locking, etc. Thumbscrews require replacement or safety wiring. Cantilevered items needing support will have brackets to available solid structure to which component is mounted.

3.3.1.6 Glass Viewing Covers

Meter faces, CRT shields, etc., must be replaced with transparent Lexan when not identified as acceptable.

3.3.1.7 Material/Hardware Replacement

Wire harness and other PVC parts will be replaced with Teflon. Cadmium-plated chassis must be stripped and replated or reformed of new non-cadmium material. All panel knobs are to be of custom-fabricated Vespel with steel shaft inserts and double shaft set screw fasteners, staked to prevent loosening in vibration. Metal knobs are also acceptable. It will be assumed that (except in unusual circumstances) no transformers or other electrical components will be replaced for chemical composition reasons. Where required, and dependent upon individual circumstances, components will be coated with a suitable material (e.g., epoxy).

3.3.1.8 Non-Rack-Mounted Instruments


Portable instruments, such as bench-top microscopes, must be secured in assigned stowage spaces when not in use. When used, provision must be made to (1) secure the instrument to a bench, and (2) ensure against injury to personnel and damage to the instrument. Securing of non-portable, non-rack-mounted instruments such as furnaces will be treated on an individual basis. Where necessary, the mount must provide for shock/vibration isolation and alignment caging. Size the structure per SEEIR safety factors to eliminate proof-testing.

3.3.1.9 Mechanical Securing of Rack-Mounted Instruments

Heavier instruments such as the 20-inch (50.8 cm) deep Sorenson power supply, with heavy transformers in the rear of the unit, and units weighing over 30 pounds (13.6 kg) will require additional mechanical support. This will be accomplished by addition of tapered-pin assemblies attached to the rear of the instrument, and modification of the rack assembly to accept the tapered pins, plus appropriate guide rails, as illustrated in Figure 3-25. Catalog off-the-shelf mounts may be used when satisfactory for the 9-g crash load and worst case vibration. Structure sizing is per SEEIR safety factors to eliminate proof-testing.

3.3.1.10 Printed Circuit Board Securing

Printed circuit assemblies and terminal boards requiring additional support to prevent oil-canning or vibration must use drilled and tapped



**Figure 3-25. Mechanical Securing
of Heavy Rack-Mounted Instruments**

standoffs which can be attached to a convenient bulkhead. In locations where a metallic standoff would short out printed circuit conductors, rectangular blocks of glass epoxy may be used as shown in Figure 3-26. The blocks are to be attached with low curing-temperature epoxy adhesive.




Figure 3-26. Mechanical Securing - Printed Circuit Boards

3.3.1.11 Protrusion Protection

The prevention of injury to personnel due to zero-g collisions with knobs, switches, etc., on the front panels of rack-mounted instruments will be accomplished by a modification of the rack cabinet consisting of rails or handles mounted horizontally across the front of the cabinets. The rails, 1.0 in. (2.5 cm) in diameter, will be secured to the exposed vertical section of the cabinets; see Figure 3-27. They will extend far enough out from the front panels of the instrument, e.g. 2-1/2 in. (6.3 cm), that in any accidental encounter by personnel the rails would be contacted before any knobs, switches, etc. The rails will be so located that they will not obscure the view or interfere with any function of the instruments. Edge protection on non-rack-mounted items will be by fastening a rolled tubing or extruded plastic along the edge by an attachment flange which is drilled, screwed, and staked for vibration.




Figure 3-27. Protrusion Protection

3.3.1.12 Venting of Chassis, Covers, or Bulkheads

Forced convection cooling and elimination of pressurization hazards dictate venting where required. Where glass components such as CRT's are involved, cover inside of vents with metallic fine-mesh screen at the most convenient enclosure that would contain any fragments.

3.3.1.13 Modification Certification

All modified equipment items require acceptance testing. Additional criteria to verify safety and performance requirements compatible with Spacelab are applied.

3.3.1.14 Chassis/Signal Return Grounds

Isolate all circuit boards and other components grounded to chassis through mounting screws by using insulating washer kits. Remove all wired grounds to the chassis. Interconnect all isolated returns to a bus wire, going through connector pins as necessary.

3.3.2 Summary of Recommended Modifications

Table 3-4 summarizes the specific modifications accomplished for the 34 items analyzed. Included in the summary are the dewar, blood cell counter, and freezer which were not priced due to functional incompatibilities with zero g. Modifications to these items would be meaningless since they would necessitate custom-design efforts.

3.3.2.1 Staking

Fasteners on most items required positive retention to prevent loosening during vibration. Self-tapping or threaded-in screws and oversized holes/slots for holding structural members can quickly loosen and come apart or alter dynamic responses. Such screws were replaced by bolts or Locktite-type fasteners. Also, control panel knobs on some equipment items were push-on type or had a single screw for holding to control shafts. Knob set screws were staked with epoxy.

Some items have up to 80 potentiometer adjustments. If one or more adjustment shifts slightly during launch, the items' usefulness would likely be severely compromised. In-flight realignment is not recommended due to relative inaccessibility and lack of support equipment. Potentiometer adjustments were epoxied into place after calibration. Mechanical adjustments were few and appear to be relatively easy to readjust for the items inspected, although exceptions may be possible.

Wiring tie-downs were recommended to keep "flapping" wires from breaking at terminations.

The use of push-on thermal dissipation fins on transistors and IC's is relatively common. Some were tight-fitting, while others were found relatively insecure and could easily work loose causing a short or a component to over-heat. Epoxy adhesive was placed on these fins to prevent them from loosening.

Table 3-4. Detailed Summary of Item Modifications

ITEMS	VIBRATION - SHOCK				EMI	MATERIALS		ZERO-G CONTAMINANTS		CRASH INTEGRITY				ALUMINUM NOSE		EXHAUST FAN	4000 COOLING	SHATTER-RESISTANT PROTECTION	PER-FLUOROPOLYMER	NOTES
	STRENGTH	CONSTRUCTION	MOUNTING	LATCHING		OTHER	FUNCTIONAL	OTHER	ITEM	ITEM MOUNTING	SOURCE TYPE	CONTINUOUS	OTHER							
														VIEW GLASS	OTHER					
1. TAPE RECORDER-ANALYZER 4070																				
2. TAPE RECORDER-HONEYWELL 1400																				
3. COMPUTER-DIGITAL EQUIPMENT																				
4. KEYBOARD SERIAL TERMINAL-1100 21 NW																				
5. BLOCKED COUNTER-COUNTDOWN 1100																				
6. THERMISTOR-3150																				
7. EVENT COUNTER-THERMISTOR																				
8. SPECTROMETER-RESEARCHER 24																				
9. SPECTROMETER-RESEARCHER 24																				
10. SPECTROMETER-RESEARCHER 24																				
11. SPECTROMETER-RESEARCHER 24																				
12. SPECTROMETER-RESEARCHER 24																				
13. SPECTROMETER-RESEARCHER 24																				
14. SPECTROMETER-RESEARCHER 24																				
15. SPECTROMETER-RESEARCHER 24																				
16. SPECTROMETER-RESEARCHER 24																				
17. SPECTROMETER-RESEARCHER 24																				
18. SPECTROMETER-RESEARCHER 24																				
19. SPECTROMETER-RESEARCHER 24																				
20. SPECTROMETER-RESEARCHER 24																				
21. SPECTROMETER-RESEARCHER 24																				
22. SPECTROMETER-RESEARCHER 24																				
23. SPECTROMETER-RESEARCHER 24																				
24. SPECTROMETER-RESEARCHER 24																				
25. SPECTROMETER-RESEARCHER 24																				
26. SPECTROMETER-RESEARCHER 24																				
27. SPECTROMETER-RESEARCHER 24																				
28. SPECTROMETER-RESEARCHER 24																				
29. SPECTROMETER-RESEARCHER 24																				
30. SPECTROMETER-RESEARCHER 24																				
31. SPECTROMETER-RESEARCHER 24																				
32. SPECTROMETER-RESEARCHER 24																				
33. SPECTROMETER-RESEARCHER 24																				
34. SPECTROMETER-RESEARCHER 24																				

ORIGINAL PAGE IS
OF POOR QUALITY

3.3.2.2 Circuit Board and Parts Stiffening

Most circuit boards were judged to be in need of added support to reduce deflections. The size of the board, mounting, mass, and mass distribution largely determines the dynamic response. Solutions included firm but resilient guide/holding rails, braces across the board, and added center post mounting standoffs for bolted-down boards. Parts mounted on the boards were also often judged in need of support. Typical military-style clips, ties and/or conformal coating provided such support. Board encapsulation is discussed further in Section 3.3.3.

3.3.2.3 Other Supports

Large parts, such as transformers, mounted on flat-plate chassis, and end-mounted parts creating cantilevered masses, required custom braces to reduce secondary resonances and damaging torques to base mounts. In a few cases, doors or modules were judged not reliably latched to assure retention during launch and/or reentry. Each design modification was customized for the situation. Additionally, some equipment items utilize connectors that do not have positive engagement. Suitable connectors and/or special brackets were installed in these cases.

3.3.2.4 Electromagnetic Interference (EMI)

Chassis and signal ground isolation and input/output circuit voltage and impedance characteristics are the main areas of concern. Several units had single-point connections of signal/power ground to chassis (usually at the power supply) which can be easily corrected. Others (particularly spectrum analyzers) have multiple point grounding at circuit board mountings, connectors, etc. These units require separate ground buses throughout the unit. This could adversely affect performance of some units and either the requirement may need waiving, with suitable noise filters added, or a different design selected.

A few items required one or more buffer amplifiers added to the output to achieve the required 5-volt signal levels. These changes could be waived with good external cabling practices, assuming that the Spacelab remote acquisition units (RAU's) can accommodate wide-range interfaces. Although IC amplifiers are envisioned in most cases, power supply voltages, space and other factors could make this type modification relatively difficult.

3.3.2.5 Material Replacement

Wire harnesses required the most frequent replacement because PVC insulated wire is commonly used on commercial units. Modification difficulty varies with the unit as some harnesses are accessible by removing all plug-in assemblies, while others are hardwired and/or require nut and bolt disassembly for changeout. Some harnesses are suitable for prefabrication while others require point-to-point rewiring.

All plastic knobs required replacement. Metal or Vespel knobs were identified as replacements.

Viewing windows, such as meter faces, made of glass were modified. Glass faces were replaced or covered with a Lexan-type plastic to contain shatterables.

Foam materials used for vibration isolation, anti-rattle packing between modules, and hold-down bars and around case and door lips were replaced with non-flammable foam. Significant foam insulation replacement on the freezer chest would also be required, however, its modification became impractical for functional reasons.

It was determined that nixie tubes generally contain a trace of mercury to prolong cathode life. Four items examined (other than the blood cell counter) contain nixies. Life of tubes without mercury was reported to be only about 500 hours, or more than an order of magnitude less than those with mercury traces. Alternatives are to decide that the mercury traces are insignificant risk, use the shorter-life tubes, or develop a direct replacement using light-emitting diodes (LED's) with integrated circuit (IC) logic drivers on a nixie-sized plug-in circuit. The Navy is reportedly working on this problem due to concern of nixie tubes being used on submarines. Units were modified by installing non-mercury type nixie tubes and replacing them after each mission.

Cadmium-plated chassis were identified in three units; however, cadmium may be more extensively used than the evaluation suggests since some overfinishes concealed the base material. Several other units have small cadmium-plated parts such as brackets. Absolute prohibition could significantly increase modification costs for some equipment unless the modifications can be made during original buildup. An effort should be made to determine whether cadmium plating can be acceptably modified by a special process.

Other specialized replacements included change of asbestos insulation in the gas chromatograph and rubber bumpers/feet on certain bench-mounted units.

3.3.2.6 Zero-G Compatibility

It was found that the blood cell counter, freezer and dewar were functionally unsuitable for zero-g operation (see Section 3.2.3). Also, there were other very extensive modifications for the counter and freezer. Therefore, new functional approaches were deemed necessary.

Several other items require less extensive modifications to operate in a zero-g environment. The microscope requires addition of a spring to replace gravity seating of the barrel. The biomedical items require redesigned liquid sample handling/container methods. Basically, the electronics associated with the sample sensors are suitable for zero g. Strip chart recorders, such as the CRT recorder and the densitometer in the electrophoresis system, need chart paper guides/rollers to constrain output paper from "flapping" around and/or getting out of the existing guides.

3.3.2.7 Item Mounting and Integrity

Except for input probe accessories and latches on some doors, the items appear well suited to withstand 9-g acceleration forces. Cathode ray tubes are generally well supported, but Lexan overlays and screening were added to eliminate hazards if implosion should occur.

It was determined that many catalog rack guide mounts may be inadequate since they are seldom intended for the vibration or acceleration anticipated. Guide rails with a rear lock-down feature were added for dynamics considerations and to facilitate installation and removal on the ground. If the item face panel fits a 19-in. (48.2 cm) rack, the rail kit was the recommended modification. If a face panel bracket or other feature is also needed (i.e., to fit a NIM bin), then a new mounting design was defined. This distinction is made in the "Item Mounting" column of Table 3-4.

3.3.2.8 Audible Noise

Noise was considered to be a problem only if continuous during item operation. Noisy fans were replaced with quieter units. These replacements may have been conservative since the units were evaluated in the open by themselves. An alternative to replacing noisy fans would be to supply equivalent cooling air, from rack forced-air sources, directly to the fan inlets. The fans could be removed. However, provisions for tight connections and flex ducting, or other design solutions, would be necessary.

Since forced air is provided by the rack to each item, few modifications for cooling were required. Units with dead air spaces were modified with additional porting, or air deflectors.

Units using fans expel air from the sides and/or top of the units. Baffles and spacing in the rack will be needed to ensure proper cooling air distribution.

3.3.2.9 Particle Containment

The main sources of possible shattered particles were CRT's. A few glass tube "auto" fuses in open clips were also identified. Exposed powdered iron core transformers and meter faces were other possible fragment sources. Containment of shattered and flaming particles (in case of breakage or ignition) was accomplished by covering all openings with clear-viewing plastic or mesh screen, as appropriate. Optimum screen mesh size is subject to further definition. The allowable size of escaping solid particle, the size versus burning time (or travel distance) of flaming droplets, and the stopping mechanism of a filter screen, bear upon this optimization. Such an analysis is beyond the scope of this study.

In general, ported cases use holes of about 3/16 in. (0.48 cm) diameter, or slots 1/8 (0.32 cm) to 3/16 in. (0.48 cm) by 1.0 in. (2.5 cm) long, for which screens were assumed to be needed. Internal screens were applied to openings in CRT shields, which generally completely enclose the CRT except for the socket end and the high-voltage cable entrance.

3.3.3 Summary of Recommended Processes

Special processes recommended for meeting SEEIR are as follows:

1. Coat circuit boards with a conformal coating to encapsulate unknown materials and provide additional support to lead-mounted components.
2. Spray painted surfaces, prone to flaking, with clear adhesive/binder.
3. Bake out unit to drive off trace contaminants.

In addition, tests are required in general for:

1. Verification of modification adequacy.
2. Tests may also be required to resolve whether a marginal design is adequate or requires modification.

3.3.3.1 Flame Retardant

Use of a flame-retardant coating is primarily intended for electronic circuit boards. Certain fiberglass-epoxy boards are non-flammable and, when so identified, may not require treatment unless undue concentrations of potential flammable plastic parts are loaded on them. In general, except for high-voltage CRT circuits, all items that were evaluated operate on low circuit voltages (+24 vdc or less) once past the input power transformer. All electrical items were fused. No items contained volatile chemicals. Therefore, it appears that the risk of ignition within almost all off-the-shelf items is very low. The presence of high voltages and bundled wiring somewhat increases the risk due to the increased chances of heat-generating leakage current paths.

An encapsulant may also be applicable to transformers. Many commercial items utilize open-frame construction with a slightly higher risk of generating smoke (if not flames) from small amounts of exposed "fish" paper and insulation varnishes. Assuming that heat dissipation is not adversely affected, these transformers could be coated in areas of exposed windings/fishpaper. The effect and suitability of this process should be among test objectives on representative equipment items purchased to verify modification feasibility and ultimate suitability for Spacelab use.

As noted previously, a conservative modification costing stance was adopted; therefore, all electronic boards are assumed to be coated with flame-retardant material, except in totally enclosed units where combustion would be limited. However, inasmuch as most boards also require encapsulation for vibration, there is no added cost delta for flame retardant.

3.3.3.2 Vibration Encapsulation

Current judgment is that most off-the-shelf electronic circuit boards cannot withstand the random-plus-sinusoidal vibration specifications placed in SEEIR without failure. Some relatively large-mass lead-mounted components may pull loose from circuit board solder joints or leads may break. A solution to supporting small parts is to coat the circuit boards with an encapsulant that lends lateral and longitudinal support to small components. Adjustments, however, should not be coated for continued usefulness (separate staking of these adjustments may be necessary).

The encapsulant should be tough, but not brittle. Pliability should compensate for differences in thermal coefficients of expansion for the encapsulant and the encompassed parts to prevent component breakage and adverse shifting of components with temperature. Shifting could affect performance in some RF circuits. Other requirements are good dielectric qualities, inertness to adjacent materials, low outgassing after curing, and non-flammability. Consultation with material experts indicates that acrylic or epoxy encapsulants with good outgassing and flammability characteristics can be provided.

To summarize, further study is recommended to:

1. Finalize vibration input specification to equipment items.
2. Develop "cook book" guides, including simple and inexpensive tests, to evaluate individual equipment modification needs.

3.3.3.3 Overspray (Painted Surfaces)

A clear overspray was assumed to be satisfactory for prevention of flaking of painted surfaces. This process avoids the necessity of reapply printed matter. In actual practice it appears that there will be little need for this process, owing to the relatively short missions and benign atmosphere that minimizes degradation of the painted surfaces. Further, it appears to be standard practice to use high-quality paints and applications on quality equipment items. These appear compatible with a Class 100,000 atmosphere or better. Most items inspected had painted front panels only. It is likely that if a flaking problem developed on a unit, it would generally be only after a number of missions. Pre-mission inspections and post-mission evaluations would avoid excessive problems on subsequent missions.

3.3.3.4 Expelling Trace Contaminants

Equipment containing organics, particularly plastics, may exude relatively large emissions of solvents, plasticizers, and other chemicals used in manufacturing processes. While some continued offgassing will occur as chemicals work to the surface of materials, it decreases significantly with time. Another source of offgassing is entrapment of solvents and cleaners that may be applied to assemblies or components before, during and after manufacture.

A recent study for NASA by Beckman on four commercial equipment items indicated that 150 hours of operation reduced outgassing to a relatively low, essentially constant level. This process was adopted as a general requirement for SEEIR.

3.3.3.5 Testing of Modified Hardware

It will be necessary to requalify/accept items that have been modified. Unless the SEEIR vibration values are reduced, it appears prudent to accomplish some type of vibration tests to verify ability to survive launch. Functional tests operating in six orientations have been successfully performed on the Apollo program to verify zero-g operational compatibility. A repeat of all normal functional tests will be required to assure no unintentional damage during the modification process. The cost analysis included necessary estimates for these activities.

3.4 COST ANALYSIS

3.4.1 Cost Analysis Approach

The approach employed to determine the cost of modified and custom-built hardware is pictorially depicted on Figure 3-28.




Figure 3-28. Cost Analysis Approach

Two separate costing activities were accomplished, leading to the selection of the cost-effective procurement approach for each item. Analytical techniques were different for both activities. Modified hardware costs were determined by a "grass roots" analysis which estimated the cost of each activity associated with the modification. Costs for custom-built hardware were estimated based on the traditional cost estimating relationship (CER) approach.

Analysis of modified hardware costs is dependent on the engineering description of the modification required to make each item suitable for installation aboard the Spacelab. Based on this information, engineering, manufacturing, test, and other costs were determined. These cost elements are summed up according to their location in the work breakdown structure defining the total modification costs. Custom costs were determined from either adjusted CER's or vendor quotes. Custom cost estimates were received from vendors for 6 of the 34 items.

3.4.1.1 Modified Equipment Costs

Modified hardware costs were determined by estimating the elemental cost incurred by the modification activity on each unit. The work breakdown structure for the costs of hardware modification is shown on Figure 3-29. Each cost category is identified along with the cost elements that make up that category. Costs were defined by cost element for engineering and manufacturing. Costs for verification testing includes testing, test fixtures, and operations. New test facilities were assumed not required. Documentation and project management estimates could not be determined at the level of their cost elements. These costs were determined at the category level, including the elements listed.




Figure 3-29. Work Breakdown Structure

In order to facilitate the detailed estimates required, engineering researched the available data on the selected hardware by reviewing sketches, drawings, photographs, and specifications; and visually examined as many units as possible. Drawings and narratives were prepared describing the necessary modifications needed to each piece of equipment to make it acceptable for use in the Spacelab. These data, along with photographs, were forwarded to Beckman Instruments for use in developing the manufacturing estimates.

Verification testing includes only tests to verify that the modifications performed on a given unit make the unit suitable for installation in the Spacelab. Characteristics analyzed to be suitable do not require verification testing.

Documentation costs include the cost of documentation not included in the engineering, manufacturing and testing activities. Documentation costs in these categories vary with modification activities. Costs in the documentation category cover general documentation requirements which tend to be more constant and less dependent upon the amount of modification undertaken.

Project management costs were considered to vary with the modification activity.

The sum of the costs by category define the cost of the modification activity. The total cost of modified hardware equals the retail cost plus this modification cost.

Engineering. Engineering hour estimates were generated by the responsible engineer, based upon his knowledge of the modification involved and the necessary engineering tasks. The estimates generated for each task were compared between items of equipment in order to reflect consistency from one to another. Where possible, estimates were based on historical data from other programs, such as an hours-per-drawing factor extracted from Saturn S-II data. A sample of the determination of engineering cost estimates is shown in Figure 3-30. The Tennelec Timer (Model TC 545), a supporting amplifier module, and a power supply module are used for the example.

A matrix of engineering tasks was employed to ensure consistency of engineering estimates. Seventy-one different tasks were identified for the 34 selected items. Baseline hours were defined for each modification task as shown in the Hours/Task column. Then, the modification actions for each package were determined. The number shown for a given modification-type represents the number of times that different engineering activities are required for the same type of modification. If, for example, more than one circuit board in a package required a conformal coating but the boards are similar, the number used for engineering estimates could be reduced. The total hours for each task were determined to indicate areas where tradeoff of relaxed specification requirements might be justified as an approach to reducing modification costs.




Figure 3-30. Engineering Cost Estimate
(Modified Equipment)

The engineering estimates for the Tennelec timer are shown. It is estimated that 60 hours would be consumed analyzing and inspecting the equipment to determine the modification required to make the equipment suitable for the Spacelab. The circuit boards in the scaler amplifier and the timer modules are very similar in mounting and configuration, indicating that engineering for both would be about the same. Therefore, only 40 hours were spent to define a rugged circuit board installation, and 20 hours were spent to assure that the conformal coating would have no adverse effect on the circuit board component. Ten hours were spent to define the bakeout procedure used to volatilize contaminants possibly present in the package. Only 6 hours were spent changing the structural fasteners to a locked type that will not shake loose in vibration; 40 hours were spent defining guards for protrusions on the package; 80 hours were spent to define a mounting reinforcement for a cantilevered trimpot on the amplifier module which could come loose in vibration. Six knobs were replaced because they were made of plastic. Since it is estimated that approximately two hours are required to define a source for new knobs of the proper configuration, a total of 12 engineering hours were estimated for this activity. Forty hours are required to provide support for two lead-mounted capacitors in the power supply.

Engineering hours directly related to the modification and process activities defined for the Tennelec unit totaled 312 hours. Preparation of the verification test procedure required an additional 40 hours, bringing the total hours for engineering to 352 hours. The engineering cost was then determined by multiplying these hours by an average engineering labor rate.

Fabrication. Beckman Instruments performed the fabrication and associated manufacturing cost estimates for all modifications. The manufacturing engineering department utilized standard data, wherever possible, based on historic information from instruments presently being manufactured. Consideration was given to the selection of fabrication methods and assembly techniques which were most appropriate and least costly for single-unit production to accommodate proper work planning and the setup time involved.

Manufacturing engineering performed a detailed review of each modification on an operation-by-operation basis, including disassembly and reassembly, where required. Estimates were generated for any and all fabrication, tooling, test, inspection and material costs involved. Hourly estimates for labor were prepared, extended by appropriate direct labor rates, burden and administrative expense and summed to a total. The same approach was used to develop the material cost associated with the modification. A separate work sheet was generated for each piece of equipment selected in the modification program.

A typical worksheet used to estimate manufacturing costs is shown in Table 3-5. Beckman Instruments performed this activity to assure that manufacturing hour estimates conform to experience in the instrumentation industry.

Table 3-5. Manufacturing Cost Estimate (Modified Equipment)

ASSEMBLY NO. <u>Tennelec TC545A</u>		ASSEMBLY NAME <u>Counter Assembly</u>		NEXT ASSY NO. _____		DATE <u>5/1/74</u>								
I	PART NO.	PART DESCRIPTION	QTY	LABOR HOURS							COST \$			
				FAB	ASS'Y	INSP.	TEST	INFL	E	MFG. TECH.	MAT'L	LABOR	TOTAL	
I	A	Protrusion protection	2	2.00			20					10.00		
	B	Fab. and install PCB stiffeners	1	1.50		50	20					7.50		
	C	Secure trimpot	1	1.00			10					3.80		
	D	Secure capacitors	2	.40	1.00	40						1.20		
	E	Epoxy on heat sinks	2			10	10					.30		
	F	Replace knobs	8		1.50	60						12.00		
	G	Replace fasteners	30		1.00	50		1.00				2.00		
	H	Conformally coat PCB	1		1.00	1.00		1.00	50			2.00		
II	A	150-hour bakeout					10.00	10.00				350.00		
	B	Test and calibration					5.00		5.00					
TOTAL				4.90	5.10	3.10	15.00	12.00	5.50			388.80	656.96	1045.76

Hours for each type of manufacturing activity are identified for all modifications. Manufacturing activities included fabrication, assembly, inspection and test activities and also include additional hours for a manufacturing engineer and a manufacturing technician support.

The Cost columns show the material costs associated with each modification activity and the extension of the labor hours into dollars. The sum of the activity costs leads to a total manufacturing cost of \$1045.76 for the Tennelec modules.

Individual costs are self-explanatory except the \$350 for the 150-hour bakeout. This process is assumed to be performed by an outside vendor and shows as a material cost. Labor hours identified support this testing activity.

Test. The cost of development testing was estimated by engineering in conjunction with test laboratory personnel. Consideration was given to the type of environment to which the equipment would be subjected as well as its basic functional requirement. Each piece of equipment was reviewed individually for such testing as shock, mechanical vibration, accoustical vibration, EMI generation, and off-gassing. Estimates were generated based upon the complexity of the unit and the amount of modification required.

Documentation. Documentation costs were estimated by factoring off the total of fabrication, engineering and test costs and proportioning this cost equally among all equipment items. The documentation effort represents the fixed programmatic documentation effort. Its cost does not vary with equipment costs. The factor used to develop the documentation estimate was based upon prior historic cost data from other in-house programs. An example of the formula as used is as follows:

Total Manufacturing/Engineering/Test	\$669,000
Documentation Factor	10 percent
$\$669,00 \times .10 = \$66,900$	
<u>\$66,900</u>	
31	= \$2153 (rounded to \$2160) cost to each unit

Program Management. The program management effort was also estimated by a factor developed from prior in-house cost experience. The factor was applied to total cost on the theory that management will pay more attention to costly items because (1) their cost reflects the number and complexity of the modification actions, and (2) a greater financial risk is incurred with the higher cost items. Tabulated below is the formula used to cost program management effort:

Selected Unit - Tennelec Timer	
Total Manufacturing/Engineering/Test/Data Management	\$12,158
Program Management Factor	5 percent
$\$12,158 \times .05 = \607	

The hour estimates generated for engineering and test effort were extended by appropriate rates which included labor, burden and administrative expense for the 1974 time frame. Any material required, such as test fixtures, was costed in total by adding procurement and administrative expenses to the basic estimate. The total cost of each modified off-the-shelf unit was determined by adding the retail cost to the cost of modification. Total modified hardware costs are reported in Volume III in the same manner as shown for the Tennelec Timer in Table 3-6.

Table 3-6. Modified Tennelec Timer Costs

Manufacturing	\$ 1,129
Engineering	4,453
Test	4,516
Documentation	2,160
Program Management	607
Total Modification	12,765
Basic Cost	1,500
Total Cost	\$14,265

3.4.1.2 Custom-Built Equipment

Parametric pricing techniques were employed to cost the new development hardware. This involved a comparison of technical characteristics and cost of the new hardware to that of existing or previously defined equipment where similar history was available. Cost projections were accomplished by the use of cost estimating relationships (CER's), knowledge of the technical characteristics of both programs, state-of-the-art position of the hardware, and a complexity analysis.

A list containing the hardware identified as representative equipment was compiled and reviewed to determine its similarity to a known piece of equipment or system. System in this sense pertains to a group of items or an assembly within which the equipment under consideration would be an integral part. A matchup was made to information contained in the data bank with regard to weight, technical characteristics and function. In some cases, it was more prudent to select multiple-source data in order to establish the more reliable hardware.

Once the source data were identified, a more in-depth review was made to establish a complexity relationship. While this complexity analysis can become subjective, the technical information contained in the data bank (if properly utilized) aids the engineer in preparation of a more objective result. To assess the complexity required in this study, the cognizant engineer was given the opportunity to apply his knowledge to an in-depth technical data source. This enabled the engineer to realistically ascertain the relative complexity of the hardware under consideration as compared to the previously built and/or designed equipment. The complexity multiplier reflects the relation of the design complexity of the selected hardware item to the design complexity of the benchmark unit used for the CER basis. In some cases the relative weight of the selected and benchmark units also reflect complexity, requiring the engineer to normalize his estimate on a complexity-per-pound basis. To finalize the CER used in this study, a state-of-the-art calculation was applied. The impact of this cost driver in regard to development cost is the state-of-the-art "know-how" status at the point in time of the new item of equipment versus that at the time of development of the benchmark unit.

The four numerical ratings as listed in MF003N¹ and an additional category called "off-the-shelf hardware" were used in making the state-of-the-art determination. They reflect the variation in degrees of definition from beyond the state of the art to off-the-shelf hardware. The designation graduations of 20 percent of the 5 categories, on a scale from 20 to 100, were applied to further adjust the CER used in the study for each individual equipment item selected. The state-of-the-art rationale and its percentage application are reflected in Table 3-7.

¹Document, Program Study Cost Estimates, revised 1 March 1973, NASA-MSFC.

Table 3-7. State-Of-The-Art Rationale

Description	Percent
OFF THE SHELF. The equipment can be placed on board in an as-is configuration.	20
MINOR MODIFICATION. Equipment requiring packaging modifications but no functional design changes are necessary. Unit retains its basic configuration. Most electrical hardware falls into this category.	40
REDESIGN REQUIRED. New configuration is necessary to be suitable for space operations. Basic operating principle remains the same; includes straightforward solutions to replacement of gravity-dependent functions. (The Beckman spectrophotometer is a unit requiring ruggedizing as well as sample containment in zero g.) The latter problem can be solved by evacuated sample containers and syringe sample injection.	60
SLIGHTLY BEYOND STATE OF THE ART. New configuration is necessary with improvements in functional capability requiring development. Zero-g operation requires new approach. Laboratory prototypes may exist. For example, a working prototype of the Coulter blood cell counter has been developed by Beckman for NASA.	80
SUBSTANTIALLY BEYOND STATE OF THE ART. Designs meeting requirements have not yet advanced beyond conceptual level. As an example, approaches for two-phase fluid containment for a zero-g cryogenic dewar have only been hypothesized on paper. Electrophoretic retention of the liquid phase material is one such suggested approach.	100

Figure 3-31 illustrates the process by which the custom cost of an item was determined. The Tennelec Timer is again used as an example. Design development, test and engineering (DDT&E) costs are coupled with the cost to manufacture the first unit to arrive at a total first article cost. The elements of the cost analysis were determined as follows.




Figure 3-31. Sample Custom Cost Estimate (Tennelec Counter/Timer)

1. The Space Station central timing package was selected as the most relevant hardware item for which cost data existed in the data bank. The cost per pound for DDT&E and the first unit of this package are 14,641 and 13,306 dollars per pound.
2. The estimated weight of the custom-built item was then determined by evaluating how much of the unit could be reconfigured into the more compact-type packaging used in the aerospace industry. The commercial modular assembly weighs 14 pounds. Engineering estimated that the weight could be reduced to 4 pounds. The existing unit has considerable structural weight that can be eliminated. Also, much of the power conditioning can be eliminated. These modules normally use ± 12 volts or ± 24 volts dc, which can easily be obtained from the 28 volts dc on-board power supply, thereby eliminating the integral power conditioning currently required to change 110 volts ac to ± 12 volts and ± 24 volts dc.
3. Since minor packaging changes are required to make the unit consistent with aerospace packaging design practices, a $2/5$ factor was assigned to the state-of-the-art multiplier. Since the CER used was for hardware at a conceptual level, a multiplier of 1 would have been used if the development status of the subject timer would also be conceptual. If the subject timer were in the prototype stage, a $4/5$ multiplier would have been used. If major modifications were required (such as development of a gravity-insensitive function), a $3/5$ multiplier would have been used. Off-the-shelf hardware would have a multiplier of $1/5$.
4. Complexity is broken down into two types. Nonrecurring complexity represents the relative functional complexity of the reference and selected item. The figure compares functional capabilities of the source and sample hardware elements. The Tennenec unit can only accomplish two of five functions accomplished by the reference unit, leading to a complexity multiplier of 0.4. The recurring cost complexity factor is also 0.4, reflecting the simpler construction and assembly required by this less complex package.

The custom cost then is determined by multiplying the above factors together. The accuracy of these estimates is verified by comparing the values to those derived on other studies and vendor estimates. Such a comparison exists for the Fluke frequency synthesizer. Fluke personnel indicated that their Mil-Spec unit cost approximately \$500K to develop. The analytical approach used in the study estimated the cost for a new frequency synthesizer to be \$452K, indicating good correlation.

3.4.2 Modification Cost Summary

The results of the detailed cost estimates for each of the modified selected hardware items is presented in this section. Table 3-8 identifies the modification cost of each item by its work breakdown structure category, and totals the costs incurred for the selected items in each category.

Table 3-8. Modified Equipment Cost

Selected Equipment	Mfg.	Eng.	Test	Doc.	Program Mgmt.	Total
Tape Recorder						
Ampex	6,034	8,979	4,416	2,160	1,616	23,205
Honeywell	14,620	15,548	4,416	2,160	1,616	38,360
Computer	20,907	17,333	5,888	2,160	2,228	48,516
Keyboard Display	4,979	19,136	5,888	2,160	1,565	33,728
Timer						
Datatron	3,029	7,452	4,416	2,160	844	17,901
Tennelec	1,129	4,453	4,416	2,160	607	12,765
Spectrophotometer	17,135	17,885	8,832	2,160	2,239	48,251
EMI-Intens. Meter	3,367	11,261	4,416	2,160	1,064	22,268
Transmitter						
RHG	835	2,392	2,944	2,160	417	8,748
Collins	1,302	3,128	2,944	2,160	477	10,011
Signal Conditioner	3,611	6,100	2,944	2,160	741	15,556
Spec. Wave Analyzer						
Hewlett-Packard	7,891	17,056	4,416	2,160	1,576	23,099
Singer	3,040	11,684	4,416	2,160	1,065	22,365
Fluke	2,941	7,250	4,416	2,160	838	17,605
Nuclear Data	5,665	15,070	4,416	2,160	1,366	28,677
Laser Assembly	4,801	10,488	4,964	2,160	1,121	23,534
Gas Chromatograph	8,163	17,811	5,988	2,160	1,661	35,783
Microscope	1,330	6,090	2,208	2,160	589	12,377
Strip Chart Recorder	6,882	11,297	4,416	2,160	995	25,750
Volt-Ohm Meter						
8200	4,652	12,935	2,944	2,160	1,056	23,747
8125	1,535	5,778	2,944	2,160	621	13,038
Oscilloscope						
Tektronix	9,352	13,119	4,416	2,160	1,328	30,375
Hewlett-Packard	8,640	7,139	2,944	2,160	965	21,848
Electrophoresis	14,280	49,000	35,872	2,160	5,016	106,328
Furnace	1,216	11,206	4,416	2,160	950	19,948
Power Supply						
Power Design	1,728	6,109	4,416	2,160	680	15,093
Sorenson	1,996	7,360	2,944	2,160	723	15,183
Centrifuge	1,014	3,956	1,472	2,160	391	8,993
Particle Counter	2,259	7,250	2,944	2,160	745	15,358
Microtome	1,751	3,662	2,208	2,160	476	10,257
pH Meter	2,197	9,586	5,888	2,160	946	20,777
Total	168,281	347,513	160,168	66,960	36,522	779,444

The equipment selected for review is grouped by category of equipment, so that where more than one of a similar item has been analyzed, a rapid comparison can be made of the difference in modification cost. The higher the cost, in a comparison of this nature, the less flightworthy the item was in its original configuration. As an example, the Tennelec timer required approximately two-thirds the modification of the Datatron timer because its basic design was more vibration-proof than was the Datatron model. Each of the differences, accounting for the cost variations, is discussed in depth in Section 2.0, Volume III.

Figure 3-32 shows the distribution of costs of modified available hardware. Darkened portions of each bar represent the retail cost for each unit, while the open part of the bar reflects the estimated cost of modification. Estimated modified hardware cost can be read from the left-hand ordinate. The cumulative quantity of modified hardware costs, below given cost plateaus, is shown on the right-hand ordinate. For example, 28 of 31 items can be purchased and modified for less than \$50,000

Statistical information is presented in the center of the figure. The average retail, modification, and modified hardware costs apply to the 31 items analyzed. The average of modification cost, divided by retail cost, was determined by calculating this ratio for each item and finding the mean of these values. This value is different from the ratio of total modification and retail costs. It was determined from the ratio of the sums of both costs. Care must be used applying these ratios to other groups of hardware. The average of all cost ratios tends to emphasize ratios for items with low retail costs. The ratio of total costs tends to emphasize the high-cost items. Together, they provide an indication of the range of values expected for modification activities.




Figure 3-32. Modification Versus Retail Cost

Equipment analyzed was divided into four classes--airborne, Mil-Spec, NIM, and commercial. Figure 3-33 indicates the relative cost to modify different types of equipment. All data shown reflect the same trend. Modification costs are the least for airborne equipment. They are somewhat higher for Mil-Spec equipment, still higher for NIM equipment, and highest for commercial hardware. These data are consistent with the general impressions developed during the suitability analysis.




Figure 3-33. Design Specification Impact on Modification Costs

A curve plot of the data is shown along with modification cost ratios previously used. The curves represent the best fit to the data which vary over broader bands than indicated by the thin lines.

Airborne equipment has the lowest modification costs because they tend to have rugged construction, are sealed preventing atmosphere interchange, and operate at altitudes up to 55,000 feet (16.5 km), requiring minimum reliance on convection for cooling. Some Mil-Spec units display similar characteristics such as the Fluke Mil-Spec voltmeter built for submarine installation, however, others meet military specifications not as important to Spacelab installation. For example, the Hewlett-Packard military oscilloscope is built for field use so it is protected from all types of hostile environments. It still has the same vibration inadequacies as its commercial counterpart. NIM units show several good characteristics such as approved material usage and a design compatible with forced-air ventilation. However, these units have poor vibration characteristics. Volumetric efficiency is also very poor with these items.

The curve displays a common characteristic among all design types. As the retail cost decreases, the ratios increase. This effect is a result of certain types of modifications or processes being common with all units such as a trace contaminant bakeout or some ruggedizing for vibration. Also, certain fixed costs such as documentation and engineering to determine the need for modification cause the modification cost to diminish more slowly than the retail cost. As a result, all ratios are higher at lower retail costs.

The major cost driver to equipment modification was the lack of a gravitational field in space. The three units which could not be modified depended on gravity to function. Two of these units, the refrigerator and dewar, require completely new design approaches to operate in zero gravity. The third unit had additional inadequacies for operation in the Spacelab. The Coulter counter uses gravity as the motivating force for sample flow. A falling mercury column reduces the pressure behind an orifice causing the sample to flow through the orifice. It also had little structural strength and used glass throughout. The refrigerator and the dewar employ gravity to separate gas and liquid fluid phases. Such separation maintains liquid in the evaporator and allows only gas to enter the compressor of the refrigerator. The same separation keeps liquid in the dewar and allows only gas to escape the vent located at the top of the dewar. Totally new developments are required in all three cases.

Modifications were identified for the units having gravity-dependent functions (see Figure 3-34). Modifications vary from development of a new batch-type electrophoresis cell to adding a few retaining screws to the turrets of a microscope. The retail and modification costs as well as the ratio of modification and retail costs are shown on the figure to indicate the range of values determined during the study.

SELECTED ANALYZED UNITS	COST (\$1000)		(MODIFICATION) RETAIL
	RETAIL	MODIFICATION	
SPECTROPHOTOMETER	2.7	47.0	17.5
ELECTROPHORESIS	6.8	113.0	16.7
pH METER	0.6	20.0	33.4
MICROSCOPE	1.6	12.4	7.9
CENTRIFUGE	0.3	8.2	25.3

• SUMMARY OF ALL EQUIPMENT

(MODIFICATION COST) RETAIL COST	GRAVITY DEPENDENT	GRAVITY INDEPENDENT
AV	20	5.4
(TOTAL MODIFICATION COST) TOTAL RETAIL COST	20.1	3.1

Figure 3-34. Gravity Dependence Impact on Modification Costs

When the tabular data are expanded to include all the units analyzed, the cost impact is quite clear. Either ratio shows a 20-to-1 cost multiple for gravity-dependent units, indicating good correlation over all retail price ranges. When this value is compared to the same ratio for units which operate independent of gravity, it is clear that equipment dependent upon gravity incur significantly higher modification costs than those that are independent.

The impact of the basic operating principle on modification costs was assessed. Figure 3-35 shows that chemical analysis hardware is the most costly to modify. It also indicates that mechanical and electrical hardware is the least costly.

Different conclusions can be drawn from the ratio data, depending upon which ratio is used for the electrical, mechanical and optical equipment. Therefore, conclusions are deferred to more in-depth study. The conclusion is clear, however, for the chemical analysis hardware. It is the most costly to modify.



$\left(\frac{\text{MODIFICATION COST}}{\text{RETAIL COST}} \right)_{AV}$	ELECTRICAL	MECHANICAL	OPTICAL	CHEMICAL
	7.0	8.4	9.0	18.0
$\left(\frac{\text{TOTAL MODIFICATION COST}}{\text{TOTAL RETAIL COST}} \right)$	3.0	2.3	5.8	11.5

• CUSTOM-BUILT HARDWARE, EITHER CHEMICAL ANALYSIS INSTRUMENTATION OR SUPPORT EQUIPMENT

CHEMICAL ANALYSIS HARDWARE	COST (\$1000)		$\left(\frac{\text{MODIFICATION}}{\text{RETAIL}} \right)$
	RETAIL	MODIFICATION	
SPECTROPHOTOMETER	2.7	47.0	17.0
GAS CHROMATOGRAPH	8.0	35.0	4.4
ELECTROPHORESIS APPARATUS	6.8	105.0	15.5
pH METER	0.6	20.0	34.8

Figure 3-35. Operating Principle Impact on Modification Costs

The high cost of modification is, for the most part, due to reliance upon gravity for many functions and construction suited only to the benign environment of the laboratory. The gas chromatograph, an exception, demonstrates that chemical analysis hardware, operating independent of gravity and of more rugged construction, can have a low modification cost ratio more typical of electrical hardware.

The dewar, refrigerator, and Coulter counter amplify the high-cost conclusion regarding chemical analysis hardware. The Coulter counter is a chemical analysis instrument. The dewar and refrigerator, in general, support chemical or biological analysis activities. None of these three could be used in a modified form.

3.4.3 Cost of Custom Equipment Summary

The custom, or new development hardware cost, has been summarized in Table 3-9. An explanation of the columnar headings is given below.

Weight. The projected weight of the equipment with the modifications taken into account.

CER. The cost estimation relationship, reflected in dollars per pound after consideration was given to weight scaling, complexity and state-of-the-art of the base CER.

Complexity Factor. The complexity of engineering, fabrication and test of the item as related to the benchmark equipment from which the CER was derived.

State-Of-The-Art. The numerical designation, as discussed in Section 3.4.1, reflects the relative state-of-the-art of the new item when compared to the benchmark hardware. The number designation is converted to a percentage for formula use.

Source Data. The source of the applicable CER for each unit is presented in this column.

Table 3-9. New Hardware Development Cost

SELECTED EQUIPMENT	WEIGHT		CER	COST	COMPLEXITY	STATE OF THE ART	SOURCE DATA
	LB	(KG)	\$/LB	MILLION \$			
TAPE RECORDER - AMPEX	80	(36.3)	2,500	0.200	0.7	1	AMPEX QUOTE
- HONEYWELL	80	(36.3)	2,500	0.200	0.7	1	AMPEX QUOTE
COMPUTER	119	(54.0)	1,782	0.212	1.0	2	SHUTTLE ORBITER
KEYBOARD DISPLAY TERMINAL	30	(13.6)	16,767	0.503	1.0	2	ORBITER AVIONICS
BLOOD CELL COUNTER	15	(6.8)	27,800	0.417	2.0	4	BECKMAN - BCC
TIMER - DATATRON	9	(4.1)	19,222	0.173	1.0	2	ORBITER - TIMING UNIT
- TENNELEC	4	(1.8)	7,750	0.031	0.4	2	
SPECTROPHOTOMETER	30	(13.6)	16,867	0.506	1.0	3	UV GRATING SPECTROMETER
EMI FIELD INTENSITY METER	25	(11.3)	7,480	0.187	1.0	2	LUNAR ORB STUDY, RADIO NOISE
TRANSMITTER - RHG	20	(9.1)	12,900	0.258	1.0	2	SHUTTLE ORBITER - COMMUNICA.
TRANSCEIVER - COLLINS	17	(7.7)	12,900	0.219	1.0	2	SHUTTLE ORBITER - COMMUNICA.
AMPLIFIER - NEFF	16	(7.3)	3,188	0.051	0.7	2	SHUTTLE ORBITER - AVIONICS
SPECTRUM ANALYZER - HEWLETT-PACKARD	40	(18.2)	12,900	0.516	1.0	2	SHUTTLE ORBITER - COMMUNICA.
- SINGER	35	(15.9)	12,914	0.452	1.0	2	SHUTTLE ORBITER - COMMUNICA.
FREQUENCY SYNTHESIZER - FLUKE	35	(15.9)	12,914	0.452	1.0	2	SHUTTLE ORBITER - COMMUNICA.
MULTICHANNEL ANALYZER - NUCLEAR	45	(20.4)	12,911	0.581	1.0	2	SHUTTLE ORBITER - COMMUNICA.
LASER	30	(13.6)	5,700	0.171	1.0	2	LASER ALT. - LUNAR ORB STUDY
GAS CHROMATOGRAPH	30	(13.6)	13,333	0.400	1.0	2	CSM - ECS
MICROSCOPE	20	(9.1)	1,800	0.036	1.0	2	SPACE STA EXP 5.13-SE
STRIP CHART RECORDER	25	(11.3)	14,520	0.363	1.0	2	CSM - DISPLAY & CONTROLS
VOLT/OHM METER - 8200	10	(4.5)	4,000	0.040	1.0	2	
- 8125A	10	(4.5)	4,000	0.040	1.0	2	SPACE STA EXP 5.22-SE
OSCILLOSCOPE - TEKTRONIX	20	(9.1)	16,300	0.326	2.0	2	SHUTTLE ORBITER - COMMUNICA.
- HEWLETT-PACKARD	10	(4.5)	14,900	0.149	1.0	2	SHUTTLE ORBITER - COMMUNICA.
ELECTROPHORESIS APPARATUS	30	(13.6)	6,200	0.186	1.5	4	SPACE STATION EXP 5.13
FURNACE	100	(45.4)	4,530	0.453	0.2	3	SPACE STATION EXP 5.16
DEWAR	5	(2.3)	262,000	1.310	1.0	5	BEECH CRYOGENIC TANKAGE
POWER SUPPLY - SORENSON	50	(22.6)	520	0.026	0.2	2	SHUTTLE ORBITER EPS
- POWER DESIGNS	10	(4.5)	2,400	0.024	0.3	2	SHUTTLE ORBITER EPS
CENTRIFUGE	25	(11.3)	960	0.024	0.2	3	SPACE STATION EXP 5.13C
REFRIGERATOR/FREEZER	50	(22.7)	2,540	0.127	1.0	4	SPACE STA EXP 5.13/5.16
PARTICLE COUNTER	6	(2.7)	34,833	0.209	0.4	2	LOW-ENERGY PARTICLES
MICROTOME	50	(22.7)	620	0.031	0.2	2	SPACE STATION EXP 5.13
pH METER	1.5	(0.7)	56,667	0.085	1.0	3	BECKMAN CSM - pH ANALYZER

3.4.4 Comparison of Modified Hardware Cost to Custom-Built Hardware Cost

This section summarizes and compares the cost of hardware being modified or custom built. It provides a visual overview at the total cost level and rapidly identifies the wide range of differences between the two concepts.

The costs in Table 3-10 are reflected in millions of dollars and the columnar headings are explained below.

Modification Cost. Total cost of item after modification, including basic or retail cost.

Development Cost. The cost of developing or custom building hardware for its specific purpose.

Delta. The net difference in cost between developing a new piece of hardware or modifying an off-the-shelf item.

Comparison of the custom costs generated to vendor information on development costs of like equipment shows good correlation. Vendor development information is shown below.

Fluke Spectrum Wave Analyzer

Present retail cost	\$ 17,200
Original development cost	\$500,000
Ratio	29:1

Ampex Tape Recorder

Present retail cost	\$ 175,000
Development cost	\$5,500,000
Ratio	31:1

Retail and development cost are based on data for Skylab tape recorder.

As a comparison to the above, the study reflects the following:

Average Custom Unit Cost	
Total development cost	\$225,000
Average retail cost	6,300
Ratio	36:1

It is concluded that the costing approach used is reasonable, and that the magnitude of the costs is representative of the costs of custom-built hardware. The slightly higher ratio of the study hardware results from the inclusion of some low-cost items which tend to have high retail-to-custom ratios.

Table 3-10. Comparison of Modification Cost to Custom-Built Cost

SELECTED EQUIPMENT	COST IN MILLION \$		
	MODIFIED EQUIPMENT	CUSTOM EQUIPMENT	DELTA
TAPE RECORDER			
AMPEX	0.058	0.200	0.142
HONEYWELL	0.050	0.200	0.150
COMPUTER	0.062	0.212	0.150
KEYBOARD DISPLAY TERMINAL	0.035	0.503	0.468
BLOOD CELL COUNTER	--	0.417	0.417
TIMER			
DATATRON	0.021	0.173	0.152
TENNELEC	0.014	0.031	0.017
SPECTROPHOTOMETER	0.050	0.506	0.456
EMI FIELD INTENSITY METER	0.038	0.187	0.149
TRANSMITTER - RHG	0.018	0.258	0.240
TRANSCIVER - COLLINS	0.015	0.219	0.204
AMPLIFIER - NEFF	0.028	0.051	0.028
MULTICHANNEL ANALYZER - NUCLEAR	0.043	0.581	0.538
SPECTRUM ANALYZER - HEWLETT-PACKARD	0.048	0.516	0.468
WAVE ANALYZER - SINGER	0.030	0.452	0.422
FREQUENCY SYNTHESIZER - FLUKE	0.037	0.452	0.415
LASER ASSEMBLY	0.035	0.171	0.136
GAS CHROMATOGRAPH	0.045	0.400	0.355
MICROSCOPE	0.014	0.036	0.022
STRIP CHART RECORDER	0.033	0.363	0.330
VOLT/OHM METER			
8200	0.025	0.040	0.015
8125A	0.015	0.040	0.025
OSCILLOSCOPE			
TEKTRONIX	0.035	0.326	0.295
HEWLETT-PACKARD	0.025	0.149	0.124
ELECTROPHORESIS APPARATUS	0.114	0.187	0.072
FURNACE	0.026	0.453	0.427
DEWAR	--	1.310	1.310
POWER SUPPLY			
SORENSEN	0.020	0.026	0.006
POWER DESIGNS	0.015	0.024	0.009
CENTRIFUGE	0.009	0.024	0.015
REFRIGERATOR/FREEZER	--	0.127	0.127
PARTICLE COUNTER	0.019	0.209	0.190
MICROTOME	0.012	0.031	0.019
pH METER	0.021	0.085	0.064
TOTAL	1.005	8.958	7.953

The bar chart shown in Figure 3-36 illustrates the distribution of custom-built and modified hardware costs. Refer to Figure 3-32 for the equipment identification key. The darkened portion of each bar represents the cost of modifying available hardware. Cost magnitudes can be read from the left-hand ordinate, while the cumulative quantity of items costing less than a given cost plateau can be determined from the right-hand ordinates. For example, 22 of 34 items cost less than \$300,000 to develop.




Figure 3-36. Modified Versus Custom Cost

Large savings are evident for many of the selected items. Savings could be achieved for all but 3 of the 34 selected items. Custom costs averaged approximately \$225,000 of which \$195,000 could be saved on the average if available equipment were modified rather than custom building the same hardware. These average savings may be understated since new hardware was not built for three units--two tape recorders and one computer. An asterisk appears above each unit. The tape recorder used on the Skylab could be made available, built to NASA specifications, for \$175K (according to an Ampex quote). Therefore, rather than building a new recorder, the price of this unit was used. The computer price is based on a preliminary quote for a unit planned for use in space.

The three units which could not be cost effectively modified are the Coulter blood cell counter, the REVCO refrigerator freezer, and the Cryogenic Associates dewar (numbers 32 through 34).

Figure 3-37 shows the results of an analysis of the cost effectiveness of modified electronic instrumentation. As a group, it presents the greatest opportunity for savings to the Spacelab program.

Data from several of the units analyzed are shown on this chart. All the units show significant savings. These savings could be expected since their

operation is independent of gravity, they are built for rack-mounting and, in general, for forced-air cooling. Modifications were limited to minor repackaging without requiring any significant developments. The summary modification ratio data indicate good correlation with either ratio.

SELECTED ANALYZED UNITS	COST (\$1000)			(CUSTOM MODIFIED)
	RETAIL	MODIFIED	CUSTOM	
EMI FIELD INTENSITY METER	13.5	38.0	187.0	4.9
SPECTRUM WAVE ANALYZER	12.9	48.0	516.0	10.8
STRIP CHART RECORDER	6.0	28.0	363.0	13.0
OSCILLOSCOPE	4.2	33.0	326.0	10.0
PARTICLE COUNTER	3.0	19.0	209.0	11.0

SUMMARY OF ALL EQUIPMENT

(CUSTOM BUILT MODIFIED) AV	ELECTRONIC INSTRUMENTATION	OTHER
	10.1	5.7
(TOTAL CUSTOM COST) (TOTAL MODIFIED COST)	10.6	6.0

Figure 3-37. Electronic Instrumentation Most Cost Effective to Modify

3.5 SELECTED EQUIPMENT ANALYSIS CONCLUSIONS

The analysis of the 34 selected equipment items generated the following conclusions. Modifications required to make the equipment suitable for operation in space were not extensive. Vibration integrity and use of allowable materials were a common modification cause. The simple modifications recommended to make these equipment items suitable for operation in the Spacelab increased the cost of the equipment from 4 to 8 times. These costs are considerably less than the same equipment custom built. Custom building this already developed hardware cost an average of 36 times more than the retail hardware.

4.0 SPECIFICATION ANALYSIS

4.1 INTRODUCTION

NASA equipment specifications are being developed consistent with the practice on previous manned spacecraft programs such as Apollo and Skylab. These specifications tend to be rigorous and preclude the use of commercially available equipment. This section presents an analysis of the impact of a typical NASA-developed specification on the selected equipment, and examines several alternative requirements which might reduce the cost impact of these specifications.

Section 4.2 presents the comparison of the NASA-developed ground equipment specification, EC006M00000A - "Part I - Performance, Design and Verification Requirements for Experiment and General Laboratory Equipment to be Selected for Project Sortie Lab," with the interface document developed during this study. (These specifications are presented in Appendixes D and E of this volume.) Considerable equipment modification is required by the NASA specification in addition to those required by the SEEIR. Additional engineering activities are also generated by this more rigorous specification. Delta costs for these additional modifications to the selected modified hardware items are presented to quantify the impact of the NASA specification in dollars.

Section 4.3 presents the results of tradeoffs of alternative approaches to reducing equipment costs. Three tradeoffs of alternative requirements for the equipment specification are discussed. An analysis was performed to determine if improved vibration isolation of the equipment would reduce the available hardware modifications necessary to survive the boost vibration environment. Reduction of material control constraints is the issue evaluated in the second tradeoff. Alternatives to degradation-type requirements such as fungus and corrosion are addressed in a third tradeoff. A fourth tradeoff, unrelated to equipment specifications, was also accomplished. This analysis addresses the possibility of reducing the power demand of available equipment by redesigning the power supplies internal to the unit and, instead, using Spacelab-conditioned power directly in each unit.

4.2 NASA SPECIFICATION VERSUS INTERFACE SPECIFICATION

The cost delta of imposing the preliminary NASA Specification EC006M00000A over the minimum requirement approach of SEEIR was determined. The result indicates EC006M and/or the SEEIR changes that might lead to better overall cost effectiveness in equipment modification. The two specifications are broadly compared in Table 4-1 in order to provide insight into the origin and philosophy used in the development of these documents. In all cases, the EC006M specification requirements are equal to or in excess of the SEEIR. The SEEIR is intended as a minimum requirement baseline and has not been approved for use outside of this study.

Table 4-1. NASA and Study Specification Origins

	SPACELAB/EXPERIMENT EQUIPMENT INTERFACE REQUIREMENTS	EC006M00000A GENERAL EQUIPMENT SPEC. FOR EXPERIMENT EQUIPMENT
<u>GENERATION</u>		
BY:	ROCKWELL	NASA
FOR:	THIS STUDY	MANNED SORTIE MISSIONS
WHEN:	4-74	11-72 (PRELIMINARY)
<u>PHILOSOPHY</u>	<ul style="list-style-type: none"> • MAXIMIZE EXPERIMENTOR RESPONSIBILITY • MAXIMIZE COMMERCIAL EQUIPMENT POTENTIAL • OPERATIONAL SPACELAB ENVIRONMENT, NOMINAL GROUND ENVIRONMENT 	<ul style="list-style-type: none"> • ASSURE EXPERIMENTER SUCCESS • COMPATIBLE WITH SPACELAB SUBSYSTEM REQUIREMENTS • SURVIVE VACUUM, WORST CASE TRANSPORTATION AND GROUND ENVIRONMENT

Figure 4-1 illustrates the logic process in determining the differences between the EC specification and the SEEIR. Initially the EC specification was evaluated to separate those requirements which were written as firm (uses the verb "shall"), and those that were not firm but advisory in nature ("shall consider", "should"). The objective of this activity is to reduce the probability of over-enforcing the NASA specification. For example, paragraph 3.2.6.2 of the EC specification states that "The equipment shall contain no ordnance devices," which is interpreted to mean that any item with ordnance devices (squib valves, etc.) would have to be replaced with alternate parts or designs. Another example is in paragraph 3.2.4.3 which states, "Accessibility




Figure 4-1. Specification Comparison - Approach



and ease of operation of latches, lockdowns, fasteners, etc., shall be considered for ease of maintenance." This is interpreted as an advisory requirement because implementation of the requirement is judgmental and generally is only accomplished if there is little or no cost impact.

Advisory, or insignificant cost impact requirements were promptly dropped from further consideration. The remaining significant cost impact requirements have a number of TBD values or reported data made obsolete by recent ESRO documents. These requirements were made consistent with the requirements defined in the Spacelab RFP so that the reported cost delta did not reflect specification differences resulting only from the date of issue of the specification.

Next, the actual comparison of specification requirements was performed. Those NASA requirements causing additional modifications or engineering activities were identified. The costs of these modifications were then determined indicating a delta cost to the previously defined modification costs.

4.2.1 Areas of Disagreement

The significant delta cost impact areas of the EC specification requirements over the SEEIR are summarized in Table 4-2. Several other requirements potentially qualify for the table, but applicability to off-the-shelf experiment equipment is negligible and would not constitute significant cost deltas. The type of delta activity incurred by an EC specification requirement is indicated by an "X" in the appropriate column. The impacts identified are generated for all equipment. Specific items may have different impacts.

As the table shows, all the listed requirements require delta design evaluation efforts to define requirements, provide technical coordination, manage program implementation efforts and other aspects of project and systems engineering. These activities lead to additional customer review coordination, documentation, design reviews, and approvals in one form or another, to verify that requirements have been met.

4.2.1.1 Requirements Increasing Engineering Costs

Reliability, maintainability and material control activities usually have special program plans to identify specific requirements for specific items. The magnitude of the latter efforts are subject to the customer's wishes and, therefore, vary from program to program. An estimated middle-of-the-road stance was adopted, assuming that demonstration of compliance via documentation and/or testing was needed on a fairly austere basis.

Experimenters, interested in obtaining useful results, will generally demand satisfactory equipment operation. Therefore, imposed assurance-type efforts will result in insignificant delta redesign. The EC specification, however, implies program efforts to verify or demonstrate adequacy. Therefore, increased system engineering costs will be incurred by the NASA specification.

Table 4-2. Delta EC006M Requirements Over SEEIR

TOPIC	EC006M PARAGRAPH	ESSENCE OF REQUIREMENT	ACTIVITY INCURRING DELTA COST			
			EVALUATION/ SYST ENGR	DATA/ APPROVALS	TESTING* NO REDESIGN	REDESIGN, + MOD & QUAL **
RELIABILITY	3.2.3	DEFINE AND CONFORM TO A RELIABILITY OBJECTIVE	X	X	X	
MAINTAINABILITY	3.2.4	PROVIDE EASE OF MAINTENANCE	X	X	X	
MATERIALS APPROVALS	3.2.6.3	NASA APPROVAL ON ALL MATERIALS APPLICATIONS	X	X	X	
NON-OPERATING TEMPERATURE	3.2.7.1	NO DAMAGE OVER -40 F TO +167F (-40 C TO 75 C) RANGE	X		X	X
NON-OPERATING PRESSURE	3.2.7.1	NO DAMAGE IN VACUUM FOR UNDEFINED TIME	X		X	X
TEST CONNECTORS	3.3.5.1.2	EXTERNAL TEST POINT ACCESS	X			X
MOISTURE SEALING	3.3.5.1.3	SEAL ALL CONNECTORS AND WIRING JUNCTIONS	X			X
SAFETY INTERLOCKS	3.3.5.1.4	DEMONSTRATE ADEQUACY	X	X		
LOOSE ACCESSORIES	3.3.6.1.3	INCORPORATE LOOSE ITEMS INTO BASIC ITEM WHEREVER POSSIBLE	X			
PVC	3.3.6.2.1	PROHIBITED--NO EXCEPTIONS	X			X
FUNGUS	3.3.8	NO UNTREATED NUTRIENTS	X	X	X	
CORROSION	3.3.9	NO METAL COUPLES PER CHART (MIL-SPEC)	X	X	X	
REPLACEABILITY/ INTERCHANGEABILITY	3.3.12.1	DEMONSTRATE ASSURANCE CONTROLS FOR INTERCHANGEABILITY	X		X	
HUMAN PERFORMANCE	3.3.15	MEET 10M33222 REQUIREMENTS 1/4-TURN FASTENERS, FRICTION HINGES, PANEL INSTRUCTIONS	X			X
*INCLUDES DEMONSTRATION TESTS, PARTS QUALIFICATION, DELTA ACCEPTANCE AND QUALIFICATION TESTS FOR AN UNMODIFIED ITEM						
**INCLUDES ALL EFFORTS ASSOCIATED WITH THE MODIFICATION DELTA EXCEPT INITIAL EVALUATION FOR MODIFICATION REQUIREMENTS.						



4.2.1.2 Requirements Causing Additional Equipment Modification

Nonoperating Temperature. The EC document specifies a -40 F (-40 C) to +167 F (75 C) nonoperating range. It is assumed that this could include worst case transportation conditions such as air cargo transport and dock area storage. The SEEIR specifies a nonoperating temperature of +40 F (4.5 C) to +130 F (54.5 C). The SEEIR implication is that the experimenter will be responsible to protect and/or ship by a nondamaging means and that Spacelab environment is compatible with commercial uses. Most equipment items examined specify or imply a nonoperating temperature range capability. Some provide acceptance tests or burn-in over a specified range. In general, the upper EC006M specified limit is not significantly higher than most equipment capabilities. Many appear to be impacted by the -40 F (-40 C) lower limit. It is estimated that 5 percent of an item's constituent parts are not capable of surviving the EC006M conditions. Additional thermal test chamber equipment and time will also be required in order to qualify and accept the redesigned item.

Nonoperating Vacuum. The EC specification requires a 0 to 15 psia range capability. The SEEIR requires 14.7 ± 0.25 psia. The slight delta overpressure by EC006M is negligible. The vacuum requirement could seriously impact parts; sealed components may explode. In addition, volatile lubricants may boil away during vacuum, making a unit inoperative. Significant additional test time is required to verify that the units can survive the vacuum environment.

For this study, 5 percent of the electronic parts of an item were assumed to be impacted by the vacuum environment. In addition, motors and other items requiring lubrication were identified so that replacement of lubricant with an acceptable type could be costed.

Test Connectors. The EC specification requires test equipment to mate with the item via test connectors. The worst case requirement would be for troubleshooting to logical replaceable assemblies. The SEEIR assumes that the usual method of troubleshooting or verifying performance is acceptable and does not specify a requirement.

In order to implement the requirement, it is assumed that one or more test connectors, a test harness, a test signal isolation circuit board and patch wiring on PCB modules through their connectors are required. Test connectors are 25 or 50-pin shell type as needed. More than one may be required. The number of test signals was estimated for each unit, ranging from 2 to 5 per replaceable assembly, depending upon item complexity and modularity. Spare PCB connector pins were assumed to be available for patch-wiring some signals from the boards.

To prevent loading, accidental short-circuit or other damage to sensitive signal sources, circuit isolation components must be installed. These may be installed on the PC boards. However, for costing it was assumed they would all go on a central special PCB to be fabricated. It was estimated that 25 percent of the signals would require isolation and that of these an average of 1 to 2 discrete electronic components would be required per signal.

It should be noted that it was assumed that the test requirement does not apply to replaceable assemblies that are tested out of the item. The EC specification does not make this clear, however.

Installation of test connectors largely meets the EC specification maintainability requirement. However, some items not modularly designed with plugable PCB's could require further redesign for ease of replacement. Redesign for ease of maintenance was identified for some units.

Moisture Sealing. The EC specification requires all connectors and wire junctions to be sealed. The SEEIR has no equivalent requirement. The rationale for the relaxation of this requirement is that normal practice should be sufficient for the currently planned Spacelab atmosphere.

Since the SEEIR requires encapsulation of PCB's in most cases for vibration and flammability reasons, it was assumed that all wiring junctions on the PCB's would be sealed, therefore, no delta would be required. Other junctions, besides connectors, are primarily at terminal strips. Therefore, the delta was assumed to consist of sealing or potting all connectors and terminal strips in each item. Some connectors are easier to seal than others. Printed-pin PCB connectors and some "auto harness" type molded mylar connectors will be very difficult to seal. The cost analysis assumes some percentage of connector replacements in addition to sealing. However, some PCB connectors may need re-sealing after each module removal/replacement. Two halves of internal connectors are sealed, while only one half of external interface connectors are sealed. Coaxial cable connectors are sealed as well as high-density pin connectors.

Polyvinyl Chloride (PVC) Replacement. The EC specification prohibits PVC without exception. The SEEIR allowed PVC in totally enclosed containers and on single loose wires. Costing delta consists of delta replacements of wiring (loose and harnesses) in sealed and totally enclosed units, unless it is TFE.

Human Performance Engineering. The EC specification requires NASA Specification 10M33322, "Man/Systems Design Requirements for Sortie Lab," to be met. Three requirements were judged to cause potentially significant cost increases:

1. All fasteners to be 1/4-turn self-contained.
2. All doors to have friction hinges.
3. Operating instructions to be silk-screened on an item's panel wherever possible.

It was assumed that each item would need at least four fasteners for rack securing. Also, any plug-in modules from the panel face would need such fasteners to replace existing thumbscrews, etc. The latter requirement affects NIM modules as well as commercial oscilloscopes, RF analyzers, etc.

All item doors were assumed to have hinge replacements.

All items were assumed to have some silk-screened instructions placed on panels except oscilloscopes and voltmeters. These two units were exempted for two reasons: (1) the panels on these are too crowded, and (2) operation of such familiar devices should not require special instructions. Instructions are assumed to be proportional to the number of panel controls. Some controls have little space available but effort to "design" the right words for this problem will offset the lack of available space.

4.2.2 Impact of EC006M00000A on Selected Equipment

This section provides a detailed summary of the delta modifications resulting from the EC specification as compared to the SEEIR. The basic modification information is discussed in Section 4.2.2.1, while detailed costs are summarized in Section 4.2.2.2.

4.2.2.1 Modifications

Table 4-3 lists the EC006M delta modifications for 31 items. The blood cell counter, freezer, and dewar were omitted from the list since these items could only be custom built. The modifications were determined by inspection and data similar to the manner in which they were previously defined. The following is an explanation of the data contained in Table 4-3.

Vacuum and Thermal Impacts. The nonoperating temperature and non-operating vacuum columns estimate the new parts and test facilities, including test time, needed to satisfy the NASA specification. As a general rule, one half of the parts replaced were assumed to be different. In the case of the Neff amplifier, a fewer number of part types is needed because the assembly consists of 16 identical modules.

As a general rule, 5 percent of the total parts count is estimated to need replacement to satisfy the additional thermal requirements if specifications indicate qualification below 0 F (-18 C). If the current manufacturing practice essentially meets the EC006M requirement, then no parts replacement is required. If the current practice lies somewhere between about -30 F (-34.4 C) and 0 F (-18 C), then 2.5 percent of the parts are assumed to need replacement. The laser was felt to be a special case where internal temperature control is needed to avoid low temperatures that could affect critical optics alignments. A heater and thermal switch were added to the unit.

A similar approach was used to estimate replaced parts for vacuum exposure. However, it was assumed that all items with significant populations of electronic parts would require the same replacement rate. That is, no reduction was assumed for parts qualified for 55,000 feet (1.68 km). The RHG transmitter was an exception because available information indicates a minimum usage of conventional discrete electronic packaging for active devices. The stripline and hybrid circuit type construction in the transmitter indicates a reduced potential for the problems often caused by vacuum; therefore, no replacements were required.

Table 4-3. EC006M Delta Modification Summary

ITEM	NON-OPERATING TEMPERATURE			NON-OPERATING VACUUM				①	ASSEMBLIES	CONNECTOR AND JUNCTION SEALING						PVC REPLACEMENTS		HUMAN FACTORS			
	PARTS		TESTING	REPLACE	REVISE SPECS	LUBRICATION	TESTING			INTERNAL			EXTERNAL			HARNESSES WIRES	LOOSE WIRES	CONTROL PANEL INSTRUCTIONS	CAPTIVE FASTENERS	FRICTION HINGES	MAINTAINABILITY REDESIGN
	REPLACE	REVISE SPECS								STANDARD	PRINTED CIRCUIT BOARD	TERMINAL BOARD	STANDARD	PRINTED CIRCUIT BOARD	TERMINAL BOARD						
KEYBOARD TERMINAL	25	*	NEW	25	*		NEW	25	5				1			50	50			X	
SPECTROPHOTOMETER	50	*	NEW	50	*	2	NEW	50	16	10	8		3			50	10		X	X	
TAPE RECORDER (AMPEX)				100	*	4	NEW	50	50		66		4			50	10	6			
TAPE RECORDER (HONEYWELL)				100	*	6	NEW	50	30	40	25	4	3			50	10	6			
COMPUTER	75	*	NEW	75	*	2	NEW	50	20	5	60		6			100	22	4		X	
MULTICHANNEL ANAL (ND 100)	100	*	NEW	100	*	2	NEW	50	13	10	7		8			100	32	4			
TIME CODE GEN (DATATRON)	20	*	DELTA	40	*		NEW	25	5		4		10			25	10	4			
POWER SUPPLY, NIM (P.D.)	7	*	NEW	8	*		NEW		1				1			50		4			
LASER	②	②		10	*		NEW	15	3			1	1			25	15	2			
VOLTMETER - 8200A				60	*		NEW	50	15		12		2	3			50				
VOLTMETER - 8125A				35	*		NEW		1								10			X	
TIMER COUNTER (TENNELEC)	25	*	NEW	30	*		NEW	25	3		4		10				10	12			
RF SPECTRUM ANAL (SINGER)				30	*		NEW	50	8	7	15		6			105	100	39	8		
EMI METER (SINGER)				100	*		NEW	100	30	32	22		7				13	4			
FREQ SYNTHESIZER - 645M	100	*	NEW	100	*		NEW	100	30				16			60	100	18	4		
MICROSCOPE ③						1		5													
CENTRIFUGE								2										1			
POWER SUPPLY (SORENSEN)				25	*		NEW	25	4	4					2		200	7	4		
GAS CHROMATOGRAPH	90	*	DELTA	90	*	1	NEW	15A	3A	5P	7P		1P				15	4	XA	XA	
ELECTROPHORESIS SYSTEM ④	60	*	NEW	60	*	2	NEW	50	13	4	10						200	10			
AMPLIFIER	100	9	NEW	100	9		NEW	75	17		17		2				64 ③	4		X	
CRT STRIP RECORDER	50	*	NEW	50	*	4	NEW	150	45	3	36		36				50	31	4		
VHF TRANSCEIVER (COLLINS)			DELTA	100	*		NEW	50	11	6			2			100	50	6			
TRANSMITTER (RHG)						1	NEW	50	10			3					50	9	4		
NIM COUNTER ASSY (ORTEC)	40	*	NEW	40	*		NEW	25	5	5			5				50	10	16		
RF SPECTRUM ANAL (H-P)				160	*	2	NEW	50	12	8	20	1	5				175	22	6		
FURNACE								2										4		⑤	
OSCILLOSCOPE (TEKTRONIX)				100	*	1	NEW	100	15	20			9				50	21		X	
OSCILLOSCOPE - MIL (H-P)				75	*		NEW	50	5	5							100				
MICROTOME						2		2										8		X	
pH METER	5	*	NEW	5	*		NEW	15	3		2		3				50	6		X	

① ESTIMATED SIGNALS NEEDED TO EXTERNALLY IDENTIFY FAILED REPLACEABLE ASSEMBLY (IN CONJUNCTION WITH AVAILABLE INPUT-OUTPUT AND DESIGN CHARACTERISTICS).

② ADD HEATER AND THERMAL SWITCH INSTEAD OF PARTS REPLACEMENT.

③ INCORPORATE LOOSE OPTICAL FILTERS AND/OR WHEELS INTO DESIGN.

④ INCORPORATE CELL AND POWER SUPPLY INTO UNIT DESIGN.

⑤ GROUP OF FOUR, REPEATED EACH AMPLIFIER MODULE.

⑥ INCORPORATE FURNACE ACCESS COVER WITH FRICTION HINGES.

* ASSUME ONE NEW SPEC FOR EACH TWO PARTS REPLACED.

A = ANALYZER ONLY

P = PROGRAMMER ONLY

- ① ESTIMATED SIGNALS NEEDED TO EXTERNALLY IDENTIFY FAILED REPLACABLE ASSEMBLY (IN CONJUNCTION WITH AVAILABLE INPUT-OUTPUT AND DESIGN CHARACTERISTICS).
- ② ADD HEATER AND THERMAL SWITCH INSTEAD OF PARTS REPLACEMENT.
- ③ INCORPORATE LOOSE OPTICAL FILTERS AND/OR WHEELS INTO DESIGN.
- ④ INCORPORATE CELL AND POWER SUPPLY INTO UNIT DESIGN.
- ⑤ GROUP OF FOUR, REPEATED EACH AMPLIFIER MODULE.
- ⑥ INCORPORATE FURNACE ACCESS COVER WITH FRICTION HINGES.
- * ASSUME ONE NEW SPEC FOR EACH TWO PARTS REPLACED.
- A = ANALYZER ONLY
- P = PROGRAMMER ONLY

Lubricants may require replacement if exposed to vacuum, depending upon sealing and duration of exposure. The items that contain significant lubricants (in terms of replacement costs) were identified by the number of lubricated assemblies such as motors, linkage assemblies, etc. The Honeywell tape recorder, for instance, contains three tape-handling motors, two fan motors, and pinch-rollers, a total of six lubricated assemblies. The CRT strip recorder has one motor, one tachometer and a paper-handling roller assembly, i.e., three lubricated assemblies. The microtome has numerous gears, ratchets, etc., making up two lubricated assemblies.

Test Connector Signals. Estimates of the number of test signals required for each item were based upon assessment of the number of logical assemblies that were functionally independent and could be repaired or replaced. As explained previously, wiring harness, connector, signal isolation circuitry and PCB rework requirements were estimated from this information.

Moisture Sealing. A count of internal and external connectors and terminal strips/boards was made to estimate sealing requirements. Connectors are divided into standard and printed-pin PCB type. These categories represent different sealing problems that could result in significantly different relative costs. Access to each connector for sealing is a function of the concerned item.

PVC Replacement. Estimates are provided for the wires in harnesses that are in totally enclosed units and therefore not replaced due to the SEEIR. No entry (-) indicates no PVC in the item (either no wiring, or TFE is reportedly used).

Human Factors. The number of controls per unit is shown since this indicates the effort required to define and apply silk-screened panel instructions. No entry (-) indicates no controls (NIM power supply) or the requirement for instructions is assumed to not be possible due to panel control density.

The number of captive 1/4-turn fasteners for each item is based on the assumption that rack-mounted units attach by their front panel and require at least four per item. The tape recorders were allotted two additional fasteners to provide more support for vertical rack mounting. Where panel access to replaceable modules are used, the mounting screws (usually thumb-screws) are replaced with the required fastener. Where new rack-mounting brackets were required by SEEIR and logically could include the fasteners with negligible added cost, a fastener delta for the EC specification was not estimated.

Finally, all doors or lids with hinges were identified and replaced to meet the EC006M and 10M33222 requirement for friction hinges.

4.2.2.2 Delta Modification Cost Analysis

The additional costs incurred by requiring the selected equipment items to conform to the EC006 specification are of two types. The first type includes additional modifications identified to satisfy the requirements of this specification which are in excess of the study-developed specification. A second cost type, particular to the EC006 specification, includes additional engineering requirements which were not associated with a modification. These costs reflect the additional analytical, testing, and reporting efforts typical of most aerospace-type specifications.

Approach. The approach to defining each cost element is similar to the approach explained for the determination of modification costs in section 3.0. The delta activities were categorized as follows: vacuum and thermal capability, test connectors, sealants, wiring, fungus, corrosion, material control, safety interlocks, interchangeability/replacement, maintainability, reliability, and parts replacement. Wherever possible a count of hardware was estimated and used as a basis for estimating, at a grass roots level, the number of engineering and manufacturing hours required to perform the required tasks.

The engineering tasks include parts identification, requirements documentation, vendor coordination, test procedures, wire definition, test signal selection, connector and harness design, circuitry and voltage analysis, materials analysis for fungus control, elimination of electrolytic couples which might corrode, maintainability design, documentation and customer approval, reliability predictions and change identification. A major input is the systems engineering effort to make these identifications and ensure consistency with the Spacelab requirements. Each of these tasks were reviewed individually by hardware item to ensure first that the effort was a requirement, and secondly, that no overlap or omission occurred. The hours were then extended by the appropriate rate which included labor, burden, and administrative expense to project the total cost.

To determine the manufacturing cost impact, man-hour standards were utilized to estimate the cost of modifications derived by engineering. These standards are the same as applied in the cost analysis of the first set of modifications and were applied to fabrication, material, quality assurance, and test procedure costs.

The test effort required was projected by engineering personnel and related to the changes required to meet the NASA specification. Caution was taken to avoid any duplication and where appropriate, such as in vibration tests, it was assumed that the test estimate for the original modification was adequate and would be performed at the conclusion of the modification program whether the study specification or EC006 specification was applied. In either case, this test is performed only once.

The same factors used for the original modification costs were applied to the documentation and project management costs; i.e., documentation equals 10 percent of engineering, manufacturing, test and program management equals 5 percent of total cost.

Results. Hardware modification costs are significantly impacted by an aerospace-type equipment specification. The average cost of hardware modified to the MSFC specification is twice the magnitude of the average cost incurred when equipment is modified to the study-developed specification. Table 4-4 summarizes the cost impact by equipment item and function. The basic cost categories are itemized including a breakdown of engineering costs into costs incurred by the modification of hardware and those incurred by systems engineering activities necessary to satisfy the EC006 specification. The total cost of additional modifications are shown for each unit. Adding this cost to the cost of the modified hardware described in Section 3 establishes the total cost of each item when modified to meet the EC006 specification.

Fabrication costs are lower relative to modification engineering costs than generally estimated for the equipment modified to SEEIR requirements. This lower multiple results from inclusion of assembly, disassembly, certification, and general manufacturing support costs with the previous modifications. The fabrication costs shown in Table 4-4 reflect only the specific manufacturing activities related to the delta modification required.

The test costs did not vary as much between equipment items as it did in the previous modification activity. Testing in support of materials control requirements of the EC006 specification and a maintainability demonstration cause this smoothing of test costs. The cost of testing materials to identify their toxicity and flammability characteristics were apportioned equally across all equipment items, since this information is applicable to more than one equipment type. (Information regarding the characteristics of a specific material would apply to all units containing that material.) Cost of a demonstration of ease of maintenance was also assessed against each unit. The cost of these tests was much higher than the cost of all other delta testing associated with the EC006 specification. Therefore, major differences in test costs did not result.

The cost of the additional modifications and system engineering amounted to \$936K for the 31 units modified. This amount is nearly the same as the cost of modifying the hardware to meet the SEEIR. The sum of these two costs equals the cost of modified hardware satisfying the EC006 specification. As Table 4-4 shows, a major fraction of the cost of meeting the EC006 was incurred by the systems engineering activities necessitated by the specification. The SEEIR minimizes reporting activities, eliminating the need for much of these costs.

Figure 4-2 compares the cost of modifying each of the selected equipment items to the EC006 specification or the study-developed equipment specification (SEEIR). The shaded area in each bar represents the cost of modification to the SEEIR. The total height of each bar is the cost of modifying the equipment to the EC006 specification. The blank area then represents the cost added by the EC006 specification.

Modified hardware costs doubled when the requirements of the MSFC specification are met. The average cost of the selected hardware is \$62,000 as compared to \$32,000 for the average cost of the hardware modified to the SEEIR. The average of the ratios of modification cost to retail cost rises from 8 to 18 when the NASA specification is imposed. The ratio of the total modification cost to the total retail cost rises from 4 to 10.

Figure 4-3 compares the hardware modified to the EC006 specification with estimated cost of the same hardware custom built. This figure should be

Table 4-4. Delta Modification Cost Summary*

UNIT	MANUFACTURER	FABRI- CATION	MODIFI- CATION ENGR	SYSTEM ENGR	TEST	DOCUMEN- TATION	PROGRAM MGMT	TOTAL MODIFICATION COSTS		
								DELTA FOR EC006	TO MEET SEEIR	TO MEET EC006
TAPE RECORDER	AMPEX	4,840	13,760	15,919	5,594	2,562	2,133	44,808	58,222	103,030
COMPUTER	HONEYWELL	1,763	13,204	15,040	5,502		1,903	39,974	49,531	89,505
	DEC	2,233	12,480	13,475	5,649		1,820	38,219	62,293	100,512
	DATATRON	690	5,620	9,008	5,262		1,157	24,299	21,598	45,897
SPECTROPHOTOMETER	TENNELEC	725	3,016	8,171	5,382		995	20,869	14,213	35,082
	BECKMAN	999	8,690	16,757	5,465		1,725	36,197	51,385	87,582
	SINGER	1,657	13,928	13,929	5,410		1,874	39,360	37,767	77,124
EMI/FIELD INTENSITY MTR										
TRANSMITTER	RHG	424	4,985	12,679	4,876		1,276	26,802	17,933	44,735
TRANSCEIVER	COLLINS	907	9,675	12,791	5,318		1,563	32,816	15,177	47,993
AMPLIFIER	NEFF	1,421	10,620	12,104	5,686		1,620	34,013	22,736	56,749
SPECTRUM ANALYZER	H/P	1,727	16,080	19,450	5,851		2,284	47,954	48,099	96,053
	SINGER	1,255	8,600	11,125	5,483		1,451	30,474	30,402	60,876
	NUCLEAR DATA	1,447	15,000	16,206	5,741		2,048	43,024	37,288	80,312
MULTICHANNEL ANALYZER										
FREQUENCY SYNTHESIZER	FLUKE	1,524	18,120	16,361	5,851		2,209	46,627	37,352	83,979
LASER	SYLVANIA	258	2,720	7,749	4,637		896	18,822	35,015	53,837
GAS CHROMATOGRAPH	BECKMAN	971	10,430	13,038	5,391		1,620	34,022	44,968	78,990
VOLT/OHM METER 8200 8125A	FLUKE	872	4,965	10,556	5,170		1,206	25,331	20,889	46,220
	FLUKE	133	2,776	8,687	4,582		937	19,677	15,156	34,833
OSCILLOSCOPE	TEKTRONIX	1,143	13,270	25,912	5,362		1,910	40,115	35,243	75,358
	H/P	698	8,555	13,157	5,005		1,499	31,476	25,407	56,883
ELECTROPHORESIS	BECKMAN	1,240	10,705	12,939	5,778		1,626	34,150	114,175	148,325
PARTICLE COUNTER	ORTEC	854	9,160	12,908	5,341		1,494	31,369	18,778	50,147
pH METER	BECKMAN	396	2,866	7,960	5,170		898	18,582	21,431	40,013
KEYBOARD DISPLAY	RESEARCH	614	5,190	14,848	5,281		1,175	24,665	35,542	60,207
CENTRIFUGE	CLAY-ADAMS	-	-	7,544	4,420		800	15,362	9,365	24,727
MICROSCOPE	AMER OPTICAL	-	-	6,164	4,420		707	13,853	19,186	33,039
FURNACE	ASTRO	-	-	11,500	4,420		1,066	19,548	25,545	45,093
MICROTOME	AMER. OPTICAL	-	-	6,311	4,420		719	14,012	11,697	25,709
POWER SUPPLY	SORENSEN	719	4,140	8,372	5,373		1,058	22,224	20,015	42,239
	POWER DESIGN	318	3,840	8,801	4,766		1,014	21,284	14,680	35,964
STRIP CHART RECORDER	HONEYWELL	2,191	19,130	14,229	5,823	2,562	2,197	46,142	32,725	78,867
TOTAL		32,079	251,525	380,685	162,479	79,422	44,880	936,070	1,003,813	1,939,883
*COSTS ARE GIVEN IN DOLLARS										




Figure 4-2. EC006M00000A Versus SEEIR Modification Costs




Figure 4-3. Modified Versus Custom Costs



compared to Figure 3-36. Previously, all 31 items modified were estimated to cost less than the cost of equivalent custom-built hardware. Enforcement of the EC006 specification caused five of these items to have modification costs in excess of their custom costs. The Tennelec timer/counter, the Power Designs power supply, the Neff amplifier, the Fluke Model 8200 volt-ohm meter, and the Sorensen power supply can be more cost effectively procured if custom built in this case. The excess of the modification cost above the custom-built cost is shown in black on the bars of these units. The average savings on all 31 modified units dropped from 5195K to \$163K. Therefore, application of either specification results in the same conclusion--use of commercially available equipment and instrumentation is cost effective when compared to custom-built hardware.

The average per-unit savings is still large even though the EC006 specification is applied. However, a potential savings of \$25 million could result to the program if the SEEIR specification were used rather than the EC006 specification. Review of the specification comparison presented in Section 4.1 indicates that the requirements of the EC006 specification could be relaxed without affecting the performance, reliability or safety of the experiment equipment. Therefore, adopting equipment requirements consistent with those defined in the study-developed specification is recommended.

4.3 ALTERNATIVE SPECIFICATION REQUIREMENTS

The study thus far has examined the requirements and costs of modifying CAM equipment to meet specified requirements for placement in the Spacelab. An alternative to modifying equipment is relaxing these specified requirements allowing available equipment to be placed on board with as little modification as possible. This section explores possible changes to the Spacelab equipment specification. Candidate requirements for tradeoff were identified by examining the characteristics of the CAM equipment in the available configuration, the modifications identified in Section 3.0 to make the equipment suitable for installation in the Spacelab, and the differing requirements identified by the specification comparison discussed in Section 4.2.

Results are presented of tradeoffs of alternative specifications for the vibration environment, material control requirements, and corrosion and fungus control, along with an analysis of the impact of changing equipment to be compatible with the 28 vdc primary Spacelab power source. Detailed tradeoff information is presented in Volume III.

4.3.1 CAM Equipment Rack Vibration Isolation

4.3.1.1 Objective

The suitability analysis determined that almost all items would require modifications to meet the ERNO Spacelab vibration spectrum used in the SEEIR specification. Sixty to 80 percent of the total CAM modification costs apply to vibration modifications. While the evaluations show that the modifications are highly cost effective, compared to the cost of custom equipment, these modifications could be significantly reduced if the items were more effectively isolated from the Spacelab vibration. Therefore, it is appropriate to



investigate the feasibility of vibration isolation and to evaluate the costs of several alternative methods.

The tradeoff is directed primarily at standard rack-mounted equipment because 25 of the 34 items studied were suited for such mounting. However, study results could just as well apply to the five items that were found to be unsuited to standard rack mounting by assuming the use of a special rack or bench with similar isolation.

4.3.1.2 Comparison of Alternatives

There are many possible system solutions to isolate the equipment items. Isolation can either be active or passive. Active isolation is not cost-effective for the smaller equipment clusters. It may be potentially effective for a number of consoles mounted in one rigid system. It may also be cost-effective at the total Spacelab level since costs would be shared by all equipment items.

Low-frequency passive isolators are feasible at the unit, rack, or group-of-racks levels. They would be less desirable at the total Spacelab level because of the resultant interface complexities from large peak displacements. The assessments made were therefore aimed at the potential of (1) passive isolation of equipment units in hard-mount racks, (2) hard-mounting of units in passive isolated racks, and (3) hard-mounting of units in racks that are hard-mounted in groups to a passively isolated framework.

Isolated Equipment Units in Hard-Mounted Racks. The most straightforward isolation approach would be to isolate each item in an equipment rack. This alternative is not cost-effective and impacts the physical dimensions of the rack. Development of an isolator effective in the low-frequency range necessary for the equipment is estimated to cost \$150 thousand. Since isolators are custom built to the characteristics of the package attenuated, this cost would be incurred for each package with dissimilar mass or resonance characteristics. The cost of the internal modifications of each package averaged \$10,000 and is therefore more cost-effective. The excursions of the vibrating equipment could be up to ± 1.5 in. (3.8 cm). These excursions would necessitate the construction of an over-sized rack, impacting the Spacelab design. Inter-unit connections would also be impacted by these excursions.

Rack and Multi-Rack Isolation. Distribution of the costs of isolator development among more equipment items by either isolation of a rack or a group of racks is more cost-effective. The most serious potential problem is that total weight variations in each rack may require custom isolators. The flexibility of using Spacelab would be affected. Either the isolators could not be designed until the rack was loaded, or extensive coordination would be needed to assign equipment items to racks and establish weight goals. Neither is compatible with a flexible user policy. This problem would be solved if a standard range of isolator designs were made available to accommodate various weight ranges. This issue requires an in-depth study to verify design feasibility and to establish the required design parameters.

4.3.1.3 Conclusions

The following conclusions were made during the tradeoff study:

1. Grouping of equipment items into multiple-rack isolated systems would be cost-effective provided that the proper isolators are made available for the spectrum of weights and mass distributions. Isolation of individual equipment items within racks is, as a general rule, impractical and not cost-effective.
2. A test program should be undertaken to understand the true non-operating vibration environment capability of CAM equipment. Two thirds of those items studied were not tested to a defined vibration spectrum. Data for the other one-third quoted mostly sea-going vessel specifications for operating conditions and many may be capable of more stringent low-frequency inputs.

4.3.2 Material Control Tradeoff

4.3.2.1 Objective

The objectives of this trade study are to (1) determine whether it is feasible to reduce the Spacelab internal experiment equipment material control requirements from the level defined by the NASA equipment specification (EC006M00000A) to the study SEEIR specification, and (2) determine whether it is feasible to eliminate all material control requirements (except for mercury, cadmium and zinc).

The first objective was evaluated by a risk analysis which considered the cost differential between imposing the study SEEIR specification requirements and the NASA specification requirements upon the equipment, compared to the cost of modifying and refllying experiments which failed because of a material control type of problem.

The second objective was evaluated by comparing the cost savings expected by further relaxing the study specification material control requirements, with the costs of adding an active method of cleaning the Spacelab environment of increased outgassing products or isolating the increased outgassing products.

4.3.2.2 Discussion of Results

Relaxation of the material control requirements of the EC006 specification to the requirements of the SEEIR does not reduce crew safety and is cost-effective. The specifications differ only in the depth of material identification required and the treatment of enclosed units. The 150-hour bakeout required for all units by the SEEIR and conformal coating of printed circuit boards eliminates the need for identification of all materials. Enclosed packaging restricts co-mingling of the package atmosphere with the habitable environment. Therefore, both specifications are equally safe.

A risk analysis was used to evaluate the cost-effectiveness of relaxing the material control requirements down to the level of the study equipment specification. Relaxation of a specification increases the risk of equipment failure.

The results of the risk analysis are shown in Figure 4-4. The horizontal line represents the cost differential between the NASA and the study level of




Figure 4-4. Risk Analysis for Relaxing Material Control Requirements from NASA Spec to Study Spec

material control requirements. This line represents the additional cost of identifying all materials in the equipment, approving the listed materials, and replacing any wire runs (internal and sealed units) from PVC to Teflon insulated wire. Since the cost is not probabilistic, because it would be expended for all items of equipment, it is represented as a horizontal line. The slanted lines represent the expectation of incurring the cost of modifying up to five units that failed due to a material failure. Since such a failure is probabilistic, the expected cost varies with the probability of occurrence. For example, the probability of five equipment failures occurring during the program may be only 1 in 10. Therefore, the expectation is only about \$700 thousand. Another way of explaining this cost is that with a probability of 1 in 10, the total cost of five failures would only be incurred once in ten experiment programs. Therefore, one-tenth the cost of these five failures is apportioned to one experiment program.

The cost of a failure is based on the cost to refly the failed experiment hardware and the cost to fix that hardware prior to reflight. Ten experiments are assumed to fly on any Spacelab; therefore, one-tenth the cost of reflighting the Spacelab is assessed against that experiment. The cost per flight used was \$10.45 million (1974 dollars) which is an average cost for 439 flights.

The intersection of the slanted and horizontal line defines the probability of occurrence where the cost of imposing the EC006 specification equals the cost of relaxing the specification to the SEEIR requirements. If only one failure should occur, it is always cost effective to relax the specification. The probability of two failures occurring during the program has to exceed 0.72 to make imposition of the EC006 specification cost effective. The probability must be greater than 0.5 for three failures to occur to make the EC006 specification cost effective and so forth.

The differences between the two specifications are small and should not result in significantly different failure rates. The major cost difference results from different material reporting requirements. An increased risk exists only for enclosed units. The SEEIR allows materials to exist in these units which would not be allowed by the EC006 specification. In the opinion of the study team, the probability of a material failure causing an equipment failure is not significantly different for the two specifications. It is believed that the probability of two failures occurring as a result of the specification differences is much less than 0.72. The probabilities for additional failures are also much less than the cross-over points shown on the figure. Therefore, relaxation of the EC006 specification to material control requirements defined by the SEEIR specification is recommended.

Elimination of all material control requirements was also examined. Total isolation of all equipment, or addition of a catalytic burner and appropriate sorbents to the Spacelab environment control subsystem (ECS), were examined. Equipment could be placed in a sealed rack charged with an inert gas to prevent fire. Sealing of the rack prevents toxins evolved from the equipment from entering into the habitable environment. Scrubbing of the atmosphere to maintain trace contaminants in the atmosphere below allowable thresholds can be accomplished by modifying the Spacelab ECS. Fire extinguishers would also be required to control a fire in this latter case.

Both alternatives are not cost-effective and were rejected. The inert gas isolation approach would cost at least as much as instituting a materials control program and would only be effective for the rack-mounted equipment. Sealing of rack-mounted hardware would require equipment modifications as extensive as those required to meet the material control specifications identified in the Spacelab/Experiment Equipment Interface Requirements (SEEIR) document. Isolation of bench-mounted equipment requires special packaging and could affect equipment operations. Most likely, a materials control program would still be necessary for these packages. Institution of a materials control program for a portion of the equipment would not achieve the objective of further relaxation of material control requirements.

The contaminant removal subassembly and fire extinguisher approach also is not cost effective. Development of a catalytic burner and associated sorbent beds is estimated to cost \$5.7 million. Each unit installed in a Spacelab would cost an additional \$0.8 million. If three Spacelab modules are required for the Sortie Experiment program, a total cost of \$8.1 million would be incurred. The cost of fire extinguishers is negligible by comparison, and has been recommended by the Phase II Spacelab study. The \$8.1 million cost does not compare favorably with the estimated programmatic cost for modification of equipment to meet the requirements of the SEEIR specification, which is \$3.7 million.

4.3.2.3 Conclusion

Reduction of material control requirements to a level defined by SEEIR is cost-effective. A risk analysis showed that the costs incurred by a failure resulting from relaxed material controls are less than the cost of imposing the NASA equipment specification. No compromise to crew safety is envisioned by implementing the study-developed material control requirements.

Further specification relaxation is not recommended. The cost of removing the increased offgassing products by means of a catalytic burner or the isolation of the offgassing products by cooling all rack-mounted equipment with a closed, gaseous nitrogen system exceeded the cost of implementing the study-developed material control requirements. More important, crew safety could very well be jeopardized by the increased amount of offgassing products if special precautions were not exercised.

4.3.3 Moisture and Fungus and Corrosion Requirements Tradeoff

4.3.3.1 Objective

The NASA specification for CAM equipment currently requires all materials used to be non-nutrient to fungus growth, or treated so the exposed surfaces will be fungus resistant. Also, the specification requires all metals to be corrosion resistant or to be processed to resist corrosion. Both of these requirements are intended to protect ground signal equipment from the deteriorating effects of climatic and service conditions encountered in military use.

The objective of this trade study then is to determine if the requirements for moisture and fungus resistance, and corrosion resistance (as stated in NASA Specification EC006M00000A), can be waived without significantly reducing the life of the equipment.

4.3.3.2 Tradeoff Conclusions

It is recommended that the EC006M00000A requirements for moisture and fungus-resistant materials be waived. The basis of this recommendation is as follows:

1. Ground operations for Spacelab, being studied for the NASA under Contract NAS1-12933, implements a study ground rule requiring a 50-percent relative humidity (or less) environment for all Spacelab ground operations including experiments. Studies have shown that relative humidities below 50 percent will not promote fungus growth nor electrolytic couples.
2. Periodic flight operations at 70-percent relative humidity (maximum) are not of sufficient duration to promote fungus growth and/or corrosion.

4.3.4 Power Distribution for Experiment Equipment Study

4.3.4.1 Objective

The issue addressed by this tradeoff is whether it is more cost-effective to modify experiment equipment items to utilize primary 28 vdc directly than to convert this power to 115 vac (60 and 400 Hz) and 220 vac (50 Hz) to be compatible with the equipment as currently designed.

Conversion of the 28 vdc Spacelab primary power to various voltages and frequencies for input to these units (with subsequent reconversions to internal secondary voltages) represents a double conversion loss. Such conversion

losses directly reduce the amount of available useful power, increase the heat rejection requirement for a given experiment load, and increase the weight of the conversion hardware. However, 28 vdc primary power may also not be a convenient power form, thereby also incurring conversion losses. This analysis compares the main elements of both power conversion systems in order to determine the relative electrical efficiencies, weights and costs involved.

4.3.4.2 Discussion of Results

The two alternative concepts shown in Figure 4-5 were compared with respect to conversion efficiency, weight and cost. No clear advantage exists for the efficiency of either alternative. Both concepts actually operate as dc-dc converters. The main difference is that the ac approach uses an extra transformer and harmonics filter. Also, the conversion frequency for the typical dc-dc converter is much higher in order to use low-cost parts and decrease weight.

Typical ranges of power losses for the key components used in the two approaches are noted in Figure 4-5. Without detailed designs, a clear-cut efficiency advantage cannot be shown for either approach, but some advantage is indicated for the ac approach.

Dc-to-ac inverters in the kw range typically operate with a loss of 10 to 15 percent. An Apollo inverter unit had a loss of 15 percent, for instance. Smaller inverters and dc-dc converters exhibit greater losses due to the greater proportion of power required by chopper and regulation electronics. However, available space in a package, left by the removal of its power supply, will probably not be adequate to accommodate efficient transformers in the dc-dc converter. The transformers used in the dc-dc converter concept would probably be less efficient than the transformers in the ac concept.


The delta weight between the two approaches is comprised of differences in component weights and support hardware weights. Table 4-5 summarizes the weights of the main elements for each system. A significant weight saving is indicated with high-voltage ac power distribution as compared with low-voltage dc power distribution.

The weight of support systems (reactants and heat rejection), incurred by the delta inefficiencies between the two approaches, can be discounted since the efficiency evaluation showed little difference between the two systems.


The ac power distribution approach requires no power supply modifications in the experiment units, but will require development of Spacelab power conversion equipment to provide the standard voltages and frequencies required by the experiments.

The dc power distribution approach requires no Spacelab power conversion equipment for experiments (some will be required for Spacelab subsystem equipment), but will require extensive redesign and modification of almost all experiment units. Due to the large number of combinations of internal

EQUIVALENT CIRCUITS NORMALIZED TO 1-Ø SOURCE REPRESENTS MOST POWER USED BY CAM



AC CONCEPT



DC CONCEPT

() Efficiency Losses

Figure 4-5. Equivalent Circuits for Power Distribution Alternatives

Table 4-5. Delta Weight Summary - Pounds (Kilograms)

Based on 3 kw Distribution; 2% Line Volt Drops; 100 ft, 2-Wire Systems

Concept	Electronics	Transformers		Distribution		Total Weight
		Units	Distribution	Filters	Wiring	
60 Hz ac	5	60 (27)	92 (42)	9-28 (4-13)	16 (7)	182-201 (83-91)
400 Hz ac	5	60 (27)	41 (19)	4-13 (2-6)	16 (7)	126-135 (57-61)
DC	29 (14)	29 (13)	--	--	256 (116)	314 (143)

NOTE: Excludes primary regulators, structures for support/heat sinks, and unit rectifiers/filters/regulator/loads (assumed equivalent for both concepts).



voltages and characteristics, and the variation in space available in the units, custom power converter designs will be required for each unit. Experiment ac power requirements will result in dc-ac inverters that can affect volume and, in some cases, thermal design. Although the cost of the chopper, transformer, rectifiers and filters for each small dc-dc converter is less than a large dc-ac inverter, the frequency control, overall design complexity and development costs are comparable. It can be reasonably expected that the requirements for experiment package power supply modifications (one for each package), will result in many times the cost of the dc-ac inverters used in the Spacelab ac power distribution.

Table 4-6 summarizes the comparative costs of the two Spacelab power distribution approaches. Estimated costs are for the total program. As shown by the table, additional costs to install the modifications into the units for the dc concept were not estimated.

4.3.4.3 Conclusions

It was concluded that there is no easy way to eliminate the double conversions that take place in both approaches, short of complete redesign of the equipment for optimum use of source power. This would mean completely redesigned custom units with the attendant cost impact.

The ac power distribution system is the lowest cost approach to accommodate the wide range of CAM equipment studied. The cost of many new dc-dc converter designs, unit modifications, and attendant problems for the 28 vdc primary power distribution system far overshadows the cost of the standard inverters used in the ac distribution system.

Table 4-6. Delta Cost Comparison of AC and DC Power Distribution

CONVERSION UNIT		ESTIMATED DEV. COST EA. UNIT (\$K) *	DEVELOP- MENTS IN PROGRAM	ESTIMATED UNIT COST (\$K)	NUMBER OF UNITS IN PROGRAM *	HARDWARE COST (\$K)	PROGRAM TOTAL (\$K)
AC POWER DISTRIBUTION	3-PHASE, 115/200 VAC, 400 HZ INVERTER: 2.2 KVA	100-150	1	6-7	3	20	145
	115 VAC, 60 HZ INVERTER: 1.0 KVA	75-100	1	5-6	3	17	92
	220 VAC, 50 HZ INVERTER: 2.0 KVA	75-100	1	5-6	3	17	92
	28 VDC, REGULATED +2%, DC-DC CONVERTER: 500 WATTS	50-70	1	4-5	3	14	74
	TOTAL						403
DC POWER DISTRIBUTION	DC-DC CONVERSION MOD, 1-3 OUTPUTS	50-70	285	2-3	788	1,970	10,430
	DC-DC CONVERSION MOD, 4-7 OUTPUTS	50-70		5-7	788	4,728	13,188
	TOTAL						23,618 + EXPERIMENT EQUIPMENT MODIFICATION COSTS
*BASED ON 3 SPACELABS, 282 UNIQUE EXPERIMENT ITEMS (1577 TOTAL) PER PROGRAM							

4-24

SD 74-SA-0047-2

5.0 MODIFICATION SOURCE TRADEOFF

The purpose of this tradeoff is to determine the most cost-effective agency to perform equipment modifications. The objectives for the tradeoff are as follows:

1. Recommend source for equipment modifications.
2. Assess the sensitivity of the above recommendation to changes in the assumptions.

5.1 GENERAL APPROACH TO THE TRADEOFF

The following assumptions were used as the basis for conducting this tradeoff:

1. The range of the recommended equipment modifications for the selected equipment items (Section 3.0) is assumed to be representative of all equipment modifications planned for the Spacelab operations.
2. The frequency of Spacelab operations (flights) will be in accordance with the latest NASA Payload Mission Model (October, 1973).

The options considered in this trade include the following:

1. Each equipment manufacturer incorporates modifications to his own equipment.
2. One centralized contractor incorporates modifications to all equipment items.
3. NASA incorporates modifications to all equipment items.

The criteria used for assessing the various options included the following:

1. Cost to do equipment modifications.
2. Time to do equipment modifications.
3. Ability of each organization to incorporate the modifications.
4. Modified equipment performance risk.



5.2 EQUIPMENT MODIFICATION PROCESS FLOWS

In order to evaluate both the cost and time to perform equipment modifications, process flow charts were constructed in a modified PERT format. These flows were constructed to a level sufficient to highlight some of the differences that could occur among the various options. They also become the basis for identifying succinct modification tasks and allocating time spans for accomplishing each task.

Figure 5-1 represents the original equipment manufacturer's modification process flow. The modified PERT format allows time values (days) to be allocated for performing the major modification tasks. These time values are in nominal days, determined from discussions with equipment manufacturers, estimates based on the total number of tasks to be performed, and the general character of each task.

Three methods of accomplishing the equipment modifications by the original manufacturer were identified. Method 1 represents equipment taken from the assembly line at an appropriate stage of buildup, modified, and put back in the assembly line for final assembly and test. Method 2 represents the construction of modified equipment from the initial piece-parts to final assembly and test (separate and distinct from all other on-going manufacturing operations). Method 3 represents the modification of completely assembled equipment taken from the manufacturer's inventory. There is virtually no difference in process duration for any of the three methods presented (considering the restraint of the overall series of tasks to be performed). Methods 1 and 3 show a nominal overall time span of 103 days to accomplish equipment modifications, while Method 2 requires 104 days. However, Method 1 is probably the most cost-effective sequence since it maximizes the use of existing production procedures.

Figure 5-2 depicts the equipment modification process flow for a centralized contractor or the NASA. The total time to accomplish a representative equipment modification is seen to take considerably longer than the original manufacturer (approximately 75 percent longer). The increase in time results mainly from the ordering of the equipment to be modified, waiting for the equipment to arrive, the inspection of the equipment construction to determine the amount of modification required, and waiting for special parts and/or materials which were ordered. The relatively long time assigned to the last item for the NASA is a result of the larger number of items to be procured. It is assumed that the centralized contractor selected to perform equipment modifications will be a broad-range instrument manufacturer who would have more of this material on hand than would the NASA. Conversely, the time required to develop new drawings may be less for the centralized contractor and the NASA because of the learning curve effect brought about by repetition of similar modifications. This advantage, however, may be offset by a lack of familiarity with the equipment.

The actual time to perform the equipment modification (an average of two days) is the same in both process flows because the actual modification task should be equally unfamiliar for all three alternatives.





Figure 5-1. Original Equipment Manufacturers' Modification Process

FOLDOUT FRAME

ORIGINAL PAGE IS
OF POOR QUALITY



Equipment Manufacturers' Modification Process Flow

DO NOT FRAME 5

5-3, 5-4

SD 74-SA-0047-2

LEGEND:

- CONTRACTOR/NASA FLOW
- NASA FLOW VARIATION
- PATH OF NOMINAL CONSTRAINT (RECOMMENDED FLOW)
- TROUBLESHOOT PATH





Figure 5-2. Process Flows for Centralized

ORIGINAL PAGE IS
OF POOR QUALITY

FOLDOUT FRAME /



Process Flows for Centralized Contractor and NASA

5.3 FACILITIES AND OTHER RESOURCE REQUIREMENTS

From the equipment modifications recommended in Volume 3, the facilities and other resources required to accomplish the various equipment modifications were determined. A key question to be answered concerned the existence of these resources at each equipment manufacturer and the availability of these facilities to perform the modification tasks. The approach taken to answer this question was a questionnaire mailed to the 25 equipment manufacturers of the items selected for suitability analysis. In addition to a query about resources, the questionnaire also addressed the method of doing the modifications, the willingness of the company to participate in such a modification program, and an indication of minimum quantities of equipment required to initiate the modification activities. Figure 5-3 is a sample of the commercial equipment questionnaire.

The results of the survey are shown in Table 5-1 and summarized as follows:

1. Sixty-eight percent of the questionnaires were answered.
2. All manufacturers who replied were in general agreement that they could undertake the typical equipment modifications (they had the resources). Several cases were identified where changes to vendor-furnished equipment would have to be negotiated.
3. Fifty percent of the companies responding would pull units off the assembly line at an appropriate stage of buildup, modify them, and place the units back in the line for final assembly; 25 percent said they would modify the completely assembled units taken out of inventory; 12-1/2 percent stated that they would special-build the equipment; and 12-1/2 percent failed to answer this question or could not answer it until specific quantities of modified equipment were identified.
4. Eighty-eight percent of all companies responding said they would be willing to perform the equipment modifications; 12 percent said they might be interested.
5. The minimum number of units which would have to be ordered for a company to undertake equipment modifications ranged from 1 to 50--the average being 10 units.

The NASA was also questioned about the facilities and resources available to perform equipment modifications. It was determined that no central facility existed, but rather the capability was widely scattered among the various NASA centers. It was also determined that certain centers probably had excellent experience with specific types of equipment, depending on the kinds of development and operational activity inherent within each center.



COMMERCIAL EQUIPMENT QUESTIONNAIRE

The analysis of commercial equipment and instrumentation for possible use on Shuttle Spacelab payloads has shown that modifications are almost always required. The intent of this questionnaire is to determine the most cost-effective agency to perform these modifications.

1. The following items are typical modifications identified for a wide range of commercial products. Does your plant/division have the capability to modify your equipment in these representative areas?

	YES	NO	REMARKS
a. Replace PVC insulated wire with Teflon insulated wire.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> _____
b. Change circuit board and chassis installations to withstand Shuttle vibration environment by adding stiffening members, special holders, clamps, etc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> _____
c. Add or enlarge air passages for adequate forced-air ventilation aboard the Spacelab.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> _____
d. Replace/protect glass with Lexan.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> _____
e. Contain all shatterables (such as CRTs) within shield housings or equipment case.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> _____
f. Apply conformal coatings to circuit boards.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> _____

2. In general, how would you plan to modify your products?

☐ Pull off production line at an appropriate stage of buildup, modify, and put back in line for final assembly steps and test?

☐ Modify the completely assembled unit in a separate area?

☐ Other (Please specify)? _____

3. Would you be willing to modify some of your products for use on Shuttle Spacelab payloads?

YES ☐ NO ☐ MAYBE ☐

4. What minimum quantity of equipment would have to be ordered for you to undertake the aforementioned modification activities (please check one)?

1	2	4	8	12	18	25	50
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Company _____

Prepared by _____

Title _____

Telephone Number _____

Figure 5-3. Sample of Commercial Equipment Questionnaire



Table 5-1. Results of Equipment Manufacturer Survey

ITEM	REPLACE PVC INSULATED WIRE WITH TEFLON INSULATED WIRE?	CHANGE CB AND CHASSIS TO WITHSTAND VIBRATION?	ADD OR ENLARGE AIR PASSAGES?	REPLACE/PROTECT GLASS WITH LEXAN?	CONTAIN ALL SHATTERABLES?	APPLY CONFORMAL COATING TO CB?	PULL OFF PRODUCTION LINE, MODIFY, PUT BACK IN LINE?	MODIFY A COMPLETELY ASSEMBLED UNIT?	OTHER?	WILLING TO MODIFY?			MINIMUM ORDER QUANTITY
										YES	NO	MAYBE	
1	YES	YES	YES	NA	NA	YES	X			X			4
2	NO	YES	YES	NO	YES	YES	X			X			25
3	YES	YES	YES	YES	YES	YES	X			X			1
4	YES	YES	YES	YES	YES	YES			SB			X	1
5	YES	YES	YES	NA	NA	YES		X		X			
6	YES	YES	YES	YES	YES	YES	X			X			1
7	YES	YES	YES	YES	YES	YES	X			X			8
8	YES	YES	YES	YES	YES	YES		X		X			1
9	NA	NA	NA	NA	NA	NA							
10	YES	YES	NO	YES	YES	NO	X			X			
11	YES	VP	YES	YES	YES	YES	X			X			1
12	YES	YES	YES	NA	NA	NA			SD	X			1
13	YES	NA	NA	YES	YES	YES	?	?				X	18
14	VP	YES	YES	YES	YES	YES		X				X	25
15	YES	YES	YES	NA	NA	NA	X						
16	YES	YES	YES	NA	NA	YES		X		X			50
17	YES	YES	YES	YES	YES	YES	?	?		X			1
NA NOT APPLICABLE VP VENDOR PRODUCTS INVOLVED SB SPECIAL BUILD													

In order to more fully explore the possibility of performing equipment modifications right on the assembly line (as a part of the original equipment buildup), two representative commercial facilities were inspected--Beckman Instruments and Hewlett-Packard. In each case, assembly techniques were studied for different types of equipment.

At the Beckman Instruments facility, spectrometers were being individually assembled on a conveyor line. Assembly records traveled with each unit, so special assembly instructions could be accommodated. The most difficult part of the operation would be ensuring proper distribution of custom subassemblies or materials to the various work stations. With sufficient lead time, these procedures could be worked out with little or no impact on the normal assembly operations.

The manufacturing flow of a gas chromatograph was much simpler because these instruments were assembled at a single work station. For this kind of assembly operation, special parts, materials and subassemblies could be kitted for the modified instruments. Accordingly, it is anticipated that any number of custom units could be assembled without significantly affecting the normal production activity.

At the Hewlett-Packard facility, timers, digital counters and receivers were observed during their manufacturing flow. Here again, the single work station per assembled instrument approach was used. This same production technique appears to be quite effective for production rates in the range of 100 to 1200 units per month.

The conclusion reached after studying the manufacturing, assembly, and test operations of several different kinds of instrumentation equipment are twofold. First, it is entirely possible to incorporate modifications as a part of the original equipment buildup. Special parts and/or materials could be distributed to the various work stations along with special instructions. Second, the substitution of Teflon-insulated wire in place of PVC-insulated wire would probably cause the most problems at an original manufacturer's facility. The basic problem is the identification of the new wire so that it would not be mistaken for PVC-insulated wire commonly used in that plant. The current PVC-insulated wire uses many different color codings. The same color combinations could also be used for Teflon-insulated wire. So, rather than assign particular color combinations exclusively for Teflon-insulated wire, a more practical method would be to utilize wire which has the word "Teflon" stamped or printed on the insulation. For this identification method, a solid color (or white) would give the best results.

Should the Teflon-insulated wire inadvertently become mixed with PVC-insulated wire, there are three quick tests (other than the printed word, "Teflon") to determine the differences. These tests are as follows.

1. Touch. Teflon feels slippery, soapy, and is hard to hold between fingers when tension is applied to wire. PVC is easy to hold with fingers.

2. Adhesive. Most adhesives will stick to PVC, but will not stick to Teflon.
3. Flame. Teflon swells and sags but does not burn; PVC burns.

For each of these tests, any nylon sleeving on the PVC-insulated wire should be removed.

5.4 EVALUATION CRITERIA ANALYSIS

In evaluating the various options in terms of the criteria previously established (cost, schedule, ability, and risk), it should be stated that actual modification costs were not available from the original equipment manufacturer. The estimated modification costs that were generated in the study actually represent activities by a centralized contractor or the NASA because they include disassembly, reassembly, recalibration and retest operations. Although these kinds of operations could be performed by an original equipment manufacturer and included in his estimated costs, he would (in all probability) not incur these costs by choosing to modify the equipment during initial build-up (either on the assembly line, adjacent to the assembly line, or at a separate assembly station).

It should also be recognized that for each different category of equipment, the centralized contractor and the NASA must first become familiar with the construction details before specific modifications can be implemented. This learning curve effect is only important when the equipment acquisition rate is very low, in which case the original equipment manufacturer has a distinct time advantage. If the acquisition rate is significantly high, both the centralized contractor and the NASA would soon become familiar with the construction details of each specific item.

5.4.1 Modification Costs

The original equipment manufacturer can incorporate the recommended modifications for less cost than a centralized contractor or the NASA. The next least-cost option appears to be the centralized contractor. The rationale for this decision is based upon two prime considerations.

1. The first consideration is the basic familiarity of a centralized contractor with many types of instrumentation being considered for Spacelab operations. This organization would not take as long as the NASA, for instance, to become familiar with basic construction details, calibration methods, test procedures, and troubleshooting techniques (should they become necessary).
2. The second consideration is the inherent flexibility of personnel to do the various tasks and the high utilization of plant and equipment (low depreciation costs) of an instrument manufacturing company, as compared to the NASA. The company can place personnel performing modifications on other instrumentation fabrication activities,

while NASA does not have parallel fabrication activities to support the excess personnel normally required for peak equipment processing periods.

3. The total number of tasks to accomplish a typical modification is 14 for the original manufacturer and 16 for the centralized contractor and the NASA (see Figures 5-1 and 5-2).
4. The original manufacturer can perform modification cheaper because he does not necessarily incur double costs for disassembly, reassembly, recalibration and retest while a centralized modifier always incurs these costs. The original manufacturer can take a partially completed unit from the assembly line and modify it, while a centralized modifier always starts with a finished product. The following cost details for the PDP-8E Computer have been extracted from Volume III to illustrate this point.

Disassemble

40 man-hours, manufacturing engineer	
40 man-hours, manufacturing technician	
	\$1280

Reassemble

20 man-hours, manufacturing engineer	
60 man-hours, manufacturing technician	
40 dollars, material costs	
	\$1283

Retest to Original Specifications

100 man-hours, test	
40 man-hours, manufacturing engineer	
	\$2095

Other Support for Above Operations

15 man-hours, assembly	
20 man-hours, inspection	
20 man-hours, manufacturing engineer	
40 man-hours, expedite	
	\$1271

Total modification costs	\$5829
--------------------------	--------

The above example reveals that approximately \$5800 in additional costs would be incurred for each modified computer if the modifications were incorporated by a centralized contractor or the NASA.

Extending this example of program proportions and assuming the \$5800 applies to each of the modified units as an upper limit, the maximum program cost for this function could be as much as \$10 million. Also, unique test setups for individual equipment types and the acquisition of special test equipment would incur additional costs for a centralized contractor or the NASA.

Conversely, a centralized contractor or the NASA would, in all probability, be able to buy certain materials in sufficient quantity that considerable cost savings could be realized. This conclusion is based on the total quantities of acquired equipment (shown on bottom line of Table 5-1).

5.4.2 Modification Time

Original Equipment Manufacturer

The nominal length of time to modify a typical item of equipment, as shown in Figure 5-1, requires approximately 103 days. The path of nominal constraint (and minimum total time) follows either Method 1 or Method 3. Method 1 would modify equipment during the buildup cycle, while Method 3 would take completed items of equipment from the manufacturer's inventory for modification and subsequent testing. The other method (assemble modified equipment completely when all parts are available) does not show an appreciable increase in time to process equipment modifications. One conclusion reached is that the specific method is not a real-time constraint and should be left to the individual manufacturer's choice.

Contractor/NASA

The nominal length of time to modify a typical item of equipment is approximately 182 days (Figure 5-2). The path of nominal constraint (and minimum total time) follows the estimated times shown for the centralized contractor. It should be noted that some of this time (60 days) is spent simply waiting for the equipment to be delivered after an order has been placed. Also, 10 days (average) is required to thoroughly inspect the equipment for construction details and determine the extent of modification required. Once these modification details have been documented for an individual equipment item, the 10 days would not be repeated; but this time period would be required for each new item to be modified.

The variation anticipated for the NASA to do the equipment modification would add about 20 days to the total processing time. The extra time stems from the concept that the NASA is not currently a producer of commercial-type instrumentation equipment and would, therefore, be unfamiliar with the design details, some manufacturing processes, shop instructions, test procedures, material usage, etc., for specific types of instrumentation equipments. The learning curve effect has been discounted to this instance because of the many different types and makes of commercial equipment to be modified.

5.4.3 Modification Ability

Original Equipment Manufacturer

The results of the questionnaire revealed that all companies who responded were equipped with the necessary facilities and resources to do the recommended modifications. The only exception to this conclusion would be in the area of certain environmental test equipment. In those cases where an original equipment manufacturer does not have the necessary environmental test equipment, the tradeoff has assumed that the necessary verification tests would be subcontracted to environmental test laboratories for virtually no additional cost to the program.

Contractor/NASA

Because of the kind of business in which a centralized contractor and/or the NASA is engaged, and their size relative to a typical original equipment manufacturer, it was assumed that both parties would have the necessary facilities and resources. There are two exceptions: (1) personnel trained to work on specific types of equipment, special test equipment and test setups, and (2) a basic stock of certain parts/materials peculiar to an original equipment manufacturer constructing a specialized product line. These shortages can all be overcome in time, but they do represent additional tasks--much more so for the NASA than for a centralized contractor. The centralized contractor would have a distinct advantage because he already has a basic knowledge of certain instrumentation details, trained personnel, and available manufacturing facilities which he can utilize.

On the other hand, the fact that the NASA's facilities are widely scattered makes the acquisition of trained instrument personnel, special test equipment, special test setups, and increased material inventories a costly venture because of the inherent duplication involved. Should a centralized instrumentation facility be justified, the additional cost involved would be a real deterrent.

A factor against both types of centralized instrumentation facilities (Contractor/NASA) is the rate of equipment acquisition. Modified equipment totals developed in Section 4, Volume III have been plotted as shown in Figure 5-4. The significance of this figure is the fluctuating nature of the experiment acquisition totals per year, which reveals a peak every five years. The resulting fluctuation of manpower requirements would present a problem of periodic hiring and firing of personnel or shifting of personnel to perform some other useful tasks. Staffing flexibility available to the centralized instrumentation manufacturer might reduce the impact of these peaks somewhat, as previously discussed.




Figure 5-4. Rate of Modified Equipment Acquisition

5.4.4 Modified Equipment Performance Risk

Original Equipment Manufacturer

The original equipment manufacturer appears to be in the best position to warrant the performance of the modified equipment. This conclusion is based upon the typical practice that equipment performance is not warranted if modifications have been performed by the buyer/user. In almost all cases, the original equipment manufacturer will guarantee performance of a modified unit providing he has the opportunity to conduct verification tests on the modified equipment.

Contractor/NASA

Centralized contractors will not normally assume responsibility for equipment that they have not designed, manufactured and tested. Therefore, in the case of modified equipment, we must assume that a centralized contractor will retain the same attitude; i.e., no warranty stated nor implied. This attitude simply points out that the performance guarantee must be assumed by the modified equipment buyer/user. But, this additional responsibility would impose no undue burden, assuming that properly trained technicians and facilities are available to make the equipment work properly and to maintain the equipment. The risk identified here, in terms of equipment performance, must be borne by the buyer/user after the equipment warranty expires. So, in a very real sense, the risk involved is limited to the time period normally covered by the original equipment manufacturer's warranty--normally one year.

5.5 RECOMMENDED APPROACH

A review of Table 5-2 indicates that the original equipment manufacturer is the best choice to accomplish the recommended equipment modifications. This conclusion is based upon the factors of cost, the shortest time to accomplish the modifications, the inherent ability to incorporate the modifications, and the least risk to guaranteed equipment performance. Even without applying weighting factors to each evaluation criterion, Table 5-2 indicates a wide margin between the original equipment manufacturer's rating and the other options considered in this tradeoff. The first recommendation, therefore, is to utilize the original equipment manufacturer for equipment modifications in all cases where practical. It should be recognized, however, that a certain amount of technical coordination will be required with each manufacturer to determine the exact changes to be implemented for each type of equipment purchased.

The above recommendations must be tempered by the willingness of the original equipment manufacturer to make the required modifications. It is conceivable that some high-production companies will not be interested in modifying their product line for just a few items (the average minimum buy is about 10 units). If this is the case, the next best source to perform the equipment modifications would be a centralized contractor. This recommendation is based primarily upon the existing capability and lower cost of a centralized contractor to perform the required equipment modifications, as compared to the NASA.

Table 5-2. Equipment Modification Trade

	Cost	Time	Ability	Responsibility
Each Vendor Does Own Modifications	GOOD High equipment utilization, low depreciation costs Test equipment already available Special test setups already available Less total tasks Flexible personnel	GOOD More time to interpret specification More time to make new drawings Shortest time to do modifications Purchasing response short (have equipment in stock)	GOOD Equipment available Skilled personnel available Procedures available Troubleshooting easy	GOOD Warranty in effect
One Centralized Vendor Does All Modifications	FAIR Additional test equipment required Buy certain materials in large quantities Special test setups required Large number of total tasks	FAIR Less time to interpret specifications (repetition) Less time to make new drawings (repetition) Medium time to do mods. Purchasing response long (number of items)	FAIR Equipment must be obtained (part) Skilled personnel not available Procedures not available Troubleshooting difficult	POOR No warranty
NASA Does All Modifications	FAIR Additional test equipment required Buy certain materials in large quantities Special test setups required Large number of total tasks	FAIR Less time to interpret specification (repetition) Less time to make new drawings (repetition) Longest time to do mods. Purchasing response long (number of items)	FAIR Equipment must be obtained (part) Skilled personnel not available Procedures not available Troubleshooting hard	POOR No warranty



Another factor to be considered, however, is the benefit to be expected from a centralized facility for certain common modification operations such as conformal coating application and verification testing (vibration, EMI, and thermal environments). These items would appear to benefit from a centralized facility because of the better quality control (conformal coatings) and a more uniform application of common environments (vibration, thermal) and EMI. Because the benefits of a centralized facility are seen to be intangible, it is difficult to precisely evaluate the degree of improvement which might be expected. Here again, the effect of the equipment acquisition rate is seen to be a negative factor because the fluctuation between high and low extremes would tend to make the operation of a centralized facility quite inefficient. Further study is required to determine whether outside vending of these activities is more cost effective than the buildup of a dedicated facility.

6.0 EQUIPMENT RECOMMENDED FOR PROCUREMENT

Adoption of a policy endorsing the purchase of CAM equipment for Spacelab experiments cannot be based solely upon analytical data. Many of the conclusions drawn during this study require physical validation by testing. NASA-MSFC currently plans to procure certain CAM items to supplement the results of this study providing information for the CAM equipment-use decision. This section presents the items recommended by the study team for procurement.

6.1 CHARACTERISTICS OF EQUIPMENT SELECTED FOR PROCUREMENT

Procurement recommendations are based upon the results of the analysis of 34 selected equipment items, and the study team's understanding of actions taken on two related programs, the Spacelab and the Beckman study, "Analysis of Multipurpose Equipment for Space Applications." In selecting specific items for recommended procurement and test, the following factors were considered:

1. Cost-Effectiveness. The equipment recommended for procurement and test can be cost-effectively modified for suitability in the Spacelab.
2. Frequency of Use. The number of equipment types required by the sortie experiment program, identified in Section 4.1 (Volume III), was reviewed to assure that the equipment recommended for procurement is that which is truly useful to the sortie experiment program, is representative of a broad range of discipline requirements, and accounts for a significant proportion of payload funding requirements.
3. Resolution of Issues for Further Analysis. Equipment items were identified which, if procured and tested, will help resolve the suitability of characteristics which could not be verified during the suitability analysis.
4. Design, Materials, and Processing Differences. Units were selected so that differences due to variations in specifications and construction practices can be assessed. Differences among standard commercial, aircraft, and military hardware must be represented by the selected hardware list.
5. Functional Diversity. It will also likely be useful to procure and test a variety of types of items: electronic (e.g., amplifiers), electro-optical (e.g., lasers, spectrophotometers), electromechanical (e.g., tape recorders), chemical (e.g., gas chromatographs, pH meters) to understand construction differences among hardware types.

6. Modification Feasibility Verification. All items recommended require some modification. The modification activity will physically illustrate and verify study findings.

6.2 TECHNICAL ASSESSMENT OF EQUIPMENT TYPES

Equipment types examined during the study are listed on Table 6-1 in order of preference for procurement on the basis of the technical merit of testing these units. This list reflects the opinion of the study team as to which equipment items would yield the most information when tested. In many cases the type of unit defines the issues that will be examined by testing. For these units, a "no preference" comment has been placed in the Manufacturer column. Obviously, the best correlation with study results will occur if the same unit analyzed is also purchased for testing. When the unit of a manufacturer displays specific characteristics meriting testing, that manufacturer has been identified. For example, the Nuclear Data spectrum analyzer is not only constructed according to NIM packaging specifications, but its wire connections are wire-wrap type instead of the typical solder connections. To assure that both wire-wrap and solder-type connections are tested, a specific selection was made.

The items listed have been grouped into three categories. The first category, Units of Significant Interest, include those items with diverse characteristics yet representative of the hardware that could be placed in the Spacelab. Units with similar characteristics are not repeated. The second group includes items which have characteristics similar to the units in the first group, or items which are expected to have limited use in the Spacelab. For example, many of the electronic packages fall into this category because the computer and the strip chart recorder were selected as more desirable units and placed in the top group. Finally, the three units that could not be cost-effectively modified for installation in the Spacelab were placed in a separate category. These units are excluded from consideration. The characteristics that make the unit interesting and/or the reasoning behind placement of the unit in a given category are stated in the Rationale column.

6.3 PRELIMINARY EQUIPMENT-TYPE ASSESSMENT

Equipment types were assessed for possible procurement based upon potential programmatic savings and the availability of equipment and information from other studies.

The analyzed equipment types were ranked by potential savings to the Spacelab experiment program. The 21 types represented by the 31 units that could be cost-effectively modified are ranked in order of highest to lowest savings to the program in Table 6-2. The magnitude of the savings was determined by multiplying the quantity of hardware estimated for the Spacelab experiment program by the estimated first-unit cost savings. The dewar, refrigerator, and Coulter counter were not included on this list because each unit could not be used in its available configuration, even if modified. It has to be custom-built and is, therefore, out of the scope of the proposed test program.

Table 6-1. Equipment Listed by Technical Importance

Type	Manufacturer/Model No.	Rationale
UNITS OF SIGNIFICANT INTEREST		
Tape Recorder	AMPEX, AR700	Airborne unit, was modified and used on Skylab
Spectrophotometer	No preference	Sensitive lab instrument, gravity-dependent junction
Laser	No preference	Glass components, exacting alignment requirements
Computer	No preference	Typical electronic equipment-type construction
Receiver/Transceiver	Collins, 618M-2	ATR packaging airborne qualification; does not overheat up to 55,000 ft (1.68 km) altitude
Spectrum Analyzer	Nuclear Data, 100	NIM packaging, maximum use of integrated circuits; wire-wrap connections
Strip Chart Recorder	Honeywell, 1856	Electromechanical; CRT containment required; design compatible with zero g--no galvanometer, no ink pen
Oscilloscope	Hewlett-Packard AN/USM-281A	Mil-Spec oscilloscope qualified for field use
Microscope	No preference	Optical unit; addresses glass containment issue for unit that can only be made of glass
Gas Chromatograph	Beckman, 6700	Process-type unit; chemical analysis hardware
Digital Voltmeter	Fluke, 8125	Ruggedized Mil-Spec unit; provides for comparison of Mil-Spec designs
Furnace	Astro, 1000 A	Heavy mechanical unit; touch-temperature problems; water cooled
Centrifuge	No preference	Electric motor with potential EMI problem
UNITS OF INTEREST IF SOME OF ABOVE ARE NOT CHOSEN		
Oscilloscope	Tektronix	Typical electronic instrumentation; sensitive circuits; unit is basic type of laboratory equipment
Digital Voltmeter	No preference	Typical electronic instrumentation; sensitive circuits; unit is basic type of laboratory equipment
Assorted NIM Equipment	Tennelec or Ortec	Typical NIM-type hardware
Electrophoresis Apparatus	No preference	Unit has unique characteristics: gravity-dependent design; high voltages; gravity-dependent batch-type process
Spectrum Analyzer	Fluke 645M or Singer 55B	Mil-Spec electronic units; characteristics already examined by previously identified Mil-Spec electronics
Receiver	No preference	Many airborne and Mil-Spec units available; characteristics examined by Collins unit
Tape Recorder	Honeywell 5600	Portable tape recorder; construction typical of high-quality electronic units
pH Meter	No preference	Simple electronic device; gravity-dependent measurement easily replaced by forced-flow technique
Spectrum Analyzer	Hewlett-Packard	Typical of commercial-type electronics
Amplifier	Neff	Typical of commercial-type electronics
Display Terminal	Research Inc.	Typical of commercial-type electronics; poor suitability for Spacelab use

**ORIGINAL PAGE IS
OF POOR QUALITY**

Table 6-1. Equipment Listed By Technical Importance (Cont)

Type	Manufacturer/Model No.	Rationale
Timer	No preference	Typical of commercial-type electronics
Microtome	Americal Optical, 820	Unit is interesting mechanical unit with major safety and operational problem; downgraded because it is no longer identified as needed for the life science sortie mission
Power Supply	Power Designs	NIM bin power supply; typical electronics; need is dependent on Spacelab power supply
Power Supply	Sorenson	Typical of commercial-type electronics
UNITS WHICH CANNOT BE COST EFFECTIVELY MODIFIED FROM AVAILABLE EQUIPMENT		
Dewar	No preference	Two-phase fluid containment requires new development
Refrigerator	No preference	Zero-g refrigeration cycle requires new development
Blood Cell Counter	Coulter	Totally gravity-dependent with hazardous design for Spacelab environment

The dewar, the refrigerator, and the Coulter counter were not included on this list because each unit could not be used in its available configuration even if it were modified. It has to be custom built and is, therefore, out of the scope of the proposed test program.

Table 6-2. Equipment Types Ranked According to Spacelab Experiment Program Savings

Spectrophotometer
Spectrum/Wave Analyzer
Receiver/Transmitter
Gas Chromatograph
Keyboard Display Terminal
Strip Chart Recorder
Tape Recorder
Furnace
Oscilloscope
Computer
Laser
Particle Counter
Timer
Electrophoresis Apparatus
Amplifier
Microscope
pH Meter
Digital Voltmeter
Power Supply
Centrifuge
Microtome

Table 6-2 was screened; items were eliminated from consideration for the following reasons.

1. MSFC has already obtained the unit from the Beckman study and need not procure an additional item.
2. Current Spacelab documentation indicates that the unit will be included as part of the Spacelab. Therefore, the Europeans will be testing this hardware and MSFC need not expend their funds for duplicate testing.
3. The unit is similar to previously selected hardware. Test data for this item would add little information to the previously selected unit.
4. Few of the given type of equipment are required by Spacelab experiment payloads.

Table 6-3 lists each item and indicates which of the rejection reasons justified elimination of that item.

Table 6-3. Equipment-Type Rejection Rationale

Item	Rationale
Spectrophotometer	Beckman study equipment
Spectrum/Wave Analyzer	
Receiver/Transmitter	
Gas Chromatograph	
Keyboard Display Terminal	Spacelab core hardware
Strip Chart Recorder	
Tape Recorder	Spacelab core hardware
Furnace	
Oscilloscope	Beckman study equipment
Computer	Spacelab core hardware
Laser	
Particle Counter	Represented by spectrum analyzer
Timer	Spacelab core hardware
Electrophoresis Apparatus	Major zero-g development
Amplifier	Spacelab core hardware
Microscope	Issues represented by laser
pH Meter	Zero-g development, low savings
Digital Voltmeter	Beckman study equipment
Power Supply	Spacelab core hardware
Centrifuge	Low program savings
Microtome	Low program savings

6.4 PROCUREMENT RECOMMENDATIONS

The final equipment selection recommendations are presented in Table 6-4. Equipment type, manufacturer, issues represented, the predicted magnitude of the first-unit savings and the quantity of each type of hardware predicted for the Spacelab program are shown for each selection.

Table 6-4. Procurement Recommendations

TYPE	MANUFACTURER	ISSUES REPRESENTED	FIRST-UNIT SAVINGS (\$1000)	QUANTITY
SPECTRUM ANALYZER	NUCLEAR DATA	NIM PACKAGING WIRE-WRAP CONNECTIONS	415-540	58
RECEIVER	COULTER	AIRBORNE EQUIPMENT	150-240	82
GAS CHROMATOGRAPH	BECKMAN	PROCESS-TYPE INSTRUMENTATION CHEMICAL ANALYSIS HARDWARE	355	55
STRIP CHART RECORDER	HONEYWELL	FIBEROPTICS, PAPER WELL SUITED TO SPACELAB	335	58
FURNACE	ASTRO	HIGH TEMPERATURE WATER COOLING	425	20
LASER	NO PREFERENCE	FRAGILE OPTICS SENSITIVE ALIGNMENT PROBLEMS	135	39
OSCILLOSCOPE	HEWLETT-PACKARD	MIL-SPEC DESIGN (ARMY FIELD USE) COMPARISON WITH COMMERCIAL TEKTRONIX	115-290	26
DIGITAL VOLTMETER	FLUKE	MIL-SPEC DESIGN (SUBMARINE) COMPARISON WITH COMMERCIAL FLUKE	17-25	21

The list represents a composite of the cost effectiveness and technical comparisons. Units of significant interest on Table 6-1 were selected if they were not eliminated by the rationale described in Table 6-3. This approach placed emphasis on defining selections which maximized the information return for as low a cost as possible.

The top six units are the units selected by the analysis in Section 6.2. Review of the characteristics represented by these units indicated that no Mil-Spec units had been selected. Therefore, two additional units were included to provide such representation.

Two Mil-Spec units were selected rather than one for the following reasons:

1. Each unit is representative of a different class of military specification. One unit is built for army field use, while the other is built for submarine use.
2. MSFC already has commercial counterparts to these units and the results of comparative testing should be very informative.



3. These units are typically found in all laboratories so they merit having some additional samples for testing.

Manufacturers identified are the same as those examined during the study so that correlation of test and analytical results can be obtained. The laser is the only exception. In the opinion of the study team, similar conclusions can be obtained from testing of cheaper units than the \$10,000 Sylvania CO₂ laser examined during the study.

The list of issues represented shows that selected hardware satisfies the goals of the originally defined objectives. All units require some modification. It is strongly recommended that these modifications be accomplished prior to testing. A test program without hardware modification will not verify the study conclusions. The most complete test program would purchase two units. One would be modified and the other left unmodified. Both units would then be subjected to the same testing program to verify the effectiveness of the modification program. If funding is not available for two sets of equipment, then modified hardware should be tested.



7.0 PROGRAMMATIC SAVINGS

Sections 3.0 and 4.0 indicated that it is cost-effective to modify commercially available equipment and instrumentation rather than procuring custom-built hardware. This section discusses the potential savings that could accrue to the entire shuttle experiment program if the savings on the selected equipment items are extended to the program level. All extensions apply to the segment of experiment equipment installed on the Spacelab which does not require functional modification.

7.1 PROGRAMMATIC EXTENSION SCENARIOS

The extension of specific equipment costs to the program level is dependent upon the scenario assumed for experiment equipment development and replacement. Dependence on hardware developed in other manned spacecraft programs could reduce the cost of custom-built hardware significantly. If, on the other hand, Spacelab experiment hardware and supporting equipment requirements are unique relative to other space programs, new developments will be required escalating the total cost of the experiment program. Equipment replacement rates establish the number of sets of equipment required by the program. An infrequent rate of equipment replacement would reduce total equipment costs. The definition of new laboratories during the "Shuttle Era" would affect both the equipment replacement rate and number of new developments required for the program.

Figure 7-1 illustrates the scenarios analyzed during the study. Two extremes of development commonality were assessed. Maximum commonality, and therefore fewer new developments, with other Space Transportation System (STS) hardware was examined. The other extreme is represented by examining the cost of all experiment equipment if each equipment type required a new development. This latter case is unrealistic, but it indicates an upper limit on the savings that could be achieved by using CAM equipment instead of custom-built hardware.

Replacement rate alternatives are also shown on the figure. The five-year-life study guideline for replacement of equipment was applied in both cases. Exceptions were made for a few equipment items which either have not experienced rapid technological change historically, or would not provide increased experiment benefits if improved items were used. The microtome is an example of a unit that has not experienced rapid technological changes. Its current design is essentially the same as the design originally developed in 1902. Certain support equipment not directly related to experiment equipment may not need replacement if increased performance requirements do not change. Power and weight considerations may, however, still dictate change-out of these units if lighter weight, lower powered units should be developed.




Figure 7-1. Programmatic Extension Scenarios

The study is based upon the payloads currently defined for the Spacelab. New laboratories will most probably be identified within the next few years, and most certainly will occur when the sortie experiment program gathers momentum in the 1980's. The Zero-G Cloud Physics Laboratory cannot be the only earth observations sortie laboratory. Also, the equipment complements identified for the currently defined payloads will change as Spacelab operations progress. The second replacement rate scenario assumes that a new equipment complement will be identified for each discipline every five years, either in the form of an entirely new laboratory or changeover of an existing laboratory.

7.2 SAVINGS PROJECTIONS

Projections of custom and modified equipment costs were based upon the relationship of all equipment identified in the functional requirements data bank that do not require functional modification to the selected equipment items. These relationships are explained in Section 4.0 of Volume III. The quantities identified for each equipment type were grouped with the selected equipment item that it most closely resembled. The modification and custom costs associated with each selected item were assumed to apply to all items similar to it. Total costs were arrived at by multiplying the total number of similar equipment types by the custom and modification costs associated with the related, selected equipment item.

The total cost for a given equipment type included the development costs associated with each unit having different functional requirements and the recurring cost of the repeated procurements of the same hardware. Development costs for the custom and modified hardware are shown in Section 3.0 of this

volume. Recurring costs for each custom-built unit were estimated, based upon the relationships of development and recurring cost CER's (explained in Figure 3-31). These costs were adjusted by assuming manufacturing lot sizes consistent with equipment quantity needs of the program. The cost of each lot was reduced to reflect learning curve affects. Recurring costs for modified hardware included all fabrication costs, but only included 25 percent of estimated engineering and documentation costs to provide for changes occurring in available hardware over the years. Custom-built hardware designs were assumed not to change when repeat sets of equipment were required. Available hardware can take advantage of technology changes without impacting costs, while custom-built hardware would have to be redesigned, incurring additional development costs. Conversely, custom-built hardware could be designed for a longer life-time than five years providing some programmatic savings by reducing recurring costs.

All modification costs assumed application of the SEEIR specification for equipment modification. Modified hardware costs would be higher by a factor of 2 if the EC006 specification had been used. The total cost of modified hardware would be \$42.4 million if only one set of new developments is required, and \$68 million if a new development is required every five years if this specification were in force. A slight reduction in the programmatic savings would result from the use of these modified hardware costs.

Figure 7-2 shows the equipment cost expenditures by year for each of the scenarios. All costs are assumed to occur in the year prior to the year the equipment is scheduled to fly. The areas under each curve yield the total programmatic cost. All costs are in 1974 dollars. If the costs were in actual dollars, the savings would be much higher. The higher funding costs in 1979 and 1980 on Figure 7-2 (a) and (b) result from initial hardware development. These peaks could be impacted by funding needs of the Shuttle program. The peaks could be lowered if the hardware development were spread out over two or three years rather than occurring in the year prior to flight. Lower funding in subsequent years results from only incurring recurring-type costs after the initial development. No funding is required in 1982 because no equipment is needed in 1983, assuming the five-year replacement guideline.

Figure 7-2 (a) and (b) depicts the funding requirements if no new payloads or equipment complements occur throughout the program. The \$230 million saved is the lowest possible savings that could result from the use of CAM equipment rather than custom-built equipment. This savings increases to \$510 million if unique developments are required for all Spacelab experiment equipment and instrumentation.

Figure 7-2 (c) and (d) reflects the programmatic savings if new payloads or equipment complements are needed every five years. It is believed that the nearly \$400 million savings shown on (c) is closest to the real world. The savings of \$1.2 billion shown in (d) is an extreme case and very unlikely. It is shown only to provide an upper limit on the savings that could be generated. This cost indicates that every effort to use already-developed hardware should be made to lower Spacelab experiment costs.





Figure 7-2. Programmatic Cost Comparison



Space Division
Rockwell International





Figure 7-2. Programmatic Cost Comparison (Cont)

8.0 CONCLUSIONS

The following major conclusions were arrived at during this study.

1. The use of commercially available equipment in the Spacelab is both feasible and cost-effective. Sortie experiment program savings achieved by the use of CAM experiment equipment are between \$230 million and \$1230 million in 1974 dollars. These savings are dependent upon the development commonality and the equipment replacement rates assumed. The most realistic scenario indicates a savings of approximately \$400 million. Further savings may be achieved by using available equipment for pallet-mounted experiments.
2. Between 55 and 65 percent of the equipment defined by all the recent Shuttle sortie payload studies can be accommodated by CAM equipment without functional modification. A more extensive search of supplier data may increase this percentage. The cost savings reported in this study apply only to the equipment items that do not require functional modification. If this quantity of equipment were increased, greater programmatic savings would result.
3. Modifications to a specification compatible with the use of commercial equipment will add approximately 4 to 8 times the retail cost to each item. Modification to a typical aerospace specification will add from 10 to 18 times the retail cost to the item. However, custom-built hardware is an average of 36 times as expensive as the CAM unit. Therefore, no matter which specification is in effect, it is cheaper to use modified available hardware. These savings occur because the functional aspects of the hardware do not have to be redeveloped to be within a compact envelope, meet a minimum weight or power target, or be compatible with a different power form than used in its original design.
4. A specification compatible with the use of commercial equipment can be written which is cost-effective as well as technically sound and compatible with crew safety requirements. The use of commercially available hardware obviates the need for many typical specification requirements which are formulated because order lot sizes are small (such as configuration control on replacement parts). Also, CAM equipment has prior operational experience, reducing the need for performance verification testing and reporting.



5. The most cost-effective modifier of the equipment to Spacelab specifications is the original manufacturer of each item. His knowledge of his own equipment and its immediate availability from inventory, make the original supplier both a cost-effective and timely modifier of his equipment.

APPENDIX A

EXPERIMENT EQUIPMENT FUNCTIONAL REQUIREMENTS LISTED BY SCIENCE AND TECHNOLOGY DISCIPLINE

ASTRONOMY

DISPLAY CRT (CONTROL/DISPLAY CONSOLE) AS001 21111
PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS. INCLUDES REFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.

SCREEN SIZE 14 IN (0.355 M)
SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) AS002 11111
PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CRT

CHANNELS 2
FORMAT VIDEO AND DATA
VIDEO INPUT PRESENTATION RASTER SCAN
SYMBOL WRITING TECHNIQUE STACK
INTERFACE DESCRIPTION DIGITAL-12-BIT DATAWORD

FUNCTION KEYBOARD AS003 11111
ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY, FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.

KEYBOARD TYPEWRITER TYPE
INTERFACE DESCRIPTION DIGITAL-32-BIT WORD

ALPHANUMERIC KEYBOARD AS004 11111
ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE ONBOARD COMPUTER FOR EXPERIMENT CONTROL AND DATA ANALYSIS.

KEYBOARD TYPEWRITER TYPE
INTERFACE DESCRIPTION DIGITAL-32-BIT WORD

MICROFILM VIEWER AS005 11111
PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUBSYSTEM OPERATIONAL PROCEDURES, ONBOARD CHECKOUT PROCEDURES, SIMPLEX SCHEMATICS AND OTHER WRITTEN OR PICTORIAL INFORMATION.

FILM FORMAT 16 MM DUAL TRACK
FILM LOADING CASSETTE
FILM SLEWING MANUAL AND COMPUTER SELECT
INTERFACE DESCRIPTION DIGITAL 13-BIT BINARY

TIMER - EVENT AS006 11111
DIGITAL DISPLAY OF TIME REMAINING OR EXPENDED FOR A PARTICULAR EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE START AND STOP COMMANDS (MINUTES AND SECONDS).

DISPLAY 4 DIGITS
DISPLAY TYPE LED

TIMER - MISSION AS007 11111
PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECOND UPDATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE, SECOND.

DISPLAY 7 DIGITS BCD
DISPLAY TYPE LED

CONTROLLER, THREE-AXIS HAND AS008 11132
THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT POINTING AND INITIAL TARGET ACQUISITION.

CONTROLLED AXES 3

CIRCUIT BREAKER/DISTRIBUTOR PANEL AS009 11131
CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO THOSE STANDARD ITEMS FOUND IN ALL ASTRONOMY EXPERIMENTS.

CIRCUIT BREAKERS 18
POWER CAPABILITY 10 AMPS

INDICATORS-CALLION AND WARNING AS010 11132
A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF SUBSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.

INDICATORS 40
..... PASTER BLARM MEMORY
..... POWER/TEST CONTROLS

RECORDER-TAPE AS011 11111
RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND DISPLAY CONSOLE.

BANDWIDTH 60 TO 200 KHZ
CHANNELS 14
TIME MARKS EDGE TRACK
RECORDED FORMAT DIGITAL

ANNUNCIATORS-ADVISORY AS012 21132
DUAL BANKS MOUNTED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY ONBOARD EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.

ANNUNCIATOR-VISUAL 12 IN DUAL BANKS
DISPLAY TYPE RED, WHITE, GREEN
ANNUNCIATOR-AUDIO NEPA

ORIGINAL PAGE IS
OF POOR QUALITY



HIGH ENERGY PHYSICS

DISPLAY CRT (CONTROL/DISPLAY CONSOLE) ME001 21111

PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS. INCLUDES DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.

SCREEN SIZE 14 IN (0.355 M)

SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) ME002 11111

PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CRT

CHANNELS 2
FORMAT VIDEO AND DATA
VIDEO INPUT PRESENTATION PASTER SCAN
SYMBOL WRITING TECHNIQUE STROBE
INTERFACE DESCRIPTION DIGITAL-12-BIT DATA/CRG

FUNCTION KEYBOARD ME003 11111

ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY, FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.

KEYBOARD TYPEWRITER TYPE
INTERFACE DESCRIPTION DIGITAL-32-BIT WORD

ALPHANUMERIC KEYBOARD ME004 11111

ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE ONBOARD COMPUTER FOR EXPERIMENT CONTROL AND DATA ANALYSIS.

KEYBOARD TYPEWRITER TYPE
INTERFACE DESCRIPTION DIGITAL-32-BIT WORD

MONOFILM VIEWER ME005 11111

PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUBSYSTEM OPERATIONAL PROCEDURES, ONBOARD CHECKOUT PROCEDURES, SIMPLER SCENARIOS AND OTHER WRITTEN OR PICTORIAL INFORMATION.

FILM FORMAT 16 MM DUAL TRACK
FILM LOADING CASSETTE
FILM SLEWING MANUAL AND COMPUTER SELECT
INTERFACE DESCRIPTION DIGITAL 13-BIT BINARY

TIMER - EVENT ME006 11111

DIGITAL DISPLAY OF TIME REMAINING OR EXPANDED FOR A PARTICULAR EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE START AND STOP COMMANDS (MINUTES AND SECONDS).

DISPLAY 4 DIGITS
DISPLAY TYPE LED
INTERFACE DISCRETE

TIMER - MISSION ME007 11111

PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECOND UP-DATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE, SECOND.

DISPLAY 7 DIGITS DEC
DISPLAY TYPE LED

CONTROLLER, THREE-AXIS HAND ME008 11112

THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT POINTING AND INITIAL TARGET ACQUISITION.

CONTROLLED APES 3

CIRCUIT BREAKER/DISTRIBUTION PANEL ME009 11111

CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO THESE STANDARD ITEMS FOUND IN ALL ASTROPHYSICS EXPERIMENTS.

CIRCUIT BREAKERS 18
POWER CAPABILITY 10 AMPS

INDICATORS-CALLION AND WARNING ME010 11112

A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF SUBSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.

INDICATORS 40
..... MASTER ALARM MEMORY
..... POWER/TEST CONTROLS

RECORDER-TAPE ME011 11111

RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND DISPLAY CONSOLE.

BANDWIDTH DC TO 200 KHZ
CHANNELS 14
TONE MARKS EDGE TRACK
RECORDED FORMAT DIGITAL

ANNUNCIATORS-ADVISORY ME012 11112

DUAL BARS MOLATED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY ONBOARD EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.

ANNUNCIATOR-VISUAL 12 IN OVAL BARS
DISPLAY TYPE RED, WHITE, GREEN
ANNUNCIATOR-AUDIC BELL

ORIGINAL PAGE IS
OF POOR QUALITY



SCLAR PHYSICS

DISPLAY CRT (CONTROL/DISPLAY CONSOLE) SC001 21111

PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS. INCLUDES REFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.

SCREEN SIZE 14 IN (10.955 IN)

SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) SC002 11111

PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CRT

CHANNELS 2
FORMAT VIDEO AND DATA
VIDEO INPUT PRESENTATION RASTER SCAN
SYMBOL WRITING TECHNIQUE STROKE
INTERFACE DESCRIPTION DIGITAL-12-BIT CATHODE

FUNCTION KEYBOARD SC003 11111

ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY, FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.

KEYBOARD TYPEWRITER TYPE
INTERFACE DESCRIPTION DIGITAL-32-BIT WORD

ALPHANUMERIC KEYBOARD SC004 11111

ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE ONBOARD COMPUTER FOR EXPERIMENT CONTROL AND DATA ANALYSIS.

KEYBOARD TYPEWRITER TYPE
INTERFACE DESCRIPTION DIGITAL-32-BIT WORD

PICOPFILM VIEWER SC005 11111

PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUBSYSTEM OPERATIONAL PROCEDURES, ONBOARD CHECKOUT PROCEDURES, SIMPLEX SCHEMATICS AND OTHER WRITTEN OR PICTORIAL INFORMATION.

FILM FORMAT 16 MM DUAL TRACK
FILM LOADING CASSETTE
FILM SLEWING MANUAL AND COMPUTER SELECT
INTERFACE DESCRIPTION DIGITAL 13-BIT BINARY

TIMER - EVENT SC006 11111

DIGITAL DISPLAY OF TIME REMAINING OR EXPENDED FOR A PARTICULAR EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE START AND STOP COMMANDS (MINUTES AND SECONDS).

DISPLAY 4 DIGITS
DISPLAY TYPE LED

TIMER - MISSION SC007 11111

PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECOND UPDATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE, SECOND.

DISPLAY 7 DIGITS BCD
DISPLAY TYPE LED

CONTROLLER, THREE-AXIS HAND SC008 11132

THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT POINTING AND INITIAL TARGET ACQUISITION.

CONTROLLED AXES 3

CIRCUIT BREAKER/DISTRIBUTOR PANEL SC009 11131

CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO THOSE STANDARD ITEMS FOUND IN ALL ASTRONAUTIC EXPERIMENTS.

CIRCUIT BREAKERS 18 EACH
POWER CAPABILITY 10 AMPS

INDICATORS-CAUTION AND WARNING SC010 11132

A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF SUBSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.

INDICATORS 40
..... MASTER ALARM MEMORY
..... POWER/TEST CONTROLS

RECORDER-TAPE SC011 11111

RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND DISPLAY CONSOLE.

BANDWIDTH 60 TO 200 KHZ
CHANNELS 14
TAPES EDGE TRACK
RECORDED FORMAT DIGITAL

ANNUNCIATORS-ADVISORY SC012 21132

DUAL BANKS MOUNTED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY ONBOARD EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.

ANNUNCIATOR-VISUAL 12 IN DUAL BANKS
DISPLAY TYPE RED, WHITE, GREEN
ANNUNCIATOR-AUDIO MORN

ORIGINAL PAGE IS
OF POOR QUALITY



ATMOSPHERIC AND SPACE PHYSICS

DEEP POSITION TV DISPLAY	AP001	11111
DISPLAY MEDIUM RESOLUTION TELEVISION OF EXPERIMENT DEEP POSITIONS BEFORE, DURING AND AFTER EACH EXPERIMENT REPEATITION.		
RESOLUTION 525 LINES		
FRAME RATE 30 FRAMES/SEC		
BANDWIDTH VHF RANGE		
EXPERIMENT TV DISPLAY	AP002	11111
DISPLAY HIGH RESOLUTION TELEVISION PICTURES OF AREAS OF INTEREST IN REAL TIME ON PHUM VIDEO RECORDERS AT VARIOUS SPEEDS.		
RESOLUTION 525 LINES		
FRAME RATE 30 FPS		
BANDWIDTH UHF RANGE		
SPECTRUM ANALYZER	AP003	11111
MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE RATIO AND MULTIPLE RATIO OF VARIOUS FREQUENCY SOURCES.		
BANDWIDTH 20 HZ TO 300 KHZ		
SENSITIVITY 20 MICROVOLTS		
ACCURACY +/- 0.6 DB		
MULTI-CHANNEL ANALYZER	AP004	11111
MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE RATIO AND MULTIPLE RATIO OF MULTIPLE FREQUENCY SOURCES SIMULTANEOUSLY.		
BANDWIDTH 0 TO 10 KHZ		
TIME BASE 200 NHz CLOCK		
MEMORY SIZE 1048, 4096, 8192 CHANNELS		
WAVE ANALYZER	AP005	11111
DETECTION OF SIGNAL AMPLITUDE AND FREQUENCY INFORMATION		
BANDWIDTH 1C TO 600 KHZ		
FREQUENCY RESOLUTION 10 HZ		
SENSITIVITY 3 MV TO 30 V FS		
CHANNELS 2 MINIMUM		
WAVE ANALYZER	AP006	11111
DETERMINE, SEPARATE AND ANALYZE VARIOUS FREQUENCY COMPONENTS OF INPUT SIGNALS (IE FUNDAMENTAL, HARMONICS, INTERMODULATION PRODUCTS, ETC.).		
BANDWIDTH 1 KHZ TO 1.5 KHZ		
FREQUENCY RESOLUTION +/- 1.0 %		
SENSITIVITY 10 MV TO 100 V		
CHANNELS 2 MINIMUM		
PATCH PANEL, COAXIAL	AP007	11111
MULTIPLE INPUT PANEL CAPABLE OF CHANNELING HIGH FREQUENCY RF ENERGY TO APPROPRIATE INSTRUMENTS WITH MINIMUM LOSS AND INTERFERENCE.		
FREQUENCY COUNTER	AP008	11111
DIRECT MEASUREMENT OF FREQUENCY AND/OR SIGNAL REPETITION RATE.		
BANDWIDTH 0 TO 35 KHZ		
SENSITIVITY 10 MV RMS		
TIME BASE 100KHZ TO 10 MHZ		
AUTOMATIC DISPLAY GENERATOR	AP009	11111
PROVIDES 3U DISPLAY		
CAMERA STILL	AP010	11111
RECORD EXPERIMENT PHENOMENA		
FILM SIZE 35MM		
SHUTTER SPEED 1 SEC TO 1/500 SEC		
X-Y RECORDER	AP011	21111
REAL-TIME ACCURATE REPRODUCTION OF SPECTRAL DATA SOURCES PLOTTED ON CARTESIAN COORDINATES.		
CHANNELS 1		
ACCURACY +/- 0.2 % FS		
RECORDER - STRIP CHART	AP012	21111
RECORD HARD COPY DATA FROM LOW FREQUENCY ANALOG INVESTIGATIONS		
FREQUENCY RESPONSE DC TO 150 KHZ		
CHANNELS EVENT, TIME, 2- DATA		
ACCURACY +/- 0.5 PCT		
RECORDER - STRIP CHART	AP013	21111
RECORD HARD COPY DATA.		
CHANNELS 1		
TAPE RECORDER - ANALOG	AP014	11111
RECORD AND PLAYBACK OF ANALOG DATA IN SUPPORT OF ATMOSPHERIC AND SPACE PLASMA PHYSICS SCIENTIFIC LABORATORY INVESTIGATIONS.		
BANDWIDTH 2 KHZ		
CHANNELS 16		
RECORDING FORMAT PP		

OSCILLOSCOPE	AP015	21111
MONITOR, MEASURE AND MAINTAIN ELECTRONIC EQUIPMENT OPERATING		
BANDWIDTH 500 MHz		
CHANNELS 2		
DEFLECTION FACTOR 1 MV		
PHOTOGRAPHIC CAPABILITY YES		
CAMERA STILL	AP016	21111
PLANAR CAMERA TO RECORD TRACKS ON OSCILLOSCOPE		
FILM FORMAT 35MM 137 BLACK/WHITE		
LENS 3.25X4.25X102.55X162.55MM		
SHUTTER SPEED 1/100 TO 1/200		
FIELD OF VIEW 40 DEGREES		
TAPE RECORDER - DIGITAL	AP017	21111
RECORDING AND PLAYBACK OF DIGITAL DATA COMPATIBLE WITH ANY CAMBOARD COMPUTING EQUIPMENT.		
BANDWIDTH 2 MHz		
CHANNELS 8		
SIGNAL 0 TO 5 VDC DIGITAL		
COMPUTER	AP018	21111
COMPUTER CAPABILITY TO SUPPORT ATMOSPHERIC AND SPACE PLASMA PHYSICS INVESTIGATIONS.		
WORD SIZE 16 BIT		
IP/OP MEMORY TRANSFER 16K-16 BIT WORDS/SEC		
CORE MEMORY 32K-64 KWORDS		
IP/OP CHANNELS 45		
KEYBOARD DISPLAY TERMINAL	AP019	11111
GENERAL PURPOSE INPUT KEYBOARD ALLOWING FOR PROGRAMING, DATA INPUT AND UPDATE, AND DATA MANIPULATION FOR ONBOARD ANALYSIS		
STATUS PANEL	AP020	21111
DISPLAY TERMINAL WITHOUT KEYBOARD PROVIDING STATUS OF EXPERIMENT FUNCTIONS		
SPECIAL DATA ACQUISITION PANEL	AP021	11120
BATTERY CONTROL AND MONITOR	AP022	11114
CAMERA - LINE	AP023	11111
TIME CODE GENERATOR AND DISPLAY	AP024	11111
ACCURATE FREQUENCY, TIME INTERVAL AND TIME KEEPING CAPABILITIES.		
OUTPUTS 100 KHZ, 1 MHz, 5 MHz		
ACCURACY 5E-10 PARTS PER DAY		
REMOTE SENSING PLATFORM GIMBAL CONTROL	AP025	11114
ALV NUMERICAL INCIDENCE SPECTROMETER CONTROL	AP026	11114
UV-VIS-NIR SCANNING SPECTROMETER CONTROL	AP027	11114
HIGH RESOLUTION FOURIER SWIR SPECTROMETER CONTROL	AP028	11114
CIRCULAR IR FOURIER SPECTROMETER CONTROL	AP029	11114
IR RADIOMETER CONTROL	AP030	11114
FAIRBANKS INTERFEROMETER CONTROL	AP031	11114
UV-VIS OCCURRENCE CAMERA CONTROL	AP032	11114
ELECTROSTATIC ANALYZER CONTROL	AP033	11114
MAGNETIC ANALYZER CONTROL	AP034	11114
KEY-REV PARTICLE DETECTOR CONTROL	AP035	11114
TOTAL ENERGY DETECTOR CONTROL	AP036	11114
LICAR AND GIMBAL PLATFORM MONITOR	AP037	11114
TRANSMITTER - 0.2 TO 2.0 MHZ	AP038	11111
VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.		
FREQUENCY BANDWIDTH 0.2 TO 2.0 MHZ		
TRANSMITTER POWER 10 W		
TRANSMITTER - 2.0 TO 20.0 MHZ	AP039	11111
VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.		
FREQUENCY BANDWIDTH 2.0 TO 20.0 MHZ		
TRANSMITTER POWER 10 W		
TRANSMITTER - 0.3 TO 20.0 MHZ	AP040	11111
VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.		
FREQUENCY BANDWIDTH 0.3 TO 20.0 MHZ		
TRANSMITTER POWER 10 W		
ELECTROSTATIC WAVE TRANSMITTER	AP041	11114
COUPLER - ANTENNA (0.4 TO 2.0 MHZ)	AP042	11111
ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.		

FREQUENCY BANDWIDTH 0.2 TO 20.0 MHZ
COUPLER - ANTENNA (2.0 TO 20.0 MHZ) AP910 11111
ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.
FREQUENCY BANDWIDTH 2.0 TO 20.0 MHZ
COUPLER - ANTENNA (0.3 TO 200 KHZ) AP911 11111
ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.
FREQUENCY BANDWIDTH 0.3 TO 200 KHZ
BANDPASS FILTER - WAVE ANALYSIS AP912 21111
DETERMINATION OF FREQUENCY AND AMPLITUDE INFORMATION ABOUT THE CARRIER AND SIDEBANDS.
FREQUENCY RANGE 300 HZ TO 20.0 MHZ
ATTENUATION 75 DB
ACCURACY +/- 2 %
PULSE WAVE ANALYZER- C.W. MODES AP913 11111
PATCH PANEL AP914 51111
INTERFACE CONNECTIONS OF VARIOUS RF RECEIVERS WITH SELECTED TEST EQUIPMENT FOR WAVE ANALYSIS.
AMPLIFIERS - WAVE ANALYSIS AP915 21111
FREQUENCY SYNTHESIZER - WAVE ANALYSIS AP916 21111
TRANSLATION OF A STABLE FREQUENCY OF A PRECISE STANDARD TO ANY OF A SELECTED INVESTIGATION REQUIREMENT.
FREQUENCY BANDWIDTH 0.0 TO 13 MHZ IN 5 RANGES
FREQUENCY RESOLUTION 0.01 TO 10 KHZ
FREQUENCY STABILITY 1E-7 PARTS PER DAY
MAIN BOOM A CONTROL AP917 11134
PLATFORM BOOM A CONTROL AP918 11134
MAIN BOOM B CONTROL AP919 11134
ALIGNMENT TV CONTROL AP920 11134
GIMBALLED PLATFORM CONTROLS AP921 11134
SM BOOM CONTROL AP922 11134
BOOM A POWER SUPPLY & DATA SYS CONTROL AP923 11134
BOOM B TARGET CONTROL AP924 11134
ACCEL-DECEL CONTROL AP925 11134
DISCHARGE FILAMENT HEATER CONTROL AP926 11134
DISCHARGE POTENTIAL CONTROL AP927 11134
PULSE SEQUENCE & BURST LENGTH CONTROL AP928 11134
GAS SELECTION & PRESSURE CONTROL AP929 11134
NEUTRALIZER EMISSION & BIAS CONTROL AP930 11134
CHARGE EXCHANGE CHANNEL ACTUATOR CONTROL AP931 11134
BEAM CURRENT MONITOR AP932 11134
ELECTRON BEAM VOLTAGE; CURRENT HEATER CONTROL AP933 11134
ELECTRON BEAM BURST LENGTH & MAGNITUDE CONTROL AP934 11134
ELECTRON BEAM EXPANSION LENS CONTROL AP935 11134
ELECTRON BEAM CURRENT MONITOR AP936 11134
PHASE, ANGLE, CURRENT MONITOR AP937 11134
MPG ARC VOLTAGE LEVEL CONTROL AP938 11134
MPD ARC BURST CURRENT DURATION CONTROL AP939 11134
MPG ARC PULSE SEQUENCE AP940 11134
MPD ARC BEAM CURRENT MONITOR AP941 11134
GAS SELECTION & PRESSURE CONTROL AP942 11134
SPHERICAL ION PROBE CONTROL AP943 11134
CYLINDRICAL ION PROBE CONTROL AP944 11134
PLANAR ELECTRON PROBE CONTROL AP945 11134
SEGMENTED PLANAR PROBE CONTROL AP950 11134
DEP. SAT. INSTRUMENT CONTROL & HOUSEKEEPING AP955 11134
DEP. SAT EJECTION MECHANISM CONTROL AP960 11134
PHOTOMETER HIGH VOLTAGE SUPPLY CONTROL AP961 11134
PHOTOMETER AMPLIFIERS CONTROL AP962 11134
TV SYS CONTROL - IMAGE INTENSIFIER CONTROL AP963 11134
CANNISTER EJECTION CONTROL AP964 11134

PROJECTILE CAMERA CONTROL AP965 11134
CANNISTER MONITOR AP966 11134
SHAPED CHARGE EJECTION CONTROL AP967 11134
SHAPED CHARGE MONITOR AP968 11134
BALLOON EJECTION CONTROL AP969 11134
GAS CONTROL SYSTEM AP970 11134

ORIGINAL PAGE 15
OF POOR QUALITY

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

EARTH OBSERVATIONS

HEAT EXCHANGER EC001 11121
EXCHANGER TO TRANSFER HEAT FROM THE CLOUD PHYSICS LAB TO THE SPACELAB HEAT REJECTION SUBSYSTEM
HEAT TRANSFER RATE.....185 BTU/HR (54 WATTS THERMAL)
SPACELAB COOLANT INLET TEMP.....43F (4C)
CPL COOLANT OUTLET TEMP.....45F (7C)

HEAT EXCHANGER EC001 11121
PROVIDE HOT FLUID TO CPL
CPL COOLANT OUTLET TEMP.....100F (38C)
HEAT RATE (EST).....170BTU/HR (50 WATTS THERMAL)

RESERVOIR EC001 11113
STORES COOLANT.

PUMP EC001 11111
PUMP COOLANT IN CPL COOLANT SUBSYSTEM
FLOW RATE(EST).....40LBS/HR (300CC/MIN)
HEAD(EST).....1 PSI (6.9E03 N/M2)

TEMPERATURE CONTROLLER EC001 11111
CONTROLS COOLANT TEMPERATURE TO CONTROL CHAMBER WALL TEMPS
TEMPERATURE RANGE.....-76 TO 104F (-60 TO 40C)
ACCURACY.....+/-0.2F (0.1C)

CRYOGENIC COOLER ASSEMBLY EC001 11114
PROVIDE LOW TEMPERATURE HEAT SINK TO EXPERIMENT CHAMBERS
MAXIMUM TEMPERATURE.....-76F (-60C)
HEAT SINK CAPACITY(EST).....1000BTU (290W THERMAL)

TEMPERATURE SENSORS EC001 41111
MEASURE COOLANT TEMPERATURES
OPERATING TEMPERATURE RANGE.....-76 TO 104F (-60 TO 40C)
ACCURACY.....+/-0.2F (+/-0.1C)

CRT DISPLAY EC001 11111
PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERIC, DYNAMIC AND DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.

DIGITAL VISUAL DISPLAY EC001 11111
DISPLAY NUMERIC DATA

ALPHANUMERIC KEYBOARD EC001 11111
ALLOWS CREWMAN TO COMMUNICATE WITH THE ONBOARD COMPUTER FOR EXPERIMENT CONTROL AND DATA ANALYSIS
KEYBOARD.....TYPEWRITER TYPE
INTERFACE DESCRIPTION.....DIGITAL-32-BIT WORD

SOLENOID VALVES EC001 101121
CONTROL FLUID FLOW

MODULATING VALVE EC001 21121
ADJUST FLUID FLOW

PRESSURE REGULATORS EC002 51121
PRESSURE REGULATING VALVE FOR GASEOUS PRESSURIZATION SUBSYSTEM
REGULATION RANGE.....2.7 TO 14.7 PSIA (14G TO 760TORR)
UPSTREAM PRESSURE(EST).....20 PSI (1.4E05 N/M2)

PRESSURE REGULATOR EC002 11121
VALVE REGULATING SUPPLY PRESSURE
REGULATION PRESSURE(EST).....20 PSI (1.4E04 N/M2)
MAX UPSTREAM PRESSURE.....300 PSI (2.1E06 N/M2)

SOLENOID VALVES EC002 111121
GAS SHUTOFF VALVES

CONDITIONING CHAMBER EC015 11134
THE CONDITIONING CHAMBER IS MADE UP OF A PRESSURE SHELL, A POSITIVE EXPULSION REGULATOR AND A HV AEROSOL CONDENSER. THIS CHAMBER CONDITIONS THE GAS SAMPLE PRIOR TO ITS INJECTION INTO THE CLOUD CHAMBER
OPERATING PRESSURE(EST).....20PSIA (1.4E03N/M2)
OPERATING TEMPERATURE.....32 TO 86F (0 TO 30C)
VOLUME.....2 CU.FT. (0.0565 CU. METERS)

PRESSURE SENSORS EC002 51111
MEASURE PRESSURES THROUGH OUT THE SYSTEM
RANGES.....2 TO 20 PSIA (100 TO 1000 TORR)
.....10 TO 50 PSI (6.9E4-3.4E4N/M2)
.....10 TO 300 PSI (3.4E4-2.1E6N/M2)
ACCURACY.....1% OF FULL SCALE

HUMIDIFIER ASSEMBLY

EC003 11135

THIS ASSEMBLY CONSISTS OF AN EVAPORATOR, A HEATER, A WATER TRAP, A METERING PUMP AND A WATER TANK WITH AN INTEGRAL BLEEDER. THIS ASSEMBLY GENERATES THE WATER VAPOR NECESSARY FOR CLOUD CREATION.

DRY AIR FLOW RATE.....175CFM (5 LITERS/MIN)
MAXIMUM TEMPERATURE.....77F (45C)
HUMIDITY RANGE.....0 TO 100%
HEATER CAPACITY(MIN).....178BTU/HR (5W THERMAL)
HEATER FEED RATE (MAX).....0.02 LB/HR (9.1 CC/HR)
WATER TANK CAPACITY.....3LES (1.35KG)

DEW POINT SENSOR

EC003 11113

THIS SENSOR MEASURES THE DEW POINT OF THE GAS SAMPLE

MEASUREMENT RANGE.....10 TO 90%
ACCURACY.....+/-1%

LIQUID WATER CONTENT METER

EC003 11235

MONITORS LIQUID WATER CONTENT WITHIN THE EXPANSION CHAMBER FOR THE PURPOSE OF ACCURATE WATER BUDGET ACCOUNTS. CURRENT CONCEPTS USE OPTICAL TECHNIQUES TO MEASURE THE WATER DROPLETS IN AIR.

ACCURACY.....+/-0.05%

CONDENSER - WATER

EC003 11131

COLLECTS WATER VAPOR BY CONDENSING WATER CONTENT OF GAS SAMPLE

CONDENSATION RATE(MAX).....0.4LBS/HR (51CC/HR)
TEMPERATURE(MIN).....45F (7C)

STORAGE TANK - EARTH SAMPLE

EC004 11113

THIS TANK STORES SAMPLES OF EARTH GASES. THE ASSEMBLY CONSISTS OF A PRESSURE VESSEL, A POSITIVE EXPULSION REGULATOR AND A BLEEDER.

PRESSURE.....1 ATM (760 MM HG)
VOLUME.....10 CU.FT. (0.28 CU. METERS)

GAS SAMPLE STORAGE TANKS

EC004 51111

PRESSURE VESSELS FOR STORING EXPERIMENT GAS SAMPLES.

MAX PRESSURE(EST).....300 PSI (2.1E06 N/M2)
VOLUME.....1 CU.FT. (0.028 CU. METERS)

SUMP STORAGE TANK

EC004 11111

STORE GAS SAMPLES AFTER EXPERIMENT EVALUATION.

MAX PRESSURE (EST).....150PSI (1.05E06N/M2)
VOLUME.....10 CU.FT. (0.28 CU. METERS)

FLOW CONTROL ASSEMBLY

EC004 11124

THIS ASSEMBLY CONTROLS THE FLOW OF GASES IN THE CPL. IT IS MADE UP OF PLUMBING COMPONENTS, SOLENOID VALVES, A FILL ASSEMBLY, A VENT ASSEMBLY, SAFETY COMPONENTS AND A MIXING SUBASSEMBLY.

SUMP COMPRESSOR

EC004 11111

COMPRESS EXPANDED EXPERIMENT GAS SAMPLES INTO SUMP STORAGE TANK

COMPRESSION RATIO(MAX).....10 TO 1

DC-DC CONVERTER(VOLTAGE CONTROLLED DC)

EC005 11111

CONVERTS SPACELAB POWER TO VOLTAGE USABLE BY ELECTRICAL FIELD GENERATING COMPONENTS

INPLT VOLTAGE.....28V DC
OUTPLT VOLTAGE.....TBD

DC-DC CONVERTER (LOW FREQUENCY AC)

EC005 11111

CONVERT POWER INTO USABLE FORM FOR ELECTRIC FIELD GENERATING CAPACITOR

FREQUENCY.....0 TO 100KHz
WAVE SHAPE.....SQUARE
AMPLITUDE.....0 TO 3000 V/M

ELECTRIC FIELD CONTROLLER

EC005 11114

PROVIDE CONTROL FUNCTION FOR ELECTRIC FIELD GENERATOR. CONTROLLER SUPPORT BOTH AC AND DC POWER SUPPLIES

ACOUSTICAL GENERATOR

EC013 11111

VARIABLE FREQUENCY AUDIO DRIVERS USED TO POSITION AND CONTROL WATER DROPLETS OR ICE PARTICLES IN EACH OF THREE MUTUALLY PERPENDICULAR AXES. EACH DRIVER CAN BE CONTROLLED INDEPENDENTLY IN FREQUENCY AND AMPLITUDE

FREQUENCY RANGE.....0 TO 100KHz
AMPLITUDE.....0 TO TBD LB

TRANSDUCER - ACOUSTICAL

EC013 241111

MEASURE SOUND INTENSITY LEVEL

INTENSITY LEVEL.....0 TO TBD DB

AMPLIFIER - ACOUSTIC GENERATOR

EC013 11111

AMPLIFY SIGNAL TO ACOUSTIC GENERATOR

POWER.....45W

LASER

EC013 11111

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

POLARIZED, COLLIMATED LIGHT SOURCE USED FOR PARTICLE HEATING AND MOTION CONTROL AS WELL AS TO ASSESS OPTICAL SCATTERING PROPERTIES OF ICE CRYSTALS

LASER - POSITIONING SERVO CONTROL EC013 41114

PROVIDES FEEDBACK AND CONTROL SIGNAL TO LASER DRIVER FOR PROPER POSITIONING AND HEATING OF PARTICULATE SPECIMENS

LIQUID DROP GENERATOR E0004 41235

GENERATES SINGLE DROPLETS BY MANUAL OR SYSTEM CONTROL. DROPLETS ARE PROPELLED INTO THE EXPERIMENT CHAMBERS FOR EXAMINATION

FORMATION FREQUENCY(MAX).....10/SEC
DROPLET DIAMETER.....10 MICROMETERS TO 1 CM
DROPLET UNIFORMITY.....WITHIN 5 %

ELECTRONIC DRIVER ASSEMBLY EC006 41114

THE ELECTRONIC DRIVER ASSEMBLY CONTROLS THE LIQUID DROP GENERATOR. IT CONSISTS OF A GENERATOR POWER SUPPLY, A PARTICLE CHARGER POWER SUPPLY AND A DUAL PULSE GENERATOR.

ICE PARTICLE GENERATOR EC007 31235

GENERATE ICE PARTICLES FOR TEST IN ICE DIFFUSION, EXPANSION AND GENERAL CHAMBERS

ICE CRYSTAL SIZE.....50 MICROMETERS TO 1 CM

DROPLET CLOUD AEROSOL GENERATORS E0008 11135

GENERATE AEROSOLS FOR EXPERIMENTATION

PARTICLE CONCENTRATION.....10 UP TO 10000/CC
CLOUD FORM.....POLYDISPERSED AND MONODISPERSED
PARTICLE SIZE.....10E-4 TO .1 CM
PRODUCTION RATE.....UP TO 10E5 HZ

LARGE AND GIANT NUCLEI AEROSOL GENERATOR E0008 11135

SAME AS DROPLET CLOUD AEROSOL GENERATOR

ATKIN AEROSOL GENERATOR E0008 11135

SAME AS DROPLET CLOUD AEROSOL GENERATOR

DIGITAL COMPUTER EC010 11111

PROVIDES DATA PROCESSING AND CONTROL FUNCTIONS FOR THE CPL

TAPE RECORDER EC010 11111

RECORD DATA ON MAGNETIC TAPE

BIT RATE.....1000/S
STORAGE CAPACITY.....3.5E7BITS

SIGNAL CONDITIONING ELECTRONICS EC010 11111

ELECTRONICS SUPPORTING DATA MANAGEMENT SUBSYSTEM

OPTICAL PARTICLE COUNTER EC011 11235

THIS DEVICE WILL BE A SIZE ANALYZER. THE ANALYZER CONSISTS OF A PHOTOMULTIPLIER, AN AMPLIFICATION STAGE AND SUPPORTING ELECTRONICS

SENSITIVE PARAMETER.....LIGHT SCATTERED FROM PARTICLES
PARTICLE SIZE.....>0.3 MICROMETERS

NUCLEI MASS MONITOR SYSTEM EC011 11235

THIS DEVICE MEASURES THE TOTAL PARTICULATE MASS PER UNIT VOLUME. THE SYSTEM USES A CRYSTAL OSCILLATOR WHICH CHANGES ITS RESONANT FREQUENCY AS PARTICLES ARE DEPOSITED ON IT.

PARTICLE COLLECTION.....ELECTROSTATIC PRECIPITATION
PARTICLE SIZE.....0.01 TO 20 MICROMETERS

ELECTRICAL PARTICLE COUNTER EC011 11135

THIS ANALYZER MEASURES THE SIZE DISTRIBUTION OF PARTICLES. THE PARTICLES ARE IONIZED AND THEIR MOBILITY AS A FUNCTION OF ELECTRIC FIELD IS MEASURED TO GIVE AN INTEGRAL SIZE DISTRIBUTION

PARTICLE SIZE.....0.01 TO 1.0 MICROMETER
FLOW RATE.....1.8CFM (50LITERS/MIN)
INTEGRAL DISTRIBUTION CAPABILITY.....YES

ANALYZER - MULTICHANNEL EC011 11111

ALLOW SELECTIVE ANALYSIS OF COUNTER DATA

EXPANSION CLOUD CHAMBER EC014 11135

CHAMBER PROVIDING SIMULATION OF LONG TERM, NATURAL CLOUD DURING AN ADIABATIC EXPANSION

DIMENSIONS.....11.8IN.(30CM) DIA. X 17.7IN.
(45CM) LONG
VOLUME.....1.12CU.FT.(31.8L)
OPERATING TEMPERATURE.....-76F TO 104F (-60 TO 40C)
UPPER AND LOWER ENDS TEMP
DIFFERENTIAL.....0 TO 18F(0 TO 10C)
END TEMPERATURE TOLERANCE.....+/-0.1F (+/-0.05C)
OPERATING PRESSURE.....2.7 TO 14.7 PSIA(140 TO 760TORR)
PRESSURE RATE OF CHANGE.....19.4PSI/SEC (1000TORR/SEC)
PRESSURE RATE TOLERANCE.....+/-1%
MAXIMUM VOLUME EXPANSION.....V/V OF 0.5+/-0.1%
CIRCULAR ENDS CONTAIN ELECTRICALLY CONDUCTIVE PLATES TO SERVE AS EQUIPOTENTIAL SURFACES FOR THE ELECTRIC FIELD MOTION CONTROL SYSTEM

CONTINUOUS FLOW DIFFUSION CLOUD CHAMBER E0014 11135

SUPERSATURATION IS CONTROLLED BY THE TEMPERATURE OF WATER COVERING THE UPPER AND LOWER SURFACES OF THE CHAMBER. IT WILL BE USED FOR CLOUD CONDENSATION NUCLEATION EXPERIMENTS

DIMENSIONS.....11.8IN(30CM) X 11.8IN(30CM) X
2IN(5CM)
VOLUME.....0.16CU.FT.(4.5L)
INNER PLATE DIMENSIONS.....11.8IN(30CM) X 10IN(25CM) SPACED
0.5IN(1.3CM) APART
ISOTHERMAL TEMP CONTROL.....+/-0.1F(0.05C)
UPPER AND LOWER SURFACE
TEMPERATURE DIFFERENCE.....0 TO 18F (0 TO 10C)
SURFACE TEMP CONTROL.....+/-0.1F(0.05C)
OPERATING PRESSURE.....14.3 TO 14.7PSIA(740 TO 760TORR)
PRESSURE MEAS. ACCURACY.....+/-0.2%
MAX PRESSURE CHANGE RATE.....0.08PSI/SEC(0.4TORR/SEC)

STATIC DIFFUSION ICE CLOUD CHAMBER EC014 11135

THIS IS A NAKAYA TYPE CHAMBER WHICH UTILIZES ICE SURFACES TO PROVIDE CONTROLLED SUPERSATURATION RELATIVE TO ICE.

DIMENSIONS.....15.75IN(40CM) DIA X 7.9IN(20CM)
LONG
VOLUME.....0.446CU.FT.(12.5L)
OPERATING TEMPERATURE.....-40 TO 77F(-40 TO 25C)
TEMP MEAS. ACCURACY.....+/-1.8F(+/-1.0C)
UPPER AND LOWER SURFACE
TEMPERATURE DIFFERENCE.....UP TO 36F (20C)
SURFACE TEMP MEAS. ACCURACY.....+/-0.1F(+/-0.05C)
OPERATING PRESSURE.....1.5 TO 14.7PSIA(100 TO 760TORR)
PRESSURE MEAS. ACCURACY.....+/-0.19PSI(+/-10TORR)

STATIC DIFFUSION LIQUID CLOUD CHAMBER EC014 11135

A TWOMEY TYPE CHAMBER USED FOR EXPERIMENTS REQUIRING ABOVE FREEZING TEMPERATURES AND SUPERSATURATION OF THE LIQUID RELATIVE TO WATER. SUPERSATURATION IS CONTROLLED BY THE TEMPERATURES OF THE WATER COVERING CHAMBER UPPER AND LOWER SURFACES.

DIMENSIONS.....6IN.(15 CM) DIA X 0.6IN(1.5CM) LONG
VOLUME.....0.01CU.FT.(0.27L)
OPERATING TEMPERATURE.....32F TO 68F (0 TO 30C)
END SURFACE MAXIMUM
TEMPERATURE DIFFERENTIAL.....18F(10C)
SURFACE TEMP MEAS. ACCURACY.....0.1V(0.05C)
OPERATING PRESSURE.....2.84 TO 14.7PSIA(140 TO 760TORR)

GENERAL CLOUD CHAMBER EC014 1111

THIS CHAMBER WILL BE USED FOR MANY EXPERIMENTS THAT REQUIRE A RELATIVE HUMIDITY BELOW 100% AND MINIMUM TEMPERATURE CONTROL. PROVISIONS WILL BE MADE FOR GENERATING VARIOUS ELECTRIC FIELDS, POSITIONING DEVICES(SCAND, OPTICAL, ELECTRICAL), AND REMOTE DROPLET SIZING.

DIMENSIONS.....11.8IN(30CM) ALL SIDES (CUBE)
VOLUME.....0.95CU.FT.(27L)
OPERATING TEMPERATURE.....5C TO 86F(10 TO 30C)
OPERATING PRESSURE.....2.84 TO 14.7PSIA(140 TO 760TORR)
ELECTRIC FIELD PLATES LOCATED ON OPPOSITE SIDES
ACOUSTIC DRIVERS LOCATED ON THREE MUTUALLY PERPENDICULAR SIDES

AC POWER CONTROLLER EC017 11111

REGULATES AND CONTROLS AC POWER

REGULATION RANGE.....TBD

DC POWER CONTROLLER EC017 11111

REGULATES AND CONTROLS DC POWER

REGULATION RANGE.....TBD

MOTION PICTURE CAMERA (FAST FRAME) E0009 11111

FAST FRAME RATE CAMERA TO RECORD DROPLET COLLISIONS

FRAME RATE.....100/SEC
RECORDING DURATION.....20SEC
FIELD OF VIEW.....WIDE ANGLE

STILL CAMERA E0009 11111

RECORD CLOUD CHAMBER PHENOMENA

FRAME RATE.....2/SEC
SHUTTER SPEED.....1SEC TO 1/500SEC
CAPACITY.....500FRAMES
INTERCHANGABLE LENSES
FILM SIZE.....35MM

STILL CAMERA E0009 11111

RECORD CLOUD CHAMBER PHENOMENA

SHUTTER SPEED.....1SEC TO 1/500SEC
RESOLUTION.....2 MICROMETERS
DEPTH OF FIELD.....1MM
FRAME SIZE.....7CMH

MICROSCOPE EC017 11111

MAGNIFIED OBSERVATIONS

MAGNIFICATION FACTOR.....10 TO 100X
OBSERVATION MODE.....VISUAL OR PHOTOGRAPHIC

STEREO MICROSCOPE EC017 11111



MAGNIFIED OBSERVATIONS

MAGNIFICATION FACTOR.....10 TO 100X
OBSERVATION MODE.....VISUAL OR PHOTOGRAPHIC
WORKING DISTANCE.....UP TO 4 IN (10 CM)

PICKSCOPE (HIGH MAGNIFICATION) ECG17 11111

MAGNIFIED OBSERVATIONS

MAGNIFICATION FACTOR.....100 X 1000X
OBSERVATION MODE.....TRINOCULAR (CAMERA ATTACHMENT)

CSCILLOSCOPE ECG17 11111

SUPPORT MAINTENANCE ACTIVITIES AND SPECIAL DATA ACQUISITION

FREQUENCY RANGEDC TO 10.0 MHZ
OBSERVATION MODEVISUAL AND CAMERA (1 HR STORE)

TV MONITOR ECG12 11111

DISPLAY TV CAMERA OUTPUT SHOWING CHAMBER PHENOMENA

VICICCA ECG12 11131

RECORD CHAMBER PHENOMENA

RESOLUTION.....HIGH
LIGHT LEVEL.....LOW
IMAGE MAGNIFICATION.....YES
IMAGE INTENSIFICATION.....YES

IMAGE DEVICE LIGHTING ECG09 11114

PROVIDES LIGHTING FOR VARIOUS OPTICAL RECORDING ACTIVITIES.
CONSISTS OF CONTINUOUS LIGHTING, STROBSCOPE AND LASER BEAM LIGHT SOURCES

HIGH INTENSITY LIGHT DURATION.....<1 SEC

OPTICAL DETECTOR ECG05 11235

SOLID STATE DETECTOR TO OBTAIN DATA CONCERNING SCATTERING PROPERTIES OF ICE CRYSTALS AND POSSIBLY EXTREME DROPLET SIZE CHANGE DATA.

DROPLET SIZE DISTRIBUTION METER ECG05 11235

THIS UNIT WILL DETERMINE DROPLET SIZE DISTRIBUTIONS WITHIN A CHAMBER UTILIZING OPTICAL TECHNIQUES

IR MICROSCOPE ECG17 11111

DETERMINE SURFACE TEMPERATURE OF DROPLETS AND ICE CRYSTALS

SENSITIVITY.....+/-0.54F (+/-0.3C) AT 5F (-15C)

OPTICAL SENSORS ECG13 12111

SENSOR FOR CONTROL OF PARTICLE CHARACTERISTICS AND OPTICAL ENERGY SOURCES

DIGITAL TAPE RECORDER ECG06 11111

PROVIDES FOR RECORDING SENSOR DIGITAL OUTPUT DATA.

TAPE SIZE1 INCH
TRACKS9
TAPE SPEEDVARIABLE
DATA RATE1.3E+06 BPS
REEL SIZE12 INCHES DIA

LICAR TRANSMITTER ECG08 11112

PROVIDES SIGNAL SOURCE TO MEASURE CLOUD HEIGHTS AND AEROSOL DISTRIBUTION.

WAVE LENGTH6328 ANGSTROMS
ILLUMINATION MODESRUBY, DOUBLED RUBY, EYE,
CURLED EYE

RECEIVER, TELESCOPE ECG09 11112

RECEIVES TRANSMITTER RETURNED SIGNALS.

INSTANTANEOUS POV0.0018 RAD
POINTING ACCURACY0.0007 RAD
POINTING STABILITY0.0007 RAD
POINTING STABILITY RATE0.0017 RAD PER SEC
POINTING DURATION7200 SEC

DIGITAL RECORDER ECG91 11111

PROVIDES FOR RECORDING DIGITAL OUTPUT DATA.

TAPE SIZE1 INCH
TRACKS9
REEL SIZE12 INCHES
TAPE SPEEDVARIABLE
DATA RATE80 BPS

CSCILLOSCOPE ECG06 11111

PERMIT, MEASURE AND MAINTAIN ELECTRONIC EQUIPMENT OPERATION.

BANDWIDTH100 MHZ
CHANNELS2
SENSITIVITY10 MV/DIVISION
TIME BASEVARIABLE
STORAGE CAPABILITYNO

TELESCOPE - TRACKING MEC 11135

USED FOR HIGH RESOLUTION VIEWING OF SPECIFIC TARGETS

INSTANTANEOUS FIELD OF VIEW.....0.50 DEG (12X)
.....4.0 DEG (16X)
FIELD OF VIEW (CROSS-TRACK).....0.50 DEG AND 4.0 DEG
POINTING ANGLE FROM NADIR
(PITCH, ALONG TRACK).....+7C TO -40 DEG
POINTING ANGLE FROM NADIR
(ROLL, CROSS TRACK).....+/-75 DEG
TOTAL ANGULAR COVERAGE (CROSS-
TRACK) FROM NADIR (11FV/2)
POINTING ANGLE.....+/-75 DEG
POINTING ACCURACY (1-SIGMA).....0.10 DEG

VIEWER - WIDE ANGLE MEC 11235

USED FOR LARGE AREA VIEWING AND ORIENTATIONAL INSTRUMENT SIMILAR TO THE WILD H2 NAVIGATION SIGHT USED WITH THE WILD-HEERBRUGG RC-10 METRIC CAMERA

INSTANTANEOUS FIELD OF VIEW.....110 DEG X 110 DEG
.....55 DEG X 55 DEG
.....20 DEG X 20 DEG
FIELD OF VIEW.....30 DEG IN AZIMUTH
POINTING ANGLE FROM NADIR
(PITCH, ALONG TRACK).....0 TO 60 DEG
TOTAL ANGULAR COVERAGE (CROSS-
TRACK) FROM NADIR (11FV/2)
POINTING ANGLE.....+/-60 DEG
POINTING ACCURACY (1-SIGMA).....0.20 DEG

LONG WAVE INFRARED SPECTROMETER MEC 11111

USED FOR IDENTIFICATION OF TYPES OF ROCKS, SOILS AND SEDIMENTS

SPECTRAL RANGE.....0.4 TO 2.4 MICRONS
.....0.2 TO 15.5 MICRONS
RADIOMETRIC CHANNEL RANGE.....10 TO 12 MICRONS
INSTANTANEOUS FIELD OF VIEW.....1 MILLI RADIAN DIA.
POINTING ANGLE FROM NADIR
(PITCH, ALONG TRACK).....+45 DEG TO -100 DEG
POINTING ANGLE FROM NADIR
(ROLL, CROSS TRACK).....+/-20 DEG
POINTING ACCURACY.....0.30 DEG

MAGNETIC TAPE RECORDER MEC 21111

CANCARD DATA RECORDING OF SCIENTIFIC INSTRUMENT DATA

HEAD IN RATE.....100 MBPS
STORAGE CAPACITY/MISSION.....3E12 BITS

DATA BUFFERS, FORMATTERS MEC 21111

ELFFERS DATA FROM HIGH RESOLUTION WIDEBAND MULTISPECTRAL SCANNER

INPUT DATA RATE.....200 MBPS
OUTPUT DATA RATE.....100 MBPS
INPUT DATA DUTY CYCLE.....33%

MULTIPLEXER MEC 11111

MULTIFLEX DATA FOR TRANSMISSION VIA SHUTTLE COMMUNICATION LINK

CRT DISPLAY MEC 61111

TWO CRT'S ARE LOCATED IN THREE CONSOLES FOR DATA DISPLAY OF ALPHANUMERICS AND GRAPHICS.

RESOLUTION.....1000 LINES
CONTROL INPUT FORM.....DIGITAL
SCAN.....MASTER
.....CCNICAL
PHOSPHOR CHARACTERISTICS.....FULL-COLOR PENETRATION
SCREEN SIZE.....16 IN (40 CM)
CONTRAST.....8 GRAY LEVELS

MICROFILM READER MEC 11111

READER CONCURRENTLY DISPLAYS SEQUENCES OF MICROFILM RECORDS OF PAPS, REFERENCE SIGNATURES AND OTHER MATERIALS NEEDED FOR REAL TIME EVALUATION OF SENSOR FORMATS SELECTED FOR CRT DISPLAY

STORAGE CAPACITY.....73,000 PAGES
RETRIEVAL SPEED.....4 SECONDS MAX
INTERFACE.....STANDARD COMPUTER

MULTIFUNCTION KEYBOARD MEC 31111

CONSOLE CONTROL KEYBOARD WITH FUNCTION KEYS LABELLED BY COMPUTER. A SUPPORTING SOFTWARE LOGIC TREE FOR DATA ENTRY AND DISPLAY CALL UP CAN BE ENTER AT ANY LEVEL.

OPERATIONS CONTROL COMPUTER MEC 11111

THE COMPUTER SHALL PROVIDE THE FOLLOWING FUNCTIONS: COMMAND INTERPRETATION AND EXECUTION, ATTITUDE AND POINTING PROGRAM CONTROL AND COMPUTATION AND DATA MANAGEMENT SEQUENCING OPERATIONS

MEMORY SIZE.....14550 WORDS
ADD/MULT EXECUTE TIME.....5/30 MICRO SECONDS



Space Division
Rockwell International

EARTH AND OCEAN PHYSICS

FILM RECORDER

CP002 41121

PROVIDE PHOTOGRAPHIC RECORDINGS TO PROVIDE ALL WEATHER STEREO
IMAGERY FOR TOPOGRAPHIC MAPPING, 3-D STRAIN FIELDS, EROSION, VOL-
CANIC MOTION AND POST-GLACIAL UPLIFT.

FILM TYPE	70 MM	
FILM RATE	6 MM/SEC	
INSTANTANEOUS FOV	0.088 RAD	
MULTIFREQUENCY PROPAGATION RECEIVER/PROCESSOR & ANTENNA SYSTEM CONTROL UNIT	OP004	11134
MULTISPECTRAL SCANNER ELECTRONICS PACKAGE	CP004	11134
IN TEMPERATURE PROFILE RADIOMETER CONTROL UNIT	OP018	11134
CAMERA CONTROL	OP020	11134
GIMBAL CONTROL	OP020	11134
OPTICAL MONITOR CONTROL	OP021	11134
CONTAMINATION MONITOR GAUGE CONTROL	OP021	11134
MASS SPECTROMETER CONTROL	OP021	11134

ORIGINAL PAGE IS
OF POOR QUALITY



ORIGINAL PAGE IS
OF POOR QUALITY

SPACE PROCESSING APPLICATIONS

FLUID SUPPLY SYSTEM SPO01 11134

THIS UNIT CONSISTS OF THREE ELEMENTS SUPPORTING ALL MATERIAL SCIENCE AND MANUFACTURING EXPERIMENTS. IT IS PART OF THE CORE SUBELEMENT PRESENT ON ALL MS/MS PAYLOADS. THE THREE ELEMENTS ARE AN INERT GAS UNIT, AN OXIDIZING GAS UNIT, AND A REGCING GAS UNIT.

DIGITAL CLOCK SPO01 11111

DIGITAL TIME DISPLAY, PART OF DIGITAL PROCESS PROGRAMMER ASSY

DISPLAY 4 DIGITS
DISPLAY TYPE LED

SCANNER PROGRAMMER SPO01 11111

TRACKS COMPLEX FUNCTION AND TRANSLATES IT INTO A SIGNAL WHICH DRIVES EXPERIMENT CONTROLLER

ACCURACY..... ± 0.58 OF SPAN
REPEATABILITY..... ± 0.1 OF SPAN

SIGNAL CONDITIONER SPO01 11111

AMPLIFIES AND CONDITIONS LOW LEVEL ANALOG SIGNALS AS PART OF DATA ACQUISITION UNIT

CAPACITY.....UP TO 33 PAIRS OF SIGNAL WIRES
OUTPUT VOLTAGE..... $\pm 10V$
OUTPUT CURRENT..... $\pm 5mA$

DIGITAL VOLTMETER SPO01 11111

VOLTMETER FOR DATA ACQUISITION AND CONTROL UNIT

DC VOLT RANGE.....100MV TO 100V
AC VOLT RANGE.....100MV TO 10V
RESOLUTION.....0.1% FULL SCALE
ACCURACY..... $\pm 0.58 \pm 1$ DIGIT

SET POINT CONTROLLER SPO01 11111

ESTABLISHES CONTROLLER SET POINT. SET POINT CAN BE VARIED FROM PRE-PROGRAMMED INPUT

SIGNAL INPUT.....0 - 50MV
DATA OUTPUT.....1 - 50MA, 0 - 5 VDC

DIGITAL PROCESS PROGRAMMER SPO01 11111

PERFORM PREDETERMINED LOGIC ROUTINES SUCH AS EVENT SEQUENCING, PARAMETER COMMANDS OR RELAY TIME SEQUENCING.

INPUT/OUTPUT STAGE SPO01 11111

THE INPUT/OUTPUT UNIT, THROUGH THE USE OF PLUG IN CORDS, CAN ACCOMMODATE A WIDE RANGE OF PROCESS INPUTS AND OUTPUTS FOR MAXIMUM CONTROL VERSATILITY.

LOGIC LEVELS, LOGICAL "0".....CV
LOGICAL "1".....5V
SIGNAL INPUTS.....COMPATIBLE WITH PDP8 TYPE COMPUTERS
SIGNAL OUTPUTS.....0 - 5V
DATA OUTPUTS.....12 BIT PARALLEL

CONTROL UNIT - OPERATOR SPO01 11111

THE OPERATOR'S CONTROL PANEL WILL HAVE THE CAPABILITY OF CONTROLLING ABOUT 1000 PROCESS LOOPS. DATA MUST BE READOUT ON A QUICK-LOOK BASIS. IT MUST PROVIDE A WAY TO ENTER NEW PARAMETERS

CHANNEL NUMBER.....3 DIGIT NIXIE (2 EACH)
DATA LEVEL.....4 DIGIT NIXIE PLUS SIGNAL (2 EACH)
LOGIC LEVELS, LOGICAL "0".....CV
LOGICAL "1".....5V
DATA OUTPUT.....VARIABLE

PRINTER SPO01 11111

PRINT RECORDED DATA IN THE FORM OF ALPHA-NUMERIC CHARACTERS IN HARD COPY FORM UPON REQUEST BY OPERATOR.

FORMAT PICA TYPE - STANDARD CHARACTERS, NUMBERS SYMBOLS

TELEPRINTER SPO01 11111

CHARACTER ORIENTED CRT TERMINAL. DATA IN THE LINE MEMORY IS CONVERTED, BY THE CHARACTER GENERATOR, INTO THE APPROPRIATE CGT PATTERN. THE DOTS ARE SHIFTED OUT OF THE SHIFT REGISTER TO FORM A VIDEO SIGNAL TO THE MONITOR.

SCREEN FORMATS.....40, 72 OR 80 CHARACTERS PER LINE
.....12 OR 24 LINES
TRANSFER RATE.....110 TO 2400 BAUD, 10 OR 11 BIT CHARACTERS
MODES.....HALF OR FULL DUPLEX

DIGITAL TAPE RECORDER SPO01 11111

PROVIDE DATA STORAGE FOR DIGITAL PROCESS PROGRAMMER

CHANNELS.....15

ANALOG (SER) CONTROLLER SPO01 11111

SILICON CONTROL RECTIFIER EMPLOYED AS A RELAY SWITCH

LOGIC LEVELS, LOGICAL "0".....0 - 0.4V
LOGICAL "1".....3.0 - 5.0V
RANGE OF CONTROL.....0 - 100% IN 1 CYCLE INCREMENTS

MULTIPLEXER A/D CONVERTER

SPO01 11111

CONVERT INPUT SIGNAL FROM ANALOG TO DIGITAL FORMAT FOR I/O UNIT

ANALOG INPUT..... -10.24 TO $+10.235V$
FULL SCALE ACCURACY..... ± 0.025
RESOLUTION.....12 BITS BINARY
CONVERSION SPEED.....8 KHZ
DATA OUTPUT.....2'S COMPLEMENT

TAPE INPUT UNIT

SPO01 11111

RECEIVES, STORES, TRANSMITS DATA OF OTHER EQUIPMENT

TAPE CAPACITY.....12 COFT (36CM), 0.5IN (1.2CM) WIDE
TAPE SPEED.....25, 18.75 AND 12.5 IN/S
.....164.46 AND 32 CM/S
LONG TERM VARIATION..... $\pm 1\%$
DATA DENSITY.....1600 CPI (PHASE ENCODED)
.....200,000,000 CPI (NRZ)

STORAGE PERIPHERALS

SPO01 11111

THIS UNIT ACCEPTS ADDRESSED BINARY DATA FROM THE DIGITAL PROCESSOR AND TRANSLATES THIS DATA INTO EQUIVALENT ANALOG FORM.

MAXIMUM SIGNAL AND HOLD CHANNELS: 32
SIGNAL AND HOLD CHANNELS PER CARD: 4

ALTERNATE PHOTOGRAPHIC PROCESSOR

SPO01 11111

UNIT PROVIDES PROCESSED MICROGRAMS AS PART OF THE ELECTRO-OPTICAL IMAGING SYSTEM

CYCLE TIME.....BETWEEN 5 TO 20 SEC
PLATE HOLDER SIZE.....4 IN (10CM) X 5 IN (13CM)
EFFECTIVE APERTURE AT PHOTO PLATE.....3.25 IN (8.5CM) X 4.0 IN (10CM)
DATA OUTPUT.....PROCESSOR MICROGRAMS

CCTV CAMERA

SPO01 11131

HIGH RESOLUTION TELEVISION CAMERA FOR VIEWING, RECORDING OR TRANSMITTING PICTURES OF MATERIAL PROCESSING EXPERIMENTS

BANDWIDTH.....32MHZ
VERTICAL SCAN LINES.....1125 LINES
HORIZONTAL RESOLUTION.....1100 LINES
VERTICAL SWEEP RATE.....25 OR 30 FRAMES PER SECOND
SENSITIVITY.....100 IRE UNITS OUTPUT PRODUCED WITH 0.5 FOOTCANDLE ILLUMINATION

CCTV CAMERA CONTROL UNIT

SPO01 11131

THIS UNIT WILL PERFORM THE FOLLOWING CONTROL FUNCTIONS: LENS FOCUS, LENS F STOP, CONTRAST BRIGHTNESS AND MODE. THE UNIT WILL WORK IN TWO MODES: (1) CONTINUOUS OBSERVATION MODE, WHERE ILLUMINATION IS CONTINUOUS, OR (2) PULSED OBSERVATION MODE USED WHEN ILLUMINATION IS PROVIDED BY A GYE LASER OR A NAKED FLASH LAMP.

SCANNING.....INTERLACED 241, THROUGH 1125 LINES PER FRAME, 25 OR 30 FRAMES PER SECOND
HORIZONTAL SCAN FREQUENCY CONTROL.....36 TO 750 KHZ
VERTICAL SCAN FREQUENCY CONTROL.....50 TO 60HZ

FRAME STORAGE UNIT

SPO01 11111

THE FRAME STORAGE UNIT PROVIDES NONDESTRUCTIVE READOUT, ERASE, AND STORAGE ON ELECTRONIC COMMAND AND WORKS IN EITHER CONTINUOUS OR PULSED MODE. IT PROVIDES A REAL-TIME, NC PHOTOGRAPHIC PROCESSING, HIGH SPEED, AUTOMATIC LIGHT ADJUSTMENT CAMERA CAPABLE OF PROVIDING REAL-TIME MOVIES OR FROZEN STILL

TV MONITOR

SPO01 11111

THE MONITOR WILL DISPLAY CCTV IMAGES AS PART OF THE ELECTRO-OPTICAL IMAGING SYSTEM.

SCANNING FREQUENCIES.....HORIZONTAL: 15-40 KHZ
.....VERTICAL: 15-60 FIELDS/SECOND
VERTICAL SCAN LINES.....1125 LINES
LINEARITY.....1% PICTURE HEIGHT

CSCILLESCOPE

SPO01 11111

DISPLAY OF THE PRESENCE AND/OR NATURE AND FORM OF OSCILLATIONS OR IRREGULARITIES OF AN ELECTRIC CURRENT. STORAGE CAPABILITY OF A CRT DISPLAY IS DESIRED FOR PROLONGED OBSERVATIONS.

BANDWIDTH.....DC TO 10MHZ
SENSITIVITY.....INTERNAL: 0.2CM DEFLECTION TO 1 MHZ
.....EXTERNAL: 250MV P-P TO 15V P-P
VIDEO STORAGE TIME.....1 HOUR

FLUID COOLING/REFRIGERATION UNIT

SPO01 11114

THIS UNIT MAINTAINS THE ELECTROPHORETIC SEPARATION SYSTEM, THE BUFFER SUPPLY AND SAMPLE AT A CONSTANT LOW TEMPERATURE. IT CONTROLS THE TEMPERATURE OF THE GAS ELIMINATION SYSTEM, ALSO IT COOLS (OR HEATS) THE BUFFER/SAMPLE SOLUTIONS AND THE BUFFER USED IN THE ELECTRODE COMPARTMENTS OF THE COLUMNS.

CONTROLLED TEMPERATURE RANGE..... $-13F$ TO $113F$ (25 TO 45C)
TEMPERATURE ACCURACY..... $\pm 1.0F$ ($\pm 0.56C$)
SYSTEM RECOVERY TIME.....2 MINUTES
MAXIMUM HEAT REJECTION.....2KW THERMAL

LASER OPTICAL SCATTERING MONITOR

SPO01 21111

THE LASER OPTICAL SCATTERING MONITOR WILL PROVIDE A REAL TIME MEASURE OF LIGHT SCATTERING IN A SPECIMEN AT ANY DESIRED MEASUR-

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

ING ANGLE. THIS DATA WILL PROVIDE INFORMATION ON THE SIZE, SHAPE, ORIENTATION, INDEX OF REFRACTION, CONCENTRATION AND LOCATION OF OPTICAL SCATTERING CENTERS.

LASER TYPE.....LON-POWER UV GAS DISCHARGE
FILTER BANDWIDTH.....4E-9 IN. (0.1NM) AT HALF POWER
SCATTERED LIGHT INTENSITY
ACCURACY.....0.18
ANGULAR POSITION ACCURACY.....1 MIN (0.3MRAD)

UV-VIS SPECTROMETER SPOC1 21111

MEASUREMENT OF SAMPLE CONSTITUENTS BASED ON THE UV AND VISIBLE SPECTRUM.

RANGE150 TO 1000 NM
WAVELENGTH/WAVENUMBER ACCURACY..... ± 0.1 NM FS
RESOLUTION0.05 NM
RESPONSE TIME2 SECONDS
SCANNING SPEED0.02 TO 2.0 NM/SEC

EYE LASER/FLASH LAMP SPOC1 31111

THIS UNIT WILL BE USED FOR FREE RADICAL GENERATION, SURFACE DAMAGE THRESHOLD DETERMINATION, AND HOLOGRAPHIC MICROSCOPY

ACCURACY:
WAVELENGTH.....0.018
VOLTAGE OF CAPACITOR BANK.....0.18
LASER AND FLASH LAMP ENERGY
OUTPUT.....1.08
SYNC. ACCURACY.....10CMS
COHERENCE.....208
RANGE OF OPERATION:
ENERGY INPUT.....500-5000J/PULSE (1.14-1.4NM/PULSE)
1 AT LOW REPS, LARGE LAMP
50-500J/PULSE (0.14-1.4NM/PULSE)
AT HIGH REPS, SMALL LAMP
ENERGY OUTPUT.....1-30J/PULSE (3E-4 - 8E-3NM/PULSE)
1 AT LOW REPS, DYE FUND. FREQ.
250MJ TO 7J/PULSE (7E-5 - 2E-3
NM/PULSE) FREQUENCY DOUBLED.
100MJ TO 3J/PULSE (3E-5 TO 8E-4
NM/PULSE) AT HIGH REPS AT DYE
FUNDAMENTAL, 25 - 750MJ (7E-6 -
2E-4 NM) FREQUENCY DOUBLED
REPETITION RATE.....LARGE LAMP - SINGLE TO 10/MIN
SMALL LAMP - SINGLE TO 20HZ
FREQUENCY RANGE.....FUNDAMENTAL 430 - 750NM (1.7E-5
TO 3E-5 IN.) FREQUENCY DOUBLED
215 - 375NM (0.9E-5 TO 1.5E-5
IN.)
FLASH DURATION.....500 - 100NS
HOLOGRAM.....35NM (1.4 IN.)
SINGLE FRAME TO 20FRAMES/SEC.

FLASH LAMP SPOC1 31111

THE FLASH LAMP HAS TWO FUNCTIONAL REQUIREMENTS: (A) TO PRODUCE FREE RADICALS, AND (B) AS A PUMP SOURCE FOR EXCITING THE DYE LASER

TYPE.....MERCURY-XENON
OPERATING TEMPERATURE.....68 TO 77F (20 TO 25C)
WAVELENGTH.....1.2E-5 TO 2.4E-5 IN. (300-600NM)
ENERGY INPUT(MAX).....1.5NM (5KJ)
(MIN).....0.3NM (1KJ)

RETROCONSTRUCTION HIGH RESOLUTION HOLOGRAPHIC SCOPE SPOC1 31235

DIALYSIS UNIT SPOC1 11113

METERING PUMPS SPOC1 11111

BUFFER AND SPECIMEN METERING PUMPS FOR CONTINUOUS FLOW ELECTROPHORETIC COLUMN

FLOW RATE.....0.07ML/MIN TO 30ML/MIN
PRESSURE HEAD.....4E04 - 5.0 TO 14 PSI (4E04 - 10E04NM/2)
REPRODUCIBILITY..... ± 0.5

CONTINUOUS FLOW ELECTROPHORETIC COLUMN WITH PUMP SPOC1 11114

THIS UNIT SEPARATES BIOLOGICAL SAMPLES BY ELECTROPHORESIS. RELATIVELY LARGE SAMPLES, APPROXIMATELY 10 CC, WILL BE PROCESSED. SAMPLE AND BUFFER FLOW WILL BE CONSTANT AND SAMPLE FRACTION COLLECTION CONTINUOUS.

TOTAL BUFFER AND SAMPLE VOLUME.....6000 IN. (1000CC)
OPERATING TEMPERATURE.....14 TO 41 F (-10 TO 5C)
OPERATING TEMPERATURE REGULATION
SEPARATION PROFILE ACCURACY..... ± 0.04 IN (± 1 MM)
POSITION MONITOR ACCURACY..... ± 0.04 IN (± 1 MM)
FRACTION OUTLET PORT LOCATION ACCURACY..... ± 0.04 IN (± 1 MM)
INPUT VOLTAGE.....500V

STATIONARY ELECTROPHORETIC COLUMNS SPOC1 11114

THE STATIONARY ELECTROPHORETIC COLUMNS HAVE TWO ANTICIPATED USAGES: (A) TO ESTABLISH PROCESSING AND OPERATING PARAMETERS FOR THE CONTINUOUS ELECTROPHORETIC COLUMNS, AND (B) FOR SMALL BATCH PROCESSING (APPROXIMATELY 0.1CC (0.04 IN³)) OF VARIOUS BIOLOGICAL COMPONENTS

TOTAL BUFFER AND SAMPLE VOLUME.....1.500 IN. (25CC)
SEPARATION VELOCITY ACCURACY..... ± 0.04 IN (± 1 MM)
POSITION MONITOR ACCURACY..... ± 0.04 IN (± 1 MM)
OPERATING TEMPERATURE RANGE.....14 TO 41 F (-10 TO 5C)
TEMPERATURE TOLERANCE..... ± 2 F (± 1 C)

GAS ELIMINATION SYSTEM SPOC1 11114

REMOVES GASSES PRODUCED BY ELECTROLYSIS IN THE ELECTRODE COMPARTMENTS OF THE ELECTROPHORETIC COLUMNS

OUTLET OXYGEN CONCENTRATION... 31 (BY VOLUME)
OUTLET HYDROGEN CONCENTRATION... 21 (BY VOLUME)
MAX FLOW RATE.....180
FLOW RATE CONSTANCY..... ± 1
OPERATING TEMPERATURE.....68 - 77 F (20 - 25 C)

BUFFER AND ELECTROLYTE SUPPLY TANKS SPOC1 11112

PH MONITOR SPOC1 21111

MEASURE ACIDITY AND ALKALINITY OF SOLUTIONS. USED IN GP LAB ALSO

PH RANGE.....4.2 TO 12
SELECTABLE RANGES.....ANY 2
ANY 5
ANY 10
READING STABILITY.....0.02 PH
LIQUID TEMPERATURE RANGE.....32 TO 212 F (0 TO 100 C)

FRACTION COLLECTION SYSTEM SPOC1 11235

THIS APPARATUS WILL HAVE TWO FUNCTIONS: (1) TO COLLECT THE DESIRED FRACTIONS AFTER SEPARATION, AND (2) TO REMOVE THE EXCESS BUFFER SOLUTION AND REMAINING UNWANTED SAMPLES

PUMP TYPE.....POSITIVE DISPLACEMENT
1.8 CU. IN/MIN (30 CC/MIN)
PUMPING RATE ACCURACY..... ± 0.1
ALLOWABLE PUMP PULSATION.....0.01PSI (± 7 C N/M2)
TEMPERATURE REGULATION..... ± 1 OF (± 0.5 C)

FLOW METER SPOC1 11111

GAS FLOW METERING AND MONITORING DEVICE.

GAS TYPEOXYGEN, NITROGEN, HELIUM
GAS FLOW RATE0.4 TO 0.3 LB/MR
0.3 TO 3.0 LB/MR

LYOPHILIZATION UNIT SPOC1 11111

RECIRCULATING FLUID INCUBATOR SPOC1 11112

DISSOLVED OXYGEN ANALYZER SPOC1 11111

THE UNIT WILL DETERMINE THE AMOUNT OF DISSOLVED OXYGEN PRESENT IN THE ELECTROLYTE AFTER IT HAS BEEN PROCESSED THROUGH THE GAS ELIMINATION SYSTEM. IT WILL MONITOR OXYGEN DURING A REACTION EXPERIMENT. IT WILL ALSO FUNCTION AS A GASEOUS OXYGEN MEASUR

CONCENTRATION RANGE.....0 TO TOTAL OXYGEN SATURATION
ACCURACY..... ± 1
RESPONSE TIME.....501 ACTUAL IN 10SECONDS
SAMPLE TEMPERATURE RANGE.....32 TO 110F (0 TO 43C)

SPECIMEN/SAMPLE SUPPLY TANKS SPOC1 11113

HIGH VOLTAGE POWER CONDIGNER SPOC1 11112

PROVIDES DC VOLTAGE TO EITHER THE STATIONARY OR CONTINUOUS FLOW ELECTROPHORETIC COLUMNS. REGULATION MUST REMAIN CONSTANT OVER A RELATIVELY LARGE CURRENT DEMAND.

OUTPUT VOLTAGE REGULATION..... ± 0.1
VOLTAGE RANGE.....25 TO 2500 VDC
HEADOUTS.....VOLTAGE LEVEL
POWER OUTPUT TO COLUMNS.....20W SUSTAINED
200W PEAK

DARK FIELD ILLUMINATOR SPOC1 21111

REFRIGERATOR SPOC1 11111

THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES

TEMPERATURE.....32 - 68F (0 - 20C)
LENGTH.....2FT (60CM)
WIDTH.....1.8FT (54CM)
HEIGHT.....2.8FT (84CM)

FREEZER SPOC1 11111

THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES

CONTROLLED TEMPERATURE RANGE.....-13F TO 113F (-25C TO 45C)
LENGTH.....1FT (30CM)
WIDTH.....1.8FT (54CM)
HEIGHT.....2.8FT (84CM)

GENAR SPOC1 11111

THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES

TEMPERATURE.....-112F (-80C)
LENGTH.....1FT (30CM)
WIDTH.....1.8FT (54CM)
HEIGHT.....2.8FT (84CM)

WASTE LIQUID TANK SPOC1 11113

MOLECULAR SIEVE SPOC1 11113

THE MOLECULAR SIEVE IS USED TO CLEAN GASSES GENERATED IN FURNACES AND EXHAUSTED TO SPACE. THIS UNIT IS ALSO USED IN THE FURNACE, GENERAL PURPOSE AND LIFTATION SUBELEMENTS

CECLANT SUPPLY TANK SPOC1 11113

ISOELECTRIC FOCUSING UNIT SPOC1 11114

HOT WALL TUBE FURNACE SPOC2 31111

THIS UNIT IS A GENERAL PURPOSE HOT WALL HEATING DEVICE



PROVIDING ACCURATE CONTROL OVER THE HOT ZONE TEMPERATURE WITH RESPECT TO A FLAT PROFILE OR TO A SPECIFIED GRADIENT PROFILE.

VACUUM LEVEL.....2E-8 PSI (1.5E-4N/M2)
OPERATING TEMPERATURE.....2200F (1200C)
MAX EXTERNAL SURFACE TEMP.....110F (45C)
HOT ZONE DIMENSIONS.....2.5 IN (6CM) DIA. X 6IN (15CM) LG
MAXIMUM DEVIATION FROM FLAT PROFILE.....+/-4F (+/-2C)

HOT WALL FURNACE SP002 21111

THIS UNIT IS A GENERAL PURPOSE, HOT WALL HEATING DEVICE PROVIDING ACCURATE CONTROL OVER THE HOT ZONE TEMPERATURE WITH RESPECT TO A FLAT PROFILE OR TO A SPECIFIED GRADIENT PROFILE.

VACUUM LEVEL.....2E-8 PSI (1.5E-4N/M2)
OPERATING TEMPERATURE.....3300F (1800C)
EXTERNAL SURFACE TEMPERATURE.....110F (45C)
HOT ZONE DIMENSIONS(MIN).....1IN. (2.5CM) DIA. X 5IN. (12.5CM) LG
MAXIMUM DEVIATION FROM FLAT PROFILE.....+/-7F (+/-4C)
VIEWPOINTS.....2

CHEST-GENERAL PURPOSE ENCLOSURE SP002 31113

THE PRIMARY FUNCTION OF THE ENCLOSURE IS TO ISOLATE THE HEATED SPECIMEN FROM THE LABORATORY ENVIRONMENT, BOTH THERMALLY AND ATMOSPHERICALLY, AND TO PROVIDE A MEANS TO STRUCTURALLY ATTACH THE VARIOUS HEATING ELEMENTS AND ACCESSORY ITEMS

VACUUM LEVEL.....2E-10PSI (1.5E-6N/M2)
MAXIMUM INTERNAL TEMP.....5400F (3000C)
MAX EXTERNAL SURFACE TEMP.....110F (45C)
HOT ZONE DIMENSIONS(MIN).....1IN. (2.5CM) DIA. X 4IN (10CM) LG
FULL LENGTH FRONT OPENING DOOR
VIEW POINTS.....FRONT, REAR, 2 SIDES, BOTTOM, TOP
FLANGE OPENINGS.....TOP, BOTTOM, 2 EACH SIDE, 2 REAR

DIRECTIONAL SOLIDIFICATION UNIT SP002 11114

THIS UNIT CONSISTS OF A FURNACE AND A CHILL UNIT. SPECIMENS ARE HEATED, MOVED THROUGH A THERMAL GRADIENT AREA AND THEN INTO A COOLING AREA.

CHILL UNIT DIMENSIONS.....0.6 IN (1.5CM) I.D. X 8.0 IN (20CM) LG
CHILL CAPACITY.....1KM THERMAL
HEATING UNIT DIMENSIONS.....0.8 IN (2CM) I.D. X 8IN (20CM) LG
HEATER CAPACITY.....1KM THERMAL
TEMPERATURE GRADIENT ADJUST RANGE.....9 TO 900F/IN (2 TO 200C/CM)
HEATER MAX TEMPERATURE.....2900F (1600C)
CHILL UNIT MIN TEMP.....68 TO 77F (20 TO 25C)
EXTERNAL SURFACE TEMP.....60 TO 110F (15 TO 50C)

ACOUSTIC MIXING AND DISPERSAL UNIT SP002 21235

ELECTROMAGNETIC MIXING AND DISPERSAL UNIT SP002 21235

IR PYROMETER SP002 21111

THE PYROMETER OPTICALLY SENSES THE INTENSITY OF RADIANT ENERGY EMITTED FROM THE SAMPLE SURFACE

TEMPERATURE RANGE104 TO 5400 F (40 TO 3000 C)
SPECTRAL RESPONSE.....0.70 - 0.97 MICROMETERS
CALIBRATION ACCURACY.....+/-1% FULL SCALE
REPEATABILITY.....+/-3% FULL SCALE

ZONE REFINER SP002 11112

THIS UNIT WILL HEAT, MELT AND MANIPULATE A SPECIMEN FOR ZONE REFINING AND MELTEN ZONE CRYSTAL GROWTH. IT WILL CAUSE A MELTEN REGION TO TRAVERSE ALONG A LONG CYLINDRICAL ROD OF MATERIAL.

SPECIMEN SIZE.....4 TO 12 IN (10 TO 30 CM) LG
0.1 TO 0.5IN (0.25 TO 1.25CM) DIA
SCANNING DISTANCES.....UP TO 12IN (30CM)
SCANNING SPEED.....UP TO 20IN/MIN (50CM/MIN)
SPECIMEN ROTATION.....0 TO 60 RPM
OPERATING TEMPERATURE.....575 TO 2900F (300 TO 1600C)

FEED AND CRYSTAL HOLDER SP002 21113

RESIDUAL GAS ANALYZER SP002 21111

IDENTIFY TYPE AND QUANTITY OF GAS REMAINING IN A PROCESSING CHAMBER OR GENERATED DURING A PROCESS. THE ANALYZER IS ALSO USED IN THE LEVITATION SUBELEMENT. THE UNIT MUST HAVE A TUNABLE MASS UNIT SECTION FOR LEAK DETECTION, BE A QUADROUPLE TYPE INSTRUMENT CAPABLE OF SEPARATING MASSES. THE INSTRUMENT MUST ACT AS A PARTIAL AND TOTAL PRESSURE ANALYZER AS WELL AS GAS ANALYZER.

MASS RANGE RESOLUTION.....UP TO 200AMU
DETAILED SPECTRAL ANAL MASS RANGES.....1-50 AND 40-65AMU
SCOPE SCAN SPEED.....VARIABLE 3MS TO 300S
DATA OUTPUT.....VISUAL ANALOG
ANALOG
DIGITAL

USEFUL PRESSURE RANGE.....1.0E-5 PSI (7.5E-2N/M2)
TO 1.0E-11 PSI (7.5E-8N/M2)
MASS SEPARATION CAPABILITY.....1AMU
MIN DETECTABLE PARTIAL PRESSURE.....2E-12 PSI (1.5E-8N/M2)
MIN DETECTABLE TOTAL PRESSURE.....1E-10 PSI (7.5E-7N/M2)

RF INDUCTION PWR COND. (MIXING & DISPERSAL) SP002 21111

LOW VOLT/HIGH AMP POWER CONDITIONER (10KW) SP002 31111

TWO - COLOR PYROMETER SP002 21111

RESISTANCE HEATER (CONTACT) SP002 21111

HEATING ELEMENT WHICH HEATS ENCAPSULATED SAMPLES BY PASSING A CURRENT THROUGH IT

POWER.....N/A SUSTAINED, 9KW PEAK
TEMPERATURE(MAX).....2900F (1600C)
ATMOSPHERE.....INERT GAS OR VACUUM

MICROWAVE HEATER SP002 31235

PROVIDE MICROWAVE ENERGY SOURCE FOR SPECIMEN HEATING

TEMPERATURE(MAX).....1800C (16,000F)
POWER(MAX).....2KW
FREQUENCY RANGE(MAX).....CF TO 1E10 HZ
OPERATIONAL ATMOSPHERE.....ANY BUT VACUUM

GRADIENT FURNACE SP002 11113

MECHANICAL MIXING AND DISPERSAL UNIT SP002 11235

THREE AXIS MANIPULATOR SP002 21132

PIEZOELECTRIC DRIVE SP002 21114

VACUUM/PRESSURE REGULATOR SP002 31111

THE UNIT REGULATES THE PRESSURE IN THE PROCESSING CHAMBERS BY ADMITTING FLUID FROM A FLUID SUPPLY CONTAINER. THIS UNIT IS ALSO EMPLOYED ON THE GENERAL PURPOSE AND LEVITATION SUBELEMENTS

PRESSURE-ABSOLUTE0.1 TO 1000 PSIA

HIGH VACUUM PUMP SP002 21111

UNIT IS USED IN CONJUNCTION WITH MOLECULAR SIEVE. IT IS USED WITH THE SPACE VACUUM TO BLEED DOWN EVENS. THIS UNIT IS ALSO USED ON THE LEVITATION SUBELEMENT.

PRESSURE-ABSOLUTEAMBIENT TO 0.1 PSIA

VACUUM/PRESSURE MEASUREMENT UNIT SP002 31111

MEASURES PRESSURES IN PROCESS CHAMBERS. UNIT IS ALSO REQUIRED FOR GENERAL PURPOSE AND LEVITATION SUBELEMENTS.

PRESSURE-ABSOLUTEAMBIENT TO 0.01 PSIA

HIGH VOLTAGE POWER CONDITICNER (17KV) SP002 31111

RF INDUCTION PWR CONDITICNER (2KW - 2MHZ) SP002 21111

THERMAL CHEN SP003 21111

IR SPECTROMETER SP003 11111

INFRARED ANALYSIS OF CONSTITUENTS IN GASEOUS OR LIQUID STREAMS.

RANGE10 TO 100 MICRONS
RESOLUTION+/- 0.5 MICRONS
RESPONSE TIME2 SECONDS

LASER PYROMETER SP003 21235

CAS CHROMATOGRAPH SP003 21111

THE FUNCTIONAL REQUIREMENTS OF THE EQUIPMENT ARE TO ANALYZE THE ATMOSPHERIC COMPOSITION AND THE IMPURITY CONCENTRATION CONTAINED WITHIN THE EXPERIMENTS' EQUIPMENT ENCLOSURES. THIS ELEMENT IS ALSO USED ON THE LEVITATION SUBELEMENT.

FLUX RATE10 PL/MINUTE
SENSITIVITY.....0.1 PPM

TIME LAPSE HIGH SPEED CAMERA SP001 21111

CAMERA ALSO USED IN LEVITATION SUBELEMENT

NUCLEAR PARTICLE COUNTING UNIT SP001 11111

CAPABILITY TO DISCRIMINATE AND COUNT ALPHA, BETA, GAMMA, X-RAYS AND K-ELECTRONS.

COUNTING RATE3.5E+06 C/M
RESOLUTION LOSS1 1/3E+06 C/M
PLATEAU CHARACTER
ALPHA500-1200 V @ 1 1/3ICCV
BETA-GAMMA1200-1900 V @ 1 1/3ICCV
PRECISION+/- 0.1 % FS
TIMER0.05 TO 2160 MINUTES

RESISTANCE THERMOMETER SP003 11111

THERMOCOUPLES SP003 11111

RESISTANCE HEATER (NON-CONTACT) SP004 11111

TUBE TYPE HEATING UNIT FOR USE IN CHEST-GENERAL PURPOSE ENCLOSURE OR HOT WALL FURNACE

POWER.....4KW SUSTAINED, 9KW PEAK
TEMPERATURE(MAX).....2900F (1600C)
SPECIMEN SIZE.....4IN (10CM) LG X 0.75 IN (2CM) DIA.

RF INDUCTION COILS SP004 11111

PROVIDE NON-CONTACT HEATING

FREQUENCY RANGE.....2MHZ TO 2MHZ
OPERATIONAL ATMOSPHERE.....ANY

ELECTRON BEAM SOURCE SP004 11114



THE ELECTRON BEAM SOURCE WILL BE USED TO HEAT SAMPLES THAT REQUIRE NON-CONTACT HEATING. THE UNIT WILL BE USED IN THE CHEST-GENERAL PURPOSE ENCLOSURE

DELIVERABLE POWER.....3 - 4 KW
INPUT POWER.....6KW
TEMPERATURE.....4500F (2500C)
MAXIMUM SAMPLE SIZE.....0.15 CU. IN. (2CC)
OPERATING PRESSURE.....1E-10 PSI (7.5E-7N/M2)

LASER SOURCE SPOC4 11111

PROVIDE HEAT TRANSMITTED TO A SPECIMEN BY LASER BEAM. CONCEPT PROVIDES HEATING WITHOUT CONTAMINATING SPECIMEN. UNIT IS USED IN CONJUNCTION WITH CHEST-GENERAL PURPOSE ENCLOSURE

POWER.....6KW
COOLING CAPACITY.....7.5KW
LASER TYPE.....ANY OF FOLLOWING FIVE:
.....CONTINUOUS WAVE HYDROGEN-FLOUR-
.....INE CHEMICAL LASER, IR CLPUT
.....CONTINUOUS WAVE FORCED CONVEC-
.....TION CO2, IR OUTPUT
.....CONTINUOUS WAVE CO-E2-H2C-HE
.....CHEMICAL LASER, IR OUTPUT
.....PULSED TEA CO2, IR OUTPUT
.....PULSED TEA XENON, UV OUTPUT

MINIMUM B WITH RF HEATING SPOC4 11235

ELECTROMAGNETIC POSITIONING COILS AND DETECTOR SPOC4 11115

PROVIDE ALTERNATING ELECTROMAGNETIC FIELD FOR CONTACTLESS POSITION CONTROL

MAGNETIC PRESSURE 1.4 PSI (10 KN/SC M)

ELECTROSTATIC POSITIONING PROBES AND DETECTOR SPOC4 11115

ELECTROSTATIC FIELDS USED FOR POSITIONING EXPERIMENT SPECIMEN SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL

CHAMBER PRESSURE > 0.14 PSI (1 KN/SC M)
..... < 2E-06 PSI (0.01 N/SC M)

CONFIGURATION IN APEX OF TETRAHEDRON

ACOUSTIC TRANSDUCER AND DETECTOR SPOC4 11111

USE OF SOUND PRESSURE WAVES TO POSITION EXPERIMENT SPECIMEN SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL

CHAMBER PRESSURE > 0.14 PSI
..... (1 KN/SC M)
..... < 2E-06 PSI
..... (0.01 N/SC M)

CONFIGURATION IN APEX OF TETRAHEDRON

GAS JET POSITIONING PROBES AND DETECTORS SPOC4 11115

USE OF GAS JETS TO POSITION SPECIMEN DURING EXPERIMENT PRO-
CEDURES SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL

GAS FILTER 8E-06 INCHES (0.2 MICROPETERS)
GASES HELIUM, HYDROGEN, OXYGEN
CONFIGURATION JET PAIR IN APEX OF TETRAHEDRON
GAS PRESSURE MINIMUM 0.014 PSIA (100 N/SC P)

DIRECTIONAL CALORIMETER SPOC4 11113

THE DIRECTIONAL CALORIMETER WILL BE COMPOSED OF A PYROELECTRIC IR DETECTOR, DETECTOR PREAMPLIFIER, OPTICAL MODULATOR (LIGHT CHOPPER) AND DEMODULATOR (RECTIFIER). IT WILL BE USED TO DETERMINE THE BEGINNING OF SOLIDIFICATION BY OBSERVING THE ACCOMPANYING THERMAL ARREST AND WILL MEASURE THE RATE OF RELEASE OF THE HEAT OF FUSION

PYROLYTIC IR DETECTOR SPECTRAL

RANGE.....2E-5 TO 3E-3 INIC-5 TO 75
MICROMETERS)
INCIDENT FLUX RANGE.....6.5N/IN2 (1 - 10N/CM2)

MECHANICAL SAMPLE PLACEMENT AND RETRIEVAL SPOC4 11113

LIQUID SYRINGE DISPENSER SPOC4 21113

INERTIAL INJECTOR SPOC4 11235

VACUUM CATCH TUBE SPOC4 11235

SOLID SAMPLE STORAGE SPOC4 11112

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

LIFE SCIENCES

CAMERA - STILL FILM (CORE) LS988 11111

STILL PHOTOGRAPHIC CAPABILITY WHICH IS SHARED WITH OTHER EXPERIMENTS. CAMERA HAS INTERCHANGEABLE LENS AND IS COMPATIBLE WITH EQUIPMENT SUCH AS MICROSCOPES AND OSCILLOSCOPES.

FILM FORMAT 35 MM
LENS 50 MM, F1.2
VIEW ANGLE 41 DEGREES
LENS 35 MM, F1.8
VIEW ANGLE 64 DEGREES
SHUTTER SPEED 8 TO 1/500 SEC
VIEWING THRU LENS
FOCUS SPLIT IMAGE

CAMERA - STILL PCLAROID (CORE) LS002 41111

PCLAROID PHOTOGRAPHIC CAPABILITY SHARED WITH MOST EXPERIMENTS ALSO ADAPTABLE WITH MICROSCOPE AND OSCILLOSCOPE INTERFACES.

FILM FORMAT PCLAROID 107 BLACK/WHITE
..... 3.25X4.25 IN (82.55X192.95 MM)
LENS 50 MM, F3.5 TO F22
SHUTTER SPEEDS 8 TO 1/60 SEC.
RANGE FINDER SPLIT IMAGE
VIEWING THRU LENS
FIELD OF VIEW 40 DEGREES

RECORDER - STRIP CHART (CORE) LS1013 11111

RECORDS UP APPROPRIATE LIFE SCIENCES DATA.

FREQUENCY RESPONSE DC TO 150 HZ
SENSITIVITY 0.1 TO 50 MV/MM
CHANNELS 8
CHART SPEED 1 TO 500 MM/SEC
WRITING METHOD THERMAL

OSCILLOSCOPE - PERSISTENT CRT (CORE) LS976 11111

ELECTRONIC MONITOR AND DATA MEASUREMENT CAPABILITY OF LIFE SCIENCE EXPERIMENT EQUIPMENT WITH A SCREEN HAVING LONG TERM RETENTION.

BANDWIDTH DC TO 5 MHZ (MINIMUM)
CHANNELS 2
DEFLECTION FACTOR 1 MV TO 10 V PER DIVISION
PHOTOGRAPHIC CAPABILITY YES
PERSISTENCE UP TO 6 HOURS

MICROSCOPE, DISSECTING LS1004 11111

PERMIT DETAIL DISSECTING OF BIOLOGICAL SPECIMENS SUCH AS PLANTS, ANIMAL TISSUE AND ORGANS.

STEREO EYEPIECE
ADJUSTABLE BRIGHT AND DARK FIELD ILLUMINATION

MICROSCOPE, COMPOUND LS1005 11111

GENERAL PURPOSE BINOCULAR MICROSCOPE FOR MICROSCOPIC STUDIES OF TISSUES.

MAGNIFICATION 10X TO 100X
PHOTOGRAPHIC CAPABILITY PCLAROID FILM
LIFTING DARK FIELD OR LIGHT FIELD AND
PHASE CONTRAST

RECORDER - VOICE (CORE) LS1006 11111

CAPABILITY OF RECORDING ORAL COMMENTS.

CHANNELS 4
BANDWIDTH 50 TO 20000 HZ

MICRODISSECTION TOOL KIT LS998 11111

STANDARD MEDICAL MICRODISSECTION KIT

MICROTOME (CORE) LS021 11111

CAPABILITY TO THINLY SLICE EXPERIMENT SPECIMEN.

SECTIONING 1 TO 50 MICROMETERS

CRYSTAT (CORE) LS022 11111

CAPABILITY OF STORING QUANTITIES OF CRYOGENIC LIQUID

QUANTITY 5 GALLONS (10.0019 CU M)
TEMPERATURE -319 F (-195 C)

VETERINARY MEDICAL KIT (CORE) LS997 11111

INCLUDES OTOSCOPE/OPHTHALMOSCOPE, REFLEX HAMMER, HEMOGLOBINOMETER, Tourniquet, SYRINGES, NEEDLES, SCALPERS, HEMOSTATS, TWEEZERS, AND OTHER STANDARD ITEMS.

EVEN LS027 21111

CAPABILITY TO HEAT, MAINTAIN CONSTANT TEMPERATURE, AND DRY TEST EQUIPMENT, LIQUIDS AND TEST SPECIMEN.

TEMPERATURE RANGE 32 TO 450 F (0 TO 233 C)
CAPACITY 1 CU FT (0.03 CU M)

HISTOLOGY KIT/SLICE CABINET (CORE) LS994 11111

STANDARD KIT CONTAINING PROBES, KNIVES, NEEDLES, SYRINGES, COVER SLIPS, LADLES PLUS OTHER NORMALLY INCLUDED ITEMS.

CONSTANT TEMPERATURE BATH (CORE) LS1015 21111

HEAT SINK MAINTAINING A CONSTANT TEMPERATURE FOR TEST TUBES AND VIALS.

TEMPERATURE 248.7-1.0 F (124.7-0.5 C)

SPECTROPHOTOMETER IN (CORE) LS032 11111

PERFORM SPECTROPHOTOMETER ANALYSIS OF LIFE SCIENCES SPECIMEN IN THE INFRARED RANGE.

SPECTRAL RANGE 2.5 TO 12 MICROMETERS
RESOLUTION 0.1 PERCENT OF WAVELENGTH
SCAN TIME 5 SEC TO 50 FOLR
SAMPLE TEMPERATURE AMBIENT TO 462 F (250 C)

GAS ANALYZER, AUTOMATIC LS993 11111

DETERMINE THE PARTIAL PRESSURES OF OXYGEN AND CARBON DIOXIDE DISSOLVED IN BLOOD SAMPLES AND TO DETERMINE THE HYDROGEN ION CONCENTRATION.

(REQUIREMENTS UNSPECIFIED, COMMERCIAL UNIT USED AS MODEL)
PH 6.000 TO 8.000
CO2 PARTIAL PRESSURE 0 TO 200 MM HG
O2 PARTIAL PRESSURE 0 TO 200 MM HG
ACCURACY:
PH +/- 0.003P
CO2 PARTIAL PRESSURE +/- 0.5PH
O2 PARTIAL PRESSURE +/- 1PH +G AT 200 PC2
..... +/- 10PH +G AT 200 PC2

URINE ANALYZER - AUTOMATIC (CORE) LS994 11111

CAPABILITY TO AUTOMATICALLY PERFORM URINE ANALYSES.

MASS SPECTROMETER LS994 11111

PROVIDE ANALYSIS OF UNKNOWN IONIZED GAS BY ATOMIC MASS MEASUREMENT.

MASS RANGE 1 TO 400 AMU
SCANNING RANGE VARIABLE
VARIABLE SCANNING RATE 50 MILLISECOND TO 600 SEC/SCAN
MAX OPERATING PRESSURE 1E-10 TORR
SELECTABLE CENTER MASS
SELECTABLE SCANNING WIDTH

PH METER LS1016 11111

METER MEASURES HYDROGEN ION CONCENTRATION OF SOLUTIONS.

PH RANGE 0 - 14
MEASUREMENT ACCURACY 0.02PH

ULTRASONIC CLEANER (CORE) LS047 11111

HIGH FREQUENCY ENERGY FOR CLEANING EXPERIMENT AND MAINTENANCE EQUIPMENT.

OPERATING FREQUENCY 20000 HZ
CAPACITY 1 CU FT

VOLTMETER LS1012 11111

PORTABLE RELIABLE MULTI-FUNCTION METER FOR GENERAL PURPOSE LABORATORY WORK AND TROUBLE SHOOTING.

DC VOLTS 100 MV TO 10 KV
DC CURRENT 1 MA TO 10 AMPS
AC VOLTS 10 MV TO 300 VOLTS
AC CURRENT 1 MA TO 1 AMP
OHMMETER 10 OHMS TO 10 MOHMS

SIGNAL GENERATOR LS052 11111

ACCURATELY REPRODUCE SELECTED FREQUENCIES AS REFERENCES FOR PURPOSES OF MAINTENANCE AND REPAIR.

BANDWIDTH DC TO 20 KHZ
OUTPUT VOLTAGE 3 V TO 0.1 MV

PLETHYMOGRAPH, LIMB INCLUDING COUPLER LS1014 11111

MEASURE CHANGES IN BLOOD VOLUME AND VASCULAR RESPONSES.

MEASUREMENT RANGE HIGH RANGE-15 TO 500 CHMS
..... LOW RANGE-1.5 TO 50 OHMS

ACCELEROMETER - SPECIMEN MOVEMENT LS901 12111

ORGANISM CONTAINER ACTIVITY AND ACTIVITY LEVEL MONITOR.

ACCURACY +/- 0.000 G
FREQUENCY RESPONSE 0 TO 100 HZ

ALCIMETER LS902 11111

FREQUENCY SOURCE ADJUST IN BOTH FREQUENCY AND AMPLITUDE FOR EARTHQUAKE OR BOAT VIBRATORS.

FREQUENCY 125 TO 8000 ADJUSTABLE IN OCTAVES
ATTENUATOR CONTINUOUSLY ADJUSTABLE IN 100 STEPS

ANALYZER - ATOMIC ABSORPTION SPECTROPHOTOMETER LS903 11111

QUANTITATIVE DETERMINATION OF METALLIC AND SEMIMETALLIC ELEMENTS IN SOLUTIONS AND INDICATING CONCENTRATION.

WAVELENGTH RANGE 200 TO 800 MILLIMICRONS
MODES OF OPERATION ATOMIC ABSORPTION, FLAME
EMISSION, UV-VISIBLE SPECTRO-
PHOTOMETRY
COUNTER ACCURACY +/- 0.4 NANOMETERS

ORIGINAL PAGE IS
OF POOR QUALITY



CELL COUNTER - BLOOD LS9C7 11111
MEASURE BLOOD CELL PROPERTIES OF BLOOD SAMPLES. ESSENTIALLY A
COUNTER COUNTER
PARAMETERS MEASUREDHEMOGLOBIN,HEMATOCRIT,RED BLOOD
CELL COUNT,WHITE BLOOD CELL
COUNT,MEAN CELL VOLUME,MEAN
CELL HEMOGLOBIN AND MEAN CELL
HEMOGLOBIN CONCENTRATION
REPRODUCIBILITY.....+/-1%

BACTERIAL COLONY COUNTER LS9C8 11111
MANUAL COUNT OF BACTERIAL COLONIES. UNIT IS STANDARD TYPE

DIGITAL PLTTER, PRINTER LS9C9 11111
KEYBOARD LS910 11111
TYPEWRITER TYPE KEYBOARD USED FOR PROVIDING INSTRUCTIONS FOR THE
COMPUTER

ELECTROPHYSIOLOGY BACKPACK LS912 11235
BACKPACK CONTAINS THE NECESSARY ELECTRONICS FOR SENSING AND
TRANSMITTING MAN'S PHYSIOLOGICAL DATA SUCH AS ECG, EEG, EPG,
ECG, ETC. INCLUDES SENSORS, SIGNAL CONDITIONERS, MULTIPLEXERS,
A/C CONVERTERS AND TRANSMITTERS

ELECTROPHYSIOLOGY RECEIVER LS913 11111
PROVIDE SPECIAL RECEPTION OF CARDIOVASCULAR AND NEURAL ELECTRO-
PHYSIOLOGICAL EVENTS VIA BIOTELEMETRY SYSTEMS
SIGNAL TYPES RECEIVED.....ELECTROENCEPHALOGRAPH
.....ELECTROCARDIOGRAPH
.....VECTOCARDIOGRAPHY
.....BALISTOCARDIOGRAPHY
.....IMPEDANCE CARDIOGRAPHY
.....PHONOCARDIOGRAPHY

BLOOD CLOT FIBRINOMETER LS915 11111
PROVIDE AUTOMATIC MEASUREMENT OF PLASMA COAGULATION TIME
ACCURACY.....+/-0.1SECCND

FREEZER, CRYOGENIC LS916 11111
PROVIDE A MEANS FOR FREEZING AND STORING BIOLOGICAL SPECIMENS
STORAGE CAPACITY.....115 CL. IN. (1.89LITERS)
OPERATING TEMPERATURE.....-320F(178K)
MIN CRYOGEN HOLDING TIME.....3 WEEKS BETWEEN REFILLS

HEMATOCRIT, ELECTRONIC LS917 11111
AUTOMATICALLY DETERMINE THE PERCENT HEMOCRIT IN BLOOD. RAPID &
READOUT AND TEMPERATURE COMPENSATOR REQUIRED

DOSIMETER, RADIATION DETECTOR LS918 11111
PORTABLE DEVICE WHICH WILL ALERT USER TO RADIATION LEVELS IN
EXCESS OF A PREDETERMINED LEVEL AND TO PROVIDE A DIRECT READOUT
OF THE TOTAL CUMULATIVE RADIATION DOSE.
RATE MEASUREMENT RANGE 0 TO 0.1 MR/HR
C TO 50 R/HR
ADJUSTABLE RADIATION ALARM 0.05 MR/HR TO 50 R/HR
CUMULATIVE DOSE RANGES 0 TO 200 MR
C TO 500 R
ACCURACY (X- & GAMMA RAY) +/- 10 %

RADIATION DETECTOR LS919 11111
ACCURATE MEASUREMENT OF EITHER ALPHA, BETA, GAMMA, X-RAY OR
NEUTRON RADIATION AS IT IS APPLIED TO A BIOLOGICAL SPECIMEN.
RATE RANGES 0 TO 0.1 MR/HR
C TO 1000 R/HR
CALIBRATED READOUTS COUNTS/SEC, /MIN, AND TOTAL

ISOTOPE SOURCE -SELF CONTAINED LS920 11111
RADIATION STRESS SOURCE USED IN CONDUCTING RED BLOOD CELL
SURVIVAL STUDIES. TYPICAL ASSORTMENT - NITROGEN 15 COMPOUNDS,
RADIO-LABELED SUBCELLULAR PRODUCTS (RIBOSOMES), RADIO-LABELED POLY-
NUCLEOTIDES OR TITRATED THYMIDINE.
RADIATION LEVEL 1 TO 100 MICROCURIE

BIOTELEMETRY RECEIVER - COMPACT CASE MODULE LS921 121113
MICROBACKPACKS FOR SMALL VERTEBRATE ACQUIRING AND TRANSMITTING
BIO DATA.
SIGNAL CONDITIONER LS922 11111
ELECTRICAL SIGNAL TRANSFORMATION BETWEEN ANY NUMBER OF TRANS-
DUCERS AND MULTIPLEXER/ANALOG TO DIGITAL CONVERTER
COMPUTER INPUT +/- 10 VOLTS

SOUND LEVEL METER LS923 11111
SOUND LEVEL SURVEYS AND LIMITED FREQUENCY ANALYSIS IN REAL TIME.
FREQUENCY RESPONSE 125 TO 8000 HZ
SOUND LEVEL RANGE 40 TO 100 DB

STAINING SYSTEM - BACTERIOLOGICAL LS924 11111

STAIN AND FIX BIOLOGICAL SPECIMEN. FOR MICROSCOPIC EXAMINATION.

TASKBOARD - FORCE/TORQUE LS925 11131
MEASURE MAN'S CAPABILITY TO APPLY FORCES AND TORQUES IN A VARIETY
OF DIRECTIONS FROM VARIOUS HAND/BODY ORIENTATIONS AND RESTRAINT
CONDITIONS.
HAND FORCE APPLICATIONS +/- 0.5 POUNDS

RADIATION COUNTER - BIOCHEMICAL SAMPLE LS926 11111
FAST, ACCURATE AND EASY BIOCHEMICAL RADIATION COUNTING SYSTEM
BASED ON STANDARD SIZE BIOCHEMICAL SAMPLES.
STANDARD PLANCHET LP TO 2 IN DIA
CAPACITY OF PLANCHETS 120

MASS SPECTROMETER LS928 11111
PROVIDE AN ANALYSIS OF AN UNKNOWN IONIZED GAS BY ATOMIC MASS
NUMBER
MASS RANGE.....0 TO 40 AMU

AIR PARTICLE SAMPLE COLLECTOR LS929 11111
AIR PARTICLE AND MICRO-ORGANISM SAMPLES FOR AIR QUALITY INTERPRE-
TATION
AIRFLOW 1 CL FT/MIN
SAMPLING TIME - CONTINUOUS ... 1 HOUR

AUTOCALYZER - MULTIPLE LS930 11111
AUTOMATIC ANALYSIS OF BLOOD AND APPROXIMATELY 30 OTHER FLUIDS.

SPECTROPHOTOMETER - GENERAL ANALYZER LS931 11111
SPECTRAL ANALYSIS OF GASES AND LIQUIDS INCLUDING SOLIDS OR LIGHT
SOURCES.
FREQUENCY 250 TO 2500 NANOMETERS
(2500 TO 25000 ANGSTROMS)
ACCURACY +/- 0.4 NANOMETERS

ANESTHETIZER - INVERTEBRATES LS932 11111
TENDER INVERTEBRATE ORGANISMS INSENSIBLE TO FACILITATE HANDLING.
ANESTHETIC GAS CARBON DIOXIDE

ANTHRUPOMETRIC GRID LS933 11111
ANTHROPOMETRIC MEASUREMENTS OF VERTEBRATES EITHER REMOTELY OR
AUTOMATICALLY.
SMALL GRID 5 X 15 CM SPACED AT 1 MM
MEDIUM GRID 50 X 75 CM SPACED AT 2 MM
LARGE GRID 2 X 2 M SPACED AT 5 MM

BENCH, LAMINAR AIRFLOW LS934 11111
A GLOVE BOX WITH RELATIVELY HIGH AIR FLOW FOR CONTROL OF
PARTICULATE AND GASEOUS CONTAMINANTS WITHIN

BENCH, GENERAL EXPERIMENTS LS935 11111
PROVIDE WORK AREA FOR PREPARATION OF EXPERIMENTS, EXPERIMENT
MAINTENANCE OBSERVATIONS, AND PREPARATIONS FOR RETURN TO EARTH
ELECTRICAL UTILITY.....28VDC, 400HZ AC, 60HZ AC
VACUUM UTILITY.....10E-6 TORR
PRESSURE SOURCE.....50 PSIG (3.4564 N/M2)
AVAILABLE GASSES.....OXYGEN, NITROGEN, CARBON DIOXIDE

CAMERA CINE LS936 11111
PROVIDE VISUAL RECORDS SUCH AS PHOTOMICROGRAPHY, TIME LAPSE
RECORDS, PHOTOMICROSCOPY AND GENERAL EXPERIMENT DOCUMENTATION
OPERATION MODE.....VARIABLE FRAME RATE OR SINGLE
FRAME PULSE OPERATION UNDER
REMOTE CONTROL
FILM SIZE.....16 OR 35MM
LEN SELECTION.....CM MOUNT LENS SYSTEM; STANDARD
ZOOM LENS AS A MINIMUM
EXPOSURE SYSTEM.....AUTOMATIC OR ZOOM
FILM CAPACITY.....50, 100 OR 200FT (15.2,30.4 OR
60.9M) ROLLS

CAMERA CONTROLLER LS937 11111
DEVICE TO CONTROL THE OPERATION OF VIDEO CAMERAS THROUGHOUT THE
LABORATORIES

CAMERA, PLATE FILM LS938 11111
PROVIDE STILL PHOTOGRAPHIC COVERAGE
FILM TYPES.....HIGH RESOLUTION GLASS PLATE
.....STANDARD 2 1/4 X 2 1/4 IN (5.75
X 5.75 CM) ROLL
.....ZOOM
.....POLAROID
SHUTTER SPEED.....1 SEC TO 1/500 SEC

VIDEO CAMERA, B/W LS939 401131
PROVIDE MEANS OF ACTIVITY MONITORING, EXPERIMENT DATA ACQUISITION
ETC. SYSTEM WILL INTERFACE WITH A 40 INPUT VIDEO MULTIPLEXER TO
PERMIT THE MONITORING OF 40 SIMILAR VIDEO CAMERAS.
VISUAL RESPONSE.....APPROXIMATE HUMAN EYE



ELECTROPHYSIOLOGY DISPLAY LS954 11114

DEDICATED DISPLAY FOR PHYSIOLOGICAL DATA TRANSMITTED BY THE ELECTROPHYSIOLOGY BACKPACK

ELECTROMETER LS955 11111

MEASURE SMALL VOLTAGES AND CURRENTS FOUND IN NERVES AND MUSCLES.

DC VOLTAGES MEASURED.....+/-10MHCVOLTS TO +/-1VOLT
AC VOLTAGES MEASURED.....+/-1 PA TO +/-3 PA
DRIFT.....+/-5V/DAY AFTER 10MIN WARMUP
OUTPUT VOLTAGE.....+/-10V

ELECTROPHORESIS APPARATUS LS956 11111

SEPARATES PROTEIN AND AMINO ACID CONSTITUENTS IN SERUM, PLASMA, URINE OR SPINAL FLUID FOR QUANTITATIVE ANALYSIS.

VOLTAGE.....0 TO 500 V DC
SUBSTRATE.....STARCH BLOCK, PAPER OR GEL
SAMPLE VOLUME.....0.001 TO ML

FREEZER, GENERAL LS957 11111

STORAGE OF SERUM, PLASMA, SPECIMENS AND ORGANISMS

OPERATING TEMPERATURE.....-4F TO -200F
TEMPERATURE TOLERANCE.....+/-4F/+/-20F
STORAGE VOLUME.....4 CU. FT. TO 1100 CU. FT.

FREEZER, LOW TEMP LS958 11111

STORAGE OF EXPERIMENT SPECIMENS

OPERATING TEMPERATURE.....-54F(-70C)
STORAGE VOLUME.....1 CU. FT. TO 0.28 CU. FT.

REFRIGERATOR LS959 11111

STORE SERUM AND PLASMA

OPERATING TEMPERATURE.....32 TO 40F(0 TO 4C)
STORAGE VOLUME.....1 CU. FT. TO 0.28 CU. FT.

REFRIGERATOR, RADIO ISOTOPE STORAGE LS960 11114

STORE LIQUID AND SOLID RADIOACTIVE MATERIALS

OPERATING TEMPERATURE.....-41F(-45F/+/-20C)
VOLUME.....100 LIT. TO 0.28 CU. FT.
ISOTOPE STORED.....P3, C14, F259, CR51, I131, CE45

GAS ANALYZER, CO2 SPECIFIC LS961 11111

THIS DEVICE IS USED TO MONITOR ATMOSPHERIC CARBON DIOXIDE LEVELS

PARTIAL PRESSURE RANGE.....130 - 270CMHG(10-20 MM HG)

GAS CHROMATOGRAPH LS962 11111

MEASURES CONCENTRATION OF GAS, LIQUID, AND SOLID CONSTITUENTS OF BIOLOGICAL SAMPLES

OPERATING TEMPERATURE.....-40 TO 752F(-40 TO 400C)
CARRIER GAS.....HELIUM
DETECTORS.....HYDROGEN FLAME
.....THERMAL CONDUCTIVITY
FRACTION COLLECTION SYSTEM

GAS ANALYZER, WATER VAPOR LS963 11111

MEASURE WATER VAPOR CONTENT

MOISTURE RANGE.....0 TO 99% RELATIVE HUMIDITY
.....(0.001 TO 2000MICROGRAMS/LITER)
TEMPERATURE RANGE.....-100 TO 140F (-110 TO 60C)

KIT, BENCH CHEM ANAL LS965 11111

KIT CONTAINS THE TOOLS AND EQUIPMENT TO MANAGE CHEMICALS AND BIOLOGICALS DURING VARIOUS MANUAL PROCEDURES GENERALLY PERFORMED WITHIN THE GLOVE BOX.

TOOLS.....SOLID TRANSPORT TOOLS, GRAVITY INDEPENDENT PIPETTES, VIALS, BOTTLES, AND TEST TUBES, CHEMICALS, STOPPERS, FILTERS, AND SAFETY SHIELDS

KIT, HEMATOLOGY LS966 11111

PREVIEW TOOLS FOR SAMPLING, HANDLING, TRANSFERRING AND ANALYZING BLOOD

ITEMS.....HEMATOMETER KIT; 10 LAMPOA DISCIPIPETTES; COVERSLIPS; SLIDES; WBC AND RBC CILIENT; CRITSEAL; HEMOCRIT TUBES; MICROHEPARINIZED; HEMOCRIT TUBES (MICRO, PLAIN); BLOOD CULTING PIPETTES (WBC, RBC); COUNTER (URINE TEST STRIPS); ALP ADAPTERS; VACUTAINER; VACUTAINER VACUTAINER-NEEDLE UNIT; VACUTAINER TUBES; ASSETEC(2M); SYRINGE (LARGE); PEDIATRIC ALCOHOL SWABS; LANCETS; NEEDLE, 25GA, 5/8IN; NEEDLE 25GA, 1 1/2IN

KIT, LINEAR MEASUREMENT LS967 11111

DETERMINE SIZE, AMPLITUDE, DISTANCE, CIRCUMFERENCE, ETC

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

ITEMS.....RULERS; TAPES; LINEAR COMPARATORS;
GRIDS; CALIPERS, INSIDE AND
OUTSIDE; VERNIER CALIPERS;
MICROMETERS, INSIDE AND OUTSIDE

KIT, MICROBIOLOGY LS968 11111

PROVIDE TOOLS TO FACILITATE GROWING AND ANALYZING MICROBIAL
ORGANIZING

ITEMS.....INOCULATING LOOPS; INOCULATING
NEEDLES; COTTON SWABS, PACKAGES;
STERILE; SYRINGE, 5ML, STERILE;
SLIDES, MICRO; ALCOHOL TUBES,
15X75MM, STERILE; CAPPECS;
ZEPHIRAN, TINCTURE, 1:500;
BACTERIOGENATOR/STERILIZING UNIT;
THIOGLYCOLLATE, TUBES; STUART
TRANSPORT MEDIA, VIALS; TSA
SLANTS

KIT, ORGANISM HOLDING AND MANAGEMENT LS969 11111

PROVIDE TOOLS AND DEVICES USED IN THE HOLDING AND HANDLING OF THE
ORGANISMS

ITEMS FOR SMALL VERTEBRATES.....GLOVE POWDER; GAGE SHIELD,
PLASTIC; PLASTIC LINER WITH PACE
FOOD PELLET DISPENSER; ALCOHOL
SWABS; TOWELS, PAPER, DISPS;
PLASTIC BAGS (FOR EXPENSES);
PLASTIC BAGS (FOR NONEXPENSES);
ORGANISM TRANSFER CAPSULE;
ANIMAL TAGS

ITEMS FOR PLANTS.....WATERING DEVICE (SPRAY BOTTLE);
LABELS; FERTILIZER PACKETS;
STAKES (LOCH)

KIT, PLANT TOOLS LS970 11111

PROVIDE TOOLS FOR VARIOUS PLANT MANIPULATIONS

ITEMS.....SCISSORS, TWEZERS, SPATULA,
APPLICATOR, SPLINTS, TAPE
DISPENSER, HYPODERMIC NEEDLES,
SYRINGES, FORCEPS, PLIERS,
SCREWDRIER, SCALPELS, SEALERS

KIT, GENERAL TOOL LS971 11111

MECHANICAL AND ELECTRICAL TOOLS AND HARDWARE TO PROVIDE
CONVENTIONAL DIAGNOSTIC, MAINTENANCE, AND SERVICE FUNCTIONS.

TOOLS.....HAMMER, WRENCHES, PLIERS, SCREW
DRIVERS, DRILL, ADHESIVE TAPE,
WIRE CUTTERS, WIRE TIES, WIRE,
LUBRICANTS, FLASHLIGHT, SCISSORS,
MULTIMETER, FASTENERS, CLAMPS,
LAMP

LYOPHILIZER LS972 11111

PROVIDE FREEZE DRYING OF SMALL SPECIMENS

SAMPLE CAPACITY.....20-100 IN (110CC)
TEMPERATURE RANGE.....-40F TO 252F (-40 TO 122C)

MASS MEASUREMENT DEVICE, MACRO LS973 11132

MEASURE MASS OF ITEMS SUCH AS BAGS OF FOOD, BEVERAGE CONTAINERS,
CONTAINERS OF URINE AND FECES, LARGE SPECIMENS, ETC

MASS RANGE.....0.22 TO 66LB (0.1 TO 30KG)

MASS MEASUREMENT DEVICE, MICRO LS974 11132

MEASURE SMALL TEST SPECIMENS

MASS MEASUREMENT RANGE.....1PG TO 100 G

OPTISCAN-FIELD AND FIXATION POINT RECORDER LS975 11111

THIS INSTRUMENT PERMITS THE RECORDING OF THE VISUAL FIELD AS
THE SUBJECT'S FIXATION WITHIN THE FIELD WITHOUT RESTRICTION OF
HEAD MOVEMENTS

ACCURACY.....2 DEG OF ARC

CATALYTIC OXIDIZER SYSTEM LS977 11235

REMOVES UNDESIRABLE COMBUSTION PRODUCTS

OPERATING TEMPERATURE.....70CF (170C)
AIR FLOW RATE.....2CFM AT 14.7PSIA AND 7CF

BLOOD PRESSURE CUFF SYSTEM LS978 11132

PRESSURE MEASURING TRANSDUCER CUFF USED TO MEASURE BLOOD PRESSURE
OF PRIMATES. SYSTEM INCLUDES AIR PUMP FOR AUTOMATED MEASURING,
PRESSURE CUP AND TRANSDUCER, AND SIGNAL CONDITIONER.

RADIATION SOURCE STORAGE LS979 11112

STORAGE FACILITY FOR RADIOACTIVE MATERIAL.

VOLUME APPROXIMATELY 1 CU FT
RADIATION PROTECTION LIMIT LP TO 500 MICROCURIES

RECEIVER - BIOTELEMETRY LS980 11111

RECEPTION OF BIOTELEMETRY SIGNALS FROM DEEP BODY TEMPERATURE AND
ANIMAL ACTIVITY AS WELL AS FOR ELECTROMAGNETIC FIELD MONITORING.

FREQUENCY RANGE TO 5 MHZ

SENSITIVITY 0.5 MV FOR 20 DB QUIETING
TUNING VARIABLE AND CRYSTAL CONTROLLED

VISION TESTER LS981 11235

OPTICAL DEVICE (WITH HOG AND RESPONSE KEYBOARD) THAT MEASURES
A VARIETY OF VISUAL FUNCTIONS.

WASTE STORAGE SYSTEM LS982 11111

PROVIDE FOR HANDLING, STORAGE AND DISPOSAL OF SOLID AND LIQUID
WASTE MATTER FROM THE EXPERIMENTS.

CAPACITY 5 CU FT

STERILIZER - TIGL LS983 11111

STERILIZE MISCELLANEOUS SMALL METAL HAND TOOLS SUCH AS SCALPELS,
BY MEANS OF ELECTRICAL INDUCTION HEATING.

CAPACITY 0.1 CU FT

AUDIO STEREO HEADSET LS985 11111

EARPHONES FOR VARIOUS HSI HEARING TESTS

FLOW METER, WATER MANIFOLD LS989 11235

IN LINE MEASUREMENT OF WATER FLOW, GENERALLY A LOW RATE
ASSOCIATED WITH WATER CONSUMPTION BY ORGANISMS

FLOW METER, GAS LS992 11111

MEASURE AIR FLOW

IMPEDANCE PNEUMOGRAPH LS995 11111

MEASURE BREATHING CYCLE CHARACTERISTICS

KIT, CLEAN UP LS999 11111

GENERAL PURPOSE CLEAN UP - SPONGES AND WIPES

KIT, MEDICAL SURGICAL LS1000 11111

FOR VARIOUS MINOR SURGICAL PROCEDURES

KIT, PHYSIOLOGY LS1001 11111

SPONGES, SPONGE SQUEEZER, VIALS, CALORIC STIMULATOR FOR EAR
CANAL, SYRINGES, THERMISTERS, ETC.

KIT, TOOL - INSECT MANIPULATION LS1002 11111

TOOL KIT FOR COUNTING, SORTING, EXAMINATION, ETC.

MICROPHONE LS1007 11111

MICROPHONES USED FOR BIOLOGY, BIOMEDICINE AND MANNEQ SYSTEMS
INTEGRATION

MICROPHONE AMPLIFIER LS1008 11111

AMPLIFY SIGNALS FROM AMPLIFIER

MONITOR, VIDEO LS1009 11111

MONITOR TV PICTURES OF ANIMAL AND OTHER LAB ACTIVITIES

PUMP GAS, CIRCULATING LS1011 11111

PUMP FOR SEALED PLANT GROWTH CHAMBER

TIMER, EVENT LS1020 11111

GENERAL PURPOSE ELAPSEC TIME DEVICE USED FOR VARIOUS PSYCHOMOTOR
TESTS

CARDIOGRAPH - IMPEDANCE LS100 11111

BLOOD FLOW RATE FROM HEART.

BLOOD FLOW RATE 2 TO 25 LITERS/MIN +/- 5 PCT FS

ECG/VCG LS101 11131

MONITOR HEART CONDITION. ELECTRO/VECTORCARDIOGRAPH.

SENSITIVITY 0 TO 3 MV +/- 1 PCT FS
HEART RATE 40 TO 180 BEATS/MIN

EEG - ELECTROENCEPHALOGRAPH LS102 11131

MONITOR ELECTRICAL IMPLUSES GENERATED BY THE BRAIN.

SENSITIVITY 10 TO 200 MICROVOLTS +/- 1 PCT FS
FREQUENCY RESPONSE 0.5 TO 100 HZ
CHANNELS 8

DYNAMOMETER LS103 11133

MEASUREMENT OF FORCE OR POWER.

ARM FLEXION 25 TO 75 LB (10 TO 35 KG)
BACK EXTENSION 100 TO 200 LB (40 TO 80 KG)

EPC - ELECTROMYOGRAPH LS104 11131

MONITOR THE MOVEMENT OF THE BODY MUSCULAR SYSTEM.

SENSITIVITY 0.01 TO 5 MV +/- 1 PCT FS
BANDWIDTHS 0.5 TO 200 HZ, 0.5 TO 1000 HZ
CHANNELS 1

PHONOCARDIOGRAM LS106 11111
AURAL MONITOR OF HEART BEAT.
SENSITIVITY 0.1 TO 1000 HZ, +/-5 PCT FS
METABOLIC ANALYZER LS1003 11131
MEASURE O2 CONSUMPTION, CO2 PRODUCTION, MAXIMAL O2 CONSUMPTION
AND ALVEOLAR PO2 AND PCO2.
O2 VOLUME 0 TO 40 L/MIN, +/-3 PCT FS
CO2 VOLUME 0 TO 40 L/MIN, +/-3 PCT FS
ALVEOLAR O2 20 TO 120 MM HG +/-2 PCT FS
ALVEOLAR CO2 20 TO 70 MM HG +/-2 PCT FS
PULMONARY FLOWMETER LS990 11131
BLOOD FLOW RATES AND VOLUME.
FLOW RATE 0 TO 15 L/MIN, +/-3 PCT FS
VOLUME 0 TO 7 L, +/-3 PCT FS
TRANSCUTANEOUS O2 PLETH FLOWMETER LS991 11112
MEASURE PULSE WAVE LENGTH.
ULTRASONIC SIGNAL 1 AND 5 MHZ
PULSE WAVE LS112 11112
MEASURE BLOOD VELOCITY
VELOCITY 3 TO 15 M/S, +/-0.1 M/S
ERGOMETER - BICYCLE LS986 11131
WORKLOAD MEASURING DEVICE.
LOAD RANGE 0 TO 300 WATTS, +/-2 PCT FS
ROTATIONAL SPEED 40 TO 90 RPM
LITTER CHAIR - ROTATING LS114 11131
GRAVITY INDUCING DEVICE AND HORIZONTAL WORK TABLE.
ROTATION SPEED 0 TO 60 RPM
LBNP - LOWER BODY NEGATIVE PRESSURE LS115 11131
MEASURE INTERNAL BODY PRESSURE DIFFERENTIAL.
PRESSURE DIFFERENTIAL 50 MM HG BELOW CABIN AMBIENT
EXERCISER/ERGOMETER - PRIMATE LS175 11133
PRIMATE EXERCISER AND WORKLOAD MEASURING DEVICE.
CALORIMETRY MODULE - PRIMATE LS964 11135
PORTABLE AND COLLAPSABLE UNIT WHICH HOLDS, FEEDS, AND ANALYZES
THE PRIMATE SPECIMEN WITHIN.
PO2 160 +/-10 MM HG
PN2 580 +/-20 MM HG
PCO2 6 MM HG OR LESS
PH2O 12 +/-3 MM HG
TEMPERATURE 77 +/-4 F (25 +/-2 C)
MASS MEASUREMENT DEVICE - WGT LS987 11131
TOTAL BODY MASS DETERMINATION BY OSCILLATING FREQUENCY.
MASS 125 TO 230 LB, +/-0.1 PCT FS
(50 TO 100 KG)
LIGHTING SYSTEM - PLANT LS1018 21114
PROVISIONS FOR LIGHT PLANT GROWING AREAS.
TYPE LIGHT DAYLIGHT-FLUORESCENT
ILLUMINATION 102 FT-C (1100 LUX/M2)
PLANT GROWTH AND SUPPORT CONTAINERS LS1019 12111
VARIOUS SIZE CONTAINER TO ENVIRONMENTALLY HOLD SEEDLINGS.
CAPACITY UP TO 20-CM PLANTS
CELLS/TISSUE HOLDING UNIT LS1016 12114
COMBINATION HOLDING UNIT/INCUBATOR CONTAINING ELECTRONIC FLUG-IN
CAPABILITY FOR MEASUREMENT DEVICES
TEMPERATURE RANGE 41 TO 140 F (5 TO 60 C)
PCO2 160 +/-5 MM
PN2 600 +/-20 MM
PCO2 3 MM OR LESS
TOTAL PRESSURE 760 +/-20 MM
RELATIVE HUMIDITY 60 TO 100 PCT
CELLS/TISSUE HOLDING UNIT LIGHTING LS1017 21114
PROVISION FOR PROVIDING COOL DIFFUSED LIGHT
ILLUMINATION 90 +/-10 FT-C
TEMPERATURE AMBIENT

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

SPACE TECHNOLOGY

CONDENSER ST102 11131
CONDENSATE TANKS ST102 21113
DEVICE TO CONDENSE STEAM IN ZERO GRAVITY

TANKS FOR STORAGE OF CONDENSATE
PRESSURE..... 20 PSI (1.38X10⁵ N/M²)

SUPPORT ELECTRONICS ST006 21114
ELECTRONICS SUPPORT COMMUNICATIONS EXPERIMENTS INCLUDING
MULTIPLEXERS WITH CONVERTERS AND OSCILLATORS

CSCILLSCOPE ST007 21111
MONITOR, MEASURE AND MAINTAIN ELECTRONIC EQUIPMENT OPERATING IN
THE L-BAND RANGE (390 TO 1550 MHZ).
USED IN ST015, ST025, ST035, ST045, ST055, ST065, ST075, AND
ST095.

BANDWIDTH..... 390 TO 1550 MHZ
CHANNELS..... 2
DEFLECTION FACTOR..... 10 MV/DIVISION
PHOTOGRAPHIC CAPABILITY..... YES

TAPE RECORDER, DIGITAL ST008 21111
RECORD AND PLAYBACK OF DIGITAL DATA IN SUPPORT OF MICROWAVE
INTERFEROMETER AND TRACKING INVESTIGATIONS.
USED IN ST015, ST025, ST035, ST045, ST055, ST065, ST075 AND
ST095.

SIGNAL TYPE..... 0 TO 5 VDC
BANDWIDTH..... 100 TO 250 KHZ
CHANNELS..... 4 (MINIMUM)
DATA REQUIREMENT..... 2 HRS (MINIMUM)

DATA BUFFER ST010 21111
CAPABILITY TO ACCEPT VERY HIGH DATA RATES, PROVIDE SHORT DURATION
STORAGE AND LOWER RATE DUMP CAPABILITY. USED IN ST015.

WORD LENGTH..... 8 TO 32 BITS
OUTPUT RATE..... VARIABLE AT 0 TO 5 VDC
STORAGE DURATION..... INFINITE

DIGITAL PHASE METER ST011 41111
METER TO DISPLAY PHASE OF SIGNAL RECEIVED BY MICROWAVE
INTERFEROMETER

BCCM CONTROLS ST012 11114
CONTROLS FOR MICROWAVE INTERFEROMETER BCCMS

RECORDER - STRIP CHART ST015 11111
RECORD ANALOG MICROWAVE RADIOMETER DATA MONITORING SEA SURFACE
TEMPERATURE AND ROUGHNESS.
USED IN ST015.

FREQUENCY RESPONSE..... 150.0 HZ
SENSITIVITY..... 10 MV/DIVISION
CHANNELS..... 2
IMPLTS..... 2

SIGNAL GENERATOR ST016 11111
ACCURATELY REPRODUCE IN PULSE FORM THE OUTPUT OF A MICROWAVE
RADIOMETER MONITORING SEA SURFACE STATE AND TEMPERATURE.
USED IN ST015.

REPETITION RATE..... CN
OUTPUT VOLTAGE..... +/- 5 VDC

SPECTRUM ANALYZER ST017 11111
SUPPORT SETUP, CHECKOUT AND SIGNAL OUTPUT ANALYSIS OF RF RADIOMETER
UNIT AND SWEEP RECEIVER ICRT DISPLAY OF SWEEP RECEIVER AMP-
LITUDE VERSUS FREQUENCY.
USED IN ST015 AND ST025.

BANDWIDTH..... 30 TO 10000 MHZ
RESOLUTION..... 2 Hz
SENSITIVITY..... 10 MV/DIVISION

CAMERA - STILL (VISIBLE) ST020 31111
STILL PHOTOGRAPHIC COVERAGE OF THE EARTH TARGET AREA IN THE
VISIBLE REGION (REQUIRES SIMULTANEOUS PHOTOGRAPH IN THE INFRARED
RANGE). USED IN ST015, ST025 AND ST035.

VIEW ANGLE..... 14.5 DEGREES
RESOLUTION..... 100 LINES/MM
FILM FORMAT..... 5 INCHES
SHUTTER..... 8 TO 1/8000 SEC

DISPLAY TERMINAL ST020 21111
TERMINAL CONSISTING OF ALPHANUMERIC KEYBOARD AND CRT DISPLAY FOR
MONITORING EXPERIMENT ACTIVITIES

COMPUTER ST024 21111
COMPUTER CAPABILITY TO PROVIDE BASIC POINTING AND TRACKING INFER-
MATION RELATIVE TO ATL STABILITY AND FLIGHT PATH REQUIREMENTS.
USED ON ST015, ST035, ST045 AND ST055.

MONITOR - TV ST028 31111

DISPLAY HIGH RESOLUTION TELEVISION PICTURES IN REAL-TIME OR FROM
VIDEO TAPE. USED ON ST015, ST025, ST035, ST045 AND ST055.

RESOLUTION..... 625 LINES
FRAME RATE..... 30 FPS
SCREEN SIZE..... 10 IN (45.4 CM)

DATA PROCESSING EQUIPMENT ST035 11111

HIGH SPEED DIGITAL DATA CONVERSION TO ANALOG FORM COMPATIBLE FOR
REAL TIME DISPLAY. USED ON ST015, ST025 AND ST035.

CHANNELS..... 32
DIGITAL DATA..... 10E+06 BPS

STORAGE - GAS SOURCE ST045 81111

GAS SOURCE STORAGE BOTTLES CAPABLE OF STORING VARIOUS GASES AT
HIGH PRESSURES. SIMILAR TO STANDARD K-BOTTLES. USED IN ST015,
ST025 AND ST035.

VOLUME - PHYSICAL..... 5 CU FT EA (0.16 CU M)
PRESSURE MAXIMUM..... 1000 PSI

GENERATOR - AEROSOL ST046 11135

CAPABILITY TO PRODUCE FINE LIQUID PARTICLES REQUIRED FOR VARIOUS
EXPERIMENTS. USED IN ST015, ST025 AND ST035.

GENERATOR - ION ST047 11113

IONIZE AEROSOL PARTICLES SO AS TO OBSERVE ZERO GRAVITY ACTION.
USED IN ST015, ST025 AND ST035.

VAPORIZER - WATER ST048 21135

CAPABILITY TO PROVIDE COOL WATER VAPOR TO EXPERIMENT CONTAINER.
USED IN ST015, ST025 AND ST035.

LIGHT SOURCE ST049 21111

PORTABLE LIGHT SOURCE USED FOR PHOTOGRAPHING EXPERIMENT OPERA-
TIONS AND RESULTS. USED IN ST015, ST025 AND ST035.

ILLUMINATION..... 100 TO 1000 LUMENS
DISTANCE..... 10 TO 30 FEET

ENVIRONMENT CHAMBER ST050 11135

OPTICAL CHAMBER IN WHICH OPTICAL PROPERTIES OF AEROSOLS
WITH CAPABILITY TO PERTURBATE THE ENVIRONMENT WITH IONS AND
OBSERVING THE AEROSOLS THROUGH A MICROSCOPE AND PHOTOGRAPHING THE
SCENE. USED IN ST015, ST025 AND ST035.

MICROSCOPE ST051 21111

MICROSCOPIC EXAMINATION OF AEROSOL PARTICLES IN THE VISUAL RANGE.
ALSO COMPATIBLE WITH BIOLOGICAL AND METALLURGICAL EXPERIMENTS.
USED ON ST015, ST025 AND ST035.

MAGNIFICATION RANGE..... 50 TO 1000 X
PHOTOGRAPHIC CAPABILITY..... YES

CAMERA - LINE ST054 61111

PHOTOGRAPHIC RECORD OF THE ACTIVITY OF AEROSOL PARTICLES IN
NORMAL, SLOW MOTION AND TIME-LAPSE SEQUENCES
USED IN ST015, ST015, ST035, ST045, ST055 AND ST065.

VIEW ANGLE..... 40 DEGREES
RESOLUTION..... 200 LINES/MM
SHUTTER..... 1/50 TO 1/5000 SEC
FILM FORMAT..... 16 MM
CAMERA SPEED..... 1, 8, 16, 64 FPS

TAPE RECORDERS (HIGH-SPEED MULTICHANNEL) ST056 11111

MULTIPLE CHANNEL HIGH SPEED ELECTRONIC DATA RECORDING SYSTEM.
USED ON ST015, ST025 AND ST035.

BANDWIDTH..... 1 KHZ TO 2.0 KHZ
CHANNELS..... 14
SIGNAL TYPE..... ANALOG

SPECTROMETER - GAMMA RAY ST063 11111

MEASURE SPACECRAFT GAMMA RADIATION BY CRYOGENICALLY COOLED LITHI-
UM DOPED GERMANIUM DETECTOR WITH ANTI-COINCIDENCE SHIELD. USED
IN ST015 AND ST035.

ELECTROMAGNETIC ENERGY..... 100 KEV TO 6 MEV

SPECTROMETER - CHARGED PARTICLE ST064 11111

THINLY MULTI-WIRE PROPORTIONAL COUNTER WITH INTERVENING LAYERS OF
TISSUE EQUIVALENT PLASTIC. USED IN ST015 AND ST035

PROTON ENERGY RANGE..... 5 TO 230 MEV
CELESIUM ENERGY RANGE..... 10 TO 300 MEV
TRITIUM ENERGY RANGE..... 15 TO 230 MEV
HELIUM-3 NUCLEI ENERGY RANGE..... 20 TO 400 MEV
ALPHA PARTICLE ENERGY RANGE..... 25 TO 450 MEV

SPECTROMETER - NEUTRON ST065 11111

LIQUID SCINTILLATOR WITH ANTI-COINCIDENCE SHIELD PULSE SHAPE
DISCRIMINATOR. USED ON ST015 AND ST035.

NEUTRON ENERGY RANGE..... 0.5 TO 15 MEV

DATA CONVERTER - ANALOG TO DIGITAL ST067 21111

CONVERT ANALOG SIGNAL DATA TO DIGITAL DATA FOR RECORDING AND
LATER DUMPING TO GROUND STATIONS. USED ON ST015 AND ST035.



FILM CAMERA ST070 11111
PENIUM ASYMMETRIC PERSONAL RADIATION DOSE, DOSE RATE AND DOSE HISTORY. USED IN ST015, ST025 AND ST035.
DOSE RATE C - 0.1 MH/HR
C - 1000 MH/HR

ECOSINETER - THERMOLUMINESCENT ST071 11111
PASSIVE DOSIMETER USING THERMOLUMINESCENCE TO MONITOR RADIATION. USED IN ST015 AND ST025.
ENERGY RANGE DOSE 1 TO 10E+05 RADS

REFRIGERATOR ST080 11111
PROVIDE A COOL STORAGE AREA TO HOUSE SPECIMEN. USED IN ST015, ST025 AND ST035.
TEMPERATURE 2 C (278 K)
SIZE APPROXIMATELY 3 X 2 X 1 FT (76.2X50.8X25.4CM)

CELL CULTURE ENVIRONMENTAL CHAMBER ST081 11112
PLACING FOR CULTURES OF BACTERIA WHICH CAN BE OBSERVED DURING EXPERIMENT. USED IN ST015, ST025 AND ST035.
SIZE APPROXIMATELY 5 X 5 X 7 IN (13 X 13 X 18 CM)

MICROBIOLOGICAL SAMPLE STORAGE ST082 11112
STORAGE HARK FOR SPARE-TYPE PLASTIC BAG, SWABS AND CAPSULES TO DESTROY BACTERIA IN CASE OF BREAKAGE. USED IN ST015, ST025 AND ST035.
SIZE APPROXIMATELY 10X13X6 IN (40X33X15 CM)

INCUBATOR ST083 11111
CONTAINER USED TO PROVIDE VARIABLE TEMPERATURE ENVIRONMENTS FOR BIOLOGICAL SPECIMEN. USED IN ST015, ST025 AND ST035.
SIZE APPROXIMATELY 2X1.5X0.5 FT (60X25X15 CM)
TEMPERATURE TO 100 F (40 C)

CELL CULTURE PRINTER-CELLCATER ST084 11111
PROVIDE AUTOMATIC BLOOD ANALYSIS (HEMOGLOBIN, HEMATOCRIT, RED BLOOD CELL COUNT, WHITE BLOOD CELL COUNT, PLATELET COUNT, MEAN CELL VOLUME AND ITS CONCENTRATION). USED IN ST015, ST025 AND ST035.

SPECTROMETER - VOLUMETRIC ST085 11111
PROVIDE A PHOTOGRAPHIC MEANS TO MEASURE INDEPENDENT VOLUME OF CELLS AND CELL SHAPES. USED IN ST015, ST025 AND ST035.

ELECTROPHORESIS APPARATUS ST088 11111
MEASURE ELECTROPHORESIS MOBILITY, SURFACE-ZETA POTENTIAL, AND SURFACE CHARGE DENSITY OF CELL LINES OVER THEIR LIFE CYCLE. USED IN ST015 AND ST025.

CENTRIFUGE ST090 11111
PROVIDE VARIABLE GRAVITY ENVIRONMENTS TO BIOLOGICAL SPECIMEN. GENERALLY LIMITED TO TEST-TUBE SIZE CONTAINERS. USED IN ST015, ST025 AND ST035.
GRAVITY LEVEL INDUCED UP TO 10 G REQUIRED

PRESSURE TANK ST098 11111
PROVIDE PRESSURE SOURCE FOR WATER FEED SYSTEM TO STEAM GENERATOR EXPERIMENT
OPERATING PRESSURE 100 PSIG (6.89476 MPa)

WATER RESERVOIR WITH HEATER ST099 11111
STORE WATER FOR STEAM GENERATOR EXPERIMENT
OPERATING PRESSURE 30PSIG (2.06832 MPa)

PUMP ST100 11111
PUMP WATER TO STEAM GENERATOR
FLOW RATE 0.1 LB/HR (0.04536 GPM)

STEAM GENERATOR ST101 11111
DEVICE TO GENERATE STEAM IN ZERO GRAVITY

AIR SAMPLE UNIT ST105 11111
PERIODIC COLLECTION OF AIR SAMPLES FROM THE SHUTTLE AND AT PROVIDING TYPE, RATE OF CHANGE AND QUALIFICATION OF PARTICLES AND NON-VIABLE PARTICLES. USED IN ST015 AND ST035.

TIMER ST107 11111
TIMING AND VARIABLE PROGRAMMING CAPABILITY TO OPERATE AND CONTROL EXPERIMENT AND SUPPORTING EQUIPMENT. USED IN ST015 AND ST035.

WORK BENCH ST108 11111
PROVIDE A WORK AREA UPON WHICH EXPERIMENTS AND MAINTENANCE MAY BE PERFORMED. ENCLOSED STORAGE TO BE INCORPORATED IN BENCH. USED IN ST015, ST025, ST035, ST045 AND ST065.
SIZE APPROXIMATELY 2834 FT (86.1X9.9X1.2 M)

TAPE RECORDER, ANALOG ST000 11111
RECORD VIDEO, FM AND VIDEO DATA OF LANDMARK AND STAR TRACKING DATA. USED IN ST025, ST035 AND ST045.
BANDWIDTH 10 TO 20 MHZ
RECORDING FORMAT VIDEO, FM AND VIDEO CHANNELS 2

LASER ST116 11111
PULSE-LOCKED LASER RANGING SYSTEM MEASURING RANGE, LINE-OF-SIGHT ANGLES, AND RANGE RATE BETWEEN AIR/SHUTTLE AND SEPARATE SPACE VEHICLES INCLUDING RANGE OF AIR TO GROUND.
RANGE TOLERANCE +/- 1.18 IN (47.3 CM) AIR-GND
COOPERATIVE TARGET DISTANCE 300 NM (556 KM)

RECEIVER - RF ST117 11111
MEASURE ELECTROMAGNETIC INTERFERENCE AT ORBITAL ALTITUDES IN THE FREQUENCY SPECTRUM OF 400 MHZ TO 15 GHZ. USED IN ST025 AND ST035.
OVERALL SPECTRUM 400 MHZ TO 15 GHZ
SPECIFIC MULTIFATH MEASURE 1.5 TO 1.6 GHZ
2.065 TO 2.3 GHZ
12.4 TO 15.35 GHZ

TELESCOPE - POINTING ST154 11111
POINTING TELESCOPE WHICH ALLOWS EXPERIMENT PHENOMENA TO BE OBSERVED AND TRACKED IN SPACE

POWER SUPPLY ST171 11111
ELECTRICAL POWER REGULATION AND DISTRIBUTION AS REQUIRED AND AT VARIABLE VOLTAGE RANGES. USED IN ST025, ST035 AND ST045.
VOLTAGE 0 TO 24 VDC
POWER 200 WATTS
SENSITIVITY 0.05 VOLT

FREQUENCY OPTICAL ST273 11111
RECORD AND PLAYBACK DIGITALLY GENERATED DATA FROM A LASER SYSTEM. USED IN ST025 AND ST035.
BANDWIDTH 10 TO 50 MHZ
SIGNAL DIGITAL
CHANNELS 2

FREQUENCY METER ST156 11111
MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD, AVERAGE RATE AND MULTIPLE RATE OF VHF FREQUENCY SOURCES. USED IN ST035.
BANDWIDTH 10 TO 100 MHZ
SENSITIVITY +/- 2 dB

CAMERA, TV ST025 11111
REAL-TIME VIEWING OF TARGET AREAS, RANGE MARK AND STAR FIELD TRACKING, AND SENSER HOMESIGHTING. USED IN ST045, ST055 AND ST065.
VIEW ANGLE 8 TO 10 DEGREES
RESOLUTION 1025 LINES

MASS SPECTROMETER ST163 11111
ANALYSIS OF GASEOUS COMBUSTION PRODUCTS FORMED IN ZERO GRAVITY VIA MASS SPECTROSCOPY. USED IN ST045 AND ST065.
MASS RANGE 1 TO 1000 AMU
SCANNING RATE 80 PSEC TO 600 SEC
SENSITIVITY 1 N

LAS CHROMATOGRAPH ST175 11111
ANALYSIS OF COMBUSTION PRODUCTS FOR FLAME CHEMISTRY EXPERIMENTS. USED IN ST045 AND ST065.

PYROMETER ST176 11111
MEASURE TEMPERATURE OF COMBUSTION PROCESS IN FLAME CHEMISTRY EXPERIMENTS. USED IN ST045 AND ST065.
TEMPERATURE RANGE 0 TO 3000 DEG F

COMBUSTION CHAMBER ST178 11111
CHAMBER WHERE FUELS AND OXIDIZERS ARE INJECTED, ALLOWED TO BURN AND THE FLAME CAPABLE OF BEING MONITORED.
SIZE APPROXIMATELY 1.67 DIA X 3.5 FT
(0.51 DIA X 1.0 M)

TANKS - FUEL/OXIDIZER ST179 11111
SIMILAR TYPE TANKS CAPABLE OF STORING CORROSIVE FUELS AND OXIDIZERS IN THE SHUTTLE LABORATORY AT LOW PRESSURES. USED IN ST045 AND ST065.
SIZE APPROXIMATELY 0.75 DIA X 3 FT
(0.225 DIA X 0.9 M)

CHAMBER, CUBICAL TEST ST186 11111
CHAMBER FOR PERFORMANCE OF SUPERFLUID HELIUM AND LIQUID CRYSTAL PARTICLE EXPERIMENTS
DIMENSIONS 1.0 FT X 0.55 FT X 0.55 FT



ACOUSTIC SPEAKER	ST105	11111
PARTICLE POSITIONING DRIVER		
FREQUENCY RANGE..... 90 TO 1000 Hz		
MAGNET, SUPERCONDUCTING	ST190	11112
MAGNET FOR DROP/PARTICLE POSITIONING		
MOTOR, ELECTRIC	ST191	11111
MOTOR FOR CHAMBER DRIVE SYSTEM		
CRYOSTAT - HELIUM	ST192	11111
CEWAT TYPE STORAGE CONTAINER TO STORE AND DISPENSE LIQUID HELIUM.		
SIZE APPROXIMATELY 1.5 DIA X 1.7 FT (0.5 DIA X 0.05 M)		
PRESSURE MONITOR	ST193	11114
MEASURES PRESSURE IN TEST CHAMBER AND DISPLAYS DIGITALLY		
TEMPERATURE MONITOR	ST194	11114
MONITORS AND CONTROLS TEST CHAMBER TEMPERATURE		
COMBUSTION CHAMBER	ST202	11135
CHAMBER TO EXAMINE COMBUSTION IN ZERO GRAVITY		
DIMENSIONS..... 1.67 FT (0.48 M) DIA X 4.0 FT (1.22 M)		
ELECTRON SPIN RESONANCE/MASS SPECTROMETER	ST258	11112
IDENTIFY CONSTITUENTS OF COMBUSTION PRODUCTS		
TEMPERATURE SENSORS	ST260	11111
SENSOR TO MEASURE TEMPERATURE OF COMBUSTION PROCESS		
TEMPERATURE RANGE(EST)..... 1000 TO 3000 F		
PRESSURE SENSORS	ST260	11111
SENSORS TO MEASURE PRESSURE OF COMBUSTION EXPERIMENT		
PRESSURE RANGE..... TBD		
FUEL TANK	ST261	11135
STORE FUEL FOR COMBUSTION EXPERIMENT		
VOLUME..... 2 CU.FT. (0.055 CU.M.)		
OXIDIZER TANK	ST262	11135
STORE OXIDIZER FOR COMBUSTION REACTION		
VOLUME..... 2 CU.FT. (0.055 CU.M.)		
INERT GAS SUPPLY	ST263	11111
SUPPLY PRESSURE FOR EXPULSION OF FUEL AND OXIDIZER FROM STORAGE TANKS		
VOLUME 0.22 CU.FT. (0.006 CU.M.)		
STROBE LAMP	ST197	11111
VARIABLE DISCHARGING STROBE LIGHTS WHICH ALLOW THE STOPPING OF ACTION.		
STROBE DISCHARGE RATE 1 TO 1000 FLASHES/SEC		
CRITICAL STATE TEST CELL	ST207	11235
TEST CELL TO TEST CRITICAL FLAME PHENOMENA. TEST CELL INCLUDES A THERMAL SHIELD, HEAT EXCHANGER, MICROPHONE, CHAMBER, SUPPORT ELECTRONICS AND PRESSURE TRANSDUCERS		
ACOUSTIC MEASUREMENT DEVICE	ST214	11114
DEVICE CONSISTS OF AN ACOUSTIC CAVITY, AN OSCILLATOR AND A DIGITAL DISPLAY		
POOL BOILING TEST CHAMBER	ST218	21235
CHAMBER TO EVALUATE POOL BOILING PHENOMENA IN ZERO GRAVITY		
INSTRUMENTATION SENSOR SYSTEM	ST223	11114
SYSTEM INCLUDE TEMPERATURE, PRESSURE AND ACCELEROMETER SENSORS		
CRYSTALLIZER SAMPLE CHAMBER	ST225	11135
CHAMBER FOR GROWING CRYSTALS IN ZERO G		
SCHLIEREN TYPE OPTICS SYSTEM	ST227	11135
OPTICS SYSTEM EMPLOYING SCHLIEREN METHOD TO EXAMINE THE STRUCTURE OF CRYSTALS		
PHOTOGRAPHIC INTERFEROMETER	ST231	11235
DEVICE TO EXAMINE STRUCTURE OF GROWING CRYSTALS		
ACCELEROMETER, TRIAXIAL	ST239	21111

ORIGINAL PAGE IS
OF POOR QUALITY



COMMUNICATIONS AND NAVIGATION

TRACKING ANTENNA SERVO ELECTRONICS CN044 11114

ELECTRONICS FOR ANTENNA POSITIONING AND CONTROL

TYPE TYPE II SERVO LOOP
GAIN CONSTANT 3.984 PER SEC SQUARED
BANDWIDTH 100 HZ

RECEIVER, L-BAND CN045 21111

RECEIVE L-BAND SIGNALS TO DEMONSTRATE THE OPERATIONAL PERFORMANCE OF TYPICAL SPACEBORNE INTERFEROMETERS FOR AIRCRAFT AND MARINE NAVIGATION/TRAFFIC CONTROL. TRACK SIGNAL PHASE.

CARRIER FREQUENCY 1509.75 TO 1570.25 MHZ
SENSITIVITY -120 DBM
RECEIVER TYPE PHASE LOCK, SUPERHETERODYNE,
FIXED FREQUENCY, AGC
DYNAMIC RANGE 40 DB INSTANTANEOUS
100 DB TOTAL
PRE-DETECTION BW 5 MHZ
SWEEP RATE TUNE 3 PPM OF CARRIER FREQ
PHASE LOCK LOOP SNR 40 DB
RECEIVER NOISE FIGURE 5 DB
PHASE LOCK LOOP BANDWIDTH 10 HZ

RECEIVER, VHF CN046 51111

RECEIVE VHF SIGNALS FOR DETAILED EVALUATION OF PERFORMANCE CHARACTERISTICS OF TORS/GROUND/SHUTTLE DATA LINK. PROVIDES ALPHABET OF RECEIVER FRONT ENDS WHICH CAN BE CONNECTED IN DIFFERENT RECEIVE CONFIGURATIONS.

CARRIER FREQUENCY 120 TO 130 MHZ
SENSITIVITY +/- 25 MHZ
CARRIER/NOISE DENSITY RATIO 53 DBM
DATA RATE 100 TO 1000 BPS

RECEIVER, KU-BAND CN047 11235

RECEIVE KU-BAND SIGNAL FOR DETAILED EVALUATION OF PERFORMANCE CHARACTERISTICS OF TORS/GROUND/SHUTTLE DATA LINK.

CARRIER FREQUENCY 13.4 TO 14.2 GHZ
SENSITIVITY +/- 25 MHZ
CARRIER/NOISE DENSITY RATIO 56 DB/HZ
DATA RATE 100 TO 1000 BPS
DATA UP TO VIDEO

RECEIVER, S-BAND CN048 11111

RECEIVE S-BAND SIGNAL FOR DETAILED EVALUATION OF PERFORMANCE CHARACTERISTICS OF TORS/GROUND/SHUTTLE DATA LINK.

CARRIER FREQUENCY 2025 TO 2120 MHZ
SENSITIVITY +/- 25 MHZ
CARRIER/NOISE DENSITY RATIO 65 DB/HZ
DATA RATE 100 TO 1000 BPS
VIDEO

RECEIVER, LASER ELECTRONICS CN049 11235

PERFORMS DATA DEFORMATTING FROM LASER CARRIER BEAM

TRANSMITTER, VHF CN050 11111

PROVIDES MEDIUM POWER OUTPUT FOR TRANSMISSION TO TORS FOR EVALUATION OF PERFORMANCE CHARACTERISTICS.

CARRIER FREQUENCY 132 TO 144 MHZ
STABILITY +/- 20 KHZ
OUTPUT POWER 10 WATTS MAX
DATA RATE 100 TO 10,000 BPS
SIGNAL BANDWIDTH 0.1 MHZ
DUTY CYCLE CONTINUOUS

TRANSMITTER, KU-BAND CN051 11235

PROVIDE KU-BAND SIGNAL FOR EVALUATION OF PERFORMANCE CHARACTERISTICS OF TORS/GROUND/SHUTTLE DATA LINK.

CARRIER FREQUENCY 14.4 TO 15.35 GHZ
STABILITY +/- 25 MHZ
TRANSMITTER POWER 100 WATTS
EFFECTIVE RADIATED POWER 33 DBM
DATA RATE 1.0 MB
VIDEO
DUTY CYCLE CONTINUOUS
SIGNAL BANDWIDTH 500 MHZ

LASER ASSEMBLY - CO2 CN052 11111

PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFINE AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMMUNICATIONS APPLICATIONS. CONTAINS OPTICAL TRANSMITTER INCLUDING BEAMWIDTH CONTROL, MODULATOR, AND COOLING SYSTEM.

LASER CAVITY
WAVELENGTH 10.6 MICRONS
LASING MODE FUNDAMENTAL (TEM₀₀)
OPTICAL OUTPUT POWER 5 WATTS AVERAGE
EFFICIENCY 5 PERCENT
MODULATOR
TYPE DIGITAL
DATA RATE 200 MBPS
EXTINCTION RATIO 40 DB
COOLING SYSTEM
COOLING METHOD CONDUCTIVE, NO LIQUIDS ALLOWED
HEAD LOAD 200 WATTS

LASER ASSEMBLY - NEODYM CN053 11111

PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFINE AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMMUNICATIONS APPLICATIONS. CONTAINS OPTICAL TRANSMITTER INCLUDING BEAMWIDTH CONTROL, MODULATOR, AND COOLING SYSTEM.

LASER CAVITY
WAVELENGTH 1.06 MICRONS
LASING MODE FUNDAMENTAL (TEM₀₀)
OPTICAL OUTPUT POWER 1 WATT AVERAGE
EFFICIENCY 1 PERCENT
MODULATOR
TYPE DIGITAL
DATA RATE 1 MBPS
EXTINCTION RATIO 40 DB
COOLING SYSTEM
COOLING METHOD CONDUCTIVE, NO LIQUIDS ALLOWED
HEAD LOAD 200 WATTS

LASER BEAM, NEODYM CN054 11111

PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFINE AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMMUNICATIONS APPLICATIONS. CONTAINS OPTICAL TRANSMITTER INCLUDING BEAMWIDTH CONTROL, MODULATOR, MODULATED, AND COOLING SYSTEM.

LASER CAVITY
WAVELENGTH 1.06 MICRONS
LASING MODE FUNDAMENTAL (TEM₀₀)
OPTICAL OUTPUT POWER 0.1 WATTS AVERAGE
EFFICIENCY 0.1 PERCENT
MODULATOR
TYPE DIGITAL
DATA RATE 1 MBPS
EXTINCTION RATIO 40 DB
COOLING SYSTEM
COOLING METHOD CONDUCTIVE, NO LIQUIDS ALLOWED
HEAD LOAD 200 WATTS

LASER BEAM, NEODYM CN055 11111

SERVE AS A LIGHT SOURCE TO FACILITATE ACQUISITION AND/OR TRACKING FROM A DISTANT COMMUNICATION TERMINAL.

BEAM TYPE G-SWITCHED COOLED, MULTIMODE
CAVITY
WAVELENGTH 0.53 MICRONS (PREFERRED)
REP RATE 5 PPS
ENERGY/PULSE 20 MJ AT 0.53 MICRONS
PULSEWIDTH 20 PSEC
PUMP METHOD FLASH PUMP
COOLING TECHNIQUE FLUID POSITION FLUORESCENCE
BEAMWIDTH VARIABLE 1 DEG TO 1 ANG MINUTE
LIFETIME 10 MEG FLASHES OR GREATER
POWER CONSUMPTION 100 WATTS

FREQUENCY SYNTHESIZER AND DRIVER CN056 11111

RECEIVER FIRST LOCAL OSCILLATOR

FREQUENCY 0.1 TO 500 MHZ
STABILITY +/- 3E-11 PER DEG C
OUTPUT VOLTAGE 1 V RMS +/- 1.5 DB

NOISE FIGURE TEST SET CN057 11111

MEASURE SYSTEM NOISE FIGURE

NOISE FIGURE RANGE 10 TO 18 DB
ACCURACY +/- 0.5 DB, G - 30 DB
INPUT FREQUENCY 100 TO 1000 MHZ
BANDWIDTH 1 MHZ

SCENE CAMERA CN058 21111

OBTAIN FILM RECORD OF DATA AND DISPLAYS.

RECOLLECTION RATIO 110:1 TO 7
LENS F/1.9 TO F/16
SHUTTER 8 TO 5 SECONDS
FILM TYPE PELLARCID

FREQUENCY COUNTER CN059 11111

MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE, TIME INTERVAL AND RATIO OF FREQUENCY SOURCES.

BANDWIDTH 1 HZ TO 10 MHZ
SENSITIVITY -100 DBM
TIME BASE INTERNATIONAL SECOND

RF POWER METER CN060 11111

MEASURE RF POWER OF VARIOUS RF SOURCES.

FREQUENCY RANGES 120 TO 130 MHZ
130 TO 144 MHZ
2025 TO 2120 MHZ
13.4 TO 14.2 GHZ
14.4 TO 15.35 GHZ
ACCURACY 0.02 DB/10 DB

AC/DC VOLTMETER CN061 21111

MEASURE THE VOLTAGE LEVEL OF VARIOUS SIGNAL SOURCES.

VOLTAGE RANGE +/- 100 MV; 1, 10, 100, 1000 V
FREQUENCY RANGE 45 HZ TO 1 MHZ
ACCURACY +/- 0.01 % OVERALL

A/D CONVERTER CN062 21111

CONVERT INPUT SIGNAL FROM ANALOG TO DIGITAL FORMAT FOR 120 HZ

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

CAMERA, CINE CNO63 11111
PROVIDE PHOTOGRAPHIC RECORD OF TELESCOPE FIELD-OF-VIEW.

FILM FORMAT 35 MM
OPERATION REMOTE
FRAME RATE 16 FPS
FILM TYPE CCLGR
RESOLUTION 56 LINES/MM
LENS 100 MM
FIELD OF VIEW 40 DEGREES

OPTICAL ANTENNA SERVO ELECTRONICS CNO64 11114

SUPPORT ELECTRONICS FOR OPTICAL ANTENNA

SWEPTBAND RECEIVER AND DEMODULATOR CNO67 41111

RECEIVE SIGNALS FROM TERRESTRIAL NOISE AND TERRESTRIAL INTERFERENCE SOURCES AS A FUNCTION OF TIME OF DAY AND SEASON OF YEAR.

FREQUENCY RANGE 100 TO 1000 MHZ SWEEP MODE OR
TUNABLE
RECEIVER TYPE SUPERHETERODYNE, LINEAR, SWEPT
FREQUENCY, AGC
SENSITIVITY -120 DBM
PRE-DETECTION BW ADJUSTABLE
DYNAMIC RANGE 60 DB INSTANTANEOUS
120 DB TOTAL
SWEEP RATE 120 DB TOTAL
RECEIVER NOISE FACTOR 4.0 DB
BANDWIDTH 1 KHZ

ATTENUATOR CALIBRATOR CNO68 41111

DYNAMIC RANGE ADJUSTMENT, REMOTE ELECTRICAL CONTROL, DIGITAL READ-OUT OF VALUE, MANUAL OVER-RIDE.

IMPEDANCE 50 OHMS
FREQUENCY DC TO 1.5 MHZ
RANGE 0 TO 40 DB

SCAN PROGRAM GENERATOR CNO69 31111

GENERATES DIGITAL CONTROL SIGNALS TO CONTROL RECEIVER SWEEP.
SUPPLIES DIGITAL READOUT OF FREQUENCY AND LEVEL.

RANGES 24
SCAN TIME 0.03 TO 300 SECONDS
CAPABILITY ANALOG DATA COLLECT & DISPLAY

SIGNAL FORMATTER CNO70 21111

PROVIDES DATA FORMATTING, ACCUMULATES PRE/POST AMBLE, CLOCKS
START OF EACH TEST AND DATA ACQUISITION.

DATA FORMAT 8, 16, 32, BITS/WORD
DATA RATE 125000 WORDS PER SECOND

RF VARIABLE POWER SUPPLY CNO71 21111

PROVIDES SOURCE OF VARIABLE POWER.

OUTPUT 100 V PEAK
VOLTAGE GAIN 20 X
DISTORTION 0.1 & THD
FREQUENCY RANGE DC TO 30 KHZ

POWER CALIBRATION UNIT CNO72 21111

PROVIDES KNOWN INCREMENTALLY CONTROLLABLE POWER LEVEL FOR
STANDARDIZATION.

OUTPUT -100 TO +100 V
LOAD REGULATION 1 MV + 0.01 %
LINE REGULATION 10 MV + 0.01 %

TRACKER ELECTRONICS - FINE CNO73 21114

PROVIDES ACCURATE TRACKING OF AN OPTICAL BEACON, PROVIDES CONTROL
SIGNALS TO THE TELESCOPE GIMBALS AND VERNIER BEAM DEFLECTOR.

TYPE IMAGE DISSECTOR TRACKER
FOV +/- 300 MICRORAD
ACCURACY 0.5 MICRORAD
SENSITIVITY TO 1.06 MICRONS OR 10.6 MICRONS

BIT ERROR COUNTER CNO74 11111

COUNT ERROR PULSES OUTPUT BY THE BIT ERROR DETECTOR TO DETER-
MINE THE BIT ERROR RATE OF THE DIGITAL TRANSMISSION SYSTEM.

SWITCH/DUPLEX/PREAMP UNIT CNO77 11134

FOR SWITCHING ANTENNAS FOR DIFFERENT ALTITUDES, USE OF RF
CANCELLATION, FORMING ORNI PATTERN ETC

MODEM CNO78 11114

MODULATE AND DEMODULATE VOICE OR DATA TRANSMISSIONS WITH OR WITH-
OUT SPREAD SPECTRUM CAPABILITY. PROVIDES BASELINE COMMUNICATION
SYSTEM, ANALOG OR DIGITAL VOICE, AND DIGITAL DATA.

DATA 75 BS TO 50 KBS (VARIABLE)
MODEABLE SS
BANDWIDTH 2 MHZ

DEMODULATOR CNO79 31111

DEMODULATE VOICE OR DATA TRANSMISSIONS WITH OR WITHOUT SPREAD
SPECTRUM CAPABILITY. PROVIDES A NUMBER OF RECEIVERS CONNECTED TO
DIFFERENT ANTENNAS.

MODEM, WIDEBAND CNO80 21231

MODULATE AND DEMODULATE WIDEBAND FM OR DIGITAL SIGNALS. WIDEBAND
FM DETECTOR WITH THRESHOLD EXTENSION, DIGITAL DEMOD WITH BIT,
SYNCH, AND ERROR CORRECTION CODING. MODULATOR PORTION HAS COMPLI-
MENTARY CHARACTERISTICS.

D/A AND A/D CONVERTER CNO81 11111

PROVIDES FOR A/D AND D/A CONVERSION.

A/D CONVERTER
SAMPLING RATES 15 MEGASAMPLES/SEC
DATA RATES 105 MBPS
D/A CONVERTER
BANDWIDTH 6 MHZ

CONTROL UNIT - ANTENNA SCAN CNO82 11114

PROVIDES VARIABLE ANGULAR SWEEP PATTERNS (SPIRAL RECTANGULAR
SCANS, ETC.) OF VARIABLE ANGULAR COVERAGE.

DATA BIT STREAM GENERATOR CNO83 11111

PROVIDES TEST BIT STREAMS FOR TRANSMISSION.

DATA RATES 1 TO 105 MBPS

COARSE TRACKER CNO84 11114

PROVIDES INITIAL GROUND BEACON ACQUISITION AND COARSE TRACKING.

DYNAMIC (SEARCH) FCV 1 DEG BY 1 DEG
ACCURACY (1 SIGMA) +/- 200 MICRORAD
TRACKING BANDWIDTH 1 MHZ
SEARCH TIME 10 SEC MAX
SENSITIVITY TO .53 MICRONS
FALSE TARGET REJECTION MUST REJECT MOON, PLANETS,
STARS, AND SUNLET CLOUES.

LASER TRANSMITTER ELECTRONICS CNO85 11235

DATA FORMATTING FOR TRANSMISSION VIA LASER CARRIER

OPTICAL COLLIMATOR CNO86 11111

PROVIDES FOR OPTICAL ALIGNMENT OF THE LASER AND TELESCOPE. ALIGN-
MENT CHECKS TO BE BOTH CONTINUOUS AND PULSED.

LASER POWER SUPPLY CNO87 21111

CONVERTS AND REGULATES RAW SPACECRAFT POWER IN A FASHION SUITABLE
FOR PUMPING A LASER. PRIOR LASER POWER CONDITIONER AND DISTRIBUTION
TICKER FOR ALL LASERS.

VOLTAGE OUTPUT
REGULATION +/- 0.1 %
POWER 100 WATTS MAX
TRANSIENT RECOVERY 0.1 SECONDS

BEAM EXPANDER OPTICS CNO88 21111

PROVIDES VARIABLE BEAM DIVERGENCE CONTROL FOR LASER TRANSMITTER
AND THE OPTICAL BEACON TO FACILITATE ACQUISITION.

TYPE AFCCAL WITH VARIABLE ELEMENT
SEPARATION.
INPUT BEAM SIZE 0.25 INCH DIA
EXIT BEAM SIZE 1 INCH DIA
INPUT BEAM CHARACTERISTICS
LASER TRANSMITTER TEM WAVE, GAUSSIAN DIST
OPTICAL BEACON 3 MMAD DIVERGENCE, GAUSSIAN
DIST

OUTPUT BEAM DIVERGENCE
CO2 LASER TRANSMITTER DIFF LIMIT TO 200 ARC-SEC
NDYAG LASER TRANSMITTER DIFF LIMIT TO 50 ARC-SEC
OPTICAL BEACON 2 ARC-PIN TO 1 DEGREE
ADJUSTMENT RESPONSE TIME 0.5 SEC MAX

BEAM DEFLECTOR CNO89 41111

ONE TYPE PROVIDES A VERNIER POINTING CONTROL WITH FASTER RESPONSE
AND GREATER PRECISION THAN THE MAIN OPTICS GIMBALS. THE OTHER
TYPE PROVIDES A SMALL OFFSET POINTING CAPABILITY BETWEEN THE
TRACKER AXIS AND THE TRANSMITTER AXIS.

DEFLECTION CAPABILITY 160 AXIS
DEFLECTION RANGE +/- 200 MICRORAD
PRECISION +/- 0.2 MICRORAD
SPECTRAL RANGE ACCURATE LASER TRANSMITTER,
RECEIVER AND BEACON.

LASER POWER METER CNO90 11111

MEASURE OUTPUT LASER POWER.

WAVELENGTH 0.53, 1.06 AND 10.6 MICRONS
ACCURACY +/- 5 PERCENT
RANGE 300 MW

CALIBRATION SIGNAL GENERATOR CNO91 11111

CALIBRATES RECEIVER SENSITIVITY AND PHASE METER.

FREQUENCY RANGE 1565.75 TO 1570.25 MHZ
FREQUENCY ACCURACY +/- 5 MHZ
FREQUENCY STABILITY 0.005 H/DEGREE C
RF OUTPUT POWER +10 DB TO -127 DB

RFEN H.F. AMPLIFIER 1-6.1 CNO92 51111

RECEIVER - NAVIGATION CNO96 11111

RECEIVE SIGNALS FROM DIRECT OPTICAL SENSING OF STARS TO INVESTI-
GATE THE PRECISION WITH WHICH A SATELLITE CAN BE PLACED VERTIC-



Space Division
Rockwell International

ALLY USING STARS AS THE SOLE INERTIAL REFERENCE.

CARRIER FREQUENCY 400 MHZ
SENSITIVITY 200 MV FOR 1C CH S+R/A
DATA RATE 2000 BPS

SPECTRUM ANALYZER/OSCILLOSCOPE CN900 21111

SINGLE DESIGN FUNCTIONS AS POWER SPECTRAL DENSITY AND WAVEFORM
DISPLAY FOR QUICKLOOK AND AS MODULATION MONITOR.

BANDWIDTH 10 G-HZ
SENSITIVITY 2 MV/DIVISION
CHANNELS 4
PHOTOGRAPHIC CAPABILITY YES

COMPUTER, GENERAL PURPOSE CN901 21111

COMPUTER CAPABILITY TO COMMAND EXPERIMENTS ON AND OFF, PROGRAM
EXECUTION OF EXPERIMENT PROCEDURAL STEPS, COMPUTE GIMBAL
ANGLES FOR TELESCOPE AND ANTENNA POINTING, AND COMPUTATION OF
POINT-TO-POINT ANGLES.

TIMER, PRECISION CLOCK CN902 11111

PROVIDES PRECISE TIME REFERENCE. MEASURE TIME AT WHICH PHASE PEAKS
OCCUR.

DISPLAY 24 BIT DIGITAL CLOCK WITH 1 HR.
TIME PERIOD
DISPLAY TYPE DIGITAL
ACCURACY 1E-60

TELEPRINTER CN903 11111

PROVIDES CAPABILITY TO SEND OR RECEIVE TELETYPE MESSAGES.

CHARACTER RATE 30 CHARACTERS/SEC (150 EPS)
CONSISTS OF SOLID STATE KEYBOARD
SOLID STATE THERMAL PRINTER
EXPANDABLE MESSAGE MEMORY

PHASE METER, DIGITAL CN904 11111

MEASURES RELATIVE PHASE ANGLES OF RECEIVED SIGNALS FOR INTERFER-
OMETRY NAVIGATION/SURVEILLANCE EXPERIMENT.

OUTPUT GENERATES 13 BIT PHASE PEAKS
BASED ON THE 00 PHZ SINUSOIDS
AT 10 DB.
FREQUENCY RANGE 400 MHZ
RESOLUTION 25 KHZ
ACCURACY +/- 2.5 KHZ

RECORDER, DIGITAL CN905 11111

ACCURATE RECORDING OF RAW DIGITAL DATA

SIGNAL TYPE 0 TO 5 VDC DIGITAL
CHANNELS 9
TAPE SPEED VARIABLE
DATA RATE 1.9E+000BPS

POWER DIVIDER, WIDEBAND CN906 41111

DIVIDES ANTENNA POWER OUTPUT.

INPUT PORTS ONE 50 OHM PORT
OUTPUT PORTS THREE 50 OHM PORTS

X-Y PLOTTER CN907 11111

PROVIDES REAL-TIME ACCURATE REPRODUCTION OF SPECTRAL DATA SOURCES
PLOTTED ON CARTESIAN COORDINATES.

CHANNELS 4 PLUG INS
ACCELERATION 1500 IN/S/S (3000 CM/S/S)
SENSITIVITY 0.5MV/IN (0.25 MV/CM)

KEYBOARD, COMPUTER CN908 11111

PROVIDES FOR CREWMAN TO COMMUNICATE WITH COMPUTER FOR
EXPERIMENT CONTROL AND DATA ANALYSIS.

KEYBOARD TYPEWRITER TYPE

CAMERA, TELEVISION CN909 11111

PROVIDES ASTRONAUT WITH VIEW OF EARTH SCENE.

VIEW ANGLE 20 DEGREES
RESOLUTION 625 LINES
SPECTRAL BANDWIDTH 4 PHZ
SHOOT RATE 15 FRAMES PER SECOND

RECORDER, VIDEO CN910 11111

PROVIDES FOR RECORDING VIDEO DATA.

CHANNELS 2
BANDWIDTH 4 KHZ MINIMUM
DURATION 1.5 HRS MINIMUM

RECEIVER, PHASE-LOCK CN911 11111

PROVIDES FOR TRACKING SIGNAL PHASE.

FREQUENCY SPECTRUM 1570 PHZ TUNABLE
SENSITIVITY 15 KHZ

DISPLAY, CRT CN912 11111

DISPLAY REAL-TIME ON VIDEO TAPE DATA.

RESOLUTION 625 LINES
FRAME RATE 15 FPS

VIDEOCAMERA

CN914 11111

REAL-TIME RAW DATA RECORDING.

BANDWIDTH 0.1 TO 5000 KHZ
CHANNELS 16
RECORDING FORMAT CRT FRAME COPY

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

APPENDIX B

EXPERIMENT EQUIPMENT FUNCTIONAL REQUIREMENTS LISTED BY APPENDIX TYPE

001000 ACCELEROMETER - SPECIMEN MOVEMENT LS901 121111
 001000
 001000 ORGANISM CONTAINER ACTIVITY AND ACTIVITY LEVEL METER.
 001000
 001000 ACCURACY +/- 0.0004 G
 001000 FREQUENCY RESPONSE 0 TO 100 HZ
 001000
 001000 ACCELEROMETER, TRIAXIAL ST235 21111
 001000
 001110 ACOUSTICAL GENERATOR EC013 11111
 001110
 001110 VARIABLE FREQUENCY AUDIO DRIVERS USED TO POSITION AND CONTROL
 001110 WATER DROPLETS OR ICE PARTICLES IN EACH OF THREE MUTUALLY
 001110 PERPENDICULAR AXES. EACH DRIVER CAN BE CONTROLLED INDEPENDENTLY
 001110 IN FREQUENCY AND AMPLITUDE
 001110
 001110 FREQUENCY RANGE..... C TO 100KHZ
 001110 AMPLITUDE..... C TO TBD DB
 001110
 001110 ACOUSTIC MIXING AND DISPERSAL UNIT SF002 21235
 001110
 001110 ACOUSTIC SPEAKER ST189 61111
 001110
 001110 PARTICLE POSITIONING DRIVER
 001110
 001110 FREQUENCY RANGE..... 50 TO 100KHZ
 001110
 001210 TRANSDUCER - ACOUSTICAL EC013 241111
 001210
 001210 MEASURE SOUND INTENSITY LEVEL
 001210
 001210 INTENSITY LEVEL..... C TO TBD DB
 001210
 001210 ACOUSTIC TRANSDUCER AND DETECTOR. SF004 11111
 001210
 001210 USE OF SOUND PRESSURE WAVES TO POSITION EXPERIMENT SPECIMEN
 001210 SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL.
 001210
 001210 CHAMBER PRESSURE > 0.14 PSI
 001210 (1 HP/50 PSI)
 001210 < 2E-04 PSI
 001210 (0.01 N/50 PSI)
 001210 CONFIGURATION IN APEX OF TETRAHEDRON
 001210
 001210 ALTIMETER LS902 11111
 001210
 001210 FREQUENCY SOURCE ADJUST IN BOTH FREQUENCY AND AMPLITUDE FOR
 001210 EARPHONES OR PHONE VIBRATORS
 001210
 001210 FREQUENCY..... 125 TO 1000 ADJUSTABLE IN OCTAVES
 001210 ATTENUATOR..... CONTINUOUSLY ADJUSTABLE IN 100 STEPS
 001210
 001210 SOUND LEVEL METER LS923 11111
 001210
 001210 SOUND LEVEL SURVEYS AND LIMITED FREQUENCY ANALYSIS IN REAL TIME.
 001210
 001210 FREQUENCY RESPONSE 125 TO 2000 HZ
 001210 SOUND LEVEL RANGE 40 TO 100 DB
 001210
 001210 ACOUSTIC MEASUREMENT DEVICE ST214 11114
 001210
 001210 DEVICE CONSISTS OF AN ACOUSTIC CAVITY, AN OSCILLATOR AND A
 001210 DIGITAL DISPLAY
 001210
 001300 BANDPASS FILTER - WAVE ANALYSIS AF912 21111
 001300
 001300 DETERMINATION OF FREQUENCY AND AMPLITUDE INFORMATION ABOUT THE
 001300 CARRIER AND SIDEBANDS.
 001300
 001300 FREQUENCY RANGE 300 HZ TO 20.00KHZ
 001300 ATTENUATION 75 DB
 001300 ACCURACY +/- 2 %
 001300
 001400 DROPLET CLOUD AEROSOL GENERATORS EC008 11135
 001400
 001400 GENERATE AEROSOLS FOR EXPERIMENTATION
 001400
 001400 PARTICLE CONCENTRATION..... 10 UP TO 10000/CC
 001400 CLOUD FORM..... POLYDISPERSED AND MONODISPERSED
 001400 PARTICLE SIZE..... 1E-4 TO 4E-1 CM
 001400 PRODUCTION RATE..... UP TO 1005 HZ
 001400
 001400 LARGE AND GIANT NUCLEI AEROSOL GENERATOR EC008 11135
 001400
 001400 SAME AS DROPLET CLOUD AEROSOL GENERATOR
 001400
 001400 AITKEN AEROSOL GENERATOR EC008 11135
 001400
 001400 SAME AS DROPLET CLOUD AEROSOL GENERATOR
 001400
 001400 GENERATOR - AEROSOL ST046 11135
 001400
 001400 CAPABILITY TO PRODUCE FINE LIQUID PARTICLES REQUIRED FOR VARIOUS
 001400 EXPERIMENTS. USED IN STC15, STC25 AND STC35.
 001400
 001400 AIR PARTICLE SAMPLE COLLECTOR LS929 11111
 001400
 001400 AIR PARTICLE AND MICRO-ORGANISM SAMPLES FOR AIR QUALITY DETERMIN-
 001400 ATION
 001400
 001400 AIRFLOW 1 CU FT/MIN
 001400 SAMPLING TIME - CONTINUOUS ... 1 HOUR
 001400
 001400 AIR SAMPLE UNIT ST105 21111
 001400
 001400 PERIODIC COLLECTION OF AIR SAMPLES FROM THE SHUTTLE AND ATL
 001400 PROVIDING TYPE, RATE OF CHANGE AND QUALIFICATION OF MICRO-

ORIGINAL PAGE IS
OF POOR QUALITY


004100 ORGANISMS AND NON-VIABLE PARTICLES. USED IN STUDIES AND TESTS.
004100
004400 GAS ELIMINATION SYSTEM 54001 11114
004400
004400 REMOVES GASES PRODUCED BY ELECTROLYSIS IN THE ELECTROLYTIC
004400 COMPARTMENTS OF THE ELECTROCHEMICAL COLUMNS
004400
004400 OUTLET OXYGEN CONCENTRATION... 10 LPM VOLUME
004400 OUTLET HYDROGEN CONCENTRATION... 10 LPM VOLUME
004400 MAX FLOW RATE..... 100
004400 FLOW RATE CONSTANCY..... 4-12
004400 OPERATING TEMPERATURE..... 10 - 77 F (20 - 24 C)
004400
004400 MOLECULAR SIEVE 54001 11114
004400
004400 THE MOLECULAR SIEVE IS USED TO CLEAN GASES GENERATED IN FURNACES
004400 AND EXHAUSTED TO SPACE. THIS UNIT IS ALSO USED IN THE FURNACE.
004400 GENERAL PURPOSE AND LEVITATION SUPPLEMENT
004400
004400 CATALYTIC OXIDIZER SYSTEM 55977 11235
004400
004400 REMOVES UNDESIRABLE COMBUSTION PRODUCTS
004400
004400 OPERATING TEMPERATURE..... 700F (370C)
004400 AIR FLOW RATE..... 20CFM AT 14.7PSIA AND 70F
004400
006170 MICROPHONE AMPLIFIER 59100 11111
006170
006170 AMPLIFY SIGNALS FROM AMPLIFIER
006170
006180 AMPLIFIERS - WAVE ANALYSIS 59415 11111
006180
006180 AMPLIFIER - ACOUSTIC GENERATOR 59013 11111
006180
006180 AMPLIFY SIGNAL TO ACOUSTIC GENERATOR
006180
006180 POWER..... 45W
006180
006180 GREEN R.F. AMPLIFIER TANK 59002 11111
006180
008590 ANESTHETIZER - INVERTER/AMPLIFIER 59492 11114
008590
008590 RENDER INVERTEBRATE ORGANISMS INSENSIBLE TO FACILITATE HANDLING.
008590
008590 ANESTHETIC GAS CARBON DIOXIDE
008590
018700 RECEIVER, TELESCOPE 59094 11114
018700
018700 RECEIVES TRANSMITTER RETURNED SIGNALS.
018700
018700 INSTANTANEOUS FLD 0.0017 RAD
018700 POINTING ACCURACY 0.0007 RAD
018700 POINTING STABILITY 0.0007 RAD
018700 POINTING STABILITY RATE 0.0017 RAD PER SEC
018700 POINTING DURATION 7000 SEC
018700
018700 TELESCOPE - TRACKING 59111 11114
018700
018700 USED FOR HIGH RESOLUTION VIEWING OF SPECIFIC TARGETS
018700
018700 INSTANTANEOUS FIELD OF VIEW..... 0.0017 RAD
018700 0.0007 RAD
018700 FIELD OF VIEW (CROSS-TRACK)..... 0.0017 RAD AND 0.0007 RAD
018700 POINTING ANGLE FROM NADIR
018700 (PITCH, ALONG TRACK)..... 0.0017 RAD
018700 POINTING ANGLE FROM NADIR
018700 (ROLL, CROSS TRACK)..... 0.0017 RAD
018700 TOTAL ANGULAR COVERAGE (CROSS-
018700 TRACK) FROM NADIR (11-12) RAD
018700 POINTING ANGLE..... 0.0017 RAD
018700 POINTING ACCURACY (1-SIGMA)..... 0.0017 RAD
018700
018700 TELESCOPE - POINTING 59112 11114
018700
018700 POINTING TELESCOPE WHICH ALLOWS EXPERIMENT PHENOMENA TO BE
018700 OBSERVED AND TRACKED IN SPACE
018700
023500 MASS MEASUREMENT DEVICE - MICRO 59497 11131
023500
023500 TOTAL MICRO MASS DETERMINATION BY OSCILLATING FREQUENCY.
023500
023500 MASS 125 TO 250 LB (57-113 KG)
023500 (50 TO 100 KG)
023500
023500 MASS MEASUREMENT DEVICE, MACRO 59493 11132
023500
023500 MEASURE MASS OF ITEMS SUCH AS BAGS OF FOOD, RESEARCH CONTAINERS,
023500 CONTAINERS OF URINE AND FECES, LARGE SPECIMENS, ETC
023500
023500 MASS RANGE..... 0.001 TO 100 LB (0.1 TO 45 KG)
023500
023500 MASS MEASUREMENT DEVICE, MICRO 59494 11132
023500
023500 MEASURE SMALL TEST SPECIMENS
023500
023500 MASS MEASUREMENT RANGE..... 100 TO 1000 LB
023500
034920 ALDIO STEWART PERSET 59495 11111
034920
034920 EARPHONES FOR VARIOUS MSI HEARING TESTS
034920
034920 CONSOLE, BEHAVIORAL MEASUREMENTS 59496 11114
034920
034920 DEVICE THAT MEASURES MENTAL PREFERENCES, ATTITUDES, ETC. SUCH AS
034920 COGNITIVE PROFILES, INDIVIDUAL BEHAVIORAL TRAITS, GROUP
034920 BEHAVIORAL TRAITS, ETC. CONSOLE INCLUDES GENERAL FLUORESC DISPLAY

ORIGINAL PAGE IS
OF POOR QUALITY



03947C AND GENERAL PURPOSE RESPONSE KEYCARD
03947C
03947C PROCESSES MEASURED.....COMPLEX PERCEPTUAL
03947CCONCEPTUAL AND THINKING ABILITY
03947CMEMORY
03947CREACTION TIME-SINGLE AND COMPLEX
03947CINDIVIDUAL BEHAVIORAL TRAITS
03947CGROUP BEHAVIORAL TRAITS
039575 LIQUID DROP GENERATOR EEC6 41235
039575
039575 GENERATES SINGLE DROPLETS BY MANUAL OR SYSTEM CONTROL. DROPLETS
039575 ARE PROPELLED INTO THE EXPERIMENT CHAMBERS FOR EXAMINATION
039575
039575 FORMATION FREQUENCY(PA).....10/SEC
039575 DROPLET DIAMETER.....10 MICROMETERS TO 1 CM
039575 DROPLET UNIFORMITY.....WITHIN 5 %
039575
039580 ICE PARTICLE GENERATOR EEC6 31235
039580
039580 GENERATE ICE PARTICLES FOR TEST IN ICE DIFFUSION, EXPANSION AND
039580 GENERAL CHAMBERS
039580
039580 ICE CRYSTAL SIZE.....50 MICROMETERS TO 1 CM
039580
039640 VISION TESTER LS981 11235
039640
039640 OPTICAL DEVICE (WITH FOOD AND RESPONSE KEYCARD) THAT MEASURES
039640 A VARIETY OF VISUAL FUNCTIONS.
039640
042000 CELL COUNTER - BLOOD LS907 11111
042000
042000 MEASURE BLOOD CELL PROPERTIES OF BLOOD SAMPLES. ESSENTIALLY A
042000 COUNTER COUNTER
042000
042000 PARAMETERS MEASUREDHEMOGLOBIN,HEMATOCRIT,RED BLOOD
042000CELL COUNT,WHITE BLOOD CELL
042000COUNT,MEAN CELL VOLUME,MEAN
042000CELL HEMOGLOBIN AND MEAN CELL
042000HEMOGLOBIN CONCENTRATION
042000 REPRODUCIBILITY.....+/-1%
042000
042000 AUTANALYZER - MULTIPLE LS930 11111
042000
042000 AUTOMATIC ANALYSIS OF BLOOD AND APPROXIMATELY 30 OTHER FLUIDS.
042000
042000 COUNTER PRINTER-COUNTER ST084 11111
042000
042000 PROVIDE AUTOMATIC BLOOD ANALYSIS (HEMOGLOBIN, HEMATOCRIT, RED
042000 BLOOD CELL COUNT, WHITE BLOOD CELL COUNT, MEAN CELL VOLUME, MEAN
042000 CELL HEMOGLOBIN AND ITS CONCENTRATION). USED IN ST015, ST025
042000 AND ST035.
042000
043000 BLOOD CLOT FIBROMETER LS915 11111
043000
043000 PROVIDE AUTOMATIC MEASUREMENT OF PLASMA COAGULATION TIME
043000
043000 ACCURACY.....+/-0.1SECCND
043000
043500 CARDIOGRAPH - IMPEDANCE LS100 11111
043500
043500 BLOOD FLOW RATE FROM HEART.
043500
043500 BLOOD FLOW RATE 2 TO 25 LITERS/MIN +/-5 PCT FS
043500
043500 PULMONARY FLOWMETER LS990 11111
043500
043500 BLOOD FLOW RATES AND VOLUME.
043500
043500 FLOW RATE 0 TO 15 L/MIN, +/-3 PCT FS
043500 VOLUME 0 TO 7 L, +/-3 PCT FS
043500
043500 TRANSCUTANEOUS DOPPLER FLOWMETER LS991 11112
043500
043500 MEASURE PULSE WAVE COUNTS.
043500
043500 ULTRASONIC SIGNAL 1 AND 5 MHZ
043500
048000 PULSE WAVE LS112 11112
048000
048000 MEASURE BLOOD VELOCITY
048000
048000 VELOCITY 3 TO 15 M/S, +/-0.1 M/S
048000
048050 BLOOD PRESSURE CUFF SYSTEM LS978 11132
048050
048050 PRESSURE MEASURING TRANSDUCER CUFF USED TO MEASURE BLOOD PRESSURE
048050 OF PRIMATES. SYSTEM INCLUDES AIR PUMP FOR AUTOMATED MEASURING.
048050 PRESSURE CUP AND TRANSDUCER, AND SIGNAL CONDITIONER.
048050
048700 REFRIGERATOR SF001 11111
048700
048700 THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES
048700
048700 TEMPERATURE.....22 - 80F (0 - 20C)
048700 LENGTH.....2FT (60CM)
048700 WIDTH.....1.5FT (46CM)
048700 HEIGHT.....2.0FT (61CM)
048700
048700 REFRIGERATOR ST080 11111
048700
048700 PROVIDE A COLD STORAGE AREA TO HOLD SPECIMEN. USED IN ST015,
048700 ST025 AND ST035.
048700
048700 TEMPERATURE 3 C (27K K)
048700 SIZE APPROXIMATELY 3 X 2 X 1 FT (76,2X50,7X25,4CM)

ORIGINAL PAGE IS
OF POOR QUALITY



SD 44-SA-0047-2

055010
055010 FILM FORMAT 35 MM
055010 LENS 58 MM, F1.2
055010 VIEW ANGLE 41 DEGREES
055010 LENS 35 MM, F1.8
055010 VIEW ANGLE 64 DEGREES
055010 SHUTTER SPEED 8 TO 1/500 SEC
055010 VIEWING THRU LENS
055010 FOCUS SPLIT IMAGE
055010
055010 CAMERA - STILL POLAROID (CORE) LSC02 41111
055010
055010 POLAROID PHOTOGRAPHIC CAPABILITY SHARED WITH MOST EXPERIMENTS
055010 ALSO ADAPTABLE WITH MICROSCOPE AND OSCILLOSCOPE INTERFACES.
055010
055010 FILM FORMAT POLAROID 107 BLACK/WHITE
055010 2.25X4.25 IN (57.15X108.00 MM)
055010 LENS 58 MM, F3.5 TO F22
055010 SHUTTER SPEED 8 TO 1/500 SEC
055010 RANGE FINDER SPLIT IMAGE
055010 VIEWING THRU LENS
055010 FIELD OF VIEW 40 DEGREES
055010
055010 CAMERA, PLATE FILM LS938 11111
055010
055010 PROVIDE STILL PHOTOGRAPHIC COVERAGE
055010
055010 FILM TYPES.....HIGH RESOLUTION GLASS PLATE
055010 STANDARD 2 1/4 X 2 1/4 IN (57.15
055010 X 57.15 CM) ROLL
055010 70MM
055010 POLAROID
055010 SHUTTER SPEED.....1 SEC TO 1/500 SEC
055010
055010 CAMERA - STILL (VISIBLE) ST020 31111
055010
055010 STILL PHOTOGRAPHIC COVERAGE OF THE EARTH TARGET AREA IN THE
055010 VISIBLE REGION (REQUIRES SIMULTANEOUS PHOTOGRAPH IN THE INFRARED
055010 RANGE). USED IN ST015, ST025 AND ST035.
055010
055010 VIEW ANGLE 14.5 DEGREES
055010 RESOLUTION 100 LINES/MM
055010 FILM FORMAT 5 INCHES
055010 SHUTTER 8 TO 1/1000 SEC
055010
055010 SCOPE CAMERA CAC58 21111
055010
055010 OBTAIN FILM RECORD OF DATA AND DISPLAYS.
055010
055010 REDUCTION RATIO 1:1, 1:10, 7
055010 LENS F/1.8 TO F/16
055010 SHUTTER 8 TO 5 SECONDS
055010 FILM TYPE POLAROID
055010
055010 CAMERA - CINE APR29 11111
055010
055010 FILM RECORDER OP002 41121
055010
055010 PROVIDE PHOTOGRAPHIC RECORDINGS TO PROVIDE ALL WEATHER STEREO
055010 IMAGERY FOR TOPOGRAPHIC MAPPING, 3-D STRAIN FIELDS, EROSION, VOL-
055010 CANIC MOTION AND POST-GLACIAL UPLIFT.
055010
055010 FILM TYPE 70 MM
055010 FILM RATE 8 MM/SEC
055010 INSTANTANEOUS FOV 0.006 RAD
055010
055010 CAMERA CINE LS936 11111
055010
055010 PROVIDE VISUAL RECORDS SUCH AS PHOTOMICROGRAPHY, TIME LAPSE
055010 RECORDS, PHOTOMICROSCOPY AND GENERAL EXPERIMENT DOCUMENTATION
055010
055010 OPERATION MODE.....VARIABLE FRAME RATE OR SINGLE
055010 FRAME PULSE OPERATION UNDER
055010 REMOTE CONTROL
055010 FILM SIZE.....16 OR 35MM
055010 LENS SELECTION.....70 MM MOUNT LENS SYSTEM; STANDARD
055010 20CM LENS AS A MINIMUM
055010 EXPOSURE SYSTEM.....AUTOMATIC OR ZOOM
055010 FILM CAPACITY.....10, 100 OR 200FT (3.0, 30.4 OR
055010 60.96M) ROLLS
055010
055010 CAMERA - CINE ST054 61111
055010
055010 PHOTOGRAPHIC RECORD OF THE ACTIVITY OF AEROSOL PARTICLES IN
055010 AERIAL, SLOW MOTION AND TIME-LAPSE SEQUENCES
055010 USED IN ST015, ST015, ST025, ST045, ST055 AND ST065.
055010
055010 VIEW ANGLE 40 DEGREES
055010 RESOLUTION 200 LINES/MM
055010 SHUTTER 1/50 TO 1/500 SEC
055010 FILM FORMAT 16 MM
055010 CAMERA SPEED 1, 1/2, 10, 64 FPS
055010
055010 CAMERA, CINE CAC63 11111
055010
055010 PROVIDE PHOTOGRAPHIC RECORD OF TELESCOPE FIELD-OF-VIEW.
055010
055010 FILM FORMAT 35 MM
055010 OPERATION REMOTE
055010 FRAME RATE 16 FPS
055010 FILM TYPE CELCO
055010 RESOLUTION 50 LINES/MM
055010 LENS 100 MM
055010 FIELD OF VIEW 40 DEGREES
055010
055010 CENTREFUEE ST090 11111
055010

ORIGINAL PAGE IS
OF POOR QUALITY

067000 PROVIDE VARIABLE GRAVITY ENVIRONMENTS TO BIOSCIENCE SPECIMEN
067000 GENERALLY LIMITED TO TEST-TUBES SIZE CONTAINERS. USED IN ST01S,
067000 ST02S AND ST03S.
067000 GRAVITY LEVEL INDUCED UP TO 10 G REQUIRED
067000 CENTRIFUGE, MICRO LS942 11111
068000 CENTRIFUGAL SEPARATION OF BIOLOGICAL SAMPLES
068000 CAPACITY.....4 SAMPLES
068000 SPEED(RPM).....1750RPM
068000 GLASS TUBE SIZES.....0.9 TO 5 ML
068000 CENTRIFUGE, REFRIGERATED, HIGH SPEED LS941 11111
068500 SEPARATION EQUIPMENT TO SUPPORT MEDICAL AND BIOLOGICAL RESEARCH.
068500 OPERATING TEMPERATURE.....35 TO 50F (2 TO 10C)
068500 SPEED.....C TO 2000RPM
068500 GAS CHROMATOGRAPH SF003 21111
074700 THE FUNCTIONAL REQUIREMENTS OF THE EQUIPMENT ARE TO ANALYZE THE
074700 ATMOSPHERIC COMPOSITION AND THE IMPURITY CONCENTRATION CONTAINED
074700 WITHIN THE EXPERIMENT'S EQUIPMENT ENCLOSURES. THIS ELEMENT IS
074700 ALSO USED ON THE LEVITATION SUBELEMENT.
074700 FLOW RATE 10 ML/MINUTE
074700 SENSITIVITY.....0.1 PPM
074700 GAS CHROMATOGRAPH ST175 21111
074700 ANALYSIS OF COMBUSTION PRODUCTS FOR FLAME CHEMISTRY EXPERIMENTS.
074700 USED IN ST04S AND ST05S.
074700 GAS CHROMATOGRAPH LS962 11111
074700 MEASURES CONCENTRATION OF GAS, LIQUID, AND SOLID CONSTITUENTS OF
074700 BIOLOGICAL SAMPLES
074700 OPERATING TEMPERATURE.....-4C TO 752F(-4C TO 400C)
074700 CARRIER GAS.....HELIUM
074700 DETECTOR.....HYDROGEN FLAME
074700THERMAL CONDUCTIVITY
074700 FRACTION COLLECTION SYSTEM
074700 BACTERIAL COLONY COUNTER LS908 11111
074700 MANUAL COUNT OF BACTERIAL COLONIES. UNIT IS STANDARD TYPE
074700 AUTOMATIC COLONY COUNTER LS952 11111
074700 DETERMINE THE NUMBER OF VIABLE COLONIES GROWING ON AN ALGUA
074700 SUBSTRATE
074700 COMBUSTION CHAMBER ST170 11135
074700 CHAMBER WHERE FUELS AND OXIDIZERS ARE INJECTED, ALLOWED TO BURN
074700 AND THE FLAME CAPABLE OF BEING MONITORED.
074700 SIZE APPROXIMATELY 1.87 DIA X 3.5 FT
074700 (0.91 DIA X 1.0 M)
074700 COMBUSTION CHAMBER ST202 11135
074700 CHAMBER TO EXAMINE COMBUSTION IN ZERO GRAVITY
074700 DIMENSIONS..... 1.47FT(0.45M) DIA X 4.0FT(1.22M)
074700 FUNCTION KEYBOARD AS003 11111
101340 ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO
101340 DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY,
101340 FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340 ALPHANUMERIC KEYBOARD AS004 11111
101340 ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE COUNCIL COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS.
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340 FUNCTION KEYBOARD HE003 11111
101340 ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO
101340 DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY,
101340 FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340 ALPHANUMERIC KEYBOARD HE004 11111
101340 ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE COUNCIL COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS.
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340 FUNCTION KEYBOARD SC003 11111

101340
101340 ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO
101340 DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY,
101340 FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
101340 ALPHANUMERIC KEYBOARD SC004 11111
101340
101340 ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE CNRCARD COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
101340 KEYBOARD DISPLAY TERMINAL AF819 11111
101340
101340 GENERAL PURPOSE INPUT KEYBOARD ALLOWING FOR PROGRAMMING, DATA
101340 INPUT AND UPDATE, AND DATA MANIPULATION FOR CNRCARD ANALYSIS
101340
101340 ALPHANUMERIC KEYBOARD EC001 11111
101340
101340 ALLOWS CREWMAN TO COMMUNICATE WITH THE CNRCARD COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
101340 MULTIFUNCTION KEYBOARD MEC 11112
101340
101340 CEASE/F CONTROL KEYBOARD WITH FUNCTION KEYS LABELLED BY COMPUTER.
101340 A SUPPORTING SOFTWARE LOGIC TREE FOR DATA ENTRY AND DISPLAY CALL
101340 UP CAN BE ENTER AT ANY LEVEL.
101340
101340 CONTROL UNIT - OPERATOR SF001 11111
101340
101340 THE OPERATOR'S CONTROL PANEL WILL HAVE THE CAPABILITY OF
101340 CONTROLLING ABOUT 1000 PROCESS LOOPS. DATA MUST BE PRESENT ON A
101340 CLICK-LOCK BASIS. IT MUST PROVIDE A WAY TO ENTER NEW PARAMETERS
101340
101340 CHANNEL NUMBER 3 DIGIT NINE/12 EACH
101340 DATA LEVEL 4 DIGIT NINE PLUS SIGN/2 EACH
101340 LOGIC LEVELS, LOGICAL CV
101340 LOGICAL V
101340 DATA OUTPUT VARIABLE
101340
101340 KEYBOARD LS910 11111
101340
101340 TYPEWRITER TYPE KEYBOARD USED FOR PROVIDING INSTRUCTIONS FOR THE
101340 COMPUTER
101340
101340 KEYBOARD, COMPUTER CASCR 11111
101340
101340 PROVIDES FOR CREWMAN TO COMMUNICATE WITH CNRCARD COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340
101340 PHASE METER, DIGITAL CAP04 11111
101340
101340 MEASURES RELATIVE PHASE ANGLES OF RECEIVED SIGNALS FOR INTERFER-
101340 METER NAVIGATION/SURVEILLANCE EXPERIMENT.
101340
101340 OUTPUT GENERATES 13 BIT PHASE PEAK
101340 BASED ON TWO 60 MHZ SIGNALS
101340 AT 10 DPM.
101340 FREQUENCY RANGE 400 MHZ
101340 RESOLUTION 25 KHZ
101340 ACCURACY +/- 2.5 KHZ
101340
101340 DIGITAL PHASE METER ST011 41111
101340
101340 METER TO DISPLAY PHASE OF SIGNAL RECEIVED BY MICROWAVE
101340 INTERFEROMETER
101340
101340 DIGITAL COMPUTER LS947 11111
101340
101340 PROVIDE AUTOMATIC DATA ACQUISITION OF EXPERIMENTAL DATA, CONTROL
101340 LIGHTING AND ANIMAL FEEDING AND GENERAL LABORATORY ENVIRONMENT
101340 MONITORING. THE COMPUTER WILL SUPPORT ANALYTICAL TOOL SUCH AS
101340 MASS SPECTROMETER, GAS ANALYZERS AND ANALYSIS EQUIPMENT. ALSO
101340 DRIVES ONE OR MORE CRT INTERACTIVE DISPLAY CONSOLES
101340
101340 MEMORY CAPACITY 16 BIT 16,000 WORDS
101340 CYCLE TIME MICROSECOND
101340 PERIPHERAL CAPABILITIES DIRECT MEMORY ACCESS
101340 FLEXIBLE INTERRUPT STRUCTURE
101340 SINGLE BIT SET AND SENSE LINE
101340 CAPABILITY OF UP TO 64 LINES
101340 EACH
101340
101340 COMPUTER APB1A 21111
101340
101340 COMPUTER CAPABILITY TO SUPPORT ATMOSPHERIC AND SPACE PLASMA
101340 PHYSICS INVESTIGATIONS.
101340
101340 WORD SIZE 16 BIT
101340 I/O MEMORY TRANSFER 1600-16 BIT WORDS/SEC
101340 CORE MEMORY 32K*16 WORDS
101340 I/O CHANNELS 64
101340
101340 DIGITAL COMPUTER EC010 11111
101340
101340 PROVIDES DATA PROCESSING AND CONTROL FUNCTIONS FOR THE CFI
101340

ORIGINAL PAGE IS
OF POOR QUALITY

102700 DIGITAL PROCESS PROGRAMMER SFOC1 11111
102700 PERFORM PREDETERMINED LOGIC ROUTINES SUCH AS EVENT SEQUENCING.
102700 PARAMETER COMMANDS ON RELAY TIME SEQUENCING.
102700
102700 COMPUTER, GENERAL PURPOSE CA901 21111
102700
102700 COMPUTER CAPABILITY TO COMMAND EXPERIMENTS ON AND OFF, PROGRAM
102700 EXECUTION OF EXPERIMENT PROCEDURAL STEPS, COMPLETE GIMBAL
102700 ANGLES FOR TELESCOPE AND ANTENNA POINTING, AND COMPUTATION OF
102700 POINT-TO-POINT ANGLES.
102700
102700 COMPUTER ST024 21111
102700
102700 COMPUTER CAPABILITY TO PROVIDE BASIC POINTING AND TRACKING INFOR-
102700 MATION RELATIVE TO ATL STABILITY AND FLIGHT PATH REQUIREMENTS.
102700 USED ON ST01S, ST03S, ST04S AND ST05S.
102700
102750 OPERATIONS CONTROL COMPUTER MEC 11111
102750
102750 THE COMPUTER SHALL PROVIDE THE FOLLOWING FUNCTIONS: COMMAND
102750 INTERPRETATION AND EXECUTION, ATTITUDE AND POINTING PROGRAM
102750 CONTROL AND COMPUTATION AND DATA MANAGEMENT SEQUENCING OPERATIONS
102750
102750 MEMORY SIZE.....1455C WORDS
102750 ADD/MULT EXECUTE TIME.....5/3C PICOC SECONDS
102750
103500 MULTIPLEXER A/C CONVERTER SPOC1 11111
103500
103500 CONVERT INPUT SIGNAL FROM ANALOG TO DIGITAL FORMAT FOR I/C UNIT
103500
103500 ANALOG INPUT.....-1C.24 TO +1C.235V
103500 FULL SCALE ACCURACY.....+0.025
103500 RESOLUTION.....12 BITS BINARY
103500 CONVERSION SPEED.....8 KHZ
103500 DATA OUTPUT.....2'S COMPLEMENT
103500
103500 A/C CONVERTER CA062 21111
103500
103500 CONVERT INPUT SIGNAL FROM ANALOG TO DIGITAL FORMAT FOR I/C UNIT
103500
103500 C/A AND A/C CONVERTER CA081 11111
103500
103500 PROVIDES FOR A/C AND C/A CONVERSION.
103500
103500 A/C CONVERTER
103500 SAMPLING RATES15 MEGASAMPLES/SEC
103500 DATA RATES105 MBS
103500 C/A CONVERTER
103500 BANDWIDTH8 MHZ
103500
103500 DATA CONVERTER - ANALOG TO DIGITAL ST067 21111
103500
103500 CONVERT ANALOG SIGNAL DATA TO DIGITAL DATA FOR RECORDING AND
103500 LATER DUMPING TO GROUND STATIONS. USED ON ST01S AND ST03S.
103500
103500 ANALOG TO DIGITAL CONVERTER LS948 11111
103500
103500 SUPPORT DIGITAL COMPUTER. INCLUDES LC-LEVEL MULTIPLEXER. UNIT
103500 SHOULD BE A PROGRAMMABLE GAIN 12 BIT TYPE
103500
103520 INPUT/OUTPUT STAGE SPO01 11111
103520
103520 THE INPUT/OUTPUT UNIT, THROUGH THE USE OF PLUG IN CIRCUITS, CAN
103520 ACCOMMODATE A WIDE RANGE OF PROCESS INPUTS AND OUTPUTS FOR MAXIMUM
103520 CONTROL VERSATILITY.
103520
103520 LOGIC LEVELS, LOGICAL "0".....0V
103520 LOGICAL "1".....5V
103520 SIGNAL INPUTS.....COMPATIBLE WITH POP0 TYPE
103520 CCMPLTERS
103520 SIGNAL OUTPUTS.....0 - 5V
103520 DATA OUTPUTS.....12 BIT PARALLEL
103520
103520 DATA BIT STREAM GENERATOR CA083 11111
103520
103520 PROVIDES TEST BIT STREAMS FOR TRANSMISSION.
103520
103520 DATA RATES1 TO 105 MBS
103520
103520 DATA PROCESSING EQUIPMENT ST035 11111
103520
103520 HIGH SPEED DIGITAL DATA CONVERSION TO ANALOG FORM COMPATIBLE FOR
103520 PAPER TIME DISPLAY. USED ON ST01S, ST02S AND ST03S.
103520
103520 CHANNELS32
103520 DIGITAL-DATA10C+06 BPS
103520
104000 SUPPORT ELECTRONICS ST006 21114
104000
104000 ELECTRONICS SUPPORT COMMUNICATIONS EXPERIMENTS INCLUDING
104000 MULTIPLEXERS WITH CONVERTERS AND OSCILLATORS
104000
104020 SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) AS002 11111
104020
104020 PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CRT
104020
104020 CHANNELS2
104020 FORMATVIDEO AND DATA
104020 VIDEO INPUT PRESENTATIONRASTER SCAN
104020 SYMBOL WRITING TECHNIQUESTROKE
104020 INTERFACE DESCRIPTIONDIGITAL-12-BIT DATABUS
104020
104020 SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) ME002 11111
104020
104020 PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CRT

ORIGINAL PAGE IS
OF POOR QUALITY



```

104020 CHANNELS ..... 2
104020 FORMAT ..... VIDEO AND DATA
104020 VIDEO INPUT PRESENTATION ..... RASTER SCAN
104020 SYMBOL WRITING TECHNIQUE ..... STROKE
104020 INTERFACE DESCRIPTION ..... DIGITAL-12-BIT DATAWORD
104020 SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) SCOC2 11111
104020 PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CPT
104020 CHANNELS ..... 1
104020 FORMAT ..... VIDEO AND DATA
104020 VIDEO INPUT PRESENTATION ..... RASTER SCAN
104020 SYMBOL WRITING TECHNIQUE ..... STROKE
104020 INTERFACE DESCRIPTION ..... DIGITAL-12-BIT DATAWORD
104020 AUTOMATIC DISPLAY GENERATOR AFPC9 11111
104020 PROVIDES 30 DISPLAY
104160 PRINTER SPOO1 11111
104160 PRINT RECORDS DATA IN THE FORM OF ALPHA-NUMERIC CHARACTERS
104160 IN HARD COPY FORM UPON REQUEST BY OPERATOR.
104160 FORMAT ..... PICA TYPE - STANDARD
104160 CHARACTERS, ALPHABETIC SYMBOLS
104160 DIGITAL PLOTTER, PRINTER 159C9 11111
104160 TELEPRINTER CA903 11111
104160 PROVIDES CAPABILITY TO SEND OR RECEIVE TELETYPE MESSAGES.
104160 CHARACTER RATE ..... 30 CHARACTERS/SEC (150 CPS)
104160 CONSISTS OF ..... SOLID STATE KEYBOARD
104160 ..... SOLID STATE THERMAL PRINTER
104160 ..... EXPANDABLE MESSAGE MEMORY
104180 DISPLAY CRT (CONTROL/DISPLAY CONSOLE) ASCO1 21111
104180 PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND
104180 STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS. INCLUDES
104180 DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.
104180 SCREEN SIZE ..... 14 IN (0.355 M)
104180 DISPLAY CRT (CONTROL/DISPLAY CONSOLE) HEOC1 21111
104180 PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND
104180 STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS. INCLUDES
104180 DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.
104180 SCREEN SIZE ..... 14 IN (0.355 M)
104180 DISPLAY CRT (CONTROL/DISPLAY CONSOLE) SECC1 21111
104180 PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND
104180 STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS. INCLUDES
104180 DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.
104180 SCREEN SIZE ..... 14 IN (0.355 M)
104180 STATUS PANEL APB20 21111
104180 DISPLAY TERMINAL WITHOUT KEYBOARD PROVIDING STATUS OF EXPERIMENT
104180 FUNCTIONS
104180 SPECIAL DATA ACQUISITION PANEL APP21 11124
104180 CRT DISPLAY EOC01 11111
104180 PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND
104180 DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.
104180 DIGITAL VISUAL DISPLAY EOC01 11111
104180 DISPLAY NUMERIC DATA
104180 CRT DISPLAY PEC 61111
104180 TWO CRT'S ARE LOCATED IN THREE CONSOLES FOR DATA DISPLAY OF
104180 ALPHANUMERICS AND GRAPHICS.
104180 RESOLUTION.....1000 LINES
104180 CONTROL INPUT FORM.....DIGITAL
104180 SCAN.....RASTER
104180 .....CCNICAL
104180 PHOSPHOR CHARACTERISTICS.....FCUR-COLOR PENETRATION
104180 SCREEN SIZE.....14IN (40CM)
104180 CONTRAST.....6 GRAY LEVELS
104180 TELEPRINTER SF001 11111
104180 CHARACTER ORIENTED CPT TERMINAL. DATA IN THE LINE MEMORY IS
104180 CONVERTED, BY THE CHARACTER GENERATOR, INTO THE APPROPRIATE CPT
104180 PATTERN. THE DOTS ARE SHIFTED OUT OF THE SHIFT REGISTER TO FORM
104180 A VIDEO SIGNAL TO THE MONITOR.
104180 SCREEN FORMATS.....40,72 OR 96 CHARACTERS PER LINE
104180 .....12 OR 24 LINES
104180 TRANSFER RATE.....11C TO 24C0 PALE, 10 CP 11 BIT
104180 CHARACTERS
104180 MODES.....HALF OR FULL DUPLEX
104180

```

ORIGINAL PAGE IS
OF POOR QUALITY

C-3

104100 CDISPLAY, CRT CA912 11111
104100 CDISPLAY REAL-TIME OR VIDEO TAPE DATA.
104100 RESOLUTION 625 LINES
104100 FRAME RATE 15 FPS
104100 CDISPLAY TERMINAL ST020 21111
104100 TERMINAL CONSISTING OF ALPHANUMERIC KEYBOARD AND CRT DISPLAY FOR
104100 MONITORING EXPERIMENT ACTIVITIES
104100 CRT DISPLAY LS951 41111
104100 DISPLAY GRAPHIC AND ALPHANUMERIC DATA
104100 ELECTROPHYSIOLOGY DISPLAY LS954 11114
104100 INDICATE DISPLAY FOR PHYSIOLOGICAL DATA TRANSMITTED BY THE
104100 ELECTROPHYSIOLOGY BACKPACK
104200 SCANNER PROGRAMMER SF001 11111
104200 TRACKS COMPLEX FUNCTION AND TRANSLATES IT INTO A SIGNAL WHICH
104200 DRIVES EXPERIMENT CONTROLLER
104200 ACCURACY.....0-0.50 OF SPAN
104200 REPEATABILITY.....0-0.1 OF SPAN
104300 SCAN PROGRAM GENERATOR CAC09 31111
104300 GENERATES DIGITAL CONTROL SIGNALS TO CONTROL RECEIVER SWEEP.
104300 SUPPLIES DIGITAL READOUT OF FREQUENCY AND LEVEL.
104300 BANC 24
104300 SCAN TIME 0.03 TO 300 SECONDS
104300 CAPABILITY ANALOG DATA COLLECT & DISPLAY
104400 CISK STORAGE SYSTEM LS949 11111
105004 STORES DATA IN SUPPORT OF COMPUTER
105004 STORAGE CAPACITY 5,000 WORDS
105500 MULTIPLEXER MEC 11111
105500 MULTIPLEX DATA FOR TRANSMISSION VIA SHUTTLE COMMUNICATION LINK
106000 CRYOGENIC COOLER ASSEMBLY EC001 11114
106000 PROVIDE LOW TEMPERATURE HEAT SINK TO EXPERIMENT CHAMBERS
106000 MAXIMUM TEMPERATURE.....-76F (-60C)
106000 HEAT SINK CAPACITY(EST).....1000BTU (2900 THERMAL)
106000 CEMAR SF001 11111
106000 THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES
106000 TEMPERATURE.....-112F (-80C)
106000 LENGTH.....1FT (30CM)
106000 WIDTH.....1.8FT (54CM)
106000 HEIGHT.....2.0FT (60CM)
106000 CRYOSTAT (CMR) LS022 11111
106000 CAPABILITY OF STORING QUANTITIES OF CRYOGENIC LIQUID
106000 QUANTITY 5 GALLONS (19.0 LITERS)
106000 TEMPERATURE -319 F (-195 C)
106000 FREEZER, CRYOGENIC LS916 11111
106000 PROVIDE A MEANS FOR FREEZING AND STORING BIOLOGICAL SPECIMENS
106000 STORAGE CAPACITY.....11.0 CU. IN. (1.0 LITERS)
106000 OPERATING TEMPERATURE.....-32CF (170K)
106000 MIN CRYOGEN FILLING TIME.....2 WEEKS BETWEEN REFILLS
106000 CRYOSTAT - HELIUM ST192 11111
106000 CEMAR TYPE STORAGE CONTAINER TO STORE AND DISPENSE LIQUID HELIUM.
106000 SIZE APPROXIMATELY 10.5 DIA X 10.7 FT
106000 (10.5 DIA X 0.05 FT)
110300 BIT ERROR COUNTER CA074 11111
110300 COUNT ERROR PULSES OUTPUT BY THE BIT ERROR DETECTOR TO DETER-
110300 MINE THE BIT ERROR RATE OF THE DIGITAL TRANSMISSION SYSTEM.
110300 CATALYSIS UNIT SF001 11113
120000 DYNAMOMETER LS103 11133
120000 MEASUREMENT OF FORCE OR POWER.
120000 ARM FLEXION 25 TO 75 LB (10 TO 35 KG)
120000 PACK EXTENSION 100 TO 200 LB (40 TO 80 KG)
120000 REC/VCG LS101 11131
120000 MONITOR HEART CONDITION. ELECTROVETERCARDIOGRAPH.
120000 SENSITIVITY 0 TO 3 PV +/-1 PCT FS
120000 HEART RATE 40 TO 180 BEATS/PIN



143500 PHONOCARDIOGRAM LS106 11111
143500 ALPAL MONITOR OF HEART BEAT.
143500 SENSITIVITY 0.1 TO 1000 HZ, +/-5 PCT FS
144500 EEG - ELECTROENCEPHALOGRAPH LS102 11131
144500 MONITOR ELECTRICAL IMPULSES GENERATED BY THE BRAIN.
144500 SENSITIVITY 1C TO 20C MICRONS +/-1 PCT FS
144500 FREQUENCY RESPONSE 0.2 TO 100 HZ
144500 CHANNELS 6
145500 EPG - ELECTROMYOGRAPH LS104 11131
145500 MONITOR THE MOVEMENT OF THE BODY MUSCULAR SYSTEM.
145500 SENSITIVITY 0.01 TO 5 MV +/-1 PCT FS
145500 BANDWIDTHS 0.5 TO 200 HZ, 0.5 TO 1000 HZ
145500 CHANNELS 1
146915 ANTHROPEMETRIC CRIC LS923 11111
146915 ANTHROPEMETRIC MEASUREMENTS OF VERTEBRATES EITHER REMOTELY OR
146915 AUTOMATICALLY.
146915 SMALL GRID 5 X 15 CM SPACED AT 1 MM
146915 MEDIUM GRID 5C X 75 CM SPACED AT 2 MM
146915 LARGE GRID 2 X 2 M SPACED AT 5 MM
149300 ELECTROPHORESIS APPARATUS ST08R 11111
149300 MEASURE ELECTROPHORESIS MOBILITY, SURFACE-ZETA POTENTIAL, AND
149300 SURFACE CHARGE DENSITY OF CELL LINES OVER THEIR LIFE CYCLE. USEC
149300 IN ST01S AND ST02S.
149300 CONTINUOUS FLOW ELECTROPHRETIC COLUMN WITH PLMP SP001 11114
149370 THIS UNIT SEPARATES BIOLOGICAL SAMPLES BY ELECTROPHORESIS.
149370 RELATIVELY LARGE SAMPLES, APPROXIMATELY 10 CC, WILL BE PROCESSED.
149370 SIMPLE AND BUFFER FLOW WILL BE CONSTANT AND SAMPLE FRACTION
149370 COLLECTION CONTINUOUS.
149370 TOTAL BUFFER AND SAMPLE VOLUME... 400L IN. (1000CC)
149370 OPERATING TEMPERATURE..... 14 TO 41 F (-10 TO 5C)
149370 OPERATING TEMPERATURE REGULATION +/-2F (+/-1C)
149370 SEPARATION PROFILE ACCURACY... +/-0.004IN (+/-1PP)
149370 POSITION MONITOR ACCURACY... +/-0.004IN (+/-1PP)
149370 FRACTION OUTLET PORT LOCATION ACCURACY +/-0.004IN (+/-1PP)
149370 INPUT VOLTAGE..... 500V
149390 STATIONARY ELECTROPHRETIC COLUMNS SF001 11114
149390 THE STATIONARY ELECTROPHRETIC COLUMNS HAVE TWO ANTICIPATED
149390 USAGES: (A) TO ESTABLISH PROCESSING AND OPERATING PARAMETERS FOR
149390 THE CONTINUOUS ELECTRORETIC COLUMNS, AND (B) FOR SMALL BATCH
149390 PROCESSING (APPROXIMATELY 0.1CC (0.001 LITERS)) OF VARIOUS
149390 BIOLOGICAL COMPONENTS
149390 TOTAL BUFFER AND SAMPLE VOLUME1.5CL IN. (25CC)
149390 SEPARATION VELOCITY ACCURACY... +/-0.04IN (+/-1PP)
149390 POSITION MONITOR ACCURACY... +/-0.04IN (+/-1PP)
149390 OPERATING TEMPERATURE RANGE... 14 TO 41F (-10 TO 5C)
149390 TEMPERATURE TOLERANCE..... +/-2F (+/-1C)
152500 ELECTROPHORESIS APPARATUS LS956 11111
152500 SEPARATES PROTEIN AND AMINO ACID CONSTITUENTS IN SERUM, PLASMA,
152500 URINE OR SPINAL FLUID FOR QUANTITATIVE ANALYSIS.
152500 VOLTAGE..... 0 TO 500 V CC
152500 SUBSTRATE..... STARCH BLOCK, PAPER OR GEL
152500 SAMPLE VOLUME..... 0.001 TO ML
159100 CONDITIONING CHAMBER EC019 11134
159100 THE CONDITIONING CHAMBER IS MADE UP OF A PRESSURE SHELL, A
159100 POSITIVE EXPULSION REGULATOR AND A HYDROSCAL CONDITIONER. THIS
159100 CHAMBER CONDITIONS THE GAS SAMPLE PRIOR TO ITS INJECTION INTO
159100 THE CLOUD CHAMBER
159100 OPERATING PRESSURE(EST)..... 20PSIA (14603N/M2)
159100 OPERATING TEMPERATURE..... 32 TO 86F (0 TO 30C)
159100 VOLUME..... 2 CL. FT. (0.0565 CL. METERS)
159310 EXERCISER/ERGOMETER - PRIMATE LS175 11133
159310 PRIMATE EXERCISER AND WORKLOAD MEASURING DEVICE.
159400 ERGOMETER - BICYCLE LS986 11131
159400 WORKLOAD MEASURING DEVICE.
159400 LOAD RANGE 0 TO 300 WATTS, +/-2 PCT FS
159400 ROTATIONAL SPEED 4C TO 50 RPM
164000 MICROFILM VIEWER AS005 11111
164000 PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUB-
164000 SYSTEM OPERATIONAL PROCEDURES, CHECKOUT CHECKLIST PROCEDURES,
164000 SIMPLEX SCHEMATICS AND OTHER WRITTEN OR PICTORIAL INFORMATION.
164000 FILM FORMAT 16 MM 35MM TRACK
164000 FILM LOADING CASSETTE

ORIGINAL PAGE IS
OF POOR QUALITY

```

164000 FILM SLEWING ..... MANUAL AND COMPUTER SELECT
164000 INTERFACE DESCRIPTION ..... DIGITAL 13-BIT BINARY
164000
164000 MICROFILM VIEWER ..... HE005 11111
164000
164000 PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUB-
164000 SYSTEM OPERATIONAL PROCEDURES, ENDORSE CHECKOUT PROCEDURES,
164000 SIMPLEX SCHEMATICS AND OTHER WRITTEN OR PICTORIAL INFORMATION.
164000
164000 FILM FORMAT ..... 16 MM DUAL TRACK
164000 FILM LOADING ..... CASSETTE
164000 FILM SLEWING ..... MANUAL AND COMPUTER SELECT
164000 INTERFACE DESCRIPTION ..... DIGITAL 13-BIT BINARY
164000
164000 MICROFILM VIEWER ..... SC004 11111
164000
164000 PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUB-
164000 SYSTEM OPERATIONAL PROCEDURES, ENDORSE CHECKOUT PROCEDURES,
164000 SIMPLEX SCHEMATICS AND OTHER WRITTEN OR PICTORIAL INFORMATION.
164000
164000 FILM FORMAT ..... 16 MM DUAL TRACK
164000 FILM LOADING ..... CASSETTE
164000 FILM SLEWING ..... MANUAL AND COMPUTER SELECT
164000 INTERFACE DESCRIPTION ..... DIGITAL 13-BIT BINARY
164000
164000 MICROFILM READER ..... PFC 11111
164000
164000 READER CONCURRENTLY DISPLAYS SEQUENCES OF MICROFILM RECORDS OF
164000 PAPS, REFERENCE SIGNATURES AND OTHER MATERIALS NEEDED FOR
164000 REAL TIME EVALUATION OF SENSOR FORMATS SELECTED FOR CRT DISPLAY
164000
164000 STORAGE CAPACITY.....13,000 PAGES
164000 RETRIEVAL SPEED.....4 SECONDS MAX
164000 INTERFACE.....STANDARD COMPUTER
164000
171600 TFSKROARE - FORCE/TORQUE ..... LS925 11131
171600
171600 MEASURE MAN'S CAPABILITY TO APPLY FORCES AND TORQUES IN A VARIETY
171600 OF DIRECTIONS FROM VARIOUS HAND/BODY ORIENTATIONS AND RESTRAINT
171600 CONDITIONS.
171600
171600 HAND FORCE APPLICATIONS ..... +/- 6.5 POUNDS
171600
171600 FLOW METER, GAS ..... LS992 11111
171600
171600 MEASURE AIR FLOW
171600 FLOW METER ..... SF001 11111
171600
171600 GAS FLOW METERING AND MONITORING DEVICE.
171600
171600 GAS TYPE ..... OXYGEN, NITROGEN, HELIUM
171600 GAS FLOW RATE ..... 0.4 TO 0.3 LB/HR
171600
171600 FLOW METER, WATER MANIFOLD ..... LS989 11235
171600
171600 IN LINE MEASUREMENT OF WATER FLOW, GENERALLY A LOW RATE
171600 ASSOCIATED WITH WATER CONSUMPTION BY ORGANISMS
171600
171600 FRACTION COLLECTION SYSTEM ..... SF001 11235
171600
171600 THIS APPARATUS WILL HAVE TWO FUNCTIONS: (1) TO COLLECT THE
171600 DESIRED FRACTIONS AFTER SEPARATION, AND (2) TO REMOVE THE EXCESS
171600 BUFFER SOLUTION AND REMAINING UNWANTED SAMPLES
171600
171600 PUMP TYPE.....POSITIVE DISPLACEMENT
171600 PUMPING RATE ACCURACY.....+/-0.1%
171600 ALLOWABLE PUMP PULSATION.....0.1% (10-70 C/MIN)
171600 TEMPERATURE REGULATION.....+/-1.0 C (10-70 C/MIN)
171600
171600 FREEZER ..... SF001 11111
171600
171600 THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES
171600
171600 CONTROLLER TEMPERATURE RANGE.....13F TO 113F (-25C TO 45C)
171600 LENGTH.....1 FT (30CM)
171600 WIDTH.....1.5 FT (45CM)
171600 HEIGHT.....2.5 FT (76CM)
171600
171600 FREEZER, GENERAL ..... LS957 11111
171600
171600 STORAGE OF SERUM, PLASMA, SPECIMENS AND ORGANISMS
171600
171600 OPERATING TEMPERATURE.....-4F (-20C)
171600 TEMPERATURE TOLERANCE.....+/-4F (+/-2C)
171600 STORAGE VOLUME.....4 CU. FT. 10.11 CU. M.
171600
171600 FREEZER, LOW TEMP ..... LS958 11111
171600
171600 STORAGE OF EXPERIMENT SPECIMENS
171600
171600 OPERATING TEMPERATURE.....-64F (-70C)
171600 STORAGE VOLUME.....1 CU. FT. 10.02 CU. M.
171600
171600 SPECTRUM ANALYZER ..... APR03 41111
171600
171600 MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE RATIO AND
171600 MULTIPLE RATIO OF VARIOUS FREQUENCY SOURCES.
171600
171600 BANDWIDTH ..... 20 HZ TO 300 KHZ
171600 SENSITIVITY ..... 20 MICROVOLTS
171600 ACCURACY ..... +/- 0.5 DB
171600
171600 PULSE WAVE ANALYZER- C.W. MODES ..... APR13 11111

```




179500 PLTI-CHANNEL ANALYZER AFRC4 11111
179500 MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE RATIO AND
179500 MULTIPLE RATIO OF MULTIPLE FREQUENCY SOURCES SIMULTANEOUSLY.
179500 BANDWIDTH 0 TO 10 MHZ
179500 TIME BASE 200 MHZ CLOCK
179500 MEMORY SIZE 1048, 4096, 8192 CHANNELS
179500 WAVE ANALYZER AFRC5 11111
179500 DETECTION OF SIGNAL AMPLITUDE AND FREQUENCY INFORMATION
179500 BANDWIDTH 10 TO 600 KHZ
179500 FREQUENCY RESOLUTION 10 HZ
179500 SENSITIVITY 3 MV TO 30 V RMS
179500 CHANNELS 2 MINIMUM
179500 WAVE ANALYZER APPC6 11111
179500 DETERMINE, SEPARATE AND ANALYZE VARIOUS FREQUENCY COMPONENTS OF
179500 INPUT SIGNALS (IE FUNDAMENTAL, HARMONICS, INTERMODULATION PRO-
179500 DUCTS, ETC.).
179500 BANDWIDTH 1 KHZ TO 1.5 MHZ
179500 FREQUENCY RESOLUTION +/- 1.0 %
179500 SENSITIVITY 10 MV TO 100 V
179500 CHANNELS 2 MINIMUM
179500 FREQUENCY COUNTER AFRC6 11111
179500 DIRECT MEASUREMENT OF FREQUENCY AND/OR SIGNAL REPETITION RATE.
179500 BANDWIDTH 0 TO 35 MHZ
179500 SENSITIVITY 10 MV RMS
179500 TIME BASE 100MHZ TO 10 MHZ
179500 ANALYZER - MULTICHANNEL EFC11 11111
179500 ALLOW SELECTIVE ANALYSIS OF COUNTER DATA
179500 FREQUENCY COUNTER CPO59 11111
179500 MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE, TIME INTERVAL
179500 AND RATIO OF FREQUENCY SOURCES.
179500 BANDWIDTH 3 HZ TO 10 MHZ
179500 SENSITIVITY 0.01 RMS
179500 TIME BASE INTERNATIONAL SECCNE
179500 SPECTRUM ANALYZER/OSCILLOSCOPE CAS00 21111
179500 SINGLE DESIGN FUNCTIONS AS POWER SPECTRAL DENSITY AND WAVEFORM
179500 DISPLAY FOR QUICKLOOK AND AS MODULATION MONITOR.
179500 BANDWIDTH 1E GHZ
179500 SENSITIVITY 2 MV/DIVISION
179500 CHANNELS 2
179500 PHOTOGRAPHIC CAPABILITY YES
179500 SPECTRUM ANALYZER STC17 11111
179500 SUPPORT SETUP, CHECKOUT AND SIGNAL OUTPUT ANALYSIS OF RF RACIO-
179500 METER UNIT AND SWEEP RECEIVER (CRT DISPLAY OF SWEEP RECEIVER AMP-
179500 LITUDE VERSUS FREQUENCY).
179500 USED IN ST01S AND ST02S.
179500 BANDWIDTH 30 TO 10000 MHZ
179500 RESOLUTION 2 %
179500 SENSITIVITY 10 MV/DIVISION
179500 FREQUENCY METER ST156 11111
179500 MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD, AVERAGE RATIO AND
179500 MULTIPLE RATIO OF VHF FREQUENCY SOURCES. USED ON ST03S.
179500 BANDWIDTH 10 TO 100 MHZ
179500 SENSITIVITY +/- 2 %
181000 FREQUENCY SYNTHESIZER - WAVE ANALYSIS AF516 21111
181000 TRANSLATION OF A STABLE FREQUENCY OF A PRECISE STARGATE TO ANY OF
181000 A SELECTED INVESTIGATION REQUIREMENT.
181000 FREQUENCY BANDWIDTH DC TO 13 MHZ IN 5 RANGES
181000 FREQUENCY RESOLUTION 0.01 TO 10 MHZ
181000 FREQUENCY STABILITY 1E-07 PARTS PER DAY
181000 SIGNAL GENERATOR LSC52 11111
181000 ACCURATELY REPRODUCE SELECTED FREQUENCIES AS REFERENCES FOR PUR-
181000 POSES OF MAINTENANCE AND REPAIR.
181000 BANDWIDTH DC TO 20 KHZ
181000 OUTPUT VOLTAGE 3 V TO 0.1 MV
181000 FREQUENCY SYNTHESIZER AND DRIVER CAC56 11111
181000 RECEIVER FIRST LOCAL OSCILLATOR.
181000 FREQUENCY DC TO 500 MHZ
181000 STABILITY +/- 5E-11 PER DEG C
181000 OUTPUT VOLTAGE 1 V RMS +/- 1.5 DB
181000 CALIBRATION SIGNAL GENERATOR CAC91 11111

ORIGINAL PAGE IS
OF POOR QUALITY

101000 CALIBRATES RECEIVER SENSITIVITY AND PHASE METER.
101000
101000 FREQUENCY RANGE 1564.75 TC 1570.25 MHz
101000 FREQUENCY ACCURACY +/- 5 MHz
101000 FREQUENCY STABILITY C.CCS 7/CEGPEE C
101000 RF OUTPUT POWER +1C DB TC -127 DB
101000
101000 SIGNAL GENERATOR STC16 11111
101000
101000 ACCURATELY REPRODUCE IN PULSE FORM THE OUTPUT OF A MICROWAVE
101000 RADIOMETER MONITORING SEA SURFACE STATE AND TEMPERATURE.
101000 USED IN ST015.
101000
101000 REPETITION RATE Cb
101000 OUTPUT VOLTAGE +/- 5 VDC
101000
101000 RF INDUCTION COILS SPC04 11111
101000
101000 PROVIDE NON-CONTACT HEATING
101000
101000 FREQUENCY RANGE.....2KHZ TC 2MHZ
101000 OPERATIONAL ATMOSPHERE.....ANY
101000
101000 MICROWAVE HEATER SPCC2 31235
101000
101000 PROVIDE MICROWAVE ENERGY SOURCE FOR SPECIMEN HEATING.
101000
101000 TEMPERATURE(MAX).....18CCCF (1C.CCCC)
101000 POWER(MAX).....2KW
101000 FREQUENCY RANGE(MAX).....LP TC 181C HZ
101000 OPERATIONAL ATMOSPHERE.....ANY PUT VACUUM
101000
101000 RESISTANCE HEATER (CONTACT) SPOC2 21111
101000
101000 HEATING ELEMENT WHICH HEATS ENCAPSULATED SAMPLES BY PASSING A
101000 CURRENT THROUGH IT
101000
101000 POWER.....4KW SUSTAINED, 9KW PEAK
101000 TEMPERATURE(MAX).....29CCF (16CCC)
101000 ATMOSPHERE.....INERT GAS OR VACUUM
101000
101000 RESISTANCE HEATER (NON-CONTACT) SPC04 11111
101000
101000 TIRE TYPE HEATING UNIT FOR USE IN CHEST-GENERAL PURPOSE ENCLOSURE
101000 OR HOT WALL FURNACE
101000
101000 POWER.....4KW SUSTAINED, 9KW PEAK
101000 TEMPERATURE(MAX).....29CCF (16CCC)
101000 SPECIMEN SIZE.....4IN (1CCPILCNG X C.75 IN(2CP)
101000 DIA.
101000
101000 DIRECTIONAL SOLIDIFICATION UNIT SPOC2 11114
101000
101000 THIS UNIT CONSISTS OF A FURNACE AND A CHILL UNIT. SPECIMENS ARE
101000 HEATED, MOVED THROUGH A THERMAL GRADIENT AREA AND THEN INTO A
101000 COOLING AREA.
101000
101000 CHILL UNIT DIMENSIONS.....C.6 IN(1.9CHILL.D.X.0 IN(20CP)
101000 LCAG
101000 CHILL CAPACITY.....1KW THERMAL
101000 HEATING UNIT DIMENSIONS.....C.6 IN(2CP)X1.0X0IN(20CP) LCAG
101000 HEATER CAPACITY.....1KW THERMAL
101000 TEMPERATURE GRADIENT ADJUST
101000 RANGE.....4 TC 9CCF/IN (2 TO 20CC/CP)
101000 HEATER MAX TEMPERATURE.....29CCF(16CCC)
101000 CHILL UNIT MIN TEMP.....48 TC 77F (2C TC 25C)
101000 EXTERNAL SURFACE TEMP.....4C TC 11CF (15 TC 50C)
101000
101000 HOT WALL TUBE FURNACE SPOC2 31111
101000
101000 THIS UNIT IS A GENERAL PURPOSE HOT WALL HEATING DEVICE
101000 PROVIDING ACCURATE CONTROL OVER THE HOT ZONE TEMPERATURE WITH
101000 RESPECT TO A FLAT PROFILE OR TO A SPECIFIED GRADIENT PROFILE.
101000
101000 VACUUM LEVEL.....2E-8 PSI (1.9E-4A/M2)
101000 OPERATING TEMPERATURE.....22CCF (12CCF)
101000 MAX EXTERNAL SURFACE TEMP.....11CF (45C)
101000 HOT ZONE DIMENSIONS.....2.5 IN (6CP) DIA. X 6IN(15CP)LC
101000 MAXIMUM DEVIATION FROM FLAT
101000 PROFILE.....+/-4F (+/-2C)
101000
101000 HOT WALL FURNACE SPOC2 21111
101000
101000 THIS UNIT IS A GENERAL PURPOSE HOT WALL HEATING DEVICE PROVIDING
101000 ACCURATE CONTROL OVER THE HOT ZONE TEMPERATURE WITH RESPECT TO
101000 A FLAT PROFILE OR TO A SPECIFIED GRADIENT PROFILE.
101000
101000 VACUUM LEVEL.....2E-8 PSI (1.9E-4A/M2)
101000 OPERATING TEMPERATURE.....22CCF (12CCF)
101000 EXTERNAL SURFACE TEMPERATURE.....11CF (45C)
101000 HOT ZONE DIMENSIONS.....1IN (2.5CHILL.D.X.1IN(13CP)LCNG
101000 MAXIMUM DEVIATION FROM FLAT
101000 PROFILE.....+/-7F (+/-4C)
101000 VIEWPOINTS.....2
101000
101000 CHEST-GENERAL PURPOSE ENCLOSURE SPCC2 31113
101000
101000 THE PRIMARY FUNCTION OF THE ENCLOSURE IS TO ISOLATE THE HEATED
101000 SPECIMEN FROM THE LABORATORY ENVIRONMENT, BOTH THERMALLY AND
101000 ATMOSPHERICALLY, AND TO PROVIDE A MEANS TO MECHANICALLY ATTACH
101000 THE VARIOUS HEATING ELEMENTS AND ACCESSORY ITEMS
101000
101000 VACUUM LEVEL.....2E-1CPSI (1.9E-6A/M2)
101000 MAXIMUM INTERNAL TEMP.....24CCF (13CCC)
101000 MAX EXTERNAL SURFACE TEMP.....11CF (45C)



186000 FCT ZONE DIMENSIONS(MIN).....11A,12,5CM DIA, 4X4IN(10CPI) LONG
186000 FULL LENGTH FRONT OPENING DCR
186000 VIEW PORTS.....FRONT, REAR, 2 SIDES, TOP, BOTTOM, TOP
186000 FLANGE OPENINGS.....TOP, BOTTOM, 2 EACH SIDE, 2 REAR
186000 ZONE REFINER.....SF002 11112
186000 THIS UNIT WILL HEAT, FOLD AND MANIPULATE A SPECIMEN FOR ZONE
186000 REFINING AND MELTEN ZONE CRYSTAL GROWTH. IT WILL CAUSE A POLYMER
186000 REGION TO TRAVERSE ALONG A LONG CYLINDRICAL ROD OF MATERIAL.
186000 SPECIMEN SIZE.....4 TO 12 IN (10 TO 30 CPI) LONG
186000 C.1 TO 0.5IN(0.25 TO 1.25CPI) DIA
186000 SCANNING DISTANCES.....LP TO 12IN(3CPI)
186000 SCANNING SPEED.....LP TO 20IN/MIN (50CPI/MIN)
186000 SPECIMEN ROTATION.....C TO 6C RPM
186000 OPERATING TEMPERATURE.....575 TO 2500F (300 TO 1600C)
186000 GRACIENT FURNACE.....SF002 11113
186000 THERMAL CHAMBER.....SF003 21111
186000 CHAMBER.....LS027 21111
186000 CAPABILITY TO HEAT, MAINTAIN CONSTANT TEMPERATURE, AND DRY TEST
186000 EQUIPMENT, LIQUIDS AND TEST SPECIMEN.
186000 TEMPERATURE RANGE32 TO 450 C F (10 TO 233 C)
186000 CAPACITY1 CU FT (0.03 CU M)
187000 BENCH, LUMINAR AIRFLOW.....LS934 11114
187000 A CLOVE BOX WITH RELATIVELY HIGH AIR FLOW FOR CONTROL OF
187000 PARTICULATE AND GASEOUS CONTAMINANTS WITHIN
187000 BENCH, GENERAL EXPERIMENTS.....LS935 11134
187000 PROVIDE WORK AREA FOR PREPARATION OF EXPERIMENTS, EXPERIMENT
187000 MAINTENANCE OBSERVATIONS, AND PREPARATIONS FOR RETURN TO EARTH
187000 ELECTRICAL UTILITY.....28VDC, 40CHZ AL, 60HZ AC
187000 VACUUM UTILITY.....10E-6 TORR
187000 PRESSURE SOURCE.....20 PSIG (3.45E4 N/M2)
187000 AVAILABLE GASES.....OXYGEN, NITROGEN, CARBON DIOXIDE
187000 WORK BENCH.....ST182 21114
187000 PROVIDE A WORK AREA UPON WHICH EXPERIMENTS AND MAINTENANCE MAY BE
187000 PERFORMED. ENCLOSED STORAGE TO BE INCORPORATED IN BENCH. USFC
187000 CN ST01S, ST02S, ST03S, ST04S AND ST05S.
187000 SIZE APPROXIMATELY2X3X4 FT (0.61X0.91X1.2 M)
187110 LITTER CHAIR - ROTATING.....LS114 11131
187110 GRAVITY INDUCING DEVICE AND HORIZONTAL WORK TABLE.
187110 ROTATION SPEEDC TO 6C RPM
201000 ULTRASONIC CLEANER (CCRI).....LSC47 11111
201000 HIGH FREQUENCY ENERGY FOR CLEANING EXPERIMENT AND MAINTENANCE
201000 EQUIPMENT.
201000 OPERATING FREQUENCY20000 HZ
201000 CAPACITY1 CU FT
204200 HEAT EXCHANGER.....ECC01 11121
204200 EXCHANGER TO TRANSFER HEAT FROM THE CLCUC PHYSICS LAB TO THE
204200 SPACELAB HEAT REJECTION SUBSYSTEM
204200 HEAT TRANSFER RATE.....165 BTU/HR (54 WATTS THERMAL)
204200 SPACELAB COOLANT INLET TEMP.....43F (4C)
204200 CPL COOLANT OUTLET TEMP.....49F (7C)
205000 HEAT EXCHANGER.....ECC01 11121
205000 PROVIDE HOT FLUID TO CPL
205000 CPL COOLANT OUTLET TEMP.....100F (38C)
205000 HEAT RATE (EST).....170BTU/HR (50 WATTS THERMAL)
207500 HEMATOCRIT, ELECTRONIC.....LS917 11111
207500 AUTOMATICALLY DETERMINE THE PERCENT HEMOCRIT IN BLOOD. RATIO &
207500 READOUT AND TEMPERATURE COMPENSATOR REQUIRED
211000 HUMIDIFIER ASSEMBLY.....ECC03 11135
211000 THIS ASSEMBLY CONSISTS OF AN EVAPORATOR, A HEATER, A WATER TRAP,
211000 A METERING PUMP AND A WATER TANK WITH AN INTEGRAL BLACER. THIS
211000 ASSEMBLY GENERATES THE WATER VAPOR NECESSARY FOR CLCUC CREATION.
211000 DRY AIR FLOW RATE.....175CFM (5 LITERS/MIN)
211000 MAXIMUM TEMPERATURE.....77F (45C)
211000 HUMIDITY RANGE.....C TO 100%
211000 HEATER CAPACITY(MIN).....170BTU/HR (50 WATTS THERMAL)
211000 HEATER FEED RATE (MAX)0.02 L/HR (9.1 CC/HR)
211000 WATER TANK CAPACITY.....2LRS (1.35KG)
211000 VAPORIZER - WATER.....ST048 21135
211000 CAPABILITY TO PROVIDE COOL WATER VAPOR TO EXPERIMENT CONTAINER.
211000 USED IN ST01S, ST02S AND ST03S.
211000

ORIGINAL PAGE IS
OF POOR QUALITY

211000 STEAM GENERATOR ST101 11135
211000 DEVISE TO GENERATE STEAM IN ZERO GRAVITY
211000 DEW POINT SENSOR ECC03 11112
211000 THIS SENSOR MEASURES THE DEW POINT OF THE GAS SAMPLE
211000 MEASUREMENT RANGE.....1C TO 500
211000 ACCURACY.....+/-1%
211000 LIQUID WATER CONTENT METER ECC03 11235
211000 MONITORS LIQUID WATER CONTENT WITHIN THE EXPANSION CHAMBER FOR
211000 THE PURPOSE OF ACCURATE WATER BUDGET ACCOUNTS. CURRENT CONCEPTS
211000 USE OPTICAL TECHNIQUES TO MEASURE THE WATER DROPLETS IN AIR.
211000 ACCURACY.....+/-0.05%
211000 DROPLET SIZE DISTRIBUTION METER ECC09 11235
211000 THIS UNIT WILL DETERMINE DROPLET SIZE DISTRIBUTIONS WITHIN A
211000 CHAMBER UTILIZING OPTICAL TECHNIQUES
211000 GAS ANALYZER, WATER VAPOR L5663 11111
211000 MONITOR WATER VAPOR CONTENT
211000 MOISTURE RANGE.....C TO 50% RELATIVE HUMIDITY
211000 (C.CC1 TO 2000MICROGRAMS/LITER)
211000 TEMPERATURE RANGE.....-100 TO 140F (-110 TO 60C)
211000 RECIRCULATING FLUID INCUBATOR SF001 11112
211000 CELLS/TISSUE HOLDING UNIT LS1016 12114
211000 GERMINATION HOLDING UNIT/INCUBATOR CONTAINING ELECTRONIC FLUO-IN
211000 CAPABILITY FOR MEASUREMENT DEVICES
211000 TEMPERATURE RANGE 41 TO 140 F (5 TO 60 C)
211000 PCO2 10C +/-5 MM
211000 PH2 0C +/-2C MM
211000 PCO2 3 MM OR LESS
211000 TOTAL PRESSURE 70C +/-2C MM
211000 RELATIVE HUMIDITY 0C TO 100 PCY
211000 CELLULAR GROWTH ENVIRONMENTAL CHAMBER ST001 11112
211000 HOUSING FOR COLONIES OF BACTERIA WHICH CAN BE OBSERVED DURING
211000 EXPERIMENT. USED IN ST015, ST025 AND ST035.
211000 SIZE APPROXIMATELY 5 X 5 X 7 IN (13 X 13 X 18 CM)
211000 INCUBATOR ST003 11111
211000 CONTAINER USED TO PROVIDE VARIABLE TEMPERATURE ENVIRONMENTS FOR
211000 MICROSCIENCE SPECIMENS. USED IN ST015, ST025 AND ST035.
211000 SIZE APPROXIMATELY 2X1.5X0.5 FT (60X15X15 CM)
211000 TEMPERATURE TO 100 F (40 C)
211000 DATA BUFFERS, FORMATTERS PEC 21111
211000 BUFFERS DATA FROM HIGH RESOLUTION NARROWBAND MULTISPECTRAL SCANNER
211000 INPUT DATA RATE.....200MBPS
211000 OUTPUT DATA RATE.....100MBPS
211000 INPUT DATA DUTY CYCLE.....33%
211000 STORAGE PERIPHERALS SP001 11111
211000 THIS UNIT ACCEPTS ADDRESSED BINARY DATA FROM THE DIGITAL
211000 PROCESSOR AND TRANSLATES THIS DATA INTO EQUIVALENT ANALOG FORM.
211000 MAXIMUM SIGNAL AND HOLD CHANNELS,32
211000 SIGNAL AND HOLD CHANNELS PER CARD,4
211000 SIGNAL FORMATTER CN070 21111
211000 PROVIDES DATA FORMATTING, ACCUMULATES PRE/POST AMBLE, CLOCKS
211000 START OF EACH TEST AND DATA ACQUISITION.
211000 DATA FORMAT 8, 16, 32, BITS/WORD
211000 DATA RATE 125000 WORDS PER SECOND
211000 DATA BUFFER ST010 21111
211000 CAPABILITY TO ACCEPT VERY HIGH DATA RATES, PROVIDE SHORT DURATION
211000 STORAGE AND LOWER RATE DUMP CAPABILITY. USED IN ST015.
211000 WORD LENGTH 8 TO 32 BITS
211000 OUTPUT RATE VARIABLE AT 0 TO 5 MCC
211000 STORAGE DURATION INFINITE
211000 PATCH PANEL, COAXIAL AF007 01111
211000 MULTIPLE INPUT PANEL CAPABLE OF CHANNELING HIGH FREQUENCY RF
211000 ENERGY TO APPROPRIATE INSTRUMENTS WITH MINIMUM LOSS AND INTER-
211000 FERENCE.
211000 COUPLER - ANTENNA (0.2 TO 2.0 MHZ) AF009 11111
211000 ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.
211000 FREQUENCY BANDWIDTH 0.2 TO 2.0 MHZ



219900 COUPLER - ANTENNA (2.0 TO 20.0 MHZ) AF910 11111
219900 ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.
219900 FREQUENCY BANDWIDTH 2.0 TO 20.0 MHZ
219900 COUPLER - ANTENNA (0.2 TO 200 KHZ) AF911 11111
219900 ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.
219900 FREQUENCY BANDWIDTH 0.2 TO 200 KHZ
219900 PATCH PANEL AP914 51111
219900 INTERFACE CONNECTIONS OF VARIOUS RF RECEIVERS WITH SELECTED TEST EQUIPMENT FOR WAVE ANALYSIS.
219900 MCEM C4078 11114
219900 MODULATE AND DEMODULATE VOICE OR DATA TRANSMISSIONS WITH OR WITHOUT SPREAD SPECTRUM CAPABILITY. PROVIDES BASELINE COMMUNICATION SYSTEM, ANALOG OR DIGITAL VOICE, AND DIGITAL DATA.
219900 DATA 15 BS TO 50 KHS (VARIABLE)
219900 BANDWIDTH MOVEABLE SS
219900 BANDWIDTH 2 MHZ
219900 DEMODULATOR C4079 31111
219900 DEMODULATE VOICE OR DATA TRANSMISSIONS WITH OR WITHOUT SPREAD SPECTRUM CAPABILITY. PROVIDES A NUMBER OF RECEIVERS CONNECTED TO DIFFERENT ANTENNAS.
219900 MCEM, WIDEBAND C4080 21235
219900 MODULATE AND DEMODULATE WIDEBAND FM OR DIGITAL SIGNALS. WIDEBAND FM DETECTOR WITH THRESHOLD EXTENSION. DIGITAL DEMOD WITH FIT, SYNC, AND ERROR CORRECTION CODING. MODULATOR PERTICA HAS COMPLEMENTARY CHARACTERISTICS.
222000 DEIONIZER, WATER 15953 11111
222000 PROVIDES A MEANS OF REMOVAL OF CATIONS AND ANIONS FROM WATER SUPPLIES REQUIRED FOR SPECIALIZED LABORATORY ANALYSIS
222000 SOLIDS REMOVAL CAPABILITY.....<1 FPM
222000 IONS REMOVAL CAPABILITY.....<1 PPM
225700 OPTICAL COLLIMATOR C4086 11111
225700 PROVIDES FOR OPTICAL ALIGNMENT OF THE LASER AND TELESCOPE. ALIGNMENT CHECKS TO BE BOTH CIRCULAR AND PULSED.
225700 BEAM EXPANDER OPTICS C4088 21111
225700 PROVIDES VARIABLE BEAM DIVERGENCE CONTROL FOR LASER TRANSMITTER AND THE OPTICAL BEACON TO FACILITATE ACQUISITION.
225700 TYPE AFCCAL WITH VARIABLE ELEMENT SEPARATION.
225700 INPUT BEAM SIZE 0.25 INCH DIA
225700 EXIT BEAM SIZE 1 INCH DIA
225700 INPUT BEAM CHARACTERISTICS
225700 LASER TRANSMITTER TEM WAVE, GAUSSIAN DIST
225700 OPTICAL BEACON 3 RAD DIVERGENCE, GAUSSIAN DIST
225700 OUTPUT BEAM DIVERGENCE
225700 CO2 LASER TRANSMITTER DIFF LIMIT TO 200 ARC-SEC
225700 Nd:YAG LASER TRANSMITTER DIFF LIMIT TO 50 ARC-SEC
225700 OPTICAL BEACON 2 ARC-MIN TO 1 DEGREE
225700 ADJUSTMENT RESPONSE TIME 0.5 SEC MAX
225700 BEAM REFLECTOR C4089 41111
225700 THE TYPE PROVIDES A VERNIER POINTING CONTROL WITH FASTER RESPONSE AND GREATER PRECISION THAN THE MAIN OPTICS GIMBALS. THE OTHER TYPE PROVIDES A SMALL OFFSET POINTING CAPABILITY BETWEEN THE TRACKER AXIS AND THE TRANSMITTER AXIS.
225700 DEFLECTION CAPABILITY TWO AXIS
225700 DEFLECTION RANGE +/- 200 MICRORAD
225700 PRECISION +/- 0.2 MICRORAD
225700 SPECTRAL RANGE ACCUMULATE LASER TRANSMITTER, RECEIVER AND BEACON.
227000 LASER ECC13 11111
227000 POLARIZED, COLLIMATED LIGHT SOURCE USED FOR PARTICLE HEATING AND POSITION CONTROL AS WELL AS TO ASSESS OPTICAL SCATTERING PROPERTIES OF ICE CRYSTALS
227000 LASER SOURCE SP004 11111
227000 PROVIDE HEAT TRANSMITTED TO A SPECIMEN BY LASER BEAM. CONCEPT PROVIDES HEATING WITHOUT CONTAMINATING SPECIMEN. UNIT IS USED IN CONJUNCTION WITH CHEST-GENERAL PURPOSE ENCLOSURE
227000 POWER.....KVA
227000 COOLING CAPACITY.....LPM
227000 LASER TYPE.....ANY OF FOLLOWING FIVE:
227000CONTINUOUS WAVE HYDROGEN-FLUORINE CHEMICAL LASER, IR OUTPUT
227000CONTINUOUS WAVE FIBER COATED TION CO2, IR OUTPUT
227000CONTINUOUS WAVE CO-CO2-P2O5 CHEMICAL LASER, IR OUTPUT
227000PULSED TEA CO2, IR OUTPUT

ORIGINAL PAGE IS
OF POOR QUALITY

```

227000 .....PULSEC TEA REACH, LV CLPUT
227000
227000 LASER ST110 11111
227000 PCDE-LOCKED LASER RANGING SYSTEM MEASURING RANGE, LINE-OF-SIGHT
227000 ANGLES, AND RANGE RATE BETWEEN ATL/SHUTTLE AND SEPARATE SPACE
227000 VEHICLES INCLUDING RANGE OF ATL TO GROUND.
227000
227000 RANGE TOLERANCE ..... +/- 1.10 IN (+/-3 CM) ATL-GAC
227000 COOPERATIVE TARGET DISTANCE .. 3CC NM (556 KM)
227000
227000 LASER ASSEMBLY - NC1YAG CA053 11111
227000
227000 PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFINE
227000 AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMMUNICATIONS
227000 APPLICATIONS. CONTAINS OPTICAL TRANSMITTER INCLUDING BEAMWIDTH
227000 CONTROL, MODULATOR, AND COOLING SYSTEM.
227000
227000 LASER CAVITY
227000 WAVELENGTH ..... 1.06 MICRONS
227000 LASING MODE ..... FUNDAMENTAL (TEM..)
227000 OPTICAL OUTPUT POWER ..... 1 WATT AVERAGE
227000 EFFICIENCY ..... 1 PERCENT
227000 MODULATOR
227000 TYPE ..... DIGITAL
227000 DATA RATE ..... 1 GRPS
227000 EXTINCTION RATIO ..... 2CDB
227000 COOLING SYSTEM
227000 COOLING METHOD ..... CONDUCTIVE, NO LIQUIDS ALLOWED
227000 HEAT LOAD ..... 2CC WATTS
227000
227000 LASER LINK, DOUBLED NC1YAG CA054 11111
227000
227000 PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFINE
227000 AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMMUNICATIONS
227000 APPLICATIONS. CONTAINS OPTICAL TRANSMITTER INCLUDING BEAMWIDTH
227000 CONTROL, MODULATOR, DOUBLER, AND COOLING SYSTEM.
227000
227000 LASER CAVITY
227000 WAVELENGTH ..... 0.53 MICRONS
227000 LASING MODE ..... FUNDAMENTAL (TEM..)
227000 OPTICAL OUTPUT POWER ..... 0.1 WATTS AVERAGE
227000 EFFICIENCY ..... 0.1 PERCENT
227000 MODULATOR
227000 TYPE ..... DIGITAL
227000 DATA RATE ..... 1 GRPS
227000 EXTINCTION RATIO ..... 2CDB
227000 COOLING SYSTEM
227000 COOLING METHOD ..... CONDUCTIVE, NO LIQUIDS ALLOWED
227000 HEAT LOAD ..... 2CC WATTS
227000
227000 LASER REACON, NC1YAC CA055 11111
227000
227000 SERVE AS A LIGHT SOURCE TO FACILITATE ACQUISITION AND/OR TRACKING
227000 FROM A DISTANT COMMUNICATION TERMINAL.
227000
227000 REACON TYPE ..... Q-SWITCHED DOUBLED, PLTIPCE
227000 CAVITY.
227000 WAVELENGTH ..... 0.53 MICRONS (PREFERRED)
227000 REP RATE ..... 5 PPS
227000 ENERGY/PULSE ..... 2C MJ AT 0.53 MICRONS
227000 PULSEWIDTH ..... 2C NSEC
227000 PUMP METHOD ..... FLASH PLMP
227000 CODING TECHNIQUE ..... PULSE POSITIONAL MODULATION
227000 BEAMWIDTH ..... VARIABLE 1 DEG TO 1 ARC MINUTE
227000 LIFETIME ..... 1C MEG FLASHES OR GREATER
227000 POWER CONSUMPTION ..... 1CC WATTS
227000
227000 EYE LASER/FLASH LAMP SP001 31111
227000
227000 THIS UNIT WILL BE USED FOR FREE RADICAL GENERATION, SURFACE
227000 DAMAGE THRESHOLD DETERMINATION, AND HOLOGRAPHIC MICROSCOPY
227000
227000 ACCURACY
227000 WAVELENGTH ..... 0.01%
227000 VOLTAGE OF CAPACITOR BANK ..... 0.1%
227000 LASER AND FLASH LAMP ENERGY
227000 OUTPUT ..... 1.0%
227000 SYNC. ACCURACY ..... 1CCNS
227000 COHERENCE ..... 20%
227000 RANGE OF OPERATION
227000 ENERGY INPUT ..... 1CC-5CCJ/PULSE (1.14-1.44W/PULSE)
227000 ..... 1 AT LOW REPS, LARGE LAMP
227000 ..... 1C-5CCJ/PULSE (0.014-0.14W/PULSE)
227000 ..... AT HIGH REPS, SMALL LAMP
227000 ENERGY OUTPUT ..... 1-3CJ/PULSE (3E-4 - 3E-3W/PULSE)
227000 ..... 1 AT LOW REPS, EYE PUMP, FREQ.
227000 ..... 25CMJ TO 7J/PULSE (7E-5 - 2E-3
227000 ..... W/PULSE) FREQUENCY DOUBLED.
227000 ..... 1CCMJ TO 3J/PULSE (3E-5 TO 3E-4
227000 ..... W/PULSE) AT HIGH REPS AT EYE
227000 ..... FUNDAMENTAL, 25 - 750PJTE-4 -
227000 ..... 2E-4 W/PULSE FREQUENCY DOUBLED
227000 REPETITION RATE ..... LARGE LAMP - SINGLE TO 10/PIA
227000 ..... SMALL LAMP - SINGLE TO 20HZ
227000 FREQUENCY RANGE ..... FUNDAMENTAL 430 - 750NM (1.7E-5
227000 ..... TO 3E-5 IN.) FREQUENCY DOUBLED
227000 ..... 215 - 375NM (0.0E-5 TO 1.5E-5
227000 ..... IN.)
227000 FLASH DURATION ..... 1CC - 1CCNS
227000 HOLOGRAPHY ..... 250M (1.4 IN.)
227000 ..... SINGLE PRAPE TO 20FRAMES/SEC.
227000
227000 FLASH LAMP SP001 31111
227000
227000 THE FLASH LAMP HAS TWO FUNCTIONAL REQUIREMENTS: (A) TO PRECEDE

```


227800 FREE RADICALS, AND BIAS A PLMP SOURCE FOR EXCITING THE EYE
227800 LASER
227800 TYPE.....PERCUTANEOUS
227800 OPERATING TEMPERATURE.....20 TO 77F (12C TO 25C)
227800 WAVELENGTH.....1.2E-4 TO 2.4E-5 IN. (300-6000)
227800 ENERGY INPUT.....1.5W (50J)
227800 (MIN).....0.2W (10J)
228041 LASER ASSEMBLY - CC2 C4052 11111
228041
228041 PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFIN
228041 AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMPLICATIONS
228041 APPLICATIONS, CONTAINS OPTICAL TRANSMITTER INCLUDING REFRACTI
228041 CENTRAL, MODULATOR, AND COOLING SYSTEM.
228041
228041 LASER CAVITY
228041 WAVELENGTH.....1.2E-4 MICRONS
228041 LASING MODE.....FUNDAMENTAL (TEM..)
228041 OPTICAL OUTPUT POWER.....5 WATTS AVERAGE
228041 EFFICIENCY.....5 PERCENT
228041 MODULATOR
228041 TYPE.....DIGITAL
228041 DATA RATE.....200 MBPS
228041 EXTINGUISH RATIO.....20 DB
228041 COOLING SYSTEM
228041 COOLING METHOD.....CONDUCTIVE, NO LIQUIDS ALLOWED
228041 HEAT LOAD.....200 WATTS
228041
229500 DARK FIELD ILLUMINATOR SFOC1 21111
229500
229600 LIGHTING SYSTEM - PLANT LS1018 21114
229600
229600 PROVISIONS FOR LIGHT PLANT GROWING AREAS.
229600
229600 TYPE LIGHT.....DAYLIGHT-FLUORESCENT
229600 ILLUMINATION.....102 FT-C (1100 LUX/PM)
229600
229600 CELLS/TISSUE HOLDING UNIT LIGHTING LS1017 21114
229600
229600 PROVISION FOR PROVIDING COOL DIFFUSE LIGHT
229600
229600 ILLUMINATION.....50 +/-10 FT-C
229600 TEMPERATURE.....AMBIENT
229600
229600 LIGHT SOURCE ST049 21111
229600
229600 PORTABLE LIGHT SOURCE USED FOR PHOTOGRAPHING EXPERIMENT CREFA
229600 TIONS AND RESULTS. USED IN STC1, STC2 AND ST03.
229600
229600 ILLUMINATION.....100 TO 1000 LUX/PM
229600 DISTANCE.....LP TO 30 FEET
229600
234000 SPACE DEVICE LIGHTING F0009 11114
234000
234000 PROVIDES LIGHTING FOR VARIOUS OPTICAL RECORDING ACTIVITIES.
234000 CONSISTS OF CONTINUOUS LIGHTING, STROBE LIGHT AND LASER BEAM LIGHT
234000 SOURCES
234000
234000 HIGH-INTENSITY LIGHT ILLUMINATION, 415FC
234000
234000 STROBE LAMP ST157 11111
234000
234000 VARIABLE DISCHARGING STROBE LIGHTS WHICH ALLOW THE STOPPING OF
234000 ACTION.
234000
234000 STROBE DISCHARGE RATE.....1 TO 1000 FLASHES/SEC
234000
235700 LYOPHILIZATION UNIT SFOC1 11111
235700
235700 LYOPHILIZER LS972 11111
235700
235700 PROVIDE FREEZE DRYING OF SMALL SPECIMENS
235700
235700 SAMPLE CAPACITY.....20-100 IN. (1600)
235700 TEMPERATURE RANGE.....-40F TO 252F (-40 TO 122C)
235700
239500 MICRO, SUPERCONDUCTING ST150 11114
239500
239500 MAGNET FOR DROP/PARTICLE POSITIONING
239500
256000 MICROPHONE LS1007 11111
256000
256000 MICROPHONES USED FOR BIOLOGY, MICROSCOPY AND MANAED SYSTEMS
256000 INTEGRATION
256000
263000 MICROSCOPE, DISSECTING LS1004 11111
263000
263000 PERMIT DETAIL DISSECTING OF BIOLOGICAL SPECIMENS SUCH AS PLANTS,
263000 ANIMAL TISSUE AND ORGANS
263000
263000 STEREO EYEPIECE
263000 ADJUSTABLE BRIGHT AND DARK FIELD ILLUMINATION
263000
263000 MICROSCOPE, COMPOUND LS1005 11111
263000
263000 GENERAL PURPOSE BINOCULAR MICROSCOPE FOR MICROSCOPIC STUDIES OF
263000 TISSUES
263000
263000 IDENTIFICATION.....10X TO 100X
263000 PHOTOGRAPHIC CAPABILITY.....POLAROID FILM
263000 LIGHTING.....DARK FIELD OR LIGHT FIELD AND
263000 PHASE CONTRAST
263000
264650 IF MICROSCOPE E0017 11111
264650

ORIGINAL PAGE IS
OF POOR QUALITY



204690 DETERMINE SURFACE TEMPERATURE OF DROPLETS AND ICE CRYSTALS
204690 SENSITIVITY.....+/-C.34F(+/-C.3C) AT 5F (-15C)
204690
205000 MICROSCOPE STC51 21111
205000
205000 MICROSCOPIC EXAMINATION OF AEROSOL PARTICLES IN THE VISUAL RANGE.
205000 ALSO COMPATIBLE WITH BIOLOGICAL AND METALLURGICAL EXPERIMENTS.
205000 USED ON ST01, ST02S AND STC3S.
205000
205000 MAGNIFICATION RANGE 5C TO 1000 X
205000 PHOTOGRAPHIC CAPABILITY YES
205000
205000 MICROSCOPE EC017 11111
205000
205000 MAGNIFYING OBSERVATIONS
205000
205000 MAGNIFICATION FACTOR.....1C TO 1CC
205000 OBSERVATION MODE.....VISUAL OR PHOTOGRAPHIC
205000
205000 MICROSCOPE (HIGH MAGNIFICATION) EC017 11111
205000
205000 MAGNIFYING OBSERVATIONS
205000
205000 MAGNIFICATION FACTOR.....1CC X 10CC
205000 OBSERVATION MODE.....TRINOCULAR(CAMERA ATTACHMENT)
205000
205000 RETRACTION/RESTRUCTION HIGH RESOLUTION MICROSCOPE SF001 31231
205000
205000 STEREO MICROSCOPE EC017 11111
205000
205000 MAGNIFYING OBSERVATIONS
205000
205000 MAGNIFICATION FACTOR.....1C TO 1CC
205000 OBSERVATION MODE.....VISUAL OR PHOTOGRAPHIC
205000 WORKING DISTANCE.....UP TO 4IN(10CM)
205000
205000
205000 MICROSCOPE (CORE) LS021 11111
205000
205000 CAPABILITY TO THINLY SLICE EXPERIMENT SPECIMEN.
205000
205000 SECTIONING 1 TO 50 MICROMETERS
205000
205000 ATTENUATOR CALIBRATOR CNO28 41111
205000
205000 DYNAMIC RANGE ADJUSTMENT, REMOTE ELECTRICAL CONTROL, DIGITAL READ
205000 OUT OF VALUE, MANUAL OVER-RIDE.
205000
205000 IMPEDANCE 5C CMPS
205000 FREQUENCY DC TO 1.5 MHZ
205000 RANGE C TO 60 DB
205000
205000 METER, ELECTRIC ST151 11111
205000
205000 METER FOR CHAMBER DRIVE SYSTEM
205000
205000 EXPANSION CLOUD CHAMBER EC014 11135
205000
205000 CHAMBER PROVIDING SIMULATION OF LONG TERM, NATURAL CLOUD CLIPPING
205000 AND ACCELERATED EXPANSION
205000
205000 DIMENSIONS.....11.01A(30CM) DIA. X 17.71A
205000 (45CM) LONG
205000 VOLUME.....1.12CL.FT.(31.8L)
205000 OPERATING TEMPERATURE.....-76F TO 104F (-60 TO 40C)
205000 UPPER AND LOWER ENDS TEMP
205000 DIFFERENTIAL.....0 TO 10FIC TO 1CC
205000 END TEMPERATURE TOLERANCE.....+/-C.1F (+/-C.5C)
205000 OPERATING PRESSURE.....1.7 TO 14.7 PSIA(140 TO 760TCRR)
205000 PRESSURE RATE OF CHANGE.....19.4PSI/SEC (1000TCRR/SEC)
205000 PRESSURE RATE TOLERANCE.....+/-10
205000 MAXIMUM VOLUME EXPANSION.....V/V OF C.3+/-0.15
205000 CIRCULAR ENDS CONTAIN ELECTRICALLY CONDUCTIVE FLATES TO SERVE
205000 AS EQUIPOTENTIAL SURFACES FOR THE ELECTRIC FIELD ACTION
205000 CONTROL SYSTEM
205000
205000 CONTINUOUS FLOW DIFFUSION CLOUD CHAMBER EC014 11135
205000
205000 SUPERSATURATION IS CONTROLLED BY THE TEMPERATURE OF WATER
205000 COVERING THE UPPER AND LOWER SURFACES OF THE CHAMBER. IT WILL BE
205000 USED FOR CLOUD CONDENSATION NUCLEATION EXPERIMENTS
205000
205000 DIMENSIONS.....11.01A(30CM) X 11.01A(30CM) X
205000 21A(5CM)
205000 VOLUME.....C.16CL.FT.(4.5L)
205000 INNER PLATE DIMENSIONS.....11.6IN(30CM) X 10IN(25CM) SPACE
205000 C.9IN(1.3CM) APART
205000 ISOTHERMAL TEMP CONTROL.....+/-C.1F(0.05C)
205000 UPPER AND LOWER SURFACE
205000 TEMPERATURE DIFFERENCE.....0 TO 10F TO 10C
205000 SURFACE TEMP CONTROL.....+/-C.1F(0.05C)
205000 OPERATING PRESSURE.....14.3 TO 14.7PSIA(740 TO 760TCRR)
205000 PRESSURE MEAS. ACCURACY.....+/-C.28
205000 MAX PRESSURE CHANGE RATE.....C.0011/SEC(0.47CRR/SEC)
205000
205000 STATIC DIFFUSION ICE CLOUD CHAMBER EC014 11135
205000
205000 THIS IS A NAMURA TYPE CHAMBER WHICH UTILIZES ICE SURFACES TO
205000 PROVIDE CONTROLLED SUPERSATURATION RELATIVE TO ICE.
205000
205000 DIMENSIONS.....79.79IN(40CM) DIA X 7.91A(20CM)
205000 LONG
205000 VOLUME.....C.444CL.FT.(14.5L)
205000 OPERATING TEMPERATURE.....-4C TO 77F(-40 TO 25C)
205000 TEMP MEAS. ACCURACY.....+/-1.0F(+/-1.0C)



297500 UPPER AND LOWER SURFACE
297500 TEMPERATURE DIFFERENCE.....LP TO 34F (2CC)
297500 SURFACE TEMP. MEAS. ACCURACY..... $\pm 0.1F$ ($\pm 0.05C$)
297500 OPERATING PRESSURE.....1.5 TO 14.7PSIA(100 TO 760TCRR)
297500 PRESSURE MEAS. ACCURACY..... $\pm 0.15PSI$ ($\pm 1CTCRR$)
297500
297500 STATIC DIFFUSION LIQUID CLOUD CHAMBER ECC14 11134
297500
297500 A THIN TYPE CHAMBER USED FOR EXPERIMENTS REFLINING ABOVE
297500 FREEZING TEMPERATURES AND SUPERSATURATION OF THE LIQUID RELATIVE
297500 TO WATER. SUPERSATURATION IS CONTROLLED BY THE TEMPERATURES OF
297500 THE WATER COVERING CHAMBER UPPER AND LOWER SURFACES.
297500
297500 DIMENSIONS.....6IN. (15 CM) DIA. X 0.6IN. (1.5 CM) HIG
297500 VOLUME.....0.01CU. FT. (0.27L)
297500 OPERATING TEMPERATURE.....32F TO 68F (1C TO 30C)
297500 END SURFACE MAXIMUM
297500 TEMPERATURE DIFFERENTIAL.....18F (1CF)
297500 SURFACE TEMP. MEAS. ACCURACY..... $\pm 0.1V$ ($\pm 0.05C$)
297500 OPERATING PRESSURE.....2.84 TO 14.7PSIA(140 TO 760TCRR)
297500
297500 GENERAL CLOUD CHAMBER ECC14 11135
297500
297500 THIS CHAMBER WILL BE USED FOR MANY EXPERIMENTS THAT REQUIRE A
297500 RELATIVE HUMIDITY BELOW 100% AND MINIMUM TEMPERATURE CONTROL.
297500 PROVISIONS WILL BE MADE FOR GENERATING VARIOUS ELECTRIC FIELDS.
297500 POSITIONING DEVICES (SOUND, OPTICAL, ELECTRICAL), AND FIBRE
297500 DROPLET SIZING.
297500
297500 DIMENSIONS.....11.6IN. (30CM) ALL SIDES (CLBF)
297500 VOLUME.....0.55CU. FT. (15.7L)
297500 OPERATING TEMPERATURE.....1C TO 68F (1C TO 30C)
297500 OPERATING PRESSURE.....1.5 TO 14.7PSIA(140 TO 760TCRR)
297500 ELECTRIC FIELD PLATES LOCATED ON OPPOSITE SIDES
297500 ACOUSTIC DRIVERS LOCATED ON THREE MUTUALLY PERPENDICULAR SIDES
297500
297500 ENVIRONMENT CHAMBER ST050 11136
297500
297500 PROVIDE SPECIAL CHAMBER TO MONITOR OPTICAL PROPERTIES OF AEROSOLS
297500 WITH CAPABILITY TO PERTURBATE THE ENVIRONMENT WITH IONS AND
297500 OBSERVING THE AEROSOLS THROUGH A MICROSCOPE AND PHOTOGRAPHING THE
297500 SCENE. USED IN ST015, ST025 AND ST035.
297500
297500 CHAMBER, CUBICAL TEST ST188 11137
297500
297500 CHAMBER FOR PERFORMANCE OF SUPERFLUID HELIUM AND LIQUID DECF/
297500 PARTICLE EXPERIMENTS
297500
297500 DIMENSIONS.....1.8FT. (55CM) CUBE
297500
301000 RADIATION COUNTER - BIOCHEMICAL SAMPLE L5926 11111
301000
301000 FAST, ACCURATE AND EASY BIOCHEMICAL RADIATION COUNTING SYSTEM
301000 BASED ON STANDARD SIZED BIOCHEMICAL SAMPLES.
301000
301000 STANDARD PLANCHET UP TO 2 IN. DIA
301000 CAPACITY OF PLANCHETS 12C
301000
305000 NUCLEAR PARTICLE COUNTING UNIT SFG01 11111
305000
305000 CAPABILITY TO DISCRIMINATE AND COUNT ALPHA, BETA, GAMMA, X-RAYS
305000 AND K-ELECTRONS.
305000
305000 COUNTING RATE 3.5E+06 C/P
305000 RESOLUTION LOSS 1.8E+06 C/P
305000 PLATEAU CHARACTER
305000 ALPHA 5CC-12CC V & 1.8/100V
305000 BETA-GAMMA 17CC-15CC V & 1.8/100 V
305000 PRECISION ± 0.1 % FS
305000 TIMER 0.05 TO 2160 MINUTES
305000
305000 RADIATION DETECTOR L5919 11111
305000
305000 ACCURATE MEASUREMENT OF EITHER ALPHA, BETA, GAMMA, X-RAY OR
305000 NEUTRON RADIATION AS IT IS APPLIED TO A BIOLOGICAL SPECIMEN.
305000
305000 RATE RANGES 0 TO 0.1 PR/HR
305000 0 TO 1000 R/HR
305000 CALIBRATED READOUTS COUNTS/SEC, /MIN, AND TOTAL
305000
308500 DOSIMETER, RADIATION DETECTOR L5918 11111
308500
308500 PORTABLE DEVICE WHICH WILL ALERT USER TO RADIATION LEVELS IN
308500 EXCESS OF A PREDETERMINED LEVEL AND TO PROVIDE A DIRECT READOUT
308500 OF THE TOTAL CUMULATIVE RADIATION DOSE.
308500
308500 RATE MEASUREMENT RANGE 0 TO 0.1 PR/HR
308500 0 TO 50 R/HR
308500 ADJUSTABLE RADIATION ALARM 0.05 MR/HR TO 50 R/HR
308500 CUMULATIVE DOSE RANGES 0 TO 200 MR
308500 0 TO 500 R
308500 ACCURACY (X- & GAMMA RAY) ± 10 %
308500
308500 FILM RACE ST070 12111
308500
308500 MONITOR ASTRONAUT PERSONAL RADIATION DOSE, DOSE RATE AND DOSE
308500 HISTORY. USED IN ST015, ST025 AND ST035.
308500
308500 DOSE RATE 0 - 0.1 PR/HR
308500 0 - 1000 R/HR
308500
308500 DOSIMETER - THERMOLUMINESCENT ST071 41111
308500
308500 PASSIVE DOSIMETER USING THERMOLUMINESCENCE TO MONITOR RADIATION.
308500 USED IN ST015 AND ST025.
308500

ORIGINAL PAGE IS
OF POOR QUALITY



300900 ENERGY RANGE DOSE 1 TO 10E+05 RADS
300900
314000 GENERATOR - ICM ST047 11111
314000
314000 IONIZING AEROSOL PARTICLES SO AS TO OBSERVE ZERO GRAVITY PARTICLES.
314000 USED IN ST015, ST025 AND ST035.
314000
319900 ISOTOPE SOURCE - SELF CONTAINED LS920 11111
319900
319900 RADIATION STRESS SOURCE USED IN CONDUCTING RED KLECC CELL
319900 SURVIVAL STUDIES. TYPICAL ASSORTMENT - NITROGEN 15 CAPPCLACS,
319900 RADIOLABELED SUBCELLULAR PRODUCTS (RIBOSOMES), RADIOLABELED POLY-
319900 NUCLEOTIDES OR TITRATED TRYPTICINE.
319900
319900 RADIATION LEVEL 1 TO 100 MICRORAD
319900
320000 CONTROLLER, THREE-AXIS HAND AS008 11132
320000
320000 THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT
320000 POINTING AND INITIAL TARGET ACQUISITION.
320000
320000 CONTROLLED AXES 3
320000
320000 CONTROLLER, THREE-AXIS HAND ME008 11132
320000
320000 THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT
320000 POINTING AND INITIAL TARGET ACQUISITION.
320000
320000 CONTROLLED AXES 3
320000
320000 CONTROLLER, THREE-AXIS HAND SC008 11132
320000
320000 THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT
320000 POINTING AND INITIAL TARGET ACQUISITION.
320000
320000 CONTROLLED AXES 3
320000
320000 THREE AXIS MANIPULATOR SF002 21132
320000
320000 MECHANICAL SAMPLE PLACEMENT AND RETRIEVAL SF004 11111
320000
330000 COMPACTOR, WASTE SOLIDS LS945 11111
330000
330000 COMPACTS WASTE SOLIDS FOR EASY HANDLING
330000
330000 COMPACTION (MINI) 1/4 INITIAL VOLUME
330000
330000 WASTE STORAGE SYSTEM LS982 11111
330000
330000 PROVIDE FOR HANDLING, STORAGE AND DISPOSAL OF SOLID AND LIQUID
330000 WASTE MATTER FROM THE EXPERIMENTS.
330000
330000 CAPACITY 5 CL FT
330000
355900 OPTICAL DETECTOR ECC09 11235
355900
355900 SOLID STATE DETECTOR TO OBTAIN DATA CONCERNING SCATTERING
355900 PROPERTIES OF ICE CRYSTALS AND POSSIBLY EXTREME PARTICLE SIZE
355900 CHANGE DATA.
355900
355900 OPTICAL SENSORS EC013 12111
355900
355900 SENSOR FOR CONTROL OF PARTICLE CHARACTERISTICS AND OPTICAL
355900 ENERGY SOURCES
355900
355900 LASER OPTICAL SCATTERING MONITOR SF001 21111
355900
355900 THE LASER OPTICAL SCATTERING MONITOR WILL PROVIDE A REAL TIME
355900 MEASURE OF LIGHT SCATTERING IN A SPECIMEN AT ANY DESIRED MEASURE-
355900 ING ANGLE. THIS DATA WILL PROVIDE INFORMATION ON THE SIZE, SHAPE,
355900 ORIENTATION, INDEX OF REFRACTION, CONCENTRATION AND LOCATION OF
355900 OPTICAL SCATTERING CENTERS.
355900
355900 LASER TYPE LOW-POWER LV GAS DISCHARGE
355900 FILTER BANDWIDTH 4E-5 IN. (C.3AP) AT HALF POWER
355900 SCATTERED LIGHT INTENSITY
355900 ACCURACY 0.1%
355900 ANGULAR POSITION ACCURACY 1 MIN (0.3RADI)
355900
358900 OSCILLOSCOPE AF015 21111
358900
358900 MONITOR, MEASURE AND MAINTAIN ELECTRONIC EQUIPMENT OPERATING
358900
358900 BANDWIDTH 500 MHZ
358900 CHANNELS 2
358900 DEFLECTION FACTOR 2 PV
358900 PHOTOGRAPHIC CAPABILITY YES
358900
358900 OSCILLOSCOPE ECC06 11111
358900
358900 MONITOR, MEASURE AND MAINTAIN ELECTRONIC EQUIPMENT OPERATING
358900
358900 BANDWIDTH 100 MHZ
358900 CHANNELS 2
358900 SENSITIVITY 10 MV/DIVISION
358900 TIME BASE VARIABLE
358900 STORAGE CAPABILITY NC
358900
358900 OSCILLOSCOPE ST007 21111
358900
358900 MONITOR, MEASURE AND MAINTAIN ELECTRONIC EQUIPMENT OPERATING IN
358900 THE L-PANE RANGE (300 TO 1000 MHZ).
358900 USED IN ST015, ST025, ST035, ST045, ST055, ST065, ST075, AND
358900 ST095.
358900
358900 BANDWIDTH 350 TO 1500 MHZ



358500 CHANNELS 2
358500 DEFLECTION FACTOR 1C MV/DIVISION
358500 PHOTOGRAPHIC CAPABILITY YES
358500
358500 CECILLOSCOPE EFO17 11111
358500
358500 SUPPORT MAINTENANCE ACTIVITIES AND SPECIAL DATA ACQUISITION
358500
358500 FREQUENCY RANGE DC TO 10.0 MHZ
358500 OBSERVATION MODE VISUAL AND CAMERA (1 HR STORE)
358500
358500 CECILLOSCOPE SFOC1 11111
358500
358500 DISPLAY OF THE PRESENCE AND/OR NATURE AND FORM OF OSCILLATIONS
358500 OF IRREGULARITIES OF AN ELECTRIC CURRENT. STORAGE CAPABILITY OF
358500 A CRT DISPLAY IS DESIRED FOR PRELIMINARY OBSERVATIONS.
358500
358500 BANDWIDTH DC TO 10MHZ
358500 SENSITIVITY INTERNAL: 0.2CM DEFLECTION TO
358500 1 MHZ
358500 EXTERNAL: 25CMV P-P TO 15V P-P
358500 VIDEO STORAGE TIME 1 MCLR
358500
358500 CECILLOSCOPE - PERSISTENT CRT (CCRF) LS976 11111
358500
358500 ELECTRONIC MONITOR AND DATA MEASUREMENT CAPABILITY OF LIFE
358500 SCIENCE EXPERIMENT EQUIPMENT WITH A SCREEN HAVING LONG TERM
358500 RETENTION.
358500
358500 BANDWIDTH DC TO 5 MHZ (MINIPLP)
358500 CHANNELS 2
358500 DEFLECTION FACTOR 1 MV TO 1C V PER DIVISION
358500 PHOTOGRAPHIC CAPABILITY YES
358500 PERSISTENCE LP TO 6 MCLRS
358500
367000 OPTICAL PARTICLE COUNTER ECO11 11235
367000
367000 THIS DEVICE WILL BE A SIZE ANALYZER. THE ANALYZER CONSISTS OF A
367000 PHOTOMULTIPLIER, AN AMPLIFICATION STAGE AND SUPPORTING
367000 ELECTRONICS
367000
367000 SENSITIVE PARAMETER LIGHT SCATTERED FROM PARTICLES
367000 PARTICLE SIZE 0.3 MICROMETERS
367000
367000 ELECTRICAL PARTICLE COUNTER ECO11 11135
367000
367000 THIS ANALYZER MEASURES THE SIZE DISTRIBUTION OF PARTICLES. THE
367000 PARTICLES ARE IONIZED AND THEIR MOBILITY AS A FUNCTION OF
367000 ELECTRIC FIELD IS MEASURED TO GIVE AN INTEGRAL SIZE DISTRIBUTION
367000
367000 PARTICLE SIZE 0.01 TO 1.0 MICROMETER
367000 FLOW RATE 0.001 TO 1.0 CFPM (0.001 TO 1.0 LPM)
367000 INTEGRAL DISTRIBUTION CAPABILITY YES
367000
367000 NUCLEI MASS MONITOR SYSTEM ECO11 11235
367000
367000 THIS DEVICE MEASURES THE TOTAL PARTICULATE MASS PER UNIT VOLUME.
367000 THE SYSTEM USES A CRYSTAL OSCILLATOR WHICH CHANGES ITS RESONANT
367000 FREQUENCY AS PARTICLES ARE DEPOSITED ON IT.
367000
367000 PARTICLE COLLECTION ELECTROSTATIC PRECIPITATION
367000 PARTICLE SIZE 0.01 TO 20 MICROMETERS
367000
367000 PH MONITOR SFOO1 21111
367000
367000 MEASURE ACIDITY AND ALKALINITY OF SOLUTIONS. USED IN GF LAB ALSO
367000
367000 PH RANGE 4.2 TO 12
367000 SELECTABLE RANGES ANY 2
367000 ANY 5
367000 ANY 10
367000 READING STABILITY 0.02 PH
367000 LIQUID TEMPERATURE RANGE 32 TO 212 F (0 TO 100 C)
367000
367000 PH METER LS1016 11111
367000
367000 METER MEASURES HYDROGEN ION CONCENTRATION OF SOLUTIONS
367000
367000 PH RANGE 0 - 14
367000 MEASUREMENT ACCURACY 0.02PH
367000
367000 AUTOMATIC PHOTOGRAPHIC PROCESSOR SFOC1 11111
367000
367000 UNIT PROVIDES PROCESSED MICROGRAMS AS PART OF THE ELECTRO-OPTICAL
367000 IMAGING SYSTEM
367000
367000 CYCLE TIME BETWEEN 5 TO 20 SEC
367000 PLATE HOLDER SIZE 4 IN (10CM) X 5 IN (13CM)
367000 EFFECTIVE APERTURE AT PMTC
367000 PLATE 3.25 IN (8.3CM) X 4.0 IN (10CM)
367000 DATA OUTPUT PROCESSED MICROGRAMS
367000
367000 FRAME STORAGE UNIT SPOO1 11111
367000
367000 THE FRAME STORAGE UNIT PROVIDES NONDESTRUCTIVE READOUT, ERASE,
367000 AND STORAGE ON ELECTRONIC COMPAND AND WORKS IN EITHER CONTINUOUS
367000 OR PULSED MODE. IT PROVIDES A REAL-TIME, NO PHOTOGRAPHIC
367000 PROCESSING, HIGH SPEED, AUTOMATIC LIGHT ADJUSTMENT CAMERA CAPABLE
367000 OF PROVIDING REAL-TIME MOVIES OR FROZEN STILLS
367000
367000 PLANT GROWTH AND SUPPORT CONTAINERS LS1019 12111
367000
367000 VARIOUS SIZE CONTAINER TO ENVIRONMENTALLY HOUSE SEEDLINGS.
367000
367000 CAPACITY UP TO 20-CM PLANTS
367000

ORIGINAL PAGE IS
OF POOR QUALITY



396400 OPTICAM-FIELD AND FIXATION POINT RECORDER LS975 11111
396400 THIS INSTRUMENT PERMITS THE RECORDING OF THE VISUAL FIELD AS
396400 THE SUBJECT'S FIXATION WITHIN THE FIELD WITHOUT RESTRICTIONS OF
396400 HEAD MOVEMENTS
396400 ACCURACY.....2DEG OF ARC
396400 PLETHYSMOGRAPH, LIMB(INCLUDING CLOUTIER) LS1010 11111
396500 MEASURE CHANGES IN BLOOD VOLUME AND VASCULAR RESPONSES
396500 MEASUREMENT RANGE.....HIGH RANGE-15 TO 500 CMPS
396500 LOW RANGE-1.5 TO 50 CMPS
401000 IMPEDANCE PNEUMOGRAPH LS995 11111
401000 MEASURE BREATHING CYCLE CHARACTERISTICS
401000 RF POWER METER CA060 11111
405400 MEASURE RF POWER OF VARIOUS RF SOURCES.
405400 FREQUENCY RANGES 126 TO 130 MHZ
405400 136 TO 144 MHZ
405400 2025 TO 2120 MHZ
405400 13.4 TO 14.2 GHZ
405400 14.4 TO 15.95 GHZ
405400 ACCURACY 0.02 DB/10 DB
405400 LASER POWER METER CA090 11111
405400 MEASURE OUTPUT LASER POWER.
405400 WAVELENGTH 0.53, 1.06 AND 10.6 MICRONS
405400 ACCURACY..... +/- 5 PERCENT
405400 RANGE 3CC AM
405500 CC-CC CONVERTER(VOLTAGE CONTROLLED CC) ECC05 11111
405500 CONVERTS SPACELAB POWER TO VOLTAGE USABLE BY ELECTRICAL FIELD
405500 GENERATING COMPONENTS
405500 INPUT VOLTAGE.....28V DC
405500 OUTPUT VOLTAGE.....TBD
405500 CC-CC CONVERTER (LCM FREQUENCY AC) ECC05 11111
405500 CONVERT POWER INTO USABLE FORM FOR ELECTRIC FIELD GENERATING
405500 CAPACITOR
405500 FREQUENCY.....0 TO 10000
405500 WAVE SHAPE.....SQUARE
405500 AMPLITUDE.....0 TO 3000 V/P
405900 AC POWER CONTROLLER ECC17 11111
405900 REGULATES AND CONTROLS AC POWER
405900 REGULATION RANGE.....TBD
405900 DC POWER CONTROLLER ECC17 11111
405900 REGULATES AND CONTROLS DC POWER
405900 REGULATION RANGE.....TBD
405900 RF VARIABLE POWER SUPPLY CA071 21111
405900 PROVIDES SOURCE OF VARIABLE POWER.
405900 OUTPUT 100 V PEAK
405900 VOLTAGE GAIN 20 X
405900 DISTORTION 0.1 % THD
405900 FREQUENCY RANGE DC TO 30 MHZ
406900 RF INDUCTION PWR COND. (FIXING & DISPERSAL) SP002 21111
406900 RF INDUCTION POWER CONDITIONER (2000 - 20000) SP002 21111
406900 HIGH VOLTAGE POWER CONDITIONER SP001 11112
407000 PROVIDES CC VOLTAGE TO EITHER THE STATIONARY OR CENTRALEX FLOW
407000 ELECTROPHORETIC COLUMNS. REGULATION PLST REPEATS CONSTANT OVER
407000 A RELATIVELY LARGE CURRENT DEMAND.
407000 OUTPUT VOLTAGE REGULATION.....-0.1%
407000 VOLTAGE RANGE.....25 TO 2500 VDC
407000 REACTORS.....VOLTAGE LEVEL
407000 POWER OUTPUT TO COLUMNS.....20W SUSTAINED
407000 2000 PEAK
407000 POWER CALIBRATION UNIT CA072 21111
407000 PROVIDES KNOWN INCREMENTALLY CONTROLLABLE POWER LEVEL FOR
407000 STANCAIZATION.
407000 OUTPUT -100 TO +100 V
407000 LINE REGULATION 1 MV + 0.01 %
407000 LINE REGULATION 10 MV + 0.01 %
407000 LASER POWER SUPPLY CA007 21111
407000 CONVERTS AND REGULATES RAB SPACECRAFT POWER IN A MANNER SUITABLE
407000 FOR PUMPING A LASER. PRIOR LASER POWER CONDITIONER AND DISTRIBUTION
407000 FOR ALL LASERS.

407000 VOLTAGE OUTPUT
407000 REGULATION +/- C.1 %
407000 POWER 100 WATTS MAX
407000 TRANSIENT RECOVERY C.1 SECONDS
407000
407000 POWER SUPPLY ST171 21111
407000
407000 ELECTRICAL POWER REGULATION AND DISTRIBUTION AS REQUIRED AND AT
407000 VARIABLE VOLTAGE RANGES. USED IN STC25, ST035 AND ST045.
407000
407000 VOLTAGE C TO 24 VDC
407000 POWER 200 WATTS
407000 SENSITIVITY C.05 VOLT
407000
407500 ELECTRONIC DRIVER ASSEMBLY EC006 41114
407500
407500 THE ELECTRONIC DRIVER ASSEMBLY CONTROLS THE LIQUID DRIP
407500 GENERATOR. IT CONSISTS OF A GENERATOR POWER SUPPLY, A PARTICLE
407500 CHARGER POWER SUPPLY AND A DUAL PULSE GENERATOR.
407500
407500 HIGH VOLTAGE POWER CONDITIONER (17KV) SP002 31111
407500
408000 LOW VOLT/HIGH AMP POWER CONDITIONER (1CKW) SF002 31111
408000
408000 PRESSURE SENSORS EC002 51111
408000
408000 MEASURE PRESSURES THROUGH OUT THE SYSTEM
408000
408000 RANGES..... 2 TO 20 PSIA (100 TO 1000 TCRR)
408000 10 TO 50 PSI (6.9E4-34.5E4/P2)
408000 10 TO 300 PSI (34.5E4-2.1E6/P2)
408000 ACCURACY..... 1% OF FULL SCALE
408000
409500 SLMP COMPRESSOR EC004 11111
409500
409500 COMPRESS EXPERIMENT GAS SAMPLES INTO SLMP STORAGE TANK
409500
409500 COMPRESSION RATIO(PAX)..... 10 TO 1
409500
410000 RESERVOIR EC001 11113
410000
410000 STORES COOLANT.
410000
410000 STORAGE TANK - EARTH SAMPLE EC004 11113
410000
410000 THIS TANK STORES SAMPLES OF EARTH GASES. THE ASSEMBLY CONSISTS
410000 OF A PRESSURE VESSEL, A POSITIVE EXPLUSION REGULATOR AND A
410000 PLATE.
410000
410000 PRESSURE..... 1 ATU (760 MM HG)
410000 VOLUME..... 10 CU.FT. (0.28 CU.METERS)
410000
410000 GAS SAMPLE STORAGE TANKS EC004 51111
410000
410000 PRESSURE VESSELS FOR STORING EXPERIMENT GAS SAMPLES.
410000
410000 MAX PRESSURE(EST)..... 300 PSI (2.1E6 N/P2)
410000 VOLUME 1 CU.FT. (0.028 CU.METERS)
410000
410000 SLMP STORAGE TANK EC004 11111
410000
410000 STORE GAS SAMPLES AFTER EXPERIMENT EVALUATION.
410000
410000 MAX PRESSURE (EST)..... 150 PSI (1.05E6 N/P2)
410000 VOLUME..... 10 CU.FT. (0.28 CU.METERS)
410000
410000 BUFFER AND ELECTROLYTE SUPPLY TANKS SF001 11113
410000
410000 SPECIMEN/SAMPLE SUPPLY TANKS SP001 11113
410000
410000 WASTE LIQUID TANK SF001 11113
410000
410000 COOLANT SUPPLY TANK SP001 11113
410000
410000 VACUUM CATCH TUBE SP004 11235
410000
410000 FUEL TANK ST261 11135
410000
410000 STORE FUEL FOR COMBUSTION EXPERIMENT
410000
410000 VOLUME..... 2 CU.FT. (0.055 CU.M)
410000
410000 CRIDIZER TANK ST262 11135
410000
410000 STORE OXIDIZER FOR COMBUSTION REACTION
410000
410000 VOLUME..... 2 CU.FT. (0.055 CU.M)
410000
410000 INERT GAS SUPPLY ST263 11111
410000
410000 SUPPLY PRESSURE FOR EXPULSION OF FUEL AND CRIDIZER FROM STORAGE
410000 TANKS
410000
410000 VOLUME 0.22 CU.FT. (0.006 CU.M)
410000
410000 CONDENSATE TANKS ST102 21113
410000
410000 TANKS FOR STORAGE OF CONDENSATE
410000
410000 PRESSURE..... 20 PSI (1.38E5 N/P42)
410000
410000 STORAGE - GAS SOURCE ST045 41111
410000
410000 GAS SOURCE STORAGE BOTTLES CAPABLE OF STORING VARIOUS GASES AT
410000 HIGH PRESSURES. SIMILAR TO STANDARD R-BOTTLES. USED IN ST015.

ORIGINAL PAGE IS
OF POOR QUALITY



410000 STORES AND STORES.
410000
410000 VOLUME - PHYSICAL 9 CU FT EA (C.14 CL P)
410000 PRESSURE MAXIMUM 1000 PSI
410000
410000 PRESSURE TANK STC58 41111
410000
410000 PROVIDE PRESSURE SOURCE FOR WATER FEED SYSTEM TO STEAM GENERATOR
410000 EXPERIMENT
410000
410000 OPERATING PRESSURE..... 100 PSI (6.9X10E5 N/M²)
410000
410000 WATER RESERVOIR WITH HEATER STC59 11113
410000
410000 STORE WATER FOR STEAM GENERATOR EXPERIMENT
410000
410000 OPERATING PRESSURE..... 30PSI (2.07X10E5 N/M²)
410000
410000 TANKS - FUEL/OXIDIZER ST170 41135
410000
410000 SIMILAR TYPE TANKS CAPABLE OF STORING CORROSIVE FUELS AND
410000 OXIDIZERS IN THE SORTIE LABORATORY AT LOW PRESSURES. USED IN
410000 STORES AND STORES.
410000
410000 SIZE APPROXIMATELY 0.75 DIA X 3 FT
410000 (C.225 DIA X 0.9 M)
410000
410000 TEMPERATURE SENSORS ECCC1 41121
410000
410000 MEASURE COOLANT TEMPERATURES
410000
410000 OPERATING TEMPERATURE RANGE...-76 TO 104F (-60 TO 40C)
410000 ACCURACY..... +/-0.2F (+/-0.1C)
410000
410000 RESISTANCE THERMOMETER SP003 41111
410000
410000 THERMOCOUPLES SP003 41111
410000
410000 INSTRUMENTATION SENSOR SYSTEM ST223 11114
410000
410000 SYSTEM INCLUDE TEMPERATURE, PRESSURE AND ACCELEROMETER SENSORS
410000
410000 TEMPERATURE SENSORS ST260 11111
410000
410000 SENSOR TO MEASURE TEMPERATURE OF COMPLETION PROCESS
410000
410000 TEMPERATURE RANGE(EST)..... 100C TO 3000C
410000
410000 PLMP EC001 11111
410000
410000 PLMP COOLANT IN CPL COOLANT SUBSYSTEM
410000
410000 FLOW RATE(EST)..... 40LBS/HR (10000G/MIN)
410000 HEAD(EST)..... 1 PSI (6.9003 N/M²)
410000
410000 METERING PUMPS SP001 11111
410000
410000 BUFFER AND SPECIMEN METERING VALVE FOR CONTINUOUS FLOW
410000 ELECTROPHORETIC COLUMN
410000
410000 FLOW RATE..... 0.07ML/MIN TO 30ML/MIN
410000 PRESSURE HEAD..... 4PSI - 5.6 TO 14 PSI (4004 - 10004N/M²)
410000 REPRODUCIBILITY..... +/-0.50
410000
410000 PLMP ST100 11111
410000
410000 PLMP WATER TO STEAM GENERATOR
410000
410000 FLOW RATE..... 0.1 LB/HR (C.045KG/MIN)
410000
410000 PLMP GAS, CIRCULATING LS1C11 11111
410000
410000 PLMP FOR SEALEC PLANT GROWTH CHAMBER
410000
410000 HIGH VACUUM PUMP SP002 21111
410000
410000 UNIT IS USED IN CONCERT WITH MOLECULAR SIEVE. IT IS USED WITH
410000 THE SPACE VACUUM TO BLEED DOWN EVENS. THIS UNIT IS ALSO USED ON
410000 THE LEVITATION SUBSYSTEM.
410000
410000 PRESSURE-ABSOLUTE AMBIENT TO 0.1 PSIA
410000
410000 PYROMETER SP002 21111
410000
410000 THE PYROMETER OPTICALLY SENSES THE INTENSITY OF RADIANT ENERGY
410000 EMITTED FROM THE SAMPLE SURFACE
410000
410000 TEMPERATURE RANGE 104 TO 5400 F (40 TO 3000 C)
410000 SPECTRAL RESPONSE..... 0.7C - 0.97 MICROMETERS
410000 CALIBRATION ACCURACY..... +/-19 FULL SCALE
410000 REPEATABILITY..... +/-19 FULL SCALE
410000
410000 TWO - COLOR PYROMETER SP002 21111
410000
410000 LASER PYROMETER SP003 21239
410000
410000 PYROMETER ST176 11111
410000
410000 MEASURE TEMPERATURE OF COMPLETION PROCESS IN FLAME CHEMISTRY
410000 EXPERIMENTS. USED IN STORES AND STORES
410000
410000 TEMPERATURE RANGE 0 TO 3000 DEC F
410000
410000 DIGITAL TAPE RECORDER SP001 11111
410000
410000 PROVIDE DATA STORAGE FOR DIGITAL PROCESS PROGRAMMER



423000 CHANNELS.....15
423000
423000
424000 RECORDER-TAPE AS011 11111
424000
424000 RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND
424000 DISPLAY CONSOLE.
424000
424000 BANDWIDTH DC TO 200 KHZ
424000 CHANNELS 14
424000 TIME MARKS EDGE TRACK
424000 RECORDED FORMAT DIGITAL
424000
424000 RECORDER-TAPE ME011 11111
424000
424000 RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND
424000 DISPLAY CONSOLE.
424000
424000 BANDWIDTH DC TO 200 KHZ
424000 CHANNELS 14
424000 TIME MARKS EDGE TRACK
424000 RECORDED FORMAT DIGITAL
424000
424000 RECORDER-TAPE SCC11 11111
424000
424000 RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND
424000 DISPLAY CONSOLE.
424000
424000 BANDWIDTH DC TO 200 KHZ
424000 CHANNELS 14
424000 TIME MARKS EDGE TRACK
424000 RECORDED FORMAT DIGITAL
424000
424000 TAPE RECORDER - ANALOG APR14 11111
424000
424000 RECORD AND PLAYBACK OF ANALOG DATA IN SUPPORT OF ATMOSPHERIC AND
424000 SPACE PLASMA PHYSICS SORTIE LABORATORY INVESTIGATIONS.
424000
424000 BANDWIDTH 2 MHZ
424000 CHANNELS 14
424000 RECORDED FORMAT FM
424000
424000 TAPE RECORDER - DIGITAL APB17 21111
424000
424000 RECORDING AND PLAYBACK OF DIGITAL DATA COMPATIBLE WITH ANY CA-
424000 PEARC COMPUTING EQUIPMENT.
424000
424000 BANDWIDTH 2 MHZ
424000 CHANNELS 8
424000 SIGNAL C TO 5 VCC DIGITAL
424000
424000 TAPE RECORDER EC010 11111
424000
424000 RECORD DATA ON MAGNETIC TAPE
424000
424000 BIT RATE.....1000/5
424000 STORAGE CAPACITY.....3.5FC7HITS
424000
424000 DIGITAL TAPE RECORDER ECC06 11111
424000
424000 PROVIDES FOR RECORDING SENSOR DIGITAL OUTPUT DATA.
424000
424000 TAPE SIZE 1 INCH
424000 TRACKS 5
424000 TAPE SPEED VARIABLE
424000 DATA RATE 1.3E+06 BPS
424000 REEL SIZE 12 INCHES CIA
424000
424000 DIGITAL RECORDER EC091 11111
424000
424000 PROVIDES FOR RECORDING DIGITAL OUTPUT DATA.
424000
424000 TAPE SIZE 1 INCH
424000 TRACKS 5
424000 REEL SIZE 12 INCHES
424000 TAPE SPEED VARIABLE
424000 DATA RATE..... 1E+06 BPS
424000
424000 MAGNETIC TAPE RECORDER PEC 21111
424000
424000 ENCARD DATA RECORDING OF SCIENTIFIC INSTRUMENT DATA
424000
424000 READ IN RATE.....1000RPS
424000 STORAGE CAPACITY/MISSION.....3E12 BITS
424000
424000 TAPE INPUT UNIT SFC01 11111
424000
424000 RECEIVES, STORES, TRANSMITS DATA OF OTHER EQUIPMENT
424000
424000 TAPE CAPACITY.....1200FT(36CPI, 0.51A, 11.3CPI WIDE
424000 TAPE SPEED..... 25, 18.75 AND 12.5 IN/S
424000 164.4E AND 32 CM/SI
424000 LONG TERM VARIATION.....+10
424000 DATA DENSITY.....1000 CPI (PHASE ENCODED)
424000 200, 556, 800 CPI (NRZII)
424000
424000 MAGNETIC TAPE RECORDER LS950 11111
424000
424000 RECORD DATA
424000
424000 DATA RATE 200 KPS
424000 TOTAL DATA STORED 1E+06 BITS
424000
424000 TAPE RECORDERS (HIGH-SPEED MULTICHANNEL) ST056 11111
424000
424000 MULTIPLE CHANNEL HIGH SPEED ELECTRONIC DATA RECORDING SYSTEM.

ORIGINAL PAGE IS
OF POOR QUALITY

424000 USED ON ST01S, ST02S AND ST03S.
424000
424000 BANDWIDTH 1 KHZ TO 2.0 PHZ
424000 CHANNELS 14
424000 SIGNAL TYPE ANALOG
424000
424000 TYPE RECORDER, ANALOG ST009 31111
424000
424000 RECORD VOICE, PH AND VIDEO DATA OF LANDMARK AND STAR-TRACKING
424000 PLOTS, USED IN ST02S, ST03S AND ST05S.
424000
424000 BANDWIDTH TO 2.0 PHZ
424000 RECORDING FORMAT VOICE, PH AND VIDEO
424000 CHANNELS 7
424000
424000 RECORDER, DIGITAL CA905 11111
424000
424000 ACCURATE RECORDING OF RAW DIGITAL DATA
424000
424000 SIGNAL TYPE 0 TO 5 VDC DIGITAL
424000 CHANNELS 9
424000 TAPE SPEED VARIABLE
424000 DATA RATE 1.5E+6BPS
424000
424000 TYPE RECORDER, DIGITAL ST008 21111
424000
424000 RECORD AND PLAYBACK OF DIGITAL DATA IN SUPPORT OF MICROWAVE
424000 INTERFEROMETER AND TRACKING INVESTIGATIONS.
424000 USED IN ST01S, ST02S, ST03S, ST04S, ST05S, ST06S, ST07S AND
424000 ST09S.
424000
424000 SIGNAL TYPE 0 TO 5 VDC
424000 BANDWIDTH 100 TO 250 KHZ
424000 CHANNELS 4 (MINIPLM)
424000 DATA REQUIREMENT 2 HRS (MINIPLM)
424000
424000 RECORDER, OPTICAL ST273 11235
424000
424000 RECORD AND PLAYBACK DIGITALLY GENERATED DATA FROM A LASER SYSTEM,
424000 USED ON ST02S AND ST03S.
424000
424000 BANDWIDTH TO 50 PHZ
424000 SIGNAL DIGITAL
424000 CHANNELS 3
424000
424000 RECORDER, VIDEO CA910 11111
424000
424000 PROVIDES FOR RECORDING VIDEO DATA.
424000
424000 CHANNELS 2
424000 BANDWIDTH 4 PHZ MINIPLM
424000 DURATION 1.5 HRS MINIPLM
424000
424000 RECORDER - VOICE (CORE) LS1004 11111
424000
424000 CAPABILITY OF RECORDING ORAL COMMENTS.
424000
424000 CHANNELS 4
424000 BANDWIDTH 50 TO 20000 HZ
424000
424000 RECORDER - STRIP CHART AP012 21111
424000
424000 RECORD WARE COPY DATA FROM LOW FREQUENCY ANALOG INVESTIGATIONS
424000
424000 FREQUENCY RESPONSE DC TO 150 HZ
424000 CHANNELS EVENT, TIME, 2- DATA
424000 ACCURACY +/- 0.5 PCT
424000
424000 RECORDER - STRIP CHART AP013 21111
424000
424000 RECORD WARE COPY DATA.
424000
424000 CHANNELS 8
424000
424000 RECORDER - STRIP CHART (CORE) LS1013 11111
424000
424000 WARE RECORD OF APPROPRIATE LIFE SCIENCES DATA.
424000
424000 FREQUENCY RESPONSE DC TO 150 HZ
424000 SENSITIVITY 0.1 TO 50 MV/PP
424000 CHANNELS 8
424000 CHART SPEED 1 TO 500 PP/SEC
424000 WRITING METHOD THERMAL
424000
424000 VISICORDER CA914 11111
424000
424000 REAL-TIME RAW DATA RECORDING.
424000
424000 BANDWIDTH DC TO 5000 HZ
424000 CHANNELS 18
424000 RECORDING FORMAT CRT WARE COPY
424000
424000 RECORDER - STRIP CHART ST015 11111
424000
424000 RECORD ANALOG MICROWAVE RADICETER DATA MONITORING SEA SURFACE
424000 TEMPERATURE AND ROUGHNESS.
424000 USED IN ST01S.
424000
424000 FREQUENCY RESPONSE 150.0 HZ
424000 SENSITIVITY 10 PP/DIVISION
424000 CHANNELS 2
424000 INPUTS 2
424000
424000 X-Y RECORDER AP011 21111
424000
424000 REAL-TIME ACCURATE REPRODUCTION OF SPECTRAL DATA SOURCES PLECTEC

430000 CA CARTESIAN COORDINATES.
430000 CHANNELS 1
430000 ACCURACY +/- 0.2 % FS
430000 X-Y PLOTTER CANC7 11111
430000 PROVIDES REAL-TIME ACCURATE REPRESENTATION OF SPECTRAL DATA SOURCES
430000 PLOTTED ON CARTESIAN COORDINATES.
430000 CHANNELS 4 PLLG IN'S
430000 ACCELERATION 1500 IN/S/S (3800 CM/S/S)
430000 SENSITIVITY 0.5 MV/IN (0.25 MV/CM)
430000 DISSOLVED OXYGEN ANALYZER SP001 11111
430000 THE UNIT WILL DETERMINE THE AMOUNT OF DISSOLVED OXYGEN PRESENT IN
430000 THE ELECTROLYTE AFTER IT HAS BEEN PROCESSED THROUGH THE GAS
430000 ELIMINATION SYSTEM. IT WILL MEASURE OXYGEN DURING A REACTION
430000 EXPERIMENT. IT WILL ALSO FUNCTION AS A GASEOUS OXYGEN MEASURER
430000 CONCENTRATION RANGE..... 0 TO TOTAL OXYGEN SATURATION
430000 ACCURACY..... +/- 1%
430000 RESPONSE TIME..... 500 ACTUAL IN 10 SECONDS
430000 SAMPLE TEMPERATURE RANGE..... 32 TO 1100 C TO 430 C
430000 RESIDUAL GAS ANALYZER SP002 21111
430000 IDENTIFY TYPE AND QUANTITY OF GAS REMAINING IN A PROCESSING
430000 CHAMBER OR GENERATED DURING A PROCESS. THE ANALYZER IS ALSO USED
430000 IN THE LEVITATION EXPERIMENT. THE UNIT MUST HAVE A TUNABLE PASS
430000 UNIT SECTION FOR LEAK DETECTION. BE A CATHODIC TYPE INSTRUMENT
430000 CAPABLE OF SEPARATING PHASES. THE INSTRUMENT MUST ACT AS A
430000 PARTIAL AND TOTAL PRESSURE ANALYZER AS WELL AS GAS ANALYZER.
430000 MASS RANGE RESOLUTION..... UP TO 200 AMU
430000 DETAILED SPECTRAL ANAL. MASS RANGES..... 1-50 AND 40-65 AMU
430000 SCOPE SCAN SPEED..... VARIABLE 3MS TO 300S
430000 DATA OUTPUT..... VISUAL ANALOG
430000 ANALOG
430000 DIGITAL
430000 USEFUL PRESSURE RANGE..... 1.0E-5 PSI (7.5E-2N/P2)
430000 TO 1.0E-11 PSI (7.5E-11N/P2)
430000 MASS SEPARATION CAPABILITY..... 1 AMU
430000 MIN DETECTABLE PARTIAL PRESSURE..... 2E-12 PSI (1.5E-11N/P2)
430000 MIN DETECTABLE TOTAL PRESSURE..... 1E-10 PSI (7.5E-11N/P2)
430000 GAS ANALYZER, AUTOMATIC LS953 11111
430000 DETERMINE THE PARTIAL PRESSURES OF OXYGEN AND CARBON DIOXIDE
430000 DISSOLVED IN BLOOD SAMPLES AND TO DETERMINE THE HYDROGEN ION
430000 CONCENTRATION
430000 (REQUIREMENTS UNSPECIFIED. COMMERCIAL UNIT USED AS MODEL)
430000 PH..... 0.000 TO 0.000
430000 CO2 PARTIAL PRESSURE..... 0 TO 200 MM HG
430000 O2 PARTIAL PRESSURE..... 0 TO 200 MM HG
430000 ACCURACY:
430000 PH..... +/- 0.003
430000 CO2 PARTIAL PRESSURE..... +/- 0.003
430000 O2 PARTIAL PRESSURE..... +/- 100 MM HG AT 200 FC2
430000 +/- 100 MM HG AT 2000 FC2
430000 GAS ANALYZER, CO2 SPECIFIC LS961 11111
430000 THIS DEVICE IS USED TO MEASURE ATMOSPHERIC CARBON DIOXIDE LEVELS
430000 PARTIAL PRESSURE RANGE..... 13C - 2700 PAH - 20 MM Hg
430000 MICROBIOLOGICAL SAMPLE STORAGE ST002 21112
430000 STORAGE RACK FOR SPACE-TYPE PLASTIC BAG, SNAPS AND CAPSULES TO
430000 DESTROY BACTERIA IN CASE OF BREAKAGE. USED IN ST01S, ST02S AND
430000 ST03S.
430000 SIZE APPROXIMATELY 16X13X6 IN (40X33X15 CM)
430000 MASS SPECTROMETER LS994 11111
430000 PROVIDE ANALYSIS OF UNKNOWN IONIZED GAS BY ATOMIC MASS MEASUREMENT.
430000 MASS RANGE..... 1 TO 400 AMU
430000 SCANNING RANGE..... VARIABLE
430000 VARIABLE SCANNING RATE..... 50 MILLIseconds TO 600 SEC/SCAN
430000 MAX OPERATING PRESSURE..... 1E-4 TORR
430000 SELECTABLE CENTER MASS
430000 SELECTABLE SCANNING WIDTH
430000 MASS SPECTROMETER LS928 11111
430000 PROVIDE AN ANALYSIS OF AN UNKNOWN IONIZED GAS BY ATOMIC MASS
430000 MEASUREMENT
430000 MASS RANGE..... 0 TO 60 AMU
430000 MASS SPECTROMETER ST103 11111
430000 ANALYSIS OF GASEOUS COMBUSTION PRODUCTS FORMED IN ZERO GRAVITY
430000 VIA MASS SPECTROSCOPY. USED IN ST04S AND ST05S.
430000 MASS RANGE 1 TO 1000 AMU
430000 SCANNING RATE 0.5 SEC TO 600 SEC
430000 SENSITIVITY 1 V
430000 ELECTRON SPIN RESONANCE/MASS SPECTROMETER ST250 11112
430000

ORIGINAL PAGE IS
OF POOR QUALITY

459605 IDENTIFY CONSTITUENTS OF COMBUSTION PRODUCTS
459605 IR SPECTROMETER SPG03 11111
459605 TAPPAEC ANALYSIS OF CONSTITUENTS IN GASES AND LIQUID STREAMS.
459605 RANGE 1C TO 1CC MICRONS
459605 RESOLUTION +/- C.5 MICRONS
459605 RESPONSE TIME 2 SECONDS
459605 SPECTROMETER - GAMMA RAY STG03 11111
459605 MEASURE SPACECRAFT GAMMA RADIATION BY CRYOGENICALLY COOLED LITH-
459605 IUM CRYSTAL GERMANIUM DETECTOR WITH ANTI-COINCIDENCE SHIELD. USED
459605 IN STC15 AND ST035.
459605 ELECTROMAGNETIC ENERGY..... 1CC MEV TO 6 MEV
459605 SPECTROMETER - CHARGE PARTICLE STG04 11111
459605 THIRTY MULTI-WIRE PROPORTIONAL COUNTER WITH INTERVENING LAYERS OF
459605 TISSUE EQUIVALENT PLASTIC. USED IN STC15 AND ST035
459605 PROTON ENERGY RANGE 5 TC 23C MEV
459605 DEUTERIUM ENERGY RANGE..... 1C TO 3CC MEV
459605 TRITIUM ENERGY RANGE..... 15 TC 23C MEV
459605 HELIUM-3 NUCLEI ENERGY RANGE . 2C TO 4CC MEV
459605 ALPHA PARTICLE ENERGY RANGE ... 25 TC 45C MEV
459610 SPECTROMETER - NEUTRON STG05 11111
459610 LIQUID SCINTILLATOR WITH ANTI-COINCIDENCE SHIELD PULSE SHAPE
459610 DISCRIMINATOR. USED ON ST015 AND STC35.
459610 NEUTRON ENERGY RANGE 0.5 TC 15 MEV
459610 UV-VIS SPECTROMETER SPOC1 21111
459610 MEASUREMENT OF SAMPLE CONSTITUENTS BASED ON THE UV AND VISIBLE
459610 SPECTRUM.
459610 RANGE 15C TO 1CC00 AN
459610 WAVELENGTH/WAVENUMBER ACCURACY... +/- 0.1 NM FS
459610 RESOLUTION C.05 NM
459610 RESPONSE TIME 2 SECONDS
459610 SCANNING SPEED 0.02 TO 2.0 CM/SEC
459620 SPECTROMETER - VOLUMETRIC STG05 11111
459620 PROVIDE A PHOTOGRAPHIC MEANS TO MEASURE INDEPENDENT VOLUME OF
459620 CELLS AND CELL SHAPES. USED IN STC15, ST025 AND ST035.
459620 LONG WAVE INFRARED SPECTROMETER PEC 11111
459620 USED FOR IDENTIFICATION OF TYPES OF ROCKS, SOILS AND SEDIMENTS
459620 SPECTRAL RANGE.....0.4 TC 2.4 MICRONS
459620 0.2 TC 15.5 MICRONS
459620 RADIOMETRIC CHANNEL RANGE.....1C TO 12 MICRONS
459620 INSTANTANEOUS FIELD OF VIEW.....1 MILLI RADIAN DIA.
459620 POINTING ANGLE FROM RADAR
459620 (PITCH, ALONG TRACK).....+45DEG TC -1CDEG
459620 POINTING ANGLE FROM RADAR
459620 (ROLL, CROSS TRACK).....+/-2CDEG
459620 POINTING ACCURACY.....C.3DEG
460300 SPECTROPHOTOMETER - GENERAL ANALYZER LS931 11111
460300 SPECTRAL ANALYSIS OF GASES AND LIQUIDS INCLUDING SOLIDS OF LIGHT
460300 SOURCES.
460300 FREQUENCY 25C TO 2500 NANOMETERS
460300 (2500 TO 250000 ANGSTROMS)
460300 ACCURACY +/- C.4 NANOMETERS
462500 ANALYZER - ATOMIC ABSORPTION SPECTROPHOTOMETER LS903 11111
462500 QUANTITATIVE DETERMINATION OF METALLIC AND SEMI-METALLIC ELEMENTS
462500 IN SOLUTIONS AND INDICATING CONCENTRATION.
462500 WAVELENGTH RANGE 2CC TO 8CC MILLIMICRONS
462500 MODES OF OPERATION ATOMIC ABSORPTION, FLAME
462500 EMISSION, UV-VISIBLE SPECTRO-
462500 PHOTOMETRY
462500 COUNTER ACCURACY +/- C.4 NANOMETERS
464400 SPECTROPHOTOMETER - IR (ICPE1) LSC22 11111
464400 PERFORM SPECTROPHOTOMETER ANALYSIS OF LIFE SCIENCES SPECIMEN IN
464400 THE INFRARED RANGE .
464400 SPECTRAL RANGE 2.5 TC 12 MICROMETERS
464400 RESOLUTION 0.1 PERCENT OF WAVELENGTH
464400 SCAN TIME 9 SEC TO 90 MCLR
464400 SAMPLE TEMPERATURE AMBIENT TO 402 F (250 C)
470500 STERILIZER - TOOL LS903 11111
470500 STERILIZE MISCELLANEOUS SMALL METAL HARD TOOLS SUCH AS SCALPERS,
470500 BY MEANS OF ELECTRICAL INDUCTION HEATING.
470500 CAPACITY 0.1 CU FT
480000 VIDECON ECG12 11131
480000 RECORD CHAMBER PHENOMENA

480000
480000 RESOLUTION.....HIGH
480000 LIGHT LEVEL.....LOW
480000 IMAGE MAGNIFICATION.....YES
480000 IMAGE INTENSIFICATION.....YES
480000
480000 CCTV CAMERA SFC01 11131
480000
480000 HIGH RESOLUTION TELEVISION CAMERA FOR VIEWING, RECORDING OF
480000 TRANSMITTING PICTURES OF MATERIAL PROCESSING EXPERIMENTS
480000
480000 BANDWIDTH.....32MHZ
480000 VERTICAL SCAN LINES.....1125 LINES
480000 HORIZONTAL RESOLUTION.....1100 LINES
480000 VERTICAL SWEEP RATE.....25 OR 30 FRAMES PER SECOND
480000 SENSITIVITY.....100 THE UNITS OUTPUT FREQUENCY
480000 WITH 0.5 FOOTCANDLE ILLUMINATION
480000
480000 VIDEO CAMERA, B/W LS939 401131
480000
480000 PROVIDE MEANS OF ACTIVITY MONITORING, EXPERIMENT DATA ACQUISITION
480000 ETC. SYSTEM WILL INTERFACE WITH A 4C INPUT VIDEO MULTIPLIER TO
480000 PERMIT THE MONITORING OF 4C SIMILAR VIDEO CAMERAS.
480000
480000 VISUAL RESPONSE.....APPROXIMATE HUMAN EYE
480000 VIDEO OUTPUT.....CONSTANT WITH LIGHT LEVEL
480000 CHANGES OF 10 TO 1000FT CANDLES
480000 1.4V P-P COMPOSITE, COMPOSITE TO
480000 EIA RS-170 STANDARD
480000
480000 CAMERA, TELEVISION C9009 11131
480000
480000 PROVIDES ASTRONAUT WITH VIEW OF EARTH SCENE.
480000
480000 VIEW ANGLE20 DEGREES
480000 RESOLUTION625 LINES
480000 SPECTRAL BANDWIDTH4 MHZ
480000 SWEEP RATE15 FRAMES PER SECOND
480000
480000 CAMERA, TV ST025 11131
480000
480000 REAL-TIME VIEWING OF TARGET AREAS, LAND MARK AND STAR FIELD
480000 TRACKING, AND SENSOR BORESIGHTING.
480000 USED ON ST045, ST055 AND ST065.
480000
480000 VIEW ANGLE0 TO 10 DEGREES
480000 RESOLUTION1025 LINES
480000
480000 VIDEO CAMERA, COLOR LS940 11111
480000
480000 PROVIDE COLOR VIDEO OBSERVATIONS OF EXPERIMENTAL ANIMALS, SETUPS,
480000 ETC. FOR STORAGE VIA VIDEO-TAPE FOR TRANSMISSION TO THE GROUND.
480000
480000 RESOLUTION (MIN).....250 LINES
480000 LENS.....STANDARD 16MM, C MOUNT
480000
480000 PCCP POSITION TV DISPLAY AF001 11111
480000
480000 DISPLAY MEDIUM RESOLUTION TELEVISION OF EXPERIMENT BECP POSITIONS
480000 BEFORE, DURING AND AFTER EACH EXPERIMENT REPETITION.
480000
480000 RESOLUTION525 LINES
480000 FRAME RATE30 FRAMES/SEC
480000 BANDWIDTHVHF RANGE
480000
480000 EXPERIMENT TV DISPLAY AF002 11111
480000
480000 DISPLAY HIGH RESOLUTION TELEVISION PICTURES OF AREAS OF INTEREST
480000 IN REAL TIME OR FROM VIDEO RECORDERS AT VARIOUS SPEEDS.
480000
480000 RESOLUTION.....525 LINES
480000 FRAME RATE.....30 FPS
480000 BANDWIDTH.....VHF RANGE
480000
480000 TV MONITOR EC012 11111
480000
480000 DISPLAY TV CAMERA OUTPUT SHOWING CHAMBER PHENOMENA
480000
480000 TV MONITOR SP001 11111
480000
480000 THE MONITOR WILL DISPLAY CCTV IMAGES AS PART OF THE ELECTRO-
480000 OPTICAL IMAGING SYSTEM.
480000
480000 SCANNING FREQUENCIES.....HORIZONTAL: 15-40 KHZ
480000VERTICAL: 15-60 FIELDS/SECOND
480000 VERTICAL SCAN LINES.....1125 LINES
480000 LINEARITY.....19 PICTURE HEIGHT
480000
480000 MONITOR - TV ST026 11111
480000
480000 DISPLAY HIGH RESOLUTION TELEVISION PICTURES IN REAL-TIME OF FROM
480000 VIDEO TAPE. USED ON ST015, ST025, ST035, ST045 AND ST055.
480000
480000 RESOLUTION625 LINES
480000 FRAME RATE30 FPS
480000 SCREEN SIZE10 IN (25.4 CM)
480000
480000 MONITOR, VIDEO LS1009 11111
480000
480000 MONITOR TV PICTURES OF ANIMAL AND OTHER LAB ACTIVITIES
480000
480000 CONSTANT TEMPERATURE BLOCK (CCP) LS1015 11111
480000
480000 HEAT SINK MAINTAINING A CONSTANT TEMPERATURE FOR TEST TUBES AND
480000 VIALS.
480000

ORIGINAL PAGE IS
OF POOR QUALITY

486500 TEMPERATURE 240+/-1.0 C F (120+/-0.5 C)
486500
487000 TEMPERATURE CONTROLLER 1C001 11111
487000
487000 CATHODES COOLANT TEMPERATURE TO CATHODE CHAMBER WALL TEMPS
487000
487000 TEMPERATURE RANGE.....-76 TO 104F (-60 TO 40C)
487000 ACCURACY.....+/-0.2F (0.1C)
487000
487000 TEMPERATURE MONITOR ST194 11114
487000
487000 PENITORS AND CONTROLS TEST CHAMBER TEMPERATURE
487000
488200 TIMER - EVENT AS006 11111
488200
488200 DIGITAL DISPLAY OF TIME REMAINING OR EXPENDED FOR A PARTICULAR
488200 EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE
488200 START AND STOP COMMANDS (MINUTES AND SECONDS).
488200
488200 DISPLAY 4 DIGITS
488200 DISPLAY TYPE LED
488200
488200 TIMER - EVENT ME006 11111
488200
488200 DIGITAL DISPLAY OF TIME REMAINING OR EXPENDED FOR A PARTICULAR
488200 EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE
488200 START AND STOP COMMANDS (MINUTES AND SECONDS).
488200
488200 DISPLAY 4 DIGITS
488200 DISPLAY TYPE LED
488200 INTERFACE DISCRETE
488200
488200 TIMER - EVENT SE006 11111
488200
488200 DIGITAL DISPLAY OF TIME REMAINING OR EXPENDED FOR A PARTICULAR
488200 EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE
488200 START AND STOP COMMANDS (MINUTES AND SECONDS).
488200
488200 DISPLAY 4 DIGITS
488200 DISPLAY TYPE LED
488200
488200 TIMER, EVENT LS1020 11111
488200
488200 GENERAL PURPOSE ELAPSED TIME DEVICE USED FOR VEHICLE PSYCHOMETER
488200 TESTS
488200
488200
488300 DIGITAL CLOCK SF001 11111
488300
488300 DIGITAL TIME DISPLAY, PART OF DIGITAL PROCESS PROGRAMMER ASSY
488300
488300 DISPLAY 4 DIGITS
488300 DISPLAY TYPE LED
488300
488300 TIMER, PRECISION CLOCK CA902 11111
488300
488300 PROVIDES PRECISE TIME REFERENCE, MEASURE TIME AT WHICH PHASE PEAK
488300 OCCURS.
488300
488300 DISPLAY 24 BIT DIGITAL CLOCK WITH 1 HR.
488300 TIME PERIOD
488300 DISPLAY TYPE DIGITAL
488300 ACCURACY 1E-06
488300
488500 TIMER - MISSION AS007 11111
488500
488500 PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECCAD UP-
488500 DATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE,
488500 SECCAD.
488500
488500 DISPLAY 7 DIGITS BCD
488500 DISPLAY TYPE LED
488500
488500 TIMER - MISSION ME007 11111
488500
488500 PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECCAD UP-
488500 DATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE,
488500 SECCAD.
488500
488500 DISPLAY 7 DIGITS BCD
488500 DISPLAY TYPE LED
488500
488500 TIMER - MISSION SE007 11111
488500
488500 PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECCAD UP-
488500 DATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE,
488500 SECCAD.
488500
488500 DISPLAY 7 DIGITS BCD
488500 DISPLAY TYPE LED
488500
488500 TIME CODE GENERATOR AND DISPLAY AP032 11111
488500
488500 ACCURATE FREQUENCY, TIME INTERVAL AND TIME KEEPING CAPABILITIES.
488500
488500 OUTPUTS..... 100 MHZ, 1 PHZ, 5 PHZ
488500 ACCURACY 5E-10 PARTS PER DAY
488500
488500 TIMER ST107 11111
488500
488500 TIMING AND VARIABLE PROGRAMMING CAPABILITY TO OPERATE AND CONTROL
488500 EXPERIMENT AND SUPPORTING EQUIPMENT. USED IN ST01S AND ST03S.
488500
488500 STAINING SYSTEM - BACTERIOLOGICAL LS924 11111
488500
488500 STAIN AND FIT BIOLOGICAL SPECIMENS FOR MICROSCOPIC EXAMINATION.
488500
488500

5079CC MICRODISSECTION TOOL KIT LS958 11111
5079CU
5079CO STANDARD MEDICAL MICRODISSECTION KIT
5079CO
5079CC VETERINARY MEDICAL KIT (CCRF) LS957 11111
507900
5079CC INCLUDES OTOSCOPE/OPHTHALMOSCOPE, REFLEX HAMMER, HEMOCLETHROMETER,
507900 THERMIST, SYRINGES, NEEDLES, SCALPERS, HEMOSTATS, THERMISTERS, AND
507900 OTHER STANDARD ITEMS.
5079CO
507900 HISTOLOGY KIT/SLIDE CABINET (CCRF) LS956 11111
5079CO
5079CC STANDARD KIT CONTAINING PROBES, KNIVES, NEEDLES, SYRINGES, CERV
507900 SLIPS, LACLES PLUS OTHER NORMALLY INCLUDED ITEMS.
5079CO
5079CC KIT, BENCH CHEM ANAL LS965 11111
507900
5079CC KIT CONTAINS THE TOOLS AND EQUIPMENT TO MANAGE CHEMICALS AND
507900 BIOLOGICALS DURING VARIOUS MANUAL PROCEDURES GENERALLY PERFORMED
507900 WITHIN THE GLOVE BOX.
5079CO
5079CC TOOLS..... SOLID TRANSPORT TOOLS, CAVITY
507900 INDEPENDENT PIPETTES, VIALS,
507900 PIPETTES, AND TEST TUBES.
507900 CHEMICALS, STROBERS, FILTERS,
507900 AND SAFETY SHIELDS
5079CO
5079CC KIT, HEMATOLOGY LS966 11111
5079CO
5079CC PROVIDE TOOLS FOR SAMPLING, HANDLING, TRANSFERRING AND ANALYZING
507900 PLCC
5079CO
5079CC ITEMS..... HEMATOCYTER KIT: 10 LAPED
507900 DISC PIPETTES; CONVEYERS;
507900 SLIDES; WBC AND RBC TILERS;
507900 CRITOCAL; HEMOCYT TUBES (MICRO,
507900 HEPARINIZED); HEMOCYT TUBES
507900 (MICRO, PLAIN); PLCC CILLING
507900 PIPETTES (WBC, RBC); CEMISTIX
507900 (URINE TEST STRIPS); ALER
507900 ADAPTERS, VACUTAINER; VACUTAINER
507900 VACUTAINER-NEEDLE UNIT;
507900 VACUTAINER TUBES, ASSORTED (2ML);
507900 SYRINGE (1/4CC), PEDIATRIC
507900 ALCOHOL SWABS; LANCETS; NEEDLE,
507900 25GA, 5/8IN; NEEDLE 22GA, 1 1/2IN
5079CO
5079CC KIT, LINEAR MEASUREMENT LS967 11111
5079CO
5079CC DETERMINE SIZE, AMPLITUDE, DISTANCE, CIRCUMFERENCE, ETC
5079CO
5079CC ITEMS..... RULERS; TAPES; LINEAR COMPARATOR;
507900 GRIDS; CALIPERS, INSIDE AND
507900 OUTSIDE; VERNIER CALIPERS;
507900 MICROMETERS, INSIDE AND OUTSIDE
5079CO
5079CC KIT, MICROBIOLOGY LS968 11111
5079CO
5079CC PROVIDE TOOLS TO FACILITATE GROWING AND ANALYZING MICROBIAL
507900 ORGANIZING
5079CO
5079CC ITEMS..... INOCULATING LOOP; INOCULATING
507900 NEEDLE; COTTON SWABS, PACHES,
507900 STERILE SYRINGE SPL, STERILE;
507900 SLIDES; MICRO; ALCOHOL; TUBES,
507900 15X75MM, STERILE, CAPPEC;
507900 ZEPHYRUS, TINCTURE, 1/2CC;
507900 BACTERIOTRANSFER (STERILIZING UNIT)
507900 1 THICGLYCOLATE, TUBES; STUART
507900 TRANSPORT MEDIA, VIALS; TSA
507900 SLANTS
5079CO
5079CC KIT, ORGANISM HOLDING AND MANAGEMENT LS969 11111
5079CO
5079CC PROVIDE TOOLS AND DEVICES USED IN THE HOLDING AND HANDLING OF THE
507900 ORGANISMS
5079CO
5079CC ITEMS FOR SMALL VERTEBRATES..... GLUE POWDER; CAGE SHIELD;
507900 PLASTIC; PLASTIC LINER WITH FAC;
507900 PLCC PELLET DISPENSER; ALCOHOL
507900 SWABS; TUBES; PAPER, DISC;
507900 PLASTIC PACHES (EXPANSIBLES);
507900 PLASTIC HATFOR ACNEFFRACILES;
507900 ORGANISM TRANSFER CAPSULE;
507900 ANIMAL TAGS
5079CO
5079CC ITEMS FOR PLANTS..... WATERING DEVICE (SPRAY BOTTLE);
507900 LABELS; FERTILIZER PACHETS;
507900 STARKS (LOCH)
5079CO
5079CC KIT, PLANT TOOLS LS970 11111
5079CO
5079CC PROVIDE TOOLS FOR VARIOUS PLANT MANIPULATIONS
5079CO
5079CC ITEMS..... SCISSORS, THERZERS, SPATULA,
507900 APPLICATOR, SPLINTS, TAPE
507900 DISPENSER, HYPODERMIC NEEDLES,
507900 SYRINGES, FORCEPS, PLIERS,
507900 SCREWDRIVER, SCALPERS, SEALERS
5079CO
5079CC KIT, GENERAL TOOL LS971 11111
5079CO
5079CC MECHANICAL AND ELECTRICAL TOOLS AND HARDWARE TO PROVIDE
507900 CONVENTIONAL DIAGNOSTIC, MAINTENANCE, AND SERVICE FUNCTIONS.
5079CO
5079CC TOOLS..... HAMMER, WRENCHES, PLIERS, SCREW
507900 DRIVERS, DRILL, ACMESTIVE TAPE.

ORIGINAL PAGE IS
OF POOR QUALITY



5079CC WIRE CUTTERS, WIRE TIES, WIPE,
5079CC LUBRICANTS, FLASHLIGHT, SCISSORS
5079CC + PLTIMETER, FASTENERS, CLAMPS,
5079CC LAMP
5079CC KIT, CLEAN UP LS959 11111
5079CC GENERAL PURPOSE CLEAN UP - SPONGES AND WIPES
5079CC KIT, MEDICAL SURGICAL LS1000 11111
5079CC FOR VARIOUS MINOR SURGICAL PROCEDURES
5079CC KIT, PHYSIOLOGY LS1001 11111
5079CC SPONGES, SPONGE SQUEEZER, VIALS, CALCRIC STIMULATOR FOR EAR
5079CC CANAL, SYRINGES, THERMISTERS, ETC.
5079CC KIT, TOOL - INSECT MANIPULATION LS1002 11111
5079CC TOOL KIT FOR COUNTING, SORTING, EXAMINATION, ETC.
5079CC URINE ANALYZER - AUTOMATIC (CORE) LS984 11111
5217CC CAPABILITY TO AUTOMATICALLY PERFORM URINE ANALYSES.
5217CC VACUUM/PRESSURE MEASUREMENT UNIT SPOC2 31111
5249CC MEASURES PRESSURES IN PROCESS CHAMBERS. UNIT IS ALSO SECUTHEC
5249CC FOR GENERAL PURPOSE AND LEVITATION SUBELEMENTS.
5249CC PRESSURE-ABSOLUTE ADHENT TO 0.01 PSIA
5249CC PRESSURE SENSORS ST200 11111
5249CC SENSORS TO MEASURE PRESSURE OF COMBUSTION EFFICIENT
5249CC PRESSURE RANGE..... TOC FCO04 11124
5262CC FLOW CONTROL ASSEMBLY
5262CC THIS ASSEMBLY CONTROLS THE FLOW OF GASES IN THE CPL. IT IS MADE
5262CC UP OF PLUMBING COMPONENTS, SCIENTIC VALVES, A FILL ASSEMBLY, A
5262CC VENT ASSEMBLY, SAFETY COMPONENTS AND A MIXING SUBASSEMBLY.
5262CC PRESSURE MONITOR ST192 11114
5262CC MEASURES PRESSURE IN TEST CHAMBER AND DISPLAYS DIGITALLY
5262CC SELENOID VALVES FCO01 101121
5270CC CONTROL FLUID FLOW
5270CC REGULATING VALVE ECO01 21121
5270CC ADJUST FLUID FLOW
5270CC PRESSURE REGULATORS ECO02 51121
5270CC PRESSURE REGULATING VALVE FOR GASEOUS PRESSURIZATION SUBSYSTEM
5270CC REGULATION RANGE.....2.7 TO 14.7 PSIA(140 TO 740TCFPA)
5270CC UPSTREAM PRESSURE(EST).....20 PSI (1.46CS A/P21)
5270CC PRESSURE REGULATOR ECO02 11121
5270CC VALVE REGULATING SUPPLY PRESSURE
5270CC REGULATION PRESSURE(EST).....20 PSI (1.46CS A/P21)
5270CC MAX UPSTREAM PRESSURE.....300 PSI (2.1 EOB A/P21)
5270CC SELENOID VALVES FCO02 111121
5270CC GAS SHUTOFF VALVES
5270CC VACUUM/PRESSURE REGULATOR SFO02 31111
5270CC THE UNIT REGULATES THE PRESSURE IN THE PROCESSING CHAMBERS BY
5270CC ADMITTING FLUID FROM A FLUID SUPPLY CONTAINER. THIS UNIT IS ALSO
5270CC EMPLOYED ON THE GENERAL PURPOSE AND LEVITATION SUBELEMENTS
5270CC PRESSURE-ABSOLUTE 0.1 TO 1000 PSIA
5307CC POWER DIVIDER, WIDENAND CANC6 41111
5307CC DIVIDES ANTENNA POWER OUTPUT.
5307CC INPUT PORTS ONE SC CMP PORT
5307CC OUTPUT PORTS THREE SC CMP PORTS
5307CC DIGITAL VOLTMETER SFO01 11111
5320CC VOLTMETER FOR DATA ACQUISITION AND CONTROL UNIT
5320CC DC VOLT RANGE..... 100V TO 100V
5320CC AC VOLT RANGE..... 100V TO 10V
5320CC RESOLUTION..... 0.1% FULL SCALE
5320CC ACCURACY..... $\pm 0.5\%$ 1 DIGIT
5320CC VOLTCOM METER LS1012 11111
5320CC PORTABLE RELIABLE MULTI-FUNCTION PETER FOR GENERAL PURPOSE EXPER-
5320CC IMENTAL WORK AND TROUBLE SHOOTING.
5320CC DC VOLTS 100 PV TO 10 MV
5320CC DC CURRENT 1 PA TO 10 APPS

532000 AC VOLTS 10 MV TO 300 VOLTS
532000 AC CURRENT 1 MA TO 1 AMP
532000 OHMMETER 10 OHMS TO 10 MHOHS
532000 AC/DC VOLTMMETER 100V 21111
532000 MEASURE THE VOLTAGE LEVEL OF VARIOUS SIGNAL SOURCES.
532000 VOLTAGE RANGE +/- 100 MV: 1, 10, 100, 1000 V
532000 FREQUENCY RANGE 45 HZ TO 1 MHZ
532000 ACCURACY +/- 0.01 % OVERALL
601002 ANNUNCIATORS-ADVISORY 21132
601002 CLAL BANKS MOUNTED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL
601002 ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY CREW CARD
601002 EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.
601002 ANNUNCIATOR-VISUAL 12 IN DUAL BANKS
601002 DISPLAY TYPE RED, WHITE, GREEN
601002 ANNUNCIATOR-AUDIO MCRN
601002 ANNUNCIATORS-ADVISORY 21132
601002 CLAL BANKS MOUNTED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL
601002 ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY CREW CARD
601002 EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.
601002 ANNUNCIATOR-VISUAL 12 IN DUAL BANKS
601002 DISPLAY TYPE RED, WHITE, GREEN
601002 ANNUNCIATOR-AUDIO MCRN
601002 ANNUNCIATORS-ADVISORY 21132
601002 CLAL BANKS MOUNTED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL
601002 ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY CREW CARD
601002 EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.
601002 ANNUNCIATOR-VISUAL 12 IN DUAL BANKS
601002 DISPLAY TYPE RED, WHITE, GREEN
601002 ANNUNCIATOR-AUDIO MCRN
603001 ANALOG (SCPI) CONTROLLER 11111
603001 SILICON CONTROL RECTIFIER EMPLOYED AS A RELAY SWITCH
603001 LOGIC LEVELS, LOGICAL "C" C - 0.6V
603001 LOGICAL "I" I - 5.0V
603001 RANGE OF CONTROL C - 100% IN 1 CYCLE INCREMENTS
603002 BATTERY CONTROL AND MONITOR 11114
603002 REMOTE SENSING PLATFORM CIRCULAR CONTROL 11134
603002 XUV NORMAL INCIDENCE SPECTROMETER CONTROL 11114
603002 UV-VIS-NIR SCANNING SPECTROMETER CONTROL 11114
603002 HIGH RESOLUTION FOURIER TRANSFORM SPECTROMETER CONTROL 11114
603002 CIRCULAR IR FOURIER SPECTROMETER CONTROL 11114
603002 IS RADIOMETER CONTROL 11114
603002 FARRE-PEROT INTERFEROMETER CONTROL 11114
603002 UV-VIS DOCUMENTATION CAMERA CONTROL 11114
603002 ELECTROSTATIC ANALYZER CONTROL 11134
603002 MAGNETIC ANALYZER CONTROL 11134
603002 KEV-MEV PARTICLE DETECTOR CONTROL 11134
603002 TOTAL ENERGY DETECTOR CONTROL 11134
603002 MAIN ROOM A CONTROL 11134
603002 PLATFORM ROOM A CONTROL 11134
603002 MAIN ROOM B CONTROL 11134
603002 TIME-ALLEG PLATFORM CONTROLS 11134
603002 SP ROOM CONTROL 11134
603002 ROOM A POWER SUPPLY & DATA SYS CONTROL 11134
603002 ROOM B TARGET CONTROL 11134
603002 ACCEL-DECEL CONTROL 11134
603002 DISCHARGE FILAMENT HEATER CONTROL 11134
603002 DISCHARGE POTENTIAL CONTROL 11134
603002 FLISE SEQUENCE & BURST LENGTH CONTROL 11134
603002 GAS SELECTION & PRESSURE CONTROL 11134
603002 NEUTRALIZER EMISSION & BIAS CONTROL 11134
603002 CHARGE EXCHANGE CHANNEL ACTUATOR CONTROL 11134
603002 BEAM CURRENT MONITOR 11134

ORIGINAL PAGE IS
OF POOR QUALITY



601002	ELECTRON BEAM VOLTAGE COMPART HEATER CONTROL	AF933	11134
601002	ELECTRON BEAM PULSE LENGTH & MAGNITUDE CONTROL	AF934	11134
601002	ELECTRON BEAM EXPANSION LENS CONTROL	AP935	11134
601002	ELECTRON BEAM CURRENT MONITOR	AF936	11134
601002	PHASE, ANGLE, CURRENT MONITOR	AF937	11134
601002	MFC ARC VOLTAGE LEVEL CONTROL	AF938	11134
601002	MFC ARC BURNST CURRENT DURATION CONTROL	AF939	11134
601002	MFC ARC PULSE SEQUENCE	AF940	11134
601002	MFC ARC BEAM CURRENT MONITOR	AP941	11134
601002	C/S SELECTION & PRESSURE CONTROL	AF942	11134
601002	SPHERICAL ION PROBE CONTROL	AF943	11134
601002	CYLINDRICAL ION PROBE CONTROL	AF944	11134
601002	PLANAR ELECTRON PROBE CONTROL	AF945	11134
601002	SEGMENTED PLANAR PROBE CONTROL	AF950	11134
601002	DEP. SAT. INSTRUMENT CONTROL & HOUSEKEEPING	AF959	11134
601002	LEP. SAT EJECTION MECHANISM CONTROL	AF960	11134
601002	PHOTOMETER HIGH VOLTAGE SUPPLY CONTROL	AF961	11134
601002	PHOTOMETER AMPLIFIERS CONTROL	AF962	11134
601002	TV SYS CONTROL - IMAGE INTENSIFIER CONTROL	AF963	11134
601002	CANNISTER EJECTION CONTROL	AF964	11134
601002	PROJECTILE CAMERA CONTROL	AF965	11134
601002	CANNISTER MONITOR	AF966	11134
601002	SHAPED CHARGE EJECTION CONTROL	AP967	11134
601002	SHAPED CHARGE MONITOR	AF968	11134
601002	MILICON EJECTION CONTROL	AP969	11134
601002	C/S CONTROL SYSTEM	AF970	11134
601002	MULTIFREQUENCY PROPAGATION RECEIVER/PROCESSOR	CP004	11134
601002	& ANTENNA SYSTEM CONTROL UNIT		
601002	IR TEMPERATURE PROFILE RADIOMETER CONTROL UNIT	CF018	11134
601002	CAMERA CONTROL	CF020	11134
601002	CINEMA CONTROL	CF020	11134
601002	ALIGNMENT TV CONTROL	AF920	11134
601002	CCTV CAMERA CONTROL UNIT	SF001	11131
601002	THIS UNIT WILL PERFORM THE FOLLOWING CONTROL FUNCTIONS: LENS		
601002	FOCUS, LENS F STOP, CONTRAST BRIGHTNESS AND MODE. THE UNIT WILL		
601002	WORK IN TWO MODES: (1) CONTINUOUS OBSERVATION MODE, WHERE ILLUMINATION		
601002	TICK IS CONTINUOUS, OR (2) PULSED OBSERVATION MODE USED WHEN		
601002	ILLUMINATION IS PROVIDED BY A DYE LASER OR A NARROW FLASH LAMP.		
601002	SCANNING.....INTERLACED 2:1, 1/40 SEC 1125		
601002	LINES PER FRAME, 25 OR 30 FRAMES		
601002	PER SECOND		
601002	HORIZONTAL SCAN FREQUENCY CONTROL 30 TO 750 HZ		
601002	VERTICAL SCAN FREQUENCY CONTROL 50 TO 60 HZ		
601002	LASER - POSITIONING SERVO CONTROL	EC013	11114
601002	PROVIDES FEEDBACK AND CONTROL SIGNAL TO LASER DRIVER FOR PROPER		
601002	POSITIONING AND FOCUSING OF PARTICULATE SPECIMENS		
601002	OPTICAL MONITOR CONTROL	CF021	11134
601002	CONTAMINATION MONITOR GUAGE CONTROL	CF021	11134
601002	MASS SPECTROMETER CONTROL	DF021	11134
601002	SET POINT CONTROLLER	SF001	11111
601002	ESTABLISHES CONTROLLER SET POINT. SET POINT CAN BE VARIED		
601002	FROM PRE-PROGRAMMED INPUT		
601002	SIGNAL INPUT.....C - 50V		
601002	DATA OUTPUT.....1 - 50MA, 0 - 5 VDC		
601002	RFCH CONTROL	ST012	11114
601002	CONTROLS FOR MICROWAVE INTERFEROMETER ROOMS		
601002	OPTICAL ANTENNA SERVO ELECTRONICS	CP066	11114
601002	SUPPORT ELECTRONICS FOR OPTICAL ANTENNA		
601002	LASER TRANSMITTER ELECTRONICS	CP085	11235
601002	DATA FORMATTING FOR TRANSMISSION VIA LASER CARRIER		



603002 SIGNAL CONDITIONING ELECTRONICS ECG10 11111
603003 ELECTRONICS SUPPORTING DATA MANAGEMENT SUBSYSTEM
603003 SIGNAL CONDITIONER SF001 11111
603003 APPLIES AND CONDITIONS LOW LEVEL ANALOG SIGNALS AS PART OF
603003 DATA ACQUISITION UNIT
603003 CAPACITY.....UP TO 33 PAIRS OF SIGNAL WIRES
603003 OUTPUT VOLTAGE.....+1CV
603003 OUTPUT CURRENT.....+5PA
603003 SIGNAL CONDITIONER LS922 11111
603003 ELECTRICAL SIGNAL TRANSFORMATION BETWEEN ANY ALPHER OF TRANS-
603003 CECERS AND MULTIPLEXER/ANALOG TO DIGITAL CONVERTER
603003 COMPUTER INPUT +/- 1C VOLTS
603003 SWITCH/CEPLEX/PREAMP UNIT CAG77 11134
603003 FOR SWITCHING ANTENNAS FOR DIFFERENT ALTITUDES, USE EFF FFI
603003 CANCELLATION, FORMING CPNI PATTERN ETC
603003 CONTROL UNIT - ANTENNA SCAN CAC82 11114
603004 PROVIDES VARIABLE ANGULAR SWEEP PATTERNS (SPIRAL RECTANGULAR
603004 SCANS, ETC.) OF VARIABLE ANGULAR COVERAGE.
603004 CIRCUIT BREAKER/DISTRIBUTOR PANEL ASG09 11131
603004 CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO
603004 THOSE STANDARD ITEMS FOUND IN ALL ASTRONOMY EXPERIMENTS.
603004 CIRCUIT BREAKERS 1E
603004 POWER CAPABILITY 1C APPS
603004 CIRCUIT BREAKER/DISTRIBUTOR PANEL FECC9 11131
603004 CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO
603004 THOSE STANDARD ITEMS FOUND IN ALL ASTRONOMY EXPERIMENTS.
603004 CIRCUIT BREAKERS 1E
603004 POWER CAPABILITY 1C APPS
603004 CIRCUIT BREAKER/DISTRIBUTOR PANEL SC009 11131
603004 CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO
603004 THOSE STANDARD ITEMS FOUND IN ALL ASTRONOMY EXPERIMENTS.
603004 CIRCUIT BREAKERS 1E EACH
603004 POWER CAPABILITY 1C APPS
603009 CALORIMETRY MODULE - PRIMATE LS964 11135
603009 PORTABLE AND COLLAPSABLE UNIT WHICH HOLDS, FEEDS, AND ANALYZES
603009 THE PRIMATE SPECIMEN WITHIN.
603009 P02 1EC +/-1C MM HG
603009 P02 5EC +/-2C MM HG
603009 P02 6 MM HG OR LESS
603009 P020 12 +/-3 MM HG
603009 TEMPERATURE 77 +/-4 F (25 +/-2 C)
603012 CONDENSER - WATER EC002 11131
603012 COLLECTS WATER VAPOR BY CONDENSING WATER CONTENT OF GAS SAMPLE
603012 CONDENSATION RATE(MPH).....C.2LB/HR (51CC/HR)
603012 TEMPERATURE(MIN).....45F(7C)
603012 CONDENSEP ST102 11131
603012 DEVICE TO CONDENSE STEAM IN ZERO GRAVITY
603012 CAMERA CONTROLLER LS937 11114
603014 DEVICE TO CONTROL THE OPERATION OF VIDEO CAMERAS THROUGHOUT THE
603014 LABORATORIES
603015 VACUUM CLEANER LS943 11111
603015 COLLECT GENERAL LABORATORY WASTE MATERIALS, DIRT, CLST AND DEBRIS
603015 VACUUM BAG VOLUME.....C.34CL/FT. (1C.2LITERS)
603015 CLIMOSTAT LS944 11111
603016 MOTOR DRIVEN MOUNTING STRUCTURE TO PROVIDE CONTINUOUS ROTATION
603016 OF PLANT SAMPLES
603016 ROTATION RATES.....TBD
603016 PLATFORM LIGHTING.....TBD
605001 ELECTROMAGNETIC POSITIONING COILS AND DETECTOR SF004 11115
605001 PROVIDE ALTERNATING ELECTROMAGNETIC FIELD FOR CONTACTLESS
605001 POSITION CONTROL.
605001 MAGNETIC PRESSURE 1.4 PSI (1C MM/SC P)
605001 ELECTROSTATIC POSITIONING PROBE AND DETECTOR SF004 11115
605002 ELECTROSTATIC FIELDS USED FOR POSITIONING EXPERIMENT SPECIMEN

ORIGINAL PAGE IS
OF POOR QUALITY

605002 SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL.
 605002 CHAMBER PRESSURE 2 G.14 PSI (1 MP/SC P)
 605002 2 F-06 PSI (0.01 MP/SC P)
 605002 CONFIGURATION IN APEX OF TETRAPOD
 605002 ELECTRIC BEAM SOURCE SPOC4 11114
 605002 THE ELECTRIC BEAM SOURCE WILL BE USED TO HEAT SAMPLES THAT
 605002 REQUIRE NON-CONTACT HEATING. THE UNIT WILL BE USED IN THE CHEST-
 605002 GENERAL PURPOSE ENCLOSURE
 605002 DELIVERABLE POWER..... 3 - 4 KW
 605002 INPUT POWER..... 4500W
 605002 TEMPERATURE..... 4500C (12500F)
 605002 MAXIMUM SAMPLE SIZE..... 0.15 CM DIA. (1200)
 605002 OPERATING PRESSURE..... 1F-1C PSI (17.5E-7A/P21)
 605002 ELECTRIC FIELD CONTROLLER ECC05 11114
 605002 PROVIDE CONTROL FUNCTION FOR ELECTRIC FIELD GENERATOR. CONTROLLER
 605002 SUPPORT BOTH AC AND DC POWER SUPPLIES
 605002 ELECTROMETER LSC05 11111
 605002 MEASURE SMALL VOLTAGES AND CURRENTS FOUND IN NERVE AND MUSCLES.
 605002 DC VOLTAGES MEASURED..... +/-1000V (VCLTS TO +/-10VLT
 605002 DC VOLTAGES MEASURED..... +/-1 FA TO +/-1 PA
 605002 DRIFT..... +/-5V/DAY AFTER 10MIN WARMUP
 605002 OUTPUT VOLTAGE..... +/-10V
 605002 FLUID SUPPLY SYSTEM SPOC1 11114
 605002 THIS UNIT CONSISTS OF THREE ELEMENTS SUPPORTING ALL MATERIAL
 605002 SCIENCE AND MANUFACTURING EXPERIMENTS. IT IS PART OF THE CCH
 605002 SHELVEFLY PRESSAT ON ALL MS/PS PAYLOADS. THE THREE ELEMENTS
 605002 ARE AN INERT GAS UNIT, AN OXIDIZING GAS UNIT, AND A RECLING
 605002 GAS UNIT.
 605002 GAS JET POSITIONING PROBE AND DETECTOR SPOC4 11114
 605002 USE OF GAS JETS TO POSITION SPECIMEN DURING EXPERIMENT PRE-
 605002 CIGURES SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL.
 605002 GAS FILTER EF-06 INCHES (0.2 MICROMETERS)
 605002 GASES HELIUM, HYDROGEN, OXYGEN
 605002 CONFIGURATION JET PATH IN APEX OF TETRAPOD
 605002 GAS PRESSURE MINIMUM 0.014 PSIA (100 MP/SC P)
 605002 INDICATORS-CAUTION AND WARNING ASO10 11112
 605002 A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF
 605002 SPSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.
 605002 INDICATORS AC
 605002 MASTER ALARM MEMORY
 605002 POWER/TEST CENTRLS
 605002 INDICATORS-CAUTION AND WARNING ME010 11112
 605002 A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF
 605002 SPSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.
 605002 INDICATORS AC
 605002 MASTER ALARM MEMORY
 605002 POWER/TEST CENTRLS
 605002 INDICATORS-CAUTION AND WARNING SCC10 11112
 605002 A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF
 605002 SPSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.
 605002 INDICATORS AC
 605002 MASTER ALARM MEMORY
 605002 POWER/TEST CENTRLS
 612015 LEMP - LOWER BODY NEGATIVE PRESSURE LS115 11111
 612015 MEASURE INTERNAL BODY PRESSURE DIFFERENTIAL.
 612015 PRESSURE DIFFERENTIAL 50 MP HG BELOW CARIN APERT
 612015 SWEPT-RECEIVER AND DEMODULATOR CAC07 41111
 612015 RECEIVE SIGNALS FROM TERRESTRIAL NOISE AND TERRESTRIAL INTERFER-
 612015 ENCE SOURCES AS A FUNCTION OF TIME OF DAY AND SEASON OF YEAR.
 612015 FREQUENCY RANGE 100 TO 1000 PHZ SWEEP RATE CA
 612015 TUNABLE
 612015 RECEIVER TYPE SUPERHETERODYNE, LINEAR, SWEPT
 612015 FREQUENCY, AGC
 612015 SENSITIVITY -120 DB
 612015 PRE-DETECTION BW ADJUSTABLE
 612015 DYNAMIC RANGE 60 DB INSTANTANEOUS
 612015 120 DB TOTAL
 612015 SWEPT RATE 100 RECEIVER BW PER HZ/LTE
 612015 RECEIVER NOISE FACTOR 4.0 DB
 612015 BANDWIDTH 1 MHZ
 612015 RECEIVER - RF ST127 11111
 612015 MEASURE ELECTROMAGNETIC INTERFERENCE AT ORBITAL ALTITUDES IN THE
 612015 FREQUENCY SPECTRUM OF 400 MHZ TO 19 GHZ. USED ON ST025 AND ST035.
 612015



610015 OVERALL SPECTRUM 400 MHZ TO 15 GHZ
610015 SPECIFIC MULTIPATH MEASURE..... 1.5 TO 1.6 GHZ
610015 2.025 TO 2.3 GHZ
610015 13.4 TO 15.35 GHZ
610015
610016 RECEIVER - BIOTELEMETRY LS980 11111
610016
610016 RECEPTION OF BIOTELEMETRY SIGNALS FROM DEEP BODY TEMPERATURE AND
610016 ANIMAL ACTIVITY AS WELL AS FOR ELECTROMAGNETIC FIELD MONITORING.
610016
610016 FREQUENCY RANGE TC 5 MHZ
610016 SENSITIVITY 0.5 MV FOR 20 DB GUESSING
610016 TUNING VARIABLE AND CRYSTAL CONTROLLED
610016
610020 RECEIVER, KU-BAND CAC47 11235
610020
610020 RECEIVE KU BAND SIGNAL FOR DETAILED EVALUATION OF PERFORMANCE
610020 CHARACTERISTICS OF TOPS/GROUND/SHUTTLE DATA LINK.
610020
610020 CARRIER FREQUENCY 13.4 TO 14.2 GHZ
610020 SENSITIVITY +/- 25 MHZ
610020 CARRIER/NOISE DENSITY RATIO .. 56 DB/MHZ
610020 DATA RATE 100 TO 1000 BPS
610020 DATA UP TO VIDEO
610020
610030 RECEIVER, L BAND CAC45 21111
610030
610030 RECEIVE L BAND SIGNALS TO DEMONSTRATE THE OPERATIONAL PERFORMANCE
610030 OF TYPICAL SPACEBORNE INTERFEROMETERS FOR AIRCRAFT AND MARINE
610030 NAVIGATION/TRAFFIC CONTROL, TRACK SIGNAL PHASE.
610030
610030 CARRIER FREQUENCY 1565.75 TO 1570.25 MHZ
610030 SENSITIVITY -120 DBM
610030 RECEIVER TYPE PHASE LOCK, SUPERHETERODYNE,
610030 FIXED FREQUENCY, AGC
610030 DYNAMIC RANGE 40 DB INSTANTANEOUS
610030 100 DB TOTAL
610030 PRE-DETECTION BW 5 MHZ
610030 SWEEP RATE TUNE 3 PPM OF CARRIER FREQ
610030 PHASE LOCK LOOP SNR 40 DB
610030 RECEIVER NOISE FIGURE 5 DB
610030 PHASE LOCK LOOP BANDWIDTH 10 MHZ
610030
610040 RECEIVER - NAVIGATION CAC96 11111
610040
610040 RECEIVE SIGNALS FROM DIRECT OPTICAL SENSING OF STARS TO INVESTIGATE
610040 THE PRECISION WITH WHICH A SATELLITE CAN BE ORIENTED VERTICALLY
610040 USING STARS AS THE SCALE INERTIAL REFERENCE.
610040
610040 CARRIER FREQUENCY 400 MHZ
610040 SENSITIVITY 2.0 MV FOR 10 DB SNR/A
610040 DATA RATE 2000 BPS
610040
610050 RECEIVER, VHF CAC46 51111
610050
610050 RECEIVE VHF SIGNALS FOR DETAILED EVALUATION OF PERFORMANCE CHARACTERISTICS
610050 OF TOPS/GROUND/SHUTTLE DATA LINK. PROVIDES SUPPORT OF
610050 RECEIVER FRONT ENDS WHICH CAN BE CONNECTED IN DIFFERENT RECEIVER
610050 CONFIGURATIONS.
610050
610050 CARRIER FREQUENCY 126 TO 130 MHZ
610050 SENSITIVITY +/- 25 MHZ
610050 CARRIER/NOISE DENSITY RATIO .. 53 DBM
610050 DATA RATE 100 TO 1000 BPS
610050
610060 RECEIVER, S-BAND CAC48 11111
610060
610060 RECEIVE S BAND SIGNAL FOR DETAILED EVALUATION OF PERFORMANCE
610060 CHARACTERISTICS OF TOPS/GROUND/SHUTTLE DATA LINK.
610060
610060 CARRIER FREQUENCY 2025 TO 2120 MHZ
610060 SENSITIVITY +/- 25 MHZ
610060 CARRIER/NOISE DENSITY RATIO .. 65 DB/MHZ
610060 DATA RATE 100 TO 1000 BPS
610060 VIDEO
610060
610070 RECEIVER, LASER ELECTRONICS CAC49 11235
610070
610070 PERFORMS DATA DEFORMATTING FROM LASER CARRIER REAP
610070
610080 RECEIVER, PHASE-LOCK CAC11 11111
610080
610080 PROVIDES FOR TRACKING SIGNAL PHASE.
610080
610080 FREQUENCY SPECTRUM 1570 MHZ TUNABLE
610080 SENSITIVITY 15 MHZ
610080
610090 ELECTROPHYSIOLOGY RECEIVER LS913 11111
610090
610090 PROVIDE SPECIAL RECEPTION OF CARDIOVASCULAR AND NEURAL ELECTRICAL
610090 PHYSIOLOGICAL EVENTS VIA BIOTELEMETRY SYSTEMS
610090
610090 SIGNAL TYPES RECEIVED.....ELECTROENCEPHALOGRAPHY
610090ELECTROCARDIOGRAPHY
610090VECTROCARDIOGRAPHY
610090PALISTECARDIOGRAPHY
610090IMPEDANCE CARDIOGRAPHY
610090PHONOCARDIOGRAPHY
610090
610090 BIOTELEMETRY RECEIVER - COMPACT CASE MODULE LS921 12113
610090
610090 MICROPACKS FOR SMALL VERTEBRATE ACQUIRING AND TRANSMITTING
610090 PIC DATA.
610090
610100 NOISE FIGURE TEST SET CAC57 11111
610100

ORIGINAL PAGE IS
OF POOR QUALITY



*1810 MEASURE SYSTEM NOISE FIGURE.
*1810 NOISE FIGURE RANGE 10 TO 14 DB
*1810 ACCURACY +/- 0.5 DB, C - 30 DB
*1810 INPUT FREQUENCY 100 TO 1000 MHZ
*1810 BANDWIDTH 1 MHZ
*19011 RADIATION SOURCE STORAGE LS979 11112
*19011 STORAGE FACILITY FOR RADIOACTIVE MATERIAL.
*19011 VOLUME APPROXIMATELY 1 CL FT
*19011 RADIATION PROTECTION LIMIT ... LP TO 500 MICROREM/HRS
*20001 TRACKING ANTENNA SERVO ELECTRONICS CPC44 11114
*20001 ELECTRONICS FOR ANTENNA POSITIONING AND CONTROL
*20001 TYPE TYPE II SERVO MOTOR
*20001 GAIN CONSTANT 3.564 PFR SEC SCLARED
*20001 BANDWIDTH 100 HZ
*20002 TRACKER ELECTRONICS - FINE CNO73 21114
*20002 PROVIDES ACCURATE TRACKING OF AN OPTICAL BEACON. PROVIDES CENTRAL
*20002 SIGNALS TO THE TELESCOPE CIPHERS AND VERNIER BEAM DEFLECTOR.
*20002 TYPE IMAGE DISSECTOR TRACKER
*20002 FOV +/- 300 MICRORAD
*20002 ACCURACY 0.5 MICRORAD
*20002 SENSITIVITY TO 1.00 MICRONS OR 10.0 MICRONS
*20003 COARSE TRACKER CNO84 11114
*20003 PROVIDES INITIAL GROUND BEACON ACQUISITION AND COARSE TRACKING.
*20003 DYNAMIC (SEARCH) FOV 1 DEG BY 1 DEG
*20003 ACCURACY (1 SIGMA) +/- 200 MICRORAD
*20003 TRACKING BANDWIDTH 1 HZ
*20003 SEARCH TIME 10 SEC MAX
*20003 SENSITIVITY TO 4.53 MICRONS
*20003 FALSE TARGET REJECTION MUST REJECT MOON, PLANETS,
STARS, AND SLLIT CLOUDS.
*20020 TRANSMITTER - 2.0 TO 20.0 MHZ AF907 11111
*20020 VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.
*20020 FREQUENCY BANDWIDTH 2.0 TO 20.0 MHZ
*20020 TRANSMITTER POWER 10 KW
*20021 TRANSMITTER - 0.2 TO 2.0 MHZ AF901 11111
*20021 VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.
*20021 FREQUENCY BANDWIDTH 0.2 TO 2.0 MHZ
*20021 TRANSMITTER POWER 10 KW
*20022 ELECTROSTATIC WAVE TRANSMITTER AF904 11134
*20022 TRANSMITTER, KU BAND CNO51 11235
*20025 PROVIDE KU BAND SIGNAL FOR EVALUATION OF PERFORMANCE CHARACTER-
*20025 ISTICS OF TORS/GROUND/SATLITE DATA LINK.
*20025 CARRIER FREQUENCY 14.4 TO 15.35 GHZ
*20025 STABILITY +/- 25 MHZ
*20025 TRANSMITTER POWER 100 WATTS
*20025 EFFECTIVE RADIATED POWER 33 DB
*20025 DATA RATE 1.0 MB
*20025 DUTY CYCLE VIDEO
*20025 SIGNAL BANDWIDTH 900 MHZ
*20030 TRANSMITTER, VHF CNO50 11111
*20030 PROVIDES MEDIUM POWER OUTPUT FOR TRANSMISSION TO TORS FOR EVAL-
*20030 UATION OF PERFORMANCE CHARACTERISTICS.
*20030 CARRIER FREQUENCY 130 TO 144 MHZ
*20030 STABILITY +/- 20 KHZ
*20030 OUTPUT POWER 10 WATTS MAX
*20030 DATA RATE 100 TO 10,000 BPS
*20030 SIGNAL BANDWIDTH 0.1 MHZ
*20030 DUTY CYCLE CONTINUOUS
*20050 TRANSMITTER - 0.3 TO 200 MHZ AF903 11111
*20050 VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.
*20050 FREQUENCY BANDWIDTH 0.3 TO 200.0 MHZ
*20050 TRANSMITTER POWER 1.0 KW
*20060 ELECTROPHYSIOLOGY BACKPACK LS912 11235
*20060 BACKPACK CONTAINS THE NECESSARY ELECTRONICS FOR SENSING AND
*20060 TRANSMITTING MAN'S PHYSIOLOGICAL DATA SUCH AS ECG, EEG, EMG,
*20060 ECG, ETC. INCLUDES SENSORS, SIGNAL CONDITIONERS, MULTIPLEXERS,
*20060 A/D CONVERTERS AND TRANSMITTERS
*20070 LIDAR TRANSMITTER ECCRP 11112
*20070 PROVIDES SIGNAL SOURCE TO MEASURE CLOUD HEIGHTS AND AEROSOL DIS-
*20070 TRIBUTION.
*20070 WAVE LENGTH 6328 ANGSTROMS



620070 ILLUMINATION MODES RUBY, OCLARED RUBY, CYE,
620070 OCLIPED CYE
62007C
622010 VIEWER - WIDE ANGLE MEC 11235
622010
62201C USED FOR LARGE AREA VIEWING AND ORIENTATION. INSTALMENT SIMILAR
622010 TO THE WILD NF2 NAVIGATION SIGHT USED WITH THE WILD-MERRUGG
622010 RC-10 METRIC CAMERA
62201C
622010 INSTANTANEOUS FIELD OF VIEW...11CDEG X 11CDEG
622010 ...55DEG X 55DEG
622010 ...20DEG X 20DEG
622010 FIELD OF VIEW.....36CDEG IN AZIMUTH
622010 POINTING ANGLE FROM NADIR
622010 (PITCH, ALONG-TRACK).....C TO 4CDEG
622010 TOTAL ANGULAR COVERAGE (CROSS-
622010 TRACK) FROM NADIR (1FOV/2)
622010 POINTING ANGLE.....+/-4CDEG
622010 POINTING ACCURACY (1-SIGMA)....2DEG
622010
623010 LOAD AND GIMBAL PLATFORM MONITOR AFB45 11134
623010
624010 MULTISPECTRAL SCANNER ELECTRONICS PACKAGE CFCC6 11134
624010
625010 ISOELECTRIC FOCUSING UNIT SPCC1 11114
625010
626010 PISTON PUMP WITH RF HEATING SFOC4 11235
626010
627010 FEED AND CRYSTAL HOLDER SPCC2 21113
627010
627020 PIEZOELECTRIC DRIVE SFOC2 21114
627020
628010 ELECTROMAGNETIC MIXING AND DISPERSAL UNIT SFOC2 21235
628010
628020 MECHANICAL MIXING AND DISPERSAL UNIT SPCC2 11235
628020
629010 INERTIAL INJECTOR SFCC4 11235
629010
629020 LIQUID SYRINGE DISPENSER SPOC4 21113
629020
630010 SOLID SAMPLE STORAGE SPOC4 11112
630010
631001 CRITICAL STATE TEST CELL ST207 11235
631001
631001 TEST CELL TO TEST CRITICAL POINT PHENOMENA. TEST CELL INCLUDES
631001 A THERMAL SHIELD, HEAT EXCHANGER, MICROMETER, CHAMBER, SUPPORT
631001 ELECTRONICS AND PRESSURE TRANSDUCERS
631001
631002 POOL BOILING TEST CHAMBER ST218 21235
631002
631002 CHAMBER TO EVALUATE POOL BOILING PHENOMENA IN ZERO GRAVITY
631002
631003 CRYSTALLIZER SAMPLE CHAMBER ST225 11135
631003
631003 CHAMBER FOR GROWING CRYSTALS IN ZERO G
631003
631004 SCHLIEREN TYPE OPTICS SYSTEM ST227 11135
631004
631004 OPTICS SYSTEM EMPLOYING SCHLIEREN METHOD TO EXAMINE THE STRUCTURE
631004 OF CRYSTALS
631004
631005 HOLOGRAPHIC INTERFEROMETER ST231 11235
631005
631005 DEVICE TO EXAMINE STRUCTURE OF GROWING CRYSTALS
631005

ORIGINAL PAGE IS
OF POOR QUALITY

APPENDIX C

COMMERCIAL EQUIPMENT SHOPPING LIST



ASTRONOMY

104180 DISPLAY CRT (CONTROL/DISPLAY CONSOLE) AS001 21111
104180
104180 PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND
104180 STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS, INCLUDES
104180 DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.
104180
104180 SCREEN SIZE 14 IN (10.355 IN)
104180
RESEARCH INC #3300 \$1580
FORMAT 24 LINESX72 OR 80 CHARACTERS
..... 24 LINESX40 CHARACTERS
..... 12 LINESX72 OR 80 CHARACTERS
REFRESH RATE 60 HZ
TRANSFER RATE 110 TO 2400 BAUD
CHARACTER FORM 5X7 DOT MATRIX
DIMENSIONS 15-1/2WX13-1/2HX23-1/20 INCH
WEIGHT 39LB (17.7KG)
HEWLETT PACKARD #2400A \$3580 111
RESEARCH INC #812-3301 \$1555 111
SCREEN FORMAT 12 LINES X 72 OR 80 CHARACTERS,
24 LINES X 40 CHARACTERS
CHARACTER FORMAT 5 X 7 DOT MATRIX
TRANSFER RATE 110 TO 2400 BAUD, 10 CR 11 BIT
CHARACTERS
MODES HALF OR FULL DUPLEX-SWITCHABLE
LOCAL OR REMOTE
RESEARCH INC #DPS-812 \$21,700 111
COMPLETE DATA DISPLAY SYSTEM MADE UP OF CRT ALPHANUMERIC DISPLAY
DEC-PDP8-F MINS-COMPUTER, #3301 OPERATIONS CONTROL CONSOLE,
PRINTER/CARD READER, #812-4A ANALYZER AND #812-13 UNIDRIVER
CHANNELS 32
MEMORY SIZE 8K WORDS
TEKTRONIX INC #4002A/021-0033-00/ \$10975 111
021-00XX-00/4901/
4951
ALPHANUMERIC FORMAT 39 LINES-85 NORMAL/ITALIC CHAR,
1 LINE- 84 CHAR IN SCRATCH AREA
CHARACTER SET 96 UPPER/LOWER CHAR (ASCII)
CHARACTER SIZE 70 X 90 MILS
CHARACTER GENERATION 7 X 5 DOT MATRIX
CURSOR PULSATING 7 X 9 MATRIX
GRAPHIC INPUT MODE 1024(X) X 1024(Y) ADDRESSABLE,
1024(X) X 768(Y) VIEWABLE PTS;
JOYSTICK CONTROLLED, CROSS-HAIR
CURSOR
DISPLAY MEDIUM 11 IN DIA CRT
DISPLAY AREA 8.3IN.HOR.X6.1IN.VERT.
ACCUMULATORS 4
CYCLE TIME 300 NANOSEC
CORE MEMORY 32 K
BUNKER-RAND CORP #2217-12CRT 111
CHARACTER CAPACITY (MAX) 960
CHARACTERS/LINE (MAX) 80
LINES/DISPLAY (MAX) 24
CHARACTERS REPERTOIRE 62 OR 92
VIEWING AREA 8.75 X 6.25
REFRESH AREA 54 FRAMES/SEC
104020 SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) AS002 11111
104020
104020 PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CRT
104020
104020 CHANNELS 2
104020 FORMAT VIDEO AND DATA
104020 VIDEO INPUT PRESENTATION RASTER SCAN
104020 SYMBOL WRITING TECHNIQUE STROKE
104020 INTERFACE DESCRIPTION DIGITAL-12-BIT DATAWORD
104020
TEKTRONIX #4002A/021-00XX-00 \$9550 111
DISPLAY MEDIUM 11 INCH DIRECT-VIEW, BISTABLE
STORAGE CRT WITH REFRESHED
SCRATCH PAD AREA
DISPLAY AREA 8.3 INCHES HORIZONTAL X 6.1
INCHES VERTICAL
ALPHANUMERIC MODE
FORMAT 39 LINES OF 85 NORMAL OR ITALIC
CHARACTERS IN MAIN AREA, ONE
LINE OF 84 CHARACTERS IN
CHARACTER SET 96 UPPER AND LOWER CASE PRINT-
ING CHARACTERS (ASCII CCODE)
CHARACTER SIZE 70X90 MILS (CAN BE DOUBLE SIZE)
CHARACTER GENERATION 7X9 DOT MATRIX
CURSOR PULSATING 7X9 MATRIX
GRAPHIC MODES LINEAR INTERPOLATE, INCREMENTAL
PLOT, POINT PLOT, 1024X1024
ADDRESSABLE POINTS, 1024X768
VIEWABLE POINTS
GRAPHIC INPUT MODE 1024 (X), 768 (Y) POINTS, JOY-
STICK CONTROLLED, CROSS-HAIR
CURSOR

ORIGINAL PAGE IS
OF POOR QUALITY

101340 FUNCTION KEYBOARD AS003 11111
101340
101340 ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO
101340 DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY,
101340 FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
BUKER-RAND CORP #2200 111
KEYS TYPEWRITER SET STD
PROGRAM ASSIST KEYS 16
INTERFACE CRT
EDITING KEYS FULL SET
101340
101340 ALPHANUMERIC KEYBOARD AS004 11111
101340
101340 ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE ONBOARD COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
BUKER-RAND CORP #2200 #440 111
KEYS TYPEWRITER SET STD
PROGRAM ASSIST KEYS 16
INTERFACE CRT
EDITING KEYS FULL SET
104000 MICROFILM VIEWER AS005 11111
104000
104000 PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUB-
104000 SYSTEM OPERATIONAL PROCEDURES, ONBOARD CHECKOUT PROCEDURES,
104000 SIMPLEX SCHEMATICS AND OTHER WRITTEN OR PICTORIAL INFORMATION.
104000
104000 FILM FORMAT 16 MM DUAL TRACK
104000 FILM LOADING CASSETTE
104000 FILM SLEWING MANUAL AND COMPUTER SELECT
104000 INTERFACE DESCRIPTION DIGITAL 13-BIT BINARY
104000
CALMA CO. #3C3H #28000 111
FUNCTION OFF LINE DIGITIZING SYSTEM
WITH FILM PROJECTION,
PROVIDES CAPABILITY OF
DIGITIZING FILM DATA
498200 TIMER - EVENT AS006 11111
498200
498200 DIGITAL DISPLAY OF TIME REMAINING OR EXPENDED FOR A PARTICULAR
498200 EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE
498200 START AND STOP COMMANDS (MINUTES AND SECONDS).
498200
498200 DISPLAY 4 DIGITS
498200 DISPLAY TYPE LED
498200
ROHDE AND SCHWARTZ #CAD 1100.6597.911 / #5440 111
DISPLAY 4 DIGITS: HR, MIN, SEC
OUTPUT 1-OUT-OF-N AND/OR BCD CDOF
PROGRAMMER TIME PROGRAMS (COUNT DOWN) AND
TIME SIGNALS
DATATRON #3350-506 #1170 111
INPUT FREQUENCY 60HZ
OUTPUTS VISUAL, BCD
DISPLAY HRS, MIN, SEC
TENNELEC #TC545 #550 111
COUNT RATE 20MHZ
TIMEBASE 0.1 OR 0.01 SEC
ACCURACY SAME AS LINE FREQUENCY
CONFIGURATION NIM COMPATIBLE
498500 TIMER - MISSION AS007 11111
498500
498500 PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECOND UP-
498500 DATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE,
498500 SECOND.
498500
498500 DISPLAY 7 DIGITS BCD
498500 DISPLAY TYPE LED
498500
ROHDE AND SCHWARTZ #HS 62P4110.3226.211 #12,700 111
DISPLAY 9 DIGITS: DAY, HR, MIN, SEC
DATATRON #3350 #1020 111
INPUT FREQUENCY 60HZ
OUTPUTS VISUAL, BCD
DISPLAY HRS, MIN, SEC
DATATRON #3030 #7940 111
CHANNELS 32
TIMING RANGE 5NANOSEC TO 5MICROSEC
RESOLUTION 1.0E-8
320000 CONTROLLER, THREE-AXIS HAND AS008 11132
320000
320000 THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT
320000 POINTING AND INITIAL TARGET ACQUISITION.
320000
320000 CONTROLLED AXES 3
320000



603006 CIRCUIT BREAKER/DISTRIBUTOR PANEL AS009 11131
603006
603006 CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO
603006 THOSE STANDARD ITEMS FOUND IN ALL ASTRONOMY EXPERIMENTS.
603006
603006 CIRCUIT BREAKERS 1P
603006 POWER CAPABILITY 10 AMPS
603006
609003 INDICATORS-CAUTION AND WARNING AS010 11132
609003
609003 A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF
609003 SUBSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.
609003
609003 INDICATORS 40
609003 MASTER ALARM MEMORY
609003 POWER/TEST CONTROLS
609003
424000 RECORDER-TAPE AS011 11111
424000
424000 RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND
424000 DISPLAY CONSOLE.
424000
424000 BANDWIDTH DC TO 200 KHZ
424000 CHANNELS 14
424000 TIME MARKS EDGE TRACK
424000 RECORDED FORMAT DIGITAL
424000
HEWLETT PACKARD #7670B/C \$4600 111
TAPE FORMAT 800, 556 OR 200 CPI NRZI AND
1600 CPI PHASE ENCODED
TRACKS 7 OR 9
TAPE SPEED 10 TO 45 IPS
TAPE 0.5 INCH, 1.5 MIL. IRP/ANSI
COMPATIBLE
SANGAM ELECTRIC #SAHER 111 \$26,500 111
TRACKS 14
TAPE WIDTH 1 INCH
TAPE SPEEDS A SELECTABLE SPEEDS FROM 15/16
TO 120 IPS
FREQUENCY RESPONSE 400 HZ TO 2.0 MHZ
RECORDING RATE 600 KBPS AT 120 IPS SERIAL MODE
WEIGHT 100 POUNDS
LEACH #MTR 700C 311
TRACKS 12
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
TAPE LENGTH 9200 FEET
PACKING DENSITY 16.7 KB/I/T
SIGNAL/NOISE 22 DB
DATA CAPACITY 2.2E+10
BANDWIDTH 2 MB/SEC/T
BORG WARNER #PERT 121
TRACKS 30
TAPE SPEED UP TO 1000 IPS
TAPE WIDTH 1/2 INCH
TAPE LENGTH 2400 FEET
PACKING DENSITY 15 KB/I/T
SIGNAL/NOISE 24 DB
DATA CAPACITY 1.3E+9
BANDWIDTH 6-15 MB/SEC/T
AMPEX CORP #AR 1700 \$28,000 112
TRACKS 14
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KB/I/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 6.2E+10 BITS
BANDWIDTH 2 MB/SEC/T
TAPE LENGTH 9200 FEET
AMPEX CORPORATION #AR-200 121
BANDWIDTH 100 HZ TO 3125 HZ THROUGH
300 HZ TO 250 KHZ
TAPE WIDTH 0.5 INCH
TAPE SPEED 1.75 TO 60 IPS
RECORDING TIME 8 MINUTES TO 4 HRS AND 16 MINS
FORMAT DIGITAL
TRACKS 8 DIGITAL, 7 ANALOG
AMPEX CORP #AR 70C \$29,100 511
TRACKS 14
TAPE SPEED 60 IPS
TAPE WIDTH 1 INCH
RECORD MODE DIRECT
PACKING DENSITY 20 KB/I/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 8E+10 BITS
BANDWIDTH 1 MB/SEC/T
DIGI-DATA CORPORATION #1600 \$2650 111
TAPE SPEED 25, 18.75, 12.5 IPS
TRACKS 7 OR 9 TRACK
DATA DENSITY 1600 CPI PHASE ENCODED
200, 556, 800 CPI NRZI
TAPE 0.5 INCH, 1.5 MIL, 1200 FEET

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

DIGI-DATA CORPORATION #1700/POP-11/7-9 95250 111
TRACK NR/1
TAPE SPEED 45, 37.5, 25, 18.75, 12.5 IPS
TAPE 0.5 INCH, 1.5 MIL, 18"/ANSI
COMPATIBLE, 10.5 INCH REEL
TRACKS 7 OR 9 TRACK
DATA DENSITY PHASE ENCODED

HONEYWELL #96 #17420
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 240IPS
REEL SIZE 16IN
MAXIMUM BANDWIDTH(DIRECT) 24 KHZ
TAPE WIDTH 1/2 IN

HONEYWELL #5600 #9730
PORTABLE TAPE RECORDER
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 60IPS
MAXIMUM BANDWIDTH(DIRECT) 300KHZ
PACKING DENSITY UP TO 60CPH
WEIGHT 70LBS (32KG)
INPUT VOLTAGE 28VDC

HEWLETT PACKARD #34450 #10200 111
BANDWIDTH 300 KHZ
CHANNELS 7
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #34450 #14700 111
BANDWIDTH 300 KHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3450A-011 #23800 111
BANDWIDTH 500 HZ TO 2 MHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3460A/13065B/13063A #4796 111
CONFIGURATION RACK MOUNTED
TAPE SPEED 15, 3 AND 3/4, 15/16 IPS
CHANNELS 6
RECORDING FORMAT FM
PASSBAND 5 KHZ
S/N RATIO 48 DB
TRACKS 30

601002 ANNUNCIATORS-ADVISORY AS013 21132
601002
601002 DUAL BANKS MOUNTED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL
601002 ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY CHORD
601002 EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.
601002
601002 ANNUNCIATOR-VISUAL 12 IN DUAL BANKS
601002 DISPLAY TYPE RED, WHITE, GREEN
601002 ANNUNCIATOR-AUDIO HORN
601002

HIGH ENERGY PHYSICS

104180 DISPLAY CRT (CONTROL/DISPLAY CONSOLE) HE001 21111
104180 PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERIC, DYNAMIC AND
104180 STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS. INCLUDES
104180 DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.
104180
104180 SCREEN SIZE 14 IN (0.355 M)
104180
RESEARCH INC #3300 \$1580
FORMAT 24 LINES X 72 OR 80 CHARACTERS
..... 24 LINES X 40 CHARACTERS
..... 12 LINES X 72 OR 80 CHARACTERS
REFRESH RATE 60 HZ
TRANSFER RATE 110 TO 2400 BAUD
CHARACTER FORM 5 X 7 DOT MATRIX
DIMENSIONS 15-1/2 X 13-1/2 X 23-1/2 INCH
WEIGHT 39 LB (17.7 KG)
HEWLETT PACKARD #2600A \$3580 111
RESEARCH INC #R12-3301 \$1555 111
SCREEN FORMAT 12 LINES X 72 OR 80 CHARACTERS,
24 LINES X 40 CHARACTERS
CHARACTER FORMAT 5 X 7 DOT MATRIX
TRANSFER RATE 110 TO 2400 BAUD, 10 OR 11 BIT
CHARACTERS
MODES HALF OR FULL DUPLEX-SWITCHABLE
LOCAL OR REMOTE
RESEARCH INC #DPS-R12 \$21,700 111
COMPLETE DATA DISPLAY SYSTEM MADE UP OF CRT ALPHANUMERIC DISPLAY
DEC-PDP8-F MINS-COMPUTER, #3301 OPERATIONS CONTROL CONSOLE,
PRINTER/CARD READER, #R12-4A ANALYSETER AND #R12-13 UNIT DRIVER
CHANNELS 32
MEMORY SIZE 8K WORDS
104020 SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) HE002 11111
104020 PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CRT
104020
104020 CHANNELS 2
104020 FORMAT VIDEO AND DATA
104020 VIDEO INPUT PRESENTATION RASTER SCAN
104020 SYMBOL WRITING TECHNIQUE STROKE
104020 INTERFACE DESCRIPTION DIGITAL-12-BIT DATAWORD
104020
TEKTRONIC #4002A/021-00XX-00 \$9550 111
DISPLAY MEDIUM 11 INCH DIRECT-VIEW, BISTABLE
STORAGE CRT WITH REFRESHED
SCRATCH PAD AREA
DISPLAY AREA 8.3 INCHES HORIZONTAL X 6.1
INCHES VERTICAL
ALPHANUMERIC MODE
FORMAT 39 LINES OF 85 NORMAL OR ITALIC
CHARACTERS IN MAIN AREA, ONE
LINE OF 84 CHARACTERS IN
SCRATCH PAD AREA
CHARACTER SET 96 UPPER AND LOWER CASE PRINT-
ING CHARACTERS (ASCII CODE)
CHARACTER SIZE 7 X 5 MILS (CAN BE DOUBLE SIZE)
CHARACTER GENERATION 7 X 5 DOT MATRIX
CURSOR PULSATING 7 X 9 MATRIX
GRAPHIC MODES LINEAR INTERPOLATE, INCREMENTAL
PLOT, POINT PLOT, 1024 X 1024
ADDRESSABLE POINTS, 1024 X 768
VIEWABLE POINTS
GRAPHIC INPUT MODE 1024 (X), 768 (Y) POINTS, JOY-
STICK CONTROLLED, CROSS-HAIR
CURSOR
101340 FUNCTION KEYBOARD HE003 11111
101340 ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO
101340 DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY,
101340 FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
BUKKER-RAND CORP #2200 \$440 111
KEYS TYPEWRITER SET STD
PROGRAM ASSIST KEYS 16
INTERFACE CRT
EDITING KEYS FULL SET
101340 ALPHANUMERIC KEYBOARD HE004 11111
101340 ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE ONBOARD COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
BUKKER-RAND CORP #2200 \$440 111
KEYS TYPEWRITER SET STD
PROGRAM ASSIST KEYS 16
INTERFACE CRT
EDITING KEYS FULL SET

ORIGINAL PAGE IS
OF POOR QUALITY



164000 MICROFILM VIEWER HE005 11111
164000
164000 PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUB-
164000 SYSTEM OPERATIONAL PROCEDURES, ONBOARD CHECKOUT PROCEDURES,
164000 SIMPLEX SCHEMATICS AND OTHER WRITTEN OR PICTORIAL INFORMATION.
164000
164000 FILM FORMAT 16 MM DUAL TRACK
164000 FILM LOADING CASSETTE
164000 FILM SLEWING MANUAL AND COMPUTER SELECT
164000 INTERFACE DESCRIPTION DIGITAL 13-BIT BINARY
164000
CALMA CO. #3C3M #28000 111
FUNCTION OFF-LINE DIGITIZING SYSTEM
WITH FILM PROJECTION.
PROVIDES CAPABILITY OF
DIGITIZING FILM DATA

498200 TIMER - EVENT HE006 11111
498200
498200 DIGITAL DISPLAY OF TIME REMAINING OR EXPENDED FOR A PARTICULAR
498200 EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE
498200 START AND STOP COMMANDS (MINUTES AND SECONDS).
498200
498200 DISPLAY 4 DIGITS
498200 DISPLAY TYPE LED
498200 INTERFACE DISCRETE
498200
ROMBE AND SCHWARZ #CAD (100.6597.91) #5400 111
DISPLAY 6 DIGITS: HRS, MINS, SECS
OUTPUT 1-OUT-OF-N AND/OR RCD CODE

CATATRON #339C-506 #1170 111
INPUT FREQUENCY 60MHZ
OUTPUTS VISUAL, BCD
DISPLAY HRS, MIN, SEC

TENNELEC #TC945 #550 111
COUNT RATE 20MMZ
TIMERASE 0.1 OR 0.01SEC
ACCURACY SAME AS LINE FREQUENCY
CONFIGURATION NIM COMPATIBLE

498500 TIMER - MISSION HF007 11111
498500
498500 PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECOND UP-
498500 DATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE,
498500 SECOND.
498500
498500 DISPLAY 7 DIGITS BCD
498500 DISPLAY TYPE LED
498500
ROMBE AND SCHWARZ #HS 62R4(110.3226.21) #12,700 111
DISPLAY 9 DIGITS: DAY, HR, MIN, SEC

CATATRON #3350 #1020 111
INPUT FREQUENCY 60MHZ
OUTPUTS VISUAL, BCD
DISPLAY HRS, MIN, SEC

DATATRON 3030 #7940 111
CHANNELS 32
TIMING RANGE 9NANOSEC TO 5MICROSEC
RESOLUTION 1.0E-8

320000 CONTROLLER, THREE-AXIS HAND HE008 11132
320000
320000 THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT
320000 POINTING AND INITIAL TARGET ACQUISITION.
320000
320000 CONTROLLED AXES 3
320000
603006 CIRCUIT BREAKER/DISTRIBUTOR PANEL HE009 11131
603006
603006 CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO
603006 THOSE STANDARD ITEMS FOUND IN ALL ASTRONOMY EXPERIMENTS.
603006
603006 CIRCUIT BREAKERS 18
603006 POWER CAPABILITY 10 AMPS
603006
609003 INDICATORS-CAUTION AND WARNING HE010 11132
609003
609003 A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF
609003 SUBSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.
609003
609003 INDICATORS 40
609003 MASTER ALARM MEMORY
609003 POWER/TEST CONTROLS
609003
424000 RECORDER-TAPE HE011 11111
424000
424000 RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND
424000 DISPLAY CONSOLE.
424000
424000 BANDWIDTH DC TO 200 KHZ
424000 CHANNELS 14
424000 TIME MARKS EDGE TRACK
424000 RECORDED FORMAT DIGITAL
424000



Space Division
Rockwell International

HEWLETT PACKARD	#7570B/C	\$6400	111
TAPE FORMAT	800, 556, OR 200 CPI NRZI AND 1600 CPI PHASE ENCODED		
TRACKS	7 OR 9		
TAPE SPEED	10 TO 45 IPS		
TAPE	0.5 INCH, 1.5 MIL, 1PP/ANSI COMPATIBLE		
SANGAMO ELECTRIC	#SARER 111	\$26,500	111
TRACKS	14		
TAPE WIDTH	1 INCH		
TAPE SPEEDS	8 SELECTABLE SPEEDS FROM 15/16 TO 120 IPS		
FREQUENCY RESPONSE	400 HZ TO 2.0 MHZ		
RECORDING RATE	600 KBPS AT 120 IPS SERIAL MODE		
WEIGHT	100 POUNDS		
LEACH	#MTR 700C		311
TRACKS	12		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
TAPE LENGTH	9200 FEET		
PACKING DENSITY	16.7 KB/I/T		
SIGNAL/NOISE	22 DB		
DATA CAPACITY	2.2E+10		
BANDWIDTH	2 MB/SEC/T		
BORG WARNER	#BPERT		121
TAPE SPEED	UP TO 1000 IPS		
TRACKS	30		
TAPE WIDTH	1/2 INCH		
TAPE LENGTH	2400 FEET		
PACKING DENSITY	15 KB/I/T		
SIGNAL/NOISE	24 DB		
DATA CAPACITY	1.3E+9		
BANDWIDTH	6-15 MB/SEC/T		
AMPEX CORP	#AR 70C	\$29,193	511
TRACKS	14		
TAPE SPEED	60 IPS		
TAPE WIDTH	1 INCH		
RECORD MODE	DIRECT		
PACKING DENSITY	20 KB/I/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	8E+10 BITS		
BANDWIDTH	1 MB/SEC/T		
AMPEX CORP	#AR 17C0	\$28000	112
TRACKS	28		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
TAPE LENGTH	9200 FEET		
PACKING DENSITY	20 KB/I/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	6.2E+10		
BANDWIDTH	2 MB/SEC/T		
AMPEX CORPORATION	#AR-200		121
BANDWIDTH	100 HZ TO 3125 HZ THROUGH 300 HZ TO 250 KHZ		
TAPE WIDTH	0.5 INCH		
TAPE SPEED	1.875 TO 60 IPS		
RECORDING TIME	8 MINUTES TO 4 HRS AND 15 MINS		
FORMAT	DIGITAL		
TRACKS	8 DIGITAL, 7 ANALOG		
DIGI-DATA CORPORATION	#1600	\$26,000	111
TAPE SPEED	25, 18.75, 12.5 IPS		
TRACKS	7 OR 9 TRACK		
DATA DENSITY	1600 CPI PHASE ENCODED		
TAPE	0.5 INCH, 1.5 MIL, 1000 FEET		
DIGI-DATA CORPORATION	#1700/PDP-11/7-9 TRACK NRZI	\$5250	111
TAPE SPEED	45, 37.5, 25, 18.75, 12.5 IPS		
TAPE	0.5 INCH, 1.5 MIL, 1PP/ANSI COMPATIBLE, 10.5 INCH REEL		
TRACKS	7 OR 9 TRACK		
DATA DENSITY	PHASE ENCODED		
HONEYWELL	#96	\$17420	
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 240IPS		
REEL SIZE	16IN		
MAXIMUM BANDWIDTH(DIRECT)	20 KHZ		
TAPE WIDTH	1/2 IN		
HONEYWELL	#9600	\$9730	
PORTABLE TAPE RECORDER			
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 60IPS		
MAXIMUM BANDWIDTH(DIRECT)	300KHZ		
PACKING DENSITY	UP TO 600PI		
WEIGHT	70LBS (32KG)		
INPUT VOLTAGE	28VDC		

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

HEWLETT PACKARD #34950 810200 111
BANDWIDTH 300 KHZ
CHANNELS 7
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #34950C 814700 111
BANDWIDTH 300 KHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3490A-C11 823800 111
BANDWIDTH 900 HZ TO 2 KHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3460A/13065A/13063A 84786 111
CONFIGURATION RACK MOUNTED
TAPE SPEED 15, 3 AND 3/4, 15/16 IPS
CHANNELS 4
RECORDING FORMAT FM
PASSBAND 5 KHZ
S/N RATIO 48 DB
TRACKS 30

001002 ANNUNCIATORS-ADVISORY HE012 21132
001002
001002 DUAL BANKS MOUNTED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL
001002 ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY OBCARD
001002 EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.
001002
001002 ANNUNCIATOR-VISUAL 12 IN DUAL BANKS
001002 DISPLAY TYPE RED, WHITE, GREEN
001002 ANNUNCIATOR-AUDIO HORN
001002



SCALAR PHYSICS

104180 DISPLAY CRT (CONTROL/DISPLAY CONSOLE) SC001 21111
104180
104180 PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND
104180 STATIC GRAPHICS, VECTORS, CIRCLES, AND SPECIAL SYMBOLS. INCLUDES
104180 DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.
104180
104180 SCREEN SIZE 14 IN (0.355 M)
104180
104180 HEWLETT PACKARD #2800A \$3580 111
RESEARCH INC #3300 \$1580
FORMAT.....24 LINESX72 OR 80 CHARACTERS
.....24 LINESX40 CHARACTERS
.....12 LINESX72 OR 80 CHARACTERS
REFRESH RATE.....60 HZ
TRANSFER RATE.....110 TO 2400 BAUD
CHARACTER FORM.....5X7 DOT MATRIX
DIMENSIONS.....15-1/2MX13-1/2HX23-1/20 INCH
WEIGHT.....39LB (17.7KG)
RESEARCH INC #P12-3301 \$1555 111
SCREEN FORMAT 12 LINES X 72 OR 80 CHARACTERS,
24 LINES X 40 CHARACTERS
CHARACTER FORMAT 5 X 7 DOT MATRIX
TRANSFER RATE 110 TO 2400 BAUD, 10 OR 11 BIT
CHARACTERS
MODES HALF OR FULL DUPLEX-SWITCHABLE
LOCAL OR REMOTE
RESEARCH INC #DPS-P12 \$21,700 111
COMPLETE DATA DISPLAY SYSTEM MADE UP OF CRT ALPHANUMERIC DISPLAY
DEC-PDP-11 F MINS-COMPUTER, #3301 OPERATIONS CONTROL CONSOLE,
PRINTER/CARD READER, #P12-4A ANALYSER AND #P12-13 UNIDIVER
CHANNELS 22
MEMORY SIZE 8K WORDS
104020 SYMBOL GENERATOR UNIT (CONTROL/DISPLAY CONSOLE) SC002 11111
104020
104020 PROVIDES VIDEO AND COMPUTER DATA INTERFACE TO THE CRT
104020
104020 CHANNELS 2
104020 FORMAT VIDEO AND DATA
104020 VIDEO INPUT PRESENTATION RASTER SCAN
104020 SYMBOL WRITING TECHNIQUE STROKE
104020 INTERFACE DESCRIPTION DIGITAL-12-BIT DATAWORD
104020
TENTRONIC #4002A/C21-COAX-00 \$9550 111
DISPLAY MEDIUM 11 INCH DIRECT-VIEW, RISTABLE
STOPAGE CRT WITH REFRESHED
SCRATCH PAD AREA
DISPLAY AREA 8.2 INCHES HORIZONTAL X 6.1
INCHES VERTICAL
ALPHANUMERIC MODE
FORMAT 39 LINES OF 85 NORMAL OR ITALIC
CHARACTERS IN MAIN AREA, ONE
LINE OF 84 CHARACTERS IN
SCRATCH PAD AREA
CHARACTER SET 96 UPPER AND LOWER CASE PRINT-
ING CHARACTERS (ASCIT CODE)
CHARACTER SIZE 7CX90 MILS (CAN BE DOUBLE SIZE)
CHARACTER GENERATION 7X9 DOT MATRIX
CURSOR PULSATING 7X9 MATRIX
GRAPHIC MODES LINEAR INTERPOLATE, INCREMENTAL
PLOT, POINT PLOT, 1024X1024
ADDRESSABLE POINTS, 1024X761
VIEWABLE POINTS
GRAPHIC INPUT MODE 1024 (X), 761 (Y) POINTS, JOY-
STICK CONTROLLED, CROSS-HAIR
CURSOR
101340 FUNCTION KEYBOARD SD003 11111
101340
101340 ALLOWANCE BY CREWMAN TO CONFIGURE EXPERIMENTS AND SUBSYSTEMS INTO
101340 DESIRED OPERATING MODES. THESE INCLUDE SELECTION OF CATEGORY,
101340 FUNCTION, MODE, STATUS AND COMMON KEYBOARD FUNCTIONS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
BUNKER-RAMO CORP #2200 \$440 111
KEYS TYPEWRITER SET STD
PROGRAM ASSIST KEYS 16
INTERFACE CRT
EDITING KEYS FULL SET
101340 ALPHANUMERIC KEYBOARD SD004 11111
101340
101340 ALLOWANCE BY CREWMAN TO COMMUNICATE WITH THE ONBOARD COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION DIGITAL-32-BIT WORD
101340
BUNKER-RAMO CORP #2200 \$440 111
KEYS TYPEWRITER SET STD
PROGRAM ASSIST KEYS 16
INTERFACE CRT
EDITING KEYS FULL SET

ORIGINAL PAGE IS
OF POOR QUALITY



164000 MICROFILM VIEWER 50005 11111
164000
164000 PROVIDES READ-ONLY, PROCEDURAL-TYPE DATA FOR EXPERIMENT AND SUB-
164000 SYSTEM OPERATIONAL PROCEDURES, ONBOARD CHECKOUT PROCEDURES,
164000 SIMPLEX SCHEMATICS AND OTHER WRITTEN OR PICTORIAL INFORMATION.
164000
164000 FILM FORMAT 16 MM DUAL TRACK
164000 FILM LOADING CASSETTE
164000 FILM SLEWING MANUAL AND COMPUTER SELECT
164000 INTERFACE DESCRIPTION DIGITAL 13-BIT BINARY
164000
CALHA CO. #303M 628000 111
FUNCTION OFF LINE DIGITIZING SYSTEM
WITH FILM PROJECTION.
PROVIDES CAPABILITY OF
DIGITIZING FILM DATA
498200 TIMER - EVENT 50006 11111
498200
498200 DIGITAL DISPLAY OF TIME REMAINING OR EXPENDED FOR A PARTICULAR
498200 EVENT. TIMER CAPABLE OF COUNTING UP OR DOWN AND PROVIDES DISCRETE
498200 START AND STOP COMMANDS (MINUTES AND SECONDS).
498200
498200 DISPLAY 4 DIGITS
498200 DISPLAY TYPE LED
498200
ROMBE AND SCHWARTZ #CAD (100.6597.91) 55400 111
DISPLAY 6 DIGITS: HR, MIN, SEC
OUTPUT 1-OUT-OF-N AND/OR BCD CODE
PROGRAMMER TIME PROGRAMS (CCCLAT
DOWN) AND TIME SIGNALS
DATATRON #3350-5C6 91170 111
INPUT FREQUENCY 60HZ
OUTPUTS VISUAL, BCD
DISPLAY HRS, MIN, SEC
TENNELEC #TC545 9550 111
COUNT RATE 20MHZ
TIMERBASE 0.1 OR 0.01SEC
ACCURACY SAME AS LINE FREQUENCY
CONFIGURATION NIM COMPATIBLE
498500 TIMER - MISSION 50007 11111
498500
498500 PROVIDES TIME REFERENCE IN GREENWICH MEAN TIME WITH 1 SECOND UP-
498500 DATE MAINTAINED VIA DATA MANAGEMENT COMPUTER - DAY, HOUR, MINUTE,
498500 SECOND.
498500
498500 DISPLAY 7 DIGITS BCD
498500 DISPLAY TYPE LED
498500
ROMBE AND SCHWARTZ #MS 6284(110.3226.21) 912.700 111
DISPLAY 9 DIGITS: DAY, HR, MIN, SEC
DATATRON #3350 91020
INPUT FREQUENCY 60HZ
OUTPUTS VISUAL, BCD
DISPLAY HRS, MIN, SEC
DATATRON 3030 97940 111
CHANNELS 32
TIMING RANGE 5NANOSEC TO 5MICROSEC
RESOLUTION 1.0E-8
320000 CONTROLLER, THREE-AXIS HAND 50008 11132
320000
320000 THREE-AXIS MULTIFUNCTION HAND CONTROLLER PROVIDED FOR INSTRUMENT
320000 POINTING AND INITIAL TARGET ACQUISITION.
320000
320000 CONTROLLED AXES 3
320000
603006 CIRCUIT BREAKER/DISTRIBUTOR PANEL 50009 11131
603006
603006 CENTRAL LOCATION OF CIRCUIT BREAKERS AND POWER DISTRIBUTION TO
603006 THOSE STANDARD ITEMS FOUND IN ALL ASTRONOMY EXPERIMENTS.
603006
603006 CIRCUIT BREAKERS 18 EACH
603006 POWER CAPABILITY 10 AMPS
603006
609003 INDICATORS-CAUTION AND WARNING 50010 11132
609003
609003 A REDUNDANT ISOLATED SYSTEM WHICH MONITORS AND GIVES WARNING OF
609003 SUBSYSTEM AND EXPERIMENT EQUIPMENT FAILURE AND/OR MALFUNCTION.
609003
609003 INDICATORS 40
609003 MASTER ALARM MEMORY
609003 POWER/TEST CONTROLS
609003
424000 RECORDER-TAPE 50011 11111
424000
424000 RECORDING OF SELECTED PARAMETERS DISPLAYED ON THE CONTROL AND
424000 DISPLAY CONSOLE.
424000
424000 BANDWIDTH DC TO 200 KHZ
424000 CHANNELS 14
424000 TIME MARKS EDGE TRACK
424000 RECORDED FORMAT DIGITAL
424000



Space Division
Rockwell International

HEWLETT PACKARD #7970R/C \$ 4600 111
TAPE FORMAT 800, 556, OR 200 CPI NRZI AND
1600 CPI PHASE-ENCODED
CHANNELS 7 OR 9
TAPE SPEED 10 TO 45 IPS
TAPE 0.5 INCH, 1.5 MILS, IBM/ANSI
COMPATIBLE

SANGAMO ELECTRIC #SABER 111 \$26,500 111
TRACKS 14
TAPE WIDTH 1 INCH
TAPE SPEEDS 8 SELECTABLE SPEEDS FROM 15/16
TO 120 IPS
FREQUENCY RESPONSE 400 HZ TO 2.0 MHZ
RECORDING RATE 600 KBPS AT 120 IPS SERIAL MODE
WEIGHT 100 POUNDS

LEACH #MTR 7C00 311
TRACKS 12
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
TAPE LENGTH 9200 FEET
PACKING DENSITY 16.7 KB/I/T
SIGNAL/NOISE 22 DB
DATA CAPACITY 2.2E+10 BITS
BANDWIDTH 2 MB/SEC/T

BORG WARNER #PERT 121
TRACKS 30
TAPE SPEED UP TO 1000 IPS
TAPE WIDTH 1/2 INCH
TAPE LENGTH 2400 FEET
PACKING DENSITY 15 KB/I/T
SIGNAL/NOISE 24 DB
DATA CAPACITY 1.3E+9 BITS
BANDWIDTH 6-15 MB/SEC/T

AMPEX CORP #AR 17C0 \$28,000 112
TRACKS 14
TAPE SPEED 120 IPS
TAPE LENGTH 9200 FEET
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KB/I/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 6.2E+10
BANDWIDTH 2 MB/SEC/T

AMPEX CORP #AR 7CC \$29,100 511
TRACKS 14
TAPE SPEED 60 IPS
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KB/I/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 8E+10 BITS
BANDWIDTH 1 MB/SEC/T

AMPEX CORPORATION #AR-200 121
BANDWIDTH 100 HZ TO 3125 HZ THROUGH
300 HZ TO 250 KHZ
TAPE WIDTH 0.5 INCH
TAPE SPEED 1.875 TO 60 IPS
RECORDING TIME 8 MINUTES TO 4 HRS AND 16 MINS
FORMAT DIGITAL
TRACKS 8 DIGITAL, 7 ANALOG

DIGI-DATA CORPORATION #1600 \$7650 111
TAPE SPEED 25, 18.75, 12.5 IPS
TRACKS 7 OR 9 TRACK
DATA DENSITY 1600 CPI PHASE ENCODED
200, 556, 800 CPI NRZI
TAPE 0.5 INCH, 1.5 MIL, 1200 FEET

DIGI-DATA CORPORATION #17C0/PDP-11/7-9 \$5250 111
TRACK NRZI
TAPE SPEED 45, 37.5, 25, 18.75, 12.5 IPS
TAPE 0.5 INCH, 1.5 MIL, IBM/ANSI
COMPATIBLE, 10.5 INCH REEL
TRACKS 7 OR 9 TRACK
DATA DENSITY PHASE ENCODED

MCNEYMELL #96 \$17420
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 240 IPS
REEL SIZE 16IN
MAXIMUM BANDWIDTH(DIRECT) 2M HZ
TAPE WIDTH 1/2 IN

MCNEYMELL #5600 \$9730
PORTABLE TAPE RECORDER
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 60 IPS
MAXIMUM BANDWIDTH(DIRECT) 300KHZ
PACKING DENSITY UP TO 600PRI
WEIGHT 70LBS (32KG)
INPUT VOLTAGE 28VDC

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

HEWLETT PACKARD	#39550	610200	111
BANDWIDTH	300 KHZ		
CHANNELS	7		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3955C	614700	111
BANDWIDTH	300 KHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3950A-011	623800	111
BANDWIDTH	500 HZ TO 2 MHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3960A/13065A/13063A	64796	111
CONFIGURATION	RACK MOUNTED		
TAPE SPEED	15, 3 AND 3/4, 15/16 IPS		
CHANNELS	4		
RECORDING FORMAT	FM		
PASSBAND	5 KHZ		
S/N RATIO	48 DB		
TRACKS	30		
601002	ANNUNCIATORS-ADVISORY	SC012	21132
601002	DUAL BANKS MOUNTED ON CONTROL/DISPLAY CONSOLE PROVIDING VISUAL		
601002	ALERT CUE WHEN A LOW PRIORITY MALFUNCTION OCCURS IN ANY ONBOARD		
601002	EXPERIMENT OR OTHER DESIGNATED MODULE SUBSYSTEM.		
601002	ANNUNCIATOR-VISUAL	12 IN DUAL BANKS	
601002	DISPLAY TYPE	RED, WHITE, GREEN	
601002	ANNUNCIATOR-AUDIO	HORN	
601002			



ATMOSPHERIC AND SPACE PHYSICS

484000 BCCM POSITION TV DISPLAY A7801 11111
484000
484000 DISPLAY MEDIUM RESOLUTION TELEVISION OF EXPERIMENT BCCM POSITIONS
484000 BEFORE, DURING AND AFTER EACH EXPERIMENT REPETITION.
484000
484000 RESOLUTION 525 LINES
484000 FRAME RATE 30 FRAMES/SEC
484000 BANDWIDTH VHF RANGE
484000

COMU ELECTRONICS INC #CQF 4775 111
AP-05 801 1
BANDWIDTH 30 MHZ
HORIZONTAL LINES 1225
FIELDS 60 PER SEC

CONRAC CORPORATION BRQA 14/R \$1620 111
VIDEO AMPLIFIER
COMPOSITE 0.3 - 3.0 V P-P
NONCOMPOSITE 0.3 V P-P FOR 100 PCT CONTRAST
3.0 V P-P MAXIMUM BEFORE PREAMP
OVERLOAD
VIDEO INPUT IMPEDANCE 50 K OHM MIN SHUNTED BY MAX 10
PFD SIGNAL OF 20 V P-P; HUM
SUPPRESSION 40 DB 60-1000 HZ
VIDEO FREQUENCY RESPONSE
AT 50 PCT OF MODULATION FLAT +/- 1 DB TO 30 PHZ, -6 DB
AT 40 PHZ
AT 30 PCT OF DEPTH MOD FLAT +/- 1 DB TO 30 PHZ, -2 DB
AT 40 PHZ
PULSE RESPONSE LESS THAN 15 NS RISE/FALL TIME
AT 50 PCT DEPTH OF MODULATION
EXTERNAL SYNC COMPOSITE OR SEPARATE H/V SYNC;
0.5 V TO 10 V P-P, GREATER THAN
5K OHM INPUT IMPEDANCE
LINEARITY +/- 1 PCT PICTURE HEIGHT
RETRACE TIME
HORIZONTAL 5 MICROSEC MAX
VERTICAL 800 MICROSEC MAX
SCANNING FREQUENCY
HORIZONTAL 15 - 40 KHZ
VERTICAL 15 - 60 FIELDS/SECOND
COMU INC. ELECTRONICS DIV #9600 \$5260
PICTURE SIZE 7-3/16" X 5-3/4" INCHES
HORIZONTAL SCAN 945 LINES, 50 OR 60 FIELDS/SEC
873 LINES, 60 FIELDS/SEC
VIDEO BANDWIDTH 20 MHZ +/- 2 DB
POWER 25-28.5 VDC, 150 WATTS

484000 EXPERIMENT TV DISPLAY A7802 11111
484000
484000 DISPLAY HIGH RESOLUTION TELEVISION PICTURES OF AREAS OF INTEREST
484000 IN REAL TIME ON FROM VIDEO RECORDERS AT VARIOUS SPEEDS.
484000
484000 RESOLUTION 525 LINES
484000 FRAME RATE 30 FPS
484000 BANDWIDTH UHF RANGE
484000

COMU ELECTRONICS INC. #CQF 4775 111
BANDWIDTH 30 MHZ
HORIZONTAL LINES 1225
FIELDS 60/SEC.

CONRAC CORPORATION BRQA 14/R \$1620 111
VIDEO AMPLIFIER
COMPOSITE 0.3 - 3.0 V P-P
NONCOMPOSITE 0.3 V P-P FOR 100 PCT CONTRAST
3.0 V P-P MAXIMUM BEFORE PREAMP
OVERLOAD
VIDEO INPUT IMPEDANCE 50 K OHM MIN SHUNTED BY MAX 10
PFD SIGNAL OF 20 V P-P; HUM
SUPPRESSION 40 DB 60-1000 HZ
VIDEO FREQUENCY RESPONSE
AT 50 PCT OF MODULATION FLAT +/- 1 DB TO 30 PHZ, -6 DB
AT 40 PHZ
AT 30 PCT OF DEPTH MOD FLAT +/- 1 DB TO 30 PHZ, -2 DB
AT 40 PHZ
PULSE RESPONSE LESS THAN 15 NS RISE/FALL TIME
AT 50 PCT DEPTH OF MODULATION
EXTERNAL SYNC COMPOSITE OR SEPARATE H/V SYNC;
0.5 V TO 10 V P-P, GREATER THAN
5K OHM INPUT IMPEDANCE
LINEARITY +/- 1 PCT PICTURE HEIGHT
RETRACE TIME
HORIZONTAL 5 MICROSEC MAX
VERTICAL 800 MICROSEC MAX
SCANNING FREQUENCY
HORIZONTAL 15 - 40 KHZ
VERTICAL 15 - 60 FIELDS/SECOND

COMU INC. ELECTRONICS DIV #9600 \$5260
PICTURE SIZE 7-3/16" X 5-3/4" INCHES
HORIZONTAL SCAN 945 LINES, 50 OR 60 FIELDS/SEC
873 LINES, 60 FIELDS/SEC
VIDEO BANDWIDTH 20 MHZ +/- 2 DB
POWER 25-28.5 VDC, 150 WATTS


ORIGINAL PAGE IS
OF POOR QUALITY

179500 SPECTRUM ANALYZER 41111
179500 MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE RATIO AND
179500 MULTIPLE RATIO OF VARIOUS FREQUENCY SOURCES.
179500
179500 BANDWIDTH 20 HZ TO 300 KHZ
179500 SENSITIVITY 20 NANOVOLTS
179500 ACCURACY +/- 0.6 DB
179500
TEKTRONIX INC #15 (PLUG-IN UNIT) 61300 111
FREQUENCY RANGE 50 HZ TO 1 MHZ
DEFLECTION FACTOR 1 MV/CM TO 100 V/CM
HEWLETT PACKARD #1417/8552P/8556A 6044 111
SELECTABLE FREQUENCY RANGES 1-20 KHZ, 0-300 KHZ
FREQUENCY SPECTRUM 20 HZ TO 300 KHZ
SENSITIVITY 20 NANOVOLT (-140 DBM TO 0 DBM)
ABSOLUTE AMPLITUDE ACCURACY +/- 0.6 DB
DISPLAY 4 DIV LINEAR, 70 DB LOGARITHMIC
FIDELAL SCIENTIFIC #0410 61500 111
SINGER #55R-90-1P 60400 111
FREQUENCY RANGE 100 HZ TO 400 KHZ
SENSITIVITY +/- 100
179500 MULTI-CHANNEL ANALYZER AP804 11111
179500
179500 MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE RATIO AND
179500 MULTIPLE RATIO OF MULTIPLE FREQUENCY SOURCES SIMULTANEOUSLY.
179500
179500 BANDWIDTH 0 TO 10 MHZ
179500 TIME BASE 200 MHZ CLOCK
179500 MEMORY SIZE 104K, 407K, 8192 CHANNELS
179500
HEWLETT PACKARD #5401B 62425 111
BUILT UP OF 541-100AH OSCILLOSCOPE WITH 54316 DISPLAY PLUG-IN
4277H DIGITAL PROCESSOR, AND 5416A ANALOG TO DIGITAL CONVERTER
IN 5410A POWER SUPPLY INTERFACE
MOUNT BENCH
ANALYZER DATA # NO 100 61700
DATA ACQUISITION MODES PHA, MCS, OR LIST
BANDWIDTH 14 MHZ
STORAGE 12 BIT DATA WORDS
DISPLAY CRT, 6X10 CM
NUMBER OF CHANNELS 1024, 2048, 4096
COUNT CAPACITY PER CHANNEL 20480 AT 1 + 4 INTERNAL
FLAG BITS
DIMENSIONS 19X10.25X23.5 INCHES
WEIGHT 40 LBS
179500 WAV ANALYZER AP805 11111
179500
179500 DETECTION OF SIGNAL AMPLITUDE AND FREQUENCY INFORMATION
179500
179500 BANDWIDTH 10 TO 800 KHZ
179500 FREQUENCY RESOLUTION 10 HZ
179500 SENSITIVITY 1 MV TO 30 V FS
179500 CHANNELS 2 MINIMUM
179500
BOWNE AND SCHWARZ #ZF (100, #831.52) 618495 111
FREQUENCY RANGE 6 KHZ TO 1.3 MHZ
SENSITIVITY < 0.5 MICRIVOLT
TEKTRONIX INC #15 (PLUG-IN UNIT) 61300 111
FREQUENCY RANGE 50 HZ TO 1 MHZ
DEFLECTION FACTOR 1 MV/CM TO 100 V/CM
HEWLETT PACKARD #7590A/3594A 64920 111
FREQUENCY RANGE 20 HZ TO 620 KHZ
FREQUENCY ACCURACY +/- (1 HZ + TIME BASE ACCURACY)
FOR 20 HZ TO 62 KHZ BAND
FOR 110 HZ + TIME BASE ACCURACY
FOR 500 HZ TO 620 KHZ BAND
FREQUENCY RESOLUTION 1 HZ/MINOR DIV (20 HZ TO 620 KHZ)
10 HZ/MINOR DIV (500 HZ TO 620 KHZ)
AMPLITUDE RANGE 7 MICROVOLT TO 30 VFS, 16 RANGES
BANDWIDTH SPECIFIED 10, 100, 1000, 3100 HZ
SWEEP RATES 1, 10, 100, 1000, 3100 HZ/S
MAX SWEEP TIME 820 S +/- 15%
179500 WAVE ANALYZER AP806 11111
179500
179500 DETACHING, SEPARATE AND ANALYZE VARIOUS FREQUENCY COMPONENTS OF
179500 INPUT SIGNALS (IE FUNDAMENTAL, HARMONICS, INTERMODULATION PRO-
179500 DUCTS, ETC.).
179500
179500 BANDWIDTH 1 KHZ TO 1.5 MHZ
179500 FREQUENCY RESOLUTION +/- 1.0 %
179500 SENSITIVITY 10 MV TO 100 V
179500 CHANNELS 2 MINIMUM
179500
HEWLETT PACKARD #312A 64275 111
FREQUENCY RANGE 1 KHZ TO 1.5 MHZ IN 14 OVERLAP-
PING BANDS, 200 KHZ OVERLAP
FREQUENCY ACCURACY +/- (10 HZ + TIME BASE ACCURACY)
AMPLITUDE RANGE -67 DBM TO 23 DBM FS 50-150 OHM
-107 DBM TO 13 DBM 600 OHM



219900 PATCH PANEL, COAXIAL APR07 81111
219900
219900 MULTIPLE INPUT PANEL CAPABLE OF CHANNELING HIGH FREQUENCY RF
219900 ENERGY TO APPROPRIATE INSTRUMENTS WITH MINIMUM LOSS AND INTER-
219900 FERENCE.
219900
TROMPETER ELEC. #J51-52 \$75 111
VARIOUS SIZE PATCH PANELS, JACKS AND PLUGS ARE AVAILABLE
FROM TROMPETER. CATALOG SHOULD BE CONSULTED FOR SPECIFIC
APPLICATION. REVIEW INDICATES ITEMS COULD BE USED DIRECTLY
IN SPACE
179500 FREQUENCY COUNTER APR08 11111
179500
179500 DIRECT MEASUREMENT OF FREQUENCY AND/OR SIGNAL REPETITION RATE.
179500
179500 BANDWIDTH 0 TO 35 MHz
179500 SENSITIVITY 10 MV RMS
179500 TIME BASE 100KHZ TO 10 MHZ
179500
GENERAL RADIO #1191-B \$1495 111
FREQUENCY DC TO 35 MHz
COUNTING GATE TIMES 1 MICROSEC TO 10 SEC
ACCURACY +/- 1 COUNT +/- TIME-BASE ACCUR
SENSITIVITY 10 MV RMS SINE WAVE, 30 MV F-P
PULSE DECREASING ABOVE 20 MHZ
TO APPROX 100 MV RMS AT 30 PHZ
DISPLAY 8 DIGITS
HEWLETT PACKARD #5323A \$1650 111
JOHN FLUKE MFG CO INC #1552A \$695 111
FREQUENCY RANGE DC TO 80 MHZ
DISPLAY 7 DIGIT LED (8 OR 9 OPTIONAL)
SENSITIVITY 50 MV RMS DC TO 50 MHZ INCREAS-
ING TO 75 MV AT 80 MHZ
PERIOD RANGE DC TO 10 MHZ (DC COUPLED)
5 HZ TO 10 MHZ (AC COUPLED)
TIME INTERVAL RANGE 0.1 MICROSEC TO 1E+6 SEC
104020 AUTOMATIC DISPLAY GENERATOR APR09 11111
104020
104020 PROVIDES 3D DISPLAY
OPTICS TECHNOLOGY INC #215 \$1995 111
POWER LEVEL 0.3 TO 1.0MW; 1.0 TO 2.0MW
EXPOSURE TIME 45 TO 180 SEC; 15 TO 60 SEC
MOUNT RENCH
FILM PLATE SIZE 4 X 5IN
HOLOGRAM TYPES GABOR AND FRESNEL
JEDON ENGR ASSOC #HS-2 \$8700
ACTIVE VIBRATION ISOLATION
LASER 20MW HE NE
SYSTEM COMPONENTS REAL TIME PLATE HOLDER, SHUTTER
BEAM STEEPER, REFERENCE MIRRORS
SPATIAL FILTERS, VARIABLE BEAM
SPLITTER, ENCLOSURE
JEDON ENGR ASSOC #AP-100 AUTO \$6500
PHOTO PROCESSOR
PROVIDES CHEMICAL FILLING AND CIRCULATING, AND DRAINING TO
PROCESS PHOTOGRAPHIC PLATE
104020
055050 CAMERA STILL APR10 11111
055050
055050 RECORD EXPERIMENT PHENOMENA
055050
055050 FILM SIZE 35MM
055050 SHUTTER SPEED 1 SEC TO 1/500 SEC
055050
HONEYWELL INC #SPOTOMATIC F +7191 \$639 111
LENS 55 MM, F/1.8 TAKUMAP
35 MM, F/2.0 TO F/16 TAKUMAP
SHUTTER SPEEDS 1/1000 SEC TO TIME
VIEWING THRU LENS
FILM SIZE 35 MM
430000 X-Y RECORDER APR11 21111
430000
430000 REAL-TIME ACCURATE REPRODUCTION OF SPECTRAL DATA SOURCES PLOTTED
430000 ON CARTESIAN COORDINATES.
430000
430000 CHANNELS 1
430000 ACCURACY +/- 0.2 % FS
430000
RHOE AND SCHWARTZ #25K (100,1950,02) \$10500 111
WRITING SPEED 1.50 M/SEC
DEFLECTION FACTOR 10 MICROVOLT/CM TO 10 V/CM
WRITING METHOD ELECTRIC PEN
HEWLETT PACKARD #7035A \$985 111
INPUT RANGE 1, 10, 100 MV/IN; 1 AND 10 V/IN
WRITING METHOD ELECTRIC PEN
ACCURACY +/- 0.2 % FS
LINEARITY +/- 0.1 % FS
425500 RECORDER - STRIP CHART APR12 21111
425500
425500 RECORD HARD COPY DATA FROM LOW FREQUENCY ANALOG INVESTIGATIONS
425500

ORIGINAL PAGE IS
OF POOR QUALITY



BD 74-BA-0047-2

HEWLETT PACKARD	#7100RA WITH OPTION 019 AND 023. 17500A PLUG-IN	\$2295	111
CHANNELS	2 CHANNELS PLUS EVENT MARKER ON BOTH SIDES		
WRITING METHOD	ELECTRIC		
ACCURACY	+/- 0.2 % FS		
VOLTAGE SPAN	5 MV TO 100 V		
HEWLETT PACKARD	#1702A/8801A	\$2750	111
FREQUENCY RESPONSE	DC TO 125 HZ		
WRITING METHOD	THERMAL		
SENSITIVITY	5 MV/DIV		
CHANNELS	2 CHANNELS TIMER/MARKER		
VOLTAGE RANGE	5 MV TO 5 V PER DIV		
ACCURACY	+/- 1 %		
HONEYWELL	#150RA	\$3300	111
DATA CHANNELS	24		
EVENT CHANNELS	4		
FREQUENCY RESPONSE	DC TO 25KHZ		
PAPER SPEED	0.1 TO 120 IPS		
MOUNT	RACK		
PAPER WIDTH	11IN		
HEWLETT PACKARD	#770RA/RH02A	\$8975	
CONFIGURATION	RACK MOUNTED		
CHART SPEED	0.25, 0.5, 1, 2.5, 5, 10, 25, 50, 100 MM/SEC		
FREQUENCY RESPONSE	150 HZ		
CHANNELS	8		
WRITING METHOD	THERMAL		
INPUT RANGES	1, 2, 5, 10, 20, 50, 100, 200, 500, 1000 MV/DIV		
ACCURACY	+/- 1 %		
HONEYWELL	#1856	\$	
RESOLUTION	200 ELEMENTS/INCH (1000/SWEEP)		
RECORDING RATE	RE+06 BITS/SEC		
CRT	FIBER OPTICS		
SPOT SIZE	0.005 INCH		
SWEEP SPEED	20,000 LINES/SEC		
HONEYWELL	#1858	\$	
SIZE-RACK MOUNTED	4.75H X 19.1W X 21D		
WEIGHT (W/O ACCESSORIES)	65 LBS.		
ACCESSORIES	6 LBS.		
CHANNELS	18		
RECORDING STYLUS	INERTIALESS-NO MOVING ELEMENT		
RESPONSE	DC TO 5000 HZ		
INPUT SENSITIVITY	100 MICROVOLT - 300 VOLTS		
ACCURACY	0.1 %		
PAPER SPEED	0.1 TO 120 IPS		
POWER	120/240 V, 50-60 HZ		
HONEYWELL	#1912	\$7700	111
PAPER WIDTH	12IN		
CHANNELS	24		
FREQUENCY RESPONSE	DC TO 25KHZ		
PAPER SPEED	200IPS		
SPEED VARIATION	15 FORWARD, 12 REVERSE SPEEDS		
MOUNT	RACK		
BECKMAN INSTRUMENTS	BR	\$10759	111
BANDWIDTH	DC TO 130HZ		
CHANNELS	5000 HZ WITH BR-500TAPE REC		
REC FORMAT	1 TO 8 UP TO 24 SPEC CROEP		
CHART SPEED	THERMAL RECTILINEAR		
EVENT MARKERS	0.1 TO 25 CM/SEC OTHERS AVAIL		
DIMENSIONS	2 (OPTION)		
POWER	69.5X21.5X22.5 INCHES		
	120V 50/60 HZ 250W		
424000 TAPE RECORDER - ANALOG		APR14	11111
424000 RECORD AND PLAYBACK OF ANALOG DATA IN SUPPORT OF ATMOSPHERIC AND SPACE PLASMA PHYSICS SORTIE LABORATORY INVESTIGATIONS.			
424000 BANDWIDTH	2 MHZ		
424000 CHANNELS	14		
424000 RECORDING FORMAT	FM		
424000 AMPLEX CORPORATION	#AR-20C		121
BANDWIDTH	100 HZ TO 3125 HZ THROUGH 300 HZ TO 250 KHZ		
TAPE WIDTH	0.5 INCH		
TAPE SPEED	1.875 TO 60 IPS		
RECORDING TIME	8 MINUTES TO 4 HRS AND 16 PINS		
FORMAT	DIGITAL		
TRACKS	8 DIGITAL, 7 ANALOG		
HEWLETT PACKARD	#757CB/C	\$ 4600	111
TAPE FORMAT	800, 556, OR 200 CPI NRZI AND 1600 CPI PHASE-ENCODED		
CHANNELS	7 OR 9		
TAPE SPEED	10 TO 45 IPS		
TAPE	0.5 INCH, 1.5 MILS, IBM/ANSI COMPATIBLE		

ORIGINAL PAGE IS
OF POOR QUALITY



DIGI-DATA CORPORATION		#1700/PDP-11/7-9	95290	111
T-ACKS		4		
TAPE SPEED		45, 37.5, 24, 18.75, 12.5 IPS		
TAPE		0.5 INCH, 1.5 MIL, INH/ANSI		
TRACKS		7 OR 9		
DATA DENSITY		PHASE ENCODED		
RONG WANNER		APERT		121
TRACKS		30		
SIGNAL/NOISE		24 DB		
TAPE SPEED		UP TO 1000 IPS		
TAPE WIDTH		1/2 INCH		
TAPE LENGTH		2400 FEET		
PACKING DENSITY		15 KB/1/1		
DATA CAPACITY		1.3E+9		
BANDWIDTH		5-15 MB/SEC/T		
LIALI		RMTH TCOL		311
T-ACKS		12		
TAPE SPEED		120 IPS		
TAPE WIDTH		1 INCH		
TAPE LENGTH		9200 FEET		
PACKING DENSITY		16.7 KB/1/1		
SIGNAL/NOISE		22 DB		
DATA CAPACITY		2.2E+10		
BANDWIDTH		2 MB/SEC/T		
KANDAM ELECTRIC		#SARF	926,900	111
T-ACKS		14		
TAPE WIDTH		1 INCH		
TAPE SPEED		A SELECTABLE SPEEDS FROM 15/16 TO 120 IPS		
FREQUENCY RESPONSE		400 HZ TO 2.0 MHZ		
RECORDING RATE		400 KBPS AT 120 IPS SERIAL MODE		
WEIGHT		100 POUNDS		
HEWLETT PACKARD		#3C95D	910200	111
BANDWIDTH		300 KHZ		
CHANNELS		7		
RECORDING FORMAT		DIRECT OR FM		
HEWLETT PACKARD		#3C95C	914700	111
BANDWIDTH		300 KHZ		
CHANNELS		14		
RECORDING FORMAT		DIRECT OR FM		
HEWLETT PACKARD		#3C90A-D11	923900	111
BANDWIDTH		500 HZ TO 2 MHZ		
CHANNELS		14		
RECORDING FORMAT		DIRECT OR FM		
HEWLETT PACKARD		#3560A/13065A/13065A	94704	111
CONFIGURATION		RACK MOUNTED		
TAPE SPEED		15, 3 AND 3/4, 15/16 IPS		
CHANNELS		4		
RECORDING FORMAT		FM		
PASSBAND		9 KHZ		
S/N RATIO		48 DB		
TRACKS		30		
MCNEYHILL		#5800	94730	
PORTABLE TAPE RECORDER				
CHANNELS		7		
SELECTABLE TAPE SPEED RANGE		15/16 TO 40IPS		
MAXIMUM BANDWIDTH(DIRECT)		300KHZ		
PACKING DENSITY		UP TO 400CPI		
WEIGHT		70LBS (32KG)		
INPUT VOLTAGE		28VDC		
MCNEYHILL		#CE	917420	
CHANNELS		7		
SELECTABLE TAPE SPEED RANGE		15/16 TO 240IPS		
WEIGHT		16IN		
MAXIMUM BANDWIDTH(DIRECT)		2M HZ		
TAPE WIDTH		1/2 IN		
DIGI-DATA CORPORATION		#1600	92650	111
TAPE SPEED		24, 18.75, 12.5 IPS		
TRACKS		7 OR 9 TRACK		
DATA DENSITY		1600 CPI PHASE ENCODED		
TAPE		200, 556, 800 CPI NRZI		
TAPE		0.5 INCH, 1.5 MIL, 1200 FEET		
AMPEX CORP		BAR 70C	929,180	511
TRACKS		14		
TAPE SPEED		60 IPS		
TAPE WIDTH		1 INCH		
RECORDING MODE		DIRECT		
PACKING DENSITY		20 KB/1/1		
SIGNAL/NOISE		20 DB		
DATA CAPACITY		4E+10 BITS		
BANDWIDTH		1 MB/SEC/T		
949500 CIRCLE KSCIPF			APR15	21111
949500				
MONITOR MEASURE AND MAINTAIN OPERATING ELECTRONIC EQUIPMENT				
949500				
949500		BANDWIDTH	500 KHZ	
949500		CHANNELS	2	
949500		DEFLECTION FACTOR	2 MV	
949500		PHOTOGRAPHIC CAPABILITY	YES	
949500				



HEWLETT-PACKARD	#172A	\$1500	111
BANDWIDTH	DC TO 500KHZ		
SENSITIVITY	100 MICROVOLTS/CM		
INDEPENDENT BEAMS	2		
WEIGHT	43 POUNDS		
TEKTRONIX INC	#7904/7A19 PLUG-IN	\$6050	111
	/7B70 PLUG-IN		
BANDWIDTH	DC TO 500 MHZ		
CHANNELS	2		
DEFLECTION FACTOR	10 MV/DIV TO 1 V/DIV		
TIMEBASE	2 NS/DIV TO 5 S/DIV		
HEWLETT/PACKARD	BAN/USH-281A	\$3100	111
MEETS MIL-SPEC REQUIREMENTS			
BANDWIDTH	DC TO 5MHZ		
SENSITIVITY	0.1V/DIV TO 1.0V/DIV		
CHANNELS	2		
HEWLETT PACKARD	#1800/1811A	\$2750	111
BANDWIDTH	4 OR 18 GHZ		
DEFLECTION FACTOR	2 MV		
CHANNELS	2		
PHOTOGRAPHIC CAPABILITY	YES		
TEKTRONIX	#4R5	\$4200	111
RANGE, CALIBRATED	5MV/DIV TO 5V/DIV IN 10 STEPS		
UNCALIBRATED	CONT VAR TO 12.5V/DIV		
CHANNELS	1, ALTERNATE, CHOPPED, ADDED, X-Y		
	CHANNEL 2 (UP OR INVERTED)		
DIMENSIONS	16.7 X 52.4 X 30.5 CM		
WEIGHT	6.5 KG		
POWER	115/230V 4R-440HZ 60W		
RESPONSE TIME	< 1NS		
BANDWIDTH	350MHZ		
TEKTRONIX INC	#P4A5	\$1800	111
BANDWIDTH	DC TO 100 MHZ		
CHANNELS	2		
DEFLECTION FACTOR	5 MV/DIV TO 5 V/DIV		
TIME BASE	0.01 MICROSEC/DIV TO 0.5 S/DIV		
TEKTRONIX INC	#R561R/3A6 PLUG-14/	\$1840	111
	3A4 PLUG-IN		
BANDWIDTH	DC TO 10 MHZ		
CHANNELS	2		
DEFLECTION FACTOR	10 MV/DIV TO 10 V/DIV		
TIME BASE	50 NS/DIV TO 5 S/DIV		
RISETIME	35 NS		
055050 CAMERA, STILL		APB16	21111
055050			
055050	POLAROID CAMERA TO RECORD TRACES ON OSCILLOSCOPE		
055050			
055050	FILM FORMAT	POLAROID 107 BLACK/WHITE	
055050		3.25X4.25IN (82.55X102.95MM)	
055050	LENS	50MM F3.5 TO F22	
055050	SHUTTER SPEEDS	8 TO 1/60 SEC	
055050	FIELD OF VIEW	40 DEGREES	
055050			
TEKTRONIX INC	#C-12	\$590	111
LENS	75 MM		
STOP	F/1.9 TO F/16		
MAGNIFICATION	0.95		
LENS SPEED	1 TO 1/100 SEC MECH		
	4 TO 1/60 SEC ELECT		
FILM TYPE	POLAROID		
TEKTRONIX	#C30	\$525	
MAXIMUM RELATIVE APERTURE	F/1.9		
MAGNIFICATION	0.7 TO 1.5		
RELATIVE SPEED	1.0		
FIELD OF VIEW	3.15X3.93IN (8X10CM)		
FILM TYPE	POLAROID		
HEWLETT-PACKARD	#1C4	\$420	
MAXIMUM RELATIVE APERTURE	F/3.5		
MAGNIFICATION	1:0.85		
SPEED	8 TO 1/60 SEC		
LENS	75MM		
FILM TYPE	POLAROID		
HEWLETT-PACKARD	#1C5A	\$1025	
MAXIMUM RELATIVE APERTURE	F/1.3		
MAGNIFICATION	1:0.5		
SPEED	8 TO 1/30 SEC		
LENS	80MM		
FILM TYPE	POLAROID		
424000 TAPE RECORDER - DIGITAL		APP17	21111
424000			
424000	RECORDING AND PLAYBACK OF DIGITAL DATA COMPATIBLE WITH ANY ON-		
424000	BOARD COMPUTING EQUIPMENT.		
424000			
424000	BANDWIDTH	2 MHZ	
424000	CHANNELS	8	
424000	SIGNAL	0 TO 5 VDC DIGITAL	
424000			

ORIGINAL PAGE IS
OF POOR QUALITY

ANALOG CORPORATION	AN-200		121
PANWIDTH	100 HZ TO 1120 HZ THROUGH		
	100 HZ TO 150 HZ		
TAPE WIDTH	0.5 INCH		
TAPE SPEED	1.875 TO 60 IPS		
RECORDING TIME	8 MINUTES TO 4 HRS AND 16 MINS		
FORMAT	DIGITAL		
TRACKS	8 DIGITAL, 7 ANALOG		
DIGI-DATA CORPORATION	#1700/DP-11/7-9	85250	111
	TRACK 14/21		
TAPE SPEED	45, 27.5, 25, 18.75, 12.5 IPS		
TAPE	0.5 INCH, 1.4 MIL, 18P/ANSI		
TRACKS	COMPATIBLE, 10.5 INCH REEL		
DATA DENSITY	7.89 TDS/IN		
	PHASE ENCODED		
ANALOG CORP	AN-1700	84,000	111
TRACKS	14		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
RECORDING MODE	DIRECT		
PACKING DENSITY	20 KB/1/1		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	8.2E+10		
PANWIDTH	2 MB/SEC/1		
TAPE LENGTH	4200 FEET		
SANGAMU ELECTRIC	ASABER 111	870,500	111
TRACKS	14		
TAPE WIDTH	1 INCH		
TAPE SPEEDS	8 SELECTABLE SPEEDS FROM 15/16		
	TO 120 IPS		
FREQUENCY RESPONSE	400 HZ TO 17.0 KHZ		
RECORDING RATE	800 FRPS AT 120 IPS SERIAL MODE		
WEIGHT	100 POUNDS		
HEWLETT PACKARD	#75700	9,400	111
TAPE FORMAT	800, 55A, UP 200 CPI NRZI AND		
	1400 CPI PHASE-ENCODED		
CHANNELS	7 OR 9		
TAPE SPEED	10 TO 45 IPS		
TAPE	0.5 INCH, 1.5 MILS, 18P/ANSI		
	COMPATIBLE		
ROG WARNER	WHERT		121
TRACKS	30		
TAPE SPEED	UP TO 1000 IPS		
TAPE WIDTH	1/2 INCH		
TAPE LENGTH	2400 FEET		
PACKING DENSITY	15 KB/1/1		
SIGNAL/NOISE	24 DB		
DATA CAPACITY	1.3E+9		
PANWIDTH	1-15 MB/SEC/1		
LEACH	WTR 7000		311
TRACKS	12		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
TAPE LENGTH	9200 FEET		
PACKING DENSITY	16.7 KB/1/1		
SIGNAL/NOISE	22 DB		
DATA CAPACITY	2.2E+10		
PANWIDTH	2 MB/SEC/1		
HEWLETT PACKARD	#35950	810200	111
PANWIDTH	300 KHZ		
CHANNELS	7		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3555C	814700	111
PANWIDTH	300 KHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3401-011	823000	111
PANWIDTH	500 HZ TO 2 KHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3600A/13056/13063A	84796	111
CONFIGURATION	15, 3 AND 1/4, 15/16 IPS		
TAPE SPEED	4		
CHANNELS	FM		
RECORDING FORMAT	5 KHZ		
PASSBAND	48 DB		
S/N RATIO	30		
HEWLETT	#5800	89730	
PORTABLE TAPE RECORDER			
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 40IPS		
MAXIMUM PANWIDTH (DIRECT)	300KHZ		
PACKING DENSITY	UP TO 800P/1		
WEIGHT	70LBS (32KG)		
INPUT VOLTAGE	28VDC		
HEWLETT	#96	817420	
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 240IPS		
REEL SIZE	16IN		
MAXIMUM PANWIDTH (DIRECT)	24 KHZ		
TAPE WIDTH	1/2 IN		



DIGI-DATA CORPORATION	#160C	\$2653	111
TAPE SPEED	25, 18-75, 12.5 IPS		
TRACKS	7 OR 9 TRACK		
DATA DENSITY	1600 CPI PHASE ENCODED		
	200, 956, 800 CPI NRZI		
TAPE	0.5 INCH, 1.5 MIL, 1200 FEET		
APPEX CORP	#AP 70C	\$29,183	511
TRACKS	14		
TAPE SPEED	60 IPS		
TAPE WIDTH	1 INCH		
RECORD MODE	DIRECT		
PACKING DENSITY	20 KB/I/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	8F+1C BITS		
BANDWIDTH	1 MB/SEC/T		
102700 COMPUTER		4PE12	21111
102700			
102700	COMPUTER CAPABILITY TO SUPPORT ATMOSPHERIC AND SPACE PLASMA		
102700	PHYSICS INVESTIGATIONS.		
102700			
102700	WORD SIZE	16 BIT	
102700	I/O MEMORY TRANSFER	1E+06--16 BIT WORDS/SEC	
102700	CORE MEMORY	32E+03 WORDS	
102700	I/O CHANNELS	45	
102700			
DATA GENERAL CORP	#NOVA 800	\$18,450	111
WORD LENGTH	16 BIT		
CYCLE TIME	800 NANOSEC.		
CORE MEMORY	32K MAX		
DATA GENERAL CORP	#SUPERNOVA SC	\$5600	111
WORD LENGTH	16 BIT		
VARIAN DATA MACHINES	#520/I		111
MEMORY CYCLE TIME	1.5 MICROSEC.		
MEMORY	EXPANDABLE, 4096 BYTES (8 BITS) TO 32,768 BYTES		
REGISTERS	12		
OPERAND PRECISION	UP TO 32 BITS		
ADDRESS REFERENCE	8, 16, 24, OR 32 BIT LEVEL		
NOTE:	HAS FULL RANGE INTERFACE HARDWARE		
VARIAN DATA MACHINES	#P-620/I		111
BASIC COMMANDS	OVER 100		
ADDRESSING MODES	6		
MAX. WORDS	32,768		
WORD LENGTH	16 OR 18 BIT		
REGISTERS	9		
NOTE:	HAS FULL RANGE INTERFACE HARDWARE		
CLARY DATACOMP SYSTEMS INC	#404		111
SIMULTANEOUS TERMINALS	16		
INPUT	TELETYPEWRITER KEYBOARD		
ACCUMULATOR	64 BIT		
REGISTERS	16 BIT INDEX (2 EACH)		
WORD LENGTH	16, 32, 48, 64 BITS		
MEMORY CAPACITY	1024 16 BIT WORDS, 4096 16 BIT WORDS, OR ADDITIONAL 4096 16 BIT WORDS. TOTAL 65536 WORDS		
GENERAL AUTOMATION INC.	#1R-30		111
FUNCTION	SUPERVISE SMALLER COMPUTERS		
MEMORY	8K CORE		
CYCLE TIME	1.2 MICROSEC.		
GENERAL AUTOMATION INC.	#SPC-16		111
MEMORY	16K		
WORD LENGTH	16 BITS		
READ/WRITE CYCLE MEMORY TIME	800 TO 1440 NANOSEC.		
READ ONLY MEMORY	400 TO 720 NANOSEC.		
INPUT/OUTPUT TIME	1.6 TO 2.8 MICROSEC.		
DMA TRANSFER RATE	0.694E+06 TO 2.5E+06		
GENERAL AUTOMATION INC.	#SPC-12		111
MEMORY	4K TO 16K		
WORD LENGTH	8 BITS		
CYCLE TIME	2.1A MICROSEC		
STORED PROGRAM EXECUTION RATE	0.23E+06/SEC		
I/O TRANSFER RATE	0.44E+06/SEC		
REGISTERS	6 12 BIT REGISTERS		
ACCUMULATORS	4 12 BIT		
DATA GENERAL CORP	#NOVA 1210	\$11900	111
WORD LENGTH	16 BIT		
CYCLE TIME	1200 NANOSEC.		
CORE MEMORY	32K MAX.		
DIGITAL EQUIPMENT CORP	#PDP-8/E	\$22000	111
MEMORY	4096 CORE, EXPANDABLE TO 32,768 WORDS		
WORD	12 BIT		
CYCLE TIME	1.5 MICROSEC.		

ORIGINAL PAGE IS
OF POOR QUALITY

101340 KEYBOARD DISPLAY TERMINAL APR19 11111
101340 GENERAL PURPOSE INPUT KEYBOARD ALLOWING FOR PROGRAMING, DATA
101340 INPUT AND UPDATE, AND DATA MANIPULATION FOR ONBOARD ANALYSIS
101340

RESEARCH INC #0PS-012 121.700 111
COMPLETE DATA DISPLAY SYSTEM MADE UP OF CRT ALPHANUMERIC DISPLAY
CPC-POPR-F HINS-COMPUTER, #0301 OPERATIONS CONTROL CASCLE,
PRINTER/CARD READER, #012-4A ANALYSER AND #012-13 UNIDRIVER
CHANNELS 32
MEMORY SIZE 8K WORDS

RESEARCH INC #012-1301 11555 111
SCREEN FORMAT 12 LINES X 72 OR 80 CHARACTERS,
24 LINES X 40 CHARACTERS
CHARACTER FORMAT 5 X 7 DOT MATRIX
TRANSFER RATE 110 TO 2400 AUD, 13 OR 11 PIT
CHARACTERS
MODES HALF OR FULL DUPLEX-SWITCHABLE
LOCAL OR REMOTE

HENLETT PACKARD #240CA 13343 111

TEKTRONIC #4002A/021-00XX-00 15550 111
DISPLAY MEDIUM 11 INCH DIRECT-VIEW, BISTABLE
STORAGE CRT WITH REFRESHED
SCRATCH PAD AREA
DISPLAY AREA 0.3 INCHES HORIZONTAL X 0.1
INCHES VERTICAL
ALPHANUMERIC MODE
FORMAT 39 LINES OF 85 NORMAL OR ITALIC
CHARACTERS IN MAIN AREA, ONE
LINE OF 84 CHARACTERS IN
SCRATCH PAD AREA
CHARACTER SET 96 UPPER AND LOWER CASE PRINT-
ING CHARACTERS (ASCII CODE)
CHARACTER SIZE 70X90 MILS (CAN BE DOUBLE SIZE)
CHARACTER GENERATION 7X9 DOT MATRIX
CURSOR PULSATING 7X9 MATRIX
GRAPHIC MODES LINEAR INTERPOLATE, INCREMENTAL
PLOT, POINT PLOT, 1024X1024
ADDRESSABLE POINTS, 1024X761
VIEWABLE POINTS
GRAPHIC INPUT MODE 1024 (X), 761 (Y) POINTS, JOY-
STICK CONTROLLED, CROSS-HAIR
CURSOR

104180 STATUS PANEL APR20 21111
104180
104180 DISPLAY TERMINAL WITHOUT KEYBOARD PROVIDING STATUS OF EXPERIMENT
104180 FUNCTIONS
104180

HENLETT PACKARD #240CA 13343 111

TEKTRONIC #4002A/021-00XX-00 15550 111
DISPLAY MEDIUM 11 INCH DIRECT-VIEW, BISTABLE
STORAGE CRT WITH REFRESHED
SCRATCH PAD AREA
DISPLAY AREA 0.3 INCHES HORIZONTAL X 0.1
INCHES VERTICAL
ALPHANUMERIC MODE
FORMAT 39 LINES OF 85 NORMAL OR ITALIC
CHARACTERS IN MAIN AREA, ONE
LINE OF 84 CHARACTERS IN
SCRATCH PAD AREA
CHARACTER SET 96 UPPER AND LOWER CASE PRINT-
ING CHARACTERS (ASCII CODE)
CHARACTER SIZE 70X90 MILS (CAN BE DOUBLE SIZE)
CHARACTER GENERATION 7X9 DOT MATRIX
CURSOR PULSATING 7X9 MATRIX
GRAPHIC MODES LINEAR INTERPOLATE, INCREMENTAL
PLOT, POINT PLOT, 1024X1024
ADDRESSABLE POINTS, 1024X761
VIEWABLE POINTS
GRAPHIC INPUT MODE 1024 (X), 761 (Y) POINTS, JOY-
STICK CONTROLLED, CROSS-HAIR
CURSOR

104180 SPECIAL DATA ACQUISITION PANEL APR21 11124
104180
603002 BATTERY CONTROL AND MONITOR APR27 11114
603002
095034 CAMERA - CINE APR29 11111
095034

J & MAUER INC #017-01 14500 311
MODES 9, 16, 32, 64 FPS ELECT CCNT
SHUTTER ROTARY
SHUTTER SPEEDS 1/60.5 TO 1/500 SEC
FILM 16 MM MAGAZINE
LEN 25 MM, F/1.4
RESOLUTION 23A L/PM 74 AXIS, 150 L/PM AWAP
FILM CAPACITY 150 FT THIN BASE
CINEMA REULTIEU CORP #0-10R 112
CONFIGURATION BENCH/LOOSE EQUIP
OPERATION MODE VARIABLE FRAME RATE OR SINGLE
FRAME PULSE OPERATION UNDER
REMOTE CONTROL
SPEED 2 TO 64 FRAMES/SEC (VARIABLE)
FILM SIZE 16 MM
LENS C MOUNT LENS SYSTEM
FILM CAPACITY VARIOUS FILM MAGAZINES



498500 TIME CODE GENERATOR AND DISPLAY APR32 11111
498500
498500 ACCURATE FREQUENCY, TIME INTERVAL AND TIME KEEPING CAPABILITIES.
498500
498500 OUTPUTS..... 100 KHZ, 1 MHZ, 5 MHZ
498500 ACCURACY 5E-10 PARTS PER DAY
498500
CATATRON #3150 \$3220
SERIAL CODE CAPABILITY..... UP TO 5
PULSE RATES..... 1000,100,10 AND 1 PPS
ACCURACY..... 50 MICROSECS PER DAY
DATATRON #3030 \$7940 111
CHANNELS..... 32
TIMING RANGE..... 5 NANOSECS TO 5 MICROSECS
ACCURACY..... 5E-7 PARTS PER DAY
DATATRON #3350 \$1020 111
INPUT FREQUENCY..... 60KHZ
OUTPUTS..... VISUAL, BCD
DISPLAY..... HRS,MIN, SEC
603002 REMOTE SENSING PLATFORM GIMBAL CONTROL AP833 11114
603002
603002 XUV NORMAL INCIDENCE SPECTROMETER CONTROL AP834 11114
603002
603002 UV-VIS-NIR SCANNING SPECTROMETER CONTROL AP835 11114
603002
603002 HIGH RESOLUTION FOURIER SWIR SPECTROMETER CONTROL AP836 11114
603002
603002 CRYOGENIC IR FOURIER SPECTROMETER CONTROL AP837 11114
603002
603002 IR RADIOMETER CONTROL APR38 11114
603002
603002 FABRE-PEROT INTERFEROMETER CONTROL AP839 11114
603002
603002 UV,VIS DOCUMENTATION CAMERA CONTROL AP840 11114
603002
603002 ELECTROSTATIC ANALYZER CONTROL APR41 11134
603002
603002 MAGNETIC ANALYZER CONTROL APR42 11134
603002
603002 KEV-MEV PARTICLE DETECTOR CONTROL AP843 11134
603002
603002 TOTAL ENERGY DETECTOR CONTROL AP844 11134
603002
623010 LIDAR AND GIMBAL PLATFORM MONITOR APR45 11134
623010
620021 TRANSMITTER - 0.2 TO 2.0 MHZ AP901 11111
620021
620021 VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.
620021
620021 FREQUENCY BANDWIDTH 0.2 TO 2.0 MHZ
620021 TRANSMITTER POWER 10 KW
620021
AVCIN CORP/VECTOR DIV #T-102-TV \$3680 121
TYPE UHF VIDEO TRANSMITTER
FREQUENCY RESPONSE 10 HZ TO 6 MHZ +/- 1.5 DB
DEVIATION SENSITIVITY +/- 6 MHZ/VOLT RMS
RF POWER OUTPUT 2 WATTS (50 OHM LOAD)
620020 TRANSMITTER - 2.0 TO 20.0 MHZ AP902 11111
620020
620020 VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.
620020
620020 FREQUENCY BANDWIDTH 2.0 TO 20.0 MHZ
620020 TRANSMITTER POWER 10 KW
620020
620050 TRANSMITTER - 0.3 TO 200 KHZ AP903 11111
620050
620050 VARIABLE RADIO FREQUENCY GENERATION AND TRANSMISSION.
620050
620050 FREQUENCY BANDWIDTH 0.3 TO 200.0 KHZ
620050 TRANSMITTER POWER 1.0 KW
620050
AVCIN CO. T1C85 \$6100 121
WATTS
POWER OUTPUT RES BAND1
MODULATION TYPE 1W/50 OHM LOAD
FREQ RESPONSE FM
DC TO 500KHZ +/-1.5DB
620022 ELECTROSTATIC WAVE TRANSMITTER AP904 11134
620022
219900 COUPLER - ANTENNA (0.2 TO 2.0 MHZ) AP909 11111
219900
219900 ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.
219900
219900 FREQUENCY BANDWIDTH 0.2 TO 2.0 MHZ

ORIGINAL PAGE IS
OF POOR QUALITY

219900
219900 COUPLER - ANTENNA (2.0 TO 20.0 MHZ) AP910 11111
219900 ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.
219900 FREQUENCY BANDWIDTH 2.0 TO 20.0 MHZ
219900
219900 COUPLER - ANTENNA (0.3 TO 200 KHZ) AP911 11111
219900 ANTENNA TO TRANSMITTER INTERFACE ALLOWING FOR MINIMUM RF LOSS.
219900 FREQUENCY BANDWIDTH 0.3 TO 200 KHZ
219900
003300 BANDPASS FILTER - WAVE ANALYSIS AP912 21111
003300
003300 DETERMINATION OF FREQUENCY AND AMPLITUDE INFORMATION ABOUT THE
003300 CARRIER AND SIDEBANDS.
003300
003300 FREQUENCY RANGE 300 KHZ TO 20.0 MHZ
003300 ATTENUATION 75 DB
003300 ACCURACY +/- 2 %
003300
ROCKLAND SYSTEMS CORP #1042F-01 82074 111
FREQUENCY RANGE 0.01 HZ TO 111 KHZ
FREQUENCY ACCURACY +/- 2 %
PASS BAND GAIN 0 DB OR 20 DB
ATTENUATION SLOPE 48 DB/OCTAVE
KRCN-MITE #3342R 82100 111
FREQUENCY RANGE 0.001 TO 99.9 KHZ
FREQUENCY ACCURACY +/- 2 %
PASS BAND GAIN 0 DB OR 20 DB
ATTENUATION SLOPE 48 DB/OCT
MAX ATTENUATION 90 DB
KRCN-MITE #3343R 81950 111
FREQUENCY RESPONSE 0.01 TO 99.9 KHZ
FREQUENCY ACCURACY +/- 2 %
PASS BAND GAIN 0 DB OR 20 DB
ATTENUATION SLOPE 48 DB/OCT
MAX ATTENUATION 90 DB
GENERAL RADIO #1652 81095 111
FREQUENCY RANGE 4 HZ TO 60 KHZ
ATTENUATION RATE 30 DB/OCT
GAIN 0 DB OR -20 DB
KRCN-MITE #3103 (R) 8060
FREQUENCY RANGE 10 HZ TO 3 MHZ
ATTENUATION SLOPE 24 DB/OCT
179500 PULSE WAVE ANALYZER- C.W. MODES AP913 11111
179500
RMOE AND SCHWARTZ #EZF (10C.8831.52) 819.4P5 111
FREQUENCY RANGE 6 KHZ TO 1.3 MHZ
SENSITIVITY < 0.5 MICROVOLT
TEKTRONIX INC #11E (PLUG-IN UNIT) 81300 111
FREQUENCY RANGE 50 HZ TO 1 MHZ
DEFLECTION FACTOR 1 MV/CM TO 100 V/CM
HEWLETT PACKARD #3590A/3594A 84920 111
FREQUENCY RANGE 20 HZ TO 620 KHZ
FREQUENCY ACCURACY +/- (1 HZ + TIME BASE ACCURACY)
FOR 20 HZ TO 62 KHZ BAND
+/- (10 HZ + TIME BASE ACCURACY)
FOR 500 HZ TO 620 KHZ BAND
FREQUENCY RESOLUTION 1 HZ/MINOR DIV (20 HZ TC 620KHZ)
10 HZ/MINOR DIV (500 HZ-620KHZ)
AMPLITUDE RANGE 3 MICROVOLT TO 30 VFS.16 RANGES
BANDWIDTH SPECIFIED 10, 100, 1000, 3100 HZ
SWEEP RATES 1, 10, 100, 1000, 3100 HZ/S
MAX SWEEP TIME 620 S +/- 15 %
HEWLETT PACKARD #312A 84275 111
FREQUENCY RANGE 1 KHZ TO 19 MHZ IN 18 OVERLAP-
PING RANGES, 200 KHZ OVERLAP
FREQUENCY ACCURACY +/- (10 HZ + TIME-BASE ACCURACY)
AMPLITUDE RANGE -97 DBM TO 23 DBM FS 50-150 OHM
-107 DBM TO 13 DBM 600 OHM
219900 PATCH PANEL AP914 51111
219900
219900 INTERFACE CONNECTIONS OF VARIOUS RF RECEIVERS WITH SELECTED TEST
219900 EQUIPMENT FOR WAVE ANALYSIS.
219900
TROMPETER ELEC. #J51-52 875 111
VARIOUS SIZE PATCH PANELS, JACKS AND PLUGS ARE AVAILABLE
FROM TROMPETER. CATALOG SHOULD BE CONSULTED FOR SPECIFIC
APPLICATION. REVIEW INDICATES ITEMS COULD BE USED DIRECTLY
IN SPACE
006180 APPLIFIERS - WAVE ANALYSIS AP915 21111
006180
HEWLETT PACKARD #461A 9380 111
FREQUENCY RANGE 1 KHZ TO 150 MHZ
FREQUENCY RESPONSE +/- 1 DB INTO 50 OHM LOAD
GAIN AT 500 KHZ 40 DB +/- 0.5 DB OR 20 +/- 1 DB
MAX INPUT 1 V RMS OR 2 V P-P PLS
TENNELEC #TC 2C2/2C3 8345/495 111
GAIN 2.5 TO 1000X
ACN LINEARITY LESS THAN 0.05 %
CONFIGURATION NIM COMPATIBLE



ROWDE AND SCHWARZ	#ATN (100.0899.02)		111
FREQUENCY RANGE	30 HZ TO 20 KHZ		
INPUT VOLTAGE	0.25 - 2.5 VOLTS		
MAX POWER	50 WATTS		
RESEARCH INC	#R12-11-17	\$3200	111
CAPACITY	TO 33 PAIRS OF SIG WIRES (16GA)		
REFERENCE TEMP	HEAT SINK CONTROLLED TO +/- .08 DEG C AT 65 DEG C		
AMBIENT TEMP	0 TO 49 DEG C		
LEAKAGE RESISTANCE	1000 MEG OHM		
HEAT-UP TIME FROM 22 DEG C	50 MINUTES		
SETTLING TIME	8 MS TO FINAL VALUE		
HONEYWELL	#ACCU DATA 117	\$1340	111
MULTICHANNEL WIDEBAND AMPLIFIER			
CHANNELS	7		
MOUNT	BENCH		
HONEYWELL	#ACCU DATA 115	\$5370	111
HIGH VOLTAGE DC AMPLIFIER			
INPUT SIGNAL VOLTAGE	+/-50MV TO 1500V		
OUTPUT SIGNAL VOLTAGE	UP TO 2000V		
CHANNELS	6		
NEFF	#122-222	\$610	111
GAIN	1 TO 2500		
COMMON MODE REJECTION	140DB		
DRIFT	0.5MICROVOLT/0C		
STEP GAIN ACCURACY	+/-0.01%		
SWITCH SELECTABLE FILTER			
BANDWIDTH	100KHZ		
RACK/MODULE MOUNT			
NEFF	#126-621	\$320	111
AMPLIFIER MULTIPLEXER			
EMPLOYES FET SWITCH			
GAIN RANGE	0.2 TO 2500		
COMMON MODE REJECTION	120DB		
OUTPUT	+/-10V AT 10 MILLIAMPERES		
NEFF	#127-102		
PROGRAMMABLE GAIN WITH SWITCH SELECTABLE OUTPUT FILTER			
GAIN STEPS	1, 2, 5, 10, 20, 50, 100, 200, 500, 1000		
COMMON MODE REJECTION	140DB		
FET SWITCHING			
RACK/MODULE CONFIGURATION			
TENNELEC	#TC133	\$295	111
FET PREAMPLIFIER			
CHARGE SENSITIVITY	10E12 VOLTS/COULOMB		
INPUT	NEGATIVE INPUT PULSE POLARITY WITH INPUT CAPACITANCE >2000PF		
OUTPUT	POSITIVE OUTPUT PULSE POLARITY WITH OUTPUT IMPEDANCE 50 OHM		
STABILITY	0.5%/V		
101000	FREQUENCY SYNTHESIZER - WAVE ANALYSIS	AP916	21111
101000	TRANSLATION OF A STABLE FREQUENCY OF A PRECISE STANDARD TO ANY OF A SELECTED INVESTIGATION REQUIREMENT.		
101000	FREQUENCY BANDWIDTH	DC TO 13 MHZ IN 5 RANGES	
101000	FREQUENCY RESOLUTION	0.01 TO 10 KHZ	
101000	FREQUENCY STABILITY	1E-07 PARTS PER DAY	
GENERAL RADIO	#1163-AR7C	\$4990	111
OUTPUT FREQUENCY RANGE	30HZ TO 12MHZ		
SMALLEST STEP	1HZ		
ROWDE AND SCHWARZ	#SMNH (100.4471.42)	\$21585	111
FREQUENCY RANGE	0 TO 50 MHZ		
OUTPUT LEVEL	0.1 MICROVOLT TO 2.5 VOLTS		
HEWLETT PACKARD	#3320R WITH OPTION	\$3145	111
FREQUENCY RANGE	0.01 HZ TO 13 MHZ IN 7 RANGES		
RANGES	10 MHZ, 1000 KHZ, 100 KHZ, 10 KHZ, 1000 HZ, 100 HZ, 10 HZ		
FREQUENCY STABILITY	+/- 10 IN 1E+6 CF SETTING/YEAR		
FREQUENCY RESPONSE	+/- 0.05 DB 10 HZ TO 13 MHZ		
AMPLITUDE RANGE	+26.99 DBM TO -69.99 DBM INTO 50 OHMS		
NAVETEK	#162	\$845	111
BANDWIDTH	0.00003 HZ TO 30 MHZ		
OUTPUT	80 DB ATTENUATION/10 DB STEPS 10 V P-P INTO 50 OHMS		
DC OFFSET	+/- 5 VOLTS INTO 50 OHMS		
STABILITY	+/- 0.05%/1C MIN. +/- 0.25%/24 HOURS		
SINGER	#8600	\$4000-5000	111
FREQUENCY RANGE	100KHZ TO 10GHZ		
POWER OUTPUT	UP TO 100MW		
PLUG-INS	13 FOR VARIOUS OUTPUT RANGES PLUG-INS COST \$2 TO 3000 EACH		

ORIGINAL PAGE IS
OF POOR QUALITY

JOHN F. HART, JR. 00456-156H 913,450 311
 FREQUENCY RANGE DC TO 50 MHz
 INCREMENTS 0.01 MHz
 SPURIOUS OUTPUTS > 100 DB (NON-HARMONICS)
 > 30 DB (HARMONICS)
 SIGNAL-TO-PHASE NOISE RATIO > 60 DB
 FREQUENCY STABILITY 2E-9/24 HOURS AVAILABLE

603002	MAIN ROOM A CONTROL	AP017	11134
603002			
603002	PLATFORM ROOM A CONTROL	AP018	11134
603002			
603002	MAIN ROOM B CONTROL	AP019	11134
603002			
603002	ALIGNMENT TV CONTROL	AP020	11134
603002			
603002	GEMHALLED PLATFORM CONTROLS	AP021	11134
603002			
603002	SM ROOM CONTROL	AP022	11134
603002			
603002	ROOM A POWER SUPPLY & DATA SYS CONTROL	AP023	11134
603002			
603002	ROOM B TARGET CONTROL	AP024	11134
603002			
603002	ACCEL-DECEL CONTROL	AP025	11134
603002			
603002	DISCHARGE FILAMENT HEATER CONTROL	AP026	11134
603002			
603002	DISCHARGE POTENTIAL CONTROL	AP027	11134
603002			
603002	PULSE SEQUENCE & PULSE LENGTH CONTROL	AP028	11134
603002			
603002	GAS SELECTION & PRESSURE CONTROL	AP029	11134
603002			
603002	NEUTRALIZER EMISSION & BIAS CONTROL	AP030	11134
603002			
603002	CHARGE EXCHANGE CHANNEL ACTUATOR CONTROL	AP031	11134
603002			
603002	PEAK CURRENT MONITOR	AP032	11134
603002			
603002	ELECTRON BEAM VOLTAGE & CURRENT HEATER CONTROL	AP033	11134
603002			
603002	ELECTRON BEAM PULSE LENGTH & MAGNITUDE CONTROL	AP034	11134
603002			
603002	ELECTRON BEAM EXPANSION LENS CONTROL	AP035	11134
603002			
603002	ELECTRON BEAM CURRENT MONITOR	AP036	11134
603002			
603002	PHASE, ANGLE, CURRENT MONITOR	AP037	11134
603002			
603002	WFO ARC VOLTAGE LEVEL CONTROL	AP038	11134
603002			
603002	WFO ARC PULSE CURRENT DURATION CONTROL	AP039	11134
603002			
603002	WFO ARC PULSE SEQUENCE	AP040	11134
603002			
603002	WFO ARC BEAM CURRENT MONITOR	AP041	11134
603002			
603002	GAS SELECTION & PRESSURE CONTROL	AP042	11134
603002			
603002	SPHERICAL ION PROBE CONTROL	AP043	11134
603002			
603002	CYLINDRICAL ION PROBE CONTROL	AP044	11134
603002			
603002	PLANAR ELECTRON PROBE CONTROL	AP045	11134
603002			
603002	SEGMENTED PLANAR PROBE CONTROL	AP050	11134
603002			
603002	DEP. SAT. INSTRUMENT CONTROL & HOUSEKEEPING	AP050	11134
603002			
603002	DEP. SAT. EJECTION MECHANISM CONTROL	AP060	11134
603002			
603002	PHOTOMETER HIGH VOLTAGE SUPPLY CONTROL	AP061	11134
603002			
603002	PHOTOMETER AMPLIFIERS CONTROL	AP062	11134
603002			
603002	TV SYS CONTROL - IMAGE INTENSIFIER CONTROL	AP063	11134
603002			
603002	CANNISTER EJECTION CONTROL	AP064	11134
603002			
603002	PROJECTILE CAMERA CONTROL	AP065	11134
603002			
603002	CANNISTER MONITOR	AP066	11134
603002			
603002	SHAPED CHARGE EJECTION CONTROL	AP067	11134
603002			
603002	SHAPED CHARGE MONITOR	AP068	11134
603002			
603002	BALLOON EJECTION CONTROL	AP069	11134
603002			
603002	GAS CONTROL SYSTEM	AP070	11134
603002			

EARTH OBSERVATIONS

204200 HEAT EXCHANGER EC001 11121
204200 EXCHANGER TO TRANSFER HEAT FROM THE CLOUD PHYSICS LAB TO THE
204200 SPACELAB HEAT REJECTION SUBSYSTEM
204200 HEAT TRANSFER RATE.....195 BTU/HR (54 WATTS THERMAL)
204200 SPACELAB COOLANT INLET TEMP...43F (6C)
204200 CPL COOLANT OUTLET TEMP.....45F (7C)
205000 HEAT EXCHANGER EC001 11121
205000 PROVIDE HOT FLUID TO CPL
205000 CPL COOLANT OUTLET TEMP.....100F (38C)
205000 HEAT RATE (EST).....170BTU/HR (50 WATTS THERMAL)
410000 RESERVOIR EC001 11113
410000 STORES COOLANT.
415500 PLMP EC001 11111
415500 PLMP COOLANT IN CPL COOLANT SUBSYSTEM
415500 FLOW RATE(EST).....40LRS/HR (300CC/MIN)
415500 HEAD(EST).....1 PSI (6.9E03 N/M2)
487000 TEMPERATURE CONTROLLER EC001 11111
487000 CONTROLS COOLANT TEMPERATURE TO CONTROL CHAMBER WALL TEMPS
487000 TEMPERATURE RANGE.....-76 TO 104F (-60 TO 40C)
487000 ACCURACY.....+/-0.2F (0.1C)
108600 CRYOGENIC COOLER ASSEMBLY EC001 11114
108600 PROVIDE LOW TEMPERATURE HEAT SINK TO EXPERIMENT CHAMBERS
108600 MAXIMUM TEMPERATURE.....-76F (-60C)
108600 HEAT SINK CAPACITY(EST).....1000BTU (290W THERMAL)
414500 TEMPERATURE SENSORS EC001 41111
414500 MEASURE COOLANT TEMPERATURES
414500 OPERATING TEMPERATURE RANGE...-76 TO 104F (-60 TO 40C)
414500 ACCURACY.....+/-0.02F (+/-0.01C)
TRCCN #SERIES RD \$1200
RESOLUTION0.1 DEG F
T/C TYPEIRON/CONSTANTAN - TYPE J
.....CHROMEL-ALUMEL - TYPE K
.....COPPER-CONSTANTAN - TYPE T
RANGE (RESPECTIVELY)-30 TO 1600 DEG F
.....-30 TO 1999 DEGF
.....-100 TO 600 DEG F
BLH ELECTRONICS #BULLETIN 209 OR 207
OPERATING RANGE.....-100 TO 200F
OMEGA ENGINEERING, INC. MODEL XT-10C \$295 11
TEMP. RANGE -100 TO +175C
NO. RANGES 4
ACCURACY 1C
RESPONSE TIME 1 SEC.
LENGTH EXTENT. LEADS UP TO 1000 FT.
104180 CRT DISPLAY EC001 11111
104180 PROVIDES DISPLAY CAPABILITY INCLUDING ALPHANUMERICS, DYNAMIC AND
104180 DEFLECTION AND VIDEO AMPLIFIERS AND REQUIRED POWER SUPPLIES.
104180 RESEARCH INC #330C \$1580
FORMAT.....24 LINESX72 OR 80 CHARACTERS
.....24 LINESX40 CHARACTERS
.....12 LINESX72 OR 80 CHARACTERS
REFRESH RATE.....60 HZ
TRANSFER RATE.....110 TO 2400 BAUD
CHARACTER FORM.....5X7 DOT MATRIX
DIMENSIONS.....15-1/2WX13-1/2HX23-1/2D INCH
WEIGHT.....39LB (17.7KG)
BUKNER-RAND CORP #2217-12CPT \$950
CHARACTER CAPACITY (MAX)960
CHARACTERS/LINE (MAX)80
LINES/DISPLAY (MAX)24
CHARACTERS REPERTOIRE42 OR 92
VIEWING AREA8.75 X 6.25
REFRESH AREA54 FRAMES/SEC
CHARACTER GENERATING METHOD ... 5X7 DOT MATRIX

ORIGINAL PAGE 13
OF POOR QUALITY



Space Division
Rockwell International

TEKTRONIX INC #40024/021-0013-00/ 613774 111
071-00XX-CG/4901/
-51

ALPHANUMERIC FORMAT 39 LINES-65 NORMAL/ITALIC CHAR,
1 LINE- 84 CHAR IN SCRATCH AREA
CHARACTER SET 96 UPPER/LOWER CHAR (ASCII)
CHARACTER SIZE 70 X 90 MILS
CHARACTER GENERATION 7 X 6 DOT MATRIX
CURSOR PULSATING 7 X 6 MATRIX
GRAPHIC INPUT MODE 1024(X) X 1024(Y) ADDRESSABLE,
1024(X) X 768(Y) VIEWABLE PTS;
JOYSTICK CONTROLLER, CROSS-HAIR
CURSOR.

1041R0 DIGITAL VISUAL DISPLAY 61001 11111
1041R0
1041R0 DISPLAY NUMERIC DATA
1041R0
101140 ALPHANUMERIC KEYBOARD 61001 11111
101140
101140 ALLOWS CREWMAN TO COMMUNICATE WITH THE UNMANNED CONTROLLER FOR
101140 EXPERIMENT CONTROL AND DATA ANALYSIS
101340
101340 KEYBOARD.....TYPEWRITER TYPE
101340 INTERFACE DESCRIPTION.....DIGITAL-12-BIT WORD
101140

PURKER-RAND CORP #2200 6440 111
KEYS TYPEWRITER SET STD
PROGRAM ASSIST KEYS 16
INTERFACE CRT
FOOTING KEYS FULL SET

927000 SCLENDID VALVES 61001 101121
927000
ALTEX SCIENTIFIC #200 111
VARIOUS CHEMICALLY INERT FITTINGS AND VALVES
TUBING 64 - 612/10FT
VALVES 645 - 655
INJECTION VALVES 6115 - 6140

927000 CENTRAL FLUID FLOW
927000
927000 MODULATING VALVE 61001 21121
927000
927000 ADJUST FLUID FLOW
927000
ALTEX SCIENTIFIC #200 111
VARIOUS CHEMICALLY INERT FITTINGS AND VALVES
TUBING 64 - 612/10FT
VALVES 645 - 655
INJECTION VALVES 6115 - 6140

927000 PRESSURE REGULATORS 61001 51121
927000
927000 PRESSURE REGULATING VALVE FOR GASEOLS PRESSURIZATION SUBSYSTEM
927000
927000 REGULATION RANGE.....2.7 TO 14.7 PSIA(140 TO 760TORR)
927000 UPSTREAM PRESSURE(EST).....20 PSI (1.49CS N/M2)
927000

CARLE INSTRUMENTS INC #6636 672 111
REGULATED PRESSURE 0 TO 60 PSIG
SOURCE PRESSURE UP TO 3000 PSIG
REGULATION +/- 0.1%

927000 PRESSURE REGULATOR 61002 11121
927000
927000 VALVE REGULATING SUPPLY PRESSURE
927000
927000 REGULATION PRESSURE(EST).....20 PSI (1.4604 N/M2)
927000 MAX UPSTREAM PRESSURE.....300 PSI (2.1 E06 N/M2)
927000

CARLE INSTRUMENTS INC #6636 672 111
REGULATED PRESSURE 0 TO 60 PSIG
SOURCE PRESSURE UP TO 3000 PSIG
REGULATION +/- 0.1%

ALTEX SCIENTIFIC #200 111
VARIOUS CHEMICALLY INERT FITTINGS AND VALVES
TUBING 64 - 612/10FT
VALVES 645 - 655
INJECTION VALVES 6115 - 6140

927000 SCLENDID VALVES 61002 11121
927000
ALTEX SCIENTIFIC #200 111
VARIOUS CHEMICALLY INERT FITTINGS AND VALVES
TUBING 64 - 612/10FT
VALVES 645 - 655
INJECTION VALVES 6115 - 6140

927000 GAS SHUTOFF VALVES
927000
ALTEX SCIENTIFIC #200 111
VARIOUS CHEMICALLY INERT FITTINGS AND VALVES
TUBING 64 - 612/10FT
VALVES 645 - 655
INJECTION VALVES 6115 - 6140



159100 CONDITIONING CHAMBER E0015 11134
159100
159100 THE CONDITIONING CHAMBER IS MADE UP OF A PRESSURE SHELL, A
159100 POSITIVE EXPULSION REGULATOR AND A WV AEROSOL CONDITICNER, THIS
159100 CHAMBER CONDITIONS THE GAS SAMPLE PRIOR TO ITS INJECTION INTO
159100 THE CLOUD CHAMBER
159100
159100 OPERATING PRESSURE(EST).....20PSIA (14E03N/M2)
159100 OPERATING TEMPERATURE.....32 TO 86F 10 TO 30C
159100 VOLUME.....2 CU.FT. (0.0565 CU. METERS)
400800 PRESSURE SENSORS E0002 51111
400800
400800 MEASURE PRESSURES THROUGH OUT THE SYSTEM
400800
400800 RANGES.....2 TO 20 PSIA (100 TO 1000 TORR)
40080010 TO 50 PSI (6.9E4-34.5E4N/M2)
40080010 TO 300 PSI (34.5E4-2.1E6N/M2)
400800 ACCURACY.....1% OF FULL SCALE
400800
BENDIX ENVIRONMENTAL SCI.DIV. MODEL 655 111
TEMP.RANGE -65 TO +150C
-45 TO +300C
PRESS.RANGE -2 TO +18PSI
PROOF PRESS. 0 TO 15 PSIG
ALH ELECTRONICS RDHF 5275 - 350 111
FULL SCALE RANGES.....FROM 20PSI
RATED OUTPUT.....3.0 MV/V
CALIBRATION ACCURACY.....0.25% FS
211000 HUMIDIFIER ASSEMBLY F0003 11135
211000
211000 THIS ASSEMBLY CONSISTS OF AN EVAPORATOR, A HEATER, A WATER TRAP,
211000 A METERING PUMP AND A WATER TANK WITH AN INTEGRAL BLADDER, THIS
211000 ASSEMBLY GENERATES THE WATER VAPOR NECESSARY FOR CLOUD CREATION.
211000
211000 DRY AIR FLOW RATE.....175CFM (5 LITERS/MIN)
211000 MAXIMUM TEMPERATURE.....77F (45C)
211000 HUMIDITY RANGE.....C TO 100%
211000 HEATER CAPACITY(MIN).....17BTU/HR (5% THERMAL)
211000 HEATER FEED RATE (MAX)0.02 LB/HR (9.1 CC/HR)
211000 WATER TANK CAPACITY.....9LBS (1.35KG)
211000
211000 DEW POINT SENSOR E0003 11113
211000
211000 THIS SENSOR MEASURES THE DEW POINT OF THE GAS SAMPLE
211000
211000 MEASUREMENT RANGE.....10 TO 98
211000 ACCURACY.....+/-1%
211000
211000 LIQUID WATER CONTENT METER E0003 11235
211000
211000 MONITORS LIQUID WATER CONTENT WITHIN THE EXPANSION CHAMBER FOR
211000 THE PURPOSE OF ACCURATE WATER BUDGET ACCOUNTS. CURRENT CONCEPTS
211000 USE OPTICAL TECHNIQUES TO MEASURE THE WATER DROPLETS IN AIR.
211000
211000 ACCURACY.....+/-0.05%
211000
603013 CONDENSER - WATER E0003 11131
603013
603013 COLLECTS WATER VAPOR BY CONDENSING WATER CONTENT OF GAS SAMPLE
603013
603013 CONDENSATION RATE(MAX).....0.2LB/HR (91CC/HR)
603013 TEMPERATURE(MIN).....45F(7C)
603013
410000 STORAGE TANK - EARTH SAMPLE E0004 11113
410000
410000 THIS TANK STORES SAMPLES OF EARTH GASSES. THE ASSEMBLY CONSISTS
410000 OF A PRESSURE VESSEL, A POSITIVE EXPULSION REGULATOR AND A
410000 BLADDER.
410000
410000 PRESSURE.....1 ATM (760 MM HG)
410000 VOLUME.....10 CU.FT. (0.28 CU. METERS)
410000
410000 GAS SAMPLE STORAGE TANKS E0004 51111
410000
410000 PRESSURE VESSELS FOR STORING EXPERIMENT GAS SAMPLES.
410000
410000 MAX PRESSURE(EST).....300 PSI (2.1E06 N/M2)
410000 VOLUME 1 CU.FT. (0.028 CU. METERS)
410000
410000 SUMP STORAGE TANK E0004 11111
410000
410000 STORE GAS SAMPLES AFTER EXPERIMENT EVALUATION.
410000
410000 MAX PRESSURE (EST).....150PSI(1.05E6N/M2)
410000 VOLUME.....10 CU.FT. (0.28 CU. METERS)
410000
526200 FLOW CONTROL ASSEMBLY E0004 11124
526200
526200 THIS ASSEMBLY CONTROLS THE FLOW OF GASSES IN THE CPL. IT IS MADE
526200 UP OF PLUMBING COMPONENTS, SOLENOID VALVES, A FILL ASSEMBLY, A
526200 VENT ASSEMBLY, SAFETY COMPONENTS AND A MIXING SUBASSEMBLY.
526200
600500 SLMP COMPRESSOR E0004 11111
600500
600500 COMPRESS EXPENDED EXPERIMENT GAS SAMPLES INTO SUMP STORAGE TANK
600500
600500 COMPRESSION RATIO(MAX).....10 TO 1

ORIGINAL PAGE IS
OF POOR QUALITY



409500
409500 CC-DC CONVERTER (VOLTAGE CONTROLLED DC) EC005 11111
409500
409500 CONVERTS SPACELAB POWER TO VOLTAGE USABLE BY ELECTRICAL FIELD
409500 GENERATING COMPONENTS
409500
409500 INPUT VOLTAGE.....28V DC
409500 OUTPUT VOLTAGE.....780
409500
409500 CC-DC CONVERTER (LOW FREQUENCY AC) EC005 11111
409500
409500 CONVERT POWER INTO USABLE FORM FOR ELECTRIC FIELD GENERATING
409500 CAPACITOR
409500
409500 FREQUENCY.....0 TO 100KHZ
409500 WAVE SHAPE.....SQUARE
409500 AMPLITUDE.....0 TO 3000 V/IN
409500
409500 ELECTRIC FIELD CONTROLLER EC005 11114
409500
409500 PROVIDE CONTROL FUNCTION FOR ELECTRIC FIELD GENERATOR, CONTROLLER
409500 SUPPORT BOTH AC AND DC POWER SUPPLIES
409500
409500 ACOUSTICAL GENERATOR E0013 11111
409500
409500 VARIABLE FREQUENCY AUDIO DRIVERS USED TO POSITION AND CONTROL
409500 WATER DROPLETS OR ICE PARTICLES IN EACH OF THREE MUTUALLY
409500 PERPENDICULAR AXES. EACH DRIVER CAN BE CONTROLLED INDEPENDENTLY
409500 IN FREQUENCY AND AMPLITUDE
409500
409500 FREQUENCY RANGE.....0 TO 10KHZ
409500 AMPLITUDE.....0 TO 780 DR
409500
409500
409500 WESTERN CORP #415 618860 113
409500 FREQUENCY 7 OR 3.5 KHZ
409500 ACOUSTIC LEVEL 107 DB
409500 BEAM WIDTH 30 DEGREES
409500
409500 TRANSDUCER - ACOUSTICAL E7013 24111
409500
409500 MEASURE SOUND INTENSITY LEVEL
409500
409500 INTENSITY LEVEL..... 0 TO 780 DR
409500
409500
409500 KAMAN CORP #KX-1800 61040 111
409500 RANGE 100-180 DB
409500 FREQUENCY RESPONSE 0-10 KHZ
409500
409500 AMPLIFIER - ACOUSTIC GENERATOR E0013 11111
409500
409500 AMPLIFY SIGNAL TO ACOUSTIC GENERATOR
409500
409500 POWER.....45W
409500
409500
409500 LASER E0013 11111
409500
409500 POLARIZED, COLLIMATED LIGHT SOURCE USED FOR PARTICLE HEATING AND
409500 MOTION CONTROL AS WELL AS TO ASSESS OPTICAL SCATTERING PROPERTIES
409500 OF ICE CRYSTALS
409500
409500
409500 PENNSYLVANIA ELECTRONICS #948 610000
409500 OUTPUT..... 5 WATTS AT 10.6 MICRONS
409500 FREQUENCY STABILITY..... < 15 MHZ AFTER 1/2 HR WARM-
409500 UP, OVER SEVERAL HOURS
409500 AMPLITUDE STABILITY..... < 5 %
409500 WEIGHT..... 27 LBS (HEAD), 40 LBS PWR SPLY
409500
409500 HOLCRAFT INC #256 613000 115
409500 WAVELENGTH..... 0.93 OR 1.06 MICRONS
409500 MODE MULTI OR TEM..
409500 CW POWER TO 50 WATTS
409500
409500 INTERNATIONAL LASER SYSTEMS #NRP-210 #36050 112
409500 TYPE ROTATING PRISM Q-SWITCHED OR
409500 POKELS CELL Q-SWITCHED
409500 WAVELENGTH 1.06 MICRONS
409500 PULSE ENERGY 100 OR 150 MJ MIN
409500 REPETITION RATE UP TO 20 PPS
409500
409500 PENNSYLVANIA ELECTRO-OPTICS ORG. #948 610000 111
409500 OUTPUT WAVELENGTH..... 10.6 MICRONS
409500 POWER OUTPUT..... 5 WATTS MINIMUM
409500 MODE PURITY..... TEM-000, SINGLE FREQUENCY
409500 BEAM DIAMETER..... 8 MM
409500 BEAM DIVERGENCE..... 5 MRAD FULL ANGLE
409500 BEAM POLARIZATION..... VERTICALLY
409500 SIZE & WEIGHT
409500 POWER SUPPLY..... 17 X 9.25 X 14 IN, 40 LBS
409500 POWER INPUT..... 115V, 60HZ AT 600VA
409500 COOLING..... 0.25 GPM, WATER
409500
409500 HOLCRAFT INC #255 610000 112
409500 WAVELENGTH..... 0.93 OR 1.06 MICRONS
409500 MODE MULTI OR TEM..
409500 CW POWER TO 25 WATTS

003002 LASER - POSITIONING SERVO CONTROL EC013 41114
 003002 PROVIDES FEEDBACK AND CONTROL SIGNAL TO LASER DRIVER FOR PROPER
 003002 POSITIONING AND HEATING OF PARTICULATE SPECIMENS
 003002
 039575 LIQUID DROP GENERATOR E0006 41235
 039575
 039575 GENERATES SINGLE DROPLETS BY MANUAL OR SYSTEM CONTROL. DROPLETS
 039575 ARE PROPELLED INTO THE EXPERIMENT CHAMBERS FOR EXAMINATION
 039575
 039575 FORMATION FREQUENCY(MAX).....10/SEC
 039575 DROPLET DIAMETER.....10 MICROMETERS TO 1 CM
 039575 DROPLET UNIFORMITY.....WITHIN 5 %
 039575
 407500 ELECTRONIC DRIVER ASSEMBLY E0006 41114
 407500
 407500 THE ELECTRONIC DRIVER ASSEMBLY CONTROLS THE LIQUID DROP
 407500 GENERATOR. IT CONSISTS OF A GENERATOR POWER SUPPLY, A PARTICLE
 407500 CHARGER POWER SUPPLY AND A DUAL PULSE GENERATOR.
 407500
 039580 ICE PARTICLE GENERATOR FC007 31235
 039580
 039580 GENERATE ICE PARTICLES FOR TEST IN ICE DIFFUSION, EXPANSION AND
 039580 GENERAL CHAMBERS
 039580
 039580 ICE CRYSTAL SIZE.....50 MICROMETERS TO 1 CM
 039580
 003400 DROPLET CLOUD AEROSOL GENERATORS E0008 11135
 003400
 003400 GENERATE AEROSOLS FOR EXPERIMENTATION
 003400
 003400 PARTICLE CONCENTRATION.....10 UP TO 10000/CC
 003400 CLOUD FORM.....POLYDISPERSED AND MONODISPERSED
 003400 PARTICLE SIZE.....10E-4 TO .1 CM
 003400 PRODUCTION RATE.....UP TO 10E5 HZ
 003400
 METRENICS ASSOCIATES #BULLETIN 24-70 111
 003400
 003400 LARGE AND GIANT NUCLET AEROSOL GENERATOR E0008 11135
 003400
 003400 SAME AS DROPLET CLOUD AEROSOL GENERATOR
 003400
 003400 AITKEN AEROSOL GENERATOR E0008 11135
 003400
 003400 SAME AS DROPLET CLOUD AEROSOL GENERATOR
 003400
 102700 DIGITAL COMPUTER E0010 11111
 102700
 102700 PROVIDES DATA PROCESSING AND CONTROL FUNCTIONS FOR THE CPL
 102700
 DATA GENERAL CORP #NOVA ROC \$18450 111
 WORD LENGTH 16 BIT
 CYCLE TIME 800 NANOSEC.
 CORE MEMORY 32K MAX
 DATA GENERAL CORP #SUPERNOVA SC \$5600
 WORD LENGTH 16 BIT
 VARIAN DATA MACHINES #520/I 111
 MEMORY CYCL TIME 1.5MICROSEC.
 MEMORY EXPANDABLE,4096 BYTES(8BITS)TO
 32,768 BYTES
 REGISTERS 12
 OPERAND PRECISION UP TO 32 BITS
 ADDRESS REFERENCE R,16,24,OR 32 BIT LEVEL
 NOTE: HAS FULL RANGE INTERFACE
 HARDWARE
 VARIAN DATA MACHINES #R-620/I 111
 BASIC COMMANDS OVER 100
 ADDRESSING MODES 6
 MAX.WORDS 32,768
 WORD LENGTH 16 OR 18 BIT
 REGISTERS 9
 NOTE: HAS FULL RANGE INTERFACE
 HARDWARE
 CLARY DATACOMP SYSTEMS INC #404 111
 SIMULTANEOUS TERMINALS 16
 INPUT TELETYPEWRITER KEYBOARD
 ACCUMULATOR 64 BIT
 REGISTERS 16 BIT INDEX (2 EACH)
 WORD LENGTH 16, 32, 48, 64 BITS
 MEMORY CAPACITY 1024 16 BIT WORDS, 4096 16 BIT
 WORDS,OR ADDITIONAL 4096 16 BIT
 WORDS, TOTAL 65536 WORDS
 GENERAL AUTOMATION INC. #18-30 111
 FUNCTION SUPERVISE SMALLER COMPUTERS
 MEMORY 8K CORE
 CYCLE TIME 1.2 MICROSEC.
 GENERAL AUTOMATION INC. #SPC-14 111
 MEMORY 16K
 WORD LENGTH 16 BITS
 READ/WRITE CYCLE MEMORY TIME ... 400 TO 1440 NANOSEC.
 READ ONLY MEMORY 400 TO 720 NANOSEC.
 INPUT/OUTPUT TIME 1.0 TO 2.8 MICROSEC.
 DMA TRANSFER RATE 0.694E+06 TO 2.5E+06

ORIGINAL PAGE 25
OF POOR QUALITY



GENERAL AUTOMATION INC.	#SVC-12		111
MEMORY	4K TO 16K		
WORD LENGTH	8 BITS		
CYCLE TIME	2.16 MICROSEC		
STORED PROGRAM EXECUTION RATE ..	0.23E+06/SEC		
I/O TRANSFER RATE	0.46E+06/SEC		
REGISTERS	6 12 BIT REGISTERS		
ACCUMULATORS	4 12 BIT		
DATA GENERAL CORP	#NOVA 1200	\$11400	
WORD LENGTH	16 BIT		
CYCLE TIME	1200 NANOSEC.		
CORE MEMORY	32K MAX.		
DIGITAL EQUIPMENT CORP	#PDP-9/E	\$22000	111
MEMORY	4096(CORE)EXPANDABLE TO 32,768 WORDS		
WORD	12 BIT		
CYCLE TIME	1.5 MICROSEC.		
OPTIONS:.....	DEC-DISC(131,072 WORDS STORAGE) MAGNETIC TAPE, TAPE READER, TAPE PUNCH, CARD EQUIP., LINE PRINTER, A/D & D/A CONVERT., CRT DISPLAY		
424000 TAPE RECORDER			ECO10 11111
424000			
424000 RECORD DATA ON MAGNETIC TAPE			
424000			
424000 BIT RATE.....	100B/S		
424000 STORAGE CAPACITY.....	3.5E7BITS		
424000			
AMPEX CORPORATION	#AR-200		121
BANDWIDTH	100 HZ TO 9125 HZ THROUGH 300 HZ TO 250 KHZ		
TAPE WIDTH	0.5 INCH		
TAPE SPEED	1.875 TO 40 IPS		
RECORDING TIME	8 MINUTES TO 4 HRS AND 16 MINS		
FORMAT	DIGITAL		
TRACKS	8 DIGITAL, 7 ANALOG		
AMPEX CORP	#AR 1700	\$20,000	117
TRACKS	14		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
RECORDING MODE	DIRECT		
PACKING DENSITY	20 KB/I/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	6.2E+10 BITS		
BANDWIDTH	2 MB/SEC/T		
TAPE LENGTH.....	9200 FEET		
SANGAM ELECTRIC	#SABER III	\$26,500	111
TRACKS	14		
TAPE WIDTH	1 IN		
TAPE SPEEDS	8 SELECTABLE SPEEDS FROM 15/16 TO 120 IPS		
FREQUENCY RESPONSE	400 HZ TO 2.0 MHZ		
RECORDING RATE	600 KBPS AT 120IPS SERIAL MODE		
WEIGHT	100 POUNDS		
DIGI-DATA CORPORATION	#1600	\$2650	111
TAPE SPEED	25, 10.75, 12.5 IPS		
TRACKS	7 OR 9 TRACK		
DATA DENSITY	1600 CPI PHASE ENCODED 200, 556, 800 CPI NRZI		
TAPE	0.5 INCH, 1.5 MIL, 1200 FEET		
DIGI-DATA CORPORATION	#1700/PDP-11/T-9	\$9250	111
TRACK NRZI			
TAPE SPEED	45, 37.5, 25, 10.75, 12.5 IPS		
TAPE	0.5 INCH, 1.5 MIL, 18P/ANSI COMPATIBLE, 10.5 INCH REEL		
TRACKS	7 OR 9		
DATA DENSITY	PHASE ENCODED COMPATIBLE		
HENLETT PACKARD	#7570B/C	\$4600	111
TAPE FORMAT ...	800, 556, OR 200 CPI NRZI AND 1600 CPI PHASE-ENCODED		
CHANNELS	7 OR 9		
TAPE SPEED	10 TO 45 IPS		
TAPE	0.5 INCH, 1.5 MILS, 18P/ANSI COMPATIBLE		
BORG WARNER	#PERT		121
TRACKS	30		
TAPE SPEED	UP TO 1000 IPS		
TAPE WIDTH	1/2 INCH		
TAPE LENGTH	2400 FEET		
PACKING DENSITY	15 KB/I/T		
SIGNAL/NOISE	24 DB		
DATA CAPACITY	1.3E+9		
BANDWIDTH	0-15 MB/SEC/T		
LEACH	#HT 7000		311
TRACKS	12		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
TAPE LENGTH.....	9200 FEET		
PACKING DENSITY	16.7 KB/I/T		
SIGNAL/NOISE	22 DB		
DATA CAPACITY	3.2E+10		
BANDWIDTH	2 MB/SEC/T		

HEWLETT PACKARD	#355D	\$10200	111
BANDWIDTH	300 KHZ		
CHANNELS	7		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#355C	\$14700	111
BANDWIDTH	300 KHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3950A-011	\$23800	111
BANDWIDTH	500 MZ TO 2 MHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3960A/13065A/13063A	\$4796	111
CONFIGURATION	PACK MOUNTED		
TAPE SPEED	15, 3 AND 3/4, 15/16 IPS		
CHANNELS	4		
RECORDING FORMAT	FM		
PASSBAND	5 KHZ		
S/N RATIO	48 DB		
TRACKS	30		
HONEYWELL	#560C	\$9730	
PORTABLE TAPE RECORDER			
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 60IPS		
MAXIMUM BANDWIDTH(DIRECT)	300KHZ		
PACKING DENSITY	UP TO 600PBI		
WEIGHT	70LBS (32KG)		
INPUT VOLTAGE	28VDC		
HONEYWELL	#96	\$17420	
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 240IPS		
REEL SIZE	16IN		
MAXIMUM BANDWIDTH(DIRECT)	2M HZ		
TAPE WIDTH	1/2 IN		
AMPEX CORP	#AP TCC	\$29,180	511
TRACKS	14		
TAPE SPEED	60 IPS		
TAPE WIDTH	1 INCH		
RECORD MODE	DIRECT		
PACKING DENSITY	20 KB/1/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	RE+10 BITS		
BANDWIDTH	1 MB/SEC/T		
603003 SIGNAL CONDITIONING ELECTRONICS		E0010	11111
603003			
603003 ELECTRONICS SUPPORTING DATA MANAGEMENT SUBSYSTEM			
603003			
NLS	#SERIES X-1	\$4000	111
ACCURACY	+/-0.00088FS		
FUNCTIONS	DC VOLTS, AC VOLTS, MILLIVOLTS, K OHMS AND RATIO		
COMMON MODE REJECTION	>130DB		
HEWLETT PACKARD	#461A	\$340	111
FREQUENCY RANGE	1 KHZ TO 150 MHZ		
FREQUENCY RESPONSE	+/- 1 DB INTO 50 OHM LOAD		
GAIN AT 500 KHZ	40 DB +/- 0.5 DB OR 20 +/- 1 DB		
MAX INPUT	1 V RMS OR 2 V P-P PULSE		
TENNELEC	#TC 202/203	\$355/495	111
GAIN	2.5 TO 3000X		
NON LINEARITY	LESS THAN 0.05%		
CONFIGURATION	NIM COMPATIBLE		
ROHDE AND SCHWARZ	#ATN (100.0899.02)		111
FREQUENCY RANGE	30 HZ TO 20 KHZ		
INPUT VOLTAGE	0.25 - 25 VOLTS		
MAX POWER	50 WATTS		
RESEARCH INC	#R12-11-17	\$3200	111
CAPACITY	TO 33 PAIRS OF SIG WIRES (160GA)		
REFERENCE TEMP	HEAT SINK CONTROLLED TO +/- .08 DEG C AT 65 DEG C		
AMBIENT TEMP	0 TO 49 DEG C		
LEAKAGE RESISTANCE	1000 MEG OHM		
HEAT-UP TIME FROM 22 DEG C	50 MINUTES		
SETTLING TIME	8 MS TO FINAL VALUE		
HONEYWELL	#ACCUOTA 117	\$1340	111
MULTICHANNEL WIDEBAND AMPLIFIER			
CHANNELS	7		
POINT	BENCH		
HONEYWELL	#ACCUOTA 115	\$5370	111
HIGH VOLTAGE DC AMPLIFIER			
INPUT SIGNAL VOLTAGE	+/-50MV TO 1500V		
OUTPUT SIGNAL VOLTAGE	UP TO 2000V		
CHANNELS	6		
NEFF	#122-222	\$610	111
GAIN	1 TO 2500		
COMMON MODE REJECTION	140DB		
DRIFT	0.5MICROVOLT/DC		
STEP GAIN ACCURACY	+/-0.01%		
SWITCH SELECTABLE FILTER			
BANDWIDTH	100KHZ		
PACK/MODULE MOUNT			

ORIGINAL PAGE 13
OF POOR QUALITY



ACFF #126-621 8320 111
AMPLIFIER MULTIPLEXER
EMPLOYES PET SWITCH
GAIN RANGE..... 0.2 TO 2500
COMMON MODE REJECTION..... 120DB
OUTPUT..... +/-10V AT 10 MILLIAMPERES

NEFF #127-1C2 #1190
PROGRAMMABLE GAIN WITH SWITCH SELECTABLE OUTPUT FILTER
GAIN STEPS..... 1, 2, 5, 10, 20, 50, 100, 200,
100, 1000
COMMON MODE REJECTION..... 140DB
PET SWITCHING
RACL/MODULE CONFIGURATION

TENNELEC #TC133 8295 111
PET PREAMPLIFIER
CHARGE SENSITIVITY..... 10E12 VOLTS/COULOMB
INPUT..... NEGATIVE INPUT PULSE POLARITY
WITH INPUT CAPACITANCE >2000PF
OUTPUT..... POSITIVE OUTPUT PULSE POLARITY
WITH OUTPUT IMPEDANCE 50 OHM
STABILITY..... 0.5%/V

367000 OPTICAL PARTICLE COUNTER E0011 11235
367000
367000 THIS DEVICE WILL BE A SIZE ANALYZER. THE ANALYZER CONSISTS OF A
367000 PHOTOMULTIPLIER, AN AMPLIFICATION STAGE AND SUPPORTING
367000 ELECTRONICS
367000
367000 SENSITIVE PARAMETER.....LIGHT SCATTERED FROM PARTICLES
367000 PARTICLE SIZE.....>0.3 MICROMETERS
367000 THERMO-SYSTEMS INC #3C30 810700 111
367000 PARTICLE SIZING RANGE..... .0032 TO 1 MICRON
367000 CLASSIFICATION RANGE..... .01 TO .5 MICRON
367000 SAMPLING RATE..... 5 LITERS/MIN
367000 COUNTING CLASSIFICATIONS..... 10 SIZE CLASSES
367000 READOUT..... COMPATIBLE WITH AUT DATA AC SYS

367500 NUCLEI MASS MONITOR SYSTEM E0011 11235
367500
367500 THIS DEVICE MEASURES THE TOTAL PARTICULATE MASS PER UNIT VOLUME.
367500 THE SYSTEM USES A CRYSTAL OSCILLATOR WHICH CHANGES ITS RESONANT
367500 FREQUENCY AS PARTICLES ARE DEPOSITED ON IT.
367500
367500 PARTICLE COLLECTION.....ELECTROSTATIC PRECIPITATION
367500 PARTICLE SIZE.....0.01 TO 20 MICROMETERS
367500
367000 ELECTRICAL PARTICLE COUNTER E0011 11135
367000
367000 THIS ANALYZER MEASURES THE SIZE DISTRIBUTION OF PARTICLES. THE
367000 PARTICLES ARE IONIZED AND THEIR MOBILITY AS A FUNCTION OF
367000 ELECTRIC FIELD IS MEASURED TO GIVE AN INTEGRAL SIZE DISTRIBUTION
367000
367000 PARTICLE SIZE.....0.01 TO 1.0 MICROMETER
367000 FLOW RATE.....1.5CFM (50LITERS/MIN)
367000 INTEGRAL DISTRIBUTION CAPABILITY.....YES
367000

ORTEC #401A, 1C9A, 410, 83000
#39, 415
367000 RESOLUTION.....<2.5KEV
367000 LINE SHAPING.....SINGLE OR DOUBLE DELAY, CP
SINGLE OR DOUBLE DIFFEREN-
TIATION WITH TIME CONSTANTS
OF 0.1 TO 10 MSEC
367000 COUNT RATE STABILITY.....<0.2% FOR 50000 COUNTS/SEC
367000 GAIN SHIFT<0.2% FOR 50000 COUNTS/SEC

PAUSCH AND LOEB INC #40-1A
367000 SAMPLING RATE 0.1 FT/0.3/100 SEC (1E+7 PARTICLES
PER FT/0.3)
367000 MAX CONCENTRATION MEAS 1E+7 PARTICLES/FT/0.3
35E+7 PARTICLES/LITER
367000 PARTICLE SIZE RANGE 0.3, 0.5, 1.0, 2.0, 3.0, 5.0, 10
MICRONS

170900 ANALYZER - MULTICHANNEL E0011 11111
170900
170900 ALLCH SELECTIVE ANALYSIS OF COUNTER DATA
170900
170900 NUCLEAR DATA # NO 100 87900
170900 DATA ACQUISITION MODES PHA, MCS, OR LIST
170900 BANDWIDTH 15 MHZ
170900 STORAGE 12 BIT DATA WORDS
170900 DISPLAY CRT, 6X10 CM
170900 NUMBER OF CHANNELS 1024, 2048, 4096
170900 COUNT CAPACITY PER CHANNEL 20020 +0 1 + 4 INTERNAL
FLAG BITS
170900 DIMENSIONS..... 19X10.25X23.5 INCHES
170900 WEIGHT 60 LBS

207500 EXPANSION CLOUD CHAMBER E0014 11134
207500
207500 CHAMBER PROVIDING SIMULATION OF LONG TERM, NATURAL CLOUD DURING
207500 AN ADIABATIC EXPANSION
207500



297500 DIMENSIONS.....11.8IN(30CM) DIA. X 17.7IN
297500 149CM LONG
297500 VOLUME.....1.12CU.FT.(31.8L)
297500 OPERATING TEMPERATURE.....-76F TO 104F (-60 TO 40C)
297500 UPPER AND LOWER ENDS TEMP
297500 DIFFERENTIAL.....0 TO 18F(0 TO 10C)
297500 END TEMPERATURE TOLERANCE.....+/-0.1F (+/-0.5C)
297500 OPERATING PRESSURE.....2.7 TO 14.7 PSIA(140 TO 760TORR)
297500 PRESSURE RATE OF CHANGE.....19.4PSI/SEC (1000TORR/SEC)
297500 PRESSURE RATE TOLERANCE.....+/-1%
297500 MAXIMUM VOLUME EXPANSION.....V/V OF 0.5+/-0.1%
297500 CIRCULAR ENDS CONTAIN ELECTRICALLY CONDUCTIVE PLATES TO SERVE
297500 AS EQUIPOTENTIAL SURFACES FOR THE ELECTRIC FIELD MOTION
297500 CONTROL SYSTEM
297500
297500 CONTINUOUS FLOW DIFFUSION CLOUD CHAMBER E0014 11135
297500
297500 SUPERSATURATION IS CONTROLLED BY THE TEMPERATURE OF WATER
297500 COVERING THE UPPER AND LOWER SURFACES OF THE CHAMBER. IT WILL BE
297500 USED FOR CLOUD CONDENSATION NUCLEATION EXPERIMENTS
297500
297500 DIMENSIONS.....11.8IN(30CM) X 11.8IN(30CM) X
297500 2IN(5CM)
297500 VOLUME.....0.16CU.FT.(4.5L)
297500 INNER PLATE DIMENSIONS.....11.8IN(30CM) X 10IN(25CM) SPACED
297500 0.5IN(1.3CM) APART
297500 ISOTHERMAL TEMP CONTROL.....+/-0.1F(0.05C)
297500 UPPER AND LOWER SURFACE
297500 TEMPERATURE DIFFERENCE.....0 TO 18F(0 TO 10C)
297500 SURFACE TEMP CONTROL.....+/-0.1F(0.05C)
297500 OPERATING PRESSURE.....14.3 TO 14.7PSIA(740 TO 760TORR)
297500 PRESSURE MEAS. ACCURACY.....+/-0.2%
297500 MAX PRESSURE CHANGE RATE.....0.005PSI/SEC(0.4TORR/SEC)
297500
297500 STATIC DIFFUSION ICE CLOUD CHAMBER E0014 11135
297500
297500 THIS IS A NAKAYA TYPE CHAMBER WHICH UTILIZES ICE SURFACES TO
297500 PROVIDE CONTROLLED SUPERSATURATION RELATIVE TO ICE.
297500
297500 DIMENSIONS.....15.75IN(40CM) DIA X 7.9IN(20CM)
297500 LONG
297500 VOLUME.....0.444CU.FT.(12.5L)
297500 OPERATING TEMPERATURE.....-40 TO 77F(-40 TO 25C)
297500 TEMP MEAS. ACCURACY.....+/-1.0F(+/-0.5C)
297500 UPPER AND LOWER SURFACE
297500 TEMPERATURE DIFFERENCE.....UP TO 36F(20C)
297500 SURFACE TEMP MEAS. ACCURACY.....+/-0.1F(+/-0.05C)
297500 OPERATING PRESSURE.....1.9 TO 14.7PSIA(100 TO 760TORR)
297500 PRESSURE MEAS. ACCURACY.....+/-0.19PSI(+/-10TORR)
297500
297500 STATIC DIFFUSION LIQUID CLOUD CHAMBER E0014 11135
297500
297500 A THOMY TYPE CHAMBER USED FOR EXPERIMENTS REQUIRING ABOVE
297500 FREEZING TEMPERATURES AND SUPERSATURATION OF THE LIQUID RELATIVE
297500 TO WATER. SUPERSATURATION IS CONTROLLED BY THE TEMPERATURES OF
297500 THE WATER COVERING CHAMBER UPPER AND LOWER SURFACES.
297500
297500 DIMENSIONS.....6IN.(15 CM) DIA. X 0.6IN(1.5CM) LONG
297500 VOLUME.....0.01CU.FT.(0.27L)
297500 OPERATING TEMPERATURE.....32F TO 68F(0 TO 30C)
297500 END SURFACE MAXIMUM
297500 TEMPERATURE DIFFERENTIAL.....18F(10C)
297500 SURFACE TEMP MEAS. ACCURACY.....0.1V(0.05C)
297500 OPERATING PRESSURE.....2.84 TO 14.7PSIA(140 TO 760TORR)
297500
297500 GENERAL CLOUD CHAMBER E0014 11135
297500
297500 THIS CHAMBER WILL BE USED FOR MANY EXPERIMENTS THAT REQUIRE A
297500 RELATIVE HUMIDITY BELOW 100% AND MINIMUM TEMPERATURE CONTROL.
297500 PROVISIONS WILL BE MADE FOR GENERATING VARIOUS ELECTRIC FIELDS.
297500 POSITIONING DEVICES(SOUND, OPTICAL, ELECTRICAL), AND REMOTE
297500 DROPLET SIZING.
297500
297500 DIMENSIONS.....11.8IN(30CM) ALL SIDES (CUBE)
297500 VOLUME.....0.95CU.FT.(27L)
297500 OPERATING TEMPERATURE.....50 TO 86F(10 TO 30C)
297500 OPERATING PRESSURE.....2.84 TO 14.7PSIA(140 TO 760TORR)
297500 ELECTRIC FIELD PLATES LOCATED ON OPPOSITE SIDES
297500 ACOUSTIC DRIVERS LOCATED ON THREE MUTUALLY PERPENDICULAR SIDES
297500
409500 AC POWER CONTROLLER E0017 11111
409500
409500 REGULATES AND CONTROLS AC POWER
409500
409500 REGULATION RANGE.....TBD
409500
JOHN FLUKE MFG CO INC #4210A 51445 111
PROG CODING.....BCD (BINARY AVAILABLE)
OUTPUT VOLTS.....0 TO +/- 10 VOLTS DC/PEAK AC
OUTPUT CURRENT.....100 MA
FREQUENCY RANGE.....DC TO 100 KHZ
RESOLUTION.....1 MV
REGULATION.....0.001%
JOHN FLUKE MFG CO INC #4250A 51795 111
PROG CODING.....BCD (BINARY AVAILABLE)
OUTPUT VOLTS.....0 TO +/- 65 VOLTS DC/PEAK AC
OUTPUT CURRENT.....1 AMP
FREQUENCY RANGE.....DC TO 30 KHZ
RESOLUTION.....1 MV
REGULATION.....0.001%

ORIGINAL PAGE IS
OF POOR QUALITY



JCHN FLUKE MFG CO INC #427CA #1495 111
PROG CODING BCD (BINARY AVAILABLE)
OUTPUT VOLTS 0 TO +/- 100 VOLTS DC/PEAK AC
OUTPUT CURRENT 0.5 AMPS
FREQUENCY RANGE DC TO 30 KHZ
RESOLUTION 1 MV
REGULATION 0.001%

405500 DC POWER CONTROLLER EC017 11111
405500
405500 REGULATES AND CONTROLS DC POWER
405500
405500 REGULATION RANGE.....TAD
405500

JCHN FLUKE MFG CO INC #421GA #1445 111
PROG CODING BCD (BINARY AVAILABLE)
OUTPUT VOLTS 0 TO +/- 10 VOLTS DC/PEAK AC
OUTPUT CURRENT 100 MA
FREQUENCY RANGE DC TO 100 KHZ
RESOLUTION 1 MV
REGULATION 0.001%

JCHN FLUKE MFG CO INC #425UA #1795 111
PROG CODING BCD (BINARY AVAILABLE)
OUTPUT VOLTS 0 TO +/- 65 VOLTS DC/PEAK AC
OUTPUT CURRENT 1 AMP
FREQUENCY RANGE DC TO 30 KHZ
RESOLUTION 1 MV
REGULATION 0.001%

JCHN FLUKE MFG CO INC #427CA #1899 111
PROG CODING BCD (BINARY AVAILABLE)
OUTPUT VOLTS 0 TO +/- 100 VOLTS DC/PEAK AC
OUTPUT CURRENT 0.5 AMPS
FREQUENCY RANGE DC TO 30 KHZ
RESOLUTION 1 MV
REGULATION 0.001%

054000 ACTION PICTURE CAMERA (FAST FRAME) EC009 11111
054000
054000 FAST FRAME RATE CAMERA TO RECORD DROPLET COLLISIONS
054000
054000 FRAME RATE.....100/SEC
054000 RECORDING DURATION.....20SEC
054000 FIELD OF VIEW.....WIDE ANGLE
054000
055090 STILL CAMERA EC009 11111
055090
055090 RECORD CLOUD CHAMBER PHENOMENA
055090
055090 FRAME RATE.....2/SEC
055090 SHUTTER SPEED.....1SEC TO 1/500SEC
055090 CAPACITY.....500FRAMES
055090 INTERCHANGABLE LENSES
055090 FILM SIZE.....35MM
055090

MCNEWMILL INC #SPOTOMATIC F #63A 111
LENS 55 MM, F/1.8 TAKUMAR
SHUTTER SPEEDS 1/1000 SEC TO TIMER
VIEWING THRU LENS

055090 STILL CAMERA EC009 11111
055090
055090 RECORD CLOUD CHAMBER PHENOMENA
055090
055090 SHUTTER SPEED.....1SEC TO 1/500SEC
055090 RESOLUTION.....2 MICROMETERS
055090 DEPTH OF FIELD.....1MM
055090 FRAME SIZE.....70MM
055090

GRAPLEX #XL #305 111
CONFIGURATION HAND HELD OR STAND MOUNTED
FILM FORMAT 3.25 X 4.25 INCH
LENS INTERCHANGABLE

HASSELHAG #500 EL #1411 111
TYPE SINGLE LENS REFLEX
DRIVE ELECTRICAL OR MANUAL
LENS 80 MM PLANAR F/2.8 ZEISS
CLOSE-UP LENS 90 MM DISTACON F/4.0
PROXAR 0.5 FOCUS RANGE
17.25 TO 24.25 IN
PROXAR 1.0 FOCUS RANGE
22.25 TO 42.5 IN
FILM SIZE 2.25 X 2.25 IN
MAGAZINE CASSETTE AND PLATE

268500 MICROSCOPE EC017 11111
268500
268500 MAGNIFIED OBSERVATIONS
268500
268500 MAGNIFICATION FACTOR.....10 TO 100X
268500 OBSERVATION MODE.....VISUAL OR PHOTOGRAPHIC
268500

FAUSCH AND LOMA INC #R57-63W #876 111
OBJECTIVES 4 FLAT FIELD
TYPE OBJECTIVE 5X, 0.8 NA; 10X, 0.25 NA; 40X,
0.65 NA; 100X, 1.25 NA
ILLUMINATOR 5 STEP, LOW VOLTAGE TRANSFORMER
EYEPIECES PAIRED 10X



Space Division
Rockwell International

LEITZ #LZ051 \$2101 111
CONFIGURATION..... BENCH MOUNT
EYEPiece..... BINOCULAR
MAGNIFICATION..... 6.3X TO 100X

AMERICAN OPTICAL #XH20T6-DW \$1576 111
CONFIGURATION..... BENCH MOUNT
EYEPiece..... BINOCULAR
MAGNIFICATION..... 4X TO 100X

OLYMPUS OPTICAL CO LTD #JM \$900 111
CONFIGURATION..... LOOSE EQUIPMENT
TYPE..... DARKFIELD ZOOM-STEREO
MAIN BODY..... BINOCULAR TUBE, INCLINED 45 DEG
60 DEG, ROTATABLE 360 DEG
MAGNIFICATION RANGE..... 3.5X TO 160X
WORKING DISTANCE..... 30, 45, 88, 105, 159 MM
APERTURE IRIS DIAPHRAGM..... 2 TO 40 MM (ADJUSTABLE)

BAUCH & LOMB #PB-252 \$1252 111
CONFIGURATION..... LOOSE EQUIPMENT
EYEPiece..... BINOCULAR
MAGNIFICATION..... 4X TO 100X

269000 STEREO MICROSCOPE EC017 11111
269000
269000 MAGNIFIED OBSERVATIONS
269000
269000 MAGNIFICATION FACTOR..... 10 TO 100X
269000 OBSERVATION MODE..... VISUAL OR PHOTOGRAPHIC
269000 WORKING DISTANCE..... UP TO 4 IN (10 CM)
269000

BAUCH & LOMB INC. #B57-63W \$976 111
OBJECTIVES..... 4 FLAT FIELD
TYPE OBJECTIVES..... 5X, 0.08NA
10X, 0.25NA
40X, 0.65NA
100X, 1.25NA
ILLUMINATOR..... 5-STEP, LOW VOLT. XFORMER
EYEPieces..... PAIRED 10X

LEITZ #LZ051 \$2101 111
CONFIGURATION..... BENCH MOUNT
EYEPiece..... BINOCULAR
MAGNIFICATION..... 6.3X TO 100X

AMERICAN OPTICAL #XH20T6-DW \$1576 111
CONFIGURATION..... BENCH MOUNT
EYEPiece..... BINOCULAR
MAGNIFICATION..... 4X TO 100X

OLYMPUS OPTICAL CO LTD #JM \$900 111
CONFIGURATION..... LOOSE EQUIPMENT
TYPE..... DARKFIELD ZOOM-STEREO
MAIN BODY..... BINOCULAR TUBE, INCLINED 45 DEG
60 DEG, ROTATABLE 360 DEG
MAGNIFICATION RANGE..... 3.5X TO 160X
WORKING DISTANCE..... 30, 45, 88, 105, 159 MM
APERTURE IRIS DIAPHRAGM..... 2 TO 40 MM (ADJUSTABLE)

BAUCH & LOMB #PB-252 \$1252 111
CONFIGURATION..... LOOSE EQUIPMENT
EYEPiece..... BINOCULAR
MAGNIFICATION..... 4X TO 100X

268500 MICROSCOPE (HIGH MAGNIFICATION) EC017 11111
268500
268500 MAGNIFIED OBSERVATIONS
268500
268500 MAGNIFICATION FACTOR..... 100 X 1000X
268500 OBSERVATION MODE..... TRINOCULAR (CAMERA ATTACHMENT)
268500
269000

CARL ZEISS CO #ULTRAPHOT II \$8000 111
CONFIGURATION..... LOOSE EQUIP
MAGNIFICATION..... 2.5X TO 2500X
CAMERA ATTACHMENT..... POLAROID, 4X5 FILM, 35 MM ROLL

359500 OSCILLOSCOPE EC017 11111
359500
359500 SUPPORT MAINTENANCE ACTIVITIES AND SPECIAL DATA ACQUISITION
359500
359500 FREQUENCY RANGE..... DC TO 10.0 MHZ
359500 OBSERVATION MODE..... VISUAL AND CAMERA (1 HR STORE)
359500

TEKTRONIX INC #R564B/3A6 PLUG-IN/ \$2340 111
3A6 PLUG-IN
BANDWIDTH..... DC TO 10 MHZ
CHANNELS..... 2
DEFLECTION FACTOR..... 10 MV/DIV TO 10 V/DIV
TIME BASE..... 50 NS/DIV TO 5 S/DIV
MAX VIEWING TIME..... 1 HOUR
MAX STORED WRITING SPEED..... 500 DIV/MS
TYPE OF STORAGE..... SPLIT SCREEN BISTABLE

TEKTRONIX #DM64 \$1095 111
BANDWIDTH..... DC TO 10 MHZ
RISE TIME..... 35 NS
SPEED WRITING..... 25 CM/MS (NOMINAL)
250 CM/MS (ENHANCED MODE)
STORAGE VIDEO TIME..... 1 HR
ERASE TIME..... 0.25 S
DEFLECTION FACTOR..... 10 MV/CM TO 50 V/CM IN 12 STEPS
TIME BASE..... 100 NS/CM TO 2 S/CM IN 23 STEPS
RESPONSE TIME..... LESS THAN 1 NS

ORIGINAL PAGE IS
OF POOR QUALITY



```

PERCUT-PAC-40      0122A      1150J      111
BANDWIDTH ..... DC TO 400KHZ
SENSITIVITY ..... 100 MICROVOLTS/CM
INDEPENDENT ..... 2
WEIGHT ..... 43 POUNDS

TEKTRONIX INC      #R465      11800      111
BANDWIDTH ..... DC TO 100 MHZ
CHANNELS ..... 2
DEFLECTION FACTOR ..... 5 MV/DIV TO 5 V/DIV
TIME BASE ..... 0.01 MICROSEC/DIV TO 0.5 S/DIV

TEKTRONIX      #R45      14200      111
RANGE, CALIBRATED ..... 5MV/DIV TO 5V/DIV IN 10 STEPS
UNCALIBRATED ..... CONT VAR TO 17.5V/DIV
RESPONSE TIME ..... LESS THAN 1 NS
CHANNELS ..... 1, ALTERNATE, CHOPPED, AND CD, X-Y
CHANNEL 2 (+UP OR INVERTED)
DIMENSIONS ..... 14.7 X 52.4 X 30.5 CM
WEIGHT ..... 9.5 KG
POWER ..... 115/230V 48-440W 60W

TEKTRONIX INC      #R41R/3A6 PLUG-IN/ 11840      111
3A4 PLUG-IN
FREQUENCY BANDWIDTH ..... DC TO 10 MHZ
CHANNELS ..... 2
DEFLECTION FACTOR ..... 10 MV/DIV TO 10 V/DIV
TIME BASE ..... 50 NS/DIV TO 5 S/DIV
RISETIME ..... 35 NS

484000 TV MONITOR      ECO12      11111
48400J
484000 DISPLAY TV CAMERA OUTPUT SHOWING CHAMBER PHENOMENA
484000
COMU ELECTRONICS INC      #CDF      1775      111
BANDWIDTH ..... 30 MHZ
HORIZONTAL LINES ..... 1725
FIELDS ..... 60 PER SEC

COMU INC, ELECTRONICS DIV      #DFFTC-TV      111
PICTURE SIZE ..... 7-3/16 X 5-3/16 INCHES
HORIZONTAL SCAN ..... 945 LINES, 50 OR 60 FIELDS/SEC
..... 873 LINES, 60 FIELDS/SEC
VIDEO BANDWIDTH ..... 20 MHZ +/- 2 DB
POWER ..... 25-28.5 VDC, 150 WATTS

480000 VIDICON      ECO12      11111
480000
480000 RECORD CHAMBER PHENOMENA
480000
480000 RESOLUTION ..... HIGH
480000 LIGHT LEVEL ..... LOW
480000 IMAGE MAGNIFICATION ..... YES
480000 IMAGE INTENSIFICATION ..... YES
480000
480000
480000 RCA ELECTRONICS COMPONENTS      #R*21      1520      111
TYPE ..... VIDICON-SULFIDE
PHOTOCONDUCTIVE ..... 11
IMAGE DIAGONAL ..... 25 MM (1 IN)
FOCUS ..... MANUAL
LIMITING RESOLUTION ..... 1500 LINES

COMU INC, ELECTRONICS DIV.      #1220      18430      111
IMAGE CONVERTER ..... 4846B VIDICON - COLOR
RESOLUTION-HORIZONTAL LIMIT ..... 300 LINES MIN.
GEOMETRIC DISTORTION ..... < 2 % OF PIC HEIGHT
LENS MOUNT ..... 16 MM - C-MOUNT

COMU INC, ELECTRONICS DIV.      #2000      13000      311
LENS-BUILT IN ..... 4:1 (20-80 MM) F/2.5 ZOOM
LENS ATTACHMENT ..... 4:1 (12.5-80 MM)
LENS INTERCHANGEABLE ..... 10:1 (15-150 MM) F/2.8 ZOOM
..... 16 MM C-MOUNT FFL
10 MHZ BANDWIDTH (MOD 2004) ..... 525 OR 725 LINE (VIDICON 7263A)
20 MHZ BANDWIDTH (MOD 2004) ..... 873 OR 945 LINE (VIDICON 8573)
PIL SPECIFICATIONS ..... MIL-F-5272C, MIL-STD-A10

COMU INC, ELECTRONICS DIVISION      #4500      12225      311
LENS ..... 4:1 ZOOM
..... 10:1 ZOOM
VERTICAL SWEEP RATE ..... 60 FIELDS PER SEC
HORIZONTAL SWEEP RATE ..... 525 LINES PER FRAME
IMAGE TUBE TYPE ..... 8541A STD
LENS MOUNT ..... 16 MM C-MOUNT
PIL SPECIFICATIONS ..... MIL-E-5400M, MIL-E-5300M

234000 IMAGE DEVICE LIGHTING      ECO09      11114
234000
234000 PROVIDES LIGHTING FOR VARIOUS OPTICAL RECORDING ACTIVITIES.
234000 CONSISTS OF CONTINUOUS LIGHTING, STROBOSCOPIC AND LASER BEAM LIGHT
23400J SOURCES
234000
234000 HIGH INTENSITY LIGHT DURATION...15SEC
234000
234000
355400 OPTICAL DETECTOR      ECO09      11235
355500
355500 SOLID STATE DETECTOR TO OBTAIN DATA CONCERNING SCATTERING
355500 PROPERTIES OF ICE CRYSTALS AND POSSIBLY EXTREME DROPLET SIZE
355500 CHANGE DATA.
355500

```

211300 DROPLEY SIZE DISTRIBUTION METER E0009 11239
211300
211300 THIS UNIT WILL DETERMINE DROPLEY SIZE DISTRIBUTIONS WITHIN A
211300 CHAMBER UTILIZING OPTICAL TECHNIQUES
211300
264650 IR MICROSCOPE E0017 11111
264650
264650 DETERMINE SURFACE TEMPERATURE OF DROPLETS AND ICE CRYSTALS
264650
264650 SENSITIVITY.....+/-0.54F(+/-0.3C) AT 5F (-15C)
264650
BAUSCH AND LOMB INC. #B57-634 \$876 111
OBJECTIVES 4 FLAT FIELD
TYPE OBJECTIVES 5X, 0.8 NA; 10X, 0.25NA; 40X,
0.65 NA; 100X, 1.25 NA
ILLUMINATOR 5 STEP, LOW VOLTAGE TRANSFORMER
EYEPieces PAIRED 10X
BAUSCH-LOMB INC. #MICROSCOPE TYPE II 111
NUMERICAL APERTURE 0.75
MEASUREMENT MAGNIFICATION 6X
OBSERVATION MAGNIFICATION 50X
EFFECTIVE SLIT LENGTH 12MM
EFFECTIVE SAMPLE AREA 0.3MMX2MM
BARNES ENGINEERING CO. #BM-24 \$7950 111
READOUT DIRECT
TEMP. RESOLUTION 0.5C AT ROOM TEMP.
DETECTOR INDIUM ANTIMONIDE
RESPONSE 8 MICROSEC.
VIEWING SIMULTANEOUS W/MEASUREMENT
SPECTRAL BANDWIDTH 1.8 TO 5.5 MICRONS
TEMP. RANGE SCALE A: 15C TO 65C
SCALE B: 15C TO 165C
SCALE C: 0C TO 350C
OPTIONAL SCALE: TO 1500C
OBJECTIVE 0.05X TO 74X
355500 OPTICAL SENSORS E0013 121111
355500
355500 SENSOR FOR CONTROL OF PARTICLE CHARACTERISTICS AND OPTICAL
355500 ENERGY SOURCES
355500
424000 DIGITAL TAPE RECORDER E0085 11111
424000
424000 PROVIDES FOR RECORDING SENSOR DIGITAL OUTPUT DATA.
424000
424000 TAPE SIZE 1 INCH
424000 TRACKS 9
424000 TAPE SPEED VARIABLE
424000 DATA RATE 1.3E+06 BPS
424000 REEL SIZE 12 INCHES DIA
424000
APPEX CORPORATION #AR-200 121
BANDWIDTH 100 HZ TO 3125 HZ THROUGH
300 HZ TO 250 KHZ
TAPE WIDTH 0.5 INCH
TAPE SPEED 1.875 TO 60 IPS
RECORDING TIME 8 MINUTES TO 4 HRS AND 16 MINS
FORMAT DIGITAL
TRACKS 9 DIGITAL, 7 ANALOG
APPEX CORP #AR 1700 \$28,000 112
TRACKS 14
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KB/1/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 6.2E+10
BANDWIDTH 2 MB/SEC/T
TAPE LENGTH 9200 FEET
SARGAMO ELECTRIC #SABER 111 \$26,500 111
TRACKS 14
TAPE WIDTH 1 INCH
TAPE SPEEDS 8 SELECTABLE SPEEDS FROM 15/16
TO 120 IPS
FREQUENCY RESPONSE 400 HZ TO 2.0 MHZ
RECORDING RATE 600 KBPS AT 120 IPS SERIAL MODE
WEIGHT 100 POUNDS
DIGI-DATA CORPORATION #160C \$2650 111
TAPE SPEED 25, 18.75, 12.5 IPS
TRACKS 7 OR 9 TRACK
DATA DENSITY 1600 CPI PHASE ENCODED
200, 556, 800 CPI NRZI
TAPE 0.5 INCH, 1.5 MIL, 1200 FEET
DIGI-DATA CORPORATION #1700/POP-11/7-9 \$5250 111
TRACK NRZI
TAPE SPEED 45, 37.5, 25, 18.75, 12.5 IPS
TAPE 0.5 INCH, 1.5 MIL, 18P/ANSI
COMPATIBLE, 10.5 INCH REEL
TRACKS 7 OR 9
DATA DENSITY PHASE ENCODED
COMPATIBLE

ORIGINAL PAGE IS
OF POOR QUALITY

HEWLETT PACKARD	#7570H/7	\$ 4500	111
TAPE FORMAT	PCO, 556, 14 200 CPI APTI AND 1400 CPI PHASE-ENCODER		
CHANNELS	7 OR 9		
TAPE SPEED	10 TO 45 IPS		
TAPE	0.5 INCH, 1.5 MILS, 144/INCH COMPATIBLE		
BORG WARNER	#PERT		121
TRACKS	30		
TAPE SPEED	UP TO 1000 IPS		
TAPE WIDTH	1/2 INCH		
TAPE LENGTH	2400 FEET		
PACKING DENSITY	15 KPI/T		
SIGNAL/NOISE	24 DB		
DATA CAPACITY	1.3E+9		
BANDWIDTH	6-14 MR/SEC/T		
LEACH	#MTF 7000		311
TRACKS	12		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
TAPE LENGTH	9200 FEET		
PACKING DENSITY	16.7 KPI/T		
SIGNAL/NOISE	22 DB		
DATA CAPACITY	2.2E+10		
BANDWIDTH	2 MR/SEC/T		
HEWLETT PACKARD	#3555D	\$10200	111
BANDWIDTH	300 KHZ		
CHANNELS	7		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3555C	\$14700	111
BANDWIDTH	300 KHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3555A-011	\$23900	111
BANDWIDTH	500 HZ TO 2 MHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3555A/13065A/13063A	\$4795	111
CONFIGURATION	RACK MOUNTED		
TAPE SPEED	15, 3 AND 3/4, 15/16 IPS		
CHANNELS	4		
RECORDING FORMAT	FM		
PASSBAND	5 KHZ		
S/N RATIO	44 DB		
TRACKS	30		
HONEYWELL	#5600	\$9730	
PORTABLE TAPE RECORDER			
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 60IPS		
MAXIMUM BANDWIDTH(DIRECT)	100KHZ		
PACKING DENSITY	UP TO 600PI		
WEIGHT	70LBS (32KG)		
INPUT VOLTAGE	28VDC		
HONEYWELL	#56	\$17420	
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 240IPS		
REEL SIZE	16IN		
MAXIMUM BANDWIDTH(DIRECT)	2M HZ		
TAPE WIDTH	1/2 IN		
ANPEX CORP	#AR 700	\$29,180	511
TRACKS	14		
TAPE SPEED	60 IPS		
TAPE WIDTH	1 INCH		
RECORD MODE	DIRECT		
PACKING DENSITY	20 KPI/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	RE+10 BITS		
BANDWIDTH	1 MB/SEC/T		
620070 LIDAR TRANSMITTER		ECORR	11112
620070 PROVIDES SIGNAL SOURCE TO MEASURE CLOUD HEIGHTS AND AEROSOL DIS-			
620070 TRIBUTION.			
620070 WAVE LENGTH	6328 ANGSTROMS		
620070 ILLUMINATION MODES	RUBY, DOUBLED RUBY, DYE, DOUBLED DYE		
620070			
620070			
018700 RECEIVER, TELESCOPE		ECORR	11112
018700 RECEIVES TRANSMITTER RETURNED SIGNALS.			
018700 INSTANTANEOUS RV	0.0018 RAD		
018700 PRINTING ACCURACY	0.0007 RAD		
018700 PRINTING STABILITY	0.0007 RAD		
018700 PRINTING STABILITY RATE	0.0017 RAD PER SEC		
018700 PRINTING DURATION	7200 SEC		
018700			

424000 DIGITAL RECORDER EC091 11111
424000
424000 PROVIDES FOR RECORDING DIGITAL OUTPUT DATA.
424000
424000 TAPE SIZE 1 INCH
424000 TRACKS 9
424000 PEEL SIZE 12 INCHES
424000 TAPE SPEED VARIABLE
424000 DATA RATE 80 MBPS
424000

AMPEX CORPORATION #AF-200 121
BANDWIDTH 100 HZ TO 3125 HZ THROUGH
300 HZ TO 250 KHZ
TAPE WIDTH 0.5 INCH
TAPE SPEED 1.875 TO 60 IPS
RECORDING TIME 4 MINUTES TO 4 HRS AND 16 MINS
FORMAT DIGITAL
TRACKS 8 DIGITAL, 7 ANALOG

AMPEX CORP #AR 1700 \$28,000 112
TRACKS 14
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KB/1/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 6.2E+10
BANDWIDTH 2 MB/SEC/T
TAPE LENGTH 9200 FEET

SANGAMO ELECTRIC #SARER 111 \$26,500 111
TRACKS 14
TAPE WIDTH 1 INCH
TAPE SPEEDS 8 SELECTABLE SPEEDS FROM 15/16
TO 120 IPS
FREQUENCY RESPONSE 400 HZ TO 2.0 MHZ
RECORDING RATE 600 KBPS AT 120 IPS SERIAL MODE
WEIGHT 100 POUNDS

DIGI-DATA CORPORATION #1600 \$2650 111
TAPE SPEED 25, 18.75, 12.5 IPS
TRACKS 7 OR 9 TRACK
DATA DENSITY 1600 CPI PHASE ENCODED
200, 556, 800 CPI NRZI
TAPE 0.5 INCH, 1.5 MIL, 1200 FEET

DIGI-DATA CORPORATION #1700/PDP-11/7-9 \$5250 111
TRACK NRZI
TAPE SPEED 45, 37.5, 25, 18.75, 12.5 IPS
TAPE 0.5 INCH, 1.5 MIL, IBM/ANSI
COMPATIBLE, 10.5 INCH REEL
TRACKS 7 OR 9
DATA DENSITY PHASE ENCODED
COMPATIBLE

HEWLETT PACKARD #79708/C \$ 4600 111
TAPE FORMAT 800, 556, OR 200 CPI NRZI AND
1600 CPI PHASE-ENCODED
CHANNELS 7 OR 9
TAPE SPEED 10 TO 45 IPS
TAPE 0.5 INCH, 1.5 MILS, IBM/ANSI
COMPATIBLE

BORG WARNER #BERT 121
TRACKS 30
TAPE SPEED UP TO 1000 IPS
TAPE WIDTH 1/2 INCH
TAPE LENGTH 2400 FEET
PACKING DENSITY 15 KB/1/T
SIGNAL/NOISE 24 DB
DATA CAPACITY 1.3E+9
BANDWIDTH 6-15 MB/SEC/T

LEACH #MTR 7000 311
TRACKS 12
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
TAPE LENGTH 9200 FEET
PACKING DENSITY 16.7 KB/1/T
SIGNAL/NOISE 22 DB
DATA CAPACITY 2.2E+10
BANDWIDTH 2 MB/SEC/T

HEWLETT PACKARD #39550 \$10200 111
BANDWIDTH 300 KHZ
CHANNELS 7
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3955C \$14700 111
BANDWIDTH 300 KHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3950A-011 \$23800 111
BANDWIDTH 500 HZ TO 2 MHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

ORIGINAL PAGE IS
OF POOR QUALITY



```

HEWLETT PACKARD      #3560A/13065A/13063A 84796      111
CONFIGURATION ..... RACK MOUNTED
TAPE SPEED ..... 15, 3 AND 3/4, 15/16 IPS
CHANNELS ..... 4
RECORDING FORMAT ..... FM
PASSBAND ..... 5 KHZ
S/N RATIO ..... 48 DB
TRACKS ..... 30

HONEYWELL            #5600      89730
PORTABLE TAPE RECORDER
CHANNELS ..... 7
SELECTABLE TAPE SPEED RANGE ..... 15/16 TO 40IPS
MAXIMUM BANDWIDTH(DIRECT) ..... 300KHZ
PACKING DENSITY ..... UP TO 600PAI
WEIGHT ..... 70LBS (32KG)
INPUT VOLTAGE ..... 2RVDC

HONEYWELL            #96      817420
CHANNELS ..... 7
SELECTABLE TAPE SPEED RANGE ..... 15/16 TO 240IPS
REEL SIZE ..... 16IN
MAXIMUM BANDWIDTH(DIRECT) ..... 2H HZ
TAPE WIDTH ..... 1/2 IN

APPEX CORP           #AR 70C      829,180      511
TRACKS ..... 14
TAPE SPEED ..... 60 IPS
TAPE WIDTH ..... 1 INCH
RECORD MODE ..... DIRECT
PACKING DENSITY ..... 20 KB/1/7
SIGNAL/NOISE ..... 20 DB
DATA CAPACITY ..... 8E+10 BITS
BANDWIDTH ..... 1 MB/SEC/7

350500 CSCOILSCOPE      EC096      11111
350500
350500 MONITOR, MEASURE AND MAINTAIN ELECTRONIC EQUIPMENT OPERATION.
350500
350500 BANDWIDTH ..... 100 MHZ
350500 CHANNELS ..... 2
350500 SENSITIVITY ..... 10 MV/DIVISION
350500 TIME BASE ..... VARIABLE
350500 STORAGE CAPABILITY ..... NO
350500

TEKTRONIX            #485      84200      111
BANDWIDTH ..... 350 MHZ
RANGE, CALIBRATED ..... 5MV/DIV TO 5V/DIV IN 10 STEPS
UNCALIBRATED ..... CONT VAR TO 12.5V/DIV
RESPONSE TIME ..... LESS THAN 1 NS
CHANNELS ..... 1, ALTERNATE, CHOPPED, ADDED, X-Y
..... CHANNEL 2 (UP OR INVERTED)
DIMENSIONS ..... 16.7 X 52.4 X 30.5 CM
WEIGHT ..... 9.5 KG
POWER ..... 115/230V 48-640HZ 60W

TEKTRONIX INC        #7313/7418 PLUG-IN/ 84300      111
787C PLUG-IN
BANDWIDTH ..... DC TO 80 MHZ
CHANNELS ..... 2
DEFLECTION FACTOR ..... 5 MV/DIV TO 5 V/DIV
TIMEBASE ..... 2 NS/DIV TO 5 S/DIV
MAX VIEW TIME ..... 4 HOURS
MAX STORED WRITING SPEED ..... 5000 DIV/NS
TYPE OF STORAGE ..... SPLIT SCREEN BISTABLE

TEKTRONIX INC        #7904/7419 PLUG-IN 86050      111
787D PLUG-IN
BANDWIDTH ..... DC TO 500 MHZ
CHANNELS ..... 2
DEFLECTION FACTOR ..... 10 MV/DIV TO 1 V/DIV
TIMEBASE ..... 2 NS/DIV TO 5 S/DIV

TEKTRONIX INC        #8465      81800      111
BANDWIDTH ..... DC TO 100 MHZ
CHANNELS ..... 2
DEFLECTION FACTOR ..... 5 MV/DIV TO 5 V/DIV
TIME BASE ..... 0.01 MICROSEC/DIV TO 0.5 S/DIV

018700 TELESCOPE - TRACKING      REC      11135
018700
018700 USED FOR HIGH RESOLUTION VIEWING OF SPECIFIC TARGETS
018700
018700 INSTANTANEOUS FIELD OF VIEW...0.90DEG(12X)
.....0.00EG(16X)
018700 FIELD OF VIEW(CROSS-TRACK)...0.90FG AND 4.00EG
018700 POINTING ANGLE FROM NADIR
.....+70 TO -40DEG
018700 PITCH, ALONG TRACK).....+70 TO -40DEG
018700 POINTING ANGLE FROM NADIR
.....+/-75DEG
018700 TOTAL ANGULAR COVERAGE (CROSS-
.....+/-75DEG
018700 TRACK) FROM NADIR (160V/2+
018700 POINTING ANGLE).....+/-75DEG
018700 POINTING ACCURACY (1-SIGNAL)...0.10EG
018700

```



622010 VIEWER - WIDE ANGLE MED 11235
622010
622010 USED FOR LARGE AREA VIEWING AND ORIENTATION. INSTRUMENT SIMILAR
622010 TO THE WILD NF2 NAVIGATION SIGHT USED WITH THE WILD-HEERHUGG
622010 RC-10 METRIC CAMERA
622010
622010 INSTANTANEOUS FIELD OF VIEW...1100EG X 1100EG
622010 ...550EG X 550EG
622010 ...280EG X 280EG
622010 FIELD OF VIEW.....3600EG IN AZIMUTH
622010 POINTING ANGLE FROM NADIR
622010 (PITCH, ALONG-TRACK).....0 TO 600EG
622010 TOTAL ANGULAR COVERAGE (CROSS-
622010 TRACK) FROM NADIR (1FOV/2)*
622010 POINTING ANGLE).....*/-600EG
622010 POINTING ACCURACY (1-SIGMA)....20EG
622010
459700 LONG WAVE INFRARED SPECTROMETER MEC 11111
459700
459700 USED FOR IDENTIFICATION OF TYPES OF ROCKS, SOILS AND SEDIMENTS
459700
459700 SPECTRAL RANGE.....0.4 TO 2.4 MICRONS
4597006.2 TO 15.5 MICRONS
459700 RADIOMETRIC CHANNEL RANGE.....10 TO 12 MICRONS
459700 INSTANTANEOUS FIELD OF VIEW...1 MILLI RADIAN DIA.
459700 POINTING ANGLE FROM NADIR
459700 (PITCH, ALONG TRACK).....+450EG TO -100EG
459700 POINTING ANGLE FROM NADIR
459700 (ROLL, CROSS TRACK).....*/-200EG
459700 POINTING ACCURACY.....0.30EG
459700
424000 MAGNETIC TAPE RECORDER MED 21111
424000
424000 ONBOARD DATA RECORDING OF SCIENTIFIC INSTRUMENT DATA
424000
424000 READ IN RATE.....100MBPS
424000 STORAGE CAPACITY/MISSION.....3E12 BITS
424000
AMPEX CORPORATION #AR-200 121
BANDWIDTH100 HZ TO 3125 HZ THROUGH
300 HZ TO 250 KHZ
TAPE WIDTH0.5 INCH
TAPE SPEED1.875 TO 60 IPS
RECORDING TIME8 MINUTES TO 4 HRS AND 15 MINS
FORMATDIGITAL
TRACKS8 DIGITAL, 7 ANALOG
AMPEX CORP #AR 1700 \$28,000 112
TRACKS14
TAPE SPEED120 IPS
TAPE WIDTH1 INCH
RECORDING MODE.....DIRECT
PACKING DENSITY20 KB/1/T
SIGNAL/NOISE20 DB
DATA CAPACITY6.2E+10
BANDWIDTH2 MB/SEC/T
TAPE LENGTH.....9200 FEET
SANCAMO ELECTRIC #SABER 111 \$26,570 111
TRACKS14
TAPE WIDTH1 INCH
TAPE SPEEDS.....8 SELECTABLE SPEEDS FROM 15/16
TO 120 IPS
FREQUENCY RESPONSE400 HZ TO 2.0 MHZ
RECORDING RATE600 KBPS AT 120 IPS SERIAL MODE
WEIGHT100 POUNDS
DIGI-DATA CORPORATION #1600 \$2650 111
TAPE SPEED25, 18.75, 12.5 IPS
TRACKS7 OR 9 TRACK
DATA DENSITY1600 CPI PHASE ENCODED
200, 556, 400 CPI NRZI
TAPE0.5 INCH, 1.5 MIL, 1200 FEET
DIGI-DATA CORPORATION #1700/POP-11/T-9 \$5250 111
TRACK NRZI
TAPE SPEED45, 37.5, 25, 18.75, 12.5 IPS
TAPE0.5 INCH, 1.5 MIL, 18M/ANSI
COMPATIBLE, 10.5 INCH REEL
TRACKS7 OR 9
DATA DENSITYPHASE ENCODED
COMPATIBLE
HEWLETT PACKARD #7570B/C \$ 4600 111
TAPE FORMAT800, 556, OR 200 CPI NRZI AND
1600 CPI PHASE-ENCODED
CHANNELS7 OR 9
TAPE SPEED10 TO 45 IPS
TAPE0.5 INCH, 1.5 MILS, 18M/ANSI
COMPATIBLE
BORG WARNER #PERT 121
TRACKS30
TAPE SPEEDUP TO 1000 IPS
TAPE WIDTH1/2 INCH
TAPE LENGTH2400 FEET
PACKING DENSITY15 KB/1/T
SIGNAL/NOISE24 DB
DATA CAPACITY1.3E+9
BANDWIDTH6-15 MB/SEC/T

ORIGINAL PAGE IS
OF POOR QUALITY

LEACH	WMR TC00		511
TRACKS	12		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
TAPE LENGTH	9200 FEET		
PACKING DENSITY	16.7 KB/1/T		
SIGNAL/NOISE	22 DB		
DATA CAPACITY	2.2E+10		
BANDWIDTH	2 MC/SEC/T		
HEWLETT PACKARD	#35550	\$10200	111
BANDWIDTH	300 KHZ		
CHANNELS	7		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#35550	\$14700	111
BANDWIDTH	300 KHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3550A-011	\$23900	111
BANDWIDTH	500 HZ TO 2 MHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#360A/13065A/13063A	\$4796	111
CONFIGURATION	RACK MOUNTED		
TAPE SPEED	15, 3 AND 3/4, 15/16 IPS		
CHANNELS	4		
RECORDING FORMAT	FM		
PASSBAND	5 KHZ		
S/N RATIO	48 DB		
TRACKS	30		
HONEYWELL	#5600	\$9730	
PORTABLE TAPE RECORDER			
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 20IPS		
MAXIMUM BANDWIDTH(DIRECT)	300KHZ		
PACKING DENSITY	UP TO 600PBT		
WEIGHT	70LBS (32KG)		
INPUT VOLTAGE	28VDC		
HONEYWELL	#96	\$17420	
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 240IPS		
REEL SIZE	16IN		
MAXIMUM BANDWIDTH(DIRECT)	2M HZ		
TAPE WIDTH	1/2 IN		
AMPEX CORP	#46 700	\$29,180	511
TRACKS	14		
TAPE SPEED	60 IPS		
TAPE WIDTH	1 INCH		
RECORD MODE	DIRECT		
PACKING DENSITY	20 KB/1/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	RE+10 BITS		
BANDWIDTH	1 MB/SEC/T		
219700 DATA BUFFERS, FORMATTERS			MEQ 21111
219700			
219700	DATA FROM HIGH RESOLUTION WIDEBAND MULTISPECTRAL SCANNER		
219700			
219700	INPUT DATA RATE.....200MRPS		
219700	OUTPUT DATA RATE.....100MRPS		
219700	INPUT DATA DUTY CYCLE.....73%		
219700			
RESEARCH INC	#R12-3	\$3500	111
OUTPUT VOLTAGE	5 V FS		
OUTPUT CURRENT	1-5, 4-20, 10-50 MAMP		
OUTPUT IMPEDANCE	1.2 K OHM ON 5 V SCALE		
CONVERSION RESOLUTION	10 BITS		
PERCENT OUTPUT/UPDATE	99.7		
DROOP RATE	0.1 %/HR		
SETTLING TIME	50 MICROSEC		
UPDATE TIME	2.2 MSEC		
MAX SIG AND HOLD CHANNELS	32		
SIG AND HOLD CHANS PER CARD	4		
105550 MULTIPLEXER			MEQ 11111
105550			
105550	MULTIPLEX DATA FOR TRANSMISSION VIA SHUTTLE COMMUNICATION LINK		
105550			
TENNELEC	#TC500	\$6000	111
CONVERSION TIME	3MICRO SEC		
CAPACITY	2048 CHANNELS		
CONFIGURATION	COMPATIBLE WITH AIM		
TENNELEC	#TC #20	\$750	111
CHANNELS	8		
CHANNEL DWELL TIME	300NSEC		
CONFIGURATION	NIM COMPATIBLE		
RESEARCH INC	#R12-4A	\$2630	111
CODING LINEARITY	+/- 0.5 LSR		
ANALOG INPUT	-10.240 TO +10.235 V		
INPUT IMPEDANCE	100 MEG OHM		
INPUT OVERVOLTAGE PROTECTION	+/- 100 V		
FULL SCALE ACCURACY	+/- .025 PERCENT		
RESOLUTION	12 BITS BINARY		
CONVERSION SPEED	5 KHZ		
LOGIC LEVELS			
LOGIC 0	5 +/- 1.5 V		
LOGIC 1	0 V + .4V, - .0V		

ANALOGIC CORP #AN471C-RA \$275 111
SINGLE ENDED CHANNELS 1A
DIFFERENTIAL CHANNELS B
CHRR NO DR AT 1000 MHZ

RCHDE AND SCHWARZ #UCH 210/1 111
(221.5603.02)

ANALOG INPUT RANGE +12 TO -12 VOLTS
SAMPLING RATE > 10000 MEAS/SEC (POSITIVE)
..... > 17000 MEAS/SEC (NEGATIVE)
DATA OUTPUT 1-1 DECADES (BCD CODE PARALLEL)
..... 0-1200 PULSES (SERIAL)

104180 CRT DISPLAY #ED 61111
104180
104180 TWO CRT'S ARE LOCATED IN THREE CONSOLES FOR DATA DISPLAY OF
104180 ALPHANUMERICS AND GRAPHICS.
104180
104180 RESOLUTION.....1000 LINES
104180 CONTROL INPUT FORM.....DIGITAL
104180 SCAN.....FASTER
104180CONICAL
104180 PHOSPHOR CHARACTERISTICS.....FOUR-COLOR PENETRATION
104180 SCREEN SIZE.....16IN (40CM)
104180 CONTRAST.....8 GRAY LEVELS
104180

RESEARCH INC #3300 \$1540
FORMAT.....24 LINESX72 OR 80 CHARACTERS
.....24 LINESX40 CHARACTERS
.....12 LINESX72 OR 80 CHARACTERS
REFRESH RATE.....40 HZ
TRANSFER RATE.....110 TO 2400 BAUD
CHARACTER FORM.....5X7 DOT MATRIX
DIMENSIONS.....15-1/2WX13-1/2HX23-1/2D 1ACH
WEIGHT.....39LB (17.7KG)

TEKTRONIX INC #4002A \$10975 111
DISPLAY MEDIUM 11 IN DIA CRT
DISPLAY AREA 8.3IN.HDR,X6.1IN.VERT.
ALPHANUMERIC MODE:
FORMAT 39LINES OF 85 NORMAL CHARACT.
CHARACTER SET 96 UPPER & LOWER CASE
CHARACTER SIZE 70X90MILS(CAN BE TWICE SIZE)
CHARACTER GENERATION 7X9 DOT MATRIX
CURSOR PULSATING 7X9 MATRIX
GRAPHIC MODES LINEAR INTERPOLATE,INCREMENTAL
PLOT,POINT PLOT,1024X1024 AD-
RESSABLE POINTS,1024X768 VIEW-
ABLE POINTS
GRAPHIC INPUT MODE 1024(X),768(Y),JOYSTICK CONTROLLED,
CROSSHAIR CURSOR
ACCUMULATORS 4
CYCLE TIME 300 NANOSEC
CORE MEMORY 32 K

166000 MICROFICHE READER #ED 11111
166000
166000 READER CONCURRENTLY DISPLAYS SEQUENCES OF MICROFILM RECORDS OF
166000 MAPS, REFERENCE SIGNATURES AND OTHER MATERIALS NEEDED FOR
166000 REAL TIME EVALUATION OF SENSOR FORMATS SELECTED FOR CRT DISPLAY
166000
166000 STORAGE CAPACITY.....71,000PAGES
166000 RETRIEVAL SPEED.....4 SECONDS MAX
166000 INTERFACE.....STANDARD COMPUTER
166000

CALMA CO. #303M \$28000 111
FUNCTION OFF LINE DIGITIZING SYSTEM
..... WITH FILM PROJECTION,
..... PROVIDES CAPABILITY OF
..... DIGITIZING FILM DATA

101340 MULTIFUNCTION KEYBOARD #ED 31111
101340
101340 CONSOLE CONTROL KEYBOARD WITH FUNCTION KEYS LABELLED BY COMPUTER.
101340 A SUPPORTING SOFTWARE LOGIC TREE FOR DATA ENTRY AND DISPLAY CALL
101340 UP CAN BE ENTER AT ANY LEVEL.
101340

BUAKER-RAND CORP #2200 \$440 111
KEYS TYPEWRITER SET STD
PROGRAM ASSIST KEYS 16
INTERFACE CRT
EDITING KEYS FULL SET

102750 OPERATIONS CONTROL COMPUTER #ED 11111
102750
102750 THE COMPUTER SHALL PROVIDE THE FOLLOWING FUNCTIONS: COMMAND
102750 INTERPRETATION AND EXECUTION, ATTITUDE AND POINTING PROGRAM
102750 CONTROL AND COMPUTATION AND DATA MANAGEMENT SEQUENCING OPERATIONS
102750
102750 MEMORY SIZE.....14550 WORDS
102750 ADD/MULT EXECUTE TIME.....5/30 MICRO SECONDS
102750

VARIAN DATA MACHINES #520/1 111
MEMORY CYCLE TIME 1.5MICROSEC.
MEMORY EXPANDABLE,4096BYTES(BITS)TO
..... 32,768BYTES
REGISTERS 12
OPERAND PRECISION UP TO 32 BITS
ADDRESS REFERENCE 8,16,24,OR 32 BIT LEVEL
NOTE: HAS FULL RANGE INTERFACE
..... HARDWARE

ORIGINAL PAGE IS
OF POOR QUALITY

VARIAN DATA MACHINES	#9-A20/1	111
BASIC COMMANDS	OVER 100	
ADDRESSING MODES	6	
MAX. WORDS	32,768	
WORD LENGTH	16 OR 18 BIT	
REGISTERS	9	
NOTE:	HAS FULL RANGE INTERFACE HARDWARE	
CLARY DATACOMP SYSTEMS INC	#404	111
SIMULTANEOUS TERMINALS	16	
INPUT	TELETYPEWRITER KEYBOARD	
ACCUMULATOR	64 BIT	
REGISTERS	16 BIT INDEX (2 EACH)	
WORD LENGTH	16, 32, 48, 64 BITS	
MEMORY CAPACITY	1024 16 BIT WORDS, 4096 16 BIT WORDS, OR ADDITIONAL 4096 16 BIT WORDS, TOTAL 85536 WORDS	
GENERAL AUTOMATION INC.	#18-30	111
FUNCTION	SUPERVISE SMALLER COMPUTERS	
MEMORY	8K CORE	
CYCLE TIME	1.2 MICROSEC.	
GENERAL AUTOMATION INC.	#5PC-14	111
MEMORY	16K	
WORD LENGTH	16 BITS	
READ/WRITE CYCLE MEMORY TIME	800 TO 1440 NANOSEC.	
READ ONLY MEMORY	400 TO 720 NANOSEC.	
INPUT/OUTPUT TIME	1.6 TO 2.8 MICROSEC.	
DATA TRANSFER RATE	0.694E+04 TO 2.5E+04	
GENERAL AUTOMATION INC.	#5PC-12	111
MEMORY	4K TO 16K	
WORD LENGTH	8 BITS	
CYCLE TIME	2.16 MICROSEC	
STATED PROGRAM EXECUTION RATE	0.23E+06/SEC	
I/O TRANSFER RATE	0.4E+06/SEC	
REGISTERS	8 12 BIT REGISTERS	
ACCUMULATORS	4 12 BIT	
DIGITAL EQUIPMENT CORP	#DDP-8/8	111
DIGITAL EQUIPMENT CORP	#DDP-P/1	111
MEMORY	4096 CORE EXPANDABLE TO 32,768 WORDS	
WORD	12 BIT	
CYCLE TIME	1.5 MICROSEC.	
LATA GENERAL CORP	#NOVA 400	111
WORD LENGTH	16 BIT	
CYCLE TIME	400 NANOSEC.	
CORE MEMORY	32K MAX	
DATA GENERAL CORP	#SUPERNOVA 50	111
WORD LENGTH	16 BIT	
DATA GENERAL CORP	#NOVA 1200	111
WORD LENGTH	16 BIT	
CYCLE TIME	1700 NANOSEC.	
CORE MEMORY	32K MAX.	



EARTH AND OCEAN PHYSICS

055055 FILM RECORDER OP002 51121
055055
055055 PROVIDE PHOTOGRAPHIC RECORDINGS TO PROVIDE ALL WEATHER STEREO
055055 IMAGERY FOR TOPOGRAPHIC MAPPING, 3-D STRAIN FIELDS, EROSION, VOL-
055055 CANIC MOTION AND POST-GLACIAL UPLIFT.
055055
055055 FILM TYPE 70 MM
055055 FILM RATE 6 MM/SEC
055055 INSTANTANEOUS FOV 0.088 RAD
055055
WESTINGHOUSE DEFENSE/SPACE 8APQ-97CLASS.EQUIP. 322
NOTE: SUBSYSTEMS PROCESSOR
6 FILM RECORDER ONLY
603002 MULTIFREQUENCY PROPAGATION RECEIVER/PROCESSOR OP004 11134
603002 L ANTENNA SYSTEM CONTROL UNIT
603002
624010 MULTISPECTRAL SCANNER ELECTRONICS PACKAGE OP006 11134
624010
603002 IR TEMPERATURE PROFILE RADIOMETER CONTROL UNIT OP018 11134
603002
603002 CAMERA CONTROL OP020 11134
603002
603002 GIMPAL CONTROL OP020 11134
603002
603002 OPTICAL MONITOR CONTROL OP021 11134
603002
603002 CONTAMINATION MONITOR GUAGE CONTROL OP021 11134
603002
603002 MASS SPECTROMETER CONTROL OP021 11134
603002

ORIGINAL PAGE IS
OF POOR QUALITY



RESEARCHING APPLICATIONS

600001 FLUID SUPPLY SYSTEM SP001 11134
600001 THIS UNIT CONSISTS OF THREE ELEMENTS SUPPORTING ALL MATERIAL
600001 SCIENCE AND MANUFACTURING EXPERIMENTS. IT IS PART OF THE CORE
600001 SUPPLEMENT PRESENT ON ALL MS/MS PAYLOADS. THE THREE ELEMENTS
600001 ARE AN INERT GAS UNIT, AN OXIDIZING GAS UNIT, AND A REDUCING
600001 GAS UNIT.
600001
400300 DIGITAL CLOCK SP001 11111
400300
400300 DIGITAL TIME DISPLAY, PART OF DIGITAL PROCESS PROGRAMMER ASSY
400300
400300 DISPLAY 4 DIGITS
400300 DISPLAY TYPE LED
400300
ROHDE AND SCHWARTZ #CAD 1100.6597.911 \$5400 111
DISPLAY 6 DIGITS: HR, MIN, SEC
OUTPUT 1-OUT-OF-N AND/OR HCD CODE
ROHDE AND SCHWARTZ #HS 4284(110.3226.21) \$17,700 111
DISPLAY 9 DIGITS: DAY, HR, MIN, SEC
DATATRON #3350 \$1020 111
INPUT FREQUENCY 40MHz
OUTPUTS VISUAL, BCD
DISPLAY HRS, MIN, SEC
DATATRON 3030 \$7940 111
CHANNELS 32
TIMING RANGE 5NANOSEC TO 5MICROSEC
RESOLUTION 10E-8
104200 TESTER PROGRAMMER SP001 11111
104200
104200 TAKES COMPLEX FUNCTION AND TRANSLATES IT INTO A SIGNAL WHICH
104200 DRIVES EXPERIMENT CONTROLLER
104200
104200 ACCURACY $\pm 0.5\%$ OF SPAN
104200 REPEATABILITY $\pm 0.1\%$ OF SPAN
104200
RESEARCH INC #5300 \$495 111
FOLLOWING ACCURACY
STATIC $\pm 0.25\%$ PERCENT OF SCAN
DYNAMIC $\pm 0.2\%$ PERCENT STATIC TO $\pm 0.75\%$ PERCENT AT MAX SLEW VEL
FOLLOWING REPEATABILITY $\pm 0.05\%$ PERCENT OF SCAN
PROB DEATHAND $\pm 0.025\%$ MM
FOLLOWING RESOLUTION WILL RESOLVE SQUARE WAVES WITH
PLOTTED TIME-DURATION EQUAL TO
0.7 PERCENT OF TIME AXIS
PROGRAM CARD KXP IN MYLAR BASE CONDUCTIVE
CLATED CARD
603003 SIGNAL CONDITIONER SP001 11111
603003
603003 AMPLIFIES AND CONDITIONS LOW LEVEL ANALOG SIGNALS AS PART OF
603003 DATA ACQUISITION UNIT
603003
603003 CAPACITY UP TO 33 PAIRS OF SIGNAL WIRES
603003 OUTPUT VOLTAGE $\pm 10V$
603003 OUTPUT CURRENT 5mA
603003
HENLETT PACKARD #441A \$390 111
FREQUENCY RANGE 1 KHZ TO 150 MHZ
FREQUENCY RESPONSE ± 1 DB INTO 50 OHM LOAD
GAIN AT 500 KHZ 40 DB ± 0.5 DB OR 20 ± 1 DB
MAX INPUT 1 V RMS OR 2 V P-P PULSE
TEKTRONIC #TC 2027201 \$145/495 111
GAIN 2.5 TO 1000X
NON LINEARITY LESS THAN 0.05%
CONFIGURATION NIM COMPATIBLE
ROHDE AND SCHWARTZ #RATN 1100.0899.021 111
FREQUENCY RANGE 30 HZ TO 20 KHZ
INPUT VOLTAGE 0.25 - 25 VOLTS
MAX POWER 50 WATTS
RESEARCH INC #812-11-17 \$3200 111
CAPACITY TO 33 PAIRS OF SIG WIRES (16GA)
REFERENCE TEMP HEAT SINK CONTROLLED TO ± 0.1 DB
DEF C AT 65 DEG C
AMBIENT TEMP 0 TO 49 DEG C
LEAKAGE RESISTANCE 1000 MEG OHM
HEAT-UP TIME FROM 22 DEG C 50 MINUTES
SETTLING TIME 4 MS TO FINAL VALUE
HONEYWELL #ACCU DATA 117 \$1340 111
MULTICHANNEL WIDEBAND AMPLIFIER
CHANNELS 7
MOUNT BENCH
HONEYWELL #ACCU DATA 115 \$5370 111
HIGH VOLTAGE DC AMPLIFIER
INPUT SIGNAL VOLTAGE $\pm 50MV$ TO 1500V
OUTPUT SIGNAL VOLTAGE UP TO 2000V
CHANNELS 6



NEFF #122-222 \$610 111
GAIN..... 1 TO 2500
COMMON MODE REJECTION..... 140DB
DRIFT..... 0.5MICROVOLT/OC
STEP GAIN ACCURACY..... +/-0.01%
SWITCH SELECTABLE FILTER
BANDWIDTH..... 100KHZ
BACK/MODULE MOUNT

NEFF #126-621 \$320 111
AMPLIFIER MULTIPLEXER
EMPLOYES FET SWITCH
GAIN RANGE..... 0.2 TO 2500
COMMON MODE REJECTION..... 120DB
OUTPUT..... +/-10V AT 1G MILLIAMPERES

NEFF #127-1C2 \$1150
PROGRAMMABLE GAIN WITH SWITCH SELECTABLE OUTPUT FILTER
GAIN STEPS..... 1, 2, 5, 10, 20, 50, 100, 200,
500,1000
COMMON MODE REJECTION..... 140DB
FET SWITCHING
PACK/MODULE CONFIGURATION

TENNELEC #1C133 \$295 111
FET PREAMPLIFIER
CHARGE SENSITIVITY..... 10E12 VOLTS/COULOMB
INPUT..... NEGATIVE INPUT PULSE POLARITY
WITH INPUT CAPACITANCE >200PF
OUTPUT..... POSITIVE OUTPUT PULSE POLARITY
WITH OUTPUT IMPEDANCE 50 OHM
STABILITY..... 0.5%/V

532000 DIGITAL VOLTMETER SP001 11111
532000
532000 VOLTMETER FOR DATA ACQUISITION AND CONTROL UNIT
532000
532000 DC VOLT RANGE..... 100MV TO 100V
532000 AC VOLT RANGE..... 100MV TO 10V
532000 RESOLUTION..... 0.1% FULL SCALE
532000 ACCURACY..... +/-0.5%+1 DIGIT
532000

JOHN FLUKE MFG CO INC #B1255 \$1845 311
VOLTAGE RANGE +/- 1 V TO +/- 1000 V DC
1 V TO 1000 V AC
RESOLUTION +/- 0.01%
FREQUENCY RANGE DC, 30 HZ TO 20 KHZ

JOHN FLUKE MFG CO INC #B425A-01 \$5995 311
VOLTAGE RANGE +/- 0.1 TO +/- 1000 V DC
1 TO 1000 V AC
RESOLUTION 0.001%
FREQUENCY RANGE DC, 10 HZ TO 100 KHZ

HONEYWELL #V1-100-4 \$320 111
DC RANGE 100 MV, 1 V, 10 V, 100 V
AC RANGE 100 MV, 1 V, 10 V
RESOLUTION 0.1 PCT FS
ACCURACY +/- 0.2 PCT +/- 1 DIGIT (10 V,
100 V DC SCALES); +/- 0.5 PCT
+/- 1 DIGIT (ALL OTHERS)
OVERRANGE 100 PCT

603002 SET POINT CONTROLLER SP001 11111
603002
603002 ESTABLISHES CONTROLLER SET POINT. SET POINT CAN BE VARIED
603002 FROM PRE-PROGRAMMED INPUT
603002
603002 SIGNAL INPUT..... C - 50MV
603002 DATA OUTPUT..... 1 - 50MA, 0 - 5 VDC
603002

RESEARCH INC #R12-6A-24 \$2500 131
SIGNAL INPUT 0 - 50 MV
BACKUP AMPLIFIER CARDS 16 CARDS, 4 CHANNELS/CARD
TOTAL BACKUP CHANNELS/MODULE 64
SET POINT ADJUST 1.0% PCT PER STEP
WATCHDOG TRANSFER SIGNAL 3-5 V (HOLDING VOLTAGE)
0 V (TRANSFER VOLTAGE)

102700 DIGITAL PROCESS PROGRAMMER SP001 11111
102700
102700 PERFORM PREDETERMINED LOGIC ROUTINES SUCH AS EVENT SEQUENCING,
102700 PARAMETER COMMANDS OR REAL TIME SEQUENCING.
102700
102700 INPUT/OUTPUT STAGE SP001 11111
102700
102700 THE INPUT/OUTPUT UNIT, THROUGH THE USE OF PLUG IN CARDS, CAN
102700 ACCOMMODATE A WIDE RANGE OF PROCESS INPUTS AND OUTPUTS FOR MAXIMUM
102700 CONTROL VERSATILITY.
102700
102700 LOGIC LEVELS, LOGICAL "0"..... 0V
102700 LOGICAL "1"..... 5V
102700 SIGNAL INPUTS..... COMPATIBLE WITH PDS TYPE
COMPUTERS
102700 SIGNAL OUTPUTS..... 0 - 5V
102700 DATA OUTPUTS..... 12 BIT PARALLEL

ORIGINAL PAGE IS
OF POOR QUALITY



10130 ANALOGIC CORP #MP/9125/14A, \$1040 111
/14W
DIGITAL OUTPUT COMP DTL/RTL PARALLEL AND SERIAL
MONOTONICITY GUARANTEED
QUANTIZING ERROR +/- 0.5 LSR
RESOLUTION 8 TO 16 BITS

10130 ANALOGIC CORP #ANT200 \$1500 111
#PCAL IRRATION 6 MONTHS
TEMPERATURE RANGE 0 TO 70 C
CONVERSION ACCURACY 0.01% FS (+/- 10 V)
RESOLUTION 8 THRU 14 BINARY BITS
DATA INPUTS 2, 3, OR 4 ACD DIGITS
DATA DTL/RTL POSITIVE TRUE LOGIC ONE
ADDRESS LOGIC/LINE
UP TO 16 LINES (14 BINARY BITS
OR 4 ACD DIGITS)
ADDRESS 8 CHANNELS - 3 BITS
24 CHANNELS - 5 BITS

10130 CONTROL UNIT - OPERATOR SP001 11111
10130 THE OPERATOR'S CONTROL PANEL WILL HAVE THE CAPABILITY OF
10130 CONTROLLING ABOUT 1000 PROCESS LOOPS. DATA MUST BE READOUT ON A
10130 QUICK-LOOK BASIS. IT MUST PROVIDE A WAY TO ENTER NEW PARAMETERS
10130
10130 CHANNEL NUMBER 3 DIGIT NIXIE (2 EACH)
10130 DATA LEVEL 4 DIGIT NIXIE PLUS SIGN (2 EACH)
10130 LOGIC LEVELS, LOGICAL "0" 0V
10130 LOGIC LEVELS, LOGICAL "1" 5V
10130 DATA OUTPUT VARIABLE
10130

10130 RESEARCH INC #DPS-#12 \$21,700 111
COMPLETE DATA DISPLAY SYSTEM MADE UP OF CRT ALPHANUMERIC DISPLAY
REF-POPE-F MINK-COMPUTER, #3301 OPERATOR'S CONTROL CONSOLE,
PRINTER/CARD READER, #R12-4A ANALYZER AND #R12-13 INTERVIEWER
CHANNELS 32
MEMORY SIZE 8K WORDS

10130 RESEARCH INC #R12-1 \$7400 111
LOGIC LEVELS 0 V (LOGIC 0)
5 V (LOGIC 1)
SIGNAL INPUTS COMPATIBLE POPE TYPE COMPUTER
SIGNAL OUTPUTS STANDARD ZERO AND 5 V OUTPUT
INTERUPTS 12 LEVEL HARDWARE
INTERVAL TIMES 1 MICROSEC TO 255 SEC ELAPSED
CLOCK CRYSTAL FREQUENCY TIME COUNTER
100 KHZ
BUFFER CARD LOCATIONS 11
CHANNELS 23

10130 RESEARCH INC #R12-2 \$2400 111
CHANNEL NUMBER 3 DIGIT NIXIE (2 EACH)
DATA LEVEL 4 DIGIT NIXIE PLUS SIGN (2 EACH)
SYSTEM STATUS 10 DUAL LAMP ANNUNCIATORS
ENGINEERING UNITS 10 DUAL LAMP ANNUNCIATORS
CONTROL PARAMETERS 10 DUAL LAMP ANNUNCIATORS
LOGIC LEVELS 0 V (LOGIC 0)
5 V (LOGIC 1)

10410 PRINTER SP001 11111
10410 PRINT RECORDED DATA IN THE FORM OF ALPHANUMERIC CHARACTERS
10410 IN HARD COPY FORM UPON REQUEST BY OPERATOR.
10410
10410 FORMAT PICA TYPE - STANDARD
CHARACTERS, NUMBERS SYMBOLS
10410

10410 TEK TRONIX INC #A601 \$3750 111
COPY SIZE 8.5 X 11 INCH
COPY TIME 18 SECONDS
RESOLUTION 4000 CHARACTER DISPLAY BASED ON
A 90 X 70 MIL MATRIX

10410 VERSATEC #1150P 112
FONT 7X9 DOT MATRIX (GOTHIC)
CHARACTERS/LINE 132
CHARACTER SPACING 12.5/IN
LINE SPACING 8.5 PER IN
SPEED 500 LINES/MINUTE
WRITING METHOD ELECTROSTATIC

10410 TELEPRINTER SP001 11111
10410 CHARACTER ORIENTED CRT TERMINAL. DATA IN THE LINE MEMORY IS
10410 CONVERTED, BY THE CHARACTER GENERATOR, INTO THE APPROPRIATE DOT
10410 PATTERN. THE DOTS ARE SHIFTED OUT OF THE SHIFT REGISTER TO FORM
10410 A VIDEO SIGNAL TO THE MONITOR.
10410
10410 SCREEN FORMAT 40, 72 OR 80 CHARACTERS PER LINE
10410 TRANSFER RATE 12 OR 24 LINES
10410 110 TO 2400 BAUD, 10 OR 11 BIT
CHARACTERS
10410 MODERATOR HALF OR FULL DUPLEX
10410

10410 VARIAN DATA MACHINES #ASH 33 6 ASH 75 111
SPEED 10 CHARACTERS/SEC.
PAPER CONTINUOUS ROLL 8.5 IN WIDE
TYPING LINE 85 CHARACTERS



RESEARCH INC	#R12-3301	\$1555	111
SCREEN FORMATS	22 LINES X 72 OR 80 CHARACTERS		
	24 LINES X 40 CHARACTERS		
CHARACTER FORM	5 X 7 DOT MATRIX		
TRANSFER RATE	110 TO 2400 BAUD, 10 OR 11 BIT CHARACTERS		
MODES	HALF OR FULL DUPLEX-SWITCHABLE LOCAL OR REMOTE		
BUNKER-RAND CORP	#2217-12CPT	\$950	111
CHARACTER RATE	240 PER SEC		
CHARACTER CAPACITY (MAX)	940		
CHARACTERS/LINE (MAX)	80		
LINES/DISPLAY (MAX)	24		
CHARACTERS REPERTOIRE	62 OR 92		
VIEWING AREA	8.75 X 6.25		
REFRESH AREA	54 FRAMES/SEC		
CHARACTER GENERATING METHOD	5X7 DOT MATRIX		
AN/GCC	#46		112
423000 DIGITAL TAPE RECORDER		SP001	11111
423000			
423000 PROVIDE DATA STORAGE FOR DIGITAL PROCESS PROGRAMMER			
423000			
423000 CHANNELS.....	15		
423000			
AMPEX CORPORATION	#AR-200		121
BANDWIDTH	100 HZ TO 3125 HZ THROUGH 300 HZ TO 250 KHZ		
TAPE WIDTH	0.5 INCH		
TAPE SPEED	1.875 TO 60 IPS		
RECORDING TIME	4 MINUTES TO 4 HRS AND 16 MINS		
FORMAT	DIGITAL		
TRACKS	8 DIGITAL, 7 ANALOG		
AMPEX CORP	#AR 17C0	\$28,000	112
TRACKS	14		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
RECORDING MODE	DIRECT		
PACKING DENSITY	20 KB/I/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	6.2E+10		
BANDWIDTH	2 MB/SEC/T		
TAPE LENGTH.....	9200 FEET		
SANGAMO ELECTRIC	#SABER 111	\$26,500	111
TRACKS	14		
TAPE WIDTH	1 INCH		
TAPE SPEEDS.....	8 SELECTABLE SPEEDS FROM 15/16 TO 120 IPS		
FREQUENCY RESPONSE	400 HZ TO 2.0 MHZ		
RECORDING RATE	60C KBPS AT 120 IPS SERIAL MODE		
WEIGHT	100 POUNDS		
DIGI-DATA CORPORATION	#1600	\$2650	111
TAPE SPEED	25, 18.75, 12.5 IPS		
TRACKS	7 OR 9 TRACK		
DATA DENSITY	1600 CPI PHASE ENCODED		
	200, 556, 800 CPI NRZI		
TAPE	0.5 INCH, 1.5 MIL, 1200 FEET		
DIGI-DATA CORPORATION	#17C0/PDP-11/7-9	\$5250	111
	TRACK NRZI		
TAPE SPEED	45, 37.5, 25, 18.75, 12.5 IPS		
TAPE	0.5 INCH, 1.5 MIL, 18M/ANST COMPATIBLE, 10.5 INCH REEL		
TRACKS	7 OR 9		
DATA DENSITY	PHASE ENCODED COMPATIBLE		
HEWLETT PACKARD	#7570B/C	\$ 4600	111
TAPE FORMAT	800, 556, OR 200 CPI NRZI AND 1600 CPI PHASE-ENCODED		
CHANNELS	7 OR 9		
TAPE SPEED	10 TO 45 IPS		
TAPE	0.5 INCH, 1.5 MILS, 18M/ANST COMPATIBLE		
BORG WARNER	#PERT		121
TRACKS	30		
TAPE SPEED	UP TO 1000 IPS		
TAPE WIDTH	1/2 INCH		
TAPE LENGTH	2400 FEET		
PACKING DENSITY	15 KB/I/T		
SIGNAL/NOISE	24 DB		
DATA CAPACITY	1.3E+9		
BANDWIDTH	4-15 MB/SEC/T		
LEACH	#MTR 7C00		311
TRACKS	12		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
TAPE LENGTH	9200 FEET		
PACKING DENSITY	16.7 KB/I/T		
SIGNAL/NOISE	22 DB		
DATA CAPACITY	2.2E+10		
BANDWIDTH	2 MB/SEC/T		

ORIGINAL PAGE IS
OF POOR QUALITY



HEWLETT PACKARD #39500 \$10200 111
BANDWIDTH 300 KHZ
CHANNELS 7
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #39550 \$14700 111
BANDWIDTH 300 KHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3950A-011 \$23800 111
BANDWIDTH 500 HZ TO 2 MHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3960A/13065A/13063A \$4794 111
CONFIGURATION RACK MOUNTED
TAPE SPEED 15, 3 AND 3/4, 15/16 IPS
CHANNELS 4
RECORDING FORMAT FM
PASSRND 4 KHZ
S/N RATIO 48 DB
TRACKS 30

HONEYWELL #560C \$9730
PORTABLE TAPE RECORDER
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 40 IPS
MAXIMUM BANDWIDTH(DIRECT) 300KHZ
PACKING DENSITY UP TO 600PPI
WEIGHT 70LBS (32KG)
INPUT VOLTAGE 28VDC

HONEYWELL #96 \$17420
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 240IPS
REFL SIZE 16IN
MAXIMUM BANDWIDTH(DIRECT) 2M HZ
TAPE WIDTH 1/2 IN

AMPEX CORP #AR 70C \$29,180 511
TRACKS 14
TAPE SPEED 60 IPS
TAPE WIDTH 1 INCH
RECORD MODE DIRECT
PACKING DENSITY 20 KB/1/7
SIGNAL/NOISE 20 DB
DATA CAPACITY 8E+10 BITS
BANDWIDTH 1 MB/SEC/T

603001 ANALOG (SCR) CONTROLLER SP001 11111
603001
603001 SILICON CONTROL RECTIFIER EMPLOYED AS A RELAY SWITCH
603001
603001 LOGIC LEVELS, LOGICAL"0" 0-0.6V
603001 LOGICAL"1" 3.0-5.0V
603001 RANGE OF CONTROL 0-100% IN 1 CYCLE INCREMENTS
603001

RESEARCH INC. #150SWITCH \$100 111
ONE REQUIRED PER CHANNEL

RESEARCH INC. #N1PC \$2800 111
LOGIC LEVELS 0-0.6 V (LOGIC 0)
3.0-5.0 V (LOGIC 1)
AMBIENT TEMP RANGE 4-43 DEG C
RANGE OF CONTROL 0-100 PCT IN 1 CYC INCREMENTS

103500 MULTIPLEXER A/D CONVERTER SP001 11111
103500
103500 CONVERT INPUT SIGNAL FROM ANALOG TO DIGITAL FORMAT FOR I/O UNIT
103500
103500 ANALOG INPUT -10.24 TO +10.235V
103500 FULL SCALE ACCURACY +/-0.025
103500 RESOLUTION 12 BITS BINARY
103500 CONVERSION SPEED 4 KHZ
103500 DATA OUTPUT 2'S COMPLEMENT

ANALOGIC CORP #MP2913A/14A \$3080 111
/15M
DIGITAL OUTPUT COMP DTL/TTL PARALLEL AND SERIAL
MONOTONICITY GUARANTEED
QUANTIZING ERROR +/- 0.5 LSB
RESOLUTION 8 TO 16 BITS

ANALOGIC CORP #AN720C \$1500 111
RECALIBRATION 6 MONTHS
TEMPERATURE RANGE 0 TO 70 C
CONVERSION ACCURACY 0.01% FS (+/- 10 V)
RESOLUTION 8 THRU 14 BINARY BITS
2, 3, OR 4 BCD DIGITS
DATA INPUTS DTL/TTL POSITIVE TRUE LOGIC ONE
LO-0/LINE
DATA UP TO 16 LINES (14 BINARY BITS
OR 4 BCD DIGITS)
ADDRESS 6 CHANNELS - 3 BITS
24 CHANNELS - 5 BITS

424000 TAPE INPUT UNIT SP001 11111
424000
424000 RECEIVES, STORES, TRANSMITS DATA OF OTHER EQUIPMENT
424000

424000 TAPE CAPACITY.....1200FT(360M), 0.5IN.(1.3CM) WIDE
424000 TAPE SPEED..... 25,18.75 AND 12.5 IN/S
424000(64,48 AND 32 CM/S)
424000 LONG TERM VARIATION..... \pm 1%
424000 DATA DENSITY.....1600 CPI (PHASE ENCODED)
424000200,556,800 CPI (NRZ1)

424000 STORAGE PERIPHERALS SP001 11111

219700 THIS UNIT ACCEPTS ADDRESSED BINARY DATA FROM THE DIGITAL
219700 PROCESSOR AND TRANSLATES THIS DATA INTO EQUIVALENT ANALOG FORM.

219700 MAXIMUM SIGNAL AND HOLD CHANNELS.32
219700 SIGNAL AND HOLD CHANNELS PER CARD..4

376400 AUTOMATIC PHOTOGRAPHIC PROCESSOR SPEC1 11111

376400 UNIT PROVIDES PROCESSED HOLOGRAMS AS PART OF THE ELECTRO-OPTICAL
376400 IMAGING SYSTEM

376400 CYCLE TIME.....BETWEEN 5 TO 20 SEC
376400 PLATE HOLDER SIZE.....4 IN(10CM) X 5 IN(13CM)
376400 EFFECTIVE APERTURE AT PHOTO
376400 PLATE.....3.25 IN.(8.5CM) X 4.0 IN (10CM)
376400 DATA OUTPUT.....PROCESSED HOLOGRAMS

OPTICS TECHNOLOGY INC #215 \$1995 111
POWER LEVEL..... 0.2 TO 1.0MW; 1.0 TO 2.0MW
EXPOSURE TIME..... 45 TO 180 SEC; 15 TO 60 SEC
MOUNT..... PENCH
FILM PLATE SIZE..... 4 X 5IN
HOLOGRAM TYPES..... GABOR AND FRESNEL

JCCON ENGR ASSOC #HS-2 \$8700 111

ACTIVE VIBRATION ISOLATION
LASER..... 20MW HE NE
SYSTEM COMPONENTS..... REAL TIME PLATE HOLDER, SHUTTER
BEAM STEERER, REFERENCE MIRRORS
SPATIAL FILTERS, VARIABLE BEAM
SPLITTER, ENCLOSURE

HRB-SINGER INC. RP-31 121

FILM SIZE..... 35MM
FILM LENGTH..... 250 FT.
PROCESSING RATE..... 0.5 TO 2 IN./SEC.

JCCON ENGINEERING ASSOC INC #AP-100 \$5250 111

FLUID HANDLING SYSTEM
MAXIMUM FLUID HEAD 3 M
HOSE FITTINGS 9.5 M
CYCLE TIME 5 - 20 SEC
TIMING 10 - 1000 SEC
HOSE LENGTHS UP TO 7.5 M
PLATE HOLDER
HOLDER SIZE 10 X 13 CM
PLATE THICKNESS SIZE 1 X 6.4 MM
EFFECTIVE APERTURE AT
PHOTO-PLATE 8.5 MM X 10 CM
OPTICAL CENTERLINE 11.4 CM

480000 CCTV CAMERA SP001 11131

480000 HIGH RESOLUTION TELEVISION CAMERA FOR VIEWING, RECORDING OR
480000 TRANSMITTING PICTURES OF MATERIAL PROCESSING EXPERIMENTS

480000 BANDWIDTH.....32MHZ
480000 VERTICAL SCAN LINES.....1125 LINES
480000 HORIZONTAL RESOLUTION.....1100 LINES
480000 VERTICAL SWEEP RATE.....25 OR 30 FRAMES PER SECOND
480000 SENSITIVITY.....100 IRE UNITS OUTPUT PRODUCED
480000 WITH 0.5 FOOTCANDLE ILLUMINATION

COMU ELECTRONICS INC. #6150 \$1943 111

SCANNING 1225 LINES/FRAME
INTERLACE 2:1
FRAME RATE 25 OR 30 FRAMES/SEC.
RESOLUTION(MAX) 1100 LINES,32MHZ
GEOMETRICAL DISTORTION 2 PCT MAX
SIGNAL TRANSMISSION DISTANCE ... 500 FEET
VERTICAL SWEEP RATES 50 OR 60 HZ
HORIZONTAL SWEEP RATES UP TO 75,750 HZ
SENSITIVITY 100 IRE UNITS OUTPUT PRODUCED
ON THE VIDICON FACEPLATE
WITH 0.5 FT-CANDLE HIGHLIGHT

603002 CCTV CAMERA CONTROL UNIT SP001 11131

603002 THIS UNIT WILL PERFORM THE FOLLOWING CONTROL FUNCTIONS: LENS
603002 FOCUS, LENS F STOP, CONTRAST BRIGHTNESS AND MODE. THE UNIT WILL
603002 WORK IN TWO MODES:(1)CONTINUOUS OBSERVATION MODE,WHERE ILLUMINA-
603002 TION IS CONTINUOUS, OR (2) PULSED OBSERVATION MODE USED WHEN
603002 ILLUMINATION IS PROVIDED BY A DYE LASER OR A NAKED FLASH LAMP.

603002 SCANNING.....INTERLACED 2:1, THROUGH 1125
603002 LINES PER FRAME, 25 OR 30 FRAMES
603002 PER SECOND
603002 HORIZONTAL SCAN FREQUENCY CONTROL 36 TO 750 HZ
603002 VERTICAL SCAN FREQUENCY CONTROL 50 TO 60HZ

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

LOMU ELECTRONICS INC. #6576-000 \$3140 111
SCANNING INTERLACED 2 TO 1, THROUGH 1225
LINES/FRAME, 25 OR 30 FR/SEC
SCAN FAILURE PROTECTION BEAM AUTOMATICALLY TURNED OFF
CONTROLS HORIZONTAL SCAN FREQ .. TO 36,750 HZ
CONTROLS VERTICAL SCAN FREQ 50 OR 60 HZ

376500 FRAME STORAGE UNIT SP001 11111
376500
376500 THE FRAME STORAGE UNIT PROVIDES NONDESTRUCTIVE READOUT, ERASE,
376500 AND STORAGE ON ELECTRONIC COMMAND AND WORKS IN EITHER CONTINUOUS
376500 OR PULSED MODE. IT PROVIDES A REAL-TIME, NO PHOTOGRAPHIC
376500 PROCESSING, HIGH SPEED, AUTOMATIC LIGHT ADJUSTMENT CAMERA CAPABLE
376500 OF PROVIDING REAL-TIME MOVIES OR FROZEN STILLS
376500
VIDEO STORAGE AND DISPLAY SYSTEM PRODUCED BY HUGHES AIRCRAFT,
CONSISTS OF MODEL 639 SCAN CONVERSION MEMORY, MODEL MSC-1 SCAN
CONVERTER UNIT AND INTERFACES WITH TV MONITOR, HARD COPY RECORDER
AND A VIDEOCON.

HUGHES AIRCRAFT COMPANY #639 \$5900 111
RESOLUTION 1350 TV LINES/TARGET DIAMETER
2100 TV LINES, LIMITING RES.
CAPACITY TO 750,000 BITS
WRITING SPEED 30 MICROSEC/SCAN LINE OR
25 NANOSEC/DOT
FRAME 33 MS TO 5 PCT FOR A FULLY
INSCRIBED SQUARE
STORAGE 25 MINUTES CONTINUOUSLY READ;
1-10 S STORAGE WITH CONTROLLED
FADE; OPTIONAL LONG STORAGE IN
BISTABLE OPERATION
GRAY SCALE LEVEL 10
ZOOM MAGNIFICATION TO 6 DIAMETERS
OVER ANY SELECTED PORTION OF
THE STORED IMAGE

HUGHES AIRCRAFT #MSC-1 \$2500 111
INPUT COMMANDS RASTER ERASE, SELECTIVE ERASE,
RASTER WRITE, GRAPHICS WRITE,
PRIME
VIDEO OUTPUT NORMAL AND INVERTED
OUTPUT BANDWIDTH 20HZ TO 10MHZ (FOR 525 LINE TV)
SIGNAL-TO-NOISE RATIO 35DB, TYPICAL
AMPLITUDE OV OFF, +4V ON

484000 TV MONITOR SP001 11111
484000
484000 THE MONITOR WILL DISPLAY CCTV IMAGES AS PART OF THE ELECTO-
484000 OPTICAL IMAGING SYSTEM.
484000
484000 SCANNING FREQUENCIES HORIZONTAL: 15-40 KHZ
..... VERTICAL: 15-60 FIELDS/SECOND
484000 VERTICAL SCAN LINES 1128 LINES
484000 LINEARITY 1X PICTURE HEIGHT
484000

CDHU INC. ELECTRONICS DIV #9600 \$5260
PICTURE SIZE 7-3/16 X 5-3/8 INCHES
HORIZONTAL SCAN 945 LINES, 50 OR 60 FIELDS/SEC
..... 873 LINES, 60 FIELDS/SEC
VIDEO BANDWIDTH 20 MHZ +/- 2 DB
POWER 25-28.5 VDC, 150 WATTS
359500 OSCILLOSCOPE SP001 11111
359500
359500 DISPLAY OF THE PRESENCE AND/OR NATURE AND FORM OF OSCILLATIONS
359500 OR IRREGULARITIES OF AN ELECTRIC CURRENT. STORAGE CAPABILITY OF
359500 A CRT DISPLAY IS DESIRED FOR PROLONGED OBSERVATIONS.
359500
359500 BANDWIDTH DC TO 10MHZ
359500 SENSITIVITY INTERNAL: 0.2CM DEFLECTION TO
1 MHZ
359500 EXTERNAL: 250MV P-P TO 15V P-P
359500 VIDEO STORAGE TIME 1 HOUR
359500

TEKTRONIX INC #7904/7A19 PLUG-IN \$6050 111
#7B70 PLUG-IN
BANDWIDTH DC TO 500 MHZ
CHANNELS 2
DEFLECTION FACTOR 10 MV/DIV TO 1 V/DIV
TIMEBASE 2 NS/DIV TO 5 S/DIV

HEWLETT/PACKARD #AN/USM-281A \$3100 111
MEETS MIL-SPEC REQUIREMENTS
BANDWIDTH DC TO 3MHZ
SENSITIVITY 0.1V/DIV TO 1.0V/DIV
CHANNELS 2

HEWLETT/PACKARD #1800/1A11A \$2750 111
BANDWIDTH 4 OR 18 GHZ
DEFLECTION FACTOR 2 MV
CHANNELS 2
PHOTOGRAPHIC CAPABILITY YES

TEKTRONIX INC #9465 \$1800 111
BANDWIDTH DC TO 100 MHZ
CHANNELS 2
DEFLECTION FACTOR 5 MV/DIV TO 5 V/DIV
TIME BASE 0.01 MICROSEC/DIV TO 0.5 S/DIV



TEKTRONIX PARS 84200 111
RANGE, CALIBRATED..... 5MV/DIV TO 5V/DIV IN 10 STEPS
UNCALIBRATED..... CONT VAR TO 12.5V/DIV
RESPONSE TIME..... LESS THAN 1 NS
CHANNELS..... 1, ALTERNATE, CHOPPED, ADDED, X-Y
..... CHANNEL 2 (UP OR INVERTED)
WEIGHT..... 9.5 KG
DIMENSIONS..... 16.7 X 52.4 X 30.5 CM
POWER..... 115/230V 48-60HZ 60W
BANDWIDTH..... 350 MHZ

TEKTRONIX INC 8R561R/38A PLUG-IN/ 81840 111
38A PLUG-IN
BANDWIDTH DC TO 10 MHZ
CHANNELS 2
DEFLECTION FACTOR 10 MV/DIV TO 10 V/DIV
TIMERASE 50 NS/DIV TO 5 S/DIV
RISETIME 35 NS

HEWLETT-PACKARD #132A 81500 111
BANDWIDTH DC TO 500KHZ
SENSITIVITY 100 MICROVOLTS/CM
INDEPENDENT BEAMS 2
WEIGHT 43 POUNDS

048800 FLUID COOLING/REFRIGERATION UNIT SP001 11114
048800
048800 THIS UNIT MAINTAINS THE ELECTROPHORETIC SEPARATION SYSTEM, THE
048800 BUFFER SUPPLY AND SAMPLE AT A CONSTANT LOW TEMPERATURE. IT
048800 CONTROLS THE TEMPERATURE OF THE GAS ELIMINATION SYSTEM, ALSO
048800 IT COOLS (OR HEATS) THE BUFFER/SAMPLE SOLUTION(S) AND THE BUFFER
048800 USED IN THE ELECTRODE COMPARTMENTS OF THE COLUMNS.
048800
048800 CONTROLLED TEMPERATURE RANGE...-13F TO 113F (25 TO 45C)
048800 TEMPERATURE ACCURACY..... +/-1.0F (+/-0.56C)
048800 SYSTEM RECOVERY TIME..... 2 MINUTES
048800 MAXIMUM HEAT REJECTION..... 2KW THERMAL
048800

355500 LASER OPTICAL SCATTERING MONITOR SP001 21111
355500
355500 THE LASER OPTICAL SCATTERING MONITOR WILL PROVIDE A REAL TIME
355500 MEASURE OF LIGHT SCATTERING IN A SPECIMEN AT ANY DESIRED MEASUR-
355500 IAG ANGLE. THIS DATA WILL PROVIDE INFORMATION ON THE SIZE, SHAPE,
355500 ORIENTATION, INDEX OF REFRACTION, CONCENTRATION AND LOCATION OF
355500 OPTICAL SCATTERING CENTERS.
355500
355500 LASER TYPE.....LOW-POWER UV GAS DISCHARGE
355500 FILTER BANDWIDTH.....4E-5 IN.(0.1NM) AT HALF POWER
355500 SCATTERED LIGHT INTENSITY
355500 ACCURACY.....C.1%
355500 ANGULAR POSITION ACCURACY.....1 MIN EQ.3MRAD)
355500

459619 UV-VIS SPECTROMETER SP001 21111
459619
459619 MEASUREMENT OF SAMPLE CONSTITUENTS BASED ON THE UV AND VISIBLE
459619 SPECTRUM.
459619
459619 RANGE 150 TO 1000 NM
459619 WAVELENGTH/WAVENUMBER ACCURACY... +/- 0.1 NM FS
459619 RESOLUTION 0.05 NM
459619 RESPONSE TIME 2 SECONDS
459619 SCANNING SPEED 0.02 TO 2.0 NM/SEC
459619

VARIAN INSTRUMENTS #CLARY 17 828300 111
WAVELENGTH RANGE 186-2650 NM

VARIAN INSTRUMENTS #CLARY 118 816645 111
WAVELENGTH RANGE 185-800 NM

COLEMAN #46 84800 111
FREQUENCY RANGE 195 TO 800 MILLIMICRONS
RESOLUTION..... 0.2 MILLIMICRONS
READOUT DIGITAL
WAVELENGTH ACCURACY +/-0.5NM

COLEMAN-HITACHI #EPS-3T 811500 111
FREQUENCY RANGE 170 TO 2600 MILLIMICRONS
WAVELENGTH ACCURACY..... 0.4 MILLIMICRONS IN UV
1.5 MILLIMICRONS AT 300 TO 600
MILLIMICRONS
2.0 MILLIMICRONS AT 600 TO 700
MILLIMICRONS
8.0 MILLIMICRONS IN NIR

PERKIN ELMER INSTRUMENT DIV #323 815200
WAVELENGTH RANGE 185 TO 2500 NANOMETERS
MODES X-Y RECORDER
AUTOMATIC REPETITIVE SCAN

G K TURNER #350 81609 111
WAVELENGTH RANGE 210 TO 1000 NM
ACCESSORIES 330-030 IR ACCESSORIES
370-040 UV ACCESSORIES
OPTICS REFLECTION TYPE
WAVELENGTH ACCURACY +/- 2 NM
WAVELENGTH READOUT DIGITAL IN NM
PHOTOMETRIC METER ACCURACY +/- 0.5%

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

Wavelengths 13300 ACTA CV 13300 111
Wavelengths 160-1000 NM
Resolution 140-800 NM
Accuracy ± 1 NM
Resolution 0.5 NM
Response Time 5 OR 2 SECONDS
Scanning Speed 0.025 TO 4 NM/SEC (7 SPEEDS)
Power 120V 50/60 HZ 3A
Dimensions 61X21X24 INCHES

227800 DYE LASER/FLASH LAMP SP001 31111
227800
227800 THIS UNIT WILL BE USED FOR FREE RADICAL GENERATION, SURFACE
227800 DAMAGE THRESHOLD DETERMINATION, AND HOLOGRAPHIC MICROSCOPY
227800
227800 ACCURACY 0.01%
227800 WAVELENGTH 160-1000 NM
227800 VOLTAGE 0.1%
227800 CAPACITY 1.0%
227800 OUTPUT 1.0%
227800 SYNC. ACCURACY 100NS
227800 LUMINANCE 20%
227800 RANGE OF OPERATIONS
227800 ENERGY INPUT 500-5000J/PULSE (1.1-1.4WH/PULSE)
227800 1 AT LOW REPS, LARGE LAMP
227800 50-500J/PULSE (1.014-1.4WH/PULSE)
227800 AT HIGH REPS, SMALL LAMP
227800 ENERGY OUTPUT 1-10J/PULSE (3E-4 - 3E-3WH/PULSE)
227800 1 AT LOW REPS, DYE PUMP, FREQ.
227800 250MJ TO 7J/PULSE (1E-5 - 2E-3
227800 WH/PULSE) FREQUENCY DOUBLED
227800 100MJ TO 3J/PULSE (1E-5 TO 3E-4
227800 WH/PULSE) AT HIGH REPS AT DYE
227800 FUNDAMENTAL, 25 - 750MHZ
227800 2E-4 WH FREQUENCY DOUBLED
227800 REPLICATION RATE LARGE LAMP - SINGLE TO 10/MIN
227800 SMALL LAMP - SINGLE TO 20/MIN
227800 FREQUENCY RANGE FUNDAMENTAL 430 - 750MHZ
227800 TO 3E-5 (IN.) FREQUENCY DOUBLED
227800 215 - 375NM (0.9E-5 TO 1.5E-5
227800 IN.)
227800 FLASH DURATION 500 - 100NS
227800 MICROCAMERA 35MM (1.4 IN.)
227800 SINGLE FRAME TO 20FRAMES/SEC.
227800
227800 FLASH LAMP SP001 31111
227800
227800 THE FLASH LAMP HAS TWO FUNCTIONAL REQUIREMENTS: (A) TO PRODUCE
227800 FREE RADICALS, AND (B) AS A PUMP SOURCE FOR EXCITING THE DYE
227800 LASER
227800
227800 TYPE MERCURY-XENON
227800 OPERATING TEMPERATURE 25 TO 75F (20 TO 25C)
227800 WAVELENGTH 1.2E-5 TO 1.4E-5 (IN.) (100-600NM)
227800 ENERGY (PULSED) 1.5WH (15KJ)
227800 (MIN) 0.3WH (1KJ)

227800
267500 RETROCONSTRUCTION HIGH RESOLUTION HOLOGRAPHIC MICROSCOPE SP001 31235
267500
267500
130000 DIALYSIS UNIT SP001 11113
130000
130000
CXFORD #10 159 111
GLASS CONTAINER WITH MOTORIZED AGITATOR
SPECIMEN CAPACITY 4

415500 METERING PUMPS SP001 11111
415500
415500 REFERR AND SPECIMEN METERING PUMPS FOR CONTINUOUS FLOW
415500 ELECTROPHORETIC COLUMN
415500
415500 FLOW RATE 0.001ML/MIN TO 30ML/MIN
415500 PRESSURE RANGE 5.0 TO 14 PSI (400 - 1000MM H2O)
415500 REPRODUCIBILITY $\pm 0.5\%$

415500
CHROMATONIX INC 111
FLOW RATE 2.4 TO 600 ML/HOUR
REPRODUCIBILITY 0.02%
ABSOLUTE ACCURACY 0.2%

VWR SCIENTIFIC #220 15445-2021 5995 111
FLOW RATE 7 MICRO-LITER/MIN TO 30 MILLI-
LITER/MIN
REPRODUCIBILITY ± 0.5 PCT
BACK PRESSURE RANGE 4E+4 TO 10E+4 N/50 M

SAGE #355 6475 111
DISCRETE RATES INFINITE
MINIMUM VOLUME/MIN 0.0004ML
MAXIMUM VOLUME/MIN 140ML

SAGE #375 6495 111
DISCRETE RATES INFINITE
MINIMUM VOLUME/MIN 1ML/MIN
MAXIMUM VOLUME/MIN 1200ML/MIN

149374 CONTINUOUS FLOW ELECTROPHORETIC COLUMN WITH PUMP SP001 11114
149374

ORIGINAL PAGE IS
OF POOR QUALITY



140378 THIS UNIT SEPARATES BIOLOGICAL SAMPLES BY ELECTROPHORESIS.
140378 RELATIVELY LARGE SAMPLES, APPROXIMATELY 10 CC, WILL BE PROCESSED.
140378 SAMPLE AND BUFFER FLOW WILL BE CONSTANT AND SAMPLE FRACTION
140378 COLLECTION CONTINUOUS.
140378
140378 TOTAL BUFFER AND SAMPLE VOLUME.. 60CU.IN. (1000CC)
140378 OPERATING TEMPERATURE..... 14 TO 41 F (-10 TO 5C)
140378 OPERATING TEMPERATURE REGULATION +/- 2F (+/- 1C)
140378 SEPARATION PROFILE ACCURACY... +/- .004IN (+/- .1PP)
140378 POSITION MONITOR ACCURACY..... +/- .004IN (+/- .1MM)
140378 FRACTION OUTLET PORT LOCATION ACCURACY +/- .004IN (+/- .1PP)
140378 INPUT VOLTAGE..... 5000V
140378
BUCHLER #VFR 30134-C56 431 111
INTERNAL DIAMETER..... 13MM
LENGTH..... 140MM
140390 STATIONARY ELECTROPHORETIC COLUMNS SP001 11114
140390
140390 THE STATIONARY ELECTROPHORETIC COLUMNS HAVE TWO ANTICIPATED
140390 USAGES: (1) TO ESTABLISH PROCESSING AND OPERATING PARAMETERS FOR
140390 THE CONTINUOUS ELECTRORETIC COLUMNS, AND (2) FOR SMALL BATCH
140390 PROCESSING (APPROXIMATELY 0.1CC (1.006 IN³)) OF VARIOUS
140390 BIOLOGICAL COMPONENTS
140390
140390 TOTAL BUFFER AND SAMPLE VOLUME 1.5CU.IN. (25CC)
140390 SEPARATION VELOCITY ACCURACY.. +/- 0.04IN (+/- 1MM)
140390 POSITION MONITOR ACCURACY..... +/- 0.04IN (+/- 1MM)
140390 OPERATING TEMPERATURE RANGE... 14 TO 41F (-10 TO 5C)
140390 TEMPERATURE TOLERANCE..... +/- 2F (+/- 1C)
140390
004400 GAS ELIMINATION SYSTEM SF001 11114
004400
004400 REMOVES GASSES PRODUCED BY ELECTROLYSIS IN THE ELECTRODE
004400 COMPARTMENTS OF THE ELECTROPHORETIC COLUMNS
004400
004400 OUTLET OXYGEN CONCENTRATION... 3% (BY VOLUME)
004400 OUTLET HYDROGEN CONCENTRATION.. 2% (BY VOLUME)
004400 MAX FLOW RATE..... TBD
004400 FLOW RATE CONSTANCY..... +/- 1%
004400 OPERATING TEMPERATURE..... 68 - 77 F (20 - 25 C)
004400
410000 BUFFER AND ELECTROLYTE SUPPLY TANKS SP001 11113
410000
374000 PH MONITOR SP001 21111
374000
374000 MEASURE ACIDITY AND ALKALINITY OF SOLUTIONS. USED IN GP LAB ALSO
374000
374000 PH RANGE..... 4.2 TO 12
374000 SELECTABLE RANGES..... ANY 2
374000 ANY 5
374000 ANY 10
374000 READING STABILITY..... 0.02 PH
374000 LIQUID TEMPERATURE RANGE..... 32 TO 212 F (0 TO 100 C)
374000
BECKMAN INSTRUMENTS #PHASAR 1 \$570
PH RANGE..... 0.00 TO 14.00
REPEATABILITY..... +/- 1 MV, +/- 0.01 PH
AUTO TEMP COMPENSATION..... 0 TO 100 C
DIMENSIONS 5-3/4X12-3/4X18-5/8 (INCHES)
WEIGHT 5 LB
BECKMAN INSTRUMENTS INC #0940 \$1000 111
RANGE 0 - 12 PH UNITS
SPANS 2, 5, OR 10 PH UNITS
STABILITY +/- 0.02 PH UNITS
AMBIENT TEMP RANGE -7 TO 50 DEG C
LIQUID TEMP RANGE 0 TO 100 DEG C
AMBIENT TEMP COEFFICIENT +/- 0.002 PH UNITS/DEG C
PERKIN ELMER CORP COLEMAN 28C \$265
RESOLUTION (METRIC 1V)
ACCURACY 0-14PH +/- 0.05PH
2% FULL SCALE
BECKMAN INSTRUMENTS #760C7 \$550 111
RANGE 0 TO 14PH
ACCURACY +/- .05PH (+/- .01PH (EXPANDED))
REPEATABILITY +/- .02PH (+/- .01PH (EXPANDED))
TEMP COMPENSATION 0 TO 100 C
RECORDER OUTPUT UP TO 200 MV
POWER 115/230V 50/60 HZ 15W
DIMENSIONS 11X11X9 INCHES
WEIGHT 9 LBS
BECKMAN INSTRUMENTS #100900 \$420 111
RANGE 0 TO 14PH
ACCURACY +/- .05PH FULL SCALE
REPEATABILITY +/- .02PH
TEMP COMPENSATION 0 TO 100 C
POWER 115V 50/60 HZ 2W
DIMENSIONS 11.3X9X5.38 INCHES
WEIGHT 6.25 LBS
RECORDER OUTPUT 10 TO 100 MV
177000 FRACTION COLLECTION SYSTEM SP001 11235
177000
177000 THIS APPARATUS WILL HAVE TWO FUNCTIONS: (1) TO COLLECT THE
177000 DESIRED FRACTIONS AFTER SEPARATION, AND (2) TO REMOVE THE EXCESS
177000 BUFFER SOLUTION AND REMAINING UNWANTED SAMPLES
177000

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

177000 PUMP TYPE.....POSITIVE DISPLACEMENT
177000 1.8 CU.IN/MIN (30 CC/MIN)
177000 PUMPING RATE ACCURACY.....+/-0.1%
177000 ALLOWABLE PUMP PULSATION.....C.CIPSI (+/-70 N/M2)
177000 TEMPERATURE REGULATION.....+/-1.0F (+/-0.5C)
177000
173000 FLOW METER SP001 11111
173000
173000 GAS FLOW METERING AND MONITORING DEVICE.
173000
173000 GAS TYPE OXYGEN, NITROGEN, HELIUM
173000 GAS FLOW RATE 0.4 TO 0.3 LB/HR
173000 0.3 TO 3.0 LB/HR
173000
TECHNOLOGY/VERSATRONICS INC #MFG-3 8970 111
OUTPUT(SINGLE ENDED) 0-5VDC LINEAR IN MASS FLOW
OUTPUT IMPEDANCE <10 OHM
ACCURACY(ABSOLUTE) +0R-2% OF READING 20 TO 100% FS
+0R-20MV FROM 0 TO 20% FS
+0R-0.4% FULL SCALE (FS)
LINEARITY +0R-0.1% OF READING
REPEATABILITY 1 ATMOSPHERE
CALIBRATION PRESSURE 55 TO 155F
STD CALIBRATION TEMP RANGE AIR
MEDIUM LESS THAN 0.07 SEC
TIME CONSTANT LESS THAN 1 IN H2O
PRESSURE DROP 1% TAUT BAND MIRROR SCALE
READOUT 0-100% FLOW RATE
AVAILABLE RANGE 0.12 TO 0.6 CU FT PER MIN
SIGNAL CONDITIONER 0 TO 40F
TRANSDUCER +45 TO 150F
285700 LYOPHILIZATION UNIT SP001 11111
285700
SCIENTIFIC PRODUCTS DIV #D7C2T-5 61390 111
CONFIGURATION BENCH EQUIP
TYPE FRONT LOADING, SINGLE SHELF
TEMPERATURE RANGE -40 TO 127 C (-40 TO 250 F)
VIRTIS #UNITRAP 61300 111
TEMPERATURE(MIN).....-65F(-55C)
DRUM SIZE.....8IN(21CM) DIA.X10IN(25.4CM) HIGH
CONFIGURATION..... FLOOR MOUNT, PORTABLE
VACUUM PUMP NOT INCLUDED
SARGENT WELCH SCIENTIFIC CO. #S-20081-80 63510 111
LOW TEMPERATURE -65F
CAPACITY 12 LITERS
VACUUM 5 MICRONS
PUMP CAPACITY 140 LITERS/MIN.
217400 RECIRCULATING FLUID INCUBATOR SP001 11112
217400
435300 DISSOLVED OXYGEN ANALYZER SP001 11111
435000
435000 THE UNIT WILL DETERMINE THE AMOUNT OF DISSOLVED OXYGEN PRESENT IN
435000 THE ELECTROLYTE AFTER IT HAS BEEN PROCESSED THROUGH THE GAS
435000 ELIMINATION SYSTEM. IT WILL MONITOR OXYGEN DURING A REACTION
435000 EXPERIMENT IT WILL ALSO FUNCTION AS A GASEOUS OXYGEN MONITOR
435000
435000 CONCENTRATION RANGE..... 0 TO TOTAL OXYGEN SATURATION
435000 ACCURACY..... +/- 1%
435000 RESPONSE TIME..... 90% ACTUAL IN 10SECONDS
435000 SAMPLE TEMPERATURE RANGE..... 32 TO 110F (0 TO 43C)
435000
BECKMAN INSTRUMENTS #1C0R00 FIELDAR 6475 111
RANGES:DISSOLVED O2 0-110-1010-25PPM
RANGES:GASEOUS O2 0-2.5E10-10E10-25%:0-100%
ACCURACY +/- 1.0% FS
REPEATABILITY5% FS
RESPONSE TIME 90% IN 10SEC GASEOUS O2
..... 90% IN 30SEC DISSOLVED O2
TEMPERATURE RANGE 0 TO 45 C
POWER 115V 50/60 HZ
DIMENSIONS 11.375X9X5.375INCHES
..... 28.8X22.9X13.6 CM
WEIGHT 7.5 LBS (3.3KG)
BECKMAN INSTRUMENTS INC #776 8700 111
CONCENTRATION RANGE ZERO TO TOTAL O2 SATURATION
UNITS READ PPM OR MG OF O2 PER LITER
TEMPERATURE RANGE 0 - 40 DEG C
ACCURACY +/- 1 PCT AT MEASURING TEMP
RESPONSE TIME 90 PCT OF ACTUAL CONCENTRATION
WITHIN 10 S
INC #G190R-25 9622.50 111
CONCENTRATION RANGE..... 0 TO 100%
ACCURACY..... +/-2%
TEMPERATURE RANGE..... 32 TO 122F (0 - 50C)
RESPONSE TIME..... 5% IN 15 SEC
410000 SPECIMEN/SAMPLE SUPPLY TANKS SP001 11113
410000
407000 HIGH VOLTAGE POWER CONDITIONER SP001 11112
407000
407000 PROVIDES DC VOLTAGE TO EITHER THE STATIONARY OR CONTINUOUS FLOW
407000 ELECTROPHORETIC COLUMNS. REGULATION MUST REMAIN CONSTANT OVER
407000 A RELATIVELY LARGE CURRENT DEMAND.



407000 OUTPUT VOLTAGE REGULATION.....+/- 0.1%
407000 VOLTAGE RANGE.....25 TO 2500 VDC
407000 READOUTS.....VOLTAGE LEVEL
407000 POWER OUTPUT TO COLUMNS.....20W SUSTAINED
407000 200W PEAK
407000
JOHN FLUKORFEG CO INC #4158 \$575 111
OUTPUT VOLTAGE0 TO 3100 VDC
CURRENT30 MA
REGULATION0.0005%
STABILITY+/- 0.01% PER DAY
TRANSIDYNE GENERAL CORP #2815 \$1475
VOLTAGE OUTPUT0 - 5 KV DC
CURRENT OUTPUT200 MA
RIPPLELESS THAN 1 MV RMS
ACCURACY+/- 0.1 PCT
SORENSEN #DCR60C0-425A \$1485 111
OUTPUT VOLTAGE.....60 TO 6000VOLTS
OUTPUT CURRENT (MAX).....0.287AMPS
REGULATION.....+/-0.075% OR +/-1.2V
CONFIGURATION.....RACK MOUNT
SIZE.....19X121/4X18IN
WEIGHT.....170 LB(BOXG)
229500 DARK FIELD ILLUMINATOR SP001 21111
229500
048700 REFRIGERATOR SP001 11111
048700
048700 THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES
048700
048700 TEMPERATURE.....23 - 68F (0 - 20C)
048700 LENGTH.....2FT (60CM)
048700 WIDTH.....1.8FT (54CM)
048700 HEIGHT.....2.8FT (84CM)
048700
MATHESON SCIENTIFIC #32C35-25 \$685 112
CAPACITY12.8 CU FT
TEMPERATURE-15 TO 27.8 C,
SHIELDED AREANONE
DIMENSIONS32X32X60 INCHES
EXPLOSION PROOF.....YES
COMPRESSOR0.7HP
POWER115V 50-60 HZ 3.5A
WEIGHT275 LBS
178200 FREEZER SP001 11111
178200
178200 THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES
178200
178200 CONTROLLED TEMPERATURE RANGE...-13F TO 113F (-25C TO 45C)
178200 LENGTH.....1FT (30CM)
178200 WIDTH.....1.8FT (54CM)
178200 HEIGHT.....2.8FT (84CM)
178200
REVCC BULT-1F5A \$1060 111
TEMP RANGE-7 TO -85 C
CAPACITY.....1.5 CU FT
DIMENSIONS.....34X24X37 INCHES
WEIGHT.....270 LBS (SHIPPING)
POWER115V 60 HZ 7A
CAPACITIES AVAILABLE.....1.5 TO 17 CU FT (VAR.PRICES)
108600 DEWAR SP001 11111
108600
108600 THIS UNIT WILL BE USED FOR COLD STORAGE OF BIOLOGICAL SAMPLES
108600
108600 TEMPERATURE.....-112F (-80C)
108600 LENGTH.....1FT (30CM)
108600 WIDTH.....1.8FT (54CM)
108600 HEIGHT.....2.8FT (84CM)
108600
CRYOGENIC ASSOCIATES #IR-90 \$1000
CAPACITY.....0.6L LIQ. N2
HOLDING TIME.....24 HOURS
DIMENSIONS.....4.25IN DIAM, 8.5IN LGNG
WEIGHT.....1.5LB(EMPTY),2.5LB(FULL)
410000 WASTE LIQUID TANK SP001 11113
410000
004400 MOLECULAR SIEVE SP001 41133
004400
004400 THE MOLECULAR SIEVE IS USED TO CLEAN GASSES GENERATED IN FURNACES
004400 AND EXHAUSTED TO SPACE. THIS UNIT IS ALSO USED IN THE FURNACE,
004400 GENERAL PURPOSE AND LEVITATION SUBELEMENTS
004400
UNION CARBIDE CORPORATION #TYPE 13X (POWDER) \$1.90/LB 114
(PELLETS) \$1.59/LB 114
TYPES AVAILABLEPOWDER
1/16 AND 1/8 INCH PELLETS
NOMINAL PORE DIAMETER10 ANGSTROM
POURED DENSITY(AVERAGE)38 LB/CU FT (PELLETS)
SETTLED BULK DENSITY42 LB/CU FT (PELLETS)
HYDRATED WET DENSITY1.95 G/CC (POWDER)
ACTIVATED DRY DENSITY1.53 G/CC (POWDER)
CRUSH STRENGTH9 LB (1/16 IN PELLETS)
24 LB (1/8 IN PELLETS)
CRYSTAL FORMOCTAHEDRAL (POWDER)
HEAT OF ABSORPTION1800 BTU/LB H2O

ORIGINAL PAGE IS
OF POOR QUALITY



410000 COOLANT SUPPLY TANK SP001 11113
410000
025010 ISOELECTRIC FOCUSING UNIT SF301 11114
025010
100000 HOT WALL TUBE FURNACE SP002 31111
100000
100000 THIS UNIT IS A GENERAL PURPOSE, HOT WALL HEATING DEVICE
100000 PROVIDING ACCURATE CONTROL OVER THE HOT ZONE TEMPERATURE WITH
100000 RESPECT TO A FLAT PROFILE OR TO A SPECIFIED GRADIENT PROFILE.
100000
100000 VACUUM LEVEL.....2E-8 PSI (1.9E-6N/M2)
100000 OPERATING TEMPERATURE.....2200F (1200C)
100000 MAX EXTERNAL SURFACE TEMP.....110F (45C)
100000 HOT ZONE DIMENSIONS.....2.5 IN (6CM) DIA. X 6IN (15CM) LONG
100000 MAXIMUM DEVIATION FROM FLAT
100000 PROFILE.....+/-4F (+/-2C)
100000
ASTEC INDUSTRIES INC #A222 81000 113
HOT ZONE DIMENSIONS 5 CM DIA. BY 28 CM LONG
MAX SUSTAINED TEMP 1200 DEG C
LINDBERG HEVI-DUTY #F8417 9355 111
MAX OPERATING TEMPERATURE 2200F (1200C)
CHAMBER SIZE 11N I.D. X 12IN LONG
POWER 57CW
LECO #521-3C0 81370 111
INDUCTION FURNACE
TEMPERATURE(MAX).....>3000F (>1700C)
POWER.....1.5KW
100000 HOT WALL FURNACE SP002 21111
100000
100000 THIS UNIT IS A GENERAL PURPOSE, HOT WALL HEATING DEVICE PROVIDING
100000 ACCURATE CONTROL OVER THE HOT ZONE TEMPERATURE WITH RESPECT TO
100000 A FLAT PROFILE OR TO A SPECIFIED GRADIENT PROFILE.
100000
100000 VACUUM LEVEL.....2E-8 PSI (1.9E-6N/M2)
100000 OPERATING TEMPERATURE.....3300F (1800C)
100000 EXTERNAL SURFACE TEMPERATURE.....110F (45C)
100000 HOT ZONE DIMENSIONS(MINI).....11N. (2.5CM) DIA. X 5IN. (13CM) LONG
100000 MAXIMUM DEVIATION FROM FLAT
100000 PROFILE.....+/-7F (+/-4C)
100000 VIEWPORTS.....2
100000
ASTEC INDUSTRIES INC #1100V - 1080-M1 83100 113
HOT ZONE DIMENSIONS 2.5 CM DIA. BY 15 CM LONG
MAX SUSTAINED TEMP 1800 DEG C
TIME TO TEMPERATURE 10 MINUTES MAX
HEATING ELEMENT MOLYBDENUM MESH WITH MOLYBDENUM
SHEET RADIATION SHIELDS
OPERATION ENVIRONMENTS VACUUM, INERT GAS, AND REDUCING
ATMOSPHERES
OBSERVATION SHUTTERED RADIAL VIEWPORT
ACCESS TOP AND/OR BOTTOM
COOLING WATER COOLED
LECO #VFR 30646-008 82065 111
TEMPERATURE(MAX)..... 3000F (1700C)
CHAMBER DIMENSIONS..... 4X4X5 1/4IN (10X17.8X21CM)
POWER..... 3KW
LINDBERG #5100C 81350 111
TEMPERATURE(MAX)..... 2700F (1500C)
CHAMBER DIMENSIONS..... 6X5X12IN (15.4X12.7X30.4CM)
POWER..... 5KW
100000 CHEST-GENERAL PURPOSE ENCLOSURE SP002 31113
100000
100000 THE PRIMARY FUNCTION OF THE ENCLOSURE IS TO ISOLATE THE HEATED
100000 SPECIMEN FROM THE LABORATORY ENVIRONMENT, BOTH THERMALLY AND
100000 ATMOSPHERICALLY, AND TO PROVIDE A MEANS TO STRUCTURALLY ATTACH
100000 THE VARIOUS HEATING ELEMENTS AND ACCESSORY ITEMS
100000
100000 VACUUM LEVEL.....2E-10PSI (1.5E-6N/M2)
100000 MAXIMUM INTERNAL TEMP.....5400F (3000C)
100000 MAX EXTERNAL SURFACE TEMP.....110F (45C)
100000 HOT ZONE DIMENSIONS(MINI).....11N.(2.5CM) DIA. X 4IN (10CM) LONG
100000 FULL LENGTH FRONT OPENING COOR
100000 VIEW PORTS.....FRONT, REAR, 2 SIDES, BOTTOM, TOP
100000 FLANGE OPENINGS.....TOP, BOTTOM, 2 EACH SIDE, 2 REAR
100000
SATREC SYSTEMS INC #VC-208 113
VACUUM CAPABILITY 1.3E-6 N/50 M
ENCLOSURE CONSTRUCTION DOUBLE-WALL WITH WATER COOLING
OF ALL SURFACES
COOR FULL LENGTH FRONT OPENING
VACUUM PUMPING CONNECTION 15 CM INSIDE DIA ON LOWER CHEST
INSTRUMENTATION LEAD-IN 10 CM INSIDE DIA ON UPPER CHEST
VACUUM SEALING METAL-TO-METAL SEALS
ASTEC (FURNACE) #1000A 84875 113
TEMPERATURE..... 5400F (3000C)
PRESSURE CAPABILITY..... VACUUM TO 15PSIG
POWER SUPPLY..... 12KVA
HOT ZONE DIMENSIONS..... 2.44IN DIA X 6 IN.
COOLING..... WATER



MATERIALS RESEARCH #PHZ-5A \$30,000 113
OPERATING TEMPERATURE..... UP TO 3000C
COMPATIBLE WITH RF INDUCTION OR ELECTRON BEAM
POWER..... 3KVA
VACUUM..... 10⁻⁶ TORR
PROCESS CHAMBER..... 12X12 HINGED PYREX BELL JAR

185450 DIRECTIONAL SOLIDIFICATION UNIT SP002 11114
185450
185450 THIS UNIT CONSISTS OF A FURNACE AND A CHILL UNIT. SPECIMENS A-E
185450 HEATED, MOVED THROUGH A THERMAL GRADIENT AREA AND THEN INTO A
185450 COOLING AREA.
185450
185450 CHILL UNIT DIMENSIONS.....0.6 IN(1.5CM) D. X 0.0 IN(20CM)
185450 LONG
185450 CHILL CAPACITY.....1KW THERMAL
185450 HEATING UNIT DIMENSIONS.....0.8 IN(2CM) D. X 1 IN(20CM) LONG
185450 HEATER CAPACITY.....1KW THERMAL
185450 TEMPERATURE GRADIENT ADJUST
185450 RANGE.....5 TO 500 F/IN (2 TO 200 C/CM)
185450 HEATER MAX TEMPERATURE.....2200 F(1600C)
185450 CHILL UNIT MIN TEMP.....-78 TO 77 F (-60 TO 25C)
185450 EXTERNAL SURFACE TEMP.....-60 TO 110 F (-45 TO 50C)
185450

003110 ACOUSTIC MIXING AND DISPERSAL UNIT SP002 21235
003110
628010 ELECTROMAGNETIC MIXING AND DISPERSAL UNIT SP002 21235
628010
419500 IR PYROMETER SP002 21111
419500
419500 THE PYROMETER OPTICALLY SENSES THE INTENSITY OF RADIANT ENERGY
419500 EMITTED FROM THE SAMPLE SURFACE
419500
419500 TEMPERATURE RANGE 104 TO 5400 F (-40 TO 3000 C)
419500 SPECTRAL RESPONSE.....0.70 - 0.97 MICROMETERS
419500 CALIBRATION ACCURACY.....+/- 1% FULL SCALE
419500 REPEATABILITY.....+/- 3% FULL SCALE

419500
PYROMETER INSTRUMENT CD #PHOTO: 11 \$3450 111
RANGES..... 1400 TO 2500 F (760C TO 1400C)
2400 TO 3400 F (1300 TO 1900C)
3200 TO 5400 F (1800 TO 3200C)
CAN OBTAIN 9000 F (5000C)
ACCURACY..... 1/4%
REPEATABILITY..... 0.005C/HR DRIFT AT 1045C

ANALOG DEVICES INC. #ADC-95 \$79 111
RESOLUTION 8 BITS
ERROR +/-0.2%
CONVERSION TIME 1 MILLISEC.

ANALOG DEVICES INC. #ADC-1207 \$305 111
RESOLUTION 12 BITS
ERROR +/-0.0125%
CONVERSION TIME 25 MICROSEC.

ANALOG DEVICES INC. #ADC-160 \$1350 111
RESOLUTION 16 BITS
ERROR +/-0.0015%
CONVERSION TIME 40 MICROSEC.

TENNELEC #TE900 \$6000 111
CONVERSION TIME..... 3 MICRO SEC
CAPACITY..... 2048 CHANNELS
CONFIGURATION..... COMPATIBLE WITH HIM

TENNELEC #TE 520 \$750 111
CHANNELS..... 8
CHANNEL DWELL TIME..... 300 NSEC
CONFIGURATION..... HIM COMPATIBLE

RCHDE AND SCHWARTZ #UCM 210/1 111
(221.4603.021)
ANALOG INPUT RANGE +12 TO -12 VOLTS
SAMPLING RATE > 10000 MEAS/SEC (POSITIVE)
> 17000 MEAS/SEC (NEGATIVE)
DATA OUTPUT 3+1 DECADES (A/D CODE PARALLEL)
0-1200 PULSES (SERIAL)

PYRO OPTICAL #87 \$495 111
TEMPERATURE RANGES..... 1400 TO 2200 F (760 TO 1210C)
1800 TO 3400 F (982 TO 1875C)
3200 TO 5200 F (1760 TO 2800C)

HY-CAL #P-R403-B \$650 111
ENERGY RANGE..... 1 TO 100 SOLAR CONSTANTS
VIEW ANGLES..... 5.7 TO 180 DEG
OUTPUT..... 5 OR 10 MILLIVOLTS/ SOLAR CONST.

186000 ZONE REFINER SP002 11112
186000
186000 THIS UNIT WILL HEAT, HOLD AND MANIPULATE A SPECIMEN FOR ZONE
186000 REFINING AND MOLTEN ZONE CRYSTAL GROWTH. IT WILL CAUSE A MOLTEN
186000 REGION TO TRAVERSE ALONG A LONG CYLINDRICAL ROD OF MATERIAL.
186000
186000 SPECIMEN SIZE..... 4 TO 12 IN (10 TO 30 CM) LONG
186000 0.1 TO 0.5 IN (0.25 TO 1.25 CM) DIA

ORIGINAL PAGE 13
OF POOR QUALITY



100000 SCANNING DISTANCES.....UP TO 12IN(30CM)
100000 SCANNING SPEED.....UP TO 20IN/MIN (50CM/MIN)
100000 SPECIMEN ROTATION.....C TO 40 RPM
100000 OPERATING TEMPERATURE.....575 TO 2900F (300 TO 1600C)
100000

MATERIALS RESEARCH CORP. #EBZ-6000 62900 112
USABLE SPECIMEN SIZE 10-30 CM LONG, 0.23-1.0 CM DIA
BEAM SOURCE WORK ACCELERATED, ANNULAR FILA-
MENT WITH BEAM FOCUSING SUPER-
STRUCTURE
POWER SUPPLY 5 KV, 400 MA DC
ZONE TRAVEL 23 CM DISPLACEMENT
2.5 - 25 CM/HR WITH AUTOMATIC
SHUT-OFF AT END OF SCAN

LINDBERG #54000 8675 111
THREE ZONE CONTROL
TEMPERATURE(MAX)..... 2200F (1200C)
ZONE DIMENSIONS..... 2-6IN(15.7CM)ENDS, 12IN(30.4CM)
CENTER SECTION

627010 FEED AND CRYSTAL HOLDER SPO02 21113
627010
435000 RESIDUAL GAS ANALYZER SPO02 21111
435000
435000 IDENTIFY TYPE AND QUANTITY OF GAS REMAINING IN A PROCESSING
435000 CHAMBER OR GENERATED DURING A PROCESS. THE ANALYZER IS ALSO USED
435000 IN THE LEVITATION SUBELEMENT. THE UNIT MUST HAVE A TUNABLE MASS
435000 UNIT SECTION FOR LEAK DETECTION. BE A QUADRUPOLE TYPE INSTRUMENT
435000 CAPABLE OF SEPARATING MASSES. THE INSTRUMENT MUST ACT AS A
435000 PARTIAL AND TOTAL PRESSURE ANALYZER AS WELL AS GAS ANALYZER.
435000
435000 MASS RANGE RESOLUTION..... UP TO 200AMU
435000 DETAILED SPECTRAL ANAL MASS RANGES..... 1-50 AND 40-65AMU
435000 SCOPE SCAN SPEED..... VARIABLE 3MS TO 300S
435000 DATA OUTPUT..... VISUAL ANALOG
435000 ANALOG
435000 DIGITAL
435000 USEFUL PRESSURE RANGE..... 1.0E-5 PSI (7.5E-2N/M2)
435000 TO 1.0E-11 PSI (7.5E-8N/M2)
435000 MASS SEPARATION CAPABILITY..... 1AMU
435000 MIN DETECTABLE PARTIAL PRESSURE..... 2E-12 PSI (1.5E-8N/M2)
435000 MIN DETECTABLE TOTAL PRESSURE..... 1E-10 PSI (7.5E-7N/M2)
435000

PEARLESS INSTRUMENT CO. INC. #210 85000 111
STD. GASES CO, CO2, H2O, H2S, HC, NH3,
NO, N2O, NO2, O3, SO2,
(ANY GAS MAY BE DEFINED
BY KNOWN ABSORPTION BAND
IN RANGE: 200 TO 10000 NM
0-100 OR 0-1000 PPM
FLOW 0.5 TO 1 CU.FT./MIN.
TEMPERATURE TO 200F
RESPONSE TIME 3 SEC. AT 90% READING
NOTE: EACH UNIT HANDLES 2 GASES
OPERATORS CONSOLE ALSO
REQD

VEECO #SPT-10 84000 111
SPECTRAL ANALYSIS 1-50 AND 40-65 AMU (DETAILED)
1-200 AMU (GENERAL)
RESOLUTION VALLEY < 10% OF THE SUM OF ADJ
PEAKS FOR 1 AMU APART THROUGH-
OUT 1-50 AND 40-65 AMU RANGES
MINIMUM DETECTABLE PRESSURE 1.3E-8 N/SQ M (PARTIAL)
6.5E-7 N/SQ M (TOTAL)
USEFUL PRESSURE RANGE 6.5E-8 N/SQ M
SCOPE SCAN SPEED 3 MS TO 300 SEC (VARIABLE)

BALZERS QMG 101 810000 111
TYPE..... QUADRUPOLE
PRESSURE/GAS ANALYZER..... YES
MASS RANGE..... 1 TO 400 AMU
MIN DETECTABLE PRESSURE..... E-14 TORR
MAX PERMISSIBLE PRESSURE..... E-4 TORR
RESOLUTION..... 10 PER CENT VALLEY
SCAN SPEED..... 25-2500MS PER MASS UNIT
WEIGHT CONTROL UNIT..... 25.8KG
WEIGHT ANALYZER..... 6.8KG

406500 RF INDUCTION PWR COND. (MIXING & DISPERSAL) SPO02 21111
406500
408000 LCW VOLT/HIGH AMP POWER CONDITIONER (10KW) SPO02 31111
408000

SOREASEN #DCR4C-500A 84200 311
OUTPUT VOLTAGE..... 0 - 40 VOLTS
INPUT VOLTAGE 414 TO 504 VOLTS
OUTPUT CURRENT..... 550 AMPS DC AT 30C
INPUT CURRENT..... 60 AMPS AC-3PHASE
TRANSIENT RESPONSE.....
CONSTANT VOLTAGE REGULATION..... +/-0.10% OR +/-15MV2
CONSTANT CURRENT REGULATION..... +/-1000MA
WEIGHT..... 809LBS (320KG)
SIZE..... 25X42X23IN
MOUNT..... FLOOR

419500 TWO - COLOR PYROMETER SPO02 21111
419500

104500	RESISTANCE HEATER (CONTACT)	SP002	21111
104500	HEATING ELEMENT WHICH HEATS ENCAPSULATED SAMPLES BY PASSING A		
104500	CURRENT THROUGH IT		
104500	POWER.....4KW SUSTAINED, 9KW PEAK		
104500	TEMPERATURE(MAX).....2900F (1600C)		
104500	ATMOSPHERE.....INERT GAS OR VACUUM		
104000	MICROWAVE HEATER	SP002	31235
104000	PROVIDE MICROWAVE ENERGY SOURCE FOR SPECIMEN HEATING		
104000	TEMPERATURE(MAX).....1400F (11000C)		
104000	POWER(MAX).....2KW		
104000	FREQUENCY RANGE(MAX).....UP TO 1010 HZ		
104000	OPERATIONAL ATMOSPHERE.....ANY BUT VACUUM		
106000	GRADIENT FURNACE	SP002	11113
106000	MECHANICAL MIXING AND DISPERSAL UNIT	SP002	11235
628020	THREE AXIS MANIPULATOR	SP002	21132
320000	PIEZOELECTRIC DRIVE	SP002	21114
627020	VACUUM/PRESSURE REGULATOR	SP002	31111
527000	THE UNIT REGULATES THE PRESSURE IN THE PROCESSING CHAMBERS BY		
527000	ADMITTING FLUID FROM A FLUID SUPPLY CONTAINER. THIS UNIT IS ALSO		
527000	EMPLOYED ON THE GENERAL PURPOSE AND LEVITATION SUBELEMENTS		
527000	PRESSURE-ABSOLUTE 0.1 TO 1000 PSIA		
417500	HIGH VACUUM PUMP	SP002	21111
417500	UNIT IS USED IN CONCERT WITH MOLECULAR SIEVE. IT IS USED WITH		
417500	THE SPACE VACUUM TO BLEED DOWN OVENS. THIS UNIT IS ALSO USED ON		
417500	THE LEVITATION SUBELEMENT.		
417500	PRESSURE-ABSOLUTE AMBIENT TO 0.1 PSIA		
417500	DUO-SEAL #1400	1136	111
	VACUUM LEVEL..... 0.05MICRON(0.0005TORR)		111
	CAPACITY..... 0.9CFM (25LITERS/MIN)		
	PRECISION #35	1135	111
	VACUUM LEVEL..... 15 MICRONS		
	CAPACITY..... 1.2CFM (35 LITERS/MIN)		
524500	VACUUM/PRESSURE MEASUREMENT UNIT	SP002	31111
524500	MEASURES PRESSURES IN PROCESS CHAMBERS. UNIT IS ALSO REQUIRED		
524500	FOR GENERAL PURPOSE AND LEVITATION SUBELEMENTS.		
524500	PRESSURE-ABSOLUTE AMBIENT TO 0.01 PSIA		
407500	HIGH VOLTAGE POWER CONDITIONER (17KV)	SP002	31111
407500	RF INDUCTION POWER CONDITIONER (2KHZ - 2MHZ)	SP002	21111
408500	SORENSEN #1020-30-L2	11220	111
	INCLUDES OUTPUT RIPLE ATTENUATOR FROM 2.0GRMS TO 0.01GRMS		
	OUTPUT VOLTAGE..... 0 TO 20 KV		
	INPUT VOLTAGE..... 105 TO 125 VAC		
	OUTPUT CURRENT..... 20 MILLIAMPS		
	INPUT CURRENT..... 10 AMPS AC, 1 PHASE		
	FREQUENCY..... 57 TO 63 HZ		
	WEIGHT..... 80LBS(32KG)		
	SIZE..... 19X 8 3/4X15 5/8IN		
	MOUNT..... RACK		
	SORENSEN #200	1900	111
	OUTPUT VOLTAGE..... 0 TO 30, 0 TO 15, 0 TO 7.5KV		
	INPUT VOLTAGE..... 105 TO 125 VAC 1PHASE		
	OUTPUT CURRENT..... 0-3, 0-6, 0-12MILLIAMPS		
	INPUT CURRENT..... 1.75AMPS AC, 1 PHASE		
	WEIGHT..... 55		
	SIZE..... 9X11X15 3/4 IN		
	CONFIGURATION..... PORTABLE		
186000	THERMAL OVEN	SP003	21111
186000	TENNY ENGINEERING INC #BLT TO SPEC		115
	PRESSURE VACUUM TO 1000 PSIG		
	TEMPERATURE -100 TO 600 F		
	FEAS PRECISION #R05A	1725	111
	TEMPERATURE(MAX) 660F (350C)		
	POWER 3100W		
	VOLUME 1.35CU.FT(0.04CU.M.)		
	MOUNT BENCH		
	NAPCO #620	1495	111
	TEMPERATURE(MAX) 340F(200C)		
	MOUNT BENCH		
	VOLUME 1.4CU.FT(0.1CU.M.)		

ORIGINAL PAGE IS
OF POOR QUALITY

BLUE M	#POM-143-A-1	8865	111
TEMPERATURE(MAX)	450F(230C)		
MOUNT	BENCH		
TR-ELCO	818	8395	111
TEMPERATURE(MAX)	440F(225C)		
POWER	1050W		
MOUNT	BENCH		
VOLUME	1.35CU.FT(0.04CU.M)		
499608 IR SPECTROMETER		SP003	11111
499608	INFRARED ANALYSIS OF CONSTITUENTS IN GASEOUS OR LIQUID STRFAMS.		
499608	RANGE	10 TO 100 MICRONS	
499608	RESOLUTION	+/- 0.5 MICRONS	
499608	RESPONSE TIME	2 SECONDS	
499608	BECKMAN INSTRUMENTS	84250	811500 111
	RANGE	4000-200CM-1(2.5-50MICRONS)	
	RESOLUTION	.5-1.4CM-1	
	ACCURACY	+/-1.0%	
	SCAN SPEED	1000 TO 2CM-1/MIN	
	SIZE	48X24X22 INCHES 122X61X56 CM	
	WEIGHT	210 LBS (95KG)	
419500 LASER PYROMETER		SP003	21235
419500	LASER PRECISION COPP	81120	8970
	ELEMENT DIAMETER	2.0MM	
	CURRENT RESPONSIBILITY	0.125 MICROAMPS/WATT	
	SATURATION POWER	800 MILLIWATTS	
	SATURATION ENERGY	6.4 MILLIJOULES	
074700 GAS CHROMATOGRAPH		SP003	21111
074700	THE FUNCTIONAL REQUIREMENTS OF THE EQUIPMENT ARE TO ANALYZE THE		
074700	ATMOSPHERIC COMPOSITION AND THE IMPURITY CONCENTRATION CONTAINED		
074700	WITHIN THE EXPERIMENTS' EQUIPMENT ENCLOSURES. THIS ELEMENT IS		
074700	ALSO USED ON THE LEVITATION SUBELEMENT.		
074700	FLOW RATE	10 ML/MINUTE	
074700	SENSITIVITY	0.1 PPM	
074700	BECKMAN INSTRUMENTS	86700	88000
	SENSITIVITY	1 +0.5 PPM	
	OPERATING TEMP.	55 +0.225C	
	TEMPERATURE CONTROL	+/- 0.05C	
BECKMAN	8GC-2A	82050	111
DETECTOR	THERMAL CONDUCTIVITY		
TEMPERATURE CONTROL RANGE	40C TO 240C		
CONTROL ACCURACY	+/-0.1		
CARLE INSTRUMENTS, INC.	8850C	81185	
TEMPERATURE RANGE	AMBIENT TO 200 DEG C		
DETECTOR	DUEL FEED SINGLE HYDROGEN FLAME		
SENSITIVITY	5E-12 AMPERES		
POWER REQUIREMENTS	115 VAC, 50 TO 60 HZ		
CARLE INSTRUMENTS INC	88001	88675	111
ACCURACY	+/- 0.5 PERCENT		
GASES DETECTED	CO, CO2, N2, AR, O2, CH4		
SAMPLE SIZE	1 MM		
REPETITIVE SAMPLING TIME	5 TO 60 MIN.		
CARLE INSTRUMENTS, INC.	88GC-311	82675	
CONTROL SENSITIVITY	< 0.001 DEG C		
TEMPERATURE RANGE	TO 399 DEG C (1 DEG INCREMENTS)		
HEAT CAPACITY-TOTAL	600 WATTS		
INLETS	DUAL ON-COLUMN		
IONIZATION DETECTORS	HYDROGEN FLAME		
VALVING	MICRO OR MINI VOLUME VALVES		
WEIGHT	53 POUNDS		
POWER REQUIREMENTS	115 VAC		
VARIAN AEROCGRAPH	81520	89200	112
COLUMN OVEN TEMP.	-99 TO 399 C		
TEMP STABILITY	BETTER THAN .1 PER CENT		
PLANE DET SENSITIVITY	.019 COULOMBS/GRAM		
MIN DETECTABILITY	5 E-12 GRAM/SEC		
NOISE	LESS THAN 3XE-14 A		
AMPLIFIER	SOLID STATE		
DIMENSIONS	20X18X22 INCHES		
BECKMAN INSTRUMENTS	86800	87350	111
CARRIER	H2		
DETECTOR	FLAME		
RANGE	0-1 PPM FULL SCALE		
ANALYSIS RATE	4 TO 6/MR (SELECTABLE)		
REPEATABILITY	.5% OF FS		
LINEARITY	1% OF FS		
ZERO DRIFT	ALOMATIC ZERO COMPENSATION		
DIMENSIONS	17X40X20 INCHES		
WEIGHT	250 LBS (SHIPPING)		
POWER	107-127V 50/60 HZ 500W		
884000 TIME LAPSE HIGH SPEED CAMERA		SP001	21111
884000	CAMERA ALSO USED IN LEVITATION SUBELEMENT		
884000			



J A PAUPER #500 1M - 9000 111
FILM SIZE 5 IN MS 33523
FORMAT 4.5 X 4.5 IN
SHUTTER 1/8000 TO 1/200 SEC
MODE PULSE OR AUTO CYCLE
FRAME RATE UP TO 5 FR/SEC
FOCAL LENGTH 24 IN
NO REO (LENSES) 2
RESOLUTION 104 L/MM AWAR
MAGAZINE 250 FT 5 IN FILM
MIL-SPEC / AIRBORNE

305000 NUCLEAR PARTICLE COUNTING UNIT SPO01 11111
305000
305000 CAPABILITY TO DISCRIMINATE AND COUNT ALPHA, BETA, GAMMA, X-RAYS
305000 AND K-ELECTRONS.
305000
305000 COUNTING RATE 3.5E+06 C/M
305000 RESOLUTION LOSS 1 E/3E+06 C/M
305000 PLATEAU CHARACTER
305000 ALPHA 900-1200 V @ 1 E/100V
305000 BETA-GAMMA 1700-1900 V @ 1 E/100 V
305000 PRECISION +/- 0.1 EFS
305000 TIMER 0.05 TO 2160 MINUTES
305000

SCINTIBLOC #1275E102 111
ENERGY RESOLUTION FOR CS137..... 9.5E OR BETTER
GAMMA RAY ENERGY..... 10KEV TO SEVERAL MEV

JCHASTON LABORATORIES #MM-21PAD:HV-4R 111
SYSTEM CONSIST OF FOCUSED MESH ELECTRON MULTIPLIER(#MM-21)
PREAMPLIFIER-AMPLIFIER-DISCRIMINATOR(#PAD)AND REGULATED HIGH
VOLTAGE POWER SUPPLY(HV-4R)
ELECTRON GAIN..... 10EA AJ 3KV
MAX COUNT RATE..... 10EP PER SEC
DETECTION EFFICIENCY..... ELECTRON- 90% AT 200E=V.
ICM- 90% AT 300GE.V.

HEWLETT PACKARD #5800A/5503A/ #1150/8650 112
#590A/5554 #1075/8360
/5502A #750
CONFIGURATION MULTI-UNIT BENCH ECLIF
RADIATION DETECTED ALPHA, BETA, X-RAY, GAMMA OR
NEUTRON
RATE MEASUREMENT RANGE 0 TO 1000 R/HR

414500 RESISTANCE THERMOMETER SPO03 11111
414500
414500 THERMOCOUPLES SPO03 11111
414500
104500 RESISTANCE HEATER (NON-CONTACT) SPO04 11111
104500
104500 TUBE TYPE HEATING UNIT FOR USE IN CHEST-GENERAL PURPOSE ENCLOSURE
104500 OR HOT WALL FURNACE
104500
104500 POWER.....4KW SUSTAINED, 9KW PEAK
104500 TEMPERATURE(MAX).....2900F (1600C)
104500 SPECIMEN SIZE.....4IN (10CM) LONG X 0.75 IN(2CM)
104500 DIA.
104500

103500 RF INDUCTION COILS SPO04 11111
103500
103500 PROVIDE NON-CONTACT HEATING
103500
103500 FREQUENCY RANGE.....2KHZ TO 2MHZ
103500 OPERATIONAL ATMOSPHERE.....ANY
103500

605003 ELECTRON BEAM SOURCE SPO04 11114
605003
605003 THE ELECTRON BEAM SOURCE WILL BE USED TO HEAT SAMPLES THAT
605003 REQUIRE NON-CONTACT HEATING. THE UNIT WILL BE USED IN THE CHEST-
605003 GENERAL PURPOSE ENCLOSURE
605003
605003 DELIVERABLE POWER.....3 - 4 KW
605003 INPUT POWER.....6KW
605003 TEMPERATURE.....4500F (2500C)
605003 MAXIMUM SAMPLE SIZE.....0.15 CU.IN. (2CC)
605003 OPERATING PRESSURE.....1E-10 PSI (7.5E-7A/P2)
605003

227000 LASER SOURCE SPO04 11111
227000
227000 PROVIDE HEAT TRANSMITTED TO A SPECIMEN BY LASER BEAM. CONCEPT
227000 PROVIDES HEATING WITHOUT CONTAMINATING SPECIMEN. LAIT IS USED
227000 IN CONJUNCTION WITH CHEST-GENERAL PURPOSE ENCLOSURE
227000
227000 POWER.....4FW
227000 COOLING CAPACITY.....7.5KW
227000 LASER TYPE.....ANY OF FOLLOWING FIVE:
227000CONTINUOUS WAVE HYDROGEN-FLUOR-
227000INE CHEMICAL LASER, IR OUTPUT
227000CONTINUOUS WAVE FORCED CONVEC-
227000TION CO2, IR OUTPUT
227000CONTINUOUS WAVE CO-C2-H2O-HE
227000CHEMICAL LASER, IR OUTPUT
227000PULSED TEA CO2, IR OUTPUT
227000PULSED TEA XENON, UV OUTPUT
227000

AMERICAN LASER CORP #90 #17950 111

ORIGINAL PAGE IS
OF POOR QUALITY



TYPE..... GAS (CO₂)
WAVELENGTH..... 10.6 MICRONS
POWER..... 300 TEM
BEAM DIVERGENCE..... 1.7 MRAD
RESONANT CAVITY DIMENSIONS..... 110 INCHES
Q SWITCH TYPE..... MECH
COOLING..... WATER

026010 MINIMUM B WITH RF HEATING SP004 11235
026010
009001 ELECTROMAGNETIC POSITIONING COILS AND DETECTOR SP004 11115
009001
009001 PROVIDE ALTERNATING ELECTROMAGNETIC FIELD FOR CONTACTLESS
009001 POSITION CONTROL.
009001
009001 MAGNETIC PRESSURE 1.4 PSI (10 KN/SQ M)
009001
009002 ELECTROSTATIC POSITIONING PROBES AND DETECTOR SP004 11115
009002
009002 ELECTROSTATIC FIELDS USED FOR POSITIONING EXPERIMENT SPECIMEN
009002 SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL.
009002
009002 CHAMBER PRESSURE > 0.14 PSI (1 KN/SQ M)
009002 < 2E-06 PSI (0.01 N/SQ M)
009002
009002 CONFIGURATION IN APEX OF TETRAHEDRON
009002
003210 ACOUSTIC TRANSDUCER AND DETECTOR. SP004 11111
003210
003210 USE OF SOUND PRESSURE WAVES TO POSITION EXPERIMENT SPECIMEN
003210 SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL.
003210
003210 CHAMBER PRESSURE > 0.14 PSI
003210 (1 KN/SQ M)
003210 < 2E-06 PSI
003210 (0.01 N/SQ M)
003210
003210 CONFIGURATION IN APEX OF TETRAHEDRON
003210
007001 GAS JET POSITIONING PROBES AND DETECTORS SP004 11115
007001
007001 USE OF GAS JETS TO POSITION SPECIMEN DURING EXPERIMENT FRC-
007001 CEDURES SUCH AS TO MAINTAIN CONTACTLESS POSITION CONTROL.
007001
007001 GAS FILTER 0E-06 INCHES (0.2 MICROMETERS)
007001 GASES HELIUM, HYDROGEN, OXYGEN
007001 CONFIGURATION JET PAIR IN APEX OF TETRAHEDRON
007001 GAS PRESSURE MINIMUM 0.014 PSIA (100 N/SQ M)
007001
052600 DIRECTIONAL CALORIMETER SP004 11113
052600
052600 THE DIRECTIONAL CALORIMETER WILL BE COMPOSED OF A PYROELECTRIC
052600 IR DETECTOR, DETECTOR PREAMPLIFIER, OPTICAL MODULATOR (LIGHT
052600 CHOPPER) AND OSMODULATOR (RECTIFIER). IT WILL BE USED TO
052600 DETERMINE THE BEGINNING OF SOLIDIFICATION BY OBSERVING THE
052600 ACCOMPANYING THERMAL ARREST AND WILL MEASURE THE RATE OF RELEASE
052600 OF THE HEAT OF FUSION.
052600

HY-CAL ENGINEERING 0C-119-A-10-003-003 0250 113
FLUX RANGE 0-100 W/SQ CM (SOURCE FLUX)
SPECTRAL RESPONSE 0.5-40 MICROM
RESPONSE RATE 0.25 SEC
FIELD-OF-VIEW 1.5 CM DIA, 0.25 RAD DIVERGENCE
ACCURACY +/- 3 PERCENT OF READING

LASER PRECISION CORP 0KT-2000/KT1-200 0190-900 113
SPECTRAL RESPONSE5 TO 75 MICROM
FREQUENCY RESPONSE 2 TO 500 HZ
FLUX RANGE 0 TO 10 W/SQ CM (INCIDENT)

052600 PYROLYTIC IR DETECTOR SPECTRAL
052600 RANGE..... 2E-5 TO 3E-3 INIO.9 TO 75
052600 MICROMETERS)
052600 INCIDENT FLUX RANGE..... 6.5W/IN2 (1 - 10W/CM2)
052600

320000 MECHANICAL SAMPLE PLACEMENT AND RETRIEVAL SF004 11113
320000
629020 LIQUID SYRINGE DISPENSER SP004 21113
629020
629010 INERTIAL INJECTOR SP004 11235
629010
410000 VACUUM CATCH TURE SP004 11235
410000
039010 SOLID SAMPLE STORAGE SP004 11112
039010

LIFE SCIENCES

099090 CAMERA - STILL FILM (CORE) LS900 11111
099090 STILL PHOTOGRAPHIC CAPABILITY WHICH IS SHARED WITH OTHER EXPERI-
099090 MENTS. CAMERA HAS INTERCHANGEABLE LENS AND IS COMPATIBLE WITH
099090 EQUIPMENT SUCH AS MICROSCOPES AND OSCILLOSCOPES.

099090 FILM FORMAT 35 MM
099090 LENS 50 MM, F1.2
099090 VIEW ANGLE 41 DEGREES
099090 LENS 35 MM, F1.4
099090 VIEW ANGLE 64 DEGREES
099090 SHUTTER SPEED 0 TO 1/500 SEC
099090 VIEWING THRU LENS
099090 FOCUS SPLIT IMAGE

HONEYWELL INC 05POTOMATIC F 063R 111
LENS 55 MM, F/1.8 TAKUMAR
35 MM, F/2.0 TO F/16 TAKUMAR
SHUTTER SPEEDS 1/1000 SEC TO TIMER
VIEWING THRU LENS

099090 CAMERA - STILL POLAROID (CORE) LS002 41111
099090 POLAROID PHOTOGRAPHIC CAPABILITY SHARED WITH MOST EXPERIMENTS
099090 ALSO ADAPTABLE WITH MICROSCOPE AND OSCILLOSCOPE INTERFACES.

099090 FILM FORMAT POLAROID 107 BLACK/WHITE
099090 3.25X4.25 IN (82.55X102.95 MM)
099090 LENS 56 MM, F3.5 TO F22
099090 SHUTTER SPEEDS 0 TO 1/60 SEC.
099090 RANGE FINDER SPLIT IMAGE
099090 VIEWING THRU LENS
099090 FIELD OF VIEW 40 DEGREES

TEKTRONIX INC 0C-12 0590 111
LENS 75 MM
STOP F/1.4 TO F/16
MAGNIFICATION 0.05
LENS SPEED 1 TO 1/100 SEC MECH
4 TO 1/60 SEC ELECT
FILM TYPE POLAROID

TEKTRONIX 0C30 0525
MAXIMUM RELATIVE APERTURE F/1.9
MAGNIFICATION 0.7 TO 1.5
RELATIVE SPEED 1.0
FIELD OF VIEW 3.15X3.93IN (80X100P)
FILM TYPE POLAROID

HEWLETT-PACKARD 0150 0420
MAXIMUM RELATIVE APERTURE F/3.5
MAGNIFICATION 1:0.05
SPEED 5 TO 1/60 SEC
LENS 75MM
FILM TYPE POLAROID

HEWLETT-PACKARD 0155A 01025
MAXIMUM RELATIVE APERTURE F/1.3
MAGNIFICATION 1:0.5
SPEED 8 TO 1/30 SEC
LENS 80MM
FILM TYPE POLAROID

425500 RECORDER - STRIP CHART (CORE) LS1013 11111
425500 HARD RECORD OF APPROPRIATE LIFE SCIENCES DATA.
425500 FREQUENCY RESPONSE DC TO 150 HZ
425500 SENSITIVITY 0.1 TO 50 MV/MM
425500 CHANNELS 8
425500 CHART SPEED 1 TO 500 MM/SEC
425500 WRITING METHOD THERMAL

HEWLETT PACKARD 071008R WITH OPTION 02295 111
019 AND 023
17500A PLUG-IN
CHANNELS 2 CHANNELS PLUS EVENT MARKER ON
BOTH SIDES
WRITING METHOD ELECTRIC
ACCURACY +/- 0.2 BPS
VOLTAGE SPAN 5 MV TO 100 V

HEWLETT PACKARD 077C20/0001A 02750 111
FREQUENCY RESPONSE DC TO 125 HZ
WRITING METHOD THERMAL
SENSITIVITY 5 MV/DIV
CHANNELS 2 CHANNELS TIMER/PARKER
VOLTAGE RANGE 5 MV TO 5 V PER DIV
ACCURACY +/- 1 %

HEWLETT PACKARD 077C20/0002A 00975

ORIGINAL PAGE 13
OF POOR QUALITY

CONFIGURATION	PACK MOUNTED		
CHART SPEED	0.25, 0.5, 1, 2.5, 5, 10, 25, 50, 100 MM/SEC		
FREQUENCY RESPONSE	150 HZ		
CHANNELS	8		
WRITING METHOD	THERMAL		
INPUT RANGES	1, 2, 5, 10, 20, 50, 100, 200, 500, 1000 MV/DIV		
ACCURACY	+/- 1%		
HONEYWELL	#1912	67700	111
PAPER WIDTH	12IN		
CHANNELS	24		
FREQUENCY RESPONSE	DC TO 25KHZ		
PAPER SPEED	200IPS		
SPEED VARIATION	15 FORWARD, 12 REVERSE SPEEDS		
MOUNT	RACK		
BECKMAN INSTRUMENTS	BR	61075*	111
BANDWIDTH	DC TO 130KHZ		
.....	5000 HZ WITH PR-500TAPE REC		
CHANNELS	1 TO 8 UP TO 24 SPEC CODEP		
REC FORMAT	THERMAL RECTILINEAR		
CHART SPEED	0.1 TO 25 CM/SEC OTHERS AVAIL		
EVENT MARKERS	2 (OPTION)		
DIMENSIONS	69.5X21.5X22.5 INCHES		
POWER	120V 50/60 HZ 250W		
HONEYWELL	#1956	67500	111
PARACOPY RECORDS FROM INFRARED, VIDEO AND SONIC SCANNERS, RADAR SYSTEMS AND REAL TIME SPECTRUM ANALYZERS			
RESOLUTION	1000 ELEMENTS PER SWEEP		
RECORDING RATE	4 MILLION OPS		
SPOT SIZE	0.005IN		
SWEEP SPEED	UP TO 20,000 LINES PER SEC		
FIBER OPTIC CRT			
HONEYWELL	#15CPA	63300	111
DATA CHANNELS	24		
EVENT CHANNELS	4		
FREQUENCY RESPONSE	DC TO 25KHZ		
PAPER SPEED	0.1 TO 120 IPS		
MOUNT	RACK		
PAPER WIDTH	12IN		
BECKMAN INSTRUMENTS	BRP	63820	111
FREQUENCY	DC TO 130 KHZ		
CHANNELS	EVENT, TIME, 2 DATA		
WRITING METHOD	THERMAL		
SENSITIVITY	1 MUVR/IN TO 5 V/IN		
RESPONSE TIME	4 NSEC FOR 10-90% OF 50 MM PEN TR		
CHART SPEED	1 TO 250 MM/SEC		
IRIG OUTPUT	OPTIONAL		
399600 OSCILLOSCOPE - PERSISTENT CRT (COPE)		LS976	11111
399600 ELECTRONIC MONITOR AND DATA MEASUREMENT CAPABILITY OF LIFE			
399600 SCIENCE EXPERIMENT EQUIPMENT WITH A SCREEN HAVING LONG TERM			
399600 RETENTION.			
399600 BANDWIDTH	DC TO 5 MHZ (MINIMUM)		
399600 CHANNELS	2		
399600 DEFLECTION FACTOR	1 MV TO 10 V PER DIVISION		
399600 PHOTOGRAPHIC CAPABILITY	YES		
399600 PERSISTENCE	UP TO 6 HOURS		
399600			
HEWLETT-PACKARD	#132A	61900	111
BANDWIDTH	DC TO 500KHZ		
SENSITIVITY	100 MICROVOLTS/CM		
INDEPENDENT BEAMS	2		
HEIGHT	65 POINTS		
HEWLETT PACKARD	#141B/1405A/1421A	62605	111
BANDWIDTH	5 MHZ		
DEFLECTION FACTOR	5 MV/DIV		
CHANNELS	2		
PERSISTENCE	0.2 SEC TO 31 MINUTE (VARIABLE)		
TERTRONIX INC	67313/741B PLUG-IN/	64300	111
7070 PLUG-IN			
BANDWIDTH	DC TO 80 MHZ		
CHANNELS	2		
DEFLECTION FACTOR	5 MV/DIV TO 5 V/DIV		
TIMEBASE	2 NS/DIV TO 5 S/DIV		
MAX VIEWING TIME	4 HOURS		
MAX STORED WRITING SPEED	5000 DIV/MS		
TYPE OF STORAGE	SPLIT SCREEN DISTABLE		
TERTRONIX INC	60504B/3A6 PLUG-IN/	62340	111
304 PLUG-IN			
BANDWIDTH	DC TO 10 MHZ		
CHANNELS	2		
DEFLECTION FACTOR	10 MV/DIV TO 10 V/DIV		
TIMEBASE	50 NS/DIV TO 5 S/DIV		
MAX VIEWING TIME	1 HOUR		
MAX STORED WRITING SPEED	500 DIV/MS		
TYPE OF STORAGE	SPLIT SCREEN DISTABLE		



263000 MICROSCOPE, DISSECTING LS1004 11111
263000
263000 PERMIT DETAIL DISSECTING OF BIOLOGICAL SPECIMENS SUCH AS PLANTS,
263000 ANIMAL TISSUE AND ORGANS
263000
263000 STEREO EYEPIECE
263000 ADJUSTABLE BRIGHT AND DARK FIELD ILLUMINATION
263000

BAUSCH - LOMB INC. #SSM-15 \$60 111
MAGNIFICATION 15X
IPO ADJUST. RANGE 45-75MM
MAX. SPECIMEN THICKNESS 0.875IN.
NUMERICAL APERTURE 0.042
PIN. RESOLVING POWER 60 LINES/MM
FIELD OF VIEW 4.0 MM
CONVERGENCE ANGLE(IPO) 10 DEG. AT 45MM
15 DEG. AT 75MM

OLYMPUS OPTICAL CO LTD #JPM-TR \$900 111
CONFIGURATION LOOSE FOOT
TYPE DARKFIELD STEPEO
MAIN BODY TRINOCULAR TUBE INCLUDING 45
DEG INCLINED BINOCULAR TUBE
PLUS VERTICAL PHOTO TUBE
MAGNIFICATION RANGE 6.3X TO 80X
WORKING DISTANCE 86 MM
APERTURE IRIS DIAPHRAGM 2 TO 40 MM (ADJUSTABLE)
CAMERA ATTACHMENT POLAROID, 4X5 SHEET FILM, 35 MM

AMERICAN OPTICAL #58M-03 \$420 111
EYEPIECES 10X, 15X WIDEFIELD
MAGNIFICATION 7, 10, 15, 20, 25, 30, 40, 50, 60, 80X
WORKING DISTANCE 1.4, 4 INCHES

263100 MICROSCOPE, COMPOUND LS1005 11111
263100
263100 GENERAL PURPOSE BINOCULAR MICROSCOPE FOR MICROSCOPIC STUDIES OF
263100 TISSUES

263100 MAGNIFICATION 10X TO 100X
263100 PHOTOGRAPHIC CAPABILITY POLAROID FILM
263100 LIGHTING DARK FIELD OR LIGHT FIELD AND
263100 PHASE CONTROL
263100

AMERICAN OPTICAL #SERIES V20 \$2035 111
EYEPIECES 10X (SUPPLIED) 20X 45X 60-1157
OBJECTIVES 10X, 20X, 45X, 100X
CAMERA ATTACHMENT YES (TRINOCULAR HEAD)
CONDENSER DARKFIELD
FLUORESCENCE TECHNIQUE UV ILLUMINATOR 50 W HG VAPOR
FILTER TURNEY WITH 3 FILTERS
AND OTHER ACCESSORIES
OBTAIN AO 11148

POLARIZER..... OBTAIN AO 11148

425000 RECORDER - VOICE (CORE) LS1006 11111
425000
425000 CAPABILITY OF RECORDING ORAL COMMENTS.

425000 CHANNELS 4
425000 BANDWIDTH 50 TO 20000 HZ
425000

SONY CORP #TC110 \$130 112
RECORDING TIME UP TO 90 MIN
RECORDING MEDIA MONAURAL CASSETTES, DUAL TRACK
VOICE INPUT MICROPHONE OR AUXILIARY INPUT

LEACH CORPORATION #LEM RECORDER 121
RECORDING TIME 10 HRS.
FUNCTION VOICE RECORD
TAPE CAPACITY 450FT
TAPE WIDTH 0.25 IN.
TRACKS 4
CHANNELS 1
FREQ. RESPONSE 300HZ TO 3.2KHZ AUDIO
SIGNAL/NOISE RATIO 35DB

507900 MICRODISSECTION TOOL KIT LS948 11111

507900
507900 STANDARD MEDICAL MICRODISSECTION KIT

507900
274100 MICROTOME (CORE) LS021 11111

274100
274100 CAPABILITY TO THINLY SLICE EXPERIMENT SPECIMEN.
274100
274100 SECTIONING 1 TO 50 MICRONS

AMERICAN OPTICAL #R20 \$1254 111
SECTION THICKNESS 1 TO 50 MICRONS
INPUT SPECIMEN BLOCK SIZE 32 X 27 MM

REICHERT #RM-12 \$5864 111
THICKNESS RANGE 1 TO 150 MICRONS
STEPS 0.2 MICRONS

LEITZ #1212-MINOT ROTARY \$1430 111
THICKNESS RANGE 1 TO 25 MICRONS
STEPS 1 TO 12 MICRONS

ORIGINAL PAGE 13
OF POOR QUALITY



109600 CRYOSTAT (CORE) LS022 11111
109600
109600 CAPABILITY OF STORING QUANTITIES OF CRYOGENIC LIQUID
109600
108600 QUANTITY 5 GALLONS (10.0019 CU M)
108600 TEMPERATURE -319 F (-195 C)
108600
CRYOGENIC ASSOCIATES #IR-90 \$1000
CAPACITY.....0.6L LIQ. N2
HOLDING TIME.....24 HOURS
DIMENSIONS.....4.25IN DIAM. 8.5IN LONG
WEIGHT.....1.5LB(EMPTY),2.5LB(FULL)
CRYOGENIC ASSOCIATES INC #SD-15 \$2852 112
CEWAR DIA.....17 IN.
LIQUID HELIUM CAP.....101 LITERS
LIQUID NITROGEN CAP.....26.5 LITERS
UNION CARBIDE CORP (LINDE) #LH-2180210-00001 \$225 111
CAPACITY.....20 LITERS
HEIGHT.....23 INCHES
DIAMETER.....15-3/16 INCHES
WEIGHT EMPTY.....23 LBS
FULL.....59 LBS
EVAPORATION RATE......22 LITER/DAY
MAX HOLDING TIME.....90 DAYS
UNION CARBIDE (LINDE DIV.) #0304-CC52(PART NO) \$1150 111
CAPACITY LIQ HE.....30 LITERS
EVAP. RATE/DAY......8%
DIMENSIONS.....44X18.75 INCHES
WEIGHT.....83 LB(MTY) 91.4LB(FULL)
507900 VETERINARY MEDICAL KIT (CORE) LS997 11111
507900
507900 INCLUDES OTOSCOPE/OPHTHALMOSCOPE, REFLEX HAMMER, HEMOGLOBINOMETER,
507900 Tourniquet, SYRINGES, NEEDLES, SCALPELS, HEMOSTATS, TWEEZERS, AND
507900 OTHER STANDARD ITEMS.
507900
186000 CVEN LS027 21111
186000
186000 CAPABILITY TO HEAT, MAINTAIN CONSTANT TEMPERATURE, AND DRY TEST
186000 EQUIPMENT, LIQUIDS AND TEST SPECIMEN.
186000
186000 TEMPERATURE RANGE 32 TO 450 F (10 TO 233 C)
186000 CAPACITY 1 CU FT (10.03 CU M)
186000
TENNY ENGINEERING INC #BLT TO SPEC 115
PRESSURE VACUUM TO 1000 PSIG
TEMPERATURE -100 TO 600 F
FREAS PRECISION #60SA \$725 111
TEMPERATURE(MAX) 660F(350C)
POWER 3100W
VOLUME 1.35CU.FT(0.04CU.M.)
MOUNT BENCH
NAPCC #620 \$495 111
TEMPERATURE(MAX) 390F(200C)
MOUNT BENCH
VOLUME 3.0CU.FT,10.1CU.M)
BLUE M #PDM-143-A-1 \$665 111
TEMPERATURE(MAX) 450F(230C)
MOUNT BENCH
TRELCO #IR \$395 111
TEMPERATURE(MAX).....440F(225C)
POWER.....1050W
MOUNT.....BENCH
VOLUME.....1.35CU.FT(0.04CU.M)
507900 HISTOLOGY KIT/SLIDE CABINET (CORE) LS996 11111
507900
507900 STAINING KIT CONTAINING PROBES, KNIVES, NEEDLES, SYRINGES, COVER
507900 SLIPS, LADLES PLUS OTHER NORMALLY INCLUDED ITEMS.
507900
486500 CONSTANT TEMPERATURE BLOCK (CORE) LS1015 21111
486500
486500 HEAT SINK MAINTAINING A CONSTANT TEMPERATURE FOR TEST TUBES AND
486500 VIALS.
486500
486500 TEMPERATURE 248+/-1.0 F (120+/-0.5 C)
486500
DCW DIAGNOSTICS #12/1A \$100 111
CONFIGURATION LOOSE EQUIP
TEMPERATURE RANGE AMBIENT TO 140C
CONTROL +/- 0.5 C
464600 SPECTROPHOTOMETER - IR (CORE) LS032 11111
464600
464600 PERFORM SPECTROPHOTOMETER ANALYSIS OF LIFE SCIENCES SPECIMEN IN
464600 THE INFRARED RANGE .
464600
464600 SPECTRAL RANGE 2.5 TO 12 MICRONS
464600 RESOLUTION 0.1 PERCENT OF WAVELENGTH
464600 SCAN TIME 5 SEC TO 50 HOUR
464600 SAMPLE TEMPERATURE AMBIENT TO 482 F (250 C)
464600
PERKIN ELMER INSTRUMENT DIV. #325 \$15200 111



Space Division
Rockwell International

WAVELENGTH RANGE 185 TO 2500 NANOMETERS
MODES X-Y RECORDER
AUTOMATIC REPETITIVE SCAN

G K TURNER ASSOCIATES #350 \$1609 111
WAVELENGTH RANGE 210NM TO 1000NM
ACCESSORIES 330-C30 IR ACCESSORY
330-040 UV ACCESSORY
OPTICS REFLECTION TYPE
WAVELENGTH ACCURACY +/- 2NM
WAVELENGTH READOUT DIGITAL IN NM
PHOTOMETRIC METER ACCURACY +/-0.5%

PERKIN-ELMER #221 \$10050 111
SPECTRAL RANGE 1.0 TO 15.5 MICRONS
RESOLUTION02 MICRONS AT 12 MICRONS
ACCURACY +/- .5 PER CENT
ORDINATE EXPANSION25X TO 40X
ABSCISSA EXPANSION25X TO 40X
POWER 105-125V 60 HZ 500W

BECKMAN INSTRUMENTS #IR-9 \$21000 111
SPECTRAL RANGE 2.5 TO 25 MICRONS
RESOLUTION25CM-1 AT 923CM-1
SCAN TIME 5MIN TO 18 DAYS
ACCURACY2% WITH CALIBRATION
HOT AND COLD SAMPLES YES

BECKMAN INSTRUMENTS #6800 IR-4 \$19740 111
SPECTRAL RANGE 1 TO 16 MICRONS
RESOLUTION01MICRON AT 10MICRONS
SCAN TIME 3MIN TO 24 HRS
ACCURACY015MICRON FROM 3.3-16MICRONS
HOT AND COLD SAMPLES YES

435000 GAS ANALYZER, AUTOMATIC LS993 11111
435000
435000 DETERMINE THE PARTIAL PRESSURES OF OXYGEN AND CARBON DIOXIDE
435000 DISSOLVED IN BLOOD SAMPLES AND TO DETERMINE THE HYDROGEN ION
435000 CONCENTRATION
435000
435000 (REQUIREMENTS UNSPECIFIED. COMMERCIAL UNIT USED AS MODEL)
435000 PH 6.000 TO 8.000
435000 CO2 PARTIAL PRESSURE 0 TO 200MM HG
435000 O2 PARTIAL PRESSURE 0 TO 2000MM HG
435000 ACCURACY:
435000 PH +/-0.003
435000 CO2 PARTIAL PRESSURE +/-0.5 MM HG
435000 O2 PARTIAL PRESSURE +/-1MM HG AT 200 MM HG O2
435000 +/-10MM HG AT 2000 MM HG O2
435000

INSTRUMENTATION LABORATORY #IL 313 \$4990
PH 6.0 TO 8.0
PCO2 0 TO 200 MMHG
PO2 0 TO 2000 MMHG
AUTOMATICALLY PROGRAMMED YES
SAMPLE TIME 45 SEC
PRECISION PH +/- .003
PCO2 +/-1MMHG
PO2 +/-10MMHG
SAMPLE SIZE 0.4ML (AUTOMATED)
POWER 115V 50/60 HZ

521700 URINE ANALYZER - AUTOMATIC (CORE) LS994 11111
521700
521700 CAPABILITY TO AUTOMATICALLY PERFORM URINE ANALYSES.
521700
459605 MASS SPECTROMETER LS994 11111
459605
459605 PROVIDE ANALYSIS OF UNKNOWN IONIZED GAS BY ATOMIC MASS NUMBER.
459605
459605 MASS RANGE 1 TO 400AMU
459605 SCANNING RANGE VARIABLE
459605 VARIABLE SCANNING RATE 50 MILLI SECOND TO 600 SEC/SCAN
459605 MAX OPERATING PRESSURE 1E-4TORR
459605 SELECTABLE CENTER MASS
459605 SELECTABLE SCANNING WIDTH

GRANVILLE-PHILLIPS CO #400 \$6975 112
CONFIGURATION BENCH EQUIP
MASS RANGE 1 TO 400 AMU
SCANNING RATE 50 MSEC TO 600 SEC (VARIABLE)
MAX OPERATING PRESSURE 1E-4 TORR

374000 PH METER LS1016 11111
374000
374000 METER MEASURES HYDROGEN ION CONCENTRATION OF SOLUTIONS
374000
374000 PH RANGE 0 - 14
374000 MEASUREMENT ACCURACY 0.02PH
374000

BECKMAN INSTRUMENTS #PHASAR 1 \$570
PH RANGE 0.00 TO 14.00
REPEATABILITY +/- 1 MV, +/- 0.01 PH
AUTO TEMP COMPENSATION 0 TO 100 C
DIMENSIONS 5-3/4HX12-3/4WX8-5/8D(11CH)
WEIGHT 5 LB

ORIGINAL PAGE IS
OF POOR QUALITY



HECKMAN INSTRUMENTS	#74	\$445	111
CONSTRUCTION	LOOSE EQUIPMENT		
PH RANGE	0 TO 14		
PH MEASUREMENT ACCURACY	+/- 0.02 PH		
PH MEASUREMENT REPEATABILITY	TO +/- 0.003 PH		
MILLIVOLT RANGE	0 TO 1400 MV		
MV MEASUREMENT ACCURACY	+/- 2MV		
MV MEASUREMENT REPEATABILITY	TO +/- 0.3 MV		
HECKMAN INSTRUMENTS	#123502	\$535	111
RANGE	0 TO 14PH		
ACCURACY	0.01		
RECORDING OUTPUT	1 TO 100 MV		
REPEATABILITY	+/-0.01PH		
DIMENSIONS	5.75X12.75X8.62 INCHES		
WEIGHT	5 LBS		
TEMP COMPENSATION	0 TO 100 C		
HECKMAN INSTRUMENTS	#74009	\$695	111
RANGE	0 TO 14 (ANY 1.4PH EXPANDED)		
ACCURACY	+/- .05PH +/- .007 (EXPANDED SCALE)		
REPEATABILITY	+/- .02PH +/- .002 (EXPANDED SCALE)		
TEMP COMPENSATION	0 TO 100 C		
RECORDING OUTPUT	0 TO 200 MV		
POWER	115/230V 50/60HZ 15W		
DIMENSIONS	11X11X9 INCHES		
WEIGHT	9 LBS (4KG)		
201000 ULTRASONIC CLEANER (CORE)		LS047	11111
201000			
201000 HIGH FREQUENCY ENERGY FOR CLEANING EXPERIMENT AND MAINTENANCE			
201000 EQUIPMENT			
201000			
201000 OPERATING FREQUENCY	20000 HZ		
201000 CAPACITY	1 CU FT		
201000			
PRAXSON	#2	\$55	111
TANK SIZE	3.5X3.5X2.5 IN (9X9X6.5 CM)		
CAPACITY	1 QUART		
PRAXSON	#32	\$225	111
TANK SIZE	6X11X6 IN (15.2X28X15.2 CM)		
CAPACITY	1 1/2 GALLON		
KAUSCH AND LOMB INC	#HALSONIC II	\$500	111
CAPACITY	1 GALLON		
TANK DIMS	7.25X9.125X3.875 IN #3		
HECKMAN INSTRUMENTS (ESC, ESS, DV1892-H)		\$885	111
FREQUENCY	50000 HZ		
CAPACITY	1 CU FT		
HEATING	UP TO 60 C		
532000		LS1012	11111
532000			
532000 PORTABLE RELIABLE MULTI-FUNCTION METER FOR GENERAL PURPOSE EXPER-			
532000 IMENTAL WORK AND TROUBLE SHOOTING			
532000			
532000 DC VOLTS	100 MV TO 10 KV		
532000 DC CURRENT	1 MA TO 10 AMPS		
532000 AC VOLTS	10 MV TO 300 VOLTS		
532000 AC CURRENT	1 MA TO 1 AMP		
532000 OHMMETER	10 OHMS TO 10 MOHMS		
532000			
ALS	#SFRIES X-1	\$4000	
ACCURACY	+/-0.00088FS		
FUNCTIONS	DC VOLTS, AC VOLTS, MILLIVOLTS, X OHMS AND RATIO		
COMMON MODE REJECTION	>136DB		
ROHDE AND SCHWARTZ	#UVV (110.4716-C2)	\$900	111
VOLTAGE RANGE	0.1 MV TO 1000 V (DC OR AC)		
FREQUENCY RANGE	0 TO 100 KHZ		
CURRENT RANGE	0.1 NA TO 1000 MA		
FEWLETT PACKARD	#4276/4288/456A	\$1150	112
CONFIGURATION	LOOSE EQUIP		
VOLTAGE RANGE	100 MV TO 1000 V (DC VOLTAGE)		
	10 MV TO 300 V (AC VOLTAGE)		
OHMMETER	10 OHM TO 10 M OHM		
CURRENT RANGE	1 MA TO 1 A (DC)		
	1 MA TO 1 A RMS (AC)		
SIMPSON	#249-3	\$98	111
CONFIGURATION	LOOSE EQUIPMENT		
VOLTAGE RANGE	AC-3, 8, 40, 160, 400, 800 VOLTS		
	DC-1, 6, 8, 40, 160, 400, 800, 1600, 4000 VOLTS		
OHMMETER	2K, 20K, 200K, 2MEG, 20MEG, AND 200MEG OHMS		
CURRENT RANGE	16, 160 MICROAMPERES		
	1.6, 16, 160 MILLIAMPERES		
	1.6, 16 AMPERES		



Space Division
Rockwell International

JOHN FLUKE MFG CO INC #B125A \$1945 311
VOLTAGE RANGE +/- 1 V TO +/- 1000 VDC
1 V TO 1000 VAC
RESOLUTION +/- 0.01%
FREQUENCY RANGE DC, 30 HZ TO 20 KHZ

JCHA FLUKE MFG CO INC #B425A-01 \$5995 311
VOLTAGE RANGE +/- 0.1 TO +/- 1000 VDC
1 TO 1000 VAC
RESOLUTION 0.001%
FREQUENCY RANGE DC, 10 HZ TO 100 KHZ

101000 SIGNAL GENERATOR LS052 11111
101000
101000 ACCURATELY REPRODUCE SELECTED FREQUENCIES AS REFERENCES FOR PUR-
101000 POSES OF MAINTENANCE AND REPAIR.
101000
101000 BANDWIDTH DC TO 20 KHZ
101000 OUTPUT VOLTAGE 3 V TO 0.1 MV
101000

INTERSTATE ELECTRONICS CORP #B32 \$345 111
FREQUENCY RANGE 0.03 HZ TO 3 MHZ
ACCURACY +/- 2%
AMPLITUDE 60 DB ATTENUATION, 20 V P-P
INTO OPEN CIRCUIT
OFFSET +/- 10 V (VARIABLE)

EXACT ELECTRONIC INC #122 \$395 111
BANDWIDTH 0.1 HZ TO 3 MHZ
FREQUENCY STABILITY 0.04%/10 MIN, 0.2%/24 HOURS
OUTPUT 10 DB ATTENUATION IN 10 DB STEP
10 V P-P INTO 50 OHMS
DC OFFSET +/- 10 VOLTS

PROHA-HITE #5300 \$720
WAVEFORMS SINE, SQUARE, TRIANGLE, PULSES,
PULSES, LIN-LOG SWEEP, TRIGGERED
AND RAMP
BANDWIDTH 0.002 HZ TO 3 MHZ
MOUNT RACK
FREQUENCY STABILITY 0.05% IN 10 MIN
SIZE 19X5-1/4X15-1/2
WEIGHT 16LB (BKG)

ROHDE AND SCHWARZ #SSN (204.P014.XX) \$5390 111
FREQUENCY RANGE 0.01 HZ - 120 KHZ
OUTPUT LEVEL 7.75 MV - 6 V

JOHN FLUKE MFG CO INC #B45A-156A \$13,950
FREQUENCY RANGE DC TO 50 MHZ
INCREMENTS 0.01 HZ
SPURIOUS OUTPUTS > 100 DB (NON-HARMONICS)
> 30 DB (HARMONICS)
SIGNAL-TO-PHASE NOISE RATIO > 66 DB
FREQUENCY STABILITY 2E-9/24 HOURS AVAILABLE

JCHA FLUKE MFG CO INC #B33A \$5965 111
CONFIGURATION RACK MOUNTED
FREQUENCY RANGE DC TO 11 MHZ
INCREMENTS 10 HZ
ACCURACY +/- 2% OF RANGE

INCRE-DATA CORP. #12P \$695 112
RANGE 0.1 TO 3 MHZ, 20 HZ TO 20 KHZ
STABILITY 0.05 % OF SETTING FOR 10 MIN,
0.25 % OF SETTING FOR 24 HRS.
ACCURACY +/- 2 % OF FS
PAMP GENERATOR 100 S TO 10 MICROSEC
WAVE FORMS SINE, SQUARE, TRIANGLE, PAMP,
PULSE, SYNC

WAVETER #130 \$395 111
FREQUENCY RANGE 0.2 HZ TO 2 MHZ
OUTPUT 40 DB ATTENUATION (VARIABLE)
10 V P-P INTO 50 OHMS
DC OFFSET +/- 5 VOLTS
RESOLUTION +/- 2% FS (1 HZ TO 2 PHZ)
STABILITY +/- 0.05%/10 MIN, +/- 0.25%/24
HOURS

399500 PLETHYSMOGRAPH, LIMB (INCLUDING COUPLER) LS1010 11111
399500
399500 MEASURE CHANGES IN BLOOD VOLUME AND VASCULAR RESPONSES
399500
399500 MEASUREMENT RANGE HIGH RANGE-1.5 TO 500 CMPS
399500 LOW RANGE-1.5 TO 50 CMPS
399500

HONEYWELL INC #ETP-4 111
TYPE IMPEDANCE PLETHYSMOGRAPH
001000 ACCELEROMETER - SPECIMEN MOVEMENT LS901 12111
001000
001000 ORGANISM CONTAINER ACTIVITY AND ACTIVITY LEVEL MONITOR.
001000
001000 ACCURACY +/- 0.0004 G
001000 FREQUENCY RESPONSE 0 TO 100 HZ

ORIGINAL PAGE IS
OF POOR QUALITY



001000 INVEYER-INERTIAL INSTRUMENTS #0A-116-16 1650 111
TYPE ULTR-SENSITIVE, LINEAR
RANGE +/- 1G
SENSITIVITY 1 VOLT/G
FREQUENCY RESPONSE DC TO 300 HZ
RESOLUTION 1E-6 G
THRESHOLD 1E-6 G

003210 AUDIOMETER LS902 11111
003210 FREQUENCY SOURCE ADJUST IN BOTH FREQUENCY AND AMPLITUDE FOR
003210 EARPHONES OR BONE VIBRATORS
003210 FREQUENCY.....125 TO 4000 ADJUSTIBLE IN OCTAVES
003210 ATTENUATOR.....CONTINUOUSLY ADJUSTIBLE IN 10P STEPS
003210

LAFAYETTE INSTRUMENT CO #3120-AC 5440 112
CONFIGURATION BENCH EQUIP
FREQUENCY OCTAVE FREQUENCIES 125 TO 4000
HZ, INTERMEDIATE FREQUENCIES
750 TO 4000 HZ
ATTENUATOR 1 OR 5 DB STEPS

442500 ANALYZER - ATOMIC ABSORPTION SPECTROPHOTOMETER LS903 11111
442500 QUANTITATIVE DETERMINATION OF METALLIC AND SEMIMETALLIC ELEMENTS
442500 IN SOLUTIONS AND INDICATING CONCENTRATION.
442500
442500 WAVELENGTH RANGE 200 TO 800 NM
442500 MODES OF OPERATION ATOMIC ABSORPTION, FLAME
442500 EMISSION, UV-VISIBLE SPECTRO-
442500 PHOTOMETRY
442500 COUNTER ACCURACY +/- 0.4 NANOMETERS
442500

BECKMAN INSTRUMENTS #24 52760
OPERATING RANGE 340 TO 700
WAVELENGTH RANGE 190 TO 1000
PRESENTATION DIGITAL COUNTER, 0.2
NM INCREMENTS
ACCURACY +/- 0.5 NM
REPEATABILITY +/- 0.2 NM
DIMENSIONS 24X16X14 INCHES
WEIGHT 70 LBS

BAUSCH AND LOMB #AC2-20 56280 111
CONFIGURATION MODULAR/BENCH
MODES OF OPERATION ATOMIC ABSORPTION, FLAME
EMISSION, UV-VISIBLE SPECTRO-
PHOTOMETRY
WAVELENGTH RANGE 190 TO 800 NM
WAVELENGTH READOUT DIGITAL COUNTER
WAVELENGTH COUNTER ACCURACY +/- 0.4 NM
WAVELENGTH SETTABLEITY +/- 0.001 NM

E.E.L. (AUTHOR H THOMAS CO) #140 53275 111
CONFIGURATION BENCH/LOOSE EQUIP
MODE OF OPERATION ATOMIC ABSORPTION
WAVELENGTH RANGE 200 TO 800 NM
ACCURACY +/- 1%

BECKMAN #579 55950
CONFIGURATION BENCH/LOOSE EQUIP
TYPE ATOMIC ABSORPTION WITH GRATING
MONOCHROMATOR
WAVELENGTH RANGE 190 TO 770 NM
ACCURACY +/- 0.5 NM
RESOLUTION 0.2 NM

042000 CELL COUNTER - BLOOD LS907 11111
042000
042000 MEASURE BLOOD CELL PROPERTIES OF BLOOD SAMPLES, ESSENTIALLY A
042000 COULTER COUNTER
042000
042000 PARAMETERS MEASUREDHEMOGLOBIN,HEMATOCRIT,RED BLOOD
042000 CELL COUNT,WHITE BLOOD CELL
042000 COUNT,MEAN CELL VOLUME,MEAN
042000 CELL HEMOGLOBIN AND MEAN CELL
042000 HEMOGLOBIN CONCENTRATION
042000
042000 REPRODUCIBILITY.....+/-1%

COULTER ELECTRONICS #FN 55000
PARTICLE SIZE0.003 TO 24 MICRONS
COUNTING RATE0.040 ML/SEC
DIMENSIONS.....49CM H X 35 CM W X 44 CM D
WEIGHT54LB (24.5KG)

COULTER ELECTRONICS INC #5 550000
PARAMETERSHEMOGLOBIN, HEMATOCRIT, RED
BLOOD CELL COUNT, WHITE BLOOD
MEAN CELL HEMOGLOBIN, AND MEAN
CELL COUNT, MEAN CELL VOLUME,
CELL HEMOGLOBIN CONCENTRATION
AUTOMATIC
TYPE
COULTER ELECTRONICS INC #2B1-6 510000 112
PARAMETERSHEMOGLOBIN, HEMATOCRIT, RED
BLOOD CELL COUNT, WHITE BLOOD
CELL COUNT, MEAN CELL VOLUME,
MEAN CELL HEMOGLOBIN, AND MEAN
CELL HEMOGLOBIN CONCENTRATION



TYPE MANUAL

096000 BACTERIAL COLONY COUNTER LS909 11111
096000
096000 MANUAL COUNT OF BACTERIAL COLONIES. UNIT IS STANDARD TYPE
096000

AMERICAN OPTICAL CORP #333C 1115 112
REGISTER..... MODEL AO 334A USED TO DATA
CUMULATIVE RECORDING OF COUNTS
TYPE QUEBEC COLONY COUNTER, DARK-
FIELD ILLUMINATION APH
DIMENSIONS..... 10X11X10 INCHES

104160 DIGITAL PLOTTER, PRINTER LS909 11111
104160
101340 KEYBOARD LS910 11111
101340
101340 TYPEWRITER TYPE KEYBOARD USED FOR PROVIDING INSTRUCTIONS FOR THE
101340 COMPUTER
101340

TEKTRONIX #46C1 \$3750 111
DISPLAY ALPHANUMERIC AND GRAPHIC
COPY SIZE VARIABLE BETWEEN 4.5X6 IN AND
8.5X14 IN
COPY TIME 18 SEC (FIRST COPY)
RESOLUTION ACTUAL SIZE OF 4000 CHARACTER
DISPLAY, BASED ON 90X70 PXL
MATRIX

VARIAN DATA MACHINES #ASR 33 & ASR 35 111
SPEED 10 CHARACTERS/SEC.
PAPER CONTINUOUS ROLL 8.5 IN WIDE
TYPING LINE 85 CHARACTERS

620060 ELECTROPHYSIOLOGY BACKPACK LS912 11235
620060
620060 BACKPACK CONTAINS THE NECESSARY ELECTRONICS FOR SENSING AND
620060 TRANSMITTING MAN'S PHYSIOLOGICAL DATA SUCH AS ECG, EFG, EMG,
620060 ECG, ETC. INCLUDES SENSORS, SIGNAL CONDITIONERS, MULTIPLEXERS,
620060 A/D CONVERTERS AND TRANSMITTERS
620060

HEWLETT PACKARD #78100A \$600 112
618090 ELECTROPHYSIOLOGY RECEIVER LS913 11111
618090
618090 PROVIDE SPECIAL PERCEPTION OF CARDIOVASCULAR AND NEURAL ELECTRO-
618090 PHYSIOLOGICAL EVENTS VIA BIOTELEMETRY SYSTEMS
618090
618090 SIGNAL TYPES RECEIVED.....ELECTROENCEPHALOGRAPH
618090ELECTROCARDIOGRAPH
618090VECTOCARDIOGRAPHY
618090PALISTOCARDIOGRAPHY
618090IMPEDANCE CARDIOGRAPHY
618090PHONOCARDIOGRAPHY
618090

043000 BLOOD CLOT FIRMOMETER LS915 11111
043000
043000 PROVIDE AUTOMATIC MEASUREMENT OF PLASMA COAGULATION TIME
043000
043000 ACCURACY..... +/- 0.1 SECOND
043000

B-E-L #B4180-1,-2,-3,-5,-6 \$974 111
CONFIGURATION MODULAR EQUIP
COAGULATION TIME ACCURACY +/- 0.1 SEC
THERMAL PREP-BLOCK HOLDS UP TO 20 FIBROTUBE DIS-
POSABLE REACTION CUPS, AND UP
TO TEN 12 MM TEST TUBES. PRE-
HEATS AND MAINTAINS PLASMA AND
REAGENTS AT 37 C
AUTOMATIC PIPET DELIVERS MEASURED 0.1 OR 0.2 ML
VOLUMES OF PLASMA OR REAGENT

108600 FREEZER, CRYOGENIC LS916 11111
108600
108600 PROVIDE A MEANS FOR FREEZING AND STORING BIOLOGICAL SPECIMENS
108600
108600 STORAGE CAPACITY.....115 CU. IN. (1.89 LITERS)
108600 OPERATING TEMPERATURE.....-320F (78K)
108600 MIN CRYOGEN HOLDING TIME.....3 WEEKS BETWEEN REFILLS
108600

CRYOGENIC ASSOCIATES #IR-60 \$1000
CAPACITY.....0.4 L LIT. N2
HOLDING TIME.....24 HOURS
DIMENSIONS.....4.25 IN DIAM, 8.5 IN LONG
WEIGHT.....1.5 LB (EMPTY), 2.5 LB (FULL)

UNION CARBIDE CORP/LINDE DIV #LR-10A-6 \$2500 112
CONFIGURATION LOOSE EQUIP
LIQUID NITROGEN CAPACITY 10.4 LITERS
STORAGE TEMPERATURE -320 F
CAPACITY 6 CANISTERS, 115 CU. IN.
LIQUID NITROGEN CONSUMPTION 0.34 L/DAY (w/0 CANISTERS)
0.26 L/DAY (w/CANISTERS)

207500 HEMATOCRIT, ELECTRONIC LS917 11111
207500
207500 AUTOMATICALLY DETERMINE THE PERCENT HEMATOCRIT IN BLOOD.
RAPID & READOUT AND TEMPERATURE COMPENSATOR REQUIRED

ORIGINAL PAGE IS
OF POOR QUALITY



207900
VHILW SPRING INSTRUMENT *3C 9275 112
CONFIGURATION LOOSE EQUIP
ANALYSIS TIME < 30 SEC
PIPET 0.02 ML WHOLE OR HEPARINIZED
BLOOD

308500 DOSIMETER, RADIATION DETECTOR LS91P 11111
308500
308500 PORTABLE DEVICE WHICH WILL ALERT USER TO RADIATION LEVELS IN
308500 EXCESS OF A PREDETERMINED LEVEL AND TO PROVIDE A DIRECT READOUT
308500 OF THE TOTAL CUMULATIVE RADIATION DOSE.
308500
308500 RATE MEASUREMENT RANGE 0 TO 0.1 MR/HR
308500 0 TO 50 R/HR
308500 ADJUSTABLE RADIATION ALARM ... 0.05 MR/HR TO 50 R/HR
308500 CUMULATIVE DOSE RANGES 0 TO 200 MR
308500 0 TO 500 R
308500 ACCURACY (X- & GAMMA RAY) +/- 10 %

308500
KAHL SCIENTIFIC INSTRUMENT *BFH 408 9130 112
CONFIGURATION PORTABLE, FULLY TRANSISTORIZED
RATE MEASUREMENT RANGES 0 TO 50 MR/HR
ACCURACY +/- 15% (FOR GAMMA)
TEMPERATURE RANGE 15 TO 130 F

305000 RADIATION DETECTOR LS91P 11111
305000
305000 ACCURATE MEASUREMENT OF EITHER ALPHA, BETA, GAMMA, X-RAY OR
305000 NEUTRON RADIATION AS IT IS APPLIED TO A BIOLOGICAL SPECIMEN.
305000
305000 RATE RANGES 0 TO 0.1 MR/HR
305000 0 TO 1000 R/HR
305000 CALIBRATED READOUTS COUNTS/SEC, /MIN, AND TOTAL

305000
HEWLETT PACKARD *5890A/5583A/ 91150/9650 112
5290A/5554 9175/9360
/5582A 9750
CONFIGURATION MULTI-UNIT BENCH EQUIP
RADIATION DETECTED ALPHA, BETA, X-RAY, GAMMA OR
NEUTRON
RATE MEASUREMENT RANGE 0 TO 1000 R/HR

315500 ISOTOPE SOURCE - SELF CONTAINED LS920 11111
315500
315500 RADIATION STRESS SOURCE USED IN CONDUCTING RED BLOOD CELL
315500 SURVIVAL STUDIES. TYPICAL ASSORTMENT - NITROGEN 15 COMPOUNDS,
315500 RADIOLABELED SURCELLULAR PRODUCTS (RIBOSOMES), RADIOLABELED POLY-
315500 NUCLEOTIDES OR TITRATED THYMIDINE.
315500
315500 RADIATION LEVEL 1 TO 100 MICROCURIE
315500

315500
MILES LABORATORIES ASSORTMENT 91-5K EA 111
TYPE CARBON-13 LABELED COMPOUNDS,
NITROGEN-15 LABELED COMPOUNDS,
TRITIUM AND CARBON-14 LABELED
RADIOACTIVE ISOTOPIIC COMPOUNDS,
HOMOPOLYMERS, RADIOLABELED POLY-
NUCLEOTIDE DUPLEXES

618090 BIOTELEMETRY RECEIVER - COMPACT CAGE MODULE LS921 121113
618090
618090 MICROBACKPACKS FOR SMALL VERTEBRATE ACQUIRING AND TRANSMITTING
618090 BIO DATA.
618090

HEWLETT PACKARD *78101A 9800 112
CONFIGURATION MICROBACKPACK
ANTENNA LOW PROFILE MONOPOLE

603003 SIGNAL CONDITIONER LS922 11111
603003
603003 ELECTRICAL SIGNAL TRANSFORMATION BETWEEN ANY NUMBER OF TRANS-
603003 DUCERS AND MULTIPLEXER/ANALOG TO DIGITAL CONVERTER
603003
603003 COMPUTER INPUT +/- 10 VOLTS
603003

HEWLETT PACKARD *8807A WITH OPTION 9775
OC1
CONFIGURATION PLUG IN UNIT FOR OSCILLOGRAPHIC
RECORDER
FREQUENCY RANGE 50 HZ TO 100 KHZ
INPUT RANGES 0.02, 0.05, 0.1, 0.2, 0.5, 1,
2, 5, 10 V RMS/DIV
ACCURACY +/- 2%

FENNELFC 9TC133 9295 111
SET PREAMPLIFIER
CHANGE SENSITIVITY 10E12 VOLTS/COULCME
INPUT NEGATIVE INPUT PULSE POLARITY
WITH INPUT CAPACITANCE >2000PF
OUTPUT POSITIVE OUTPUT PULSE POLARITY
WITH OUTPUT IMPEDANCE 50 OHM
STABILITY 0.5%/V

HEWLETT PACKARD *9461A 9380 111
FREQUENCY RANGE 1 KHZ TO 150 KHZ
FREQUENCY RESPONSE +/- 1 DB INTO 50 OHM LOAD
GAIN AT 900 KHZ 40 DB +/- 0.5 DB OR 20 +/- 1 DB
MAX INPUT 1 V RMS OR 2 V P-P PULSE



TENNELEC	BYC 202703	1395/495	
CATN.....	2.5 TO 300GR		
NON LINEARITY.....	LESS THAN 0.05R		
CONFIGURATION.....	NIM COMPATIBLE		
ROHDE AND SCHWARZ	#ATN (100.0099.02)		111
FREQUENCY RANGE.....	30 HZ TO 20 MHZ		
INPUT VOLTAGE.....	0.25 - 25 VOLTS		
MAX POWER.....	50 WATTS		
RESEARCH INC	#R12-11-17	\$3200	111
CAPACITY.....	TO 33 PAIRS OF SIG WIRES (16GA)		
REFERENCE TEMP.....	HEAT SINK CONTROLLED TO +/- .0R		
	DEG C AT 65 DEG C		
AMBIENT TEMP.....	0 TO 49 DEG C		
LEAKAGE RESISTANCE.....	1000 MEG OHM		
HEAT-UP TIME FROM 22 DEG C.....	50 MINUTES		
SETTLING TIME.....	4 MS TO FINAL VALUE		
HONEYWELL	#ACCUDATA 117	\$1360	111
MULTICHANNEL WIDEBAND AMPLIFIER			
CHANNELS.....	7		
MOUNT.....	BENCH		
HONEYWELL	#ACCUDATA 115	\$5370	111
HIGH VOLTAGE DC AMPLIFIER			
INPUT SIGNAL VOLTAGE.....	+/-50MV TO 1500V		
OUTPUT SIGNAL VOLTAGE.....	UP TO 2000V		
CHANNELS.....	0		
NEFF	#122-222	\$610	111
GAIN.....	1 TO 2500		
COMMON MODE REJECTION.....	140DB		
DRIFT.....	0.5MICROVOLT/DC		
STEP GAIN ACCURACY.....	+/-0.01%		
SWITCH SELECTABLE FILTER			
BANDWIDTH.....	160KHZ		
RACK/MODULE MOUNT			
NEFF	#126-621	\$320	111
AMPLIFIER MULTIPLEXER			
EMPLOYES FET SWITCH			
GAIN RANGE.....	0.2 TO 2500		
COMMON MODE REJECTION.....	120DB		
OUTPUT.....	+/-10V AT 10 MILLIAMPERES		
NEFF	#127-102	\$1150	
PROGRAMMABLE GAIN WITH SWITCH SELECTABLE OUTPUT FILTER			
GAIN STEPS.....	1, 2, 5, 10, 20, 50, 100, 200,		
	500,1000		
COMMON MODE REJECTION.....	140DB		
FET SWITCHING			
RACK/MODULE CONFIGURATION			
TENNELEC	#TC133	\$295	111
FET PREAMPLIFIER			
CHARGE SENSITIVITY.....	10E12 VOLTS/COULOMB		
INPUT.....	NEGATIVE INPUT PULSE POLARITY		
	WITH INPUT CAPACITANCE 22000PF		
OUTPUT.....	POSITIVE OUTPUT PULSE POLARITY		
	WITH OUTPUT IMPEDANCE 50 OHM		
STABILITY.....	0.5%/V		
STATFAM/V MUELLER AND CO	#SM102G/SM1011/	\$3000	111
	SM100B/SM100A		
CONFIGURATION.....	RACK MOUNTED		
TYPES.....	VENUS PRESSURE MODULE - DISPLAY		
	METER 1.5 TO 30 MP HGT		
	BRIDGE AMPLIFIER MODULE - SIX		
	SENSITIVITY RANGES 10-1 MV FS,		
	0-2 MV FS, 0-5 MV FS, 0-10 MV		
	FS, 0-20 MV FS, 0-50 MV FS WITH		
	7.5 V EXCITATION		
	ELECTRONIC THERMOMETER MODULE -		
	METER SCALES 93-107 F, 34-61C		
HONEYWELL	#ACCUDATA 106	\$400	111
MILLIVOLT/THERMOCOUPLE CONTROL UNIT			
CONTAINS COLD JUNCTION COMPENSATION, SUPPRESSION, CALIBRATION,			
AND ATTENUATION			
MOUNTS.....	MODULE IN 7X19IN. RACK		
HONEYWELL	#ACCUDATA 120	\$430	111
CL AMPLIFIER			
FREQUENCY RESPONSE.....	DC TO 100KHZ		
GAIN.....	X10 TO X1000 IN EIGHT STEPS		
COMMON MODE REJECTION.....	GREATER THAN 80DB		
INPUT IMPEDANCE.....	OVER 10 MEGOHMS		
MOUNT.....	MODULE IN 7X19IN. RACK		
HONEYWELL	#ACCUDATA 122/123	\$495/600	111
DIFFERENTIAL WIDEBAND AMPLIFIER			
INPUT RANGE.....	+/-2MV TO 100VFS.		
GAIN SETTINGS.....	16		
FILTERING.....	SELECTABLE		
BANDWIDTH.....	DC TO 100KHZ		

ORIGINAL PAGE IS
OF POOR QUALITY



MINIEMERL *ACCUDATA 214-2 1111 111
CASE CONTROL/AMPLIFIER
OUTPUT FILTERS..... 5
INPUT SUPPRESSION..... +/-100MV
GAIN..... UP TO X10,000

003210 SOUND LEVEL METER LS923 11111
003210
003210 SOUND LEVEL SURVEYS AND LIMITED FREQUENCY ANALYSIS IN REAL TIME,
003210
003210 FREQUENCY RESPONSE 125 TO 8000 HZ
003210 SOUND LEVEL RANGE 40 TO 100 DB
003210
BROEL AND KJAER INSTRUMENTS #3501/S \$2500 111
CONFIGURATION CONSISTS OF: 2203 SOUND LEVEL
METER, 1613 OCTAVE BAND FILTER,
4220 PISTON PHONE CALIBRATOR,
MICROPHONES, ACCELEROMETERS AND
CABLES

904500 STAINING SYSTEM - BACTERIOLOGICAL LS924 11111
904500
904500 STAIN AND FIX BIOLOGICAL SPECIMENS FOR MICROSCOPIC EXAMINATION.
904500
VWR SCIENTIFIC DIV #48420-050 \$1715 111
CONFIGURATION BENCH EQUIPMENT

BECKMAN INSTRUMENTS PROTOTYPE 125
STATUS NOW ON SKYLAR
CAPACITY..... 8 BLOOD SMEARS, 16 BACTERIA
SLIDES WITH 1 FILLING
ZERO GRAVITY USE..... YES
WEIGHT..... 1783 GRAMS
DIMENSIONS..... 4.1X5.1X9 INCHES
REAGENT STORAGE LIFE..... 2 YEARS
PREPARATION TIME..... ABOUT 15 MIN

171600 TASKHARD - FORCE/TORQUE LS925 11111
171600
171600 MEASURE MAN'S CAPABILITY TO APPLY FORCES AND TORQUES IN A VARIETY
171600 OF DIRECTIONS FROM VARIOUS HAND/BODY ORIENTATIONS AND RESTRAINT
171600 CONDITIONS.
171600
171600 HAND FORCE APPLICATIONS +/- 0.5 POUNDS
171600
171600 GENERAL ELECTRIC #NAS 9-RA40 512
CAPABILITIES MEASURE ONE/TWO HANDED FORCE
APPLICATIONS WITHIN +/- 0.5 LB.
RESTRAIN TEST SUBJECT IN THE
DESIRED POSITION.

301000 RADIATION COUNTER - BIOCHEMICAL SAMPLE LS926 11111
301000
301000 FAST, ACCURATE AND EASY BIOCHEMICAL RADIATION COUNTING SYSTEM
301000 BASED ON STANDARD SIZED BIOCHEMICAL SAMPLES.
301000
301000 STANDARD PLANCHET UP TO 2 IN DIA
301000 CAPACITY OF PLANCHETS 120
301000
HEWLETT PACKARD #5560A 111
CONFIGURATION BENCH/LOOSE EQUIP
TYPE NUCLEAR COUNTING SYSTEM, ULTRA-
LOW BACKGROUND, GAS-FLOW CR
SCINTILLATION DETECTORS

459605 MASS SPECTROMETER LS928 11111
459605
459605 PROVIDE AN ANALYSIS OF AN UNKNOWN IONIZED GAS BY ATOMIC MASS
459605 ALMBER
459605
459605 MASS RANGE..... C TO 60 AMU
459605

004100 AIR PARTICLE SAMPLE COLLECTOR LS929 11111
004100
004100 AIR PARTICLE AND MICRO-ORGANISM SAMPLES FOR AIR QUALITY DETERMIN-
004100 ATION
004100
004100 AIRFLOW 1 CU FT/MIN
004100 SAMPLING TIME - CONTINUOUS ... 1 HOUR
004100
PAUSCH - LOMB INC. #40-1A
SAMPLING RATE 0.1 CU. FT./100 SEC.
(0.17 LITERS/MIN.)
MAX. CONCENTRATION MEAS ... 1E+7 PARTICLES/CU. FT.
(135E+3 PARTICLES/LITER)
PARTICLE SIZE RANGE 0.3, 0.5, 1.0, 2.0, 3.0, 5.0, 10 MICRONS

042000 AUTOANALYZER - MULTIPLE LS930 11111
042000
042000 AUTOMATIC ANALYSIS OF BLOOD AND APPROXIMATELY 30 OTHER FLUIDS.
042000
460300 SPECTROPHOTOMETER - GENERAL ANALYZER LS931 11111
460300
460300 SPECTRAL ANALYSIS OF GASES AND LIQUIDS INCLUDING SOLIDS OR LIGHT
460300 SOURCES.
460300



440000 FREQUENCY 250 TO 2500 NANOMETERS
440000 (2500 TO 25000 ANGSTROMS)
440000 ACCURACY +/- 0.4 NANOMETERS

440000 VARIAN INSTRUMENTS #CARY 14 \$22000 112
CONFIGURATION BENCH EQUIP
WAVELENGTH RANGE 1840 A TO 2.65 MICRONS
ACCURACY +/- 4 A
REPRODUCIBILITY < 0.5 A IN UV-VIS
..... < 2.5 A IN NEAR IR

CARY #SIMULS
FREQ RANGE 200 TO 700 NM
RESOLUTION +/- 0.5 nm
PHOTOMETRIC ACCURACY +/- 0.3%
BANDWIDTHS 0.5 OR 5.0 NM

SHIMADZU #BAQV-5CD \$3500
WAVELENGTH RANGE 183 TO 1200NM
ACCURACY 0.1%
RESOLUTION 0.05NM AT 205NANOMETERS
PHOTOMETRIC ACCURACY +/-0.2%
BANDWIDTHS 10C TO 30C, 300 TO 400, 400 TO
600, 600 TO 850 AND 850 TO 1200
NANOMETERS

CORNING #24C \$5-50 111
FREQUENCY RANGE 190 TO 800NM
RESOLUTION +/- 0.045UNITS FULL SCALE

PERKIN ELMER INSTRUMENT DIV #323 \$15200
WAVELENGTH RANGE 185 TO 2500 NANOMETERS
MODES X-Y RECORDER
AUTOMATIC REPETITIVE SCAN

G K TURNER ASSOCIATES #350 \$1609 111
WAVELENGTH RANGE 210 NM TO 1000 NM
ACCESSORIES 330-030 IR ACCESSORY
..... 330-040 UV ACCESSORY
OPTICS REFLECTION TYPE
WAVELENGTH ACCURACY +/- 2 NM
WAVELENGTH REPRODUT DIGITAL IN NM
PHOTOMETRIC REPEAT ACCURACY +/- 0.5%

BECKMAN INSTRUMENTS #ACTA P-VII 133700 \$23000 112
SPECTRAL RANGE 19C -800NM (PHOTOMULT TUBE)
..... 800 -3000NM(LEAD SULFIDE CELL)
RESOLUTION BETTER THAN .05NM(19C-800NM)
..... BETTER THAN .3NM (800-3000NM)
PHOTOMETRIC ACCURACY003AT 1.0A .03 AT 3.0A
CHART SPEEDS02 TO 4NM/SEC (190-800NM)
..... .06 TO 1NM/SEC(800-3000NM)
SCAN EXPANSION 1 TO 100NM/INCH(190-800NM)
..... 4 TO 400NM/INCH(800-3000NM)
DIMENSIONS 24X62X22 INCHES(161X157X56CM)
WEIGHT 250 LBS (113KG)

000590 ANESTHETIZER - INVERTEBRATES LS932 11114
000590 BREWED INVERTEBRATE ORGANISMS INSENSIBLE TO FACILITATE HANDLING.
000590 ANESTHETIC GAS CARBON DIOXIDE

146915 ANTHROPOMETRIC GRID LS933 11111
146915 ANTHROPOMETRIC MEASUREMENTS OF VERTEBRATES EITHER REMOTELY OR
146915 AUTOMATICALLY.
146915 SMALL GRID 5 X 15 CM SPACED AT 1 MM
146915 MEDIUM GRID 50 X 75 CM SPACED AT 2 MM
146915 LARGE GRID 2 X 2 " SPACED AT 5 MM

107000 BENCH, LAMINAR AIRFLOW LS934 11114
107000 A GLOVE BOX WITH RELATIVELY HIGH AIR FLOW FOR CONTROL OF
107000 PARTICULATE AND GASEOUS CONTAMINANTS WITHIN
107000 BENCH, GENERAL EXPERIMENTS LS935 11134
107000 PROVIDE WORK AREA FOR PREPARATION OF EXPERIMENTS, EXPERIMENT
107000 MAINTENANCE OBSERVATIONS, AND PREPARATIONS FOR RETURN TO EARTH
107000 ELECTRICAL UTILITY.....28VDC, 400HZ AC, 60HZ AC
107000 VACUUM UTILITY.....10E-6 TORR
107000 PRESSURE SOURCE.....50 PSIG (3.45E4 N/M2)
107000 AVAILABLE GASSES.....OXYGEN, NITROGEN, CARBON DIOXIDE

055055 CAMERA CINE LS936 11111
055055 PROVIDE VISUAL RECORDS SUCH AS PHOTOMICROGRAPHY, TIME LAPSE
RECORDS, PHOTOMICROSCOPY AND GENERAL EXPERIMENT DOCUMENTATION
055055 OPERATION MODE.....VARIABLE FRAME RATE OR SINGLE
055055 FRAME PULSE OPERATION UNDER
055055 REMOTE CONTROL
055055 FILM SIZE.....16 OR 35MM
055055 LENS SELECTION....."C" MOUNT LENS SYSTEM; STANDARD
055055 ZOOM LENS AS A MINIMUM
055055 EXPOSURE SYSTEM.....AUTOMATIC OR ZOOM
055055 FILM CAPACITY.....90, 100 OR 200FT (15.2,30.4 OR
60.0M) ROLLS

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

```

099055
J A PAURER INC #3C8 115 - 20,000 911
FILM ..... 16 MM DOUBLE OR SINGLE PERFOR-
ATED
FORMAT ..... 0.295 X 0.404 INCH
FRAME RATE ..... 1, 6, 12 AND 24 FR/SEC. TIME
SHUTTER SPEED ..... TIME, 1/60, 1/125, 1/250, 1/500
AND 1/1000 SEC
FILM CAPACITY ..... 90 FT STD-BASE FILM, 145-165 FT
THIN-BASE FILM
FILM ON APOLLO

J A PAURER INC #ODYSSEY 16 911
FRAME SPEED ..... 8-96 FPS
LENS MOUNT ..... WIDEBASE BAYONET TYPE
MAGAZINE ..... 400 FT
ASA ..... 8 TO 500

BEATTIE-COLEMAN #KD-9 111

603014 CAMERA CONTROLLER LS937 11114
603014
603014 DEVICE TO CONTROL THE OPERATION OF VIDEO CAMERAS THROUGHOUT THE
603014 LABORATORIES
603014
099050 CAMERA, PLATE FILM LS938 11111
099050
099050 PROVIDE STILL PHOTOGRAPHIC COVERAGE
099050
099050 FILM TYPES.....HIGH RESOLUTION GLASS PLATE
099050 .....STANDARD 2 1/4 X 2 1/4 IN (5.75
099050 .....X 5.75 CM) ROLL
099050 .....TOMM
099050 .....POLAROID
099050 SHUTTER SPEED.....1 SEC TO 1/500 SEC
099050

HASSELBLAD #900 EL 91411 111
TYPE ..... SINGLE LENS REFLEX
DRIVE ..... ELECTRICAL OR MANUAL
LENS ..... 80 MM PLANAR F/2.8 ZEISS
50 MM DISTAGON F/4.0
CLOSE-UP LENS ..... PROXAR 0.5 FOCUS RANGE
17.25 TO 24.25 IN
PROXAR 1.0 FOCUS RANGE
22.25 TO 42.5 IN
FILM SIZE ..... 2.25 X 2.25 IN
MAGAZINE ..... CASSETTE AND PLATE

TENTONIX INC #C-12 9590 111
LENS ..... 75 MM
STOP ..... F/1.9 TO F/16
MAGNIFICATION ..... 0.75
LENS SPEED ..... 1 TO 1/100 SEC MECH
4 TO 1/60 SEC ELECT
FILM TYPE..... POLAROID

TENTONIX #C30 9525
MAXIMUM RELATIVE APERTURE..... F/1.9
MAGNIFICATION..... 0.7 TO 1.5
RELATIVE SPEED..... 1.0
FIELD OF VIEW..... 3.15X3.93IN (8X10CM)
FILM TYPE..... POLAROID

HEWLETT-PACKARD #158 9420
MAXIMUM RELATIVE APERTURE..... F/3.5
MAGNIFICATION..... 110.85
SPEED..... 8 TO 1/40 SEC
LENS..... 75MM
FILM TYPE..... POLAROID

HEWLETT-PACKARD #159A 91025
MAXIMUM RELATIVE APERTURE..... F/1.3
MAGNIFICATION..... 110.5
SPEED..... 8 TO 1/30 SEC
LENS..... 80MM
FILM TYPE..... POLAROID

480000 VIDEO CAMERA, B/W LS939 401131
480000
480000 PROVIDE MEANS OF ACTIVITY MONITORING, EXPERIMENT DATA ACQUISITION
480000 ETC. SYSTEM WILL INTERFACE WITH A 40 INPUT VIDEO MULTIPLEXER TO
480000 PERMIT THE MONITORING OF 40 SIMILAR VIDEO CAMERAS.
480000
480000 VISUAL RESPONSE.....APPROXIMATE HUMAN EYE
480000 VIDEO OUTPUT.....CONSTANT WITH LIGHT LEVEL
480000 CHANGES OF 10 TO 1000FT CANDLES
480000 1.4V P-P COMPOSITE, CONFORMS TO
480000 EIA RS-170 STANDARD
480000

GENERAL ELECTRODYNAMICS CORP #EC 6038 A-1 910000 411
CONFIGURATION ..... LOOSE EQUIP
TYPE ..... ELECTROSTATIC FOCUS VIDICON
TUBE, PLUG-IN PRINTED CIRCUITS,
SOLID STATE SYSTEM

```



RESOLUTION	760 LINES		
CHANNELS	1		
VIDICON	1 INCH 101339		
LENS	C MOUNT		
VIDEO OUTPUT	1.4 V P-P COMPOSITE SYNC 75 OHM		
AUTOMATIC CONSTANT VIDEO OUTPUT	10 TO 10,000 FT-CANDLES		
HORIZONTAL SCANNING FREQUENCY	15,750 HZ		
VERTICAL SCANNING FREQUENCY	60 HZ		
INTERLACE	2:1 30 FR/SEC		
ASPECT RATIO	3:4		
SWEEP LINEARITY	> 18		
RCA ELECTRONICS COMPONENTS #8521	\$520	111	
TYPE	VIDICON-SULFIDE		
PHOTOCONDUCTOR	11		
IMAGE DIAGONAL	25 MM (1 IN)		
FOCUS	MANUAL		
COMU INC. ELECTRONICS DIV. #1220	\$6430	111	
IMAGE CONVERTER	4846P VIDICON - COLOR		
RESOLUTION-HORIZONTAL LIMIT	300 LINES MIN.		
GEOMETRIC DISTORTION	< 2 % OF PTC HEIGHT		
LENS MOUNT	16 MM - C-MOUNT		
COMU INC. ELECTRONICS DIV. #2000	\$1000	311	
LENS-BUILT IN	4:1 (20-40 MM) F/2.5 ZOOM		
LENS ATTACHMENT	4:1 (12.5-40 MM)		
LENS INTERCHANGEABLE	10:1 (15-150 MM) F/2.8 ZOOM		
10 MHZ BANDWIDTH (MOD 2006)	525 OR 725 LINE (VIDICON 7263A)		
20 MHZ BANDWIDTH (MOD 2004)	A73 OR 945 LINE (VIDICON 4573)		
MIL SPECIFICATIONS	MIL-E-5272C, MIL-STD-810		
COMU INC. ELECTRONICS DIVISION #4500	\$2225	311	
LENS	4:1 ZOOM		
.....	10:1 ZOOM		
VERTICAL SWEEP RATE	60 FIELDS PER SEC		
HORIZONTAL SWEEP RATE	525 LINES PER FRAME		
IMAGE TUBE TYPE	941A STD		
LENS MOUNT	16 MM C-MOUNT		
MIL SPECIFICATIONS	MIL-E-5400M, MIL-E-5300M		
481000 VIDEO CAMERA, COLOR	L5940	11111	
481000			
481000 PROVIDE COLOR VIDEO OBSERVATIONS OF EXPERIMENTAL ANIMALS, SETUPS,			
481000 ETC. FOR STORAGE VIA VIDEO-TAPE FOR TRANSMISSION TO THE GROUND.			
481000			
481000 RESOLUTION (MIN)	350 LINES		
481000 LENS	STANDARD 16MM, C MOUNT		
481000			
COMU ELECTRONICS INC #1000/1900/9800	\$20000	111	
CONFIGURATION	BOTH LOOSE AND RACK-MOUNT EQUIP		
COLOR SYSTEM	3 VIDICON SYSTEM		
LENS	C MOUNT VIDICON-FERRAR		
RESOLUTION (GREEN CHANNEL)	LIMITING AT 500 LINES, 508 AT 400 LINES, 1008 AT 300 LINES		
SCAN LINEARITY	< 1%		
VIDEO OUTPUT	0.7 V P-P WITH 300 FOOTLAMPHERTS		
COMU INC. ELECTRONICS DIV. #1220	\$6430	111	
IMAGE CONVERTER	4846P VIDICON - COLOR		
RESOLUTION-HORIZONTAL LIMIT	300 LINES MIN.		
GEOMETRIC DISTORTION	< 2 % OF PTC HEIGHT		
LENS MOUNT	16 MM - C-MOUNT		
COMU INC. ELECTRONICS DIV. #2000	\$1000	311	
LENS-BUILT IN	4:1 (20-40 MM) F/2.5 ZOOM		
LENS ATTACHMENT	4:1 (12.5-40 MM)		
LENS INTERCHANGEABLE	10:1 (15-150 MM) F/2.8 ZOOM		
10 MHZ BANDWIDTH (MOD 2006)	525 OR 725 LINE (VIDICON 7263A)		
20 MHZ BANDWIDTH (MOD 2004)	A73 OR 945 LINE (VIDICON 4573)		
MIL SPECIFICATIONS	MIL-E-5272C, MIL-STD-810		
COMU INC. ELECTRONICS DIVISION #4500	\$2225	311	
LENS	4:1 ZOOM		
.....	10:1 ZOOM		
VERTICAL SWEEP RATE	60 FIELDS PER SEC		
HORIZONTAL SWEEP RATE	525 LINES PER FRAME		
IMAGE TUBE TYPE	941A STD		
LENS MOUNT	16 MM C-MOUNT		
MIL SPECIFICATIONS	MIL-E-5400M, MIL-E-5300M		
068500 CENTRIFUGE, REFRIGERATED, HIGH SPEED	L5941	11111	
068500			
068500 SEPARATION EQUIPMENT TO SUPPORT MEDICAL AND BIOLOGICAL RESEARCH.			
068500			
068500 OPERATING TEMPERATURE	15 TO 50F (12 TO 10C)		
068500 SPEED	0 TO 20000 RPM		
068500			
BECKMAN INSTRUMENTS #J-21B-JA-20	\$3490	111	
SPINNING RATE	0-20000 RPM		
APPLIED GRAVITY	0-4000 G		
ACCOMMODATION	4-50 ML TUBES		
ROTOR	4JA-20		
TEMPERATURE	-15C AT REDUCED SPEED		
.....	2 C TO 30 C FULL SPEED		

ORIGINAL PAGE IS
OF POOR QUALITY



INTERNATIONAL EQUIPMENT CO. #RR-2 \$3100 112
SPINNING RATE 0 TO 9100 RPM
SEPARATION FORCE TO 25000 G
TEMPERATURE -15 TO 30 C
AUTOMATIC TIMER 2 MIN TO 2 HRS (ADJUSTABLE)

INTERNATIONAL EQUIPMENT CO. #RR-1 \$3100 112
SPINNING RATE 0 TO 20000 RPM
SEPARATION FORCE TO 40000 G
TEMPERATURE -20 TO 40 C
ACCOMMODATION 8 PLACES, 50 ML
AUTOMATIC TIMER 2 TO 120 MINUTES

INTERNATIONAL EQUIPMENT CO. #RR-2C \$3100 112
SPINNING RATE 20000 RPM
SEPARATION FORCE 40000 G
TEMPERATURE -30 TO 30 C
LOAD ACCOMMODATIONS 12 X 14 ML

BECKMAN INSTRUMENTS INC. #LS-5C \$9450
SPINNING RATE 0-50000 RPM
APPLIED GRAVITY 0-54000 G
ACCOMMODATIONS 1.5 ML CAPACITY
TEMPERATURE 0 C TO AMBIENT

MSE #LR-A \$3450 111
A/C HIGH SPEED ATTACHMENT
SPINNING RATE 22000RPM
APPLIED GRAVITY 37500GM
TEMPERATURE -4 TO 104F (-20 TO 40C)

060000 CENTRIFUGE, MICRO LS942 11111
060000
060000 CENTRIFUGAL SEPARATION OF BIOLOGICAL SAMPLES
060000
060000 CAPACITY 4 SAMPLES
060000 SPEED(MIN) 1750RPM
060000 GLASS TUBE SIZES 0.5 TO 5 ML
060000

INTERNATIONAL EQUIPMENT CO. #20211/20215/20216 \$86 112
CONFIGURATION BENCH/LOOSE EQUIP
APPLICATION MICRO OR SEMI-MICRO ANALYSIS
SPINNING SPEED 1760 RPM
SPIN STABILIZATION TIME 10 SEC
ACCOMMODATIONS 0.5, 1, 2, 3, 5 ML TUBES

ACAMS #DYNAC \$267
SPINNING RATE 700 TO 2500RPM
APPLIED GRAVITY 65 TO 1080 GCF
SAMPLE CAPACITY 24

MSE #GT-4 \$410 111
SPINNING SPEED 4000RPM
APPLIED GRAVITY 2900
SAMPLE CAPACITY 4

BAUSCH AND LOMB INC. #SYSTEM CENTRIFUGE 112
CAPACITY 100 SAMPLES
TIME TO SPIN DOWN 3 MINUTES
MAX. INBALANCE 2 TUBES
SPEED UP TO 4000RPM
RCF 2000 AT 4000RPM

CLAY-ACAMS, INC. #R1004 \$324
SPEED RANGE 800 TO 3040 RPM
SAMPLES 12 (MAX), 15 ML EACH
DIMENSIONS 15X12X26 INCHES
CONFIGURATION BENCH MOUNT

BECKMAN INSTRUMENTS #R152 \$245 111
SPINNING RATE 0 TO 16000 RPM
ACCOMMODATION4ML TUBES(20)
DIMENSIONS 9X7X7 INCHES
WEIGHT 13 LBS(SHIPPING)
POWER 115V 50/60 HZ
DEVELOPED G 13000

603015 VACUUM CLEANER LS943 11111
603015
603015 COLLECT GENERAL LABORATORY WASTE MATERIALS, DIRT, DUST AND DEBRIS
603015
603015 VACUUM BAG VOLUME 0.36CU.FT. (10.2LITERS)
603015
603016 CLINOSTAT LS944 11111
603016
603016 MOTOR DRIVEN MOUNTING STRUCTURE TO PROVIDE CONTINUOUS ROTATION
603016 OF PLANT SAMPLES
603016
603016 ROTATION RATES TRD
603016 PLATFORM LIGHTING TRD
603016
338000 COMPACTOR, WASTE SOLIDS LS945 11111
338000
338000 COMPACTS WASTE SOLIDS FOR EASY HANDLING
338000
338000 COMPACTION(MIN) 1/4 INITIAL VOLUME
338000

SEARS, ROEBUCK AND CO 5200 112
VOLUME 0.2 HOURS (7.2 FT003)
COMPACTION (FINAL VOL) 1/4 INITIAL VOL

000470 CONSOLE, BEHAVIORAL MEASUREMENTS L5946 11114
000470
000470 DEVICE THAT MEASURES MENTAL PROCESSES, ATTITUDES, ETC. SUCH AS
000470 COGNITIVE PROCESSES, INDIVIDUAL BEHAVIORAL TRAITS, GROUP
000470 BEHAVIORAL TRAITS, ETC. CONSOLE INCLUDES GENERAL PURPOSE DISPLAY
000470 AND GENERAL PURPOSE RESPONSE KEYBOARD
000470
000470 PROCESSES MEASURED COMPLEX PERCEPTUAL
000470 CONCEPTUAL AND THINKING ABILITY
000470 MEMORY
000470 REACTION TIME-SIMPLE AND COMPLEX
000470 INDIVIDUAL BEHAVIORAL TRAITS
000470 GROUP BEHAVIORAL TRAITS

102500 DIGITAL COMPUTER L5947 11111
102500
102500 PROVIDE AUTOMATIC DATA ACQUISITION OF EXPERIMENTAL DATA, CONTROL
102500 LIGHTING AND ANIMAL FEEDING AND GENERAL LABORATORY ENVIRONMENT
102500 MONITORING. THE COMPUTER WILL SUPPORT ANALYTICAL TOOLS SUCH AS
102500 MASS SPECTROMETERS, GAS ANALYZERS AND CATALYSIS EQUIPMENT. ALSO
102500 DRIVES ONE OR MORE CRT INTERACTIVE DISPLAY CONSOLES
102500
102500 MEMORY CAPACITY 16BIT 16,000 WORD
102500 CYCLE TIME MICROSECOND
102500 PERIPHERAL CAPABILITIES DIRECT MEMORY ACCESS
102500 FLEXIBLE INTERRUPT STRUCTURE
102500 SINGLE BIT SET AND SENSE LINE
102500 CAPABILITY OF UP TO 64 LINES
102500 EACH

DATA GENERAL CORP 0NOVA 800 518450 111
WORD LENGTH 16 BIT
CYCLE TIME 800 NANOSEC.
CORE MEMORY 32K MAX

DATA GENERAL CORP 0SUPERNOVA 5C 95600
WORD LENGTH 16 BIT

DATA GENERAL CORP 0NOVA 120C 511900
WORD LENGTH 16 BIT
CYCLE TIME 1200 NANOSEC.
CORE MEMORY 32K MAX

VARIAN DATA MACHINES 0520/1 111
MEMORY CYCLE TIME 1.9MICROSEC.
MEMORY EXPANDABLE 4096BYTES(8BITS) TO
32,768BYTES
REGISTERS 12
OPERAND PRECISION UP TO 32 BITS
ADDRESS REFERENCE 0,16,24,OR 32 BIT LEVEL
NOTE: HAS FULL RANGE INTERFACE
HARDWARE

VARIAN DATA MACHINES 0R-620/1 111
FASIC COMMANDS OVER 100
ADDRESSING MODES 6
MAX. WORDS 32,768
WORD LENGTH 16 OR 10 BIT
REGISTERS 9
NOTE: HAS FULL RANGE INTERFACE
HARDWARE

CLARY DATACOMP SYSTEMS INC 0404 111
SIMULTANEOUS TERMINALS 16
INPUT TELETYPEWRITER KEYBOARD
ACCUMULATOR 64 BIT
REGISTERS 16 BIT INDEX (2 EACH)
WORD LENGTH 16, 32, 48, 64 BITS
MEMORY CAPACITY 1624 16 BIT WORDS, 4096 16 BIT
WORDS, OR ADDITIONAL 4096 16 BIT
WORDS. TOTAL 65536 WORDS

GENERAL AUTOMATION INC. 010-30 111
FUNCTION SUPERVISE SMALLER COMPUTERS
MEMORY OR CORE
CYCLE TIME 1.2 MICROSEC.

GENERAL AUTOMATION INC. 0SPC-16 111
MEMORY 16K
WORD LENGTH 16 BITS
READ/WRITE CYCLE MEMORY TIME 800 TO 1440 NANOSEC.
READ ONLY MEMORY 400 TO 720 NANOSEC.
INPUT/OUTPUT TIME 1.6 TO 2.0 MICROSEC.
DMA TRANSFER RATE 0.694E+06 TO 2.5E+06

GENERAL AUTOMATION INC. 0SPC-12 111
MEMORY 4K TO 16K
WORD LENGTH 8 BITS
CYCLE TIME 2.16 MICROSEC
STORED PROGRAM EXECUTION RATE 0.23E+06/SEC
I/O TRANSFER RATE 0.40E+06/SEC
REGISTERS 4 12 BIT REGISTERS
ACCUMULATORS 4 12 BIT

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

DIGITAL EQUIPMENT CORP. DMP-R/F 822000 111
MEMORY 4096 CORE EXPANDABLE TO 32,768 WORDS
WORD 12 BIT
CYCLE TIME 1.5 MICROSEC.
103500 ANALOG TO DIGITAL CONVERTER LS948 11111
103500
103500 SUPPORT DIGITAL COMPUTER. INCLUDES LO-LEVEL MULTIPLEXER. UNIT
103500 SHOULD BE A PROGRAMMABLE GAIN 12 BIT TYPE
103500
105004 DISK STORAGE SYSTEM L5949 11111
105004
105004 STORES DATA IN SUPPORT OF COMPUTER
105004
105004 STORAGE CAPACITY 5,096 WORDS
105004
ANALOG DEVICES INC. #ADC-#5 879 111
RESOLUTION 8 BITS
ERROR +/-0.2%
CONVERSION TIME 1 MILLI SEC.
ANALOG DEVICES INC. #ADC-1202 6305 111
RESOLUTION 12 BITS
ERROR +/-0.0125%
CONVERSION TIME 25 MICROSEC.
ANALOG DEVICES INC. #ADC-160 81350 111
RESOLUTION 16 BITS
ERROR +/-0.0015%
CONVERSION TIME 400 MICROSEC.
TENNELEC #TC-500 86000 111
CONVERSION TIME 3 MICRO SEC
CONFIGURATION COMPATIBLE WITH NIM
CAPACITY 2048 CHANNELS
TENNELEC #TC-520 8740 111
CHANNELS 8
CHANNEL DWELL TIME 300 NSEC
CONFIGURATION NIM COMPATIBLE
ROMBE AND SCHWARTZ #UCM 210/1 111
(221.5603.02)
ANALOG INPUT RANGE +12 TO -12 VOLTS
SAMPLING RATE > 10000 MEAS/SEC (POSITIVE)
..... > 17000 MEAS/SEC (NEGATIVE)
DATA OUTPUT 1-1 DECADES (BCD CODE PARALLEL)
..... 0-1200 PULSES (SERIAL)
424000 MAGNETIC TAPE RECORDER LS940 11111
424000
424000 RECORD DATA
424000
424000 DATA RATE 30+ RPS
424000 TOTAL DATA STORED 16,000 RITS
424000
AMPLEX CORPORATION #AR-200 121
BANDWIDTH 100 HZ TO 3125 HZ THROUGH
..... 300 HZ TO 250 KHZ
TAPE WIDTH 0.4 INCH
TAPE SPEED 1.875 TO 60 IPS
RECORDING TIME 8 MINUTES TO 4 HRS AND 16 MINS
FORMAT DIGITAL
TRACKS 8 DIGITAL, 7 ANALOG
AMPLEX CORP. #AR-1700 828,000 112
TRACKS 14
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KR/I/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 6.2F*10
BANDWIDTH 2 MR/SEC/T
SANGAM ELECTRIC #SABER 111 826,500 111
TRACKS 14
TAPE WIDTH 1 INCH
TAPE SPEEDS 8 SELECTABLE SPEEDS PREP 15/16
..... TO 120 IPS
FREQUENCY RESPONSE 400 HZ TO 2.0 MHZ
RECORDING RATE 600 ROPS AT 120 IPS SERIAL MODE
WEIGHT 100 POUNDS
DIGI-DATA CORPORATION #1600 82690 111
TAPE SPEED 25, 18.75, 12.5 IPS
TRACKS 7 OR 9 TRACK
DATA DENSITY 1600 CPI PHASE ENCODED
..... 200, 556, 800 CPI NRZI
TAPE 0.5 INCH, 1.5 MIL, 1200 FEET
DIGI-DATA CORPORATION #1700/PDP-11/7-9 85250 111
TRACK NRZI
TAPE SPEED 45, 37.5, 25, 18.75, 12.5 IPS
TAPE 0.5 INCH, 1.5 MIL, 10P/ANSI
..... COMPATIBLE, 10.5 INCH REEL
TRACKS 7 OR 9
DATA DENSITY PHASE ENCODED
..... COMPATIBLE



Space Division
Rockwell International

HEWLETT PACKARD	#7570H/C	\$ 4600	111
TAPE FORMAT 800, 556, OR 200 CPI KHZ AND 1600 CPI PHASE-ENCODED			
CHANNELS 7 OR 9			
TAPE SPEED 10 TO 45 IPS			
TAPE 0.5 INCH, 1.5 MILS, IBM/ANSI COMPATIBLE			
BCRG WARNER	#PERT		121
TRACKS 30			
TAPE SPEED UP TO 1000 IPS			
TAPE WIDTH 1/2 INCH			
TAPE LENGTH 2400 FEET			
PACKING DENSITY 15 KB/1/7			
SIGNAL/NOISE 24 DB			
DATA CAPACITY 1.3E+9			
BANDWIDTH 6-15 MH/SEC/T			
LEACH	#MTR 7000		311
TRACKS 12			
TAPE SPEED 120 IPS			
TAPE WIDTH 1 INCH			
TAPE LENGTH 9200 FEET			
PACKING DENSITY 16.7 KB/1/7			
SIGNAL/NOISE 22 DB			
DATA CAPACITY 2.2E+10			
BANDWIDTH 2 MH/SEC/T			
HEWLETT PACKARD	#35550	\$10200	111
BANDWIDTH 300 KHZ			
CHANNELS 7			
RECORDING FORMAT DIRECT OR FM			
HEWLETT PACKARD	#3555C	\$14700	111
BANDWIDTH 300 KHZ			
CHANNELS 14			
RECORDING FORMAT DIRECT OR FM			
HEWLETT PACKARD	#3550A-C11	\$21800	111
BANDWIDTH 500 HZ TO 2 MHZ			
CHANNELS 14			
RECORDING FORMAT DIRECT OR FM			
HEWLETT PACKARD	#3560A/13065A/13063A	\$4796	111
CONFIGURATION RACK MOUNTED			
TAPE SPEED 15, 3 AND 3/4, 15/16 IPS			
CHANNELS 4			
RECORDING FORMAT FM			
PASSBAND 5 KHZ			
S/N RATIO 48 DB			
TRACKS 30			
HONEYWELL	#5600	\$9730	
PORTABLE TAPE RECORDER			
CHANNELS 7			
SELECTABLE TAPE SPEED RANGE 15/16 TO 60IPS			
MAXIMUM BANDWIDTH(DIRECT) 300KHZ			
PACKING DENSITY UP TO 600PI			
WEIGHT 70LBS (32KG)			
INPUT VOLTAGE 28VDC			
HONEYWELL	#96	\$17420	
CHANNELS 7			
SELECTABLE TAPE SPEED RANGE 15/16 TO 240IPS			
REEL SIZE 16IN			
MAXIMUM BANDWIDTH(DIRECT) 2M HZ			
TAPE WIDTH 1/2 IN			
APPEX CORP	#AP 70C	\$29,180	511
TRACKS 14			
TAPE SPEED 60 IPS			
TAPE WIDTH 1 INCH			
RECORD MODE DIRECT			
PACKING DENSITY 20 KB/1/7			
SIGNAL/NOISE 20 DB			
DATA CAPACITY RE+10 BITS			
BANDWIDTH 1 MH/SEC/T			
104180 CRT DISPLAY		LS951	41111
104180			
104180 DISPLAY GRAPHIC AND ALPHANUMERIC DATA			
104180			
RESEARCH INC	#3300	\$1580	
FORMAT 24 LINE SX72 OR 80 CHARACTERS			
..... 24 LINE SX40 CHARACTERS			
..... 12 LINE SX72 OR 80 CHARACTERS			
REFRESH RATE 60 HZ			
TRANSFER RATE 110 TO 2400 BAUD			
CHARACTER FORM 5X7 DOT MATRIX			
DIMENSIONS 15-1/2HX13-1/2HX23-1/2D INCH			
WEIGHT 39LB (17.7KG)			
HEWLETT PACKARD	#2600A	\$3580	111
RESEARCH INC	#R12-33C1	\$1555	111

ORIGINAL PAGE IS
OF POOR QUALITY



SCREEN FORMAT 12 LINES X 72 OR 80 CHARACTERS,
24 LINES X 40 CHARACTERS
CHARACTER FORMAT 5 X 7 DOT MATRIX
TRANSFER RATE 110 TO 2400 BAUD, 10 OR 11 BIT
CHARACTERS
MODES HALF OR FULL DUPLEX-SWITCHABLE
LOCAL OR REMOTE

RESEARCH INC #DPS-R12 \$21,700 111
COMPLETE DATA DISPLAY SYSTEM MADE UP OF CRT ALPHANUMERIC DISPLAY
DEC-PDP-11 MINI-COMPUTER, #1301 OPERATORS CONTROL CONSOLE,
PRINTER/CARD READER, #812-4A ANALYSETER AND #812-13 UNIDRIVER
CHANNELS 32
MEMORY SIZE 64 WORDS

TEKTRONIX INC #4002A/021-0033-00/ \$10975 111
021-C0XX-00/4901/
4951

ALPHANUMERIC FORMAT 39 LINES-85 NORMAL/ITALIC CHAR,
1 LINE- 84 CHAR IN SCRATCH AREA
CHARACTER SET 96 UPPER/LOWER CHAR (ASCII)
CHARACTER SIZE 70 X 90 MILS
CHARACTER GENERATION 7 X 9 DOT MATRIX
CURSOR PULSATING 7 X 9 MATRIX
GRAPHIC INPUT MODE 1024(X) X 1024(Y) ADDRESSABLE,
1024(X) X 768(Y) VIEWABLE PTS;
JOYSTICK CONTROLLED, CROSS-HAIR
CURSOR
DISPLAY MEDIUM 11 IN DIA CRT
DISPLAY AREA 8.3IN.HOR, 26.1IN.VERT.
ACCUMULATIONS 4
CYCLE TIME 300 NANOSEC
CORE MEMORY 32 K

BUNKER-RAMC CORP #2217-12CPT 111
CHARACTER CAPACITY (MAX) 960
CHARACTERS/LINE (MAX) 80
LINES/DISPLAY (MAX) 24
CHARACTERS REPERTOIRE 82 OR 92
VIEWING AREA 8.75 X 6.25
REFRESH AREA 54 FRAMES/SEC

096000 AUTOMATIC COLONY COUNTER LS952 11111
096000
096000 DETERMINE THE NUMBER OF VIABLE COLONIES GROWING ON AN ALGUR
096000 SUBSTRATE
096000
AMERICAN OPTICAL CORP #3330 \$115 112
REGISTER MODEL AO 3348 USED TO OBTAIN
CUMULATIVE RECORDING OF COUNTS
TYPE QUERCY COLONY COUNTER, DARK-
FIELD ILLUMINATION APH
DIMENSIONS 10X11X10 INCHES

222000 DEIONIZER, WATER LS953 11111
222000
222000 PROVIDES A MEANS OF REMOVAL OF CATIONS AND ANIONS FROM WATER
222000 SUPPLIES REQUIRED FOR SPECIALIZED LABORATORY ANALYSIS
222000
222000 SOLIDS REMOVAL CAPABILITY.....<1 PPM
222000 IONS REMOVAL CAPABILITY.....<1 PPM
222000

VWR SCIENTIFIC DIV #24804-002/24805-005 \$65 112
CONFIGURATION PORTABLE
TYPE DETERMINE RESIN AND FILTERS
CRYSTALAB DEMINERIZER
PRODUCTIVE CAPACITY 5 GALLONS/HOUR
MINIMUM GRAIN REMOVAL 1000 GRAINS AS NaCl
1000 GRAINS AS CaCO3

VWR SCIENTIFIC DIV #24809-006/24809-050 \$127 112
CONFIGURATION PORTABLE
TYPE LAB FLOW QUIKPURE DEMINERALIZER
USES ION EXCHANGE RESINS
USEABLE EXCHANGE CAPACITY 1500 GRAINS TOTAL IONIC SOLIDS,
CALCULATED AS NaCl
PRODUCTIVE CAPACITY 60 GALLONS/HOUR
PURITY METER 100,000 TO 10,000,000 OHMS

CRYSTALAB #DF500 \$125 111
FLOW RATE..... 20 TO 120 GAL /HR

104180 ELECTROPHYSIOLOGY DISPLAY LS954 11114
104180
104180 DEDICATED DISPLAY FOR PHYSIOLOGICAL DATA TRANSMITTED BY THE
104180 ELECTROPHYSIOLOGY BACKPACK
104180
605006 ELECTROMETER LS955 11111
605006
605006 MEASURE SMALL VOLTAGES AND CURRENTS FOUND IN NERVES AND MUSCLES.
605006
605006 DC VOLTAGES MEASURED.....+/-10MICROVOLTS TO +/-1VOLT
605006 AC CURRENTS MEASURED.....+/-1 PA TO +/-3 MA
605006 DRIFF.....+/-5V/DAY AFTER 10MIN WARMUP
605006 OUTPUT VOLTAGE.....+/-10V
605006

HEWLETT PACKARD #425A \$695 112



CONFIGURATION BENCH EQUIP
VOLTAGE RANGE +/- 10 MICROV TO +/- 1 V FS. 11
STEPS, 1, 3, 10 SEQUENCE
CURRENT RANGE +/- 10 PA TO +/- 3 MA FS. 18
STEPS, 1, 3, 10 SEQUENCE
INPUT IMPEDANCE 1 MEGOHM +/- 3% (VOLTAGE RANGE)
1 MEGOHM TO 0.33 OHM (CURRENT
RANGE-DEPENDS ON RANGE)
ACCURACY < 3% OF RANGE
OUTPUT 0 TO 1 V FOR FS READING

RHODE AND SCHWARTZ #UIG (203-5111-021) \$1120 111
VOLTAGE RANGE 0.2 MICROVOLT TO 320 VOLTS
CURRENT RANGE 0.01 MICROAMP TO 300 MA
ACCURACY +/- 1.5%

152500 ELECTROPHORESIS APPARATUS LS956 11111
152500
152500 SEPARATES PROTEIN AND AMINO ACID CONSTITUENTS IN SERUM, PLASMA,
152500 URINE OR SPINAL FLUID FOR QUANTITATIVE ANALYSIS.
152500
152500 VOLTAGE 0 TO 500 V DC
152500 SUBSTRATE STARCH BLOCK, PAPER OR GEL
152500

ARTHUR H THOMAS CO #4306-810 \$710 112
CONFIGURATION BENCH/LOOSE EQUIP
TYPE HORIZONTAL STRIP APPARATUS FOR
PAPER, ACETATE STRIP OR GEL
BLOCK ELECTROPHORESIS
CURRENT 0 TO 100 MA
VOLTAGE 0 TO 500 VOLTS

BECKMAN/SPINCO DIV #R-100 \$3300-6700 112
CONFIGURATION BENCH/LOOSE EQUIP

BUCHLER INSTRUMENTS #27384-0CR \$505 112
CONFIGURATION BENCH/LOOSE EQUIP
CURRENT 0 TO 200 MA AND 0 TO 20 MA
VOLTAGE 0 TO 500 VOLTS

BUCHLER #VFR 30136-001 \$570 111
SAMPLE SIZE 25MG
VOLUME OF BUFFER SOLUTION 2 TO 10ML

170200 FREEZER, GENERAL LS957 11111
170200
170200 STORAGE OF SERUM, PLASMA, SPECIMENS AND ORGANISMS
170200
170200 OPERATING TEMPERATURE -4F (-20C)
170200 TEMPERATURE TOLERANCE +/- 4F (+/- 2C)
170200 STORAGE VOLUME 4 CU. FT. (0.11CU.M.)
170200

REVCO #ULT-1P5A \$1060 111
TEMP RANGE -7 TO -85 C
CAPACITY 1.5 CU FT
DIMENSIONS 34X24X37 INCHES
WEIGHT 270 LBS (SHIPPING)
POWER 115V 60 HZ 7A
CAPACITIES AVAILABLE 1.5 TO 17 CU FT (VAR. PRICES)

170200 FREEZER, LOW TEMP LS958 11111
170200
170200 STORAGE OF EXPERIMENT SPECIMENS
170200
170200 OPERATING TEMPERATURE -94F (-70C)
170200 STORAGE VOLUME 1 CU.FT. (0.028 CU.M.)
170200

REVCO #ULT-1P5A \$1060 111
TEMP RANGE -7 TO -85 C
CAPACITY 1.5 CU FT
DIMENSIONS 34X24X37 INCHES
WEIGHT 270 LBS (SHIPPING)
POWER 115V 60 HZ 7A
CAPACITIES AVAILABLE 1.5 TO 17 CU FT (VAR. PRICES)

040700 REFRIGERATOR LS959 11111
040700
040700 STORE SERUM AND PLASMA
040700
040700 OPERATING TEMPERATURE 32 TO 40F (10 TO 4C)
040700 STORAGE VOLUME 1 CU.FT. (0.028 CU.M.)
040700

MATHESON SCIENTIFIC #32035-25 \$685 112
CAPACITY 12.8 CU FT
TEMPERATURE -15 TO -7.8 C
SHIELDED AREA NONE
DIMENSIONS 32X31X60 INCHES
EXPLOSION PROOF YES
COMPRESSOR 0.7HP
POWER 115V 50-60 HZ 3.5A
WEIGHT 275 LBS

DILLON-LILLY #512-D \$975 112
CAPACITY 12.3 CU FT
TEMPERATURE 3.7 C
SHIELDED AREA NONE
DRAWERS 12
COMPRESSOR 0.25 HP
POWER 115 V 60 HZ
WEIGHT 420 LBS
DIMENSIONS 36X56.5X28.5 INCHES

ORIGINAL PAGE IS
OF POOR QUALITY



JORDAN WFT-1-AR 11177 112
TEMPERATURE 2.22-4.11C IN 1.11C INCH
CAPACITY 22.4 CU FT
COMPRESSOR 0.25 HP
SHIELDED AREA NONE
CRAIN NONE REQUIRED
POWER 115V 60 HZ 6A
DIMENSIONS 78X28.5X31.625 INCHES

BIRKEN MANUFACTURING CO #STM-HL(DC-AC) 1335 112
CAPACITY 8 FT005
VOLTAGE 115 VAC, 12, 24, 32 VDC

FRICID CAB #14 1725 112
CAPACITY 14 CL FT
TEMPERATURE -6.7 TO +15 C
SHIELDED AREA NONE
EXPLOSION-PROOF YES
SEPARATE FREEZING COMP. YES
COMPRESSOR 0.25HP
POWER 115V 60 HZ 3.6A
DIMENSIONS(MAIN CHMRP) 26.5X20X36.5 INCHES

048700 REFRIGERATOR, RADIO ISOTOPE STORAGE LS960 1113-
048700
048700 STORE LIQUID AND SOLID RADIOACTIVE MATERIALS
048700 OPERATING TEMPERATURE.....-41/-46F5+/-2C
048700 VOLUME.....1CU.FT,10.028CU.M
048700 ISOTOPE STORED.....H3,C14,FE59,CRS1,1131,CO45
048700
435000 GAS ANALYZER, CO2 SPECIFIC LS961 11111
435000
435000 THIS DEVICE IS USED TO MONITOR ATMOSPHERIC CARBON DIOXIDE LEVELS
435000
435000 PARTIAL PRESSURE RANGE.....130 - 2700PA11-20 MM HG
435000

BECKMAN INSTRUMENTS #6501 AND #601 \$2750 111
PARTIAL PRESSURE CO2 11 THROUGH 30 MM HG
PARTIAL PRESSURE O2 0 TO 1000 MM HG
POWER 12-32 VDC INT OR EXTERNAL
SIZE EACH UNIT 6.25X4.25X3.5 INCHES
WEIGHT (2 UNITS)..... 9 LBS
2 UNITS CAN BE COMBINED IN 1 PKG

INSTRUMENTATION LABORATORY #IL 313 \$6990 111
PH 0.0 TO 9.0
PCO2 0 TO 200 MMHG
PO2 0 TO 2000 MMHG
AUTOMATICALLY PROGRAMMED YES
SAMPLE TIME 45 SEC
PRECISION PH +/-0.003
PCO2 +/-1MMHG
PO2 +/-10MMHG
SAMPLE SIZE 0.4ML (AUTOMATED)
POWER 115V 50/60 HZ

BECKMAN INSTRUMENTS INC #CUSTOM ASSEMBLY 114
GASES CO, CO2, O2, N2, H2S, SO2
SENSITIVITY 1-2 PPM FS
AMBIENT TEMP LIMITS -2 TO 120 F
AIR REQUISITS 2.5 TO 4 SCFM AT 30 PSIG

BECKMAN INSTRUMENTS #GC-95(977563) \$6870 111
SENSITIVITY:
ELFOTACON CAPTURE DETECTOR..... .5 PICOGRAMS LINDANE
..... .53A ML-HG
THERMAL CONDUCTIVITY DET 1200C DIMBAT-PROTER-STROSS
FLAME EMISSION DET..... 10 PICOGRAMS PHOSPHOROUS
..... 4C PICOGRAMS SULFUR
OVEN TEMPERATURE -66 C TO +400 C
DIMENSIONS:
OVEN 27.75X25.75X26.25 INCHES
ELECTRONICS 19.75X26X23 INCHES

074700 GAS CHROMATOGRAPH LS962 11111
074700
074700 MEASURES CONCENTRATION OF GAS, LIQUID, AND SOLID CONSTITUENTS OF
074700 BIOLOGICAL SAMPLES
074700 OPERATING TEMPERATURE.....-40 TO 752F1-40 TO 400C
074700 CARRIER GAS.....HELIUM
074700 DETECTORS.....HYDROGEN FLAME
.....THERMAL CONDUCTIVITY
074700 FRACTION COLLECTION SYSTEM
074700

BECKMAN #GC-2A \$2050 111
DETECTOR THERMAL CONDUCTIVITY
TEMPERATURE CONTROL RANGE 40C TO 240C
CONTROL ACCURACY +/-0.1

PERKIN-ELMER #3920 \$3000 111

CARLE INSTRUMENTS, INC. #09500 \$1185
TEMPERATURE RANGE AMBIENT TO 200 DEG C
DETECTOR DUEL FEED SINGLE HYDROGEN FLAME
SENSITIVITY 5E-12 AMPERES
POWER REQUIREMENTS 115 VAC, 5C TO 60 HZ

CARLE INSTRUMENTS INC #0001 \$675 111



ACCURACY +/- 0.5 PERCENT
GASES DETECTED CO, CO2, N2, AR, O2, CH4
SAMPLE SIZE 1 PM
REPETITIVE SAMPLING TIME 5 TO 60 MIN.

CABLE INSTRUMENTS, INC. WAGC-311 \$2675
CONTROL SENSITIVITY < 0.001 DEG C
TEMPERATURE RANGE TO 459 DEG C (1 DEG INCREMENTS)
HEAT CAPACITY-TOTAL 500 WATTS
INLETS DUAL ON-COLUMN
IONIZATION DETECTORS HYDROGEN FLAME
VALVING MICRO DP MINI VOLUME VALVES
WEIGHT 53 POUNDS
POWER REQUIREMENTS 115 VAC

VARIAN AEROGRAPH N1520 \$5200 112
COLUMN OVEN TEMP -59 TO 399 C
TEMP STABILITY BETTER THAN .1 PER CENT
FLAME DET SENSITIVITY 0.015 COULOMBS/GRAM
MIN DETECTABILITY 5 E-12 GRAM/SEC
NOISE LESS THAN 3XE-14 A
AMPLIFIER SOLID STATE
DIMENSIONS 26X18X22 INCHES

BECKMAN INSTRUMENTS 8680G \$7350 111
CARRIER H2
DETECTOR FLAME
RANGE 0-1 PPM FULL SCALE
ANALYSIS RATE 4 TO 6/HR (SELECTABLE)
REPEATABILITY5% OF FS
LINEARITY 1% OF FS
ZERO DRIFT AUTOMATIC ZERO COMPENSATION
DIMENSIONS 17X40X20 INCHES
WEIGHT 250 LBS (SHIPPING)
POWER 107-127V 50/60 HZ 500W

211300 GAS ANALYZER, WATER VAPOR LS963 11111

211300 MONITOR WATER VAPOR CONTENT
211300 MOISTURE RANGE 0 TO 50% RELATIVE HUMIDITY
211300 (0.001 TO 20000 MICROGRAMS/LITER)
211300 TEMPERATURE RANGE -160 TO 1400 (-110 TO 600)
211300

507900 KIT, BENCH CHEM ANAL LS964 11111

507900 KIT CONTAINS THE TOOLS AND EQUIPMENT TO HANDLE CHEMICALS AND
507900 BIOLOGICALS DURING VARIOUS MANUAL PROCEDURES GENERALLY PERFORMED
507900 WITHIN THE GLOVE BOX.

507900 TOOLS SOLID TRANSPORT TOOLS, GRAVITY
507900 INDEPENDENT PIPETTES, VIALS,
507900 BOTTLES, AND TEST TUBES,
507900 CHEMICALS, STOPPERS, FILTERS,
507900 AND SAFETY SHIELDS

507900 KIT, HEMATOLOGY LS966 11111

507900 PROVIDE TOOLS FOR SAMPLING, HANDLING, TRANSFERRING AND ANALYZING
507900 BLOOD

507900 ITEMS HEMATOMETER KIT; 10 LAPROA
507900 DISPOPIPETTES; COVERSLIPS;
507900 SLIDES; WBC AND RBC DILUENTS;
507900 CRITOSEAL; HEMOCRIT TUBES (MICRO,
507900 HEPARINIZED); HEMOCRIT TUBES
507900 (MICRO, PLAIN); BLOOD DILUTING
507900 PIPETTES (WBC, RBC); COMBISTIX
507900 (URINE TEST STRIPS); LUER
507900 ADAPTERS; VACUTAINER; VACUTAINER
507900 VACUTAINER-NEEDLE UNIT;
507900 VACUTAINER TUBES, ASSORTED (2ML);
507900 SYRINGE (LARGE), PEDIATRIC
507900 ALCOHOL SWABS; LANCETS; NEEDLE,
507900 25GA, 5/8IN; NEEDLE 22GA, 1 1/2IN

507900 KIT, LINEAR MEASUREMENT LS967 11111

507900 DETERMINE SIZE, AMPLITUDE, DISTANCE, CIRCUMFERENCE, ETC

507900 ITEMS RULERS; TAPES; LINEAR COMPARATORS;
507900 GRIDS; CALIPERS, INSIDE AND
507900 OUTSIDE; VERNIER CALIPERS;
507900 MICROMETERS, INSIDE AND OUTSIDE

507900 KIT, MICROBIOLOGY LS968 11111

507900 PROVIDE TOOLS TO FACILITATE GROWING AND ANALYZING MICROBIAL
507900 ORGANISMS

507900 ITEMS INOCULATING LOOP; INOCULATING
507900 NEEDLES; COTTON SWABS, PACKAGES,
507900 STERILE; SYRINGE SML, STERILE;
507900 SLIDES, MICRO; ALCOHOL; TUBES,
507900 15X70MM, STERILE, CAPREO;
507900 ZEPHYRUS, TINCTURE, 1:500;
507900 BACTECINERATOR/STERILIZING UNIT;
507900 THIOGLYCOLLATE, TUBES; STUART
507900 TRANSPORT MEDIA, VIALS; TSA
507900 SLANTS

ORIGINAL PARTS
OF POOR QUALITY



LKB INSTRUMENTS INC. #HR15 ACCESS. KIT I 111
SPECIMEN HOLDERS 5
SYRINGE 2 ML
LKB INSTRUMENTS INC. #8816 ACCESS. KIT II 111

907900 KIT, ORGANISM HOLDING AND MANAGEMENT LS969 11111
907900
907900 PROVIDE TOOLS AND DEVICES USED IN THE HOLDING AND HANDLING OF THE
907900 ORGANISMS
907900 ITEMS FOR SMALL VERTEBRATES.....GLOVE POWDER; CAGE SHIELD;
907900 PLASTIC; PLASTIC LINER WITH PAD;
907900 FOOD PELLET DISPENSER; ALCOHOL
907900 SWABS; TOWELS; PAPER; DISC;
907900 PLASTIC BAGS FOR EXPENDIBLES;
907900 PLASTIC BAGS FOR NONEXPENDIBLES;
907900 ORGANISM TRANSFER CAPSULE;
907900 ANIMAL TAGS
907900 ITEMS FOR PLANTS.....WATERING DEVICES (SPRAY BOTTLE);
907900 LABELS; FERTILIZER PACKETS;
907900 STAKES (10CM)
907900 KIT, PLANT TOOLS LS970 11111
907900
907900 PROVIDE TOOLS FOR VARIOUS PLANT MANIPULATIONS
907900
907900 ITEMS.....SCISSORS, TWEEZERS, SPATULA,
907900 APPLICATOR, SPLINTS, TAPE
907900 DISPENSER, HYPODERMIC NEEDLES,
907900 SYRINGES, FORCEPS, PLIERS,
907900 SCREWDRIVER, SCALPELS, SEALERS
907900 KIT, GENERAL TOOL LS971 11111
907900
907900 MECHANICAL AND ELECTRICAL TOOLS AND HARDWARE TO PROVIDE
907900 CONVENTIONAL DIAGNOSTIC, MAINTENANCE, AND SERVICE FUNCTIONS.
907900
907900 TOOLS.....HAMMER, WRENCHES, PLIERS, SCREW
907900 DRIVERS, DRILL, ADHESIVE TAPE,
907900 WIRE CUTTERS, WIRE TIES, WIRE,
907900 LUBRICANTS, FLASHLIGHT, SCISSORS,
907900 MULTIMETER, FASTENERS, CLAMPS,
907900 LAMP
907900
LABORATORY SPECIALISTS LTS #1602 #140 112
TOOLS.....33 SELECTED TOOLS INCLUDING
WRENCHES, PLIERS, TAPPING TOOL,
SCREW DRIVERS, NUT DRIVERS,
RIVETS, TERMINALS, ETC
CASE.....19X14X5.5 INCHES
235700 LYOPHILIZER LS972 11111
235700
235700 PROVIDE FREEZE DRYING OF SMALL SPECIMENS
235700
235700 SAMPLE CAPACITY.....20-100 IN. (16CC)
235700 TEMPERATURE RANGE.....-40F TO 252F (-40 TO 122C)
235700
SCIENTIFIC PRODUCTS DIV #07C27-5 #1390 111
CONFIGURATION BENCH EQUIP
TYPE FRONT LOADING, SINGLE SHELF
TEMPERATURE RANGE -40 TO 122 C (-40 TO 250 F)
VIRTIS #UNITRAP #1300 111
TEMPERATURE (MIN).....-65F (-55C)
DRUM SIZE.....8 IN. (21CM) DIA. X 10 IN. (25.4CM) HIGH
CONFIGURATION..... FLOOR MOUNT, PORTABLE
VACUUM PUMP NOT INCLUDED
SARGENT WELCH SCIENTIFIC CO. #S-288B1-80 #3510 111
LOW TEMPERATURE -85F
CAPACITY 12 LITERS
VACUUM 5 MICRONS
PUMP CAPACITY 140 LITERS/MIN.
023500 MASS MEASUREMENT DEVICE, MACRO LS973 11132
023500
023500 MEASURE MASS OF ITEMS SUCH AS BAGS OF FOOD, BEVERAGE CONTAINERS,
023500 CONTAINERS OF URINE AND FECES, LARGE SPECIMENS, ETC
023500
023500 MASS RANGE.....0.22 TO 66LB (0.1 TO 30KG)
023500
SOUTHWEST RESEARCH INSTITUTE #074 (SKYLAB) 511
SPECIMEN MASS MEASURING DEVICE
ELEMENTS OF ASSEMBLY SPECIMEN TRAY SUPPORTED ON
PLATE-FLUOR SPRINGS, PERIOD/
TEMPERATURE MEASURING ELEC-
TRONIC SUBSYSTEM
023500 MASS MEASUREMENT DEVICE, MICRO LS974 11132
023500
023500 MEASURE SMALL TEST SPECIMENS
023500
023500 MASS MEASUREMENT RANGE.....1MG TO 100 G
023500
306480 OPTISCAN-FIELD AND FIXATION POINT RECORDER LS975 11111
306480
306480 THIS INSTRUMENT PERMITS THE RECORDING OF THE VISUAL FIELD AS
306480 THE SUBJECT'S FIXATION WITHIN THE FIELD WITHOUT RESTRICTION OF
306480 HEAD MOVEMENTS
306480
306480 ACCURACY.....2DEG OF ARC



396440
004400 CATALYTIC OXIDIZER SYSTEM LS977 11235
004400
004400 REMOVES UNDESIRABLE COMBUSTION PRODUCTS
004400
004400 OPERATING TEMPERATURE.....70CF(1370C)
004400 AIR FLOW RATE.....2CFM AT 14.7PSIA AND 70F
004400
048050 BLOOD PRESSURE CUFF SYSTEM LS978 11132
048050
048050 PRESSURE MEASURING TRANSDUCER CUFF USED TO MEASURE BLOOD PRESSURE
048050 OF PRIMATES. SYSTEM INCLUDES AIR PUMP FOR AUTOMATED MEASURING.
048050 PRESSURE CUP AND TRANSDUCER, AND SIGNAL CONDITIONER.
048050

MARTIN CO #M092-BPMS (SKYLAB) 511
CUFF VOLUME 20 CU.IN.
BPMS MAX OPERATING PRESSURE 210 PSIG
BPMS MIN OPERATING PRESSURE 90 PSIG
BPMS BURST PRESSURE 5000 PSIG
BPMS MAX OVERPRESSURE 235 PSIG
CUFF MAX OPERATING PRESSURE 250 MM HG
CUFF RELIEF PRESSURE 260 MM HG
GAS TEMPERATURE 40 TO 120F

619011 RADIATION SOURCE STORAGE LS979 11112
619011
619011 STORAGE FACILITY FOR RADIOACTIVE MATERIAL.
619011
619011 VOLUME APPROXIMATELY 1 CU FT
619011 RADIATION PROTECTION LIMIT ... UP TO 500 MICROCURIES
619011

618016 RECEIVER - BIOTELEMETRY LS980 11111
618016
618016 RECEPTION OF BIOTELEMETRY SIGNALS FROM DEEP BODY TEMPERATURE AND
618016 ANIMAL ACTIVITY AS WELL AS FOR ELECTROMAGNETIC FIELD MONITORING.
618016
618016 FREQUENCY RANGE TO 5 MHZ
618016 SENSITIVITY 0.5 MV FOR 20 DB QUIETING
618016 TUNING VARIABLE AND CRYSTAL CONTROLLED
618016

HEWLETT PACKARD #78101A 1800 112
CONFIGURATION BENCH EQUIP
FREQUENCY RESPONSE 0.1 TO 1000 HZ
SENSITIVITY 0.5 MICRO-V
ANTENNA 7 INCH MONOPOLE
INPUT IMPEDANCE 50 OHMS NOMINAL

039644 VISION TESTER LS981 11235
039644
039644 OPTICAL DEVICE (WITH HOOD AND RESPONSE KEYBOARD) THAT MEASURES
039644 A VARIETY OF VISUAL FUNCTIONS.
039644

338000 WASTE STORAGE SYSTEM LS982 11111
338000
338000 PROVIDE FOR HANDLING, STORAGE AND DISPOSAL OF SOLID AND LIQUID
338000 WASTE MATTER FROM THE EXPERIMENTS.
338000
338000 CAPACITY 5 CU FT
338000

470500 STERILIZER - TOOL LS983 11111
470500
470500 STERILIZE MISCELLANEOUS SMALL METAL HAND TOOLS SUCH AS SCALPES,
470500 BY MEANS OF ELECTRICAL INDUCTION HEATING.
470500
470500 CAPACITY 0.1 CU FT
470500

VWR SCIENTIFIC #5861C-CC7 \$255
OPERATING TEMPERATURE.....200C (MAX)
SHELF SIZE.....18X10X1/2 (3)
CHAMBER SIZE.....18X8X10
POWER.....115 VAC, 50/60 HZ, 7.5 A
900 W

034930 AUDIO STEREO HEADSET LS985
034930
034930 EARPHONES FOR VARIOUS HSI HEARING TESTS
034930

173000 FLOW METER, WATER MANIFOLD LS989 11235
173000
173000 IN LINE MEASUREMENT OF WATER FLOW, GENERALLY A LOW RATE
173000 ASSOCIATED WITH WATER CONSUMPTION BY ORGANISMS
173000

172500 FLOW METER, GAS LS992 11111
172500
172500 MEASURE AIR FLOW
172500

401000 IMPEDANCE PNEUMOGRAPH LS995 11111
401000
401000 MEASURE BREATHING CYCLE CHARACTERISTICS
401000

TECHNOLOGY/VERSATRONICS INC #MFG-3 \$970 111
OUTPUT(SINGLE ENDED) 0-5VDC LINEAR IN MASS FLOW
OUTPUT IMPEDANCE <10 OHM
ACCURACY(ABSOLUTE) +/- 2% OF READING 20 TC 100% FS
LINEARITY +/- 0.4% FULL SCALE (FS)
REPEATABILITY +/- 0.1% OF READING
CALIBRATION PRESSURE 1 ATMOSPHERE

ORIGINAL PANELS
OF POOR QUALITY



STD CALIBRATION TEMP RANGE 35 TO 155F
MEDIUM AIR
TIME CONSTANT LESS THAN 0.07 SEC
PRESSURE DROP LESS THAN 1 IN-H2O
READOUT 18 TAUT BAND MIRROR SCALE
0-100% FLOW RATE
AVAILABLE RANGE 0.12 TO 0.6 CU FT PER MIN
SIGNAL CONDITIONER 0 TO 40F
TRANSDUCER -45 TO 150F

907900	KIT, CLEAN UP	LS999	11111
907900	GENERAL PURPOSE CLEAN UP - SPONGES AND WIPES		
907900	KIT, MEDICAL SURGICAL	LS1000	11111
907900	FOR VARIOUS MINOR SURGICAL PROCEDURES		
907900	KIT, PHYSIOLOGY	LS1001	11111
907900	SPONGES, SPONGE SQUEEZER, VIALS, CALORIC STIMULATOR FOR EAR		
907900	CANAL, SYRINGES, THERMOMETERS, ETC.		
907900	KIT, TOOL - INSECT MANIPULATION	LS1002	11111
907900	TOOL KIT FOR COUNTING, SORTING, EXAMINATION, ETC.		
256000	MICROPHONE	LS1007	11111
256000	MICROPHONES USED FOR BIOLOGY, BIOMEDICINE AND MANNED SYSTEMS		
256000	INTEGRATION		
006170	MICROPHONE AMPLIFIER	LS1008	11111
006170	AMPLIFY SIGNALS FROM AMPLIFIER		
006170	MONITOR, VIDEO	LS1009	11111
484000	MONITOR TV PICTURES OF ANIMAL AND OTHER LAB ACTIVITIES		
484000	COMU ELECTRONICS INC. #CQF	8775	111
	BANDWIDTH 30 MHZ		
	HORIZONTAL LINES 1225		
	FIELDS 60/SEC.		
	CCNRAC CORPORATION #RQA 14/R	81620	111
	VIDEO AMPLIFIER		
	COMPOSITE 0.3 - 3.0 V P-P		
	NONCOMPOSITE 0.3 V P-P FOR 100 PCT CONTRAST		
	3.0 V P-P MAXIMUM BEFORE PREAMP OVERLOAD		
	VIDEO INPUT IMPEDANCE 50 K OHM MIN SHUNTED BY MAX 10		
	PFD SIGNAL OF 20 V P-P; HUM SUPPRESSION 40 DB 60-1000 HZ		
	VIDEO FREQUENCY RESPONSE		
	AT 50 PCT OF MODULATION FLAT +/- 1 DB TO 30 PHZ, -6 DB		
	AT 30 PCT OF DEPTH MOD FLAT +/- 1 DB TO 30 PHZ, -2 DB		
	PULSE RESPONSE LESS THAN 15 MS RISE/FALL TIME		
	AT 50 PCT DEPTH OF MODULATION		
	EXTERNAL SYNC COMPOSITE OR SEPARATE H/V SYNC;		
	0.5 V TO 10 V P-P, GREATER THAN 5K OHM INPUT IMPEDANCE		
	LINEARITY +/- 1 PCT PICTURE HEIGHT		
	RETRACE TIME		
	HORIZONTAL 5 MICROSEC MAX		
	VERTICAL 800 MICROSEC MAX		
	SCANNING FREQUENCY		
	HORIZONTAL 15 - 40 KHZ		
	VERTICAL 15 - 60 FIELDS/SECOND		
	COMU INC. ELECTRONICS DIV #960C	85260	
	PICTURE SIZE 7-3/16 X 5-3/8 INCHES		
	HORIZONTAL SCAN 945 LINES, 50 OR 60 FIELDS/SEC		
	VIDEO BANDWIDTH 875 LINES, 60 FIELDS/SEC		
	POWER 20 MHZ +/- 2 DB		
	25-28.5 VDC, 150 WATTS		
419500	PUMP GAS, CIRCULATING	LS1011	11111
419500	PUMP FOR SEALED PLANT GROWTH CHAMBER		
498200	TIMER, EVENT	LS1020	11111
498200	GENERAL PURPOSE ELAPSED TIME DEVICE USED FOR VARIOUS PSYCHOMOTOR		
498200	TESTS		
	ROWDE AND SCHWARZ #CAD (100,6597,91)	85400	111
	DISPLAY 6 DIGITS: HR, MIN, SEC		
	OUTPUT 1-CUT-OFF-N AND/OR BCD CCDE		
	PROGRAMMER TIME PROGRAMS/COUNT		
	DOWN) AND TIME SIGNALS		
	CATATRON #3350-506	81170	111
	INPUT FREQUENCY 60HZ		
	OUTPUTS VISUAL, BCD		
	DISPLAY HRS, MIN, SEC		



TENNELEC #TC545 8550 111
COUNT RATE..... 20PHZ
TIMEBASE..... 0.1 OR 0.01 SEC
ACCURACY..... SAME AS LINE FREQUENCY
CONFIGURATION..... NIM COMPATIBLE
043500 CARDIOGRAPH - IMPEDANCE LS100 11111
043500
043500 BLOOD FLOW RATE FROM HEART.
043500
043500 BLOOD FLOW RATE 2 TO 25 LITERS/MIN +/-5 PCT FS
043500
MARTIN CO #M093 (SKYLAB) 511
MEASURES ELECTROCARDIOGRAPHIC POTENTIALS OF EACH ASTRONAUT
143500 ECG/VCG LS101 11131
143500
143500 MONITOR HEART CONDITION. ELECTRO/VECTORCARDIOGRAPH.
143500
143500 SENSITIVITY 0 TO 3 MV +/-1 PCT FS
143500 HEART RATE 40 TO 180 BEATS/MIN
143500
TEKTRONIX INC #410 8975 111
BANDWIDTH 0.1 HZ TO 100 HZ (+/- 15%)
SENSITIVITY 10 MM/50 MICROVOLT (+/- 5%)
AUX MODE 2 MM/MV (+/- 5%)
DISPLAY 5 INCH CRT
144500 EEG - ELECTROENCEPHALOGRAPH LS102 11131
144500
144500 MONITOR ELECTRICAL IMPULSES GENERATED BY THE BRAIN.
144500
144500 SENSITIVITY 10 TO 200 MICROVOLTS +/-1 PCTFS
144500 FREQUENCY RESPONSE 0.2 TO 100 HZ
144500 CHANNELS 6
144500
BECKMAN INSTRUMENTS #9P52A(COUPLER ONLY) 895 111
COUPLER USED WITH 425500
FREQUENCY 10 TO 200 HZ
PULSE DURATION 0.6 MS TO 20 MS
SKIN ELECTRODES REQUIRED
TEKTRONIX INC #410-3 81025 111
BANDWIDTH 0.1 HZ TO 100 HZ (+/- 15%)
SENSITIVITY 10 MM/50 MICROVOLT (+/- 5%)
AUX MODE 2 MM/MV (+/- 5%)
DISPLAY 5 INCH CRT
BECKMAN INSTRUMENTS #ACCUTRACE 85890 112
SENSITIVITY 0.5 MICROVOLTS/MM (MAX)
SENSITIVITY MASTER CONTROL 2,3,5,7,5,10,15,20,30,50,75MMV
SENSITIVITY MULTIPLIER CONTROLS 120,4,2,1,0.5,0.25
FREQUENCY RESPONSE DC TO 70 HZ
CHANNELS 8 PLUS TWO EVENT CHANNELS
WEIGHT 150 LBS
DIMENSIONS 15X23X21 INCHES
POWER 105-130V 60 HZ 50W
BECKMAN INSTRUMENTS #TYPE T PORTABLE 86670 112
SENSITIVITY 1 TO 100 MICROVOLTS/MM
FREQUENCY RESPONSE DC TO 150 HZ
CHANNELS 8
WEIGHT 16 LBS
POWER 115V 60 HZ 50W AND 6 D CELLS
DIMENSIONS 17X19X10.5 INCHES/EACH OF 2 PCS
142000 DYNAMOMETER LS103 11133
142000
142000 MEASUREMENT OF FORCE OR POWER.
142000
142000 ARM FLEXION 25 TO 75 LB (10 TO 35 KG)
142000 BACK EXTENSION 100 TO 200 LB (40 TO 90 KG)
142000
145500 EMPG - ELECTROMYOGRAPH LS104 11131
145500
145500 MONITOR THE MOVEMENT OF THE BODY MUSCULAR SYSTEM.
145500
145500 SENSITIVITY 0.01 TO 5 MV +/-1 PCT FS
145500 BANDWIDTHS 0.5 TO 200 HZ, 0.5 TO 1000 HZ
145500 CHANNELS 1
145500
143500 PHONOCARDIOGRAM LS106 11111
143500
143500 AURAL MONITOR OF HEART BEAT.
143500
143500 SENSITIVITY 0.1 TO 1000 HZ, +/-5 PCT FS
143500
BECKMAN INSTRUMENTS #9R02 COUPLER ONLY 835 111
COUPLER USED WITH 425500
FREQUENCY 5 TO 2000 HZ
ABOVE 200 HZ REQUIRES SCOPE OR TAPE READCUT
TRANSDUCERS REQUIRED
052570 METABOLIC ANALYZER LS1003 11131
052570
052570 MEASURE O2 CONSUMPTION, CO2 PRODUCTION, MAXIMAL O2 CONSUMPTION
052570 AND ALVEOLAR PO2 AND PCO2.
052570

ORIGINAL PARTS
OF POOR QUALITY



082970 O2 VOLUME 0 TO 40 L/MIN, +/-3 PCT FS
082970 CO2 VOLUME 0 TO 40 L/MIN, +/-3 PCT FS
082970 ALVEOLAR O2 20 TO 120 MM HG +/-2 PCT FS
082970 ALVEOLAR CO2 20 TO 70 MM HG +/-2 PCT FS
082970

MARTIN CO 0M171 (SKYLAB) 511
MEASURE METABOLIC ACTIVITY
PARAMETERS MEASURED..... OXYGEN CONSUMPTION
..... CARBON DIOXIDE
..... RESPIRATORY EXCHANGE RATIO
..... MINUTE VOLUME
..... VITAL CAPACITY
..... BODY TEMPERATURE
..... ERGOMETER WORK RATE, RPP AND
..... TOTAL WORK
COMPONENTS OF ASSEMBLY MASS SPECTROMETER
..... GAS ANALYZER
..... ANALOG COMPUTER
..... INSPIRATION AND EXPIRATION
..... SPIROMETERS
..... CO2 STABILIZER
..... CALIBRATION GAS SUPPLY

PERKIN ELMER CORP 0NAS 9-2244 845000 511
TYPE SOLID STATE
CO2 RANGE 130 TO 4000 PA (1-30 MM HG)
OUTPUT VOLTAGE 0 - 5 VDC
OUTPUT IMPEDANCE < 500 OHMS

RECKMAN INSTRUMENTS 0554103 915000 111
SENSITIVITY O2..... 100% FS
CO2..... 10% FS
ACCURACY O2 +/-2 FS
CO2..... +/-1 % FS
RESPONSE (MSEC) O2 300 AT 500ML/MIN
CO2 900 AT 100ML/MIN
METER READOUT DIGITAL
RESOLUTION (B GAS)..... 0.01% (CO2), 1% (O2)
MEASURE ALVEOLAR O2..... YES
MEASURE ALVEOLAR CO2..... YES
CHART SPEED 1MM/510MM/S
WRITING METHOD..... HEATED STYLLS
WEIGHT 270 LBS (122.5 KG)

043500 PULMONARY FLOWMETER LS990 11131
043500
043500 BLOOD FLOW RATES AND VOLUME.
043500
043500 FLOW RATE 0 TO 15 L/MIN, +/-3 PCT FS
043500 VOLUME 0 TO 7 L, +/-3 PCT FS
043500
043500 TRANSCUTANEOUS DOPPLER FLOWMETER LS991 11112
043500
043500 MEASURE PULSE WAVE CONTOUR.
043500
043500 ULTRASONIC SIGNAL 1 AND 5 MHZ
043500
048000 PULSE WAVE LS112 11112
048000
048000 MEASURE BLOOD VELOCITY
048000
048000 VELOCITY 3 TO 15 M/S, +/-0.1 M/S
048000

BECKMAN INSTRUMENTS INC 0 84000 111
150400 ERGOMETER - BICYCLE LS986 11131
150400
150400 WORKLOAD MEASURING DEVICE.
150400
150400 LOAD RANGE 0 TO 300 WATTS, +/-2 PCT FS
150400 ROTATIONAL SPEED 40 TO 90 RPM
150400
107110 LITTER CHAIR - ROTATING LS114 11131
107110
107110 GRAVITY INDUCING DEVICE AND HORIZONTAL WORK TABLE.
107110
107110 ROTATION SPEED 0 TO 60 RPM
107110
612015 LBNP - LOWER BODY NEGATIVE PRESSURE LS115 11131
612015
612015 MEASURE INTERNAL BODY PRESSURE DIFFERENTIAL.
612015
612015 PRESSURE DIFFERENTIAL 50 MM HG BELOW CABIN AMBIENT
612015

MARTIN CO 0M092-LBNPD (SKYLAB) 511
VOLUME (INTERNAL) 10 CU.FT
OPERATING TEMPERATURE RANGE 0 TO 70F
MAX NEGATIVE PRES BELOW AMBIENT 50 MM HG
LEAKAGE 2.5 CU.FT./MIN MAX

190310 EXERCISER/ERGOMETER - PRIMATE LS175 11133
190310
190310 PRIMATE EXERCISER AND WORKLOAD MEASURING DEVICE.
190310
003009 CALORIMETRY MODULE - PRIMATE LS964 11135
003009
003009 PORTABLE AND COLAPSABLE UNIT WHICH HOUSES, FEEDS, AND ANALYZES
003009 THE PRIMATE SPECIMEN WITHIN.



603009
603009 PD2 160 +/-10 MM HG
603009 PN2 580 +/-20 MM HG
603009 PCO2 6 MM HG OR LESS
603009 PH2O 12 +/-3 MM HG
603009 TEMPERATURE 77 +/-4 F 125 +/-2 C
603009
023500 MASS MEASUREMENT DEVICE - BODY LS987 11131
023500
023500 TOTAL BODY MASS DETERMINATION BY OSCILLATING FREQUENCY.
023500
023500 MASS 125 TO 230 LB +/-0.1 PCT FS
023500 (50 TO 100 KG)
023500
SOUTHWEST RESEARCH INSTITUTE #M172 (SKYLAB) 511
ELEMENTS OF ASSEMBLY SEAT, SPRING LOADING SYSTEM
AN ELECTRONIC SUBSYSTEM MODULE
229600 LIGHTING SYSTEM - PLANT LS1018 21114
229600
229600 PROVISIONS FOR LIGHT PLANT GROWING AREAS.
229600
229600 TYPE LIGHT DAYLIGHT-FLUORESCENT
229600 ILLUMINATION 102 FT-C (1100 LUM/M2)
229600
304000 PLANT GROWTH AND SUPPORT CONTAINERS LS1019 121111
304000
304000 VARIOUS SIZE CONTAINER TO ENVIRONMENTALLY HOUSE SEEDLINGS.
304000
304000 CAPACITY UP TO 20-CM PLANTS
304000
217500 CELLS/TISSUE HOLDING UNIT LS1016 121114
217500
217500 COMBINATION HOLDING UNIT/INCUBATOR CONTAINING ELECTRONIC PLUG-IN
217500 CAPABILITY FOR MEASUREMENT DEVICES
217500
217500 TEMPERATURE RANGE 41 TO 140 F (5 TO 60 C)
217500 PCO2 160 +/-5 MM
217500 PN2 600 +/-20 MM
217500 PCO2 3 MM OR LESS
217500 TOTAL PRESSURE 760 +/-20 MM
217500 RELATIVE HUMIDITY 60 TO 100 PCT
217500
229600 CELLS/TISSUE HOLDING UNIT LIGHTING LS1017 21114
229600
229600 PROVISION FOR PROVIDING COOL DIFFUSED LIGHT
229600
229600 ILLUMINATION 90 +/-10 FT-C
229600 TEMPERATURE AMBIENT
229600

ORIGINAL PAGE
OF POOR QUALITY




SPACE TECHNOLOGY

603013	CONDENSER	ST102	11131
603013			
603013	DEVICE TO CONDENSE STEAM IN ZERO GRAVITY		
603013			
410000	CONDENSATE TANKS	ST102	21113
410000			
410000	TANKS FOR STORAGE OF CONDENSATE		
410000			
410000	PRESSURE..... 20PSI (1.38X10 ⁵ N/m ²)		
410000			
104000	SUPPORT ELECTRONICS	ST006	21114
104000			
104000	ELECTRONICS SUPPORT COMMUNICATIONS EXPERIMENTS INCLUDING		
104000	MULTIPLEXERS WITH CONVERTERS AND OSCILLATORS		
104000			
358500	OSCILLOSCOPE	ST007	21111
358500			
358500	MONITOR, MEASURE AND MAINTAIN ELECTRONIC EQUIPMENT OPERATING IN		
358500	THE L-BAND RANGE (390 TO 1550 MHZ).		
358500			
358500	BANDWIDTH 390 TO 1550 MHZ		
358500	CHANNELS 2		
358500	DEFLECTION FACTOR 10 MV/DIV TO 10 V/DIV		
358500	PHOTOGRAPHIC CAPABILITY YES		
358500			
	TERRONIX INC	PR564B/3A6 PLUG-IN/	92340 111
		304 PLUG-IN	
	BANDWIDTH DC TO 10 MHZ		
	CHANNELS 2		
	DEFLECTION FACTOR 10 MV/DIV TO 10 V/DIV		
	TIME BASE 50 NS/DIV TO 5 S/DIV		
	MAX VIEWING TIME 1 HOUR		
	MAX STORED WRITING SPEED 500 DIV/MS		
	TYPE OF STORAGE SPLIT SCREEN BISTABLE		
	TEKTRONIX	80M64	81095 111
	BANDWIDTH DC TO 10 MHZ		
	RISE TIME 35 NS		
	SPEED WRITING 25 CM/MS (NOMINAL)		
	STORAGE VIDEO TIME 1 HR		
	ERASE TIME 0.25 S		
	DEFLECTION FACTOR 10 MV/CM TO 50 V/CM IN 12 STEPS		
	TIME BASE 100 NS/CM TO 2 S/CM IN 23 STEPS		
	RESPONSE TIME LESS THAN 1 NS		
	HEWLETT-PACKARD	0132A	81500 111
	BANDWIDTH DC TO 500KHZ		
	SENSITIVITY 100 MICROVOLTS/CM		
	INDEPENDENT BEAMS 2		
	WEIGHT 43 POUNDS		
	TEKTRONIX INC	PR465	81800 111
	BANDWIDTH DC TO 100 MHZ		
	CHANNELS 2		
	DEFLECTION FACTOR 5 MV/DIV TO 5 V/DIV		
	TIME BASE 0.01 MICROSEC/DIV TO 0.5 S/DIV		
	TEKTRONIX	PR5	84200 111
	RANGE CALIBRATED 5MV/DIV TO 5V/DIV IN 10 STEPS		
	UNCALIBRATED CONT VAR TO 12.5V/DIV		
	RESPONSE TIME LESS THAN 1 NS		
	CHANNELS 1, ALTERNATE, CHOPPED, ADDED, X-Y		
 CHANNEL 2 (UP OR INVERTED)		
	DIMENSIONS 16.7 X 52.4 X 30.5 CM		
	WEIGHT 9.5 KG		
	POWER 115/230V 40-440HZ 60W		
	TEKTRONIX INC	PR561B/3A6 PLUG-IN/	81840 111
		304 PLUG-IN	
	FREQUENCY BANDWIDTH DC TO 10 MHZ		
	CHANNELS 2		
	DEFLECTION FACTOR 10 MV/DIV TO 10 V/DIV		
	TIME BASE 50 NS/DIV TO 5 S/DIV		
	RISE TIME 35 NS		
424000	TYPE RECORDER, DIGITAL	ST008	21111
424000			
424000	RECORD AND PLAYBACK OF DIGITAL DATA IN SUPPORT OF MICROWAVE		
424000	INTERFEROMETER AND TRACKING INVESTIGATIONS.		
424000			
424000	SIGNAL TYPE 0 TO 5 VDC		
424000	BANDWIDTH 100 TO 250 KHZ		
424000	CHANNELS 4 (MINIMUM)		
424000	DATA REQUIREMENT 2 HRS (MINIMUM)		
424000			
	AMPX CORPORATION	PAR-200	121
	BANDWIDTH 100 HZ TO 3125 HZ THROUGH		
	TAPE WIDTH 0.5 INCH		
	TAPE SPEED 1.875 TO 60 IPS		
	RECORDING TIME 8 MINUTES TO 4 HRS AND 16 MINS		
	FORMAT DIGITAL		
	TRACKS 8 DIGITAL, 7 ANALOG		
	AMPX CORP	PAR 1700	828,000 112



TRACKS	14		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
RECORDING MODE	DIRECT		
PACKING DENSITY	20 MB/1/T		
SIGNAL/NOISE	20 DB		
DATA CAPACITY	6.2E+10		
BANDWIDTH	2 MB/SEC/T		
SANGAM ELECTRIC	05ABER	111	426,500
TRACKS	14		
TAPE WIDTH	1 INCH		
TAPE SPEEDS	8 SELECTABLE SPEEDS FREQ 15/16 TO 120 IPS		
FREQUENCY RESPONSE	400 HZ TO 2.0 MHZ		
RECORDING RATE	600 KOPS AT 120 IPS SERIAL MODE		
WEIGHT	100 POUNDS		
DIGI-DATA CORPORATION	#160C	92690	111
TAPE SPEED	25, 18.75, 12.5 IPS		
TRACKS	7 OR 9 TRACK		
DATA DENSITY	1600 CPI PHASE ENCODED		
TAPE	200, 556, 800 CPI NRZI		
TAPE	0.5 INCH, 1.5 MIL, 1200 FEET		
DIGI-DATA CORPORATION	#1700/POP-11/7-9	85250	111
TRACK NRZI			
TAPE SPEED	45, 37.5, 25, 18.75, 12.5 IPS		
TAPE	0.5 INCH, 1.5 MIL, 100/ANSI		
TRACKS	7 OR 9		
DATA DENSITY	PHASE ENCODED		
COMPATIBLE			
HEWLETT PACKARD	#757C9/C	5 4600	111
TAPE FORMAT	800, 556, OR 200 CPI NRZI AND 1600 CPI PHASE-ENCODED		
CHANNELS	7 OR 9		
TAPE SPEED	10 TO 45 IPS		
TAPE	0.5 INCH, 1.5 MILS, 100/ANSI		
COMPATIBLE			
BORG WARNER	0PERT		121
TRACKS	30		
TAPE SPEED	UP TO 1000 IPS		
TAPE WIDTH	1/2 INCH		
TAPE LENGTH	2400 FEET		
PACKING DENSITY	15 MB/1/T		
SIGNAL/NOISE	24 DB		
DATA CAPACITY	1.3E+9		
BANDWIDTH	6-15 MB/SEC/T		
LEACH	#RTR TC00		311
TRACKS	12		
TAPE SPEED	120 IPS		
TAPE WIDTH	1 INCH		
TAPE LENGTH	9200 FEET		
PACKING DENSITY	16.7 MB/1/T		
SIGNAL/NOISE	22 DB		
DATA CAPACITY	2.2E+10		
BANDWIDTH	2 MB/SEC/T		
HEWLETT PACKARD	#35550	810200	111
BANDWIDTH	300 KHZ		
CHANNELS	7		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3555C	814700	111
BANDWIDTH	300 KHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3990A-011	823800	111
BANDWIDTH	500 HZ TO 2 MHZ		
CHANNELS	14		
RECORDING FORMAT	DIRECT OR FM		
HEWLETT PACKARD	#3560A/13065A/13063A 84700		111
CONFIGURATION	RACK MOUNTED		
TAPE SPEED	15, 3 AND 3/4, 15/16 IPS		
CHANNELS	4		
RECORDING FORMAT	FM		
PASSBAND	5 KHZ		
S/N RATIO	40 DB		
TRACKS	30		
HEWLETT	#5600	80730	
PORTABLE TAPE RECORDER			
CHANNELS	7		
SELECTABLE TAPE SPEED RANGE	15/16 TO 60 IPS		
MAXIMUM BANDWIDTH (DIRECT)	3000 HZ		
PACKING DENSITY	UP TO 6000 BI		
WEIGHT	70 LBS (32 KG)		
INPUT VOLTAGE	200VDC		
HEWLETT	#994	817420	

ORIGINAL PARTS
OF POOR QUALITY




SD 74-BA-0047-2



HONEYWELL		#1P56	1
RESOLUTION	200 ELEMENTS/INCH (1000/SWEEP)		
RECORDING RATE	RE+GE BITS/SEC		
CRT	PIPER OPTICS		
SPOT SIZE	0.005 INCH		
SWEEP SPEED	20,000 LINES/SEC		
HONEYWELL		#1P58	1
SIZE-RACK MOUNTED	8.75H X 19.1W X 21D		
WEIGHT (W/O ACCESSORIES)	65 LBS.		
ACCESSORIES	6 LRS.		
CHANNELS	18		
RECORDING STYLES	INERTIAL-NO MOVING ELEMENT		
RESPONSE	DC TO 5000 HZ		
INPUT SENSITIVITY	100 MICROVOLT - 300 VOLTS		
ACCURACY	0.1 %		
PAPER SPEED	0.1 TO 120 IPS		
POWER	120/240 V, 50-400 HZ		
HONEYWELL		#1S12	111
PAPER WIDTH	12IN		
CHANNELS	24		
FREQUENCY RESPONSE	DC TO 25KHZ		
PAPER SPEED	200IPS		
SPEED VARIATION	15 FORWARD, 12 REVERSE SPEEDS		
MOUNT	RACK		
BECKMAN INSTRUMENTS		#R	111
BANDWIDTH	DC TO 130KHZ		
CHANNELS	5000 HZ WITH PR-SOOTAPE REC		
REC FORMAT	1 TO 4 UP TO 24 SPEC ORDER		
CHART SPEED	0.1 TO 25 CM/SEC OTHERS AVAIL		
EVENT MARKERS	2 (OPTION)		
DIMENSIONS	65.5X21.5X22.5 INCHES		
POWER	120V 50/60 HZ 250W		
181000 SIGNAL GENERATOR		ST016	001111
181000			
181000 ACCURATELY REPRODUCE IN PULSE FORM THE OUTPUT OF A MICROWAVE			
181000 RADIOMETER MONITORING SEA SURFACE STATE AND TEMPERATURE.			
181000			
REPETITION RATE	ON		
OUTPUT VOLTAGE	+/- 5 VDC		
181000			
INTERSTATE ELECTRONICS CORP		#P25	111
REPETITION RATE	1 HZ TO 50 MHZ IN 4 RANGES		
POLARITY	1 POSITIVE AND 1 NEGATIVE CHAN		
WIDTH	<10 NSEC TO 1 SEC		
DELAY	<10 NSEC TO 1 SEC		
AMPLITUDE	10, 5, 1 V		
NAVETER		#162	111
BANDWIDTH	0.00003 HZ TO 30 MHZ		
OUTPUT	80 DB ATTENUATION/10 DB STEPS		
CC OFFSET	10 V P-P INTO 50 OHMS		
STABILITY	+/- 5 VOLTS INTO 50 OHMS		
	+/- 0.25% PER 24 HOURS		
TEKTRONIX INC		#FG 501	111
FREQUENCY RANGE	.001 HZ TO 1 MHZ		
AMPLITUDE OUTPUT	7.5 V P-P INTO 50 OHMS		
WAVEFORMS	SINE, PULSE, SQUARE, TRIANGLE, RAMP		
179500 SPECTRUM ANALYZER		ST017	11111
179500			
179500 SUPPORT SETUP, CHECKOUT AND SIGNAL OUTPUT ANALYSIS OF RF RADIO-			
179500 METER UNIT AND SWEEP RECEIVER (CRT DISPLAY OF SWEEP RECEIVER AMP-			
179500 LITUDE VERSUS FREQUENCY).			
179500			
BANDWIDTH	30 TO 10000 PHZ		
RESOLUTION	2 %		
SENSITIVITY	1C MV/DIVISION		
179500			
TEKTRONIX INC		#P491	112
FREQUENCY RANGE	10 MHZ TO 40 GHZ		
DISPERSION, FREQ SPAN	1 KHZ/DIV TO 10 MHZ/DIV		
RESOLUTION BANDWIDTH	1 KHZ TO 100 KHZ		
INCIDENTAL FM	<300 HZ AT FUNDAMENTAL		
TEKTRONIX INC		#1401A-4	111
FREQUENCY RANGE	1 TO 500 MHZ		
INTERMODULATION DISTORTION	< 40 DB FULL SCREEN		
SPECTRA-PHYSICS		#ATC	111
FOR USE IN LASER ALIGNMENT			
SPECTRAL RANGE	2 TO 8 GHZ		
INSTRUMENTAL BANDWIDTH	6.7 TO 27 MHZ		
055050 CAMERA - STILL (VISTAR)		ST020	31111
055050			
055050 STILL PHOTOGRAPHIC COVERAGE OF THE EARTH TARGET AREA IN THE			
055050 VISIBLE REGION (REQUIRES SIMULTANEOUS PHOTOGRAPH IN THE INFRARED			
055050			
VIEW ANGLE	14.5 DEGREES		
RESOLUTION	100 LINES/MM		
FILM FORMAT	5 INCHES		
SHUTTER	8 TO 1/4000 SEC		
055050			

ORIGINAL PARTS
OF POOR QUALITY



SD 74-SA-0047-2



WORD LENGTH 16 BIT
CYCLE TIME 1200 NANOSEC.
CORE MEMORY 32K MAX.

VARIAN DATA MACHINES #520/1 111
MEMORY CYCLE TIME 1.5 MICROSEC.
MEMORY EXPANDABLE 4096 BYTES (8 BITS) TO 32,768 BYTES
REGISTERS 12
OPERAND PRECISION UP TO 32 BITS

VARIAN DATA MACHINES #P-20/1 111
BASIC COMMANDS OVER 100
ADDRESSING MODES 6
MAX. WORDS 32,768
WORD LENGTH 16 OR 18 BIT
REGISTERS 6
NOTE: HAS FULL RANGE INTERFACE HARDWARE

CLARY DATA/COMP SYSTEMS INC #404 111
SIMULTANEOUS TERMINALS 16
INPUT TELETYPEWRITER KEYBOARD
ACCUMULATOR 64 BIT
REGISTERS 16 BIT INDEX 12 EACH
WORD LENGTH 16, 32, 48, 64 BITS
MEMORY CAPACITY 1024 16 BIT WORDS, 4096 16 BIT WORDS, OR ADDITIONAL 4096 16 BIT WORDS, TOTAL 65536 WORDS

GENERAL AUTOMATION INC. #IP-30 111
FUNCTION SUPERVISE SMALLER COMPUTERS
MEMORY 8K CORE
CYCLE TIME 1.2 MICROSEC.

GENERAL AUTOMATION INC. #SPE-14 111
MEMORY 16K
WORD LENGTH 16 BITS
READ/WRITE CYCLE MEMORY TIME 400 TO 1440 NANOSEC.
READ ONLY MEMORY 400 TO 720 NANOSEC.
INPUT/OUTPUT TIME 1.8 TO 2.4 MICROSEC.
DMA TRANSFER RATE 0.654E+06 TO 2.5E+06

GENERAL AUTOMATION INC. #SPE-12 111
MEMORY 4K TO 16K
WORD LENGTH 8 BITS
CYCLE TIME 2.36 MICROSEC
STORED PROGRAM EXECUTION RATE 0.23E+06/SEC
I/O TRANSFER RATE 0.4E+06/SEC
REGISTERS 6 12 BIT REGISTERS
ACCUMULATORS 4 12 BIT

DIGITAL EQUIPMENT CORP #PDP-8/E \$22000 111
MEMORY 4096 CORE, EXPANDABLE TO 32,768 WORDS
WORD 12 BIT
CYCLE TIME 1.5 MICROSEC.
ADDRESS REFERENCE 8, 16, 24, OR 32 BIT LEVEL
NOTE: HAS FULL RANGE INTERFACE HARDWARE

484000 MONITOR - TV ST026 31111
484000
484000 DISPLAY HIGH RESOLUTION TELEVISION PICTURES IN REAL-TIME OR FROM
484000 VIDEO TAPE. USED ON ST01S, ST02S, ST03S, ST04S AND ST05S.
484000
484000 RESOLUTION 425 LINES
484000 FRAME RATE 30 FPS
484000 SCREEN SIZE 10 IN (25.4 CM)
484000

COMU ELECTRONICS INC. #CCF \$775 111
BANDWIDTH 50 MHz
HORIZONTAL LINES 1225
FIELDS 60/SEC.

CENRAC CORPORATION #PQA 14/R \$1420 111
VIDEO AMPLIFIER
COMPOSITE 0.3 - 3.0 V P-P
NONCOMPOSITE 0.3 V P-P FOR 100 PCT CONTRAST
3.0 V P-P MAXIMUM AT CRE. PREAMP OVERLOAD
VIDEO INPUT IMPEDANCE 50 K OHM MIN SHUNTED BY MAX 10 PFD SIGNAL OF 20 V P-P; HUM SUPPRESSION 40 DB 60-1300 HZ
VIDEO FREQUENCY RESPONSE
AT 50 PCT OF MODULATION FLAT +/- 1 DB TO 30 MHz, -6 DB AT 40 MHz
AT 30 PCT OF DEPTH MOD FLAT +/- 1 DB TO 30 MHz, -2 DB AT 40 MHz
PULSE RESPONSE LESS THAN 15 NS RISE/FALL TIME AT 50 PCT DEPTH OF MODULATION
EXTERNAL SYNC COMPOSITE OR SEPARATE H/V SYNC; 0.5 V TO 10 V P-P, GREATER THAN 5K OHM INPUT IMPEDANCE
LINEARITY +/- 1 PCT PICTURE HEIGHT
RETRACE TIME
HORIZONTAL 5 MICROSEC MAX
VERTICAL 400 MICROSEC MAX

ORIGINAL PARTS
OF POOR QUALITY

SCANNING FREQUENCY
HORIZONTAL 15 - 40 KHZ
VERTICAL 15 - 60 FIELDS/SECOND

COMU INC. ELECTRONICS DIV 89600 85260
PICTURE SIZE 7-3/16" X 5-3/8" INCHES
HORIZONTAL SCAN 945 LINES, 50 OR 60 FIELDS/SEC
..... 873 LINES, 60 FIELDS/SEC
VIDEO BANDWIDTH 20 MHZ +/- 2 DB
POWER 25-28.5 VDC, 150 WATTS

103520 DATA PROCESSING EQUIPMENT ST035 11111
103520
103520 HIGH SPEED DIGITAL DATA CONVERSION TO ANALOG FORM COMPATIBLE FOR
103520 REAL TIME DISPLAY. USED ON ST01S, ST02S AND ST03S.
103520
103520 CHANNELS 32
103520 DIGITAL-DATA 10E+06 RPS
103520

TENNRELC 87C500 96000 111
CONVERSION TIME 3MICRO SEC
CAPACITY 2048 CHANNELS
CONFIGURATION COMPATIBLE WITH AEM
TENNRELC 87C 520 6750 111
CHANNELS 8
CHANNEL DWELL TIME 300NSEC
CONFIGURATION NIM COMPATIBLE

410000 STORAGE - GAS SOURCE ST045 01111
410000
410000 GAS SOURCE STORAGE BOTTLES CAPABLE OF STORING VARIOUS GASES AT
410000 HIGH PRESSURES. SIMILAR TO STANDARD A-BOTTLES. USED IN ST01S,
410000 ST02S AND ST03S.
410000
410000 VOLUME - PHYSICAL 5 CU FT EA (0.15 CL MI)
410000 PRESSURE MAXIMUM 1000 PSI
410000

003400 GENERATOR - AEROSOL ST046 11135
003400
003400 CAPABILITY TO PRODUCE FINE LIQUID PARTICLES REQUIRED FOR VARIOUS
003400 EXPERIMENTS. USED IN ST01S, ST02S AND ST03S.
003400

314000 GENERATOR - ION ST047 11113
314000
314000 IONIZE AEROSOL PARTICLES SO AS TO OBSERVE ZERO GRAVITY MOTION.
314000
METRONICS ASSOCIATES 8BULLETIN 24-70 111

211000 VAPORIZER - WATER ST048 21135
211000
211000 CAPABILITY TO PROVIDE COOL WATER VAPOR TO EXPERIMENT CONTAINER.
211000
229600 LIGHT SOURCE ST049 21111
229600
229600 PORTABLE LIGHT SOURCE USED FOR PHOTOGRAPHING EXPERIMENT OPERA
229600 TIONS AND RESULTS. USED IN ST01S, ST02S AND ST03S.
229600
229600 ILLUMINATION 100 TO 1000 LUMENS
229600 DISTANCE UP TO 30 FEET
229600

297500 ENVIRONMENT CHAMBER ST050 11135
297500
297500 PROVIDE SPECIAL CHAMBER TO MONITOR OPTICAL PROPERTIES OF AEROSOLS
297500 WITH CAPABILITY TO PERTURBATE THE ENVIRONMENT WITH IONS AND
297500 OBSERVING THE AEROSOLS THROUGH A MICROSCOPE AND PHOTOGRAPHING THE
297500 SCENE. USED IN ST01S, ST02S AND ST03S.
297500

269000 MICROSCOPE ST051 21111
269000
269000 MICROSCOPIC EXAMINATION OF AEROSOL PARTICLES IN THE VISUAL RANGE.
269000 ALSO COMPATIBLE WITH BIOLOGICAL AND METALLURGICAL EXPERIMENTS.
269000
269000 MAGNIFICATION RANGE 50 TO 1000 X
269000 PHOTOGRAPHIC CAPABILITY YES
269000
269000 USED ON ST01, ST02S AND ST03S.

BAUSCH AND LOMB INC 8B57-636 8876 111
OBJECTIVES 4 FLAT FIELD
TYPE OBJECTIVE 5X, 0.8 NA; 10X, 0.25 NA; 40X,
..... 0.65 NA; 100X, 1.25 NA
ILLUMINATOR 5 STEP, LOW VOLTAGE TRANSFORMER
EYEPieces PAIRED 10X

CARL ZEISS CO 8ULYRAPHOT II 88000 111
CONFIGURATION LOOSE EQUIP
MAGNIFICATION 2.5X TO 2500X
CAMERA ATTACHMENT POLAROID, 4X5 FILM, 35 MM ROLL

059055 CAMERA - CINE ST054 01111
059055
059055 PHOTOGRAPHIC RECORD OF THE ACTIVITY OF AEROSOL PARTICLES IN
059055 NORMAL, SLOW MOTION AND TIME-LAPSE SEQUENCES
059055 USED IN ST01S, ST01S, ST03S, ST04S, ST05S AND ST06S.
059055
059055 VIEW ANGLE 40 DEGREES
059055 RESOLUTION 200 LINES/MM
059055 SHUTTER 1/50 TO 1/500 SEC
059055 FILM FORMAT 16 MM
059055 CAMERA SPEED 1, 8, 16, 64 FPS
059055



Space Division
Rockwell International

CINEPA REVELIER CORP #R-128 \$3000 112
CONFIGURATION BENCH/LOOSE EQUIP
OPERATION MODE VARIABLE FRAME RATE OR SINGLE
FRAME PULSE OPERATION UNDER
REMOTE CONTROL
SPEED 2 TO 64 FRAMES/SEC (VARIABLE)
FILM SIZE 16 MM
LENS C MOUNT LENS SYSTEM
FILM CAPACITY VARIOUS FILM MAGAZINES

J A MAURER INC #317-E1 \$4500 311
MODES 8, 16, 32, 64 FPS ELECT CONT
SHUTTER ROTARY
SHUTTER SPEED 1/62.5 TO 1/500 SEC
FILM 16 MM MAGAZINE
FILM CAPACITY 150 FT THIN BASE
LENS 25 MM, F/1.4
RESOLUTION 236 L/MM ON AXIS, 150 L/MM AWAY

424000 TAPE RECORDERS (HIGH-SPEED MULTICHANNEL) ST056 11111
424000
424000 MULTIPLE CHANNEL HIGH SPEED ELECTRONIC DATA RECORDING SYSTEM.
424000 USED ON ST015, ST025 AND ST035.
424000
424000 BANDWIDTH 1 KHZ TO 2.0 MHZ
424000 CHANNELS 14
424000 SIGNAL TYPE ANALOG
424000

AMPEX CORPORATION #AR-200 121
BANDWIDTH 100 HZ TO 3125 HZ THROUGH
300 HZ TO 250 KHZ
TAPE WIDTH 0.5 INCH
TAPE SPEED 1.875 TO 60 IPS
RECORDING TIME 8 MINUTES TO 4 HRS AND 16 MINS
FORMAT DIGITAL
TRACKS 8 DIGITAL, 7 ANALOG

SANGAMO ELECTRIC #SARER 111 \$26,500 111
TRACKS 14
TAPE WIDTH 1 INCH
TAPE SPEEDS 8 SELECTABLE SPEEDS FROM 15/16
TO 120 IPS
FREQUENCY RESPONSE 400 HZ TO 2.0 MHZ
RECORDING RATE 600 KFPS AT 120 IPS SERIAL MODE
WEIGHT 100 POUNDS

DIGI-DATA CORPORATION #1600 \$2650 111
TAPE SPEED 25, 18.75, 12.5 IPS
TRACKS 7 OR 9 TRACK
DATA DENSITY 1600 CPI PHASE ENCODED
200, 556, 800 CPI NRZI
TAPE 0.5 INCH, 1.5 MIL, 1200 FEET

DIGI-DATA CORPORATION #1700/POV-11/7-9 \$5250 111
TRACK NRZI
TAPE SPEED 45, 37.5, 25, 18.75, 12.5 IPS
TAPE 0.5 INCH, 1.5 MIL, 1000/ANSI
COMPATIBLE, 10.5 INCH REEL
TRACKS 7 OR 9
DATA DENSITY PHASE ENCODED
COMPATIBLE

HEWLETT PACKARD #7570B/C \$4600 111
TAPE FORMAT 800, 556, OR 200 CPI NRZI AND
1600 CPI PHASE-ENCODED
CHANNELS 7 OR 9
TAPE SPEED 10 TO 45 IPS
TAPE 0.5 INCH, 1.5 MILS, 1000/ANSI
COMPATIBLE

BORG WARNER #PERT 121
TRACKS 30
TAPE SPEED UP TO 1000 IPS
TAPE WIDTH 1/2 INCH
TAPE LENGTH 2400 FEET
PACKING DENSITY 15 KB/1/2
SIGNAL/NOISE 24 DB
DATA CAPACITY 1.3E+9
BANDWIDTH 6-15 MH/SEC/T

LEACH #LTR 7000 311
TRACKS 12
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
TAPE LENGTH 9200 FEET
PACKING DENSITY 16.7 KB/1/2
SIGNAL/NOISE 22 DB
DATA CAPACITY 2.2E+10
BANDWIDTH 2 MH/SEC/T

HEWLETT PACKARD #35550 \$10200 111
BANDWIDTH 300 KHZ
CHANNELS 7
RECORDING FORMAT DIRECT OR FM
HEWLETT PACKARD #3555C \$14700 111
BANDWIDTH 300 KHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3550A-Q11 \$23800 111
BANDWIDTH 500 HZ TO 2 MHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

ORIGINAL FROM SET
OF POOR QUALITY



HEWLETT PACKARD #7460A/13045A/13063A 94-796 111
CONFIGURATION RACK MOUNTED
TAPE SPEED 15, 3 AND 3/4, 15/16 IPS
CHANNELS 4
RECORDING FORMAT FM
PASSBAND 5 KHZ
S/N RATIO 40 DB
TRACKS 30

HONEYWELL #5600 69730
PORTABLE TAPE RECORDER
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 60IPS
MAXIMUM BANDWIDTH(DIRECT) 300KHZ
PACKING DENSITY UP TO 600PRI
WEIGHT 70LBS (32KG)

HONEYWELL #96 617420
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 240IPS
REEL SIZE 16IN
MAXIMUM BANDWIDTH(DIRECT) 2H KZ
TAPE WIDTH 1/2 IN

APPEX CORP #AR 700 629.180 511
TRACKS 14
TAPE SPEED 60 IPS
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KB/1/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 8E+10 BITS
BANDWIDTH 1 MB/SEC/T

APPEX CORP #AR 1700 629.000 112
TRACKS 14
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KB/1/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 6.2E+10
BANDWIDTH 2 MB/SEC/T

459608 SPECTROMETER - GAMMA RAY ST063 11111
459608 MEASURE SPACECRAFT GAMMA RADIATION BY CRYOGENICALLY COOLED LITH-
459608 IUM DRIFTED GERMANIUM DETECTOR WITH ANTI-COINCIDENCE SHIELD. USED
459608 IN ST01S AND ST03S.
459608 ELECTROMAGNETIC ENERGY 100 KEV TO 6 MEV
459608 SPECTROMETER - CHARGED PARTICLE ST064 11111
459608 THIRTY MULTI-WIRE PROPORTIONAL COUNTER WITH INTERVENING LAYERS OF
459608 TISSUE EQUIVALENT PLASTIC. USED IN ST01S AND ST03S.
459608 PROTON ENERGY RANGE 5 TO 230 MEV
459608 DEUTERIUM ENERGY RANGE 10 TO 300 MEV
459608 TRITIUM ENERGY RANGE 19 TO 230 MEV
459608 HELIUM-3 NUCLEI ENERGY RANGE 20 TO 400 MEV
459608 ALPHA PARTICLE ENERGY RANGE ... 29 TO 450 MEV
459608 SPECTROMETER - NEUTRON ST065 11111
459618 LIQUID SCINTILLATOR WITH ANTI-COINCIDENCE SHIELD PULSE SHAPE
459618 DISCRIMINATOR. USED ON ST01S AND ST03S.
459618 NEUTRON ENERGY RANGE 0.5 TO 15 MEV
103500 DATA CONVERTER - ANALOG TO DIGITAL ST067 21111
103500 CONVERT ANALOG SIGNAL DATA TO DIGITAL DATA FOR RECORDING AND
103500 LATER DUMPING TO GROUND STATIONS. USED ON ST01S AND ST03S.
103500
300900 FILM BADGE ST070 121111
300900 MONITOR ASTRONAUT PERSONAL RADIATION DOSE, DOSE RATE AND DOSE
300900 HISTORY. USED IN ST01S, ST02S AND ST03S.
300900 DOSE RATE 0 - 0.1 MR/HR
300900 0 - 1000 R/HR
300900 DOSIMETER - THERMOLUMINESCENT ST071 41111
300900 PASSIVE DOSIMETER USING THERMOLUMINESCENCE TO MONITOR RADIATION.
300900 USED IN ST01S AND ST02S.
300900 ENERGY RANGE DOSE 1 TO 10E+05 RADS
300900
048700 REFRIGERATOR ST080 11111
048700 PROVIDE A COLD STORAGE AREA TO HOUSE SPECIMEN. USED IN ST01S,
048700 ST02S AND ST03S.
048700 TEMPERATURE 3 C (270 K)
048700 SIZE APPROXIMATELY 3 X 2 X 1 FT (76.2X50.7X25.4CM)



048700
DILLON-LILLY #512-D 1875 112
CAPACITY 12.3 CU FT
TEMPERATURE 3.7 C
SHIELDED AREA NONE
DRAWERS 12
COMPRESSOR 0.25 HP
POWER 115 V 60 HZ
WEIGHT 420 LBS
DIMENSIONS 36X56.5X28.5 INCHES

JORDAN #FT-1-BR \$1177 112
TEMPERATURE 2+22-6.11C IN 1.11C INCR.
CAPACITY 22.4 CU FT
COMPRESSOR 0.25 HP
SHIELDED AREA NONE
DRAIN NONE REQUIRED
POWER 115V 60 HZ 6A
DIMENSIONS 74X28.5X31.625 INCHES

BIRKEN MANUFACTURING CO #STM-KLDC-AC1 \$335 114
CAPACITY 8 FT³
VOLTAGE 115 VAC, 12, 24, 32 VDC

FRIGID CAR #14 \$725 112
CAPACITY 14 CU FT
TEMPERATURE -2.7 TO +15 C
SHIELDED AREA NONE
EXPLOSION-PROOF YES
SEPARATE FREEZING COMP. YES
COMPRESSOR 0.25HP
POWER 115V 60 HZ 3.5A
DIMENSIONS (MAIN CHAMBER) 26.5X20X36.5 INCHES

217500 COLONY GROWTH ENVIRONMENTAL CHAMBER STO 1 11112
217500
217500 HOUSING FOR COLONIES OF BACTERIA WHICH CAN BE OBSERVED DURING
217500 EXPERIMENT. USED IN ST01S, ST02S AND ST03S.
217500
217500 SIZE APPROXIMATELY 5 X 5 X 7 IN (13 X 13 X 18 CM)
217500
438500 MICROBIOLOGICAL SAMPLE STORAGE ST0R2 21112
438500
438500 STORAGE RACK FOR SPACE-TYPE PLASTIC BAG, SWABS AND CAPSULES TO
438500 DESTROY BACTERIA IN CASE OF BREAKAGE. USED IN ST01S, ST02S AND
438500 ST03S.
438500
438500 SIZE APPROXIMATELY 16X12X6 IN (40X33X15 CM)
438500
217500 INCUBATOR ST0R3 11111
217500
217500 CONTAINER USED TO PROVIDE VARIABLE TEMPERATURE ENVIRONMENTS FOR
217500 BIOSCIENCE SPECIMEN. USED IN ST01S, ST02S AND ST03S.
217500
217500 SIZE APPROXIMATELY 2X1.5X0.5 FT (66X51X15 CM)
217500 TEMPERATURE TO 100 F (40 C)
217500

SARGENT WELCH SCIENTIFIC CO. #S-4330E-B \$925 111
RANGE 35 TO 70C
MAX. TEMP 90C

042000 COUNTER PRINTER-COUNTER ST0R4 11111
042000
042000 PROVIDE AUTOMATIC BLOOD ANALYSIS (HEMOGLOBIN, HEMATOCRIT, RED
042000 BLOOD CELL COUNT, WHITE BLOOD CELL COUNT, MEAN CELL VOLUME, MEAN
042000 CELL HEMOGLOBIN AND ITS CONCENTRATION). USED IN ST01S, ST02S
042000 AND ST03S.
042000

COUNTER ELECTRONICS #FN \$5000
PARTICLE SIZE 0.003 TO 24 MICRONS
COUNTING RATE 0.040 ML/SEC
DIMENSIONS 49CM H X 35 CM W X 44 CM D
WEIGHT 54LB (24.5KG)

COUNTER ELECTRONICS INC #S \$5000
PARAMETERS HEMOGLOBIN, HEMATOCRIT, RED
BLOOD CELL COUNT, WHITE BLOOD
CELL COUNT, MEAN CELL VOLUME, AND MEAN
CELL HEMOGLOBIN CONCENTRATION
TYPE AUTOMATIC

COUNTER ELECTRONICS INC #2R1-6 \$10000 112
PARAMETERS HEMOGLOBIN, HEMATOCRIT, RED
BLOOD CELL COUNT, WHITE BLOOD
CELL COUNT, MEAN CELL VOLUME,
MEAN CELL HEMOGLOBIN, AND MEAN
CELL HEMOGLOBIN CONCENTRATION
TYPE MANUAL

459620 SPECTROMETER - VOLUMETRIC ST0R5 11111
459620
459620 PROVIDE A PHOTOGRAPHIC MEANS TO MEASURE INDEPENDENT VOLUME OF
459620 CELLS AND CELL SHAPES. USED IN ST01S, ST02S AND ST03S.
459620
149300 ELECTROPHORESIS APPARATUS ST0R6 11111
149300
149300 MEASURE ELECTROPHORESIS MOBILITY, SURFACE-ZETA POTENTIAL, AND
149300 SURFACE CHARGE DENSITY OF CELL LINES OVER THEIR LIFE CYCLE. USED
149300 IN ST01S AND ST02S.
149300

ORIGINAL PAGE IS
OF POOR QUALITY

PAUSCH - LOMB INC. #SPECTROPHOT I 111
MIGRATION VOLTAGE 0-50CV
MIGRATION CURRENT 0-300MA
UV LIGHT SOURCE DEUTERIUM
VISIBLE LIGHT SOURCE QUARTZ IODINE
MONOCHROMATOR RANGE 200-650 MILLIMICRONS
NO. SAMPLES 16
MIGRATION TIME 20MINUTES APPROX.
BUPPER GEL THICKNESS 1.7MM
GEL COMPOSITION AGAROSE

067000 CENTRIFUGE ST090 11111
067000
067000 PROVIDE VARIABLE GRAVITY ENVIRONMENTS TO BIOSCIENCE SPECIMEN
067000 GENERALLY LIMITED TO TEST-TUBE SIZE CONTAINERS. USED IN ST015,
067000 ST025 AND ST035.
067000
067000 GRAVITY LEVEL INDUCED UP TO 13 G REQUIRED
067000

CLAY-ADAMS, INC. #1CG4 8324
SPEED RANGE 400 TO 3040 RPM
SAMPLES 12 (MAX), 15 ML EACH
DIMENSIONS 19X12X26 INCHES
CONFIGURATION BENCH MOUNT

INTERNATIONAL EQUIPMENT CO #2C211/2C215/2C216 885 112
CONFIGURATION BENCH/LOOSE EQUIP
APPLICATION MICRO OR SEMI-MICRO ANALYSIS
SPINNING SPEED 1780 RPM
SPIN STABILIZATION TIME 10 SEC
ACCOMMODATIONS 0.5, 1, 2, 3, 5 ML TUBES

ACAMS #0YNAC 8267
SPINNING RATE 700 TO 2500RPM
APPLIED GRAVITY 55 TO 1080 GCF
SAMPLE CAPACITY 24

MSE #GT-4 8410
SPINNING SPEED 4000RPM
APPLIED GRAVITY 2900GCF
SAMPLE CAPACITY 4

PAUSCH AND LOMB INC. #SYSTEM CENTRIFUGE
CAPACITY 100 SAMPLES
TIME TO SPIN DOWN 3 MINUTES
MAX. IMBALANCE 2 TUBES
SPEED UP TO 4000RPM
RCP 2000 AT 4000RPM
BECKMAN INSTRUMENTS #152 8245
SPINNING RATE 0 TO 16000 RPM
ACCOMMODATION4ML TUBES(20)
DIMENSIONS 9X7X7 INCHES
WEIGHT 13 LBS(SHIPPING)
POWER 115V 50/60 HZ
DEVELOPED G 1300G

410000 PRESSURE TANK ST098 41111
410000
410000 PROVIDE PRESSURE SOURCE FOR WATER FEED SYSTEM TO STEAM GENERATOR
410000 EXPERIMENT
410000 OPERATING PRESSURE 100 PSI(6.9X10E5 N/M**2)
410000
410000 WATER RESERVOIR WITH HEATER ST049 11113
410000
410000 STORE WATER FOR STEAM GENERATOR EXPERIMENT
410000
410000 OPERATING PRESSURE 30PSI(2.07X10E5 N/M**2)
410000
415500 PLMP ST100 11111
415500
415500 PUMP WATER TO STEAM GENERATOR
415500
415500 FLOW RATE 0.1 LB/HR (0.045KG/MIN)
415500
211000 STEAM GENERATOR ST101 11135
211000
211000 DEVICE TO GENERATE STEAM IN ZERO GRAVITY
211000
004100 AIR SAMPLE UNIT ST105 21111
004100
004100 PERIODIC COLLECTION OF AIR SAMPLES FROM THE SHUTTLE AND ATL
004100 PROVIDING TYPE, RATE OF CHANGE AND QUALIFICATION OF MICRO-
004100 ORGANISMS AND NON-VIABLE PARTICLES. USED ON ST015 AND ST035.
004100

PAUSCH - LOMB INC. #40-1A 111
SAMPLING RATE 0.1 CL. FT./100SEC.
(0.17 LITERS/MIN.)
MAX. CONCENTRATION MEAS 1E+7 PARTICLES/CL. FT.
(13E+3 PARTICLES/LITER)
PARTICLE SIZE RANGE 0.3-0.5-1.0-2.0-3.0-5.0-10MICRONS
408500 TIMER ST107 11111
408500
408500 TIMING AND VARIABLE PROGRAMMING CAPABILITY TO OPERATE AND CONTROL
408500 EXPERIMENT AND SUPPORTING EQUIPMENT. USED IN ST015 AND ST035.
408500

ROHDE AND SCHWARZ #C40 (100.6597.91) 85400
DISPLAY 6 DIGITS: HR, MIN, SEC
OUTPUT I-OUT-OF-N AND/OR BCD CODE
PROGRAMMER TIME PROGRAMS(COUNT
DOWN) AND TIME SIGNALS



DATATRON #1350-50N \$117.0

INPUT FREQUENCY..... 60HZ

OUTPUTS..... VISUAL, ACO

DISPLAY..... HRS, MIN, SEC

TENNELEC #TC545 \$550

COUNT RATE..... 20MHZ

TIMERASE..... 0.1 OR 0.01SEC

ACCURACY..... SAME AS LINE FREQUENCY

CONFIGURATION..... NIM COMPATIBLE

187000 WORK BENCH \$1142 \$1114

187000

187000 PROVIDE A WORK AREA UPON WHICH EXPERIMENTS AND MAINTENANCE MAY BE

187000 PERFORMED. ENCLOSED STORAGE TO BE INCORPORATED IN BENCH. USED

187000 ON ST01S, ST02S, ST03S, ST04S AND ST05S.

187000

187000 SIZE APPROXIMATELY 2X3X4 FT 10.5X10.5X1.2 M

187000

424000 TAPE RECORDER, ANALOG \$1009 \$1111

424000

424000 RECORD VOICE, FM AND VIDEO DATA OF LANDMARK AND STAFF-TRACKING

424000 FILMS, USED IN ST02S, ST03S AND ST05S.

424000

424000 BANDWIDTH TO 2.5 MHZ

424000 RECORDING FORMAT VOICE, FM AND VIDEO

424000 CHANNELS 7

424000

AMPEX CORPORATION #A9-200 121

BANDWIDTH 100 HZ TO 1125 HZ THROUGH

300 HZ TO 250 KHZ

TAPE WIDTH 0.5 INCH

TAPE SPEED 1.875 TO 60 IPS

RECORDING TIME 4 MINUTES TO 4 HRS AND 16 MINS

FORMAT DIGITAL

TRACKS 4 DIGITAL, 7 ANALOG

AMPEX CORP #AR 1700 \$28,000 112

TRACKS 14

TAPE SPEED 120 IPS

TAPE WIDTH 1 INCH

RECORDING MODE DIRECT

PACKING DENSITY 20 KB/1/7

SIGNAL/NOISE 20 DB

DATA CAPACITY 6.2E+10

BANDWIDTH 2 MB/SEC/T

SANJANO ELECTRIC #SABER 111 \$26,500 111

TRACKS 14

TAPE WIDTH 1 INCH

TAPE SPEEDS 8 SELECTABLE SPEEDS FROM 15/16

TO 120 IPS

FREQUENCY RESPONSE 400 HZ TO 2.0 MHZ

RECORDING RATE 600 KIPS AT 120 IPS SERIAL MODE

WEIGHT 100 POUNDS

DIGI-DATA CORPORATION #1400 \$2650 111

TAPE SPEED 25, 18.75, 12.5 IPS

TRACKS 7 OR 9 TRACK

DATA DENSITY 1600 CPI PHASE ENCODED

200, 556, 800 CPI NRZI

TAPE 0.5 INCH, 1.5 MIL, 1200 FEET

DIGI-DATA CORPORATION #1700/PDR-11/7-5 \$5250 111

TRACK NRZI

TAPE SPEED 45, 37.5, 25, 18.75, 12.5 IPS

TAPE 0.5 INCH, 1.5 MIL, 1000/ANSI

COMPATIBLE, 10.5 INCH REEL

TRACKS 7 OR 9

DATA DENSITY PHASE ENCODED

COMPATIBLE

HEWLETT PACKARD #7470B/C \$4600 111

TAPE FORMAT 800, 556, OR 200 CPI NRZI AND

1600 CPI PHASE-ENCODED

CHANNELS 7 OR 9

TAPE SPEED 10 TO 45 IPS

TAPE 0.5 INCH, 1.5 MILS, 1000/ANSI

COMPATIBLE

RORG WARNER #PERT 121

TRACKS 10

TAPE SPEED UP TO 1000 IPS

TAPE WIDTH 1/2 INCH

TAPE LENGTH 2400 FEET

PACKING DENSITY 15 KB/1/7

SIGNAL/NOISE 24 DB

DATA CAPACITY 1.3E+9

BANDWIDTH 6-15 MB/SEC/T

LEACH #MTD 7000 311

TRACKS 12

TAPE SPEED 120 IPS

TAPE WIDTH 1 INCH

TAPE LENGTH 4200 FEET

PACKING DENSITY 16.7 KB/1/7

SIGNAL/NOISE 22 DB

DATA CAPACITY 2.2E+10

ORIGINAL PAGE IS
OF POOR QUALITY



BANDWIDTH 2 MB/SEC/T

HEWLETT PACKARD #3555D 810200 111
BANDWIDTH 300 KHZ
CHANNELS 7
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3555C 814700 111
BANDWIDTH 300 KHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3550A-G11 823800 111
BANDWIDTH 500 HZ TO 2 MHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3560A/13065A/13063A 84796 111
CONFIGURATION RACK MOUNTED
TAPE SPEED 15, 3 AND 3/4, 15/16 IPS
CHANNELS 4
RECORDING FORMAT FM
PASSBAND 5 KHZ
S/N RATIO 48 DB
TRACKS 30

MCNEWMELL #5800 89730
PORTABLE TAPE RECORDER
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 60IPS
MAXIMUM BANDWIDTH(DIRECT) 300KHZ
PACKING DENSITY UP TO 600PBI
WEIGHT 70LBS (32KG)
INPUT VOLTAGE 28VDC

MCNEWMELL #9A 817420
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 240IPS
PEEL SIZE 16IN
MAXIMUM BANDWIDTH(DIRECT) 24 HZ
TAPE WIDTH 1/2 IN

AMPEX CORP #AR 70C 824180 511
TRACKS 14
TAPE SPEED 60 IPS
TAPE WIDTH 1 INCH
RECORD MODE DIRECT
PACKING DENSITY 20 KB/I/T
SIGNAL/NOISE 20 DB
DATA CAPACITY 8E+10 BITS
BANDWIDTH 1 MB/SEC/T

227000 LASER 57118 11111
227000
227000 MODE-LOCKED LASER RANGING SYSTEM MEASURING RANGE, LINE-OF-SIGHT
227000 ANGLES, AND RANGE RATE BETWEEN ATL/SHUTTLE AND SEPARATE SPACE
227000 VEHICLES INCLUDING RANGE OF ATL TO GROUND.
227000
227000 RANGE TOLERANCE +/- 1.10 IN (+/- 3 CM) ATL-GND
227000 COOPERATIVE TARGET DISTANCE .. 300 NM (550 KM)
227000

618015 RECEIVER - RF 57127 11111
618015
618015 MEASURE ELECTROMAGNETIC INTERFERENCE AT ORBITAL ALTITUDES IN THE
618015 FREQUENCY SPECTRUM OF 400 MHZ TO 15 GHZ. USED ON ST02S AND ST04S.
618015
618015 OVERALL SPECTRUM 400 MHZ TO 15 GHZ
618015 SPECIFIC MULTIPATH MEASURE 1.5 TO 1.6 GHZ
618015 2.025 TO 2.3 GHZ
618015 13.4 TO 15.35 GHZ
618015

SINGER 813500
FREQ. RANGE 30 +0 1000 MHZ IN 8 BANDS
ACCURACY +/- 2 % FREQUENCY, +/- 3%
VOLTAGE MEAS. RANGE 600R METER SCALE + 20,40,
60,80 DB ATTENUATOR STEPS
TOTAL 1400B
DETECTOR FUNCTIONS FIELD INTENSITY, QUASI-PEAK,
DIRECT PEAK, SLIDEBACK PEAK,
BPO, PH DISCRIMINATOR, LINEAR

018700 TELESCOPE - POINTING 57152 11133
018700
018700 POINTING TELESCOPE WHICH ALLOWS EXPERIMENT PHENOMENA TO BE
018700 OBSERVED AND TRACKED IN SPACE
018700
407000 POWER SUPPLY 57171 21111
407000
407000 ELECTRICAL POWER REGULATION AND DISTRIBUTION AS REQUIRED AND AT
407000 VARIABLE VOLTAGE RANGES. USED IN ST02S, ST03S AND ST04S.
407000
407000 VOLTAGE 0 TO 24 VDC
407000 POWER 200 WATTS
407000 SENSITIVITY 0.05 VOLT
407000

POWER DESIGNS INC. #APC-420 8275
OUTPUT +/-12VDC &+ 0-2AMP



REGULATION	+/- 24VDC AT 0-1AMP	
IMPEDANCE	+/- 0.05 %	
CIRCUITS	0.3 OHM UP TO 150 KHZ	
WEIGHT	17X3-1/2X5-1/2 INCHES	
BOMDE AND SCHWARZ	16.5 LB	
VOLTAGE	BNGRE (100.0254.14) \$1500	111
CURRENT	0 TO 30 V	
	0 TO 10 A	
JOHN FLUKE MFG CO INC	#4250A	\$1795 111
PROG CODING	RCD (BINARY AVAILABLE)	
OUTPUT VOLTS	0 TO +/- 65 VOLTS DC/PEAK AC	
OUTPUT CURRENT	1 AMP	
FREQUENCY RANGE	DC TO 30 KHZ	
RESOLUTION	1 MV	
REGULATION	0.001%	
JOHN FLUKE MFG CO INC	#382A	\$1795 111
OUTPUT VOLTAGE	0 TO 50 VDC; 0 TO 5 VDC	
OUTPUT CURRENT	0 TO 2 AMPS	
REGULATION	0.005%	
STABILITY	+/- 0.0025% PER DAY	
SORENSEN	#SRL40-12	\$550 111
OUTPUT VOLTAGE	0 TO 40 VOLTS	
OUTPUT CURRENT (MAX)	12 AMPS	
REGULATION	+/- 0.1% OR +/- 2MV	
CONFIGURATION	RACK MOUNT	
SIZE	51/4X19X15 1/8IN	
WEIGHT	64LB (30KG)	
COOLING	NATURAL CONVECTION	
SORENSEN	#MUR 72-20	\$450 111
MIN COMPATIBLE POWER SUPPLY		
OUTPUT VOLTAGES	+/- 12V; +/- 24V	
OUTPUT CURRENT	2ADC; 1ADC	
REGULATION	+/- 0.05%	
CONFIGURATION	MOUNTS TO BACK OF NIM BIN	
SIZE	16.825X3.43X5.5IN	
WEIGHT	11.5LB (5.2KG)	
AYDIN CORP, VECTOR DIV.	MODEL TV-64	\$325 111
INPUT POWER	30MA AT 24VDC	
INPUT VOLTAGE RANGE	28+/- 4VDC	
OUTPUT REGULATION	0.2%	
OUTPUT RIPPLE	LESS THAN 5MVP-P	
OUTPUT CURRENT	+/- 22VDC AT 100MA EA.	
MODULAR POWER INC.	#T-27	\$129
OUTPUTS	4	
POWER OUTPUT	2 TO 50 VOLTS, 1 TO 10 AMPS	
AC INPUT @ 50/1000 HZ	105/125 OR 210/250 VAC	
RIPPLE OUTPUT	500 MICROWOLTS (RMS)	
WEIGHT	8.5 LBS.	
SIZE (INXHD)	4.25 X 5.5 X 3.63 INCHES	
TEMPERATURE BASE	71 DEG C	
TENNELEC	#TC911	\$275 111
INPUT VOLTAGE	103 TO 129VAC 50 TO 60HZ	
OUTPUT	20R TO 250VAC 50 TO 60HZ	
	+/- 24VDC AT 1AMP	
	+/- 12VDC AT 2AMP	
424100 RECORDER, OPTICAL		ST273 11235
424100		
424100 RECORD AND PLAYBACK DIGITALLY GENERATED DATA FROM A LASER SYSTEM.		
424100 USED ON STOPS AND STOPS.		
424100		
424100 BANDWIDTH	TO 50 MHZ	
424100 SIGNAL	DIGITAL	
424100 CHANNELS	3	
424100		
MRB SINGER INC.	#MMD-7	115
FILM SIZE	5 IN.	
FILM LENGTH	250FT.	
RESOLUTION	4500 LINES/TRACK	
	(2250 LINES/FRAME)	
FILM SPEED	17 FPS/IN.	
DATA RECORD	DATA CARD, HOOK CLOSER, COUNTER, E	
	FIDUCIAL MARK	
FILM TRANSPORT SPEED	UP TO 200 IN./MIN. FUNCTION V/M	
	RATIO	
RESALAB SCIENTIFIC	#HIDRABD DIG. REC. / REPRODUCE	212
DATA RATE	85 MEGABITS/SEC.	
CHANNEL BANDWIDTH	100MHZ/CHANNEL	
RECORD METHOD	MULTILASER	
INTERNAL DIGITAL DATA RATE	22 MEGABITS/SEC.	
FILM RATE	160FT/MIN	
FILM CAPACITY	4000 FT	
NUMBER RECORD TRACKS	4	
FILM SIZE	70MM	
CONTRAST RATIO	100:1	
179500 FREQUENCY METER		ST156 11111
179500		
179500 MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD, AVERAGE RATIO AND		
179500 MULTIPLE RATIO OF VHF FREQUENCY SOURCES. USED ON STOPS.		
179500		
179500 BANDWIDTH	10 TO 100 MHz	
179500 SENSITIVITY	+/- 2 %	
179500		

ORIGINAL PAGE IS
OF POOR QUALITY



JOHN FLUKE MFG CO INC #1452A-C7 1995 111
 FREQUENCY RANGE DC TO 515 MHZ
 DISPLAY 7 DIGIT LED (R OR 9 CRITICAL)
 SENSITIVITY 50 MV RMS
 PERIOD RANGE DC TO 10 MHZ (DC COUPLED)
 5 HZ TO 10 MHZ (AC COUPLED)
 TIME INTERVAL RANGE 0.1 MICROSEC TO 1E+6 SEC

ROMBE AND SCHWARZ #PET2(100.6039.02/ 111
 110.6125.02/
 100.6083.02/
 100.6091.02/
 100.6068.02/
 FREQUENCY RANGE 0 - 100 MHZ
 SENSITIVITY -9 TO +9 VOLTS
 PERIOD RANGE 0.5 MICROSEC TO 100 SEC
 TIME INTERVAL RANGE 0.2 MICROSEC TO 1E+6 SEC

480000 CAMERA, TV ST025 31131
 480000
 480000 REAL-TIME VIEWING OF TARGET AREAS, LAND MARK AND STAR FIELD
 480000 TRACKING, AND SENSOR BORESIGHTING.
 480000 USED ON ST04S, ST05S AND ST06S.
 480000
 480000 VIEW ANGLE P TO 10 DEGREES
 480000 RESOLUTION 1025 LINES
 480000

RCA ELECTRONICS COMPONENTS #8521 5520 111
 TYPE VIDICON-SULFIDE
 PHOTOCONDUCTOR II
 IMAGE DIAGONAL 25 MM, (1 IN)
 FOCUS MANUAL
 LIMITING RESOLUTION 1500 LINES

COMU INC. ELECTRONICS DIV. #1220 58430 111
 IMAGE CONVERTER 4848R VIDICON - COLOR
 RESOLUTION-HORIZONTAL LIMIT 300 LINES MIN.
 GEOMETRIC DISTORTION < 2 % OF PIC HEIGHT
 LENS MOUNT 16 MM - C-MOUNT

CCMU INC. ELECTRONICS DIV. #2000 53300 311
 LENS-BUILT IN 411 (20-80 MM) F/2.5 ZCCM
 LENS ATTACHMENT 411 (12.5-80 MM)
 LENS INTERCHANGEABLE 1C11 (15-150 MM) F/2.8 ZCCM
 16 MM C-MOUNT PFL
 10 MHZ BANDWIDTH (MOD 2006) 525 OR 729 LINE (VIDICON 7263A)
 20 MHZ BANDWIDTH (MOD 2004) 873 OR 945 LINE (VIDICON 8973)
 FIL SPECIFICATIONS MIL-E-5272C, MIL-STD-810

COMU INC. ELECTRONICS DIVISION #4500 52225 311
 LENS 411 700M
 1C11ZOOM
 VERTICAL SWEEP RATE 60 FIELDS PER SEC
 HORIZONTAL SWEEP RATE 525 LINES PER FRAME
 IMAGE TUBE TYPE 8541A STD
 LENS MOUNT 16 MM C-MOUNT
 FIL SPECIFICATIONS MIL-E-5400M, MIL-E-5300M

490000 MASS SPECTROMETER ST163 11111
 490000
 490000 ANALYSIS OF GASEOUS COMBUSTION PRODUCTS FORMED IN ZERO GRAVITY
 490000 VIA MASS SPECTROSCOPY. USED IN ST04S AND ST06S.
 490000
 490000 MASS RANGE 1 TO 1000 AMU
 490000 SCANNING RATE 60 HSEC TO 600 SEC
 490000 SENSITIVITY 1 #

EXTRANUCLEAR LABORATORIES #279-7 533900 111
 MASS RANGE 2 TO 1000 AMU
 RESOLUTION GREATER THAN UNIT TO MASS 1000
 SENSITIVITY 8X1 .075 NG/SEC METHYL STEARATE
 ACCESSORIES GC INTERFACES, PROBES, HEATER, ECT

VEBCC INSTRUMENTS INC #GA-4-BM 59495 111
 MASS RANGE 1 TO 760 AMU
 RESOLUTION RESOLVES TWO MASSES .025% APART
 SENSITIVITY 1X1E-13 TORR FOR N2
 RECORDER DECADES ON SINGLE CHART
 SCAN RATES VERY SLOW TO FAST

074700 GAS CHROMATOGRAPH ST175 21111
 074700
 074700 ANALYSIS OF COMBUSTION PRODUCTS FOR FLAME CHEMISTRY EXPERIMENTS.
 074700 USED IN ST04S AND ST06S.
 074700

BOCKMAN #GC-2A 52050 111
 DETECTOR THERMAL CONDUCTIVITY
 TEMPERATURE CONTROL RANGE 40C TO 240C
 CONTROL ACCURACY +/- 0.1

BOCKMAN INSTRUMENTS #6700 58000
 SENSITIVITY 1 +0.5 PPM
 OPERATING TEMP 95 +0.225C
 TEMPERATURE CONTROL +/- 0.05C

BOCKMAN INSTRUMENTS #6700 58000
 SENSITIVITY 1 +0.5 PPM
 OPERATING TEMP 95 +0.225C
 TEMPERATURE CONTROL +/- 0.05C



CARLE INSTRUMENTS, INC. #900-111 42075
CONTROL SENSITIVITY < 0.001 DEG C
TEMPERATURE RANGE TO 350 DEG C (1 DEG INCREMENTS)
HEAT CAPACITY-TOTAL 600 WATTS
INLETS DUAL ON-COLUMN
IGNITION DETECTORS HYDROGEN FLAME
VALVING MICRO OR MINI VOLUME VALVES
WEIGHT 53 POUNDS
POWER REQUIREMENTS 115 VAC

CARLE INSTRUMENTS, INC. #9500 41185
TEMPERATURE RANGE AMBIENT TO 200 DEG C
DETECTOR DUAL FED SINGLE HYDROGEN FLAME
SENSITIVITY 5E-12 AMPERES
POWER REQUIREMENTS 115 VAC, 50 TO 60 HZ

419500 PYROMETER ST176 11111
419500
419500 MEASURE TEMPERATURE OF COMBUSTION PROCESS IN FLAME CHEMISTRY
419500 EXPERIMENTS, USED IN ST04S AND ST05S
419500
419500 TEMPERATURE RANGE 0 TO 3000 DEG F
419500
IRCON #2000 42400 111
TEMPERATURE RANGE 40 - 3000 DEG C
SPECTRAL RESPONSE 0.70 TO 0.97 MICROM
CALIBRATION ACCURACY +/- 2 PCT FS TEMP OR +/- 5 DEG C
REPEATABILITY +/- 0.3 PCT FS TEMP
RESPONSE TIME AT METERS 95 PCT FS INPUT WITHIN 1 SEC
INDICATOR/CONTROLLER AMB TEMP .. 10 TO 65 DEG C
SENSING HEAD AMBIENT TEMP -32 TO 94 DEG C

098500 COMBUSTION CHAMBER ST174 11135
098500
098500 CHAMBER WHERE FUELS AND OXIDIZERS ARE INJECTED, ALLOWED TO BURN
098500 AND THE FLAME CAPABLE OF BEING MONITORED.
098500
098500 SIZE APPROXIMATELY 1.67 DIA X 3.5 FT
098500 (0.51 DIA X 1.0 M)
098500
410000 TANKS - FUEL/OXIDIZER ST179 41135
410000
410000 SIMILAR TYPE TANKS CAPABLE OF STORING CORROSIVE FUELS AND
410000 OXIDIZERS IN THE SORTIE LABORATORY AT LOW PRESSURES, USED IN
410000 ST04S AND ST05S.
410000
410000 SIZE APPROXIMATELY 0.75 DIA X 3 FT
410000 (0.225 DIA X 0.9 M)
410000
297500 CHAMBER, CUBICAL TEST ST188 11135
297500
297500 CHAMBER FOR PERFORMANCE OF SUPERFLUID HELIUM AND LIQUID DCP/
297500 PARTICLE EXPERIMENTS
297500
297500 DIMENSIONS 1.4 FT (0.55 M) CUBE
297500
003110 ACOUSTIC SPEAKER ST189 61111
003110
003110 PARTICLE POSITIONING DRIVER
003110
003110 FREQUENCY RANGE 50 TO 10KHZ
003110
238500 MAGNET, SUPERCONDUCTING ST190 11112
238500
238500 MAGNET FOR DCP/PARTICLE POSITIONING
238500
286000 MOTOR, ELECTRIC ST191 11111
286000
286000 MOTOR FOR CHAMBER DRIVE SYSTEM
286000
108600 CRYOSTAT - HELIUM ST192 11111
108600
108600 DEWAR TYPE STORAGE CONTAINER TO STORE AND DISPENSE LIQUID HELIUM,
108600
108600 SIZE APPROXIMATELY 1.5 DIA X 1.7 FT
108600 (0.5 DIA X 0.05 M)
108600
CRYOCENIC ASSOCIATES INC #SD-15 42852 112
DEWAR DIA 17 IN.
LIQUID HELIUM CAP 101 LITERS
LIQUID NITROGEN CAP 26.5 LITERS

UNION CARBIDE (LINDE DIV.) #G104-CC52(PART NO) 41150 111
(CAPACITY LIQ HE 30 LITERS
EVAP RATE/DAY08
DIMENSIONS 44X18.75 INCHES
WEIGHT 93 LB (MT) 91.4 LB (FULL)

526200 PRESSURE MONITOR ST193 11114
526200
526200 MEASURES PRESSURE IN TEST CHAMBER AND DISPLAYS DIGITALLY
526200
487000 TEMPERATURE MONITOR ST194 11114
487000
487000 MONITORS AND CONTROLS TEST CHAMBER TEMPERATURE
487000
098500 COMBUSTION CHAMBER ST202 11135
098500
098500 CHAMBER TO EXAMINE COMBUSTION IN ZERO GRAVITY
098500
098500 DIMENSIONS 1.67 FT (0.51 M) DIA X 4.0 FT (1.22 M)

ORIGINAL PAGE IS
OF POOR QUALITY

000500
409600 ELECTRON SPIN RESONANCE/MASS SPECTROMETER ST25R 11112
409600
409600 IDENTIFY CONSTITUENTS OF COMBUSTION PRODUCTS
409600
414000 TEMPERATURE SENSORS ST260 11111
414000
414000 SENSOR TO MEASURE TEMPERATURE OF COMBUSTION PROCESS
414000
414000 TEMPERATURE RANGE(EST)..... 1000 TO 3000F
414000
414000 PRESSURE SENSORS ST260 11111
414000
414000 SENSORS TO MEASURE PRESSURE OF COMBUSTION EXPERIMENT
414000
414000 PRESSURE RANGE..... TWO
414000
410000 OXIDIZER TANK ST262 11135
410000
410000 STORE OXIDIZER FOR COMBUSTION REACTION
410000
410000 VOLUME..... 2 CU.FT. (0.55 CU.M.)
410000
410000 FUEL TANK ST261 11135
410000
410000 STORE FUEL FOR COMBUSTION EXPERIMENT
410000
410000 VOLUME..... 2 CU.FT. (0.55 CU.M.)
410000
410000 INERT GAS SUPPLY ST263 11111
410000
410000 SUPPLY PRESSURE FOR EXPULSION OF FUEL AND OXIDIZER FROM STORAGE
410000 TANKS
410000
410000 VOLUME 0.22 CU.FT.(0.006 CU.M.)
410000
410000
234000 STORE LAMP ST167 11111
234000
234000 VARIABLE DISCHARGING STORE LIGHTS WHICH ALLOW THE STOPPING OF
234000 ACTION.
234000
234000 STORE DISCHARGE RATE 1 TO 1000 FLASHES/SEC
234000
234000
GRANVILLE PHILLIPS DSPECTRASCAN 400 ST000
MASS RANGE 1 TO 400AMU IN ONE RANGE
SCAN RATE VARIABLE 30MILLISECONDS TO 600
SECONDS/SCAN
SENSITIVITY 50 AMP/TOUR FOR M2 AT M/AP=62
031001 CRITICAL STATE TEST CELL ST207 11235
031001
031001 TEST CELL TO TEST CRITICAL POINT PHENOMENA, TEST CELL INCLUDES
031001 A THERMAL SHIELD, HEAT EXCHANGER, MICROMETER, CHAMBER, SUPPORT
031001 ELECTRONICS AND PRESSURE TRANSDUCERS
031001
003210 ACOUSTIC MEASUREMENT DEVICE ST214 11114
003210
003210 DEVICE CONSISTS OF AN ACOUSTIC CAVITY, AN OSCILLATOR AND A
003210 DIGITAL DISPLAY
003210
031002 POOL BOILING TEST CHAMBER ST210 21235
031002
031002 CHAMBER TO EVALUATE POOL BOILING PHENOMENA IN ZERO GRAVITY
031002
414000 INSTRUMENTATION SENSOR SYSTEM ST223 11114
414000
414000 SYSTEM INCLUDE TEMPERATURE, PRESSURE AND ACCELEROMETER SENSORS
414000
414000 CRYSTALLIZER SAMPLE CHAMBER ST225 11135
031003
031003 CHAMBER FOR GROWING CRYSTALS IN ZERO G
031003
031004 SCHLIEREN TYPE OPTICS SYSTEM ST227 11135
031004
031004 OPTICS SYSTEM EMPLOYING SCHLIEREN METHOD TO EXAMINE THE STRUCTURE
031004 OF CRYSTALS
031004
031005 HOLOGRAPHIC INTERFEROMETER ST231 11235
031005
031005 DEVICE TO EXAMINE STRUCTURE OF GROWING CRYSTALS
031005
OPTICS TECHNOLOGY INC #215 61995 111
POWER LEVEL..... 0.5 TO 1.0MW; 1.0 TO 2.0MW
EXPOSURE TIME..... 45 TO 180 SEC; 15 TO 60 SEC
PULSE..... 80NCH
FILM PLATE SIZE..... 4 X 5IN
HOLOGRAM TYPES..... GABOP AND FRESNEL
JCDOA ENGR ASSOC #MS-2 98700 111
ACTIVE VIBRATION ISOLATION
LASER..... 20MW HE NE
SYSTEM COMPONENTS..... REAL TIME PLATE HOLDER, SHUTTER
BEAM STEERER, REFERENCE PEARLS
SPATIAL FILTERS, VARIABLE BEAM
SPLITTER, ENCLOSURE
JCDOA ENGR ASSOC SAP-100 AUTO 98500
PHOTO PROCESSOR
PROVIDES CHEMICAL FILLING AND CIRCULATING, AND DRAINING TO
PROCESS PHOTOGRAPHIC PLATE
001000 ACCELEROMETER, TRIAXIAL ST235 21111
001000

COMMUNICATIONS AND NAVIGATION

620001 TRACKING ANTENNA SERVO ELECTRONICS CAC44 11114
620001
620001 ELECTRONICS FOR ANTENNA POSITIONING AND CONTROL
620001
620001 TYPE TYPE 11 SERVO LOOP
620001 GAIN CONSTANT 3.5E4 PER SEC SQUARED
620001 BANDWIDTH 100 HZ
620001
610030 RECEIVER, L BAND CAC45 21111
610030
610030 RECEIVE L BAND SIGNALS TO DEMONSTRATE THE OPERATIONAL PERFORMANCE
610030 OF TYPICAL SPACEBORNE INTERFEROMETERS FOR AIRCRAFT AND PARIAE
610030 NAVIGATION/TRAFFIC CONTROL. TRACK SIGNAL PHASE.
610030
610030 CARRIER FREQUENCY 1569.75 TO 1570.25 MHZ
610030 SENSITIVITY -120 DBM
610030 RECEIVER TYPE PHASE LOCK, SUPERHETERODYNE,
610030 FIXED FREQUENCY, AGC
610030 DYNAMIC RANGE 40 DB INSTANTANEOUS
610030 100 DB TOTAL
610030 PRE-DETECTION BW 5 MHZ
610030 SWEEP RATE TUNE 3 PPM OF CARRIER FREQ
610030 PHASE LOCK LOOP SNR 40 DB
610030 RECEIVER NOISE FIGURE 5 DB
610030 PHASE LOCK LOOP BANDWIDTH 10 HZ
610030
RMC ELECTRONICS LAB INC RFR W1900ATR 89000 211
FREQUENCY 1.4 TO 2.3 GHZ FIXED TUNED
IF/DISCRIMINATOR BW 50 MHZ
NOISE FIGURE 11 DB
610050 RECEIVER, VHF CAC46 51111
610050
610050 RECEIVE VHF SIGNALS FOR DETAILED EVALUATION OF PERFORMANCE CHAR-
610050 ACTERISTICS OF TORS/GROUND/SHUTTLE DATA LINK. PROVIDES NUMBER OF
610050 RECEIVER FRONT ENDS WHICH CAN BE CONNECTED IN DIFFERENT RECEIVE
610050 CONFIGURATIONS.
610050
610050 CARRIER FREQUENCY 124 TO 130 MHZ
610050 SENSITIVITY +/- 25 MHZ
610050 CARRIER/NOISE DENSITY RATIO .. 53 DBM
610050 DATA RATE 100 TO 1000 RPS
610050
ECI 4076 112
POMME AND SCHWARZ RESG (100.1472.02) 916.625 111
FREQUENCY RANGE 30 TO 330 MHZ
INPUT VOLTAGE 1 MICROVOLT TO 0.1 VOLT
NOISE FIGURE 10 TO 13 DB
OUTPUT 2XIF, 2XAF, RECORDER
610020 RECEIVER, KU-BAND CAC47 11235
610020
610020 RECEIVE KU BAND SIGNAL FOR DETAILED EVALUATION OF PERFORMANCE
610020 CHARACTERISTICS OF TORS/GROUND/SHUTTLE DATA LINK.
610020
610020 CARRIER FREQUENCY 13.4 TO 14.2 GHZ
610020 SENSITIVITY +/- 25 MHZ
610020 CARRIER/NOISE DENSITY RATIO .. 96 DB/MZ
610020 DATA RATE 100 TO 1000 RPS
610020 DATA UP TO VIDEO
610020
610060 RECEIVER, S-BAND CAC48 11111
610060
610060 RECEIVE S BAND SIGNAL FOR DETAILED EVALUATION OF PERFORMANCE
610060 CHARACTERISTICS OF TORS/GROUND/SHUTTLE DATA LINK.
610060
610060 CARRIER FREQUENCY 2025 TO 2120 MHZ
610060 SENSITIVITY +/- 25 MHZ
610060 CARRIER/NOISE DENSITY RATIO .. 85 DB/MZ
610060 DATA RATE 100 TO 1000 RPS
610060 VOICE
610060
POMME AND SCHWARZ RUSU3 (100.1608.02) 810.160 111
100.1614.021
FREQUENCY RANGE 2 TO 5.1 GHZ
VOLTAGE RANGE 25 MICROVOLT TO 80 MILLIVOLT
NOISE VOLTAGE <3.5 MICROVOLT
IMAGE FREQUENCY REJECTION 30 TO 40 DB
OUTPUT 2XIF, 2XAF, RECORDER
RMC ELECTRONICS LAB INC RFR W1900ATR 89000 211
FREQUENCY 1.4 TO 2.3 GHZ FIXED TUNED
IF/DISCRIMINATOR BW 50 MHZ
NOISE FIGURE 11 DB
610070 RECEIVER, LASER ELECTRONICS CAC49 11235
610070
610070 PERFORMS DATA DEFORMATTING FROM LASER CARRIER BEAM
610070
620030 TRANSMITTER, VHF CAC50 11111
620030
620030 PROVIDES MEDIUM POWER OUTPUT FOR TRANSMISSION TO TORS FOR EVAL-
620030 UATION OF PERFORMANCE CHARACTERISTICS.
620030

ORIGINAL PAGE IS
OF POOR QUALITY



620030 CARRIER FREQUENCY 136 TO 144 MHZ
620030 STABILITY +/- 20 KHZ
620030 OUTPUT POWER 10 WATTS MAX
620030 DATA RATE 100 TO 10,000 BPS
620030 SIGNAL BANDWIDTH 0.1 MHZ
620030 DUTY CYCLE CONTINUOUS
620030
ECI 0935 112
FREQUENCY 290 TO 315 MHZ
POWER OUTPUT 60 WATTS
BANDWIDTH UP TO 10 MHZ
SPURIOUS OUTPUTS > 80 DB DGM
AVCIN CO., VECTOR DIV TM-1225 0945 122
POWER OUTPUT 0.3 WATT, NOMINAL
OUTPUT IMPEDANCE 50 OHM
FREQ. RANGE 215 TO 260 MHZ
FREQ. ACCURACY +/- 0.005%
CARRIER DEVIATION +/- 125 KHZ
TYPE MODULATION PHASE
620025 TRANSMITTER, KU BAND CA051 11235
620025 PROVIDE KU BAND SIGNAL FOR EVALUATION OF PERFORMANCE CHARACTER-
620025 ISTICS OF TDPS/GROUND/SHUTTLE DATA LINK.
620025 CARRIER FREQUENCY 14.4 TO 19.35 GHZ
620025 STABILITY +/- 25 MHZ
620025 TRANSMITTER POWER 100 WATTS
620025 EFFECTIVE RADIATED POWER 33 DBM
620025 DATA RATE 1.0 MBPS
620025 DUTY CYCLE VIDEO
620025 SIGNAL BANDWIDTH CONTINUOUS
620025
WATKINS-JOHNSON BWJ-1049 112
COBER ELECTRONICS INC. 0915 122
FREQUENCY 16.0-16.4 GHZ
POWER 150 WATTS
MODE PULSE
PRF 6C TC 4000
PULSE WIDTH 0.2 TO 1.2 MICROSEC.
220041 LASER ASSEMBLY - CO2 CA052 11111
220041 PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFINE
220041 AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMMUNICATIONS
220041 APPLICATIONS. CONTAINS OPTICAL TRANSMITTER INCLUDING BEAMWIDTH
220041 CONTROL, MODULATOR, AND COOLING SYSTEM.
220041 LASER CAVITY
220041 WAVELENGTH 10.6 MICRONS
220041 LASING MODE FUNDAMENTAL (TEM..)
220041 OPTICAL OUTPUT POWER 5 WATTS AVERAGE
220041 EFFICIENCY 5 PERCENT
220041 MODULATOR
220041 TYPE DIGITAL
220041 DATA RATE 200 MBPS
220041 EXTINCTION RATIO 20 DB
220041 COOLING SYSTEM
220041 COOLING METHOD CONDUCTIVE, NO LIQUIDS ALLOWED
220041 HEAD LOAD 200 WATTS
220041
Sylvania Electronics 0940 810,000
OUTPUT 5 WATTS AT 10.6 MICRONS
FREQUENCY STABILITY < 15 MHZ AFTER 1/2 HR WARM-
UP, OVER SEVERAL HOURS
AMPLITUDE STABILITY < 5%
WEIGHT 27 LBS (HEAD), 40 LBS PWR SUPPLY
HUGHES 119
LASER HEAD 695 X 8 X 38 IN., 27 LBS
227400 LASER ASSEMBLY - Nd:YAG CA053 11111
227400 PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFINE
227400 AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMMUNICATIONS
227400 APPLICATIONS. CONTAINS OPTICAL TRANSMITTER INCLUDING BEAMWIDTH
227400 CONTROL, MODULATOR, AND COOLING SYSTEM.
227400 LASER CAVITY
227400 WAVELENGTH 1.06 MICRONS
227400 LASING MODE FUNDAMENTAL (TEM..)
227400 OPTICAL OUTPUT POWER 1 WATT AVERAGE
227400 EFFICIENCY 1 PERCENT
227400 MODULATOR
227400 TYPE DIGITAL
227400 DATA RATE 1 GBPS
227400 EXTINCTION RATIO 2000
227400 COOLING SYSTEM
227400 COOLING METHOD CONDUCTIVE, NO LIQUIDS ALLOWED
227400 HEAD LOAD 200 WATTS
227400
227400 LASER LINK, DOUBLED Nd:YAG CA054 11111
227400 PROVIDE SOURCE OF COHERENT OPTICAL ENERGY TO BE USED TO REFINE
227400 AND EXTEND KNOWLEDGE OF THE USE OF LASERS IN SPACE COMMUNICATIONS
227400 APPLICATIONS. CONTAINS OPTICAL TRANSMITTER INCLUDING BEAMWIDTH
227400 CONTROL, MODULATOR, DOUBLER, AND COOLING SYSTEM.
227400



227400 LASER CAVITY
227400 WAVELENGTH 0.53 MICRONS
227400 LASING MODE FUNDAMENTAL (TEM₀₀)
227400 OPTICAL OUTPUT POWER 0.1 WATTS AVERAGE
227400 EFFICIENCY 0.1 PERCENT
227400 MODULATOR
227400 TYPE DIGITAL
227400 DATA RATE 1 GRPS
227400 EXTINGUISH RATIO 20 DB
227400 COOLING SYSTEM
227400 COOLING METHOD CONDUCTIVE, NO LIQUIDS ALLOWED
227400 HEAD LOAD 200 WATTS
227400
HOLGREN #256 #13030 115
WAVELENGTH53 OR 1.06 MICRONS
MODE MULTI OR TEM₀₀
CW POWER TO 50 WATTS
227400 LASER BEACON, NO:YAG CN055 11111
227400
227400 SERV. AS A LIGHT SOURCE TO FACILITATE ACQUISITION AND/OR TRACKING
227400 FROM A DISTANT COMMUNICATION TERMINAL.
227400
227400 BEACON TYPE Q-SWITCHED COUPLED, MULTIMODE
227400 CAVITY
227400 WAVELENGTH 0.53 MICRONS (PREFERRED)
227400 REP RATE 5 PRS
227400 ENERGY/PULSE 20 MJ AT 0.53 MICRONS
227400 PULSEWIDTH 20 MSEC
227400 PUMP METHOD FLASH PUMP
227400 COOLING TECHNIQUE PULSE POSITION MODULATION
227400 BEAMWIDTH VARIABLE 1 DEG TO 1 ARC MINUTE
227400 LIFETIME 10 MEG FLASHES OR GREATER
227400 POWER CONSUMPTION 100 WATTS
227400
181000 FREQUENCY SYNTHESIZER AND DRIVER CN056 11111
181000
181000 RECEIVER FIRST LOCAL OSCILLATOR
181000
181000 FREQUENCY DC TO 500 MHZ
181000 STABILITY +/- SE-11 PER DEG C
181000 OUTPUT VOLTAGE 1 V RMS +/- 1.5 DB
181000
ROMBE AND SCHWARTZ #500W/5110R #13250 111
FREQUENCY RANGE C - 500 MHZ
OUTPUT LEVEL 0.2 MICROVOLT TO 2 VOLT
SPURIOUS SIGNAL SUPPRESSION > 40 DB
HEWLETT PACKARD #5100W/5110R #13250 111
FREQUENCY DC TO 500MHZ
OUTPUT VOLTAGE 1 V RMS +/- 1 DB 100 MHZ TO 50
MHZ, 1 V RMS +/- 2 DB, +/- 4 DB 50
HZ TO 100 MHZ INTO 50 OHM LOAD
DIGITAL FREQUENCY SELECTION 0.01 HZ TO 10 MHZ PER STEP
618100 NOISE FIGURE TEST SET CN057 11111
618100
618100 MEASURE SYSTEM NOISE FIGURE
618100
618100 NOISE FIGURE RANGE 10 TO 16 DB
618100 ACCURACY +/- 0.5 DB, 0 - 30 DB
618100 INPUT FREQUENCY 100 TO 1000 MHZ
618100 BANDWIDTH 1 MHZ
618100
GENERAL MICROWAVE 551 #1400 111
FREQUENCY RANGE 30, 45, 60, 70, 105, 120, 130,
140, OR 200 MHZ, PLUG-IN UNITS
EXTENDS RANGE TO 1400HZ.
NOISE FIGURE RANGE 0 TO 30 DB
ACCURACY +/- 0.25 DB FOR 0-20 DB
+/- 0.5 DB FOR 20-26 DB
+/- 1.0 DB FOR 26-30 DB
BANDWIDTH 1 MHZ
055050 SCOPE CAMERA CN058 21111
055050
055050 OBTAIN FILM RECORD OF DATA AND DISPLAYS.
055050
055050 REDUCTION RATIO 1:1, 1:0.7
055050 LENS F/1.9 TO F/16
055050 SHUTTER 8 TO 5 SECONDS
055050 FILM TYPE POLAROID
055050
TERTRONIX INC #C-27-P #590 111
SHUTTER MECHANICAL (ELECTRICAL OPTION)
FILM RACK PACK-FILM (ROLL FILM AND 4X5
GRAFLEX OPTIONAL)
SHUTTER SPEED 1 TO 1/50 SEC, BULB AND TIME
LENS F/1.9, 0.85 MAGNIFICATION
TERTRONIX INC #C-53-P #825 111
SHUTTER ELECTRIC
FILM RACK PACK-FILM (ROLL FILM AND 4X5
GRAFLEX OPTIONAL)
SHUTTER SPEED 4 TO 1/60 SEC, BULB AND TIME
LENS F/1.9, 0.85 MAGNIFICATION
HEWLETT PACKARD #147A #675 111
REDUCTION RATIO 1:1 TO 1:0.7 (ADJUSTABLE)
LENS 75 MM, F/1.9 HIGH TRANSMISSION
LENS; APERTURE RANGES F/1.9 TO
F/16

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

SHUTTER SPEEDS 1/30, 1/15, 1/8, 1/4, 1/2, 1,
2, 4, TIME, BULB
CAMERA BACK POLAROID USING TYPE 107

TEXTRONIX INC 8C-12 8590
LENS 75 MM
STOP F/1.4 TO F/16
MAGNIFICATION 0.85
LENS SPEED 1 TO 1/100 SEC MECH
4 TO 1/60 SEC ELEC
FILM TYPE POLAROID

TEXTRONIX 8C30 8525
MAXIMUM RELATIVE APERTURE F/1.4
MAGNIFICATION 0.7 TO 1.5
RELATIVE SPEED 1.0
FIELD OF VIEW 3.15X3.93 IN (8X10 CM)
FILM TYPE POLAROID

HEWLETT-PACKARD 8198 8420
MAXIMUM RELATIVE APERTURE F/3.5
MAGNIFICATION 0.85
SPEED 8 TO 1/60 SEC
LENS 75 MM
FILM TYPE POLAROID

HEWLETT-PACKARD 8195A 81025
MAXIMUM RELATIVE APERTURE F/1.3
MAGNIFICATION 0.9
SPEED 8 TO 1/30 SEC
LENS 80 MM
FILM TYPE POLAROID

179500 FREQUENCY COUNTER CA059 11111
179500
179500 MEASURE FREQUENCY, PERIOD, MULTIPLE PERIOD AVERAGE, TIME INTERVAL
179500 AND RATIO OF FREQUENCY SOURCES.
179500
179500 BANDWIDTH 3 HZ TO 10 MHZ
179500 SENSITIVITY 0.01 RMS
179500 TIME BASE INTERNATIONAL SECOND
179500

JOHN FLUKE MFG CO INC 81952A 8695 111
FREQUENCY RANGE DC TO 80 MHZ
DISPLAY 7 DIGIT LED (8 OR 9 OPTIONAL)
SENSITIVITY 50 MV RMS DC TO 30 MHZ INCREASE
ING TO 75 MV AT 80 MHZ
PERIOD RANGE DC TO 10 MHZ (DC COUPLED)
5 HZ TO 10 MHZ (AC COUPLED)
TIME INTERVAL RANGE 0.1 MICROSEC TO 1E+6 SEC

HEWLETT PACKARD 85340A 85600 111
BOMBE AND SCHWABZ 8FET21100.8C39.02/ 111
110.8125.02/
100.8089.02/
100.8051.02/
100.8088.02/
FREQUENCY RANGE 0 - 100 MHZ
SENSITIVITY -5 TO +5 VOLTS
PERIOD RANGE 0.5 MICROSEC TO 100 SEC
TIME INTERVAL RANGE 0.2 MICROSEC TO 1E+6 SEC

408400 RF POWER METER CA060 11111
408400
408400 MEASURE RF POWER OF VARIOUS RF SOURCES.
408400
408400 FREQUENCY RANGES 126 TO 130 MHZ
130 TO 144 MHZ
2025 TO 2120 MHZ
13.4 TO 14.2 GHZ
14.4 TO 15.35 GHZ
408400 ACCURACY 0.02 DB/10 DB
408400

BOMBE AND SCHWABZ 8NRS (100.2433.52/
100.2440.XX) 81945 111
FREQUENCY RANGE 0 TO 15 GHZ
POWER RANGE 0.1 - 350 MW

GENERAL MICROWAVE CORP 8471 81675 111
FREQUENCY RANGE 10 MHZ TO 40 GHZ
ACCURACY +/- 0.5% OF READING +/- 1 COUNT
POWER RANGE 10 NANOWATTS TO 3 WATTS
READOUT DIGITAL

BEINCHSEL 8PB-18 111
RANGE 0.1 TO 45 MB
ACCURACY +/- (0.1% + 1 MICRO-WATT)
BIAS POWER SELF BALANCING 1E+5 MIN (EFFECTIVE LOOP GAIN)
1.0 MICROWATT/DIV INDICATOR
SENSITIVITY
0.2 MICROWATT (INDICATOR RESOLU-
TION)

332000 AC/DC VOLTMETER CM061 21111
332000
332000 MEASURE THE VOLTAGE LEVEL OF VARIOUS SIGNAL SOURCES.
332000



932000	VOLTAGE RANGE	+/- 100 MV; 1, 10, 100, 1000 V	
932000	FREQUENCY RANGE	45 HZ TO 1 MHZ	
932000	ACCURACY	+/- 0.01 % OVERALL	
932000			
POHDE AND SCHWARZ	#URU (100.0060.02)	#1210	111
	FREQUENCY RANGE	DC, 10 HZ TO 1.5 GHZ	
	VOLTAGE RANGE	5 MV TO 1 KV (DC)	
		0.1 TO 1000 V (AC)	
	ACCURACY	+/- 2% (DC), +/- 3% (AC)	
FLUKE	#R200A	#995	111
	VOLTAGE RANGE	100MICROVOLTS TO 1200VOLTS DC	
		1 MV TO 1100V AC	
	RESISTANCE	10 MIL OHM TO 16 MEGOHMS	
	READOUT	NEON TUBE DISPLAY	
		PRINTER OUTPUT OPTICAL	
	ACCURACY	+/- .001 % OF INPUT	
	RESOLUTION01 PER CENT OF RANGE	
	POWER	115/230 V 50-60HZ 25W	
	DIMENSIONS	19.7X24.1X41.4CM	
	WEIGHT	6.75KG	
HEWLETT PACKARD	#3450B WITH OPTIONS	#4675	111
	001, C04		
	DISPLAY	6 DIGITS	
	DC VOLTAGE RANGE	+/- 100 MV, +/- 1, 10, 100 AND 1000 V	
	RANGE SELECTION	MANUAL OR AUTOMATIC	
	FREQUENCY RANGE	45 HZ TO 1 MHZ	
	AC VOLTAGE RANGE	1, 10, 100, 1000 V RMS	
103500	A/D CONVERTER	CN062	21111
103500			
103500	CONVERT INPUT SIGNAL FROM ANALOG TO DIGITAL FORMAT FOR 1/G UNIT		
103500			
ANALOG DEVICES INC.	#ADC-85	#79	111
	RESOLUTION	8 BITS	
	ERROR	+/-0.2%	
	CONVERSION TIME	1 MILLISEC.	
ANALOG DEVICES INC.	#ADC-120Z	#305	111
	RESOLUTION	12 BITS	
	ERROR	+/-0.0125%	
	CONVERSION TIME	25MICROSEC.	
ANALOG DEVICES INC.	#ADC-160	#1350	111
	RESOLUTION	16BITS	
	ERROR	+/-0.0015%	
	CONVERSION TIME	40MICROSEC.	
TENNELEC	#TC500	#6000	
	CONVERSION TIME3 MICRO SEC	
	CAPACITY	2048 CHANNELS	
	CONFIGURATION	COMPATIBLE WITH NIM	
TENNELEC	#TC 520	#750	111
	CHANNELS	8	
	CHANNEL DWELL TIME	300MSEC	
	CONFIGURATION	NIM COMPATIBLE	
POHDE AND SCHWARZ	#UCH 210/1		111
	(221.5603.02)		
	ANALOG INPUT RANGE	+12 TO -12 VOLTS	
	SAMPLING RATE	> 10000 MEAS/SEC (POSITIVE)	
		> 17000 MEAS/SEC (NEGATIVE)	
	DATA OUTPUT	3+1 DECADES (BCD CODE PARALLEL)	
		0-1200 PULSES (SERIAL)	
055055	CAMERA, CINE	CN063	11111
055055			
055055	PROVIDE PHOTOGRAPHIC RECORD OF TELESCOPE FIELD-OF-VIEW.		
055055			
055055	FILM FORMAT	35 MM	
055055	OPERATION	REMOTE	
055055	FRAME RATE	16 FPS	
055055	FILM TYPE	COLOP	
055055	RESOLUTION	50 LINES/MM	
055055	LENS	100 MM	
055055	FIELD OF VIEW	40 DEGREES	
J A MAURER INC	#317-E1	#4500	311
	MODES	8, 16, 32, 64 FPS ELECT CCNT	
	SHUTTER	ROTARY	
	SHUTTER SPEEDS	1/62.5 TO 1/500 SEC	
	FILM	16 MM MAGAZINE	
	FILM CAPACITY	150 FT THIN BASE	
	LENS	25 MM, F/1.4	
	RESOLUTION	236 L/IN ON AXIS, 150 L/PP AVG	
J A MAURER INC	#308	#15 - 20,000	511
	FILM	16 MM DOUBLE OR SINGLE PERFORATED	
	FORMAT	0.295 X 0.404 INCH	
	FRAME RATE	1, 6, 12 AND 24 FR/SEC. TIME	
	SHUTTER SPEED	TIME, 1/60, 1/125, 1/250, 1/500 AND 1/1000 SEC	
	FILM CAPACITY	90 FT STD-BASE FILM, 145-165 FT THIN-BASE FILM	

ORIGINAL PAGE IS
OF POOR QUALITY



FLCWN ON APOLLO

603000 OPTICAL ANTENNA SERVO ELECTRONICS CNO66 11114
603000
603000 SUPPORT ELECTRONICS FOR OPTICAL ANTENNA
603000
610010 SWEPTBAND RECEIVER AND DEMODULATOR CNO57 41111
610010
610010 RECEIVE SIGNALS FROM TERRESTRIAL NOISE AND TERRESTRIAL INTERFER-
610010 ENCE SOURCES AS A FUNCTION OF TIME OF DAY AND SEASON OF YEAR.
610010
610010 FREQUENCY RANGE 100 TO 1000 MHZ SWEEP MODE OR
610010 TUNABLE
610010 RECEIVER TYPE SUPERHETERODYNE, LINEAR, SWEPT
610010 FREQUENCY, AGC
610010 SENSITIVITY -120 DBM
610010 PRE-DETECTION BW ADJUSTABLE
610010 DYNAMIC RANGE 60 DB INSTANTANEOUS
610010 120 DB TOTAL
610010 SWEEP RATE TWO RECEIVER BW PER MINUTE
610010 RECEIVER NOISE FACTOR 4.0 DB
610010 BANDWIDTH 1 KHZ
610010
SINGER 04M-37/57 912,000 311
A PROGRAMMABLE PRECISION METER FOR MEASUREMENT OF CONDUCTED
OR RADIATIVE RF INTERFERENCE. MEETS MIL-STD-461/026.
FREQUENCY BANDWIDTH 30MHZ TO 1GHZ
DYNAMIC RANGE 140DB (0.1MICROVOLT-1VCLT)
PROGRAMMABLE FUNCTIONS BAND SELECTION, FREQUENCY
TUNING, DETECTOR FUNCTION (PEAK
50MS, PEAK 300MS, PEAK 3SEC,
QUASI-PEAK AND FIELD INTENSITY)
BANDWIDTH SELECTION, CONTIN-
UOUS IF GAIN, RF ATTENUATOR.
276000 ATTENUATOR CALIBRATOR CNO68 41111
276000
276000 DYNAMIC RANGE ADJUSTMENT, REMOTE ELECTRICAL CONTROL, DIGITAL READ
276000 OUT OF VALUE, MANUAL OVER-RIDE.
276000
276000 IMPEDANCE 50 OHMS
276000 FREQUENCY DC TO 1.5 MHZ
276000 RANGE 0 TO 60 DB
276000
104360 SCAN PROGRAM GENERATOR CNO69 31111
104360
104360 GENERATES DIGITAL CONTROL SIGNALS TO CONTROL RECEIVER SWEEP.
104360 SUPPLIES DIGITAL READOUT OF FREQUENCY AND LEVEL.
104360
104360 BANDS 24
104360 SCAN TIME 0.03 TO 300 SECONDS
104360 CAPABILITY ANALOG DATA COLLECT & DISPLAY
104360
SINGER 0P-7 89,000 311
SELECTS ANY RECEIVER COMBINATION
AND FOLLOWING RECEIVER FUNCS.. BAND(24 BANDS); BAND SEGMENT;
BANDWIDTH(NARROW, MEDIUM, WIDE);
1 DETECTOR MODE (9 MODES); SCAN
TIME (VARIABLE .03 TO 300 SEC).
BANDS SCANNED 1 TO 24 IN ANY COMBINATION
OUTPUT DATA CHANNELS 7-AMPLITUDE, FREQUENCY AND LOG
VIDEO
219700 SIGNAL FORMATTER CNO70 21111
219700
219700 PROVIDES DATA FORMATTING, ACCUMULATES PRE/POST AMBLE, CLOCKS
219700 START OF EACH TEST AND DATA ACQUISITION.
219700
219700 DATA FORMAT 8, 16, 32, BITS/WORD
219700 DATA RATE 129000 WORDS PER SECOND
219700
219700
TEKTRONIX INC #40024 111
DISPLAY MEDIUM 11 IN DIA CRT
DISPLAY AREA 9.9IN.HOR.X6.1IN.VERT.
ALPHANUMERIC MODE:
FORMAT 30 LINES OF 85 NORMAL CHARACT.
CHARACTER SET 96 UPPER & LOWER CASE
CHARACTER SIZE 70X85(USIAN RE TWICE SIZE)
CHARACTER GENERATION 7X9 DOT MATRIX
CURSOR PULSATING 7X9 MATRIX
GRAPHIC MODES LINEAR INTERPOLATE, INCREMENTAL
PLOT, POINT PLOT, 1024X1024 AD-
RESSABLE POINTS, 1024X768 VIEW-
ABLE POINTS
GRAPHIC INPUT MODE 1024(X), 768(Y), JOYSTICK CONTROLLED,
CROSSHAIR CURSOR
ACCUMULATORS 4
CYCLE TIME 300 NANOSEC
CORE MEMORY 32 K
409900 RF VARIABLE POWER SUPPLY CNO71 21111
409900
409900 PROVIDES SOURCE OF VARIABLE POWER.
409900
409900 OUTPUT 100 V PEAK
409900 VOLTAGE GAIN 20 X
409900 DISTORTION 0.1 B THD
409900 FREQUENCY RANGE DC TO 30 KHZ
409900



JOHN FLUKE MFG CO INC #4270A \$1895 111
PROG CODING RCD (PRIMARY AVAILABLE)
OUTPUT VOLTS 0 TO +/- 100 VOLTS DC/PEAK AC
OUTPUT CURRENT 0.5 AMPS
FREQUENCY RANGE DC TO 30 KHZ
RESOLUTION 1 MV
RESOLUTION 0.001%

407000 POWER CALIBRATION UNIT CN072 21111
407000
407000 PROVIDES KNOWN INCREMENTALLY CONTROLLABLE POWER LEVEL FOR
407000 STANDARDIZATION.
407000
407000 OUTPUT -100 TO +100 V
407000 LOAD REGULATION 1 MV + 0.01 %
407000 LINE REGULATION 10 MV + 0.01 %
407000

JOHN FLUKE MFG CO INC #427CA \$1895 111
PROG CODING RCD (PRIMARY AVAILABLE)
OUTPUT VOLTS 0 TO +/- 100 VOLTS DC/PEAK AC
OUTPUT CURRENT 0.5 AMPS
FREQUENCY RANGE DC TO 30 KHZ
RESOLUTION 1 MV
RESOLUTION 0.001%

JOHN FLUKE MFG CO INC #332R \$2445 111
OUTPUT VOLTAGE 0 TO 1111.1110 VDC
OUTPUT CURRENT 50 MA
REGULATION 0.0002%

620002 TRACKER ELECTRONICS - FINE CN073 21114
620002
620002 PROVIDES ACCURATE TRACKING OF AN OPTICAL BEACON. PROVIDES CONTROL
620002 SIGNALS TO THE TELESCOPE GIMBALS AND VERNIER BEAM DEFLECTOR.
620002
620002 TYPE IMAGE DISSECTOR TRACKER
620002 FOV +/- 300 MICRORAD
620002 ACCURACY 0.5 MICRORAD
620002 SENSITIVITY TO 1.06 MICRONS OR 10.6 MICRONS
620002
110300 BIT ERROR COUNTER CN074 11111
110300
110300 COUNT ERROR PULSES OUTPUT BY THE BIT ERROR DETECTOR TO DETER-
110300 MINE THE BIT ERROR RATE OF THE DIGITAL TRANSMISSION SYSTEM.
110300

HENLETT PACKARD #5376A \$1350 111
603003 SWITCH/DIPEX/PREAMP UNIT CN077 11134
603003
603003 FOR SWITCHING ANTENNAS FOR DIFFERENT ALTITUDES, USE FFT WFI
603003 CANCELLATION, FORMING DMMT PATTERN ETC
603003
219900 PCDEM CN078 11114
219900
219900 MODULATE AND DEMODULATE VOICE OR DATA TRANSMISSIONS WITH OR WITH-
219900 OUT SPREAD SPECTRUM CAPABILITY. PROVIDES BASELINE COMMUNICATION
219900 SYSTEM, ANALOG OR DIGITAL VOICE, AND DIGITAL DATA.
219900
219900 DATA 75 HS TO 50 KHS (VARIABLE)
219900 MOVEABLE SS
219900 BANDWIDTH 2 MHZ
219900
219900 DEMODULATOR CN079 31111
219900
219900 DEMODULATE VOICE OR DATA TRANSMISSIONS WITH OR WITHOUT SPREAD
219900 SPECTRUM CAPABILITY. PROVIDES A NUMBER OF RECEIVERS CONNECTED TO
219900 DIFFERENT ANTENNAS.
219900
219900 MODEM, WIDEBAND CN080 21235
219900
219900 MODULATE AND DEMODULATE WIDEBAND FM OR DIGITAL SIGNALS. WIDEBAND
219900 FM DETECTOR WITH THRESHOLD EXTENSION. DIGITAL DEMOD WITH BIT.
219900 SYNC, AND ERROR CORRECTION CODING. MODULATOR PORTION HAS COMPLE-
219900 MENTARY CHARACTERISTICS.
219900

103500 D/A AND A/D CONVERTER CN081 11111
103500
103500 PROVIDES FOR A/D AND D/A CONVERSION.
103500
103500 A/D CONVERTER
103500 SAMPLING RATES 15 MEGASAMPLES/SEC
103500 DATA RATES 105 MB/S
103500 D/A CONVERTER
103500 BANDWIDTH 6 MHZ
103500

TENNELEC #TC500 \$6000 111
CONVERSION TIME 3 MICRO SEC
CAPACITY 2048 CHANNELS
CONFIGURATION COMPATIBLE WITH NIM

TENNELEC #TC 520 \$750 111
CHANNELS 8
CHANNEL DWELL TIME 300NSEC
CONFIGURATION NIM COMPATIBLE

603004 CENTRAL UNIT - ANTENNA SCAN CN082 11114
603004

ORIGINAL PAGE IS
OF POOR QUALITY

603004 PROVIDES VARIABLE ANGULAR SWEEP PATTERNS (SPIRAL RECTANGULAR
603004 SCANS, ETC.) OF VARIABLE ANGULAR COVERAGE.
603004
103920 DATA BIT STREAM GENERATOR CAC83 11111
103920
103920 PROVIDES TEST BIT STREAMS FOR TRANSMISSION.
103920
103920 DATA PATES 1 TO 109 MBS
103920
RESEARCH INC #812-4A 82530 111
ANALOG INPUT -10.240 TO +10.235 V
INPUT IMPEDANCE 100 MEG OHM
INPUT OVERVOLTAGE PROTECTION ... +/- 100 V
FULL SCALE ACCURACY +/- 0.025 %
RESOLUTION 12 BITS BINARY
CONVERSION SPEED 5 KHZ
LOGIC LEVELS 5 +/- 1.5 V (LOGIC 1)
0 + .4, - .0 V (LOGIC 0)
620003 COARSE TRACKER CNO84 11114
620003
620003 PROVIDES INITIAL GROUND BEACON ACQUISITION AND COARSE TRACKING.
620003
620003 DYNAMIC (SEARCH) FOV 1 DEG BY 1 DEG
620003 ACCURACY (1 SIGMA) +/- 200 MICRORAD
620003 TRACKING BANDWIDTH 1 MZ
620003 SEARCH TIME 10 SEC MAX
620003 SENSITIVITY TO .53 MICRONS
620003 FALSE TARGET REJECTION MUST REJECT MOON, PLANETS,
620003 STARS, AND SUNLIT CLOUDS.
620003
603002 LASER TRANSMITTER ELECTRONICS CNO85 11235
603002
603002 DATA FORMATTING FOR TRANSMISSION VIA LASER CARRIER
603002
225700 OPTICAL COLLIMATOR CNO86 11111
225700
225700 PROVIDES FOR OPTICAL ALIGNMENT OF THE LASER AND TELESCOPE. ALIGN-
225700 MENT CHECKS TO BE BOTH CONTINUOUS AND PULSED.
225700
TROPEL #280 82990 111
SPECTRAL RANGE45 TO 1.10 MICRONS
OPTICAL TRANSMISSION 92 PER CENT
TROPEL #290 81350 225
SPECTRAL RANGE 10.6 MICRONS
DAVICSON #060C 81582 112
407000 LASER POWER SUPPLY CNO87 21111
407000
407000 CONVERTS AND REGULATES RAW SPACECRAFT POWER IN A FASHION SUITABLE
407000 FOR PUMPING A LASER. PRIME LASER POWER CONDITIONER AND DISTRIBU-
407000 TION FOR ALL LASERS.
407000
407000 REGULATION +/- 0.1 %
407000 POWER 100 WATTS MAX
407000 TRANSIENT RECOVERY 0.1 SECONDS
407000
225700 BEAM EXPANDER OPTICS CNO88 21111
225700
225700 PROVIDES VARIABLE BEAM DIVERGENCE CONTROL FOR LASER TRANSMITTER
225700 AND THE OPTICAL BEACON TO FACILITATE ACQUISITION.
225700
225700 TYPE AFDCAL WITH VARIABLE ELEMENT
225700 SEPARATION.
225700 INPUT BEAM SIZE 0.25 INCH DIA
225700 EXIT BEAM SIZE 1 INCH DIA
225700 INPUT BEAM CHARACTERISTICS
225700 LASER TRANSMITTER TEM WAVE, GAUSSIAN DIST
225700 OPTICAL BEACON 3 MRAD DIVERGENCE, GAUSSIAN
225700 DIST
225700 OUTPUT BEAM DIVERGENCE
225700 CO2 LASER TRANSMITTER DIFF LIMIT TO 200 ARC-SEC
225700 NDYAG LASER TRANSMITTER .. DIFF LIMIT TO 50 ARC-SEC
225700 OPTICAL BEACON 2 ARC-MIN TO 1 DEGREE
225700 ADJUSTMENT RESPONSE TIME 0.5 SEC MAX
225700
225700 BEAM DEFLECTOR CNO89 41111
225700
225700 ONE TYPE PROVIDES A VERNIER POINTING CONTROL WITH FASTER RESPONSE
225700 AND GREATER PRECISION THAN THE MAIN OPTICS GIMBALS. THE OTHER
225700 TYPE PROVIDES A SMALL OFFSET POINTING CAPABILITY BETWEEN THE
225700 TRACKER AXIS AND THE TRANSMITTER AXIS.
225700
225700 DEFLECTION CAPABILITY TWO AXIS
225700 DEFLECTION RANGE +/- 200 MICRORAD
225700 PRECISION +/- 0.2 MICRORAD
225700 SPECTRAL RANGE ACCOMMODATE LASER TRANSMITTER,
225700 RECEIVER AND BEACON.
225700
605400 LASER POWER METER CNO90 11111
605400
605400 MEASURE OUTPUT LASER POWER.
605400
605400 WAVELENGTH 0.53, 1.06 AND 10.6 MICRONS
605400 ACCURACY +/- 5 PERCENT
605400 RANGE 300 NM



409400
HADRCH #99 #1050 111
WAVELENGTH..... 0.3 TO 11 MICRONS
ACCURACY..... +/- 5 PER CENT
TYPE..... CONE
MAX POWER..... 5J
SENSITIVITY..... 450 MU V/J
RISE TIME..... 5 SEC

OPTICS TECHNOLOGY #615 #445 111
MOUNT..... BENCH
SPECTRAL RANGE..... 400 TO 1150NM
DYNAMIC RANGE..... 8 RANGES WITH FULL SCALE
READINGS OF 0.2,1,3,10,30,100,
300 AND 1000 MILLIWATTS
ACCURACY..... +/-5%

YELLOW SPRINGS INSTRUMENT #54 #395 111
RANGE..... 0 TO 20 DPM
OPERATING TEMPERATURE..... -28 TO 113F (-2 TO 45C)
BATTERY OPERATED

101000 CALIBRATION SIGNAL GENERATOR CN091 11111
101000
101000 CALIBRATES RECEIVER SENSITIVITY AND PHASE METER.
101000
101000 FREQUENCY RANGE 1500.75 TO 1570.25 MHZ
101000 FREQUENCY ACCURACY +/- 5 MHZ
101000 FREQUENCY STABILITY 0.005 %/DEGREE C
101000 RF OUTPUT POWER +10 DB TO -127 DB
101000
RHODE AND SCHWARZ #5NA1 (100.4594,021) #9675 111
FREQUENCY RANGE 0.5 - 1.8 GHZ
FREQUENCY ERROR +/- 0.5%
OUTPUT LEVEL +10 DBM TO -130 DBM

LEGIMETRICS INC #PER SPEC #7C00 11-
FREQUENCY 1500 TO 2000 MHZ
RESOLUTION 1 HZ OR 1 KHZ (OPTIONAL)
RF OUTPUT +10 DBM (LEVEL)
+/- 1 DB (LEVELING)
SPURIOUS OUTPUTS 25 DB BELOW CARRIER (HARMONICS)
80 DB BELOW CARRIER (AC-HARM)

HEWLETT PACKARD #8E14A #2570 111
FREQUENCY RANGE 800 TO 2400 MHZ
OUTPUT POWER +10 DBM TO -127 DBM INTO 50-DBM
FREQUENCY ACCURACY +/- 5 MHZ
FREQUENCY STABILITY 0.005 %/DEG C
< 0.003 % FOR +/- 10 % LINE
VOLTAGE VARIATION

006185 OPEN R.F. AMPLIFIER T.W.T CN092 51111
006185
610040 RECEIVER - NAVIGATION CN096 11111
610040
610040 RECEIVE SIGNALS FROM DIRECT OPTICAL SENSING OF STARS TO INVESTI-
610040 GATE THE PRECISION WITH WHICH A SATELLITE CAN BE POINTED VERTIC-
610040 ALLY USING STARS AS THE SOLE INERTIAL REFERENCE.
610040
610040 CARRIER FREQUENCY 400 MHZ
610040 SENSITIVITY 2.0 MV FOR 10 DB S+N/R
610040 DATA RATE 2000 BPS
610040
RHODE AND SCHWARZ #ESU (100.1195,021) #7980 111
ADD#PLUC IN-4 #6000 111
FREQUENCY RANGE 100 - 470 MHZ
INPUT VOLTAGE 1 MICROVOLT TO 1 VOLT
NOISE FIGURE <12 DB
OUTPUT 2XIF, 2XAF, RECORDER

179500 SPECTRUM ANALYZER/OSCILLOSCOPE CA900 21111
179500
179500 SINGLE DESIGN FUNCTIONS AS POWER SPECTRAL DENSITY AND WAVEFORM
179500 DISPLAY FOR QUICKLOOK AND AS MODULATION MONITOR.
179500
179500 BANDWIDTH 1R GHZ
179500 SENSITIVITY 2 MV/DIVISION
179500 CHANNELS 2
179500 PHOTOGRAPHIC CAPABILITY YES
179500
TERTRONIX INC #RT613/7A26 PLUG-IN #10150 111
/7L13 PLUG-IN
FREQUENCY RANGE 10 KHZ TO 1.8 GHZ
FREQUENCY SPAN 200 HZ/DIV TO 100 MHZ/DIV
RESOLUTION 30 HZ TO 3 MHZ

TERTRONIX #P556 #4200 111
BANDWIDTH DC TO > 50 MHZ
RISETIME < 7 NSEC
TYPE DUAL-BEAM

TERTRONIX INC #E491 #4995 112
FREQUENCY RANGE 10 MHZ TO 40 GHZ
CALIBRATED DISPERSION 1 KHZ/DIV TO 10 MHZ/DIV
RESOLUTION BANDWIDTH 1 KHZ TO 100 KHZ
INCIDENTAL FM <300 HZ AT FUNDAMENTAL

ORIGINAL PAGE IS
OF POOR QUALITY



Space Division
Rockwell International

TEKTRONIX INC	#1401A-4	93725	111
FREQUENCY RANGE	1 MHZ TO 900 MHZ		
INTERMODULATION DISTORTION	< 60 DB FULL SCREEN		
102700 COMPUTER, GENERAL PURPOSE		CN901	21111
102700			
102700 COMPUTER CAPABILITY TO COMMAND EXPERIMENTS ON AND OFF, PROGRAM			
102700 EXECUTION OF EXPERIMENT PROCEDURAL STEPS, COMPUTE GENERAL			
102700 ANGLES FOR TELESCOPE AND ANTENNA POINTING, AND COMPUTATION OF			
102700 POINT-TO-POINT ANGLES.			
102700			
CATA GENERAL CORP	#NOVA 800		111
WORD LENGTH	16 BIT		
CYCLE TIME	800 NANOSEC.		
CORE MEMORY	32K MAX		
CATA GENERAL CORP	#SUPERNOVA 3C		111
WORD LENGTH	16 BIT		
CATA GENERAL CORP	#NOVA 1200		111
WORD LENGTH	16 BIT		
CYCLE TIME	1200 NANOSEC.		
CORE MEMORY	32K MAX.		
VARIAN DATA MACHINES	#920/1		111
MEMORY CYCLE TIME	1.5 MICROSEC.		
MEMORY	EXPANDABLE 4096 BYTES (8 BITS TO 32,768 BYTES)		
REGISTERS	12		
OPERAND PRECISION	UP TO 32 BITS		
VARIAN DATA MACHINES	#R-620/1		111
BASIC COMMANDS	OVER 100		
ADDRESSING MODES	8		
MAX. WORDS	32,768		
WORD LENGTH	16 OR 18 BIT		
REGISTERS	9		
NOTE:	HAS FULL RANGE INTERFACE HARDWARE		
CLARY DATACOMP SYSTEMS INC	#404		111
SIMULTANEOUS TERMINALS	16		
INPUT	TELETYPE/WRITER KEYBOARD		
ACCUMULATOR	64 BIT		
REGISTERS	16 BIT INDEX (2 EACH)		
WORD LENGTH	16, 32, 48, 64 BITS		
MEMORY CAPACITY	1024 16 BIT WORDS, 4096 16 BIT WORDS, OR ADDITIONAL 4096 16 BIT WORDS, TOTAL 85936 WORDS		
GENERAL AUTOMATION INC.	#18-30		111
FUNCTION	SUPERVISE SMALLER COMPUTERS		
MEMORY	8K CORE		
CYCLE TIME	1.2 MICROSEC.		
GENERAL AUTOMATION INC.	#SPC-16		111
MEMORY	16K		
WORD LENGTH	16 BITS		
READ/WRITE CYCLE MEMORY TIME	800 TO 1440 NANOSEC.		
READ ONLY MEMORY	400 TO 720 NANOSEC.		
INPUT/OUTPUT TIME	1.4 TO 2.4 MICROSEC.		
DMA TRANSFER RATE	0.094E+06 TO 2.5E+06		
GENERAL AUTOMATION INC.	#SPC-12		111
MEMORY	4K TO 16K		
WORD LENGTH	8 BITS		
CYCLE TIME	2.16 MICROSEC		
STORED PROGRAM EXECUTION RATE	0.23E+06/SEC		
I/O TRANSFER RATE	0.40E+06/SEC		
REGISTERS	6 12 BIT REGISTERS		
ACCUMULATORS	4 12 BIT		
DIGITAL EQUIPMENT CORP	#PDP-8/E	922000	111
MEMORY	4096 CORE (EXPANDABLE TO 32,768 WORDS)		
WORD	12 BIT		
CYCLE TIME	1.5 MICROSEC.		
ADDRESS REFERENCE	8, 16, 24, OR 32 BIT LEVEL		
NOTE:	HAS FULL RANGE INTERFACE HARDWARE		
400300 TIMER, PRECISION CLOCK		CN902	11111
400300			
400300 PROVIDES PRECISE TIME REFERENCE, MEASURE TIME AT WHICH PHASE MEAS			
400300 OCCURS.			
400300			
400300 DISPLAY	24 BIT DIGITAL CLOCK WITH 1 HR. TIME PERIOD		
400300			
400300 DISPLAY TYPE	DIGITAL		
400300 ACCURACY	1E-6		
400300			
TENNELC	#TCS45	9550	111
COUNT RATE	20MHZ		
TIMERASE	0.1 OR 0.01 SEC		
ACCURACY	SAME AS LINE FREQUENCY		
CONFIGURATION	NIM COMPATIBLE		



ROME AND SCHWARTZ #CAD (100.6597.91) 85400 111
 DISPLAY 6 DIGITS: HR, MIN, SEC
 OUTPUT 1-OUT-OF-N AND/OR ACD CODE
 PROGRAMMER TIME PROGRAMS/COUNT
 DOWN) AND TIME SIGNALS

DATATRON #335C-506 81170 111
 INPUT FREQUENCY 60HZ
 OUTPUTS VISUAL, ACD
 DISPLAY HRS, MIN, SEC

104160 TELEPRINTER CN903 11111
 104160
 104160 PROVIDES CAPABILITY TO SEND OR RECEIVE TELETYPE MESSAGES.
 104160
 104160 CHARACTER RATE 30 CHARACTERS/SEC (150 BPS)
 104160 CONSISTS OF SOLID STATE KEYBOARD
 104160 SOLID STATE THERMAL PRINTER
 104160 EXPANDABLE MESSAGE MEMORY
 104160

VARIAN DATA MACHINES #ASR 33 L ASR 35 111
 SPEED 10 CHARACTERS/SEC
 PAPER CONTINUOUS ROLL 8.5" WIDE
 TYPING LINE 85 CHARACTERS

101600 PHASE METER, DIGITAL CA904 11111
 101600
 101600 MEASURES RELATIVE PHASE ANGLES OF RECEIVED SIGNALS FOR INTERFER-
 101600 METER NAVIGATION/SURVEILLANCE EXPERIMENT.
 101600
 101600 OUTPUT GENERATES 13 BIT PHASE MEAS
 101600 BASED ON TWO 60 MHZ SINUSOIDS
 101600 AT 10 DBM
 101600 FREQUENCY RANGE 400 MHZ
 101600 RESOLUTION 25 KHZ
 101600 ACCURACY +/- 2.5 KHZ
 101600

ROME AND SCHWARTZ #ZSK (100.1450.02) 810500 111
 WRITING SPEED 1.50 M/SEC
 DEFLECTION FACTOR 10 MICROVOLT/CM TO 10 V/CM
 WRITING METHOD ELECTRIC PEN

HEWLETT PACKARD #7035A 5985 111
 INPUT RANGE 1, 10, 100 MV/IN; 1 AND 10 V/IN
 WRITING METHOD ELECTRIC PEN
 ACCURACY +/- 0.2 % FS
 LINEARITY +/- 0.1 % FS

424000 RECORDER, DIGITAL CN905 11111
 424000
 424000 ACCURATE RECORDING OF RAW DIGITAL DATA
 424000
 424000 SIGNAL TYPE 0 TO 5 VDC DIGITAL
 424000 CHANNELS 9
 424000 TAPE SPEED VARIABLE
 424000 DATA RATE 1.5E+06BPS
 424000

AMPER CORPORATION #AR-200 121
 BANDWIDTH 100 HZ TO 3125 HZ THROUGH
 300 HZ TO 250 KHZ
 TAPE WIDTH 0.5 INCH
 TAPE SPEED 1.875 TO 60 IPS
 RECORDING TIME 8 MINUTES TO 4 HRS AND 16 MINS
 FORMAT DIGITAL
 TRACKS 8 DIGITAL, 7 ANALOG

SANGAM ELECTRIC #SABER 111 826,500 111
 TRACKS 14
 TAPE WIDTH 1 INCH
 TAPE SPEEDS 8 SELECTABLE SPEEDS FROM 15/16
 TO 120 IPS
 FREQUENCY RESPONSE 400 HZ TO 2.0 MHZ
 RECORDING RATE 600 KIPS AT 120 IPS SERIAL MODE
 WEIGHT 100 POUNDS
 DIGI-DATA CORPORATION #1600 82650 111
 TAPE SPEED 25, 18.75, 12.5 IPS
 TRACKS 7 OR 9 TRACK
 DATA DENSITY 1600 CPI PHASE ENCODED
 200, 556, 800 CPI NRZI
 TAPE 0.5 INCH, 1.5 MIL, 1200 FEET

DIGI-DATA CORPORATION #1700/POP-11/7-9 85250 111
 TRACK NRZI
 TAPE SPEED 45, 37.5, 25, 18.75, 12.5 IPS
 TAPE 0.5 INCH, 1.5 MIL, IBM/ANSI
 COMPATIBLE, 10.5 INCH REEL
 TRACKS 7 OR 9
 DATA DENSITY PHASE ENCODED
 COMPATIBLE

HEWLETT PACKARD #7670R/C 8 4500 111
 TAPE FORMAT RCO, 556, OR 200 CPI NRZI AND
 1600 CPI PHASE-ENCODED
 CHANNELS 7 OR 9
 TAPE SPEED 10 TO 45 IPS
 TAPE 0.5 INCH, 1.5 MILS, IBM/ANSI
 COMPATIBLE

ORIGINAL PAGE 13
OF POOR QUALITY



RONG WAPNEL APERT 121
TRACKS 30
TAPE SPEED UP TO 1000 IPS
TAPE WIDTH 1/2 INCH
TAPE LENGTH 2400 FEET
PACKING DENSITY 15 KB/I/T
SIGNAL/NOISE 24 DB
DATA CAPACITY $1.3E+9$
BANDWIDTH 6-15 MB/SEC/T

LEACH WMTN 7000 311
TRACKS 12
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
TAPE LENGTH 9200 FEET
PACKING DENSITY 16.7 KB/I/T
SIGNAL/NOISE 22 DB
DATA CAPACITY $2.2E+10$
BANDWIDTH 2 MB/SEC/T

HEWLETT PACKARD #3659D #10200 111
BANDWIDTH 300 KHZ
CHANNELS 7
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3655C #14700 111
BANDWIDTH 300 KHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3650A-011 #23900 111
BANDWIDTH 500 HZ TO 2 MHZ
CHANNELS 14
RECORDING FORMAT DIRECT OR FM

HEWLETT PACKARD #3960A/13065A/13063A #4796 111
CONFIGURATION RACK MOUNTED
TAPE SPEED 15, 3 AND 3/4, 15/16 IPS
CHANNELS 4
RECORDING FORMAT FM
PASSBAND 5 KHZ
S/N RATIO 48 DB
TRACKS 30

HONEYWELL #5600 #9730
PORTABLE TAPE RECORDER
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 60 IPS
MAXIMUM BANDWIDTH(DIRECT) 300 KHZ
PACKING DENSITY UP TO 600 PBT
WEIGHT 70 LBS (32 KG)
INPUT VOLTAGE 28VDC

HONEYWELL #96 #17420
CHANNELS 7
SELECTABLE TAPE SPEED RANGE 15/16 TO 240 IPS
REEL SIZE 16 IN
MAXIMUM BANDWIDTH(DIRECT) 2M HZ
TAPE WIDTH 1/2 IN

AMPEX CORP #AR 70C #29,180 311
TRACKS 14
TAPE SPEED 60 IPS
TAPE WIDTH 1 INCH
RECORD MODE DIRECT
PACKING DENSITY 20 KB/I/T
SIGNAL/NOISE 20 DB
DATA CAPACITY $8E+10$ BITS
BANDWIDTH 1 MB/SEC/T

AMPEX CORP #AR 170C #28,000 112
TRACKS 14
TAPE SPEED 120 IPS
TAPE WIDTH 1 INCH
RECORDING MODE DIRECT
PACKING DENSITY 20 KB/I/T
SIGNAL/NOISE 20 DB
DATA CAPACITY $4.2E+10$
BANDWIDTH 2 MB/SEC/T

930700 POWER DIVIDER, WIDEBAND CN904 41111
930700
930700 DIVIDES ANTENNA POWER OUTPUT.
930700
930700 INPUT PORTS ONE 50 OHM PORT
930700 OUTPUT PORTS THREE 50 OHM PORTS
930700

HERBIMAC #PDM 3C 111
FREQUENCY RANGE 2-32 MHZ, 10-100 MHZ, 100-600 MHZ
COUPLING -4.8 DB
ISOLATION 30 DB, 30 DB, 25 DB
VSWR 1.3:1
POWER 5 WATTS
CONNECTORS (4) SMA



Space Division
Rockwell International

430000 X-Y PLOTTER CA907 11111
430000
430000 PROVIDES REAL-TIME ACCURATE REPRODUCTION OF SPECTRAL DATA SCURCES
430000 PLOTTED ON CARTESIAN COORDINATES.
430000
430000 CHANNELS 4 PLUG INS
430000 ACCELERATION 1500 IN/S/S (1800 CM/S/S)
430000 SENSITIVITY 0.5MV/IN (0.25 PV/CP)
430000
101340 KEYBOARD, COMPUTER CA909 11111
101340
101340 PROVIDES FOR CREWMAN TO COMMUNICATE WITH ONBOARD COMPUTER FOR
101340 EXPERIMENT CONTROL AND DATA ANALYSIS.
101340
101340 KEYBOARD TYPEWRITER TYPE
101340
RESEARCH INC 80PS-912 821-703 111
COMPLETE DATA DISPLAY SYSTEM MADE UP OF CRT ALPHANUMERIC DISPLAY
SEC-PPR-F MINS-COMPUTER, 83301 OPERATIONS CONTROL CCRASLE,
PRINTER/CARD READER, 8812-4A ANALYSEK AND 8812-13 UNIDIVER
CHANNELS 72
MEMORY SIZE 8K WORDS
RESEARCH INC 8812-3701 81555 111
SCREEN FORMAT 12 LINES X 72 OR 80 CHARACTERS,
24 LINES X 40 CHARACTERS
CHARACTER FORMAT 5 X 7 DOT MATRIX
TRANSFER RATE 110 TO 2400 BAUD, 10 OR 11 BIT
CHARACTERS
MODES HALF OR FULL DUPLEX-SWITCHABLE
LOCAL OR REMOTE
HEWLETT PACKARD 8260CA 83543 111
TERTRONIC 840024/021-008A-00 85552 111
DISPLAY MEDIUM 11 INCH DIRECT-VIEW, BISTABLE
STORAGE CRT WITH REFRESHED
SCRATCH PAD AREA
DISPLAY AREA 8.3 INCHES HORIZONTAL X 6.1
INCHES VERTICAL
ALPHANUMERIC MODE
FORMAT 39 LINES OF 95 NORMAL OR ITALIC
CHARACTERS IN MAIN AREA, ONE
LINE OF 84 CHARACTERS IN
SCRATCH PAD AREA
CHARACTER SET 96 UPPER AND LOWER CASE PRINT-
ING CHARACTERS (ASCII CODE)
70X50 HILS (CAN BE DOUBLE SIZE)
CHARACTER SIZE 7X9 DOT MATRIX
CHARACTER GENERATION PULSATING 7X9 MATRIX
CURSOR LINEAR INTERPOLATE, INCREMENTAL
GRAPHIC MODES PLOT, POINT PLOT, 1024X1024
ADDRESSABLE POINTS, 1024X768
VIEWABLE POINTS
GRAPHIC INPUT MODE 1024 (X), 768 (Y) POINTS, JOY-
STICK CONTROLLED, CROSS-HAIR
CURSOR
480000 CAMERA, TELEVISION CA909 11131
480000
480000 PROVIDES ASTRONAUT WITH VIEW OF EARTH SCENE.
480000
480000 VIEW ANGLE 20 DEGREES
480000 RESOLUTION 825 LINES
480000 SPECTRAL BANDWIDTH 4-14
480000 SWEEP RATE 15 FRAMES PER SECOND
480000
RCA ELECTRONICS COMPONENTS 88521 85520 111
TYPE VIDICON-SULFIDE TYPE
PHOTOCONDUCTOR 11
IMAGE DIAGONAL 25MM LIN.
FOCUS MANUAL
LIMITING RESOLUTION 1500 LINES
COMU INC. ELECTRONICS DIV. 81220 88430 111
IMAGE CONVERTER 4846P VIDICON - COLOR
RESOLUTION-HORIZONTAL LIMIT 300 LINES MIN.
GEOMETRIC DISTORTION < 2 % OF PIC HEIGHT
LENS MOUNT 16 MM - C-MOUNT
COMU INC. ELECTRONICS DIV. 82000 83000 311
LENS-BUILT IN 4:1 (20-80 MM) F/2.5 ZOOM
LENS ATTACHMENT 4:1 (12.5-80 MM)
LENS INTERCHANGEABLE 10:1 (15-150 MM) F/2.8 ZOOM
16 MM C-MOUNT FFL
10 MMZ BANDWIDTH (MOD 2000) 525 OR 724 LINE (VIDICON 7263A)
20 MMZ BANDWIDTH (MOD 2004) 873 OR 945 LINE (VIDICON 8573)
PIL SPECIFICATIONS PIL-E-5272C, PIL-STD-810
COMU INC. ELECTRONICS DIVISION 84500 82225 311
LENS 4:1 ZOOM
10:1 ZOOM
VERTICLE SWEEP RATE 60 FIELDS PER SEC
HORIZONTAL SWEEP RATE 525 LINES PER FRAME
IMAGE TUBE TYPE 8541A STD
LENS MOUNT 16 MM C-MOUNT
PIL SPECIFICATIONS PIL-E-5400P, PIL-E-5300P

ORIGINAL PAGE IS
OF POOR QUALITY

424200 RECORDER, VIDEO CN910 11111
424200 PROVIDES FOR RECORDING VIDEO DATA.
424200 CHANNELS 2
424200 BANDWIDTH 4 MHZ MINIMUM
424200 DURATION 1.5 HRS MINIMUM
424200 SANGAMO ELECTRIC CO #7C-2 121
TAPE WIDTH 1 INCH
TAPE SPEED 120 IN./SEC.
FREQ. RESPONSE (DIRECT) 400HZ TO 2 MHZ
SIGNAL/NOISE 20DB
CHANNELS 14
FORMAT FM
618080 RECEIVER, PHASE-LOCK CN911 11111
618080 PROVIDES FOR TRACKING SIGNAL PHASE.
618080 FREQUENCY SPECTRUM 1570 MHZ TUNABLE
618080 SENSITIVITY 15 KHZ
618080
104180 DISPLAY, CRT CN912 11111
104180 DISPLAY REAL-TIME OR VIDEO TAPE DATA.
104180 RESOLUTION 625 LINES
104180 FRAME RATE 15 FPS
104180
RESEARCH INC #3300 #1580
FORMAT 24 LINESX72 OR 80 CHARACTERS
..... 24 LINESX40 CHARACTERS
..... 12 LINESX72 OR 80 CHARACTERS
REFRESH RATE 60 HZ
TRANSFER RATE 110 TO 2400 BAUD
CHARACTER FORM 5X7 DOT MATRIX
DIMENSIONS 19-1/2X13-1/2XHX23-1/20 INCH
HEIGHT 3.00 (17.7KG)
HUNKER-BAND CORP #221-12CPT #950
CHARACTER CAPACITY (MAX) 960
CHARACTERS/LINE (MAX) 80
LINES/DISPLAY (MAX) 24
CHARACTERS REPERTOIRE 62 OR 92
VIEWING AREA 8.75 X 6.25
REFRESH AREA 54 FRAMES/SEC
CHARACTER GENERATING METHOD ... 5X7 DOT MATRIX
787RONIX INC #4CC2A/021-C033-00/ #10979 111
021-00XX-00/4901/
4991
ALPHANUMERIC FORMAT 19 LINES-45 NORMAL/ITALIC CHAR,
1 LINE- 84 CHAR IN SCRATCH AREA
CHARACTER SET 96 UPPER/LOWER CHAR (ASCII)
CHARACTER SIZE 70 X 90 MILS
CHARACTER GENERATION 7 X 9 DOT MATRIX
CURSOR PULSATING 7 X 9 MATRIX
GRAPHIC INPUT MODE 1024(X) X 1024(Y) ADDRESSABLE,
1024(X) X 768(Y) VIEWABLE PTS;
JOYSTICK CONTROLLED, CROSS-HAIR
CURSOR.
429900 VISICORDER CN914 11111
429900 REAL-TIME RAW DATA RECORDING.
429900 BANDWIDTH DC TO 5000 HZ
429900 CHANNELS 18
429900 RECORDING FORMAT CRT HARD COPY
429900 HONEYWELL #1856 #
RESOLUTION 200 ELEMENTS/INCH (1000/SWEEP)
RECORDING RATE 88-06 BITS/SEC
CRT FIBER OPTICS
SPOT SIZE 0.005 INCH
SWEEP SPEED 20,000 LINES/SEC
HONEYWELL #1856 #
RESOLUTION 200 ELEMENTS/INCH (1000/SWEEP)
RECORDING RATE 88-06 BITS/SEC
CRT FIBER OPTICS
SPOT SIZE 0.005 INCH
SWEEP SPEED 20,000 LINES/SEC



Space Division
Rockwell International

APPENDIX D

SPACELAB/EXPERIMENT EQUIPMENT
INTERFACE REQUIREMENTS SPECIFICATION
(SEEIR)

PREPARED BY <i>J. J. Hendricks</i>	CODE IDENT. NO.: 03933	NUMBER
APPROVALS <i>R. P. Hues</i>	SPACE DIVISION	TYPE REQUIREMENTS
	ROCKWELL INTERNATIONAL CORPORATION	DATE 4/25/74
	SPECIFICATION	SUPERSEDES SPEC. DATED: 2/15/74
		REV. LTR. A PAGE 1 of 19

TITLE

SPACELAB/Experiment Equipment Interface Requirements

THIS DOCUMENT IS INTENDED AS A WORKING
TOOL FOR THE STUDY ENTITLED ANALYSIS OF
COMMERCIAL EQUIPMENT AND INSTRUMENTATION
FOR SPACELAB PAYLOADS. IT IS NOT A CONTRACTUAL
SPECIFICATION.

NUMBER	REVISION LETTER						PAGE 2
	A						

1. SCOPE

1.1 Scope. This document delineates the requirements for experiment equipments within pressurized areas of the Spacelab. Requirements include those that are mandatory for preservation of Spacelab-Shuttle safety margins and those necessary to achieve intended operation of equipments. Meeting the requirements herein does not, however, assure the ability to operate within Spacelab in an arbitrary fashion since combined equipment support demands could exceed total Spacelab capability. Integrated operation must therefore be coordinated with Spacelab operations management for specific missions in order to establish mutual operating timelines and allocation of Spacelab resources.

2. APPLICABLE DOCUMENTS

2.1 Applicability. The following documents were used in the derivation of the requirements and are referenced herein by title or by number in the text margins as appropriate. Interpretations were made where conflicts in documents existed or data were unavailable.

DOCUMENTS

- | | |
|-------------------------------|--|
| (1) ERNO Spacelab | Preliminary Design Requirements Specifications,
December 1973; Final Briefing, February 1974 |
| (2) MBB Spacelab | Preliminary Design Requirements Specifications,
December 1973 |
| (3) SL-E-0001 | Electromagnetic Compatibility Requirements, Systems for
the Space Shuttle Program, June 4, 1973 |
| (4) SL-E-0002 | EMI Interference Characteristics, Requirements for
Equipment for the Space Shuttle Program, June 4, 1973 |
| (5) SD 72-SH-0172 | Space Shuttle Orbiter Materials Control and Verification
Plan, November 1972 |
| (6) NHB 8060.1 | Flammability, Odor and Off-Gassing, Requirements for
Materials in Environments which Support Combustion,
November 1971 |
| (7) U719-10-206 | Material Code Directory by Material Description
- Nonmetals, December 15, 1973 and subsequent issues |
| (8) U719-10-207 | Material Code Directory by Material Description
- Metals, December 15, 1973 and subsequent issues |
| (9) U719-10-208 | Material Selection Rating and Criteria Report
- Metals, December 15, 1973 and subsequent issues |
| (10) U719-10-209 | Material Selection Rating and Criteria Report
- Nonmetals, December 15, 1973 and subsequent issues |
| (11) Federal Standard
209A | Clean Room and Work Station Requirements for Controlled
Environments, CSA, August 10, 1966 |

NUMBER	REVISION LETTER							PAGE
	A							

- (12) MSC00720 Solar-Powered Space Station Preliminary Design
- (13) No Number Draft, Interim Spacelab Reference Document, February 1974
- (14) ESRO RFP Request for Proposal for the Design and Development Contract
AO/600 (Phase C/D), Appendix 11, System Requirements, March 1, 1974
- (15) JSC 07700 Space Shuttle System Payload Accommodations, Revision A,
Volume XIV July 1973

3. REQUIREMENTS

3.1 General.

3.1.1 Spacelab Design Characteristics. The following Spacelab characteristics establish the basis for the experiment equipment requirements specified in paragraph 3.2, and establish the maximum total Spacelab support available for combined experiment and Spacelab subsystems hardware.

3.1.1.1 Physical Perspective. Figure 1 provides the general Spacelab configuration and defines the axes coordinate system used to interpret applicable specifications herein. The Spacelab consists of pressurized module and unpressurized pallet elements flown together or independently, depending upon the mission. Size of the Spacelab elements may vary, depending upon the nature of the mission. The pressurized module consists of a support section containing Spacelab subsystems as well as accommodations for experiment equipment, plus an experiment section that can be configured in various lengths to accommodate additional experiment equipments. This specification does not cover pallet-mounted equipment requirements. Pallet equipment, however, shares available power and data handling support with pressurized module equipments. (1)


3.1.1.2 Experimenter Support Characteristics. The support resources available to combined experimenter payloads follow.

3.1.1.2.1 Equipment Physical Accommodations. Approximately 500 ft³ (14 m³) is available above the floor in the short Spacelab, and approximately 800 ft³ (~22 m³) is available in the larger Spacelab for rack/bench/shelf-mounted experiment equipments. Refer to Figure 1 for example equipment locations. In addition, equipment may be attached to the underside of floor panels at up to 35 lb/ft² (172 kg/m²) loading. (1)

3.1.1.2.1.1 Rack Mounting Capability. Standard racks are available that can accommodate equipments up to approximately 60 lb/ft³ (300 kg/m³). Standard racks tie into Spacelab structure and are sized for MIL-STD-864 19-inch rack modules. NIM, ATR, ARINC 404 or MIL-STD-172 type electronic packages may be accommodated in appropriate substructure. Packages up to 39.3 inches (1 m) wide may be custom-mounted in double wide racks/shelves. Larger sizes require detailed design installation. (1)

3.1.1.2.1.1.1 Component Mounting. Components are normally fastened to their rack with screws through the front-mount surface. Equipment front panels will prevent mixing of rack cooling air with cabin air, minor leakage excepted. Guide rails with rigid-mate support fixtures are available for long/heavy items. (1)

3.1.1.2.1.1.2 Standard Rack Module Size. Maximum component size (less mounting surface) is 17-1/2 in. wide x 30 in. deep x height (0.44 m x 0.76 m x height). (1)



X = 0 is at forward end of Spacelab

Y = 0 is Orbiter $Y_0 = 0$

Z = 0 is Orbiter $Z_0 = 400$

- ① WORK BENCH
- ② CONTROL CENTER
- ③ MODULAR EQUIPMENT RACK

FIGURE 1. SPACELAB IN ORBITER ILLUSTRATING POTENTIAL GENERAL EXPERIMENTER OCCUPANCY REGIONS AND AXES COORDINATE SYSTEM

NUMBER

A

REVISION LETTER

PAGE 4

CODE IDENT. NO. 03953

NUMBER	REVISION LETTER						PAGE 5
	A						

3.1.1.2.1.2 Equipment Stowage Boxes. Standard stowage containers that meet required safety and crew habitability requirements will be available to stow equipments, as required, in standard 19-in. rack widths. Stowage may be at various Spacelab locations. Custom boxes up to 39.3 in. (1 m) in width may be utilized. Containers may be packed by experimenters and direct-loaded onto Spacelab provided that a certification is completed which states the absence of flammability/explosion/corrosive/toxicity hazards and/or appropriate hazard isolation and cleaning procedures have been implemented. Fitted foam materials meeting flammability, off-gassing and contamination requirements may be used as separators and environment isolators within the stowage boxes.

(1)

3.1.1.2.2 Electrical Power. Total power available to experiment equipments cannot exceed 4 kW continuous, or 9 kW peak for 15 minutes every 3 hours. Greater demands require special primary power sources to be added. The power forms available to experimenter equipments are provided below.

(1)

(2)

(14)

(13)

Basic

- . 24 to 32 vdc, 1 v P-P ripple: Orbiter 28 vdc plus booster battery and limiter to eliminate transients
- . 115/200 vac \pm 5%, 400 Hz \pm 1 Hz: Spacelab dc-ac converters, 3-phase Y-connected, 5% harmonics, 20% maximum phase-to-phase load unbalance usable as single phase

Add-on

- . 115 vac \pm 5%, 60 Hz \pm 1 Hz, 5% harmonics: Spacelab dc-ac converters, single phase
- . 220 vac \pm 5%, 50 Hz, \pm 1 Hz, 5% harmonics: Spacelab dc-ac converters, single phase
- . 28 vdc \pm 2%, 50 mv rms ripple: Spacelab dc-dc converters

Maximum continuous power from each secondary power source can be adjusted by adding converters.

3.1.1.2.2.1 Emergency Battery. An emergency battery is available to power only critical Spacelab subsystems and the caution and warning network.

(1)

(2)

3.1.1.2.2.2 Power Control and User Access. Power panels are conveniently located with circuit breakers and switchable to protect the power sources and to provide for emergency disconnect.

(1)

3.1.1.2.2.3 Launch, Reentry and Landing Power. Experimenter equipment will normally be OFF during launch, reentry and landing. A total of 1 kW continuous, 1.5 kW peak, will be available for Spacelab subsystems and experimenters with critical power needs during these phases of operation.

(1)

(2)

3.1.1.2.2.4 Ground Checkout Operation. Integration and ground test operations in the Spacelab are supported by a GSE Orbiter power source simulator that provides power to on-orbit operation specifications.

3.1.1.2.3 Data Management and Control Functions. A dedicated data management and control subsystem is provided by Spacelab for experimenter use. Provisions are made for input and control, tape storage, display, playback and transmit.

(1)

(2)

(13)

NUMBER	REVISION LETTER	PAGE 6
	A	

3.1.1.2.3.1 On-Board Data Management Capability. The basic total capability to support all experiments on a Spacelab mission follows.

- . Sensors accepted - Up to 512 digital or converter analog signals
- . Discrete signals accepted - Up to 512 binary-level signals
- . Combined sensor data rates - Up to 30 Mbps digital
- . Wideband analog inputs - 5 channels of 6 MHz bandwidths each
- . Digital data storage - 4 x 10¹⁰ bits/reel at up to 30 Mbps
- . Analog data storage - 30 minutes/reel dc to 6 MHz
- . Computer memory - 32K or 64K 16-bit words

3.1.1.2.3.2 Data System-Experimenter Interface. Current concepts provide standard remote acquisition and control units (RAU's) for direct connection of individual sensor for analog-to-digital conversion and multiplexing/data formatting to a compatible data bus. Up to 8 RAU's, each accepting up to 64 analog and 64 discrete inputs, are available. Command and control are routed through the RAU's.

3.1.1.2.3.3 Ground Communications. Downlink data transmission is normally handled via the Shuttle subsystems. The following total basic capability is provided.

Downlink

- . PCM Digital
 - 25 kbps through TDRS/STDN
 - 256 kbps through STDN
 - 50 Mbps through TDRS (KU-Band)
- . Analog
 - TV (color) through STDN or TDRS (KU-Band)
 - Combined data (6 MHz bandwidth) through TDRS (KU-Band)
 - [Alternate to TV through TDRS (KU-Band)]
- . Voice
 - 2 channels, one dedicated to payload, through STDN or TDRS/STDN

Uplink

- . Voice
 - 2 channels, one dedicated to payload, through STDN or TDRS/STDN
- . Command
 - 2 kbps through STDN or TDRS/STDN
- . TBD digital and analog data through TDRS (KU-Band)

3.1.1.2.3.4 Ground Checkout Operation in Spacelab. On-orbit operation is simulated.

3.1.1.3 Spacelab Internal Environments.

3.1.1.3.1 Normal Atmosphere. Atmospheric pressure throughout a mission, from launch to landing, is maintained in the range of 14.7 \pm 0.25 psia.

3.1.1.3.1.1 Atmospheric Composition. Launch-to-landing atmospheric composition is nominally equivalent to sea-level air (21 parts oxygen to 79 parts nitrogen, 760 \pm 13 mmHg), with maximum carbon dioxide of 7.6 mmHg partial pressure. Other elements as listed below (due to materials off-gassing, etc.) can be expected. Ground integration operations utilize ambient earth atmosphere conditioned to be compatible with

NUMBER	REVISION LETTER							PAGE 7
	A							

operational requirements (moisture, particles, etc.). Nominal air movement is 5 to 12 meters per minute.

3.1.1.3.1.2 Humidity. Maximum operational relative humidity will be 70 percent with the temperature of all points of the Spacelab interior above the dew point. Ground storage of Spacelab permits relative humidities of between 20 and 90 percent.

(13)

(14)

3.1.1.3.1.3 Contamination Levels.

3.1.1.3.1.3.1 Particles. The Spacelab environment will be maintained at a Class 100,000 environment per Federal Standard 209A up to launch. Subsequent to launch, Spacelab filtration removes particles larger than 5 microns.

(13)

3.1.1.3.1.3.2 Chemicals. The following contaminants may be present (but not worst effect additive) up to the indicated maximum concentrations in milligrams per cubic meter in flight or ground operations.

(1)

(12)

ACETONE	240	CUMENE	25
ACETALDEHYDE	36	CYCLOHEXANE	100
ACETIC ACID	2.5	CYCLOHEXENE	100
ACETYLENE	180	CYCLOHEXANOL	20
ACETONITRILE	7	CYCLOPENTANE	100
ACROLEIN	0.25	CYCLOPROPANE	100
ALLYL ALCOHOL	0.5	CYANAMIDE	45
AMMONIA	3.5	DECALIN	5.0
AMYL ACETATE	53	1,1 DIMETHYL CYCLOHEXANE	120
AMYL ALCOHOL	36	TRANS 1,2 DIMETHYL CYCLOHEX	120
BENZENE	8	2,2 DIMETHYL BUTANE	93
N-BUTANE	180	DIMETHYL SULFIDE	15
ISO-BUTANE	180	1,1 DICHLORO ETHANE	40
BUTENE-1	180	DI-ISO-BUTYL KETANE	29
CIS-BUTENE-2	180	1,4 DIOXANE	36
TRANS-BUTENE-2	180	DIMETHYL FURAN	3.0
1,3 BUTADIENE	220	DIMETHYL HYDRAZINE	0.1
ISO-BUTYLENE	180	ETHANE	180
N-BUTYL-ALCOHOL	30	ETHYL ALCOHOL	190
ISO-BUTYL-ALCOHOL	30	ETHYL ACETYLENE	180
SEC-BUTYL-ALCOHOL	30	ETHYL BENZENE	44
TERT-BUTYL-ALCOHOL	30	ETHYLENE DICHLORIDE	40
BUTYL ACETATE	71	ETHYL ETHER	120
BUTRALDEHYDES	71	ETHYL BUTYL ETHER	200
BUTYRIC ACID	14	ETHYL FORMATE	30
CARBON DISULFIDE	6	ETHYLENE	180
CARBON MONOXIDE	29	ETHYLENE GLYCOL	114
CARBON TETRACHLORIDE	6.5	TRANS 1, METHYL 3 ETHYL C	117
CARBONYL SULFIDE	25	ETHYL SULFIDE	97
CHLORINE	1.5	ETHYL MERCAPTAN	2.5
CHLOROACETONE	100	FREON 11	560
CHLOROBENZENE	35	FREON 12	500
CHLOROFLUORMETHANE	24	FREON 22	350
CHLOROFORM	24	FREON 23	12
CHLOROPROPANE	84	FREON 113	700
CAPRYLIC ACID	155	FREON 114	700

NUMBER	REVISION LETTER							PAGE 8
	A							

FREON 114 (UNSYM)	700	PROPYLENE	180
FREON 125	25	ISO-PENTANE	295
FORMALDEHYDE	0.6	N-PENTANE	395
FURAN	3	PENTENE-1	180
FURFURAL	2	PENTENE-2	180
HYDROGEN	215	PROPANE	180
HYDROGEN CHLORIDE	0.15	N-PROPYL ACETATE	84
HYDROGEN FLUORIDE	0.08	N-PROPYL ALCOHOL	75
HYDROGEN SULFIDE	1.5	ISO-PROPYL ALCOHOL	98
HEPTANE 9	200	N-PROPYL BENZENE	44
HEXENE-1	180	ISO-PROPYL CHLORIDE	260
N-HEXANE	180	ISO-PROPYL ETHER	120
HEXAMETHYLCYCLO-THRISILOHEXANE	240	PROPRIONALDEHYDE	30
INDOLE	126	PROPIONIC ACID	15
ISOPRENE	140	PROPYL MERCAPTAN	82
METHYLENE CHLORIDE	21	PROPYLENE ALDEHYDE	10
METHYL-ACETATE	61	PYRUVIC ACID	0.9
METHYL-BUTENE	1430	PHENOL	1.9
METHYL FURAN	3	SKATOL	141
METHYL ETHYL KETONE	59	SULFUR DIOXIDE	0.8
METHYL ISOBUTYL KETONE	41	STYRENE	42
METHYL ISOPROPL KETONE	70	TETRACHLOROETHYLENE	67
METHYL CYCLO HEXANE	200	TETRAFLUORETHYLENE	205
METHYL ACETYLENE	165	TETRAHYDROFURANE	59
METHYL ALCOHOL	26	TOLUENE	75
3-METHYL PENTANE	295	TRICHLOROETHYLENE	
METHYL METHACRYLATE	41	1,2,4 TRI METHYL BENZENE	49
METHANE	1720	1.1.3 TRI METHYL CYCLOHEX	140
MONO METHYL HYDRAZINE	0.035	VALERALDEHYDE	70
METHYL MERCAPTAN	2	VALERIC ACID	110
NAPHTHALENE	5.0	VINYL CHLORIDE	130
NITRIC OXIDE	32	VINYL METHYL ETHER	60
NITROGEN TETROXIDE	1.8	VINYLDENE CHLORIDE	20
NITROGEN DIOXIDE	0.9	O-XYLENE	44
NITROUS OXIDE	47	M-XYLENE	44
OCTANE	235	P-XYLENE	44

3.1.1.3.1.4 Emergency Conditions Atmosphere. The Spacelab can be depressurized in space to near vacuum within TBD minutes if necessary to combat fire, smoke, and toxics. Experiment equipments are not required to operate after a depressurization.

3.1.1.3.2 Spacelab Internal Temperatures. The following temperatures can occur inside the Spacelab:

- . Air Temperature: Adjustable 65 F to 85 F + 2 F (19 C to 26 C \pm 1 C)
- . Wall Temperature: 38 F (3 C) to 118 F (48 C)
- . Ground Storage or Nonoperating: 40 F to 130 F (4 C to 55 C)

Exposed equipment surfaces shall have touch temperatures of less than 115 F (46 C).

3.1.1.3.2.1 Equipment Cooling. Coldplates and forced-air cooling provided by standard Spacelab facilities with the following controlled characteristics can be assumed (with equipment operating or nonoperating):

NUMBER	REVISION LETTER						PAGE 9
	A						

- . Rack cooling air temperature: 79 F to 84 F (26 C to 29 C)
(Cooling air isolated from cabin air)
- . Maximum coldplate temperature: 104 F (40 C)
(Requires Spacelab interface coordination.)

3.1.1.3.3 Shock, Vibration, and Acceleration. The following describes the random and sinusoidal vibration spectrums and steady-state and/or very low frequency accelerations that equipments mounted within Spacelab may encounter during launch, reentry, and landing as transmitted from the Orbiter through the Spacelab mounting and structure dynamics. Controls on all ground handling will limit these inputs to values less than worst-case flight values. (1)

3.1.1.3.3.1 Accelerations. Equipment accelerations are the same as specified for Orbiter (see Figure 1 for axes orientation). Maximum values throughout the mission (launch to landing) are as follows:


	X-AXIS	Y-AXIS	Z-AXIS
. Maximum g's (normal launch, reentry, landing)	-3.3 +1.5	+1.5	+3.0 -1.0
. Crash landing (hold-down and restraint fixture limit)	+9.0 -1.5	+1.5	+4.5 -2.0

3.1.1.3.3.2 Random Vibration. The vibration seen by the modules within Spacelab depends upon their rigid-body clustering masses, location of mountings to Spacelab, and whether or not shockmounts are employed. Table I summarizes the G rms values over the 20 to 2000 Hz-band for three weight group categories, three Spacelab mounting locations, and for both with and without shock mounts. Figure 2 provides spectral density characteristic curves. (1) (13)

TABLE I. G RMS OVER THE RANGE OF 20 TO 2000 Hz (ANY AXIS)

LOCATION IN SPACELAB	G RMS FOR INDICATED GROUP WEIGHTS		
	< 33 LB	33 TO 67 LB	> 67 LB
HARD-MOUNTED			
CYLINDRICAL AREAS	14.0	11.7	10.6
BULKHEAD AREAS	10.6	10.6	10.6
WALLS AND FLOORS	9.7	9.5	9.1
SHOCK-MOUNTED			
CYLINDRICAL AREAS	4.4	3.7	3.3
BULKHEAD AREAS	3.3	3.3	3.3
WALLS AND FLOORS	3.1	3.0	2.9

3.1.1.3.3.3 Sinusoidal Vibration. All weight groups, whether hard- or shock-mounted, may experience sinusoidal vibrations as given by Table II (in addition to the random vibration of paragraph 3.1.1.3.3.2 and accelerations of paragraph 3.1.1.3.3.1), throughout the Spacelab during launch and reentry periods. (1)



(1)

ORIGINAL PAGE IS
OF POOR QUALITY

FIGURE 2.

SPECTRAL DENSITY CURVES FOR INDICATED
RIGID-BODY MASSES, SPACELAB MOUNTING
LOCATIONS AND SHOCK MOUNTING

NUMBER	REVISION LETTER							PAGE 11
	A							

TABLE II. SINUSOIDAL VIBRATION CHARACTERISTICS IN SPACELAB

AXIS	FREQUENCY (Hz)	AMPLITUDE
X	3 TO 8.5	0.8-INCH DOUBLE AMPLITUDE
	8.5 TO 35	+3.0 G PEAK TO PEAK
	35 TO 50	+1.0 G PEAK TO PEAK
Y AND Z	3 TO 7	0.8-INCH DOUBLE AMPLITUDE
	7 TO 35	+2.0 G PEAK TO PEAK

3.1.1.3.3.4 Shock. Shock transmitted by Spacelab structure in flight or on the ground will not exceed the vibration/accelerations anticipated for the flight environment. (1)

3.1.1.3.4 Acoustic Noise Levels within Spacelab. The noise levels within Spacelab interior during launch and reentry phases are as defined in Figure 3. These are the maximum levels over complete missions. The on-orbit working audio noise levels from all combined sources will not exceed the values of the NC-50 curve on standard noise criteria (NC) curves. This curve is approximately defined by 73 dB at 20 cps, linearly decreasing to 55 dB at 300 cps to 47 dB at 4800 cps (dB referenced to 0.0002 microbar). (14)

3.1.1.3.5 Atomic Radiation. Radiation levels, from natural or other sources, within Spacelab are low enough for safe human exposure. Therefore, equipment operation is not physically affected except where principle of operation is affected by background counts such as sensitive films, etc. Film vault storage will be provided. (1)

3.1.1.3.6 Electromagnetic Interference (EMI).

3.1.1.3.6.1 EMI Generator Control. The net effect of radiated and conducted EMI caused by Spacelab subsystems and experimenter equipments will be controlled by observing good design practices for grounding, shielding, signal and power transmission wiring, and bonding impedance. The intent is to achieve a 6-dB minimum safety margin between maximum imposed EMI levels and operative threshold levels of affected circuits in accordance with SL-E-0001 and SL-E-0002. (14)

3.1.1.3.6.2 EMI Susceptibility. Equipment will nominally be required to be able to operate in EMI environments in accordance with SL-E-0001 and SL-E-0002 with a 6-dB safety margin between the maximum environment levels and operating thresholds. (1)

3.1.1.3.7 Zero-G Effects. The resultant average acceleration of Spacelab is zero in all axes while on orbit (no artificial gravity). Minor perturbations between $10^{-2}g$ and $10^{-3}g$ may occur resulting from crew movement, manipulator activation, water pumps, etc. During attitude maneuvers, accelerations of up to -0.2 and +0.1 g along the X axis and +0.1 g along the Y and Z axes may occur. (13)

3.1.1.4 Spacelab Safety Provisions. The following, in conjunction with minimum requirements to be met by equipment brought on board, provides for crew and Orbiter-Spacelab safety. (15)

3.1.1.4.1 Caution and Warning (C&W) Devices. Identified potential flammability, toxicity and contamination hazards require warning devices to allow sufficient time for corrective action. Such warning devices operate independently from other operational circuitry and power supplies. Spacelab provides fire, smoke and toxics warnings for (1)
(14)

BAND LEVEL - DB (REF TO 2×10^{-5} N/M²)

30

100

300

1000

3000

ONE-THIRD OCTAVE CENTER FREQUENCY - Hz

(1)

FIGURE 3. ACOUSTIC NOISE INSIDE SPACELAB

NUMBER

A

REVISION LETTER

PAGE 12

CODE IDENT. NO. 03953

the general environment. Hazardous experiment equipment requires additional caution/warning sensors to be integrated into the monitoring and control subsystem to enable appropriate visual and audio alarms when critical tolerances are exceeded (caution) and when immediate corrective action is needed to prevent injury or loss of personnel/equipment (warning).

3.1.1.4.2 Fire Extinguishers. At least three Apollo-type fire extinguishers will be strategically located throughout Spacelab.

(1)
(14)

3.1.1.4.3 Depressurization. In emergencies the Spacelab may be depressurized in TBD minutes. The Orbiter provides a safe haven for Spacelab crew members. Equipments must not erupt with new hazards upon such depressurization and subsequent repressurization.

(1)
(14)

3.1.1.4.4 Oxygen Masks. Oxygen masks and 30-minute portable oxygen supplies will be available if the atmosphere becomes unsafe for personnel.

(1)
(14)

3.1.1.4.5 Shatterable Materials. All materials that could conceivably shatter for any specified Spacelab condition (including ground operations, depressurization, crash landing, launch and reentry forces) or credible accident, oversight, etc., shall be contained or constrained adequately to prevent debris and missiles from being emitted that could impair habitation quality or equipment operations, or cause injury or further damage.

(1)
(14)

3.1.1.4.6 Protrusions and Sharp Edges/Corners. All equipment surfaces exposed to bodily contact by personnel during normal operations shall have rounded (0.5-inch radius) edges and corners, or have guards or other protection to serve an equivalent function.

(1)

3.1.1.4.7 Crash-Landing Provisions. All equipments are mounted, restrained or stowed to withstand the crash-landing loads of paragraph 3.1.1.3.3.1 without erupting hazardous fluids or releasing loose missiles or debris.

(1)

3.1.1.4.8 Flammability. Spacelab structures and subsystem materials and design are equivalent to Orbiter design criteria and utilize materials rated "A" or "B" for use in air as defined by SD 72-SH-0127, "Space Shuttle Orbiter Materials Control and Verification Plan." Materials are flammability-rated "A" where noncombustible or self-extinguishing upward in air, without spark, drip or sputter of flaming particles within 6 inches of the point of ignition at the bottom of a vertically oriented specimen (per NHB 8060.1). Such "A" rated materials are suitable for unrestricted use in air. "B" rated materials for use in air have propagation downward rates in air of less than 0.3 in./sec and have flash and fire points greater than 450 F. "B" materials may be utilized providing that discrete quantities which could mutually support combustion are limited to less than 1 pound and 40 in.² of surface area, or are in containers (with or without ignition sources) that prevent propagation of flames or burning debris. Materials exposed to LGOX, HGOX, LOX, N₂H₄, or other non-air environments are not covered in this document, and special requirements and evaluations may be required for such applications.

(5)

3.1.1.4.8.1 Storage of Flammable Materials. Fire-proof storage is provided for stock and waste tissues, etc.

(14)
(1)

NUMBER	REVISION LETTER							PAGE 14
	A							

3.1.1.4.9 Off-Gassing, Toxics/Odors. Spacelab environmental protection requires controlling materials off-gassing similar to Shuttle requirements. In accordance with SD 72-SH-0172, materials rated "A" or "B" for toxicity are utilized. A material is toxicity-rated "A" if it has ≤ 250 micrograms per gram total organics, ≤ 50 micrograms per gram of carbon monoxide, and has an average odor rating ≤ 3.0 , when tested per NHB 8060.1. A "B" material does not meet one or more "A" requirements, but will meet an "A" requirement in conjunction with specific controls such as coatings, hermetic sealing, etc.

3.1.1.4.9.1 Prohibited Materials. The following shall not be brought into Spacelab in any measurable quantity in any form:

- . Mercury
- . Methyl chloroform
- . Cadmium

3.1.1.4.10 Physical Equipment Restraints. All materials, parts, and equipments shall be secured in such a way as to reliably prevent loose items from drifting or accelerating against or within other equipments or personnel in such a way as to cause damage, operational interference or injury for all mission modes and phases, including emergency modes.

3.1.1.4.11 Touch Temperatures. See paragraph 3.1.1.3.2.

3.2 Experimenter Equipment/Material Requirements.

3.2.1 Mandatory Safety Requirements. The following requirements are mandatory, unless specifically waived by Orbiter-Spacelab management, in order to ensure Orbiter-Spacelab and personnel safety.

3.2.1.1 Equipment Construction and Installation. The equipments and their constituent parts shall not shatter, come loose from their mountings, or erupt to generate flying missiles or release hazardous or corrosive substances when exposed to the crash environment of paragraph 3.1.1.3.3.1, or the random plus sinusoidal environments of paragraphs 3.1.1.3.3.2 and 3.1.1.3.3.3, for applicable mounting conditions, or during emergency depressurizations to near-vacuum and repressurization per paragraph

3.1.1.4.3. When the equipments are stowed in protective containers, the requirements apply to the stowage box and equipment combined. All materials and parts subject to accidental or intentional breakage during normal on-orbit operations that could result in missiles, abrasive or injurious particles, or release of toxics, corrosives or contaminants immediately or at a later mission time shall be constrained or guarded to preclude such possibilities. Parts and equipments shall be secured from drifting away or from unintended accelerations with positive mechanisms designed to function properly after exposure to launch environments.

3.2.1.1.1 Sharp Corners, Edges, Protrusions. All equipment surfaces exposed to personnel operations shall not protrude into travel lanes and shall have rounded edges and corners or guards with 0.5-inch minimum radius to prevent bodily contact with edges/corners of lesser radii and protrusions that could cause injury.

3.2.1.2 Flammability. All equipments to be placed on Spacelab shall be evaluated for flammability. Equipments meeting requirements of paragraph 3.1.1.4.8 shall be acceptable for Spacelab. The current issue of Rockwell Space Division computerized material analysis tracking and control system reports [U719-10-206 and U719-10-209 (for non-metals) and U719-10-207 and U719-10-208 (for metals)] may be used as the basis for evaluating flight acceptability of equipments with respect to specific constituent materials. Unidentified materials as a minimum shall be upgraded to a verified "B" rating by approved coating, replacement, or other means, except that unidentified materials in a completely enclosed unit that precludes the escape of potentially flaming particles need not be upgraded.

NUMBER

REVISION LETTER

PAGE 15

A

3.2.1.3 Toxicity. Usage of materials in equipments that conform to "A" or "B" ratings per paragraph 3.1.1.4.9 shall be acceptable to place on board Spacelab. The materials analysis tracking and control reports (U719-10-206 and U719-10-209) for exposure categories A and B shall be utilized to evaluate equipment material toxicity. Materials rated "A" may be utilized unrestricted. Materials rated "B" may be utilized without controls (coatings, etc.) to the extent that the total organics and carbon monoxide per equipment or group of equipments do not exceed the values that could result if the equipment consisted of all "A" materials containing maximum allowable organics and carbon monoxide. Before being placed in the Spacelab, items shall be operated or otherwise "baked out" for a sufficient time period to drive off excessive surface gases that may reside due to manufacturing or cleaning processes.

3.2.1.4 Electric Shock. External housings shall be grounded to preclude shocks to personnel from either electrical leakage or static charge buildup. Grounding satisfying EMI requirements will satisfy this requirement.

3.2.1.5 Caution and Warning Networks. Any potential hazard to crew or total Spacelab shall be instrumented to provide "0" or "1" state (safe or caution/warning) signals to the Spacelab caution and warning network. This network operates from an independent power supply. The Spacelab operations management will integrate the signal into the C&W network as warranted by the degree and conditions for which a hazard exists.

3.2.2 Mandatory Environmental Preservation Criteria. The following requirements are mandatory in order to maintain a nominal Spacelab environment for all experimenter payloads and to assure Spacelab subsystems performance.

3.2.2.1 Materials and Construction. Equipments meeting the safety requirements of applicable subparagraphs of paragraph 3.2.1 above will meet requirements of this section. Calculated verification of structural safety is acceptable providing that yield and ultimate safety factors are at least 2 and 3 respectively. Liquid and gas lines up to 1.5-in. diameter are satisfactory if ultimate pressure ratings are at least 4 times peak limit pressure, and if lines are unpressurized during launch and reentry. Liquid lines shall not be exposed to provide possible hand rails or obstacles to personnel/equipment movement.

3.2.2.2 Contamination Generation.

3.2.2.2.1 Particles. Equipment and packing materials shall be subjected to a cleaning operation compatible with installation in a Class 100,000 area prior to on-loading into Spacelab (whether or not in stowage containers). Exterior or non-sealed interior of equipments shall have coatings or finishes that do not show visible signs of flaking or peeling.

3.2.2.3 Acoustic Noise. Individual equipment items shall not emit audible noise greater than 40 dB in order to meet the Spacelab requirement in paragraph 3.1.1.3.4. However, occasional noise levels up to 80 dB for 3 minutes are permitted (i.e., fast tape rewind, etc.). Values are for 2 feet from accessible areas when installed for operation.

3.2.2.4 EMI Susceptibility and Radiation. Equipment units that facilitate the Spacelab design practices below are considered to meet the requirements of paragraph 3.1.1.3.6. Only those practices necessary to prevent transmission of excessive conducted and radiated interference external to equipments or to assure operability in the Spacelab environment as a consequence of Spacelab design practices are provided. Internal interference shielding practices not relevant to external interfaces are the user's responsibility.

3.2.2.4.1 Grounding and Signal Circuits. Equipment input/output connectors and input/output circuits shall facilitate the following practices:

NUMBER	REVISION LETTER							PAGE 16
	A							

- a. Spacelab, rack or equipment structures shall not be used as signal or power returns.
- b. Separate power grounds shall be returned to each Spacelab power source panel connection and shall not be switched.
- c. Signal and power returns shall not be shared.
- d. All equipment cases, enclosures and chassis shall facilitate electrical bonding to the vehicle structure with bonds not exceeding 10 milliohms dc impedance.
- e. No signal or power ground shall be connected with less than 1 megohm dc impedance to a chassis or case that is separately returned to structure ground. A chassis connected directly to signal ground shall be isolated from case ground to prevent system ground loops. Isolated chassis requires isolation from personnel contact where a potential shock or external short-circuit hazard could exist.
- f. Shield grounds shall be brought directly to the structure or case ground, using conductive epoxy, halo rings, or equivalent methods. Shields shall not be connected to a common wire and grounded at a remote point.
- g. Output impedances shall be less than 10K ohms and full-scale signal voltages at least 5 volts from ground reference.
- h. Use coaxial cable or twisted shielded pairs for outputs with frequency components greater than 150 kHz.
- i. Use separately returned shielded pairs or coax for signals where rise/fall rates of change are greater than 30 MA per microsecond and where maximum signal level is 30 MA or more.

3.2.2.4.2 Shielding (Generator or Receptor).

3.2.2.4.2.1 Signal Shields. Signal wire shields shall have the option of grounding or not grounding at a particular equipment input or output connector. Internal signal shielding shall be grounded internally and not rely upon connector feedthrough. Connector feedthrough shall be available to ground external shields internally for signals requiring them.

3.2.2.4.2.2 Metal Shields. Equipment external cases shall be of metal of at least 0.040 inch thickness that effectively encloses the unit in order to attenuate E-field and plane wave radiation and H-field radiation above 10 kHz. In addition, a thin copper shield may be required to maximize attenuation in the region of 60 Hz to 10 kHz for large current flows (~ 10 amp). External case shields shall not have openings that permit unattenuated radiation. Separate sections of the case shield shall be bonded or bolted, at intervals not exceeding 4 inches, through positive bare metal or conductive coating pressure contacts, or equivalent practices. RFI gaskets or equivalent may be required on some items.

3.2.2.4.3 Power Turn-On

- a. Turn-on in-rush current shall not exceed 2-1/2 times the steady-state value.
- b. Current rise times should not exceed 5 microseconds, and ringing should not exceed 5 cycles to reach ± 5 percent of steady state.

NUMBER

REVISION LETTER

A

PAGE 17

3.2.2.4.4 Output Signal Characteristics.

- a. Rise and fall times of digital signals shall be controlled in order to minimize unwanted EMI coupling.
- b. Rise and fall times should be larger than twice the propagation delay of the line from driving to receiving circuit if the line is not terminated in its characteristic impedance.

3.2.3 Typical Operating Requirements. The following criteria are solely the equipment user's responsibility. The criteria should be adequately considered to ensure successful completion of the experimenter's objectives. As a minimum, the equipment and/or its packaging should be designed to survive approximately a 3-minute launch, and operate over the required on-orbit duty cycle duration. Reuse on subsequent missions requires design/refurbishment for repeated launch and operational periods plus reentry and landing phases of approximately 15 minutes each (where worst-case vibration/acceleration may occur).

3.2.3.1 Equipment Size. To avoid delays and expense for custom-mounting racks and stowage containers, selected equipment items should conform to the standard sizes referenced in paragraph 3.1.1.2.1 and subparagraphs. While mass and volume are not critical, wherever practical, smaller sizes should be chosen consistent with integration with other equipments and objectives.

3.2.3.2 Electrical Power. Equipments shall be required to operate on one or more of the voltages identified in paragraph 3.1.1.2.2. The amount of power allotted to each equipment is not directly specified; however, the total power used by all users cannot exceed the value specified in paragraph 3.1.1.2.2, and the integrated loads on the 3-phase ac must remain within 20 percent between phase legs. Any one item that could use more than 1 kW peak shall consider means to reduce needs. Power during launch and reentry and landing shall only be utilized for critical purposes, to be identified and coordinated with Spacelab operations management. Only necessary caution and warning circuits will be accessible to emergency battery power.

3.2.3.3 Data Management and Command Control Compatibility. Outputs to be multiplexed or direct recorded through Spacelab remote data acquisition units (RAU's) require voltages in the range of 5 to 12 volts from ground reference for digital/discretes 1-state levels or analog full-scale levels. Scaling will be done in the data acquisition units. Output driver impedance shall be less than 10K ohms. Commands from RAU's will be 0 or 5 to 12 volts nominal at $\leq 1K$ ohm output impedance.

3.2.3.4 Atmosphere. Equipment will be exposed to an equivalent-air or air atmosphere per paragraph 3.1.1.3.1 and subparagraphs. The normal atmosphere also may contain the particle and chemical contaminants in paragraph 3.1.1.3.1.3.

3.2.3.5 Ambient Temperatures. Equipments shall withstand the operational and ground storage ambient air temperatures of paragraph 3.1.1.3.2.

3.2.3.6 Equipment Cooling. Where forced-air or coldplate operation is required, the equipment shall be capable of satisfactory operation with the cooling temperatures provided for by paragraph 3.1.1.3.2.1.

NUMBER

REVISION LETTER

A

PAGE 18

3.2.3.7 Shock, Vibration, Acceleration and Acoustic Noise. As a minimum, equipment shall operate satisfactorily for the duration of the mission after 3 minutes exposure to combined acceleration and random and sinusoidal vibrations of paragraph 3.1.1.3.3 and subparagraphs, plus the acoustic noise levels of paragraph 3.1.1.3.4.

3.2.3.8 Atomic Radiation. Atomic radiation from natural sources is well below the levels of any significance to equipment survival. Experiment equipment which emits radiation will be shielded to prevent impact to other experiment equipment. Equipments sensitive to radiation-induced noise or signals are the responsibility of the experimenter.

3.2.3.9 Contamination Susceptibility. The equipment shall be capable of operation for its total useful life in the particle and chemical environments of paragraph 3.1.1.3.1.3.

3.2.3.10 Zero-G Effects

3.2.3.10.1 Thermal Dissipation. Equipment points of power dissipation shall have adequate thermal conductance paths where forced-air convection cannot be relied upon. The effects of lack of gravity must be considered for points of power dissipation in "dead" air zones where positive thermal paths are not adequate. Equipments without provisions for internal-external air circulation that are designed to operate at 50,000 feet, relying on internal conduction to transfer heat to the external case, can be assumed to be acceptable. Equipments with self-contained fans and direct air flow over components or their heat sinks can be assumed acceptable. Components mounted on insulators in dead air spaces require cooling provisions if the dissipation on the insulator (circuit board) exceeds 2 watts per 20 square inches of circuit board area provided good contact or a hemispherical view to metallic structure exists. Density of electrical conductors that provide good conductance paths to reduce the effective thermal resistance can be used to justify increasing the foregoing dissipation rule.

3.2.3.10.2 Gravity-Dependent Functions. Equipment functions that rely upon gravity such as mechanical seating, containing liquid, material feed constraints, etc., require redesign for operation in weightless environments. Redesign should consider hysteresis effects on such as repeatability of operation.

3.2.3.10.3 Loose Accessories. Common table-top accessories or attachments such as adjusters, eyepieces, etc., must be provided with positive detents, lanyards, etc., to prevent loose items in Spacelab.

**ORIGINAL PAGE IS
OF POOR QUALITY**

NUMBER

REVISION LETTER

PAGE 19

A

4.0 PROCESS CONTROLS

4.1 Process Controls. Manufacture and test of experimenter equipment is the responsibility of the experimenter as he deems necessary to achieve his objectives and to verify meeting mandatory requirements of paragraph 3.2.1.

5.0 TEST AND VERIFICATION

5.1 Test and Verification. All experimenter equipment tests prior to integration are the responsibility of the experimenter, except as necessary to verify meeting the mandatory safety and environmental preservation requirements of paragraphs 3.2.1 and 3.2.2. Requirements of paragraphs 3.2.1 and 3.2.2 not met will result in the equipment not being allowed on Spacelab or operation being restricted to times when it will not adversely affect other equipment operations.

5.2 Off-Gassing. Where an assessment of specific materials in an equipment is not practical, an overall equipment test TBD may be performed to establish compliance. A general policy shall be implemented of operating equipment for an equivalent 150 hours normal operation after purchase or solvent cleaning to drive off the largest practical amount of out/off gases before taking into Spacelab.

5.3 EMI Tests. Equipment EMI generation and susceptibility tests may be accomplished in accordance with MIL-STD-462 to assure ability to operate unrestricted or to establish an operating plan to permit objectives to be achieved.

6.0 OTHER EQUIPMENT CONSIDERATIONS

6.1 Other Equipment Considerations. Equipment effectiveness factors such as reliability, maintainability, spares/maintenance equipment, reusability, etc., are the user's responsibility. On-orbit maintenance should be limited to simple replacement of plug-in-type modules.

APPENDIX E

NASA GENERAL EQUIPMENT SPECIFICATION
EC006M00000A

EC 006M00000 A

January 1973

GENERAL EQUIPMENT SPECIFICATION

PART I

PERFORMANCE, DESIGN AND
VERIFICATION REQUIREMENTS

FOR

EXPERIMENT AND GENERAL LABORATORY EQUIPMENT

TO BE SELECTED

FOR

PROJECT SORTIE LAB

APPROVED BY _____
(Preparing activity)

DATE _____

APPROVED BY _____

DATE _____

CONTENTS

<u>Paragraph</u>	<u>Title</u>	<u>Page</u>
1.0	SCOPE	1
2.0	APPLICABLE DOCUMENTS	1
3.0	REQUIREMENTS	3
3.1	DEFINITION	3
3.1.1	GENERAL DESCRIPTION	3
3.1.2	MISSIONS	3
3.1.3	OPERATIONAL CONCEPTS	3
3.2.	CHARACTERISTICS	4
3.2.1	PERFORMANCE	4
3.2.2	PHYSICAL	4
3.2.3	RELIABILITY	4
3.2.4	MAINTAINABILITY	4
3.2.5	OPERATIONAL AVAILABILITY	5
3.2.6	SAFETY	5
3.2.7	ENVIRONMENT	6
3.2.8	TRANSPORTABILITY/TRANSPORTATION	8
3.3	DESIGN AND CONSTRUCTION STANDARDS	8
3.3.1	SELECTION OF SPECIFICATIONS AND STANDARDS	8
3.3.2	GENERAL	8
3.3.3	AERONAUTICAL	8
3.3.4	CIVIL	8
3.3.5	ELECTRICAL	8
3.3.6	MECHANICAL	9
3.3.7	NUCLEAR	10
3.3.8	MOISTURE AND FUNGUS RESISTANCE	10
3.3.9	CORROSION OF METAL PARTS	10
3.3.10	CONTAMINATION CONTROL	10
3.3.11	COORDINATE SYSTEMS	10
3.3.12	INTERCHANGEABILITY AND REPLACE- ABILITY	12

CONTENTS (Cont'd)

<u>Paragraph</u>	<u>Title</u>	<u>Page</u>
3.3.13	IDENTIFICATION AND MARKING	12
3.3.14	WORKMANSHIP	12
3.3.15	HUMAN PERFORMANCE/HUMAN ENGINEERING	12
3.4	LOGISTICS	12
3.5	PERSONNEL AND TRAINING	13
4.0	VERIFICATION	13
4.1	GENERAL	13
5.0	PREPARATION FOR DELIVERY	13
6.0	NOTES	13
6.1	INTENDED USE	13

1.0 SCOPE

This specification establishes the requirements for performance, design and verification of existing commercially available or off-the-shelf Experiment and General Laboratory Equipment, for orbital operation in a pressurized compartment in Manned Sortie Missions. Equipment which meets the quantitative performance requirements of this specification and which has particularly desirable attributes for a specific use, shall not be rejected from further consideration due to failure to meet a qualitative requirement contained herein. Candidate equipment may be selected for modification or by limited waiver of requirements. Typical equipment which may be procured to these requirements is listed in Section 6.0.

2.0 APPLICABLE DOCUMENTS

The following documents, of the exact issue shown, form a part of this specification to the extent specified herein. In the event of conflict between documents referenced, and other detailed contents of this specification, the detailed requirements herein shall be considered superseding.

SPECIFICATIONS

Marshall Space Flight Center

MSFC-SPEC-101B
March 15, 1971

Flammability, Odor, and Off-gassing requirements and Test Procedures for Materials in Environments which Support combustion.

Referenced
Paragraph

3.2.6.3

Military

MIL-F-14072A
Sept. 11, 1968

Finishes for Ground Signal Equipment

3.3.9
Figure 1
Footnote 3

STANDARDS

Military

MIL-STD-461A(Not. 4)
Feb. 9, 1971

Electromagnetic Interference Characteristics Requirements for Equipment

3.3.5.1

MIL-STD-462(Not. 3)
Feb. 9, 1971

Electromagnetic Interference Characteristics, Measurement of

3.3.5.1

OTHER PUBLICATIONS

Atomic Energy Commission

Rules and Regulations

Title 10, Part 71
Dec. 31, 1968

Packaging of Radioactive
Material for Transit

3.2.7.3

MISCELLANEOUS

Marshall Space Flight Center

10M 33222
August 15, 1972

Man/Systems Design
Requirements for
Sortie Lab

3.3.15

3.0 REQUIREMENTS

3.1 DEFINITION

3.1.1 GENERAL DESCRIPTION

The Experiment and General Laboratory Equipment shall consist of selected hardware and equipment to support experiments and experiment related activities and procedures to be employed in missions such as described in 3.1.2.

3.1.2 MISSIONS

Typical Sortie Lab missions consist of support for and conduct of experiments in the following general areas:

- a. Astronomy: UV, IR, High Energy Stellar.
- b. Physics: Space Physics (Orbital Environment) - Plasma Physics - Molecular and Chemistry Phenomena.
- c. Earth Observations: Weather - Resources - Pollution.
- d. Communications/Navigation: Optical Frequency Demonstration - Surveillance, Search, Rescue - Laser Ranging - Satellite Navigation for Terrestrial Users.
- e. Materials Sciences/Manufacturing: Materials Composites, Casting, Crystal Growth - Supercooling - Biological Separation.
- f. Technology: Contamination - EVA - Advanced Spacecraft Systems Tests - Teleoperations.
- g. Life Sciences: Medical Research - Plant Research - Life Support - Manned System Integration.

3.1.3 OPERATIONAL CONCEPTS

3.1.3.1 Mission Duration: The equipment shall be required to operate during a mission of 7 to 30 days duration.

3.1.3.2 Equipment design shall be compatible with prelaunch and launch operations which dictate minimal access while on the launch pad. Any requirement for on-pad loading of equipment consumables or perishables shall be identified.

3.2 CHARACTERISTICS

3.2.1 PERFORMANCE

3.2.1.1 GENERAL PERFORMANCE

3.2.1.1.1 Commercially available or off-the-shelf equipment and hardware shall be utilized to the fullest extent possible, consistent with the requirements stated herein.

3.2.1.1.2 Equipment covered by this specification shall not be required to be operating during launch, reentry or landing.

3.2.1.2 ELECTRICAL FUNCTIONAL AREA

Type of Power: 28 VDC \pm 4.0 VDC and 115/200 VAC
 \pm TBD, 3 ϕ , 400 Hz.

NOTE: If 60 Hz power is required for an experiment or instrument, individual inverters will be provided as part of the experiment.

3.2.2 PHYSICAL

3.2.2.1 MASS PROPERTIES

Any limitations and restrictions are to be determined.

NOTE: A maximum weight of the equipment will be specified (specification weight).

3.2.2.2 DIMENSIONS AND VOLUMES

Any general limitations and restrictions are to be determined.

3.2.3 RELIABILITY

Non-catastrophic and non-propagating failures shall be permitted when compatible with individual experiment reliability goal.

3.2.4 MAINTAINABILITY

3.2.4.1 The design shall provide for accessibility, rapid fault isolation, ease of remove/replace activities, and the use of standard tools and test equipment. Ground maintenance shall be the normal mode of maintenance.

3.2.4.2 Consideration shall be given to location of test points and adjustments for ease of maintenance without disconnecting electrical connectors.

3.2.4.3 Accessibility and ease of operation of latches, lockdowns, fasteners, etc., shall be considered for ease of maintenance.

3.2.4.4 For units or components requiring frequent adjustments or for installation and removal of cassettes, rolls, samples, film, etc., front panel access is desirable.

3.2.5 OPERATIONAL AVAILABILITY

Not Applicable

3.2.6 SAFETY

3.2.6.1 The equipment shall not degrade the safety of the Sortie Lab, Space Shuttle, Ground or Flight Personnel or violate safety levels and requirements stated herein either during normal operations, maintenance or failure modes.

3.2.6.2 The equipment shall contain no ordnance devices.

3.2.6.3 All materials flammability, odor, and offgassing characteristics shall be identified for the application usage including relative locations, quantities and configurations as defined in MSFC-SPEC-101. The atmosphere pressures and gas mixtures shall be as specified in 3.2.7.2 herein for on-orbit operations.

Any materials which do not meet the criteria of MSFC-SPEC-101 for Group I must be approved by MSFC for use on an individual basis.

In lieu of performing any testing required by MSFC-SPEC-101, including batch/lot testing, materials lists with usage configurations may be submitted to MSFC for initial determination of acceptability. Results of any previous tests and supportive analyses may be submitted to MSFC for determination of compliance with the intent of MSFC-SPEC-101. Final evaluation and approval shall be conditioned on results of individual materials and/or end item configuration tests performed by or required by MSFC on an individual basis.

The requirement of MSFC-SPEC-101, paragraph 3.1, Materials Control Plan, is deleted from equipment covered by this specification. The requirement of MSFC-SPEC-101, paragraph 3.2.2.1.1 for testing metallic materials for LOX/GOX compatibility shall only apply to equipment used in high pressure LOX/GOX systems (Type D usage defined in MSFC-SPEC-101, Table 1).

3.2.6.4 Catastrophic and critical hazards whereby environment, personnel error, design characteristics, procedural deficiencies, or subsystem malfunction may result in loss of personnel capability or loss of system shall be eliminated or controlled by means of safety devices, warning devices or procedures. Residual hazards for which controls cannot be developed or provided shall be identified.

3.2.7 ENVIRONMENT

3.2.7.1 Non-Operational

Temperature: -40°C to $+75^{\circ}\text{C}$ (-40°F to $+167^{\circ}\text{F}$).

Pressure: 0-15 psia/Rate of Change Approx. 2 Min. from maximum P to minimum P

Atmosphere composition: 80% N_2 ; 20% O_2 , \pm TBD.

Relative Humidity: TBD

Acceleration: ± 4.5 g (all axes)

Launch Status: Passive

Contaminants: filter air to class 100,000 as defined in 3.3.10.1

Emergency Landing Shock: Longitudinal Lateral Vertical
8g 4.5g 4.5g

Equipment need not operate thereafter, but shall retain package and mounting integrity.

Vibration: TBD

Acoustic: TBD

Shock: TBD

ORIGINAL PAGE 13
OF POOR QUALITY

3.2.7.2 Operational

A nominal "shirt-sleeve" environment shall exist for equipment operation.

Temperature: 18°C to 29°C (65°F to 85°F)

Humidity: 45 ± 5% R. H. at 21°C (70°F)

Pressure: 14.7 psia ± 0.2 psi

Atmosphere Composition:

Total Pressure	14.7 psia
O ₂ Partial Pressure	3.1 psia
CO ₂ Partial Pressure	<7.6 mm Hg
H ₂ O Partial Pressure	6-11 mm Hg

Atmosphere Contamination:

Cleanliness: Class 100,000 as defined in 3.3.10.1

Maximum Methane Concentration < $5 \times 10^{14} \frac{\text{Molecules}}{\text{CM}^3}$
(10 PPM by weight)

Maximum Hydrocarbons (Molecular Weight > 200)
< $10^{14} \frac{\text{Molecules}}{\text{CM}^3}$

Acceleration: < 10^{-2} g's

Vibration: TBD

Touch Temperature: 3°C to 46°C (38°F to 115°F)
Equipment surface temperatures below 16°C (60°F) should be avoided to prevent atmospheric moisture from condensing on its surface.

Noise Level: Equipment and main equipment interfaces shall be compatible with an overall Sortie Module noise level less than the following:

Frequency Range (Hz)	Sound Pressure Levels
63	70 db
125	70 db
250	60 db
500	60 db
1000	55 db
2000	55 db
4000	55 db
8000	60 db
16000	60 db

Pure tone components shall be not more than 10 db above the octave band level containing the pure tone frequency.

3.2.7.3 **Radioactive Sources.** The quantities and/or concentrations of radioactive sources such as calibration sources and display dials shall not exceed the exemptions specified in USAEC Rules and Regulations, Title 10, Part 71, Paragraph 7.1.5. Also any equipment susceptible to radiation degradation shall be defined.

3.2.8 **TRANSPORTABILITY/TRANSPORTATION**

Equipment shall be compatible with commercial transportation systems.

Vibration and shock criteria - TBD

3.3. **DESIGN AND CONSTRUCTION STANDARDS**

3.3.1 **SELECTION OF SPECIFICATIONS AND STANDARDS**

Not Applicable

3.3.2 **GENERAL**

Not Applicable

3.3.3 **AERONAUTICAL**

TBD

3.3.4 **CIVIL**

Not Applicable

3.3.5 **ELECTRICAL**

3.3.5.1 Electromagnetic Compatibility (EMC) characteristics and measurements shall be in accordance with MIL-STD-461 and MIL-STD-462. The following frequency ranges shall be considered in establishing equipment emission and susceptibility levels. (Frequencies TBD)

3.3.5.1.2 Shorting clips or springs shall not be used in electrical or electronic connectors. Test equipment used shall be fitted with cable connectors that mate with the experiment equipment connectors or sockets. In no case shall socket connectors be connected by holding meter probes against the pins or by attachment of alligator clips.

3.3.5.1.3 All electrical connectors, plugs and receptacles shall be positively keyed to prevent incorrect connection with other accessible connectors, plugs and receptacles. Connectors and wiring junctions shall be sealed from moisture condensation.

3.3.5.1.4 There shall be safety interlocks where necessary to insure personnel safety.

3.3.6 MECHANICAL.

3.3.6.1 Mechanical design shall generally meet the following requirements:

3.3.6.1.1 Material which can shatter, such as glass, shall not be used unless positive protection is incorporated to prevent fragments or dust from entering the habitable environment.

3.3.6.1.2 Surfaces which are expected to be exposed to continuous or extensive abrasion and rubbing by the crew shall not be painted or coated with materials that are subject to flaking or peeling.

3.3.6.1.3 Where possible, actuating devices shall be made an integral part of the equipment to be operated. Detachable items, such as handles, pins, and ratchets shall be secured by lanyards or similar devices which will not compromise crew safety.

3.3.6.1.4 Factors of Safety

3.3.6.1.4.1 Package integrity and structural mounting provisions load carrying capability shall be based on the following minimum factors of safety in lieu of performing static load structural testing:

Yield Factor of Safety	=2.0
Ultimate Factor of Safety	=3.0

3.3.6.1.4.2 Hydraulic and pneumatic systems, where used, shall meet the following minimum requirements:

(1) Lines and Fittings, less than 1.5 inch diameter

Proof Pressure	=2.0 x limit pressure
Ultimate Pressure	=4.0 x limit pressure

(2) Lines and Fittings, 1.5 inch diameter or greater

Proof Pressure	=1.2 x limit pressure
Ultimate Pressure	=1.5 x limit pressure

(3) Hydraulic and Pneumatic Tanks & High Pressure Vessels

Proof Pressure	=1.5 x limit pressure
Ultimate Pressure	=2.0 x limit pressure

(4) Actuating Cylinder, Valves, Filters, Switches

Proof Pressure	=1.5 x limit pressure
Ultimate Pressure	=2.0 x limit pressure

ORIGINAL PAGE IS
OF POOR QUALITY

3.3.6.2 MATERIALS

The following material requirements shall be in addition to the flammability, odor and offgassing requirements of paragraph 3.2.6.3.

3.3.6.2.1 Equipment covered by this specification shall not contain mercury, liquid or gaseous halogenated hydrocarbons, or polyvinyl chloride (PVC).

3.3.6.2.2 Use of cadmium plating should be avoided in equipment containers subject to elevated temperatures (above 450°F) or where exposed cadmium in contact with breathing gas could reach temperatures that would generate toxic fumes such as might result from electrical short or fire.

3.3.7 NUCLEAR

Not Applicable.

3.3.8 MOISTURE AND FUNGUS RESISTANCE

All materials used shall be non-nutrient to fungus growth, or shall be treated so the exposed surfaces will be fungus resistant, except that nutrients may be used inside of enclosures that will be hermetically sealed. Methods of treatment for fungus resistance shall be compatible with the treated and adjacent materials.

3.3.9 CORROSION OF METAL PARTS

Metals shall be corrosion resistant or shall be processed to resist corrosion. Dissimilar metals shall not be used in intimate contact unless suitably protected against electrolytic corrosion. Compatible metal couples are defined in Figure 1. Any protective coating used shall meet the requirements of paragraph 3.2.6.3.

3.3.10 CONTAMINATION CONTROL

3.3.10.1 Equipment material or equipment finish shall not flake off, generate dust or contain releasable particles that could degrade the class 100,000 spacecraft environment. All surfaces shall be capable of being cleaned with suitable solvents to maintain surfaces visibly clean. Class 100,000 environment is defined as not more than 100,000 particles, 1/2 micron in size or larger per cubic foot and not more than 700 particles, 5 microns in size or larger per cubic foot.

3.3.10.2 Contamination protection devices such as aperture doors, window covers, electrical heaters, and dust free storage containers shall be used as necessary.

3.3.11 COORDINATE SYSTEMS

Not Applicable

Group No.	Metallurgical category <u>3/</u>	Compatible couples <u>1/</u>
1	Gold, solid and plated; gold-platinum alloys; wrought platinum	○
2	Rhodium plated on silver-plated copper	● ○
3	Silver, solid or plated; high silver alloys	● ● ○
4	Nickel, solid or plated; monel metal, high nickel-copper alloys, titanium	● ● ○
5	Austenitic stainless steels; copper, solid or plated; low brasses or bronzes; silver solder; German silver; high copper-nickel alloys; nickel chromium alloys	● ● ● ○
6	Commercial yellow brasses and bronzes	● ● ● ○
7	High brasses and bronzes; naval brass; muntz metal	● ● ● ○
8	18% chromium-type steels	● ● ● ● ○
9	Chromium, plated; tin, plated; 12% chromium type steels	● ● ● ● ○
10	Tin-plate; terneplate; tin-lead solders	● ● ● ● ○
11	Lead, solid or plated; high lead alloys	● ● ● ● ○
12	Aluminum, wrought alloys of the Duralumin type	● ● ● ● ○
13	Iron, wrought, gray, or malleable; plain carbon and low alloy steels, armco iron	○ ● ● ● ●
14	Aluminum, wrought alloys other than Duralumin type; aluminum, cast alloys of the silicon type	● ○ ● ● ●
15	Aluminum, cast alloys other than silicon type; cadmium, plated and chromated	● ● ○ ● ●
16	Hot-dip-zinc plate; galvanized steel	● ○
17	Zinc, wrought; zinc-base die-casting alloys; zinc, plated	●
18	Magnesium and magnesium-base alloys cast or wrought	● <u>2/</u>

FIGURE 1. COMPATIBLE COUPLES

Footnotes 1/, 2/, and 3/ (see next page)

Footnotes for Figure 1.

- 1/ Members of groups connected by lines will form permissible couples. O indicates the most cathodic member of the series, ● an anodic member, and the arrows anodic direction.
- 2/ Aluminum alloys 5052, 5056, 5356, 6061 and 6063 are considered to be compatible with magnesium alloys.
- 3/ A more detailed listing of metals in each group is given in Specification MIL-F-14072.

3.3.12 INTERCHANGEABILITY AND REPLACEABILITY

- 3.3.12.1 Interchangeability shall exist between identical replaceable parts, assemblies, subassemblies and supplies, regardless of manufacturer or supplier. All parts having the same part numbers, regardless of source, shall be functionally and dimensionally interchangeable.

3.3.13 IDENTIFICATION AND MARKING

Major assemblies and replaceable units shall be adequately identified and marked.

3.3.14 WORKMANSHIP

Workmanship on all equipment (electronic, electrical and mechanical) shall be in accordance with good commercial practices and parts shall be free of burrs, sharp edges, or any other damage or defect that could make the part (or equipment) unsatisfactory for the purpose intended.

3.3.15 HUMAN PERFORMANCE/HUMAN ENGINEERING

Human performance/human engineering requirements shall be as specified in MSFC Drawing 10M33222, "Man/Systems Design Requirements for Sortie Lab".

3.4 LOGISTICS

Not Applicable.

3.5 PERSONNEL AND TRAINING

Personnel usage and skills requirements shall be identified by the contractor for all items of equipment. Familiarization and operational training requirements shall also be identified.

4.0 VERIFICATION**4.1 GENERAL**

Unless otherwise specified in the contract or purchase order, the contractor is responsible for the verification of all performance and design requirements contained herein and shall submit a written statement verifying that all requirements herein have been met. This verification statement is subject to approval and the procuring activity reserves the right to perform any inspection, operation or test deemed necessary to assure that material and services conform to the prescribed requirements.

5.0 PREPARATION FOR DELIVERY

The Contract End Items shall be prepared for delivery in accordance with the terms of the contract. Marking shall be adequate for proper identification. Preservation, packaging, packing, handling and shipping requirements shall be compatible with commercial transportation unless otherwise specified by the contractor. Any special unpacking or handling procedural requirements shall be defined by the contractor.

6.0 NOTES**6.1 INTENDED USE**

The intended use of this specification is to procure Experiment and General Laboratory Equipment and Instrumentation for support of the Sortie Lab.

Typical equipment to be procured by this specification includes but is not limited to the following:

Receivers and Transmitters - Microwave thru Laser

Optical Devices

Controls/Displays

Test, Checkout and Monitor Equipment

Calibration Devices

Data Recorders

Clocks (event timing)

Cameras - TV and Film

**Cryogenic Devices - Solid Cryogenic Refrigerators, Cryostats,
Cryogenic Liquid Storage**

Equipment Racks and Consoles

Spark Chamber

Spectrographic Devices